Rejûvaskin® | Scar Heal
Advanced Skin Therapy

Scar Fx ®

Silicone Sheeting

Application and Indications

Scar Fx® silicone sheeting is available in 12 surgery specific sizes, four (4) of which are specifically designed for the breast. Scar Fx® should overlap the scar area by at least 1cm on all sides.

The therapeutic effect of silicone gel formulations on scars is well documented in the medical literature. The actual mechanism of action remains unclear, although several hypotheses exist. Most researchers agree that silicone acts by creating a hydrated occluded environment decreases capillary activity by reducing fibroblast induced collagen deposition and scar hypertrophy. Silicone sheeting and silicone gel decreases hyperemia minimizes fibroblast production of collagen and promotes wound flattening.

Apply to clean, dry, fully healed skin by removing the protective sheeting from both sides of the silicone piece. Place the sticky side onto the scar area. Scar Fx® is self-adhesive. When Scar Fx® becomes less sticky, you can secure it in place with Hypo-allergenic tape or a suitable type wrap.

Hygiene

Scar Fx® should be washed with mild soap and water and air dried. Do not use talcum powder, moisturizer or any ointments under the silicone sheet. The skin being treated must be kept clean and dry. When not using Scar Fx®, store in a zip lock bag. Scar Fx® is reusable and one-piece Scar Fx® will last between 6-8 weeks when properly cared for.

Duration of use

For optimal results use for a minimum of 8-12 hours per day. Use until scar maturity and no further improvement is noted.

Component Description Chemical Name Composition

Non-Tacky Silicone Backing/ Liquid Silicone: Reinforced Dimethyl Methylvinyl Siloxanes Tacky Silicone Gel Surface/ High Strength Firm Silicone Gel: Dimethyl Methylhydrogen, Methylvinyl Silicone Gel

FUNCTIONS

Creates occlusive barrier for ideal healing environment
Improves appearance of surgical and traumatic scars
Flattens and softens raised scars
Reduces redness and discoloration
Relieves itchiness and discomfort

FEATURES

Best option for non-invasive scar management
100% healthcare grade silicone
Available in 13 uniquely shaped self-adhesive sheets
Soft, pliable, and contours easily to the body
Used by dermatologists and plastic surgeons worldwide

Affirmed By: _____ Date: 2/4/19
Thomas Christenberry, President, Scar Heal

SCAR HEAL LIMITED LIABILITY SEAL 1988 FLORIDA COMPANY

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