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MYOSTAR™ Injection Catheter

- STERILE. Sterilized using ethylene oxide.
- For single use only. Do not resterilize.
- Do not use if the package is open or damaged.

CATHETER DESCRIPTION

The MYOSTAR™ Injection Catheter is a multi-electrode, percutaneous catheter with a deflectable tip and injection needle designed to inject agents into the myocardium.

The tip of the Injection Catheter is equipped with a Biosense Webster location sensor and a retractable, hollow 27 gauge needle for fluid delivery. Tip deflection is controlled at the proximal end by a tubular handpiece in which a piston slides. When the piston is pushed forward with the thumbknob, the tip is deflected (curved). When the piston is pulled back, the tip straightens. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. Needle movement is controlled by moving a piston on a second handle, the injection handle, located proximal to the deflection handle. Push and hold the thumbknob forward on the injection handle, to retract the needle, release the thumbknob. This handle has a standard Luer fitting for connection to a syringe. The catheter interfaces with the Biosense Webster NOGA™ System, which is standard recording equipment. The dead space of the catheter is approximately 0.1 cc.

CONSULT THE LOCAL DISTRIBUTOR OR MANUFACTURER FOR THE APPROPRIATE INTERFACE CABLES.

INDICATIONS FOR USE

The MYOSTAR™ Injection Catheter, along with related accessory devices, is intended for targeted delivery of diagnostic and therapeutic agents indicated for the myocardium. The MYOSTAR™ Injection Catheter should only be used under the guidance of a NOGA™ Cardiac Navigation System. Injection without the cardiac navigation system is contraindicated.
All agents must be used in accordance with the specific agent manufacturer's Instructions for Use.

CONTRAINDICATIONS

- Use of the catheter may not be appropriate for patients with prosthetic valves.
- Do not use this device in patients with active systemic infection.
- Use of the catheter is contraindicated for patients with left ventricular thrombus.
- The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch. Further, the catheter is not an over the wire model and does not utilize a transseptal sheath system.
- The retrograde trans-aortic approach is contraindicated in patients with aortic valve replacement.
- Patients who cannot receive anticoagulation therapy.

WARNINGS

Cardiac catheterization procedures present the potential for significant X-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the X-ray beam intensity and duration of the fluoroscopic imaging. Catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must therefore be given for the use of the device in pregnant women.

The MYOSTAR™ Injection Catheter should only be positioned using the guidance of a NOGA™ Cardiac Navigation System. Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
Animal studies have shown that liquids injected into the myocardium using the MYOSTAR™ Injection Catheter may be injected more deeply into tissue than anticipated due to the forward jet of injected material. The injection depth may vary for various injected materials, as depth is dependent on the nature of the injected material. Viscosity of the injected material, needle extension length, and cardiac wall thickness should be considered prior to injection.

PRECAUTIONS

- The MYOSTAR™ Injection Catheter should only be used by physicians trained in interventional procedures and cardiac mapping.
- Do not attempt to operate the MYOSTAR™ Injection Catheter, or the Biosense Webster NOGA™ Cardiac Navigation System prior to completely reading and understanding the applicable Instructions for Use.
- All agents delivered must be used in accordance with the specific agent manufacturer's Instructions for Use.
- Verify compatibility of the material to be injected prior to first clinical use with the MYOSTAR™ Injection Catheter.
- MYOSTAR™ Injection Catheter procedures should be performed by appropriately trained personnel in a fully equipped laboratory.
- Do not use the injection catheter for mapping. NOGA™ mapping should be performed with the NOGASTAR® Catheter. The NOGASTAR® Catheter shall be used according to its Instructions for Use.
- The MYOSTAR™ Injection catheter must be flushed with agent to be injected prior to use in the patient.
- Store in a cool, dry, dark place.
- The sterile packaging and catheter should be inspected prior to use. Do not use if package is open or damaged.
- The MYOSTAR™ Injection Catheter is intended for single patient use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Do not expose catheter to organic solvents such as alcohol.
- Do not immerse proximal handle or cable connector in fluids; electrical performance could be affected.
- Regularly inspect and test re-usable cables and accessories.
- Make sure the small vent hole(s) at the connector end of the handpiece is patent; if clogged, the piston could force air into the catheter shaft and possibly into the bloodstream.
- Target cardiac wall thickness should be determined prior to insertion of the MYOSTAR™ Injection Catheter.
- The patient should be in a stable sinus rhythm or a similarly paced rhythm, prior to initiation of an injection procedure.

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COMPLICATIONS AND ADVERSE REACTIONS

Procedures requiring the introduction of a catheter into the heart should not be attempted by physicians unfamiliar with the possible complications or adverse events.

A number of serious adverse reactions have been documented for intracardiac catheter procedures:

- Pulmonary embolism
- Myocardial infarction
- Stroke
- Cardiac tamponade
- Death

Possible Complications include, but are not limited to:

- Vascular bleeding/local hematomas
- Thrombosis
- Arteriovenous fistula
- Pseudoaneurysm formation
- Thromboembolism
- Vasovagal reactions
- Cardiac perforation
- Cardiac tamponade
- Thrombus
- Air embolism
- Arrhythmias
- Valvular damage
- Pneumothorax
- Hemothorax
- Infection

STERILIZATION/EXPIRATION DATE

This catheter has been sterilized with ethylene oxide gas. Do not use the catheter if the package is open or damaged. Use the catheter prior to the expiration date shown on the package label.

DISPOSAL

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

SUGGESTED INSTRUCTIONS FOR USE

- Remove catheter from package and place in a sterile work area.
- Connect the injection catheter to the Biosense Webster NOGA™ Cardiac Navigation System junction box using the appropriate cables.
- Prior to beginning the procedure, measure the needle as shown in Figure 1 to ensure proper needle extension.


Figure 1

1. Remove Aortic Arch Simulator (AAS) & Ruler packaged with the MYOSTAR™ Injection Catheter and place on a sterile field.
2. Place the extended injection catheter on a flat sterile surface.
3. Insert the catheter into the AAS.
4. Push the thumbknob on the injection handle forward to extend the needle. Hold the thumbknob and measure the needle length in millimeters.
5. With the catheter in the AAS, flex the tip to 90 degrees. Measure the needle length. **The needle length should measure no more than 50% of the target tissue depth**.
6. With the catheter in the AAS, measure the needle length with the catheter tip straight. **The needle length should measure no more than 50% of the target tissue depth**.
7. Use the thumbscrew on the injection handle thumbknob to make adjustments to the needle length.
8. A counterclockwise rotation will lengthen the extended needle. Adjustments to the needle length should be made with the catheter tip straight and the needle extended. The catheter remains in the AAS [illegible] tire needle measurement and adjustment process.
9. Remove the AAS from the MYOSTAR™ Injection Catheter. Place on sterile field.
10. Document the needle length in two positions while on a flat sterile surface: 1) while straight and 2) while deflected at 90 degrees (see Figure 1).

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- Connect the agent-filled 1 cc Luer hub syringe to the injection port.
- With the needle extended, flush the catheter with the agent to be injected.
- Release the thumbknob on the injector handle to retract the needle.
- Inspect the distal tip of the catheter to insure that the distal needle tip is fully retracted.
 WARNING: The distal needle tip should not protrude from the distal tip of the catheter. Extension of the needle during catheter manipulation could result in tissue damage.
- Insert the catheter.
 WARNING: Always pull the thumbknob back to straighten the catheter tip prior to insertion of the catheter. Failure to do so may result in tissue damage.
- Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and the NOGA™ Cardiac Navigation System to aid in proper positioning.
 WARNING: To avoid thromboemboli, intravenous heparin is used when entering the left side of the heart during the injection catheterization procedure.
- The catheter tip can be deflected to facilitate positioning by using the thumbknob to vary tip curvature. Pushing the thumbknob forward causes the catheter tip to bend; when the knob is pulled back, the tip straightens.
- Prior to needle extension, determine that the following criteria have been met:
 1. Target tissue is not near the mitral valve annulus or apex of the heart.
 2. Catheter tip is perpendicular to the myocardial wall using NOGA™ Cardiac Navigation System mapping and fluoroscopy.
 3. Loop Stability is < 3 mm.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Released Date: 8/11/2017

4. Unipolar voltage of ≥ 7 mV is observed if the target tissue is viable myocardium or < 7 mV for infracted myocardium.

- When it has been determined that these criteria have been met, extend the needle from the catheter distal tip by moving the piston on the injection handle and holding it in that position. After extending the needle into the myocardium, ensure that:
 1. Premature ventricular contraction(s) has been observed.
 2. No catheter pushback has been observed under fluoroscopy.
 3. The recommended dosage is delivered via a slow and steady injection by actuating the syringe.

 WARNING: Do not inject near the mitral valve annulus or apex of the heart. Injecting in these areas may result in effusion or perforation.

 CAUTION: The catheter tip should be perpendicular to the application site prior to injecting into the myocardium.
- Before attempting to reposition or maneuver the catheter, retract the needle into the distal catheter tip, by releasing the thumbknob on the injector handle.

 WARNING: Always retract the needle before attempting to reposition the catheter. Failure to do so may result in tissue damage.
- AAS & Ruler are for use with a single catheter only. Do not re-sterilize or re-use MYOSTAR™ Injection Catheter, AAS, or Ruler.
- Additional injections may be delivered to other sites using the same catheter.
- Injection catheter must be flushed with agent to be injected prior to use in the patient.
- Prior to removal of the catheter, for additional safety, confirm that the thumbknob has been pulled back completely and the needle has been retracted. Remove the catheter and dispose of it in an appropriate manner.

 WARNING: Always pull the thumbknob back to straighten the catheter tip and retract the needle before catheter withdrawal. Failure to do so may result in tissue damage.

If any of the following events / symptoms occur, the procedure should be temporarily halted and the patient should be reevaluated for suitability to continue with the treatment under investigation:

- Persistent complaints of chest pain;
- Complains of cardiac pain associated with injections;
- Persistent hypotension;
- Complaints of shortness of breath;
- ICD shocks to stop ventricular tachycardia (VT);
- DC cardioversion or defibrillation for VT;
- There is any question as to the location of the catheter tip in relation to vasculature or the left ventricle (LV); and
- Any unanticipated change in level of consciousness or neurological status.

The procedure should be aborted in the event that:

- Cardiac tamponade is strongly suspected or confirmed;
- Patient experiences recurrent, hemodynamically unstable sustained runs of VT > 150 beats / minute
- The patient experiences ventricular fibrillation (VF);
- Identification of thrombus in the LV or the aorta that was not previously present on echo or left ventriculogram;
- An aortic dissection is suspected or confirmed.

After the injections are completed, remove the MYOSTAR™ Injection Catheter at the conclusion of the endomyocardial injection portion of the procedure.

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DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Released Date: 8/11/2017

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M-5276-120L © Biosense Webster, Inc. 2005-2017 Printed in USA
BIOSENSE WEBSTER, the Biosense Webster logo, MYOSTAR, NOGA, AND NOGASTAR are trademarks of Biosense Webster, Inc.
The third-party trademarks used herein are trademarks of their respective owners.