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
DIANATAL®
obstetric gel

Product information

Made in Switzerland
SWISS INNOVATION

Dianatal® obstetric gel Set: product and application

Dianatal® obstetric gel does not contain any pharmaceutically active ingredients. Each sterile single-use tube contains 5 ml Dianatal® obstetric gel. Dianatal® obstetric gel is sterile and filled in sterile single-use tubes and packed in single use sterile blisters. Dianatal® obstetric gel is an odourless, almost colourless gel. The gel is highly viscous, with a slightly acidic pH of 5.5 – 6.5 and electrical conductivity of 5.0 – 9.0 mS . cm-1.

Dianatal® obstetric gel application

Dianatal® obstetric gel is applied after each vaginal examination during the entire delivery process after the onset of regular contractions, starting optimally with the vaginal examination around 3-4 cm dilatation of the cervix.

Dianatal® obstetric gel should be used throughout the entire labour. Approx. 3-5 ml of Dianatal® obstetric gel per application is applied with the tube or manually to the presenting part. The application of Dianatal® obstetric gel is consecutively performed in a 1 to 4 hours interval to ensure an optimal lubrication film on the birth canal to make delivery safer & easier for mother and baby. Application is best performed directly after vaginal examination. In very dry vaginal conditions ~ 2 ml fluid should be added after the gel is applied. An additional application of gel should be made 15-30 minutes after rupture of the amniotic sac. Dianatal® can be applied manually or by inserting the tube into the vagina. Do use a dry towel to deliver the baby to prevent slipping.

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DIANATAL®
obstetric gel

Important: english is masterfile for content

CE signs explanation: 2
GB: English - Product information – Master file for content 4
DE: German - Produktinformation 5
FR: French - Information concernant le produit 7
ES: Spanish - Información sobre el producto 9
IT: Italien - Informazioni sul prodotto 10
SLO: Slovenian - Navodilo za uporabo 12
PL: Polish - Informacja o produkcie 14
SK: Slovak - Informácie o výrobku 16
RO: Romanian - Instrucţiuni de utilizare 17
HU: Hungarian - Termékinformáció 19
CZ: Czech - Příbalová informace 21
LV: Latvian - Lietošanas instrukcija 22
LT: Lithuanian - Naudojimo Instrukcija 24
EST: Estonian - Kasutusjuhend 26
NL: Dutch - Productinformatie 27
RU: Russian – Информация о продукте 29
TR: Turkish - Ürün bilgileri 31
PT: Portuguese - Informação do produto 33
GR: Greec - Πληροφορίες προϊόντος 35
BG: Bulgarian - ИНСТРУКЦИИ ЗА УПОТРЕБА 37
RS: Serbian - Uputstvo za upotrebu 38
HR: Croatian - Uputa za upotrebu 40
HEB: Hebrew - 44
ARB: Arabian 47
FAR: Farsi 49

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CE signs explanation:

GB: English - Product information – Master file for content
not for reuse - sterile by aseptic process techniques – do not re-sterilize - do not use if packaging is damaged – consult instructions for use - temperature limitation – store in dry place - store away from sunlight - Latex free - DEHP / PHT free

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GB: English - Product information – Master file for content

Dianatal® obstetric gel - For an easier birth

For use by qualified obstetric professionals (midwifes and physicians) only.
Please read this information before using Dianatal® obstetric gel.

Application form and contents:

Dianatal® obstetric gel Set: sterile single-use tubes containing 5 ml obstetric gel:

Properties, galenic form and ingredients

Dianatal® obstetric gel is packed in single-use tubes. Dianatal® is non-allergenic, bioadhesive, electroconductive, slightly acidic, free of preservatives and latex and does not cause irritation to mucous membranes and or eyes.
Dianatal® obstetric gel does not contain any pharmaceutically active ingredients.
Dianatal® Obstetric Gel is an odourless, almost colourless gel packed in sterile single-use packaging. The gel is highly viscous, with a slightly acidic pH of 5.5 – 6.5 and an electrical conductivity of ≥ 5 mS . cm-1.

Ingredients:

One single-use tube of 5 ml Dianatal® obstetric gel contains following ingredients:
Propylene glycol, carbomer, hydroxyethylcellulose, purified water

Indication and uses of Dianatal® obstetric gel

Dianatal® obstetric gel has been developed specifically to facilitate vaginal delivery. Dianatal® is bioadhesive and has a high water binding capacity. Dianatal® is used during vaginal delivery to make birth easier for both mother and child and to reduce perineal trauma. Dianatal® forms a low-friction bioadhesive film on the vaginal birth canal, reducing the friction forces between the child and the vaginal birth canal to facilitate vaginal delivery.
Dianatal® obstetric gel is sterile, non-allergenic, electroconductive, slightly acidic and does not cause irritation to mucosal membranes or eyes.
Dianatal® obstetric gel can be used in first-time mothers (randomized control trial (RCT) results) as well as in multiparas and in women with a history of cesarean section (clinical experience, investigator trial results).
Dianatal® obstetric gel eases and shortens both the first and second stages of labour and offers protection to mother and child; in particular protecting the vagina, the pelvic floor and the perineum.
Dianatal® obstetric gel is used specifically:

- to facilitate vaginal delivery and to protect both mother and child by reducing the opposing friction force in first-time mothers (RCT results)
- to facilitate vaginal delivery and to protect both mother and child by reducing the opposing friction force in multiparas, water deliveries and women with history of cesarean (clinical experience, investigator trial results)
- to prevent or treat prolonged or obstructed delivery (clinical experience, investigator trial results)
- to facilitate vaginal operative procedures (vacuum/forceps) (clinical experience)
- to protect the perineum and the vagina during vaginal delivery (RCT results)
- to reduce the vaginal-urethral burn feeling post partum (clinical experience)

Dosage and application

Dianatal® obstetric gel is applied after vaginal examination during delivery after onset of regular contractions, beginning with the vaginal examination around 3-4 cm dilatation of the cervix (active phase of labor). Approx. 3-5 ml birth gel is applied ahead of the presenting part with the tube or manually. In order to achieve an optimum low-friction film, care should be taken to ensure that sufficient fluid is added in dry vaginal conditions. Standard solutions such as NaCl or aseptic solutions can be used for this. An additional application of gel should be made 15-30 minutes after rupture of the amniotic sac.
Dianatal® obstetric gel shortens both the first and second stages of labour. It should therefore be used throughout the entire delivery. The dosage is, on average, between 15 ml and 25 ml Dianatal® obstetric gel per delivery. Experience has shown that the maximum amount required is 30 ml Dianatal® obstetric gel.
The application should only be performed by qualified obstetric professionals.

Method of use:

For local application in the vaginal birth canal only.

Precautions:

Content is sterile as long as the blister pack has not been opened. You should not use Dianatal® obstetric gel if the packaging is damaged, if the packaging has been opened or if the expiry date has been exceeded. Once the sterile blister pack has been opened the tubes must

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be handled in accordance with sterile guidelines to maintain sterility. After delivery of the baby's head, the face area (eyes, mouth and nose area) should be wiped with a dry cloth. If appropriate, superficial cleaning of the mouth/nose area should also be carried out. It is advisable to use a dry towel to deliver the baby to prevent the baby from slipping out of the hands.
Due to the electro conductive properties of Dianatal® obstetric gel, electrosurgery can be employed for haemostasis at any time after the use of Dianatal® obstetric gel; for example when performing a cesarean section or a revision of the birth canal using electrosurgery. Dianatal® obstetric gel should therefore not be used in combination with other gels as this may impair the electrical conductivity of Dianatal® obstetric gel.

Biocompatibility and use during pregnancy:
The biocompatibility of Dianatal® obstetric gel has been tested according to international guidelines. There is no evidence for any harmful effect on the mother or the baby.

Adverse effects
No adverse effects have been reported to date.

Interaction
No interaction events have been reported to date.

Contraindications:
Dianatal® obstetric gel should not be used if:
- there is the possibility of an allergic reaction to any of the ingredients
- vaginal delivery is contra-indicated

Storage and shelf-life:
Dianatal® obstetric gel should be stored at room temperature.
The expiry date is printed on the outer wrapping and packaging. Dianatal® obstetric gel should not be used after the expiry date.

Date of revision of the text: 4.2017

CE-Mark: CE 0481

Medical device

Manufacturer:
Tentan AG, Dellenbodenweg 8, CH-4452 Itingen, Switzerland
Dianatal® is a registered Trademark of Tentan AG, Switzerland

Please refer to the instructions for use Must not be reused

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Dianatal® obstetric gel - key benefit:
reducing trauma & saving time and pain

- Developed to ease childbirth for mother and baby
- Clinically proven for safety, effectiveness, and tolerability for mother & baby (RCT results)
- Reduces duration of labour Stage II (RCT results)
- Reduces incidence of perineal damage and protects the vagina (RCT results)

Dianatal® obstetric gel special properties: safe for mother & baby

- Sterile gel in sterile packaging, single use
- Bioadhesive
- Bacterizide in vitro (B Strept, A Strept., Gardnerella, E. coli)
- Electroconductive
- Non-allergenic, non pyrogenic
- Latex-free, free of phthalate, free of preservatives
- Kind on vagina and on the eyes
- And biocompatible:

Dianatal® obstetric gel was examined in accordance with international and FDA guidelines in regards to its biocompatibility (Toxikon Europe, Interleuvenlaan 3/3, B-3001 Leuven, Belgium): Agar diffusion test, Buehler sensitization test, Embryotoxicity test, Human blood hemolysis test, Intracutaneous injection test, Oral irritation test, Primary ocular irritation test, Primary skin irritation test, Primary vaginal exposure test, Rabbit pyrogen test, Systemic injection test.
Dianatal® obstetric gel was further examined in accordance with international guidelines in regards to its antibacterial activity (Labor L&S AG, Germany). Following bacteria have been tested that can be relevant in obstetrics: Streptococcus pyogenes (Strept. A); Streptococcus agalactiae (Strept. B); Gardnerella vaginalis; Escherichia coli.
All bacteria types incubated in Dianatal® obstetric gel have been eliminated after 24 hours.

Dianatal® obstetric gel is a global brand of Tentan AG, Switzerland

CE 0481
Manufacturer: Tentan AG, Dellenbodenweg 8, CH-4452 Itingen, Switzerland

not for reuse - sterile by aseptic process techniques – do not re-sterilize - do not use if packaging is damaged – consult instructions for use - temperature limitation – store in dry place - store away from sunlight - Latex free - DEHP / PHT free

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