Biosense Webster®
PART OF THE Johnson & Johnson FAMILY OF COMPANIES

CE 2797

Electrophysiology Catheter: Fixed Curve; Fixed Curve Lumen; Fixed Orthogonal

INSTRUCTIONS FOR USE **Electrophysiology Catheter: Fixed Curve; Fixed Curve Lumen; Fixed Orthogonal**	English	Page	8
MODE D'EMPLOI **Cathéter d'électrophysiologie – Courbure fixe ; avec lumière, à courbure fixe ; orthogonal à courbure fixe**	Français	Page	10
GEBRAUCHSANWEISUNG **Elektrophysiologie-Katheter – Fixed Shape; Fixed Shape Lumen; Fixed Orthogonal**	Deutsch	Seite	12
ISTRUZIONI PER L'USO **Catetere per elettrofisiologia – punta fissa; punta fissa con lume; punta fissa ortogonale**	Italiano	Pagina	14
INSTRUCCIONES DE USO **Catéteres de electrofisiología de curva fija; curva fija con lumen; punta fija ortogonal**	Español	Página	16
INSTRUÇÕES DE UTILIZAÇÃO **Cateter Electrofisiológico – Curvatura Fixa; Lúmen de Curvatura Fixa; Ortogonal Fixa**	Português	Página	18
GEBRUIKSAANWIJZING **Elektrofysiologiekatheter – Vaste curve; Lumen met vaste curve; Vast orthogonaal**	Nederlands	Pagina	20
BRUGSVEJLEDNING **Elektrofysiologikateter – Fast bue; Fast bue lumen; Fast ortogonalt**	Dansk	Side	22
KÄYTTÖOHJEET **Kiinteäkaarinen elektrofysiologinen katetri; Kiinteäkaarinen luumenillinen; Kiinteäkaarinen suorakulmainen**	Suomi	Sivu	24
BRUKSANVISNING **Elektrofysiologisk kateter med fixerad Böjning, fixerad böjning och lumen, fixerad ortogonal**	Svensk	Sida	26
BRUKSANVISNING **Elektrofysiologisk kateter Fast kurve; Fast kurve lumen; fast ortogonal**	Norsk	Side	28
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ **Καθετήρας Ηλεκτροφυσιολογίας – Σταθερής καμπύλης, Σταθερής Καμπύλης με Αυλό, Ορθογώνιος Σταθερής Καμπύλης**	Ελληνικά	Σελίδα	30
KULLANMA TALİMATI **Elektrofizyoloji Kateteri: Sabit Kývrýmlý; Sabit Kývrýmlý Lümen; Sabit Ortogonal**	Türkçe	Sayfa	32
ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ **Электрофизиологический катетер: фиксированной кривизны; с просветом фиксированной кривизны; ортогональный**	Русский	Страниц	34
INSTRUKCJA OBSŁUGI **Cewnik elektrofizjologiczny o stałej krzywiźnie; o stałej krzywiźnie z kanałem; stały ortogonalny**	Polski	Strona	36
NÁVOD K POUŽITÍ **Elektrofyziologický katétr – pevné zakřivení; lumen s pevnou křivkou; pevný ortogonální**	Česky	Strana	38
HASZNÁLATI UTASÍTÁSOK **Elektrofiziológiás katéter – rögzített görbület; rögzített görbületű lumen; rögzített ortogonális**	Magyar	Oldal	40

Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru
WARNING: This is a controlled, proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Release Date: 11/1/2019

POKYNY NA POUŽITIE **Elektrofyziologický katéter – pevné zakrivenie; lúmen s pevným zakrivením; pevný ortogonálny**	Slovensky	Strana	42
NAVODILA ZA UPORABO **Elektrofiziološki kateter; fiksna krivina; fiksirana krivina svetline; fiksen ortogonalen**	Slovensko	Stran	44
УКАЗАНИЯ ЗА УПОТРЕБА **Електрофизиологичен катетър: Фиксирана извивка; Фиксирана извивка с лумен; Фиксиран ортогонален**	Български	Страница	46
INSTRUCŢIUNI DE UTILIZARE **Cateter de electrofiziologie: Curbură fixă; Lumen de curbură fixă; Ortogonal fix**	Română	Pagina	48
KASUTUSJUHEND **Elektrofüsioloogiline kateeter – fikseeritud kaarega; fikseeritud kaare valendikuga; fikseeritud ortogonaalne**	Eesti	Lehekülg	50
LIETOŠANAS NORĀDĪJUMI **Elektrofizioloģijas katetrs – fiksēts izliekums; fiksēta izliekuma lūmens; fiksēts ortogonāls**	Latviski	Lappuse	52
NAUDOJIMO INSTRUKCIJOS **Elektrofiziologinis kateteris – pastovus linkis; stabilus lenktasis liumeninis; stabilus statmenasis**	Lietuvių k.	Puslapis	54
使用说明 **电生理导管：固定弯曲、固定弯曲内腔、固定正交**	简体中文	页码	56
사용 설명서 **전기생리 카테터: 고정 커브; 고정 커브 관내강; 고정 직각**	한국어	페이지	57
使用說明 **電生理導管： 固定彎曲； 固定彎曲內腔；固定直角**	繁體中文	頁碼	59
INSTRUKSI PENGGUNAAN **Kateter Elektrofisiologi: Kurva Tetap; Lumen Kurva Tetap; Ortogonal Tetap**	Bahasa Indonesia	Halaman	60
UPUTE ZA UPORABU **Elektrofiziološki kateter: fiksna krivina; fiksna krivina lumena; fiksan ortogonalan**	Hrvatski	Strana	62
УПУТСТВО ЗА УПОТРЕБУ **Електрофизиолошки катетер: Фиксирана крива; лумен са фиксираном кривом; фиксирани ортогонални**	Српски	Страна	64
УПАТСТВО ЗА УПОТРЕБА **Електрофизиолошки катетер: Фиксна крива; Лумен со фиксна крива; Фиксен ортогонален**	Македонски јазик	Страница	66

Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru
WARNING: This is a controlled, proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Release Date: 11/1/2019

Symbol Definitions / Définitions des symboles / Symboldefinitionen / Definizioni dei simboli / Definiciones de símbolos / Definições dos símbolos / Definitie van symbolen / Symboldefinitioner / Symbolien selitykset / Symbolförklaringar / Symboldefinisjoner / Ορισμοί συμβόλων / Sembol Tanımları / Определение символов / Definicje symboli / Definice symbolů / A szimbólumok magyarázata / Definicije simbolov / Definície symbolov / Определения на символите / Definițiile simbolurilor / Sümbolite tähendused / Simbolu definīcijas / Simbolių apibrėžtys / 标识定义 / 기호 정의 / 符號定義 / Definisi Simbol / Opis simbola / Дефиниције симбола / Дефиниции за символите

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
Electrophysiology Catheter: Fixed Curve; Fixed Curve Lumen; Fixed Orthogonal

- **STERILE. Sterilized using ethylene oxide.**
- **For single use only.**
- **Do not use if the package is open or damaged.**
- **Federal (USA) law restricts this device to sale by or on the order of a physician.**

CATHETER DESCRIPTION
The BIOSENSE WEBSTER® Fixed Curve Electrophysiology (EP) Catheter has been designed to facilitate electrophysiological mapping of cardiac structures. The catheter has a high-torque shaft with an array of platinum electrodes at the distal tip that can be used for stimulation and recording. The tip section has a preformed curve with a nylon tip protector. (Standard curve types are identified by a color band on the shaft near the proximal end - see catalog).

The Fixed Curve Lumen Catheter has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque through-lumen shaft with a fixed curve section containing platinum electrodes at the distal tip that can be used for stimulation and recording. The tip section has a preformed curve and is packaged with a nylon tip protector. The Fixed Curve Lumen Catheter terminates with a luer hub that may be used with a guidewire or for injection purposes.

The Fixed Orthogonal Catheter has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with an array of platinum plaque electrodes at the distal tip that can be used for stimulation and recording. The tip section has a preformed curve and is packaged with a nylon tip protector.

The high-torque shaft allows the plane of the distal curve to be manually rotated to assist in positioning the catheter tip at the desired site.

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INDICATIONS
The catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. In addition, the Fixed Orthogonal Catheter is designed for use in the coronary sinus.

CONTRAINDICATIONS
- The BIOSENSE WEBSTER® Fixed Curve EP Catheter has not been shown to be safe and effective for electrical ablation or for use in the coronary vasculature other than the coronary sinus ostium.
- Use of the catheter may not be appropriate for patients with prosthetic valves.
- A relative contraindication for cardiac catheter procedures is active systemic infection.

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PRECAUTIONS
- Do not attempt to operate the BIOSENSE WEBSTER® Fixed Curve EP Catheter prior to completely reading and understanding these directions for use.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Store in a cool, dry, dark place.
- The sterile packaging and catheter should be inspected prior to use.
- BIOSENSE WEBSTER® EP Catheters are intended for single patient use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

“USE BY” DATE
Use the device prior to the “Use By” date on the package label.

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SUGGESTED DIRECTIONS FOR USE
- Remove the catheter from its package and place it in a sterile work area. Remove nylon tip protector.
- Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the lead pins or connectors to the appropriate recording equipment.
- Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. The electrode catheter shaft may be gently hand rotated to facilitate positioning at the desired mapping site.
- Remove the catheter and dispose it in an appropriate manner. Do not resterilize and reuse.

IF THERE ARE ANY QUESTIONS REGARDING THE USE OR PERFORMANCE OF THIS PRODUCT, PLEASE CONSULT WITH THE LOCAL DISTRIBUTOR OR THE MANUFACTURER.

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Descriptions and specifications appearing in Biosense Webster Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product at the time of manufacture and are not made or given as a warranty of the prescribed product in any way.

Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru
WARNING: This is a controlled, proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Release Date: 11/1/2019

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- **Sterile. Sterilized using ethylene oxide.**
- **For single use only.**
- **Do not use if the package is open or damaged.**
- **Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.**

CATHETER DESCRIPTION
The Biosense Webster® Deflectable Tip Electrode Catheter has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing an array of platinum electrodes that can be used for stimulation and recording. Several tip configurations are available. (Standard curve types are identified by a color band on the shaft near the proximal end - see catalog.)

The Deflectable Tip Orthogonal Catheter has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing platinum plaque electrodes that can be used for stimulation and recording.

The Deflectable Braided Tip Catheter is a multi-electrode catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing platinum electrodes. The catheter is designed to facilitate stable local electrograms due to the internal braid. All of the electrodes may be used for recording and stimulation purposes.

The Deflectable Tip CRISTACATH™ and ISMUS® CATH catheters have been designed for electrophysiological mapping of cardiac structures. The catheters have a high-torque shaft with a tip section containing ten pairs of platinum electrodes that can easily be seen under fluoroscopy. The Deflectable Tip CRISTACATH™ and ISMUS® CATH catheters facilitate simultaneous local electrograms due to the greater number of electrodes and their deflection capabilities.

The Deflectable Tip HALO® XP and HALO® catheters have been designed for electrophysiological mapping of the heart, in particular, the tricuspid annulus. The catheters have a hightorque shaft with a halo-shaped tip section containing ten pairs of platinum electrodes that can easily be seen under fluoroscopy. The tip section also contains a radiopaque mark in the center of the electrode array. The halo-shaped tip section has a preformed curve which can be positioned around the atrial aspect of the tricuspid annulus. The tip section is packaged with a plastic tip protector.

A piston in the handpiece is attached to an internal puller wire which changes the radius of curvature. When the piston is pushed forward, the radius of curvature of the preformed loop is reduced; when the thumbknob is pulled back, the radius of the curvature is increased until the tip section returns to the preformed shape. The high torque shaft allows the plane of the curve to be maneuvered in order to facilitate accurate positioning.

The Deflectable Tip HALO® XP and HALO® catheters facilitate simultaneous local electrograms spanning the tricuspid annulus, from the midseptal to anterior to lateral to posterolateral. Recordings of the entire annulus can be obtained without repositioning the catheter tip.

Tip deflection is controlled at the proximal end by a tubular handpiece in which a piston slides. When the piston is pushed forward with the thumbknob, the tip is deflected (curved). When the piston is pulled back, the tip straightens. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

The catheter interfaces with standard recording equipment via interface cables with the appropriate connectors.

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INDICATIONS
The catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. In addition, the Deflectable Tip Orthogonal Catheter is designed for use in the coronary sinus.In addition the Deflectable Braided Tip Catheter is designed for use in the right side of the heart. In addition, the CRISTACATH™ and ISMUS® CATH catheters are designed to facilitate electrograms in the atrial region of the heart. In addition, the HALO® XP and HALO® catheters' preformed halo shape of the tip section is designed specifically for the tricuspid annulus.

CONTRAINDICATIONS
The Biosense Webster® Deflectable Tip Electrode Catheter has not been shown to be safe and effective for electrical ablation or for use in the coronary vasculature other than the coronary sinus ostium. In addition, the Braided Tip Catheter is contraindicated for use in the left side of the heart. In addition, the transseptal approach is contraindicated.

Use of the catheter may not be appropriate for patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.

WARNINGS
Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Cardiac catheterization should only be performed after adequate attention has been given to potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must therefore be given for the use of this catheter in pregnant women.

- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the catheter to organic solvents such as alcohol.
- Do not autoclave the catheter.

PRECAUTIONS
- Do not attempt to operate the Biosense Webster deflectable tip electrophysiology catheter prior to completely reading and understanding these directions for use.
- Store in a cool, dry, dark place.
- The sterile packaging and catheter should be inspected prior to use.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.

"USE BY" DATE
Use the device prior to the "Use By" date on the package label.

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SUGGESTED INSTRUCTIONS FOR USE
- Remove the catheter from its package and place it in a sterile work area. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the interface connectors to the appropriate recording equipment.
- Confirm that the thumbknob is pulled back completely before insertion. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter and dispose of it in an appropriate manner. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

ADVERSE REACTIONS
A number of serious adverse reactions have been documented for cardiac catheterization procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death. The following complications associated with cardiac catheterization have also been reported in the literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, cardiac perforation, air embolism, arrhythmias, valvular damage, pneumothorax and hemothorax.

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THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.

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Descriptions and specifications appearing in Biosense Webster Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru
WARNING: This is a controlled, proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Released Date: 08/11/2017

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Electrophysiology Catheter: Deflectable Tip; Deflectable Tip Orthogonal; Deflectable Braided Tip; Deflectable Tip HALO® XP Tricuspid Mapping, HALO® & Deflectable Tip CRISTACATH™ and ISMUS® CATH Mapping Catheters protected under one or more of the following U.S. Patents: 5,827,278; 7,371,232 and other patents pending in the U.S. and other countries.

BIOSENSE WEBSTER, the Biosense Webster logo, CRISTACATH, ISMUS, and HALO are trademarks of Biosense Webster, Inc.

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CE 0086

HALO® Catheter with Auto ID Technology

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HALO® Catheter with Auto ID Technology

Caution: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

- **Sterile. Sterilized with ethylene oxide gas.**
- **For single use only.**
- **Do not use if the package is open or damaged.**

CATHETER DESCRIPTION
The Deflectable Tip HALO® Catheters have been designed for electrophysiological mapping of the heart, in particular, the tricuspid annulus. The catheters have a high torque shaft with a HALO-shaped tip section containing ten pairs of platinum electrodes that can easily be seen under fluoroscopy. The tip section also contains a radiopaque marker in the center of the electrode array. The HALO-shaped tip section has a preformed curve which can be positioned around the atrial aspect of the tricuspid annulus. The tip section is packaged with a plastic tip protector. The catheter is equipped with Electronically Erasable Programable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO™ EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information.

A piston in the handpiece is attached to an internal puller wire which changes the radius of the curvature. When the piston is pushed forward with the thumbknob, the radius of curvature of the preformed loop is reduced; when the thumbknob is pulled back, the radius of the curvature is increased until the tip section returns to the preformed shape. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

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INDICATIONS
The HALO® Catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. The preformed HALO shape of the tip section is designed specifically for the tricuspid annulus.

CONTRAINDICATIONS
The Deflectable Tip HALO® Catheter has not been shown to be safe and effective for electrical ablation or for use in the coronary vasculature other than the coronary sinus ostium. The transseptal approach is contraindicated. Use of the catheter may not be appropriate for patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.

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PRECAUTIONS
- Do not attempt to operate the HALO® Catheter prior to completely reading and understanding these directions for use.
- Store in a cool, dry place.
- The sterile packaging and catheter should be inspected prior to use.
- Do not resterilize and reuse.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.

This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

“USE BY” DATE
Use the device prior to the "Use By" date on the package label.

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SUGGESTED INSTRUCTIONS FOR USE
- Remove the catheter from its package and place it in a sterile work area. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the interface cable to the Patient Interface Unit of the appropriate CARTO™ EP Navigation System and connect the catheter to the interface cable.
- Confirm that the thumbknob is pulled back completely before insertion. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter and dispose of it in an appropriate manner. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

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Interface Cable

Federal (USA) law restricts this device to sale by or on the order of a physician.

INDICATIONS AND USE
This cable provides a means to interface a BIOSENSE WEBSTER® catheter to the appropriate equipment. The cable may also be used as an extension cable to facilitate catheter connection to equipment that is located out of the immediate area of use. This cable may be reused subject to the sterilization restrictions herein.

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PRODUCT CONDITION AND STERILIZATION
This cable may be provided sterile or non-sterile. If non-sterile, it may be sterilized using a standard ethylene oxide (EtO) gas process. Do not use either radiation or steam sterilization as these will damage the cable.

The cable may be re-sterilized for reuse. The cable and connector surfaces should be wiped clean with soap and water. The connector should not be immersed in fluids. The recommended maximum number of sterilization cycles is ten (10)*. Lead-to-lead leakage and end-to-end resistance should be checked after each use. Leakage should not exceed 1.0 microamp at 250 volts and resistance should not exceed 10 ohms. Exceeding ten uses may impair the proper functioning of the cable.

* Sterilization test data on file.

CAUTIONS
- Do not immerse the cable connector(s) in fluids.
- Patient or operator injury can result from improper handling of the cable.
- If used in the presence of electrical equipment, noise could be induced into the cable.
- The cable should be tested for leakage (lead-to-lead) and resistance (end-to-end) before reuse.
- Do not expose cable to organic solvents.

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NOTICE: This information applies only to the BIOSENSE WEBSTER® interface cable used with the BIOSENSE WEBSTER® Diagnostic/Ablation temperature sensing (thermocouple) catheter and a compatible RF generator (refer to the specifications below).

1. Use this interface cable to connect the BIOSENSE WEBSTER® Diagnostic/ Ablation Deflectable Tip catheter to a compatible RF generator.
2. Refer to the RF generator manual for detailed generator operating instructions for RF catheter ablation.

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PARAMETER	SPECIFICATION
Thermometry	Thermocouple
Temperature Limit, maximum	100°C
Modes: (must operate in all 3 modes)	1. Temperature Control 2. Temperature Monitoring 3. Power Control
Maximum Output Power	50 Watts
RF output frequency	450 kHz - 550 kHz
Impedance cut-off	high: 250 Ω low: 40 Ω

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.

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INDICATIONS FOR USE
This cable provides a means to interface a Biosense Webster electrophysiology catheter or sheath to the appropriate equipment. This cable may be reused subject to the cleaning and sterilization restrictions herein.

CONTRAINDICATIONS
There are no known contraindications to this cable.

PRODUCT CONDITION AND STERILIZATION
This cable is provided sterile and may be resterilized for up to 20 uses per the instructions provided.

PRODUCT USE
Ensure each cable end is connected to the proper equipment. The “System” end of the cable is labeled “To System” and must be connected to the System. The “Catheter” end of the cable is labeled with the cable catalog number and must be connected to the Catheter or Sheath. For information on how to connect this cable, refer to the applicable catheter/sheath and equipment instructions for use.

The cable should be tested for leakage (lead to lead) and resistance (end-to-end) before each reuse (leakage should not exceed 1.0 µA at 250 V and resistance should not exceed 10 Ω).

CAUTIONS
- Do not immerse the cable connectors in fluids.
- Do not reuse the cable after twenty (20) uses.
- If used in the presence of electrical equipment, noise may be induced into signals conducted by the cable.
- The cable should be visually inspected for any damage (i.e. cuts, kinks, nicks, crushed or elongated sections) prior to use or reuse.
- Do not attempt to repair any damage. In case of doubt, discard the cable and do not use.
- Do not expose cable to organic solvents.
- Patient or operator injury can result from improper handling or storage of the cable.To prevent damage, store the cable in a clean and secure area.

WARNINGS
- If a break occurs in the cable wire or the cable becomes otherwise electrically discontinuous, arcing may occur in the patient-return or active circuit and may burn the patient or create a fire.

RESTERILIZATION
At the point of use, remove excess contaminants with a disposable cloth or paper wipe. Cables should be resterilized as soon as is practical after use.

It is recommended that the cables are packaged in a breathable pouch prior to resterilization. Cables should be coiled into a loop when placing into the pouch. Use a pouch large enough to prevent stress to the pouch seams or excessive bending of the cable.

The cable may be resterilized using the methods specified in Table 1. Sterilization cycles should be performed using techniques and parameters compliant with the equipment manufacturer’s instructions and facility protocols.

The resterilization parameters provided in Table 2 have been validated by the manufacturer. Tables 3 and 4 provide common re-sterilization parameters for reference and have not been validated by the manufacturer. For more information on common practices and validation, refer to ANSI/AAMI TIR12:2010[1,6] and ISO 14937:2009[5].

NOTES
- This cable meets the requirements of ANSI/AAMI EC53:1998[2] and ANSI/AAMI EC53A:1998/(R) 2008[2].
- Recommended parameters are validated for device functionality. Efficacy of cleaning and sterilization is the responsibility of the user.
- Methods described here are for reference only and no guarantee is made that effective sterilization will result.
- Chamber dimensions, temperature distribution profile and the configuration of products placed in the chamber are variables that must be considered in verifying the effectiveness of a specific sterilization cycle.

CAUTIONS
- Radiation resterilization must not be used, as this method will damage the cable.
- Exceeding twenty (20) resterilization cycles may impair the proper functioning of the cable. The cable should be discarded after twenty uses.

WARNINGS
- Failure to properly clean cables could lead to inadequate resterilization.
- Do not immerse the ends of the cable in liquid.
- Automated cleaning of the cable is not recommended, and may damage the cable.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Release Date: 4/27/2018
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

Table 1: Validated Resterilization Methods per Cable Model

Catalog Number	Twenty (x20) EtO Resterilization Cycles	Twenty (x20) Steam Resterilization Cycles
CB3434CT	X	X
CR3434CT	X	X
CR3425CT	X	X
CB3410CT	X	X
CB3412CT	X	X
CY1210CT	X	X
CY1212CT	X	X
CG2025CT	X	X
D128709		X
D128622		X
D128623	X	X
D128624		X
D128627	X	X
D128635		X

Table 2: Validated Resterilization Parameters

	Exposure Temperature	Exposure Time	Gas Concentration	Humidity(RH)	Accommodation/ Aeration time
Gravity Displacement Steam Resterilization	135°C ± 1°C	10 minutes	N/A	N/A	3hr
	123°C ± 1°C	30 minutes	N/A	N/A	24hr
100% EtO Gas Resterilization	52-57°C	120-150 minutes	600mg/L	40-70%	24hr

The cycle parameters above were used to evaluate product functionality after resterilization per Table 1.
Cycle parameters falling under these conditions are also considered validated for product functionality.
Sterilization efficacy using these parameters is the responsibility of the end user as specified under ANSI/AAMI ST24:1999/R 2009[3,6], ANSI/AAMI ST79:2010[4,6], and ISO 14937:2009[5].

NON-VALIDATED REFERENCE INFORMATION

The information and resterilization parameters provided in tables 3 and 4 have not been evaluated by Biosense Webster, and are provided for reference only. Biosense Webster cannot guarantee the product functionality after resterilization. For further information and validation procedures, please see applicable standards referenced.

*FOR REFERENCE ONLY: THE RESTERILIZATION PARAMETERS IN TABLES 3 AND 4 HAVE NOT BEEN VERIFIED/VALIDATED BY BIOSENSE WEBSTER, INC.

Table 3: Common Steam Resterilization Parameters Provided for Reference*

	Exposure Temperature	Exposure Time	Drying Time
Gravity Steam Resterilization	121°C	30 minutes	15-30 minutes
	132°C	15 minutes	15-30 minutes
	135°C	10 minutes	30 minutes
Dynamic Air Removal (Pre-Vac) Steam Resterilization	132°C	4 minutes	20 minutes
	135°C	3 minutes	16 minutes

Cycle parameters are provided for reference, as described in ANSI/AAMI TIR 12:2010[1,6].
Sterilization efficacy using these parameters is the responsibility of the end user as specified under ANSI/AAMI ST79:2010[4,6] and ISO 14937:2009[5].

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Release Date: 4/27/2018
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

FOR REFERENCE ONLY

Table 4: Common EtO Resterilization Parameters Provided for Reference*

	Exposure Concentration	Exposure Time	Relative Humidity
100% EtO Gas Resterilization	880 mg/L	1 hr	50% to 80%
	740 mg/L	1 hr	80%
	735 mg/L	1 hr	50% to 80%
	725 mg/L	1 hr	30% to 80%

Cycle parameters are provided for reference, as described in ANSI/AAMI TIR 12:2010[1,6].
Sterilization efficacy using these parameters are the responsibility of the end user as specified under ANSI/AAMI ST24:1999/R 2009[3,6] and ISO 14937:2009[5].

*FOR REFERENCE ONLY: THE RESTERILIZATION PARAMETERS SPECIFIED IN TABLES 3 AND 4 HAVE NOT BEEN VERIFIED/VALIDATED BY BIOSENSE WEBSTER, INC.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]This is an American Standard or Technical Report.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Norme ou rapport technique américain.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Dies ist eine US-amerikanische Norm oder ein US-amerikanischer technischer Bericht.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Release Date: 4/27/2018
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Questa è una normativa americana o relazione tecnica.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Este es un estándar norteamericano o Informe técnico.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Trata-se de uma norma americana ou relatório técnico.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Dit is eenAmerikaanse standaard of technisch rapport.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Dette er en amerikansk standard eller teknisk rapport.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Tämä on amerikkalainen standardi tai tekninen raportti.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Det här är en amerikansk standard eller teknisk rapport.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 og ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 og A2:2011 (samlet tekst).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6] Dette er en amerikansk standard eller teknisk rapport.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Πρόκειται για αμερικάνικο πρότυπο ή τεχνική έκθεση.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Release Date: 4/27/2018
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Bu, bir Amerikan Standart veya Teknik Raporudur.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6] Это стандарт США или технический отчет.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Release Date: 4/27/2018
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Jest to amerykańska norma lub raport techniczny.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Release Date: 4/27/2018
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Toto je americká norma nebo technická zpráva.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Ez egy amerikai szabvány vagy egy műszaki jelentés.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Ide o americkú normu alebo technickú správu.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]To je ameriški standard ali tehnično poročilo.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Това е американски стандарт или технически доклад.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Release Date: 4/27/2018
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Aceasta este o normă americană sau un raport tehnic american.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]See on Ameerika standardi või tehniline aruanne.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Šis ir Amerikas standarts vai tehniskais ziņojums.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Tai JAV standartas arba techninė ataskaita.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]这是美国标准或技术报告。

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6] 이는 미국 표준 또는 기술 보고서입니다.

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6] 這是一種美國標準或一份技術報告。

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[1]Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers, 3ed. AAMI Technical Information Report 12:2010.
[2]ECG cables and leadwires, 1ed and 1998 Amendment. ANSI/AAMI EC53:1995/(R)2008 and ANSI/AAMI EC53A:1998/(R)2008.
[3]Automatic, general-purpose ethylene oxide sterilizers and ethylene oxide sterilant sources intended for use in health care facilities, 3ed. ANSI/AAMI ST24:1999/(R)2009.
[4]Comprehensive guide to steam sterilization and sterility assurance in health care facilities, 2ed and 2010/2011 amendments. ANSI/AAMI ST79:2010, A1:2010 and A2:2011 (consolidated text).
[5]Sterilization of health care products -- General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices. ISO 14937:2009
[6]Ini adalah Standar Amerika atau Laporan Teknis.

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Pertinta data MD is Protected by Copyright
www.neveacert.info@neveacert.eu
Not for Commercial Use

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Released Date: 08/11/17

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ISMUS® CATH Mapping Catheter with Auto ID Technology

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CATHETER DESCRIPTION
The Deflectable Tip ISMUS® CATH Catheters have been designed for electrophysiological mapping of cardiac structures. The catheters have a high-torque shaft with a tip section containing ten pairs of platinum electrodes that can easily be seen under fluoroscopy. The tip section also contains a radiopaque marker in the center of the electrode array. The Deflectable Tip ISMUS® CATH Catheters facilitate simultaneous local electrograms due to the greater number of electrodes and their deflection capabilities. The catheter is equipped with Electronically Erasable Programmable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO™ EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information.

A piston in the handpiece is attached to an internal puller wire which changes the radius of the curvature. When the piston is pushed forward with the thumbknob, the radius of the curvature is reduced; when the thumbknob is pulled back, the radius of the curvature is increased until the tip section straightens. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

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INDICATIONS
The catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. The ISMUS® CATH Catheters are designed to facilitate electrograms in the atrial regions of the heart.

CONTRAINDICATIONS
The Deflectable Tip ISMUS® CATH Electrode Catheter has not been shown to be safe and effective for electrical ablation. Use of the catheter may not be appropriate for patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.

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PRECAUTIONS

- Do not attempt to operate the ISMUS® CATH Catheter prior to completely reading and understanding these directions for use.
- Store in a cool, dry place.
- The sterile packaging and catheter should be inspected prior to use.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.

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SUGGESTED INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the interface cable to the Patient Interface Unit of the appropriate CARTO™ EP Navigation System and connect the catheter to the interface cable.
- Confirm that the thumbknob is pulled back completely before insertion. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter and dispose of it in an ner. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

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WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. AND ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

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WARNING: To reduce the risk of entrapping cardiac structures in the mapping electrode portion of the catheter, place the LASSO™ 2515 NAV Variable Catheter by torquing (or rotating) the shaft in a clockwise motion only.

WARNING: When not in regions intended for mapping, manipulate the catheter with the loop in the fully expanded (i.e. 25 mm diameter) position to further reduce the risk of entrapping cardiac structures.


- **STERILE. Sterilized using ethylene oxide.**
- **For single use only.**
- **Do not autoclave.**
- **Do not use if the package is open or damaged.**

CATHETER DESCRIPTION

The Biosense Webster LASSO™ 2515 NAV Variable Catheter has been designed to facilitate electrophysiological mapping of the atria of the heart with the CARTO™ 3 EP Navigation System and a reference device. It is deployed in the right or left atrium through an 8 F guiding sheath. This deflectable catheter consists of a 4 F circular spine on its distal tip, with platinum electrodes that can be used for stimulation and recording.

The LASSO™ 2515 NAV Variable Catheter features a Nitinol loop design that allows the expansion and contraction of the loop to custom-fit veins with different sizes, ranging from 25 mm to 15 mm diameter (± 15%).

With the catheter pointing away from you, rotate the handle clockwise to reduce the diameter of the loop. Rotating the handle counter-clockwise increases the loop diameter (see Figure 1). When the handle is fully rotated clockwise, the minimum 15 mm diameter is attained. Likewise, when the handle is fully rotated counter-clockwise, the maximum diameter of 25 mm is attained. Pushing forward on the catheter thumbknob deflects the tip of the catheter. When the thumbknob is pulled back, the catheter tip straightens.

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1: LASSO™ 2515 NAV Variable Catheter

The catheter interfaces with standard recording equipment and CARTO™ 3 EP Navigation System equipment via interface cables with the appropriate connectors. Consult the local distributor or manufacturer for the appropriate interface cables.

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INDICATIONS

The LASSO™ 2515 NAV Variable Catheter is indicated for multiple electrode electrophysiological mapping of the cardiac structures of the heart, i.e. recording or stimulation only. The Catheter is designed to obtain electrograms in the atrial regions of the heart.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Released Date: 08/11/2017

The LASSO™ 2515 NAV Variable Catheter provides location information when used with compatible CARTO™ EP Navigation Systems version 1.0 or higher.

CONTRAINDICTIONS

- The LASSO™ 2515 NAV Variable Catheter has not been shown to be safe and effective for radio frequency (RF) ablation.
- Use of the catheter may not be appropriate for use in patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.
- The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.
- The retrograde approach is contraindicated because of risk of entrapping the LASSO™ 2515 NAV Variable Catheter in the left ventricle or valvular apparatus. The LASSO™ 2515 NAV Variable Catheter is not recommended for use in the ventricles.

WARNINGS

1. Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury, as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
2. Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
3. Careful consideration must be given for the use of this catheter in pregnant women.
4. Do not immerse the proximal handle, pig tail or cable connector in fluids; electrical performance could be affected.
5. Do not expose the catheter to organic solvents such as alcohol.
6. Do not autoclave the catheter.
7. To avoid potential damage to anatomical structures, do not attempt to pull the catheter, or withdraw it into the sheath, with the loop in a contracted position. The loop should be fully relaxed (handle grip rotated fully to the left) to minimize tension applied to the Nitinol structure.
8. Do not introduce the LASSO™ 2515 NAV Variable Catheter tip folded into the guiding sheath.
9. The LASSO™ 2515 NAV Variable Catheter is recommended for use with an 8 F guiding sheath. Note: Do not use the LASSO™ 2515 NAV Variable Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25 mm in diameter.

PRECAUTIONS

1. The standard transeptal procedure should be followed during mapping when moving the catheter from the right atrium to the left atrium.
2. To prevent entanglement of the catheter with the valves and to prevent slippage of the catheter into the ventricles, care should be taken when using the catheter in or around the atrio-ventricular valve region.
3. Prior to mapping, the physician should ensure that the intracardiac signals are recorded by all of the electrodes on the loop of the catheter.
4. Do not use in the proximity of magnetic resonance imaging (MRI) equipment because the MRI equipment may induce movement of the catheter, resulting in perforation.
5. The catheter should be stored in its original packaging and in a cool, dry, dark place until it is used.
6. Read all accessory operating instructions prior to connection of the LASSO™ 2515 NAV Variable Catheter.
7. Do not attempt to operate the LASSO™ 2515 NAV Variable Catheter prior to completely reading and understanding these Instructions for Use.
8. Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
9. Careful manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance through a guiding sheath. Do not use excessive force to advance or withdraw the catheter through the guiding sheath, when resistance is encountered. In addition, extra care should be taken while inserting, aspirating and manipulating the guiding sheath.
10. The sterile packaging and catheter should be inspected prior to use.
11. The LASSO™ 2515 NAV Variable Catheter is intended for single patient use only.
12. This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
13. Do not resterilize and reuse.
14. Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.
15. Do not apply RF energy when the ablation catheter is in contact with one or more of the LASSO™ 2515 NAV Variable Catheter electrodes.

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“USE BY” DATE

- This catheter is sterilized with ethylene oxide gas.
- Do not use if the package is open or damaged.
- Use the device prior to the “Use By” date on the package label.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Released Date: 08/11/2017

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SUGGESTED INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area.
- Follow standard practice for vessel puncture, guidewire insertion and guiding sheath use and aspiration per its Instructions for Use.
- Connect the interface connectors to the appropriate recording equipment.
- **NOTE: Read pacing and recording equipment operator manual for proper set up and operation.**
- Confirm that the thumbknob is pulled back completely before insertion and that the loop-contraction mechanism is not activated, ensuring minimal tension to the Nitinol loop.
- Adjust the loop diameter with the handle grip. By rotating the handle grip to the right, the loop contracts; by rotating the grip to the left, the loop expands (see Figure 1).
- Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely and that the loop is in a fully relaxed position (handle grip rotated fully to the left). Remove the catheter through the guiding sheath and dispose of it in an appropriate manner. Remove the guiding sheath, vessel dilator and guidewire as a unit per its Instructions for Use. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

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The LASSO™ 2515 NAV Variable Catheter and accessories are protected under one or more of the following U.S. Patents: 5,391,199; 5,443,489; 5,480,422; 5,546,951; 5,558,091; 5,568,809; 5,626,136; 5,713,946; 5,718,241; 5,738,096; 6,484,118; 6,987,995; 7,142,903; 7,274,957; 7,593,760; 7,853,302; 7,720,517; 8,000,765; 8,275,440 and other patents pending in the U.S. and other countries.

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BIOSENSE WEBSTER, the Biosense Webster logo, CARTO, and LASSO are trademarks of Biosense Webster, Inc.

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Rev.: 2017-07

M-5276-562J © Biosense Webster, Inc. 2008-2017 Printed in USA

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Released Date: 08/11/2017

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CE 0086

LASSO™ Deflectable Circular Mapping Catheter with Auto ID Technology

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LASSO™ Deflectable Circular Mapping Catheter with Auto ID Technology

- **STERILE. Sterilized using ethylene oxide.**
- **For one single use only.**
- **Do not autoclave.**
- **Do not use if the package is open or damaged.**

CATHETER DESCRIPTION
The Biosense Webster LASSO™ Deflectable Circular Mapping Catheter has been designed to facilitate electrophysiological mapping of the atria of the heart. It is deployed in the right or left atrium through an 8F guiding sheath, similar to Biosense Webster's PREFACE™ Guiding Sheaths. This deflectable catheter consists of a 3F or 4F circular spine on its distal tip, with platinum electrodes that can be used for stimulation and recording. Models are available with different circular spine diameters. By pushing forward on the catheter thumbknob, the tip is deflected; when the thumbknob is pulled back, the catheter tip straightens. The catheter is equipped with Electronically Erasable Programable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO™ EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information.

The catheter interfaces with CARTO™ EP Navigation Systems equipped with Auto ID Technology via interface cables with the appropriate connectors.

CONSULT THE LOCAL DISTRIBUTOR OR MANUFACTURER FOR THE APPROPRIATE INTERFACE CABLES. REFER TO THE APPROPRIATE CARTO™ EP NAVIGATION SYSTEM USER MANUAL FOR DETAILS ON SYSTEM SETTINGS RELATED TO THE CATHETER.

INDICATIONS
The catheter is indicated for multiple electrode electrophysiological mapping the cardiac structures of the heart, i.e. recording or stimulation only. The Biosense Webster LASSO™ Circular Mapping Catheter is designed to obtain electrograms in the atrial region of the heart.

CONTRAINDICATIONS
- The Biosense Webster LASSO™ Circular Mapping Catheter has not been shown to be safe and effective for radio frequency (RF) ablation.
- Use of the catheter may not be appropriate for use in patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.
- The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.
- The retrograde approach is contraindicated because of risk of entrapping the LASSO™ Catheter in the left ventricle or valvular apparatus. The LASSO™ Catheter is not recommended for use in the ventricles.

WARNINGS
- Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
- Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
- Careful consideration must be given for the use of this catheter in pregnant women.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the catheter to organic solvents such as alcohol.
- Do not autoclave the catheter.
- Do not introduce the LASSO™ Catheter tip folded into the guiding sheath.
- The LASSO™ Deflectable Circular Mapping Catheter is recommended for use with the Biosense Webster PREFACE™ Braided Guiding Sheath.

 Note: Do not use the LASSO™ Deflectable Circular Mapping Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25mm in diameter.

PRECAUTIONS
- The catheter should be stored in its original packaging and in a cool, dry place until it is used.
- Do not attempt to operate the Biosense Webster LASSO™ Circular Mapping Catheter prior to completely reading and understanding these Instructions for Use.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance through a guiding sheath. Do not use excessive force to advance or withdraw the catheter through the guiding sheath when resistance is encountered. In addition, extra care should be taken while inserting, aspirating and manipulating the guiding sheath.
- The sterile packaging and catheter should be inspected prior to use.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.
- To place LASSO™ Catheter, torque (or rotate) shaft in a clockwise motion only.

“USE BY” DATE
- Use the device prior to the "Use By" date on the package label.

DISPOSAL
- Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

SUGGESTED INSTRUCTIONS FOR USE
- Remove the catheter from its package and place it in a sterile work area.
- Follow standard practice for vessel puncture, guidewire insertion and guiding sheath use and aspiration per its Instructions for Use.
- Connect the interface cable to the Patient Interface Unit of the appropriate CARTO™ EP Navigation System and connect the catheter to the interface cable.
- Confirm that the thumbknob is pulled back completely before insertion. Advance the catheter through the guiding sheath to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
- To place LASSO™ Catheter, (or rotate) shaft in a clockwise motion only.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter through the guiding sheath and dispose of it in an appropriate manner. Remove the guiding sheath, vessel dilator and guidewire as a unit per its Instructions for Use. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru
WARNING: This is a confidential document. Verify revision is current prior to use.
AVISO: Este es un documento confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.

ADVERSE REACTIONS
A number of serious adverse reactions have been documented for cardiac catheterization procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death.

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DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY
THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.

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Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

WARNING: This is a controlled confidential document. Verify revision is current prior to use.
AVISO: Este es un documento confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

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LASSO® Deflectable Circular Mapping Catheter

- **STERILE. Sterilized using ethylene oxide.**
- **For one single use only.**
- **Do not autoclave.**
- **Do not use if the package is open or damaged.**
- **CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.**

CATHETER DESCRIPTION
The Biosense Webster LASSO® Deflectable Circular Mapping Catheter has been designed to facilitate electrophysiological mapping of the atria of the heart. It is deployed in the right or left atrium through an 8 F guiding sheath, similar to Biosense Webster 's PREFACE® Guiding Sheaths. This deflectable catheter consists of a 3 F circular spine on its distal tip, with platinum electrodes that can be used for stimulation and recording. Models are available with different circular spine diameters. By pushing forward on the catheter thumbknob, the tip is deflected; when the thumbknob is pulled back, the catheter tip straightens.

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INDICATIONS
The catheter is indicated for multiple electrode electrophysiological mapping the cardiac structures of the heart, i.e. recording or stimulation only. The Biosense Webster LASSO® Circular Mapping Catheter is designed to obtain electrograms in the atrial region of the heart.

CONTRAINDICATIONS
- The Biosense Webster LASSO® Circular Mapping Catheter has not been shown to be safe and effective for radio frequency (RF) ablation.
- Use of the catheter may not be appropriate for use in patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.
- The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.
- The retrograde approach is contraindicated because of risk of entrapping the LASSO® in the left ventricle or valvular apparatus. The LASSO® is not recommended for use in the ventricles.

WARNINGS
- Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
- Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
- Careful consideration must be given for the use of this catheter in pregnant women.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the catheter to organic solvents such as alcohol.
- Do not autoclave the catheter.
- Do not introduce LASSO® Catheter tip folded into the guiding sheath.
- The LASSO® Deflectable Circular Mapping Catheter is recommended for use with the Biosense Webster PREFACE® Braided Guiding Sheath.

 Note: Do not use the LASSO® Deflectable Circular Mapping Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25 mm in diameter.

PRECAUTIONS
- Store in a cool, dry, dark place.
- Do not attempt to operate the Biosense Webster LASSO® Circular Mapping Catheter prior to completely reading and understanding these Instructions for Use.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance through a guiding sheath. Do not use excessive force to advance or withdraw the catheter through the guiding sheath when resistance is encountered. In addition, extra care should be taken while inserting, aspirating and manipulating the guiding sheath.
- The sterile packaging and catheter should be inspected prior to use.
- The Biosense Webster LASSO® Circular Mapping Catheter is intended for single patient use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru
WARNING: This is a controlled, confidential, and proprietary document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Released Date: 08/11/2017

- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.
- To place LASSO® catheter, torque (or rotate) shaft in a clockwise motion only.

“USE BY” DATE
- This catheter is sterilized with ethylene oxide gas.
- Do not use if the package is open or damaged.
- Use the device prior to the "Use By" date on the package label.

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SUGGESTED INSTRUCTIONS FOR USE
- Remove the catheter from its package and place it in a sterile work area.
- Follow standard practice for vessel puncture, guidewire insertion and guiding sheath use and aspiration per its Instructions for Use.
- Connect the interface connectors to the appropriate recording equipment.
- Confirm that the thumbknob is pulled back completely before insertion. Advance the catheter through the guiding sheath to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
- To place LASSO® catheter, (or rotate) shaft in a clockwise motion only.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter through the guiding sheath and dispose of it in an appropriate manner. Remove the guiding sheath, vessel dilator and guidewire as a unit per its Instructions for Use. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

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Descriptions and specifications appearing in Biosense Webster , Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru
WARNING: This is a controlled, proprietary, and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
Released Date: 08/11/2017

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BIOSENSE WEBSTER, the Biosense Webster logo, LASSO, and PREFACE are trademarks of Biosense Webster, Inc.

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Biosense Webster
PART OF THE Johnson & Johnson FAMILY OF COMPANIES

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
LASSO™ NAV Catheter

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNING: To reduce the risk of entrapping cardiac structures in the mapping electrode portion of the catheter, place the LASSO™ NAV Catheter by torquing (or rotating) the shaft in a clockwise motion only.

• STERILE. Sterilized using ethylene oxide.
• For single use only.
• Do not autoclave.
• Do not use if the package is open or damaged.

CATHETER DESCRIPTION
The Biosense Webster LASSO™ NAV Catheter has been designed to facilitate electrophysiological mapping of the atria of the heart. It is deployed in the right or left atrium through an 8 Fr guiding sheath. This deflectable catheter consists of a 4.5 Fr circular spine on its distal tip, with platinum electrodes that can be used for stimulation and recording. Models are available with different circular spine diameters.

The LASSO™ NAV Catheter provides location information when used with a compatible CARTO™ EP Navigation System.

The catheter interfaces with standard recording equipment and CARTO™ EP Navigation System equipment via interface cables with the appropriate connectors. Consult local distributor or manufacturer for the appropriate interface cables. For use in mapping with a CARTO™ EP Navigation System, refer to the instructions for the CARTO™ EP Navigation System.

For further description of the operation of the CARTO™ EP Navigation System, refer to the operating instructions for this instrument.

INDICATIONS
The LASSO™ NAV Catheter is indicated for multiple electrode electrophysiological mapping of the cardiac structures of the heart, i.e. recording or stimulation only. This Catheter is designed to obtain electrograms in the atrial regions of the heart.

The LASSO™ NAV Catheter provides location information when used with a compatible CARTO™ EP Navigation System.

CONTRAINDICTIONS
- The LASSO™ NAV Catheter has not been shown to be safe and effective for radiofrequency (RF) ablation.
- Use of the catheter may not be appropriate for use in patients with prosthetic valves.
- A relative contraindication for cardiac catheter procedures is active systemic infection.
- The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.
- The retrograde approach is contraindicated because of risk of entrapping the LASSO™ NAV Catheter in the left ventricle or valvular apparatus.
- The LASSO™ NAV Catheter is contraindicated for use in the ventricles.

WARNINGS
- Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury, as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
- Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
- Careful consideration must be given for the use of this catheter in pregnant women.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the catheter to organic solvents such as alcohol.
- Do not autoclave the catheter.
- Do not introduce the LASSO™ NAV Catheter tip folded into the guiding sheath.
- The LASSO™ NAV Catheter is recommended for use with an 8 Fr atraumatic soft tipped guiding sheath.
 Note: Do not use the LASSO™ NAV Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25 mm in diameter.

PRECAUTIONS
- The standard transseptal procedure should be followed during mapping when moving the catheter from the right atrium to the left atrium.
- To prevent entanglement of the catheter with the valves and to prevent slippage of the catheter into the ventricles, care should be taken when using the catheter in or around the atrial-ventricle valve region.
- Prior to use the physician should ensure that the intracardiac signals are recorded by all of the electrodes on the loop of the catheter.
- Do not use in the proximity of magnetic resonance imaging (MRI) equipment because the MRI equipment may induce movement of the catheter that could result in perforation.
- The catheter should be stored in its original packaging and in a cool, dry, dark place until it is used.
- Read all accessory operating instructions prior to connection of the LASSO™ NAV Catheter.
- Do not attempt to operate the LASSO™ NAV Catheter prior to completely reading and understanding these Instructions for Use.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance through a guiding sheath. Do not use excessive force to advance or withdraw the catheter through the guiding sheath, when resistance is encountered. In addition, extra care should be taken while inserting, aspirating and manipulating the guiding sheath.
- The sterile packaging and catheter should be inspected prior to use.
- The LASSO™ NAV Catheter is intended for single use only.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Do not resterilize and reuse.
- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its straight position.
- To place LASSO™ NAV Catheter, torque (or rotate) shaft in a clockwise motion only.
- Do not apply RF energy when the ablation catheter is in contact with one or more of the LASSO™ NAV Catheter electrodes.

“USE BY” DATE
- This catheter is sterilized with ethylene oxide gas.
- Do not use if the package is open or damaged.
- Use the device prior to the “Use By” date on the package label.

DISPOSAL
Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

SUGGESTED INSTRUCTIONS FOR USE
- Remove the catheter from its package and place it in a sterile work area.
- The catheter is recommended for use with an 8 Fr atraumatic soft tipped sheath. Confirm compatibility of the catheter with the guiding sheath by fully inserting and withdrawing the catheter through the irrigated sheath before clinical use. If excessive force is required or interference between the catheter and sheath is observed, use an alternate guiding sheath to avoid damaging the catheter or sheath.
- Follow standard practice for vessel puncture, guidewire insertion and guiding sheath use and aspiration per its Instructions for Use.
- Connect the interface connectors to the appropriate recording equipment.
 NOTE: Read pacing and recording equipment operator manual for proper set up and operation.
- Confirm that the thumbknob is pulled back completely before insertion ensuring that the tip is not deflected (straight position).
- To place LASSO™ NAV Catheter, torque (or rotate) shaft in a clockwise motion only.
- Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely (tip straightened). Remove the catheter through the guiding sheath and dispose of it in an appropriate manner. Remove the guiding sheath, vessel dilator and guidewire as a unit per its Instructions for Use. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

WARNING: This is a controlled, proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, propietario, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Released Date: 8/11/2017
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

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The following complications associated with cardiac catheterization have also been reported in the literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, cardiac perforation, air embolism, arrhythmias, valvular damage, pneumothorax and hemothorax.

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THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW.

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Release Date: 7/6/2017

PENTARAY™ NAV *eco* High-Density Mapping Catheter

- **CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.**
- **STERILE. Sterilized with ethylene oxide.**
- **For single use only.**
- **Do not autoclave.**
- **Do not use if the package is open or damaged.**

CATHETER DESCRIPTION

The Biosense Webster PENTARAY™ NAV *eco* High-Density Mapping Catheter is designed to facilitate electrophysiological mapping of the heart with the CARTO™ 3 EP Navigation System. It is designed for deployment in a heart chamber through an 8 F guiding sheath. This deflectable catheter consists of multiple 3 F spines on its distal tip, each spine having multiple platinum electrodes that are used for stimulation and recording. A magnetic location sensor embedded in the deflectable tip transmits location information to the CARTO™ 3 EP Navigation System. The catheter has two electrodes on the deflectable tip to provide for visualization of the tip when used with the CARTO™ 3 EP Navigation System. Pushing forward on the catheter thumbknob deflects the tip; pulling back on the thumbknob straightens the tip. This device includes an irrigation lumen for connection to a source of continuous drip anticoagulant fluid.

This catheter interfaces with standard recording equipment and the CARTO™ 3 EP Navigation System via interface cables with the appropriate connectors.

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For use in mapping procedures, refer to the instructions for the CARTO™ 3 EP Navigation System.

For further description of the operation of the CARTO™ 3 EP Navigation System, refer to the operating instructions for this instrument.

INDICATIONS FOR USE

The Biosense Webster PENTARAY™ NAV *eco* High-Density Mapping Catheter is indicated for multiple electrode electrophysiological mapping of cardiac structures in the heart, i.e., recording or stimulation only. This catheter is intended to obtain electrograms in the atrial and ventricular regions of the heart. The PENTARAY™ NAV *eco* High-Density Mapping Catheter provides location information when used with compatible CARTO™ 3 EP Navigation Systems. (This catheter is not compatible with CARTO™ 3 EP Navigation Systems prior to Version 3.x).

CONTRAINDICATIONS

- The Biosense Webster PENTARAY™ NAV *eco* High-Density Mapping Catheter has not been shown to be safe and effective for radio frequency (RF) ablation.
- Do not use PENTARAY™ Catheters in patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.
- The transseptal approach is contraindicated in patients with intracardiac thrombus or myxoma, or interatrial baffle or patch.

WARNINGS

- Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
- Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
- Careful consideration must be given for the use of this catheter in pregnant women.
- Do not immerse the proximal handle, pig tail or cable connector in fluids; electrical performance could be affected.
- Do not expose this catheter to organic solvents such as alcohol.
- Do not autoclave this catheter.
- Flush the catheter with heparinized saline prior to insertion into the body.
- Do not introduce the PENTARAY™ NAV *eco* Catheter into a guiding sheath with its distal spines folded backwards (i.e., toward the handle). Collapse the spines together using the insertion tube prior to insertion.
- The PENTARAY™ NAV *eco* Catheter is recommended for use with an 8F guiding sheath.

 Note: Do not use the PENTARAY™ NAV *eco* Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25 mm in diameter.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento propiedad, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.

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PRECAUTIONS

- **Do not use excessive force to advance or withdraw the catheter through the guiding sheath when resistance is encountered.**
- This catheter is recommended for use with an 8 F Guiding Sheath as the distal spines may be damaged if used with a sheath that is not compatible.
- The PENTARAY™ NAV *eco* Catheter should be stored in its original packaging and in a cool, dry place until used.
- Do not attempt to operate the PENTARAY™ NAV *eco* Catheter prior to completely reading and understanding these Instructions for Use.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance through a guiding sheath. In addition, extra care should be taken while inserting, aspirating, and manipulating the guiding sheath.
- The sterile packaging and catheter should be inspected prior to use.
- The PENTARAY™ NAV *eco* Catheter is intended for single patient use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.
- Flush the catheter with heparinized saline prior to insertion into the body. Always follow standard practices of using a continuous drip of anticoagulant fluid under pressure through the proximal Luer connector when the device is in the body.

“USE BY” DATE


- This catheter is sterilized with ethylene oxide (EO) gas.
- Do not use if the package is open or damaged.
- Use the device prior to the “Use By” date on the package label.

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INSTRUCTIONS FOR USE

- Remove the PENTARAY™ NAV *eco* High-Density Mapping Catheter from its package and place it in a sterile work area.
- Follow standard practice for vessel puncture, guidewire insertion, and guiding sheath use and aspiration per its Instructions for Use.
- Connect the interface connectors to the CARTO™ 3 EP Navigation System.
- Flush the catheter with heparinized saline prior to insertion into the body.
- Always follow standard practices of using a continuous drip of anticoagulant fluid under pressure through the proximal Luer connector when the device is in the body.
- Confirm that the thumbknob is pulled back completely before insertion.
- Advance the insertion tube along the catheter shaft to collapse the spines together prior to insertion into the sheath (see illustration). Do not bend the spines backward. After insertion, slide the insertion tube back toward the handle.

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PENTARAY Catheter is protected under one or more of the following U.S. Patents: 7,089,045; 7,228,164; 7,155,270; 7,302,285 and other patents pending in the U.S. and other countries.

BIOSENSE WEBSTER, the Biosense Webster logo, PENTARAY, and CARTO are trademarks of Biosense Webster, Inc.

The third-party trademarks used herein are trademarks of their respective owners.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento propiedad y confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.
© Biosense Webster, Inc. 2012-2017

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Biosense Webster
a Johnson & Johnson company

CE 0086

PENTARAY® NAV
High-Density Mapping Catheter

INSTRUCTIONS FOR USE	English	page	6
MODE D'EMPLOI	Français	page	9
GEBRAUCHSANWEISUNG	Deutsch	Seite	12
ISTRUZIONI PER L'USO	Italiano	pagina	15
INSTRUCCIONES DE USO	Español	página	18
INSTRUÇÕES DE UTILIZAÇÃO	Português	página	21
GEBRUIKSAANWIJZING	Nederlands	pagina	24
BETJENINGSVEJLEDNING	Dansk	side	27
KÄYTTÖOHJEET	Suomi	sivu	30
BRUKSANVISNING	Svensk	sida	33
BRUKSANVISNING	Norsk	side	36
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	Ελληνικά	σελίδα	39
KULLANIM TALİMATLARI	Türkçe	sayfa	42
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ	Русский	страница	45
INSTRUKCJA OBSŁUGI	Polski	strona	48
NÁVOD K POUŽITÍ	Česky	strana	51
HASZNÁLATI UTASÍTÁS	Magyar	oldal	54
NÁVOD NA POUŽÍVANIE	Slovensky	strana	57
NAVODILA ZA UPORABO	Slovensko	stran	60
ИНСТРУКЦИИ ЗА УПОТРЕБА	Български	страница	63
INSTRUCȚIUNI DE UTILIZARE	Română	pagina	66
KASUTUSJUHEND	Eesti	lehekülg	69
LIETOŠANAS PAMĀCĪBA	Latviski	lappuse	72
NAUDOJIMO INSTRUKCIJA	Lietuvių k.	puslapis	75
使用说明	简体中文	页码	78
사용 설명서	한국어	페이지	81
使用說明	繁體中文	頁碼	84
INSTRUKSI PENGGUNAAN	Bahasa Indonesia	halaman	87

WARNING: This is a controlled proprietary and confidential document.
AVISO: Este es M-5376-720[illegible] controlado, confidencial, y con derechos reservados.
Release Date: 6/13/2013
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

Sterilized with ethylene oxide gas / Stérilisé à l'oxyde d'éthylène / Sterilisiert mit Ethylenoxidgas / Sterilizzato con gas ossido di etilene / Esterilizado con gas óxido de etileno / Esterilizado com gás de óxido de etileno / Gesteriliseerd met ethyleenoxidegas / Steriliseret med ætylenoxidgas / Sterilisoitu etyleenioksidikaasulla / Steriliserad med etylenoxidgas / Sterilisert med etylenoksidgass / Αποστειρώνεται με αέριο οξείδιο του αιθυλενίου / Etilen oksit gazıyla sterilize edilmiştir / Стерилизовано этиленоксидом / Sterylizowane gazowym tlenkiem etylenu / Sterilizováno plynným etylénoxidem / Etilénoxid gázban sterilizálva / Sterilizované etylenoxidovým plynom / Sterilizirano z etilenoksidom / Стерилизиран с газ етиленов оксид / Sterilizat cu oxid de etilenă / Steriliseeritud etüleenoksiidgaasiga / Sterilizēts ar gāzveida etilēnoksīdu / Sterilizuota etileno oksido dujomis / 已采用环氧乙烷气体灭菌 / 산화에틸렌 가스로 멸균 처리 / 已採用環氧乙烷氣體滅菌 / Disterilkan dengan menggunakan etilen oksida

For single use only / À usage unique / Nur zum Einmalgebrauch / Esclusivamente monouso / Para un solo uso / Destina-se apenas a uma única utilização / Uitsluitend voor eenmalig gebruik / Kun til engangsbrug / Vain kertakäyttöön / Endast för engångsbruk / Kun til engangsbruk / Για μία χρήση μόνο / Tek kullanımlıktır / Только для одноразового использования / Instrument przeznaczony wyłącznie do jednorazowego użytku / Určeno k jednorázovému použití / Kizárólag egyetlen használatra / Len na jednorazové použitie / Samo za enkratno uporabo / Само за еднократна употреба / Exclusiv de unică folosință / Üksnes ühekordseks kasutamiseks / Tikai vienreizējai lietošanai / Tik vienkartinio naudojimo / 只能使用一次 / 일회용 / 只能使用一次 / Hanya untuk digunakan satu kali

Refer to accompanying Instructions for Use / Consultez le mode d'emploi joint / Siehe zugehörige Gebrauchsanleitung / Fare riferimento alle Istruzioni per l'uso allegate / Consulte las Instrucciones de uso adjuntas / Consultar as Instruções para Utilização inclusas / Raadpleeg de bijbehorende gebruiksaanwijzing / Se medfølgende brugervejledning / Katso käyttöohjeita / Se medföljande bruksanvisning / Se medfølgende bruksinstruksjoner / Ανατρέξτε στις συνοδευτικές οδηγίες χρήσης / Beraberindeki Kullanma Talimatı'na bakınız / См. Инструкцию по применению, входящую в комплект поставки / Sprawdzić w załączonej Instrukcji stosowania / Viz přiložený návod k použití / Használja útmutatóul a mellékelt Használati utasítást / Odkazujeme na sprievodný Návod na použitie / Upoštevajte priložena navodila za uporabo / Консултирайте се с указанията за употреба / Consultați Instrucțiunile de utilizare care însoțesc produsul / Tutvuge kaasasoleva kasutusjuhendiga / Iepazīstieties ar pievienotajiem lietošanas norādījumiem / Žr. pridedamą Naudojimo instrukciją / 请参阅随附的使用说明 / 동봉한 사용 설명서를 참조하십시오 / 請參閱隨附的使用說明 / Lihat instruksi penggunaan yang terlampir

Do not use open or damaged packages / Ne pas utiliser si l'emballage est ouvert ou endommagé / Geöffnete oder beschädigte Verpackungen nicht verwenden / Non utilizzare confezioni aperte o danneggiate / No utilice envases abiertos o estropeados / Não utilizar embalagens abertas ou danificadas / Geopende of beschadigde verpakkingen niet gebruiken / Bør ikke bruges, hvis emballagen har været åbnet eller er beskadiget / Älä käytä avattuja tai vaurioituneita pakkauksia / Använd inte öppnade eller skadade förpackningar / Bruk ikke åpne eller skadde pakninger / Μη χρησιμοποιείτε αν η συσκευασία έχει ανοιχτεί ή υποστεί βλάβη / Использование содержимого открытых или поврежденных упаковок запрещено / Açık veya hasarlı ambalajları kullanmayınız / Nie używać produktu, którego opakowanie było otwarte lub uszkodzone / Nepoužívejte, jsou-li obaly otevřené či poškozené / Ne használjon nyitott vagy sérült csomagokat / Nepoužívajte otvorené alebo poškodené balenia / Ne uporabljajte, če je ovojnina odprta ali poškodovana / Не използвайте отворени или повредени опаковки / A nu se utiliza produsul dacă ambalajul este deschis sau deteriorat / Ärge kasutage, kui pakend on avatud või katki / Neizmantot, ja iepakojumi ir bojāti vai atvērti / Nenaudokite atidarytų arba sugadintų pakuočių / 切勿使用已打开或破损的包装 / 포장이 개방되었거나 손상된 경우 사용하지 마십시오 / 請勿使用已開啟或破損的包裝 / Jangan menggunakan kemasan yang terbuka atau yang rusak

Store in a cool, dry place / Entreposer dans un endroit sec et frais / An einem kühlen, trockenen Ort aufzubewahren / Conservare in un luogo fresco e asciutto / Almacénese en un lugar fresco y seco / Armazenar num local fresco e seco / Bewaren op een koele, droge plaats / Opbevares et køligt, tørt sted / Säilytä viileässä ja kuivassa / Förvaras på en sval, torr plats / Oppbevar på et kjølig og tørt sted / Αποθηκεύετε σε δροσερό, ξηρό μέρος / Serin ve kuru bir yerde saklayınız / Хранить в прохладном, сухом месте / Przechowywać w chłodnym, suchym miejscu / Uchovávejte na chladném, suchém místě / Száraz, hűvös helyen tárolandó / Skladujte na chladnom a suchom mieste / Hranite v hladnem in suhem prostoru / Да се съхранява на хладно и сухо място / A se păstra la loc răcoros şi uscat / Hoida jahedas ning kuivas kohas / Uzglabāt vēsā, sausā vietā / Laikykite vėsioje ir sausoje vietoje / 贮存在阴凉干燥处 / 서늘하고 건조한 곳에 보관하십시오 / 貯存在陰涼乾燥處 / Simpan di tempat yang sejuk dan kering

Use by / Utiliser avant le / Verwendbar bis / Scadenza / Fecha de caducidad / Utilizar até / Uiterste gebruiksdatum / Anvendes inden / Viimeinen käyttöpäivämäärä / Använd före / Brukes innen / Ημερομηνία λήξης / Son kullanma tarihi / Использовать до / Użyć przed / Použít do / Felhasználható: / Použitie do / Uporabno do / Срок на годност / Utilizat de: / Kõlblik kuni / Izlietot līdz / Naudoti iki / 使用期限 / 유효기간 / 使用期限 / Gunakan Sebelum

WARNING: This is a controlled proprietary and confidential document.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados.

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Release Date: 6/13/2013

PENTARAY® NAV High-Density Mapping Catheter

- **STERILE. Sterilized with ethylene oxide (EO) gas.**
- **For single use only.**
- **Do not autoclave.**
- **Do not use if the package is open or damaged.**
- **CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician.**

CATHETER DESCRIPTION

The Biosense Webster PENTARAY® NAV High-Density Mapping Catheter is designed to facilitate electrophysiological mapping of the heart with the CARTO® 3 EP Navigation System. It is designed for deployment in a heart chamber through an 8F guiding sheath. This deflectable catheter consists of multiple 3F spines on its distal tip, each spine having multiple platinum electrodes that are used for stimulation and recording. A magnetic location sensor embedded in the deflectable tip transmits location information to the CARTO® 3 EP Navigation System. The catheter has two electrodes on the deflectable tip to provide for visualization of the tip when used with the CARTO® 3 EP Navigation system. Pushing forward on the catheter thumbknob deflects the tip; pulling back on the thumbknob straightens the tip. This device includes an irrigation lumen for connection to a source of continuous drip anticoagulant fluid.

This catheter interfaces with standard recording equipment and the CARTO® 3 EP Navigation System via interface cables with the appropriate connectors.

CONSULT THE LOCAL DISTRIBUTOR OR MANUFACTURER FOR THE APPROPRIATE INTERFACE CABLES.

For use in mapping procedures, refer to the instructions for the CARTO® 3 EP Navigation System.

For further description of the operation of the CARTO® 3 EP Navigation System, refer to the operating instructions for this instrument.

INDICATIONS FOR USE

The Biosense Webster PENTARAY® NAV High-Density Mapping Catheter is indicated for multiple electrode electrophysiological mapping of cardiac structures in the heart, i.e., recording or stimulation only. This catheter is intended to obtain electrograms in the atrial and ventricular regions of the heart. The PENTARAY® NAV High-Density Mapping Catheter provides location information when used with compatible CARTO® 3 EP Navigation Systems. (This catheter is not compatible with CARTO® 3 EP Navigation Systems prior to Version 3.0).

CONTRAINDICATIONS

- The Biosense Webster PENTARAY® NAV High-Density Mapping Catheter has not been shown to be safe and effective for radio frequency (RF) ablation.
- Use of this catheter may not be appropriate for use in patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.
- The transseptal approach is contraindicated in patients with intracardiac thrombus or myxoma, or interatrial baffle or patch.

WARNINGS

- Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
- Cardiac catheterization should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure.
- Careful consideration must be given for the use of this catheter in pregnant women.
- Do not immerse the proximal handle, pig tail or cable connector in fluids; electrical performance could be affected.
- Do not expose this catheter to organic solvents such as alcohol.
- Do not autoclave this catheter.
- Flush the catheter with heparinized saline prior to insertion into the body.
- Do not introduce the PENTARAY® NAV Catheter into a guiding sheath with its distal spines folded backwards (i.e., toward the handle). Collapse the spines together using the insertion tube prior to insertion.
- The PENTARAY® NAV Catheter is recommended for use with an 8F guiding sheath.

 Note: Do not use the PENTARAY® NAV Catheter in conjunction with transseptal sheaths featuring side holes larger than 1.25 mm in diameter.

WARNING: This is a controlled proprietary and confidential document.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados.
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

PRECAUTIONS

- **Do not use excessive force to advance or withdraw the catheter through the guiding sheath when resistance is encountered.**
- This catheter is recommended for use with an 8F Guiding Sheath as the distal spines may be damaged if used with a sheath that is not compatible.
- The PENTARAY® NAV Catheter should be stored in its original packaging and in a cool, dry place until used.
- Do not attempt to operate the PENTARAY® NAV Catheter prior to completely reading and understanding these Instructions for Use.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance through a guiding sheath. In addition, extra care should be taken while inserting, aspirating, and manipulating the guiding sheath.
- The sterile packaging and catheter should be inspected prior to use.
- The PENTARAY® NAV Catheter is intended for single patient use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.
- Flush the catheter with heparinized saline prior to insertion into the body. Always follow standard practices of using a continuous drip of anticoagulant fluid under pressure through the proximal Luer connector when the device is in the body.

"USE BY" DATE


- This catheter is sterilized with ethylene oxide (EO) gas.
- Do not use if the package is open or damaged.
- Use the device prior to the "Use By" date on the package label.

DISPOSAL

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

INSTRUCTIONS FOR USE

- Remove the PENTARAY® NAV High-Density Mapping Catheter from its package and place it in a sterile work area.
- Follow standard practice for vessel puncture, guidewire insertion, and guiding sheath use and aspiration per its Instructions for Use.
- Connect the interface connectors to the CARTO® 3 EP Navigation System.
- Flush the catheter with heparinized saline prior to insertion into the body.
- Always follow standard practices of using a continuous drip of anticoagulant fluid under pressure through the proximal Luer connector when the device is in the body.
- Confirm that the thumbknob is pulled back completely before insertion.
- Advance the insertion tube along the catheter shaft to collapse the spines together prior to insertion into the sheath (see illustration). Do not bend the spines backward. After insertion, slide the insertion tube back toward the handle.


- Advance the catheter through the guiding sheath to the area of the endocardium under evaluation. The spines will naturally return to their expanded pattern as they emerge from the tip of the guiding sheath. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to curve; pulling the thumbknob back straightens the catheter tip.

WARNING: This is a controlled proprietary and confidential document.
AVISO: Este es M-5376-7205 es controlado, confidencial, y con derechos reservados.
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

Release Date: 6/13/2013

- Under fluoroscopy, catheter spine ‘A’ can be identified by the marker band near ring 2. Catheter spine ‘B’ can be identified by the marker band near ring 8. Once spines A and B have been identified, all other spines and electrodes may be identified as follows:

- For questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely to straighten the tip. Remove the catheter from the guiding sheath and dispose of it in an appropriate manner.
- Do not resterilize and reuse this device.

ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterization procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death.

The following complications associated with cardiac catheterization have also been reported in the literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, cardiac perforation, air embolism, arrhythmias, valvular damage, pneumothorax, and hemothorax.

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW.

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product at the time of manufacture and are not made or given as a warranty of the prescribed product in any way.

WARNING: This is a controlled proprietary and confidential document.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados.
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

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European Community Representative
Biosense Webster
A Division of Johnson & Johnson Medical NV/SA
Leonardo da Vincilaan 15
1831 Diegem, Belgium
Tel: +32-2-7463-401
Fax: +32-2-7463-403

Protected under one or more of the following U.S. Patent Nos.: 5,391,199; 5,443,489; 5,480,422; 5,546,951; 5,558,091; 5,568,809; 5,694,945; 5,713,946; 5,718,241; 5,738,096; 7,089,045; 7,228,164; 7,155,270; 7,302,285 and other patents pending in the U.S. and other countries.

BIOSENSE WEBSTER, the Biosense Webster logo, PENTARAY, and CARTO® are trademarks of Biosense Webster, Inc.

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CE 0086

WEBSTER™ Deflectable Tip Electrophysiology Catheter with Auto ID Technology

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WEBSTER™ Deflectable Tip Electrophysiology Catheter with Auto ID Technology

Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

- **Sterile. Sterilized Using Ethylene Oxide.**
- **For single use only.**
- **Do not use if the package is open or damaged.**

CATHETER DESCRIPTION

The WEBSTER™ Deflectable Tip Electrode Catheter has been designed to facilitate electrophysiological mapping of the heart. The catheter has a high-torque shaft with a deflectable tip section containing an array of platinum electrodes that can be used for stimulation and recording. Several tip configurations are available. (Standard curve types are identified by a color band on the shaft near the proximal end - see catalog.) The catheter is equipped with Electronically Erasable Programable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO™ EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information.

A piston in the handpiece is attached to an internal puller wire which changes the radius of curvature. When the piston is pushed forward with the thumbknob, the radius of curvature is reduced; when the thumbknob is pulled back, the radius of the curvature is increased until the tip section straightens. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

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CONSULT THE LOCAL DISTRIBUTOR OR MANUFACTURER FOR THE APPROPRIATE INTERFACE CABLES. REFER TO THE CARTO™ EP NAVIGATION SYSTEM USER MANUAL FOR DETAILS ON SYSTEM SETTINGS RELATED TO THE CATHETER.

INDICATIONS

The catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only.

CONTRAINDICATIONS

The WEBSTER™ Deflectable Tip Electrode Catheter has not been shown to be safe and effective for electrical ablation or for use in the coronary vasculature other than the coronary sinus ostium. Use of the catheter may not be appropriate for patients with prosthetic valves. A relative contraindication for cardiac catheter procedures is active systemic infection.

WARNINGS

Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Cardiac catheterization should only be performed after adequate attention has been given to potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must therefore be given for the use of this catheter in pregnant women.

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PRECAUTIONS

- Do not attempt to operate the WEBSTER™ Deflectable Tip Electrophysiology Catheter prior to completely reading and understanding these directions for use.
- Store in a cool, dry place.
- The sterile packaging and catheter should be inspected prior to use.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Always pull the thumbknob of the catheter back before insertion or withdrawal to assure that the catheter tip assumes its original shape.

"USE BY" DATE

Use the device prior to the "Use By" date on the package label.

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SUGGESTED INSTRUCTIONS FOR USE

- Remove the catheter from its package and place it in a sterile work area. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the interface cable to the Patient Interface Unit of the appropriate CARTO™ EP Navigation System and connect the catheter to the interface cable.
- Confirm that the thumbknob is pulled back completely before insertion. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. Adjust the radius of curvature as necessary by manipulating the thumbknob. Pushing the thumbknob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
- Prior to removal of the catheter, confirm that the thumbknob has been pulled back completely. Remove the catheter and dispose of it in an appropriate manner. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterization procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death. The following complications associated with cardiac catheterization have also been reported in the literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, cardiac perforation, air embolism, arrhythmias, valvular damage, pneumothorax and hemothorax.

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WEBSTER™ CS Catheter with Auto ID Technology

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DEVICE DESCRIPTION
The WEBSTER™ Coronary Sinus Catheter is a steerable, multi-electrode catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart. The device is a 6 FR catheter with a usable length of 115 cm. The catheter has a high-torque shaft with a deflectable tip section containing an array of platinum electrodes that can be used for stimulation and recording.

Standard features of this catheter include a braided 6 FR deflectable tip section with an array of platinum electrodes that includes a 2 mm tip dome. The braided tip is controlled by a proximal hand piece that features a thumb operated sliding piston, and is offered in various curve types. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

The catheter is equipped with Electronically Erasable Programmable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO™ EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information.

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INDICATIONS AND USAGE
The WEBSTER™ CS Catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. The catheter is designed for use in the coronary sinus.

CONTRAINDICATIONS
The WEBSTER™ CS Catheter has not been shown to be safe and effective for electrical ablation.

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WARNINGS AND PRECAUTIONS
- Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff. Significant x-ray exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects due to the x-ray beam intensity and duration of the fluoroscopic imaging. The long-term risks of protracted fluoroscopy have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children and pregnant women.
- Do not attempt to operate the WEBSTER™ CS Catheter prior to completely reading and understanding the applicable Instructions for Use.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Always ensure that the thumb knob is pulled back completely to straighten the catheter tip before insertion or withdrawal of the catheter.
- Do not use excessive force to advance or withdraw the catheter. Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.
- Use both fluoroscopy and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid vascular or cardiac damage.

HANDLING AND STERILIZATION
- **For single use only.**
- **Do not sterilize.**
- **This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.**

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Catheter damage may occur due to:
- Autoclaving
- Resterilizing
- exposure to organic solvents
- immersing proximal handle or cable connector in fluids

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HOW SUPPLIED

- The available curves for WEBSTER™ CS Catheters include D, F and J.
- Appropriate interface cables are supplied separately.

Packaging

The WEBSTER™ CS Catheter is supplied STERILE (EO). The catheter is secured in a mounting card and placed into a Tyvek/Nylon film pouch, sealed, and placed in a box. Both the pouch and mounting card are sterile unless the package is damaged or opened.

Storage

The WEBSTER™ CS Catheter must be stored in a cool, dry, dark place. Storage temperature should be between 5º and 25ºC (41º and 77ºF).

Shelf-Life

Product and package testing have been conducted to support the "Use By" date printed on the product labels. DO NOT USE after the "Use By" date.

DIRECTIONS FOR USE

Physician Training

Physicians must be familiar with the techniques and appropriately trained for cardiac mapping procedures. All mapping procedures must be performed in a fully equipped electrophysiology laboratory.

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Handling and Preparation

Before use, inspect the packaging. Do not use if open or damaged. Using aseptic technique, remove the catheter from its package and place it in a sterile working area. Inspect the catheter carefully for electrode integrity and overall condition.

Operating Instructions

1. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
2. Connect the interface cable to the Patient Interface Unit of the appropriate CARTO™ EP Navigation System and connect the catheter to the interface cable.
3. Confirm that the thumb knob is pulled back completely before insertion (catheter in straight position).
4. Advance the catheter into the coronary sinus. Use both fluoroscopy and electrograms to aid in proper positioning.
5. Adjust the radius of curvature as necessary by manipulating the thumb knob. Pushing the thumb knob forward causes the catheter tip to bend (curve); when the knob is pulled back, the tip straightens.
6. Prior to removal of the catheter, confirm that the thumb knob has been pulled back completely (straightened). Remove the catheter and dispose of it in an appropriate manner.

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The following complications associated with cardiac catheterization have also been reported in literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, caridiac perforation, tamponade, thrombi, air embolism, arrhythmias, valvular damage, pneumothorax and hemothorax.

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WEBSTER™ CS Catheter with EZ STEER™and Auto ID Technology

Caution: Federal (USA) law restricts this device to sale by or on the order of the physician.

FOR SINGLE USE ONLY. DO NOT RESTERILIZE.

DEVICE DESCRIPTION

The WEBSTER™ Coronary Sinus Catheter is a steerable, multi-electrode catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart. The device is a 7 FR catheter with a usable length of 115 cm. The catheter has a high-torque shaft with a bi-directional deflectable tip section containing an array of platinum electrodes that can be used for stimulation and recording.

Standard features of this catheter include a braided bi-directional deflectable tip section with an array of platinum electrodes that includes a 2 mm tip dome. The braided bi-directional tip provides the user with two 180° opposed single plane curves. Currently, the available curves for WEBSTER™ CS Catheter include DF and FJ. A Rocker Lever is used to deflect the tip. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

The catheter is equipped with Electronically Erasable Programable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO™ EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information.

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INDICATIONS AND USAGE

The WEBSTER™ CS Catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. The catheter is designed for use in the coronary sinus.

CONTRAINDICATIONS

The WEBSTER™ CS Catheter has not been shown to be safe and effective for electrical ablation.

- Use of the catheter may not be appropriate for patients with prosthetic valves.
- A relative contraindication for cardiac catheter procedures is active systemic infection.

Do not use this device:

- Via the transseptal approach in patients with left atrial thrombus or myxoma or interatrial baffle or patch;
- Via the retrograde approach in patients with aortic valve replacement.

WARNINGS AND PRECAUTIONS

- Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff. Significant x-ray exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects due to the x-ray beam intensity and duration of the fluoroscopic imaging. The long-term risks of protracted fluoroscopy have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children and pregnant women.
- Do not attempt to operate the WEBSTER™ CS Catheter prior to completely reading and understanding the applicable Instructions for Use.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Always place the Rocker Lever in the neutral position to straighten the catheter tip before insertion or withdrawal of the catheter.
- Do not use excessive force to advance or withdraw the catheter. Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.
- Use both fluoroscopy and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid vascular or cardiac damage.

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- **This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.**

Catheter damage may occur due to:

- Autoclaving
- Resterilizing
- exposure to organic solvents
- immersing proximal handle or cable connector in fluids

Observe “Use By” Date. Sterilized with ethylene oxide gas.

The sterile packaging and catheter should be inspected prior to use. If the package or the catheter appears damaged, do not use. Contact your local Biosense Webster representative.

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HOW SUPPLIED

- The available curves for WEBSTER™ CS Catheters include DF and FJ.
- Appropriate interface cables are supplied separately.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Released Date: 08/11/2017

Packaging
The WEBSTER™ CS Catheter is supplied STERILE (EO). The catheter is secured in a two-piece thermoform tray and placed into a Tyvek/Nylon film pouch, sealed, and placed in a box. Both the pouch and thermoform tray are sterile unless the package is damaged or opened.

Storage
The WEBSTER™ CS Catheter must be stored in a cool, dry place. Storage temperature should be between 5º and 25ºC (41º and 77ºF).

Shelf-Life
Product and package testing have been conducted to support the "Use By" date printed on the product labels. DO NOT USE after the "Use By" date.

DIRECTIONS FOR USE

Physician Training
Physicians must be familiar with the techniques and appropriately trained for cardiac mapping procedures. All mapping procedures must be performed in a fully equipped electrophysiology laboratory.

Compatible Accessories: Use appropriate Biosense Webster accessory cables to connect the WEBSTER™ CS Catheter to the Patient Interface Unit of the appropriate CARTO™ EP Navigation System and connect the catheter to the interface cable.

Handling and Preparation
Before use, inspect the packaging. Do not use if open or damaged. Using aseptic technique, remove the catheter from its package and place it in a sterile working area. Inspect the catheter carefully for electrode integrity and overall condition.

Operating Instructions

1. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
2. Connect the interface cable to the Patient Interface Unit of the appropriate CARTO™ EP Navigation System and connect the catheter to the interface cable.
3. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.
4. Use the Rocker Lever to deflect the catheter tip (Figure 1). When the lever is pulled back from neutral, the tip will deflect relative to the direction of rotation. The amount of deflection is relative to the amount of lever rotation. When the lever is pushed forward, the tip will deflect in the opposite direction. To straighten the tip, return the Rocker Lever to neutral position.



Figure 1

The handle has an adjustable friction control that allows the operator to use the Rocker Lever and deflecting tip in a “free” state or adjust the friction to where the Rocker Lever and tip curve are “locked” in place (Figure 2). This knob is located on the opposite side of the Rocker Lever. Out of the package, the knob will be in the “off“ position which allows free movement for the lever and deflecting tip. The amount of friction increases as the Friction Control Knob is rotated clockwise until it reaches the fully “on” position.

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Figure 2

5. Determine that the electrodes are in stable contact with the intended mapping site.

ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterization procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death.
The following complications associated with cardiac catheterization have also been reported in literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, caridiac perforation, tamponade, thrombi, air embolism, arrhythmias, valvular damage, pneumothorax and hemothorax.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Released Date: 08/11/2017

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THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN.

UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW.

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BIOSENSE WEBSTER, the Biosense Webster logo, EZ STEER, CARTO, and WEBSTER are trademarks of Biosense Webster, Inc.

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WEBSTER™ CS Catheter with EZ STEER™ and Auto ID Technology protected under one or more of the following U.S. Patents: 7,377,906; 7,591,799; 8,021,327; 8,348,888 and other patents pending in the U.S. and other countries.

Rev: 2017-07

M-5276-581F © Biosense Webster, Inc. 2000-2017 Printed in USA

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada.

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CE
0086

WEBSTER™CS
Catheter with EZ STEER™ Technology

Instructions for Use English Page 9
WEBSTER™**CS Catheter with** EZ STEER™ **Technology**

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WEBSTER™CS Catheter with EZ STEER™ Technology

Caution: **Federal (USA) law restricts this device to sale by or on the order of a physician.**

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DEVICE DESCRIPTION
The WEBSTER™Coronary Sinus Catheter is a steerable, multi-electrode catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart. The device is a 7 FR catheter with a usable length of 115 cm. The catheter has a high-torque shaft with a bi-directional deflectable tip section containing an array of platinum electrodes that can be used for stimulation and recording.

Standard features of this catheter include a braided bi-directional deflectable tip section with an array of ten platinum electrodes that includes a 2 mm tip dome. Additionally, two asymmetric curve types are available providing two 180° opposed, single plane curves. Currently, the available curves for WEBSTER™CS include DF and FJ. A Rocker Lever is used to deflect the tip. The high-torque shaft allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site.

INDICATIONS AND USAGE
The WEBSTER™CS Catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only. The catheter is designed for use in the coronary sinus.

CONTRAINDICATIONS
The WEBSTER™CS Catheter has not been shown to be safe and effective for electrical ablation.
- Use of the catheter may not be appropriate for patients with prosthetic valves.
- A relative contraindication for cardiac catheter procedures is active systemic infection.

Do not use this device:
- Via the transseptal approach in patients with left atrial thrombus or myxoma or interatrial baffle or patch;
- Via the retrograde approach in patients with aortic valve replacement.

WARNINGS AND PRECAUTIONS
- Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff. Significant x-ray exposure can result in acute radiation injury as well as increased risk for somatic and genetic effects due to the x-ray beam intensity and duration of the fluoroscopic imaging. The long-term risks of protracted fluoroscopy have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children and pregnant women.
- Do not attempt to operate the WEBSTER™CS Catheter prior to completely reading and understanding the applicable Instructions for Use.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Always place the Rocker Lever in the neutral position to straighten the catheter tip before insertion or withdrawal of the catheter.
- Do not use excessive force to advance or withdraw the catheter. Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.
- Use both fluoroscopy and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid vascular or cardiac damage.

Handling and Sterilization

FOR SINGLE USE ONLY. DO NOT RESTERILIZE.
Catheter damage may occur due to:
- autoclaving
- resterilizing
- exposure to organic solvents
- immersing proximal handle or cable connector in fluids
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Observe “Use By” Date.
- Sterilized with ethylene oxide gas.
- The sterile packaging and catheter should be inspected prior to use. If the package or the catheter appears damaged, do not use. Contact your local Biosense Webster representative.

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HOW SUPPLIED
- The available curves for WEBSTER™CS Catheters include DF and FJ.
- Appropriate interface cables are supplied separately.

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Released Date: 08/11/2017

Packaging
The WEBSTER™CS Catheter is supplied STERILE (EO). The catheter is secured in a two-piece thermoform tray and placed into a Tyvek/Nylon film pouch, sealed, and placed in a box. Both the pouch and thermoform tray are sterile unless the package is damaged or opened.

Storage
The WEBSTER™CS Catheter must be stored in a cool, dry and dark place. Storage temperature should be between 5º and 25ºC (41º and 77ºF).

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Compatible Accessories: Use appropriate Biosense Webster accessory cables to connect the WEBSTER™CS Catheter to standard recording equipment.

Handling and Preparation
Before use, inspect the packaging. Do not use if open or damaged.
Using aseptic technique, remove the catheter from its package and place it in a sterile working area. Inspect the catheter carefully for electrode integrity and overall condition.

Operating Instructions
1. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
2. Connect the catheter to the interface cables and standard recording equipment using the appropriate interface cables.
3. Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.
4. Use the Rocker Lever to deflect the catheter tip (Figure 1). When the lever is pulled back from neutral, the tip will deflect relative to the direction of rotation. The amount of deflection is relative to the amount of lever rotation. When the lever is pushed forward, the tip will deflect in the opposite direction. To straighten the tip, return the Rocker Lever to neutral position.



Figure 1

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Clockwise rotation from the “off” position increases the friction within the Deflection Mechanism.

Counter-Clockwise rotation from the “on” position decreases the friction within the Deflection Mechanism.

Figure 2

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ADVERSE REACTIONS

A number of serious adverse reactions have been documented for cardiac catheterization procedures including pulmonary embolism, myocardial infarction, stroke, cardiac tamponade, and death. The following complications associated with cardiac catheterization have also been reported in literature: vascular bleeding, local hematomas, thrombosis, AV fistula, pseudoaneurysm, thromboembolism, vasovagal reactions, caridiac perforation, tamponade, thrombi, air embolism, arrhythmias, valvular damage, pneumothorax and hemothorax.

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RANTY AND LIMITATION OF LIABILITY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER **, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.**

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER **, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.**

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WEBSTER® Fixed Curve Electrophysiology Catheter with Auto ID Technology

Federal law (USA) restricts this device to sale by or on the order of a physician.

- **STERILE. Sterilized with ethylene oxide gas.**
- **For single use only.**
- **Do not use if the package is open or damaged.**

CATHETER DESCRIPTION
The WEBSTER® Fixed Curve Electrophysiology (EP) Catheter has been designed to facilitate electrophysiological mapping of cardiac structures. The catheter has a high-torque shaft with an array of platinum electrodes at the distal tip that can be used for stimulation and recording. The tip section has a preformed curve with a nylon tip protector. (Standard curve types are identified by a color band on the shaft near the proximal end - see catalog). The catheter is equipped with Electronically Erasable Programable Read Only Memory (EEPROM) which is used to store unique catheter identification information. CARTO™ EP Navigation Systems equipped with Auto ID Technology can access the stored information and automatically recognize the catheter information.

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INDICATIONS
The catheter is indicated for electrophysiological mapping of cardiac structures; i.e., stimulation and recording only.

CONTRAINDICATIONS
- The WEBSTER® Fixed Curve EP Catheter has not been shown to be safe and effective for electrical ablation or for use in the coronary vasculature other than the coronary sinus ostium.
- Use of the catheter may not be appropriate for patients with prosthetic valves.
- A relative contraindication for cardiac catheter procedures is active systemic infection.

WARNINGS
Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Cardiac catheterization should only be performed after adequate attention has been given to potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must therefore be given for the use of this catheter in pregnant patients.

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PRECAUTIONS
- Do not attempt to operate the WEBSTER® Fixed Curve EP Catheter prior to completely reading and understanding these directions for use.
- Cardiac catheterization procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement and placement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Store in a cool, dry place.
- The sterile packaging and catheter should be inspected prior to use.
- WEBSTER® Fixed Curve EP Catheter is intended for single patient use only.
- Do not resterilize and reuse.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

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SUGGESTED DIRECTIONS FOR USE
- Remove the catheter from its package and place it in a sterile work area. Remove nylon tip protector.
- Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
- Connect the interface cable to the Patient Interface Unit of the appropriate CARTO™ EP Navigation System and connect the catheter to the interface cable.
- Advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning. The electrode catheter shaft may be

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.
AVISO: Este es un documento controlado, confidencial, y con derechos reservados. Revisar si es la revision mas actualizada. Released Date: 08/11/2017

gently hand rotated to facilitate positioning at the desired mapping site.

- Remove the catheter and dispose it in an appropriate manner. Do not resterilize and reuse.
- If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer

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Descriptions and specifications appearing in Biosense webster printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

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BIOSENSE WEBSTER, the Biosense Webster logo, CARTO, and HALO are trademarks of Biosense Webster, Inc.
BIOSENSE WEBSTER and the Biosense Webster logo are trademarks of Biosense Webster, Inc.
BIOSENSE WEBSTER, the Biosense Webster logo, CARTO and ISMUS are trademarks of Biosense Webster, Inc.
© 2018 Biosense Webster, Inc. All rights reserved. Confidential, y con derechos reservados.
BIOSENSE WEBSTER, the Biosense Webster logo, CARTO and LASSO are trademarks of Biosense Webster, Inc.
Copyright © 2017 Biosense Webster, Inc. All rights reserved. Todos los derechos reservados. Revisar si es la revision mas actualizada.