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GB. Reamers type K manual

Instruction for use

Reamers type K are designed for enlarging and smoothing of root.

Application field - deontology.

Reamers type are designed for application and storage at the temperature + 10- 35° C with relative humidity below 80% and temperature + 25° C.

The air in the storage room shall be free from corrosion agents. Reamers type K may be used 5 times or more. Supplied non sterile.

SAFETY MEASURES

- Reamers type K may be used by qualified medical personnel with sufficient experience and good manual skills only;
- Do not use instruments if they are corroded or show signs of mechanic damage;
- Do not use instruments if its nominal diameter is smaller than that of root canal;
-Do not touch operative part of the instrument in order to avoid cuts over sharp edges of reamers;

Use the instrument according to its intended purpose.

Reamers type K material is non toxic. Used instruments shall be disposed under general method established for disposal of such instruments.

PREPARATION FOR OPERATION

Operative part material - stainless steel, handle material – aluminum alloy

Reamers type K shall be disinfected, cleaned and sterilised before application. Use neutral pH solutions for this purpose.

Disinfection: Prior to first and repeated application of the instrument sink it in the disinfection solution. Observe the soak time and follow recommendations specified by the manufacturer of the disinfection solution. When finished, rinse the instruments.

Cleaning: reamers type K are subject to ultrasonic cleaning within 10 min. After that rinse them in distilled water for 0.5 min. Dry at the temperature 85° C until water completely evaporates.

Sterilisation: prior to sterilisation place the instruments on sterilisation stands or trays. Sterilise by

autoclave: at the temperature 134°C and pressure 0.21 MPa, hold time 5 min, or at the tempera-ture

131°C and pressure 0.11 MPa, hold time 20 min.

When sterilising several instruments at once make sure that autoclave is not overloaded. After sterilisation and cooling place the instruments in sealed containers that provide sterile storage conditions and eliminate accidental mechanic damage. Remove the instrument by the handle using sterile forceps.

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