H-517 Ver.4
2012-08-29

CELLO™
Balloon Guide Catheter
Intravascular Catheter
with Cuff

English
Instructions for Use

Device Description
The Cello™ Balloon Guide Catheter is a coaxial-lumen, braid-reinforced, variable stiffness catheter with two radiopaque markers on both the distal and proximal ends of the balloon and a bifurcated Luer hub on the proximal end. A compliant silicone balloon is mounted on the distal end. Balloon Guide Catheter shaft material and recommended balloon inflation volumes are indicated on the product label. If indicated on product label, a dilator is provided.

Indications for Use
This catheter is intended to temporarily block blood flow by expanding a balloon inside blood vessels during surgery. It can be used for percutaneous, transluminal, therapeutic surgery, perfusion of blood to peripheral vessel, and arterial injection for chemotherapy.

Contraindications
Use of the Cello™ Balloon Guide Catheter is contraindicated under these circumstances:
- Patients having acute cardiac failure, irregular heartbeat or abnormal serum electrolyte.
- Patients having severe disorder of liver or kidney function.
- Patients having a blood coagulation disorder or severe changes in coagulation for extravasation.
- Patients who cannot be supine on the angiography bed because of congestive heart failure or respiratory problems.
- Pregnant women.
- Patients who are uncooperative during examination.

Complications
Possible complications include, but are not limited to, the following: infection, hematoma, distal embolization, vessel thrombosis, dissection, false aneurysm formation, acute occlusion, clot formation, hemorrhage at the puncture site, intracranial hemorrhage, arterial rupture, arterial injury, arterial perforation, Arteriovenous fistula, infection at the puncture site, nausea and vomiting, blood pressure drop, arrhythmia, bile bleeding, vasospasm, obstruction of bile duct, nervous system injury, stroke and death.

Compatibility
Polyethylene sheath size and maximum guidewire diameter are indicated on the product label.

Warnings
- Single use only. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning.
- Never advance or torque catheter against resistance without careful assessment of cause of resistance and necessary fluoroscopy. If cause cannot be determined, withdraw catheter. Movement against resistance may result in damage to vessel or catheter.
- To reduce risk of complications due to slow balloon deflation, adhere to the following recommendations:
 - Wet the distal shaft with saline before passing into the introducer sheath.
 - Minimize pushing forces on the shaft during advancement. These forces can cause wrinkles in the balloon, which may impair deflation.
 - Do not use the device if the shaft is damaged during use.
 - Prepare balloon according to the recommended procedure
- To reduce the risk of complications due to an emboli, remove air from the balloon according to the recommended procedure.
- Withdrawing the balloon deeper than the introducer sheath may damage the balloon. Do not use the catheter again after withdrawing the balloon through the introducer sheath.
- To avoid balloon leakage, do not allow the balloon to contact calcified or stented arteries and do not allow the balloon to move during inflation.
- Do not use a device that has been damaged. Use of damaged devices may result in complications.
- Do not exceed maximum recommended balloon inflation volume as indicated on the label. Excess inflation volume may rupture the balloon.
- For through-lumen, do not exceed 2.07 MPa (300 psi) maximum recommended infusion pressure. Excess pressure may result in catheter rupture or tip detachment.
- If flow through the catheter becomes restricted do not attempt to clear catheter lumen by infusion. Doing so may cause the catheter to rupture, resulting in vessel trauma. Remove and replace the catheter.
- Do not steam shape the guide catheter.
- Do not pull the catheter at a force exceeding 3.5 N (0.35 kgf).

Precautions
- Store in a cool, dry dark place.
- Do not use open or damaged packages.
- Use by "Use By" date.
- Exposure to temperatures above 54°C (130°F) may damage the device and autoclave. Do not autoclave.
- Upon removal from package, inspect the device to ensure it is not damaged.
- Do not expose the device to solvents.
- Do not use clips, control media, consisting of ethyl ester of iodinated poppy-seed oil fatty acids, such as Lipiodol.
- The guidewire used with this system should incorporate a hydrophilic coating.
- Use the device in conjunction with fluoroscopic visualization and proper anti-coagulation agents.
- Do not conduct MRI inspection while the product is in place.
- Torquing the guide catheter while linked may cause damage that could result in separation of the catheter shaft.
- If a device becomes lodged in the guide catheter, or if guide catheter becomes severely kinked, withdraw the entire system (guide catheter, guidewire and introducer) as a unit.
- To prevent thrombus formation and contrast media crystal formation maintain a constant infusion of appropriate flush solution through the guide catheter lumen.

Recommended Procedure
1. Gently remove the catheter from the pouch by grasping the hub and slowly rubbing the catheter straight out of the protective tubing without pulling on the shaft. Inspect the catheter before use to verify that its size, shape and condition are suitable for the specific procedure.
2. Dilator Instructions (if applicable):
 - Remove dilator from protective tube.
 - Gently insert dilator into guide catheter.
 - Flush dilator with heparinized saline.
3. Prepare balloon inflation media by mixing contrast with saline (50% by volume is equivalent to 150 mg/mL iodine).
4. Fill a 20 mL syringe with 1 mL of balloon inflation media. Attach a 3-way stopcock to the balloon hub. Attach a 20 mL syringe to the 3-way stopcock.
5. Attach a 1 mL syringe to the balloon's 3-way stopcock. Turn the stopcock towards the 1 mL syringe.
6. With the 20 mL syringe pushed downward:
 - Pull back on the syringe plunger to aspirate the balloon lumen. Maintain negative pressure until the air bubbling from priming in the syringe
 - Release the syringe plunger to allow media to be drawn into the balloon lumen. Do not relase media.
 - Again pull back on the syringe plunger to aspirate the balloon lumen. Maintain negative pressure until air bubbles stop forming in the syringe.
 - Release the syringe plunger to allow media to be drawn into the balloon lumen. Do not relase media.
7. Turn stopcock towards the balloon hub. Transfer maximum recommended balloon inflation volume from the 20 mL syringe to the 1 mL syringe.
8. Inflate the balloon up to the recommended inflation volume. Refer to the product label. Turn the stopcock off towards the balloon hub.
9. Inspect the balloon for leakage. Wait for balloon inflated until the air bubbles diffuse from the balloon. If air remains in the balloon, it should dissipate if left in a sterile field for several minutes.
10. Deflate the balloon by turning stopcock towards the 1 mL syringe and applying vacuum with the 20 mL syringe.
11. After ensuring that the balloon is fully deflated, wet the distal shaft with saline and insert the balloon portion of the catheter into the inserter.
12. Insert the tip of the guidewire and the inserter/guide catheter assembly through the proximal valve of the introducer sheath.
13. During insertion of the catheter, when the balloon passes the sheath introducer, add flush solution to the inserter/guide catheter as needed.
14. Place the guide catheter in the selected vessel using fluoroscopy.
15. Remove dilator (if applicable) and guidewire.
16. Attach a rotating hemostasic valve (RHV) to the guide catheter through-lumen hub. Flush the through-lumen with heparinized saline. Attach a 3-way stopcock to the RHV and flush the guide catheter.
17. To inflate the balloon, transfer up to the maximum recommended balloon inflation volume from the 20 mL syringe to the 1 mL syringe and gently inflate the balloon inflation media with the 1 mL syringe until desired balloon diameter is attained.
18. Ensure the balloon is completely deflated before withdrawing the guide catheter.

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