Boston Scientific

FilterWire EZ™

Embolic Protection System

Directions for Use	**2**
Instrucciones de uso	**8**
Mode d'emploi	**15**
Gebrauchsanweisung	**22**
Istruzioni per l'uso	**29**
Gebruiksaanwijzing	**36**
Instruções de Utilização	**43**

Boston Scientific (Master Brand DFUTemplate 8in x 8in Global, 90106041AS), DFU, MB, FilterWire EZ, Global, 50607386-01A

Black (K) ΔE ≤5.0

FilterWire EZ™

Embolic Protection System

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Refer to the instructions supplied with any interventional devices to be used in conjunction with the FilterWire EZ Embolic Protection System for their intended uses, contraindications, warnings and precautions.

WARNING

Contents supplied STERILE using a Radiation process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Boston Scientific FilterWire EZ Embolic Protection System is a temporary intravascular 0.014 in (0.36 mm) guidewire filtration system that is placed distal to the target lesion to be treated by interventional procedures. The system consists of a protection wire, an EZ Delivery Sheath, an EZ Retrieval Sheath and accessories. When deployed, the protection wire's filter bag is designed to contain and remove embolic material that may be liberated during the interventional procedure. The protection wire is used as a standard 0.014 in (0.36 mm) steerable guidewire. The spring coil tip of the protection wire and the filter loop are radiopaque to enable visual guidance during placement. At the completion of the procedure, the filter is captured using the EZ Retrieval Sheath and then removed from the patient.


Figure 1. Protection Wire

Table 1. Contents

Quantity	Description
1 each	Protection Wire
1 each	EZ Delivery Sheath
1 each	EZ Retrieval Sheath
1 each	Accessory Tool Kit (One Wire Torquer, One Peel-Away Introducer, One Valve Dilator Tool)

Table 2. Available Sizes

Vessel Diameter	Protection Wire Lengths
3.5 mm - 5.5 mm	190 cm or 300 cm

Note: The 190 cm protection wire is compatible with the Boston Scientific AddWire™ Extension Wire (catalog number 22150-01). The EZ Bent Tip Retrieval Sheath (catalog number 50100-150), an alternate retrieval sheath, is available separately.

INTENDED USE/INDICATIONS FOR USE

The Boston Scientific FilterWire EZ Embolic Protection System is indicated for use as a guidewire and embolic protection system to contain and remove embolic material (thrombus/debris) while performing angioplasty and stenting procedures in the peripheral vascular system, carotid arteries, coronary arteries and saphenous vein grafts. The diameter of the vessel at the site of filter loop placement should be between 3.5 mm and 5.5 mm.

CONTRAINDICATIONS

- Patients with severe allergy to heparin.
- Patients with bleeding diathesis or other disorders that limit the use of anticoagulant therapy.

WARNINGS

- Only physicians experienced in percutaneous intravascular techniques and procedures should use the FilterWire EZ Embolic Protection System.
- The appropriate antiplatelet and anticoagulation therapy should be administered pre- and post-procedure to minimize the risk of embolism and thrombosis.
- Avoid using power injection in the cerebral circulation.
- Filter safety and effectiveness has not been tested with power injection.
- Failure to follow recommended device preparation and delivery instructions may result in air embolism or other adverse reaction.
- Introduce and advance devices slowly to prevent air embolism or trauma to the vasculature.
- Do not attempt to move the protection wire without observing the resultant tip response.
- All distal wire tips have the potential to cause vessel injury. Confirm that the wire tip is free within the vessel.
- Do not use excessive force when attempting to cross the lesion with the FilterWire EZ System.
- Observe all protection wire movement in the vessels under fluoroscopic imaging.

Boston Scientific (Master Brand DFUTemplate 8in x 8in Global, 90106041AS), DFU, MB, FilterWire EZ, Global, 50607386-01A

Black (K) ΔE ≤5.0

- Always keep the open filter loop distal to the deployed stent. Pulling the filter loop into the stent area may lead to entanglement with the stent and possible filter loop detachment.
- Ensure that the protection wire is stabilized throughout the procedure. Failure to stabilize the protection wire could lead to inadvertent movement of the filter resulting in protection wire entanglement and/or delay in the procedure.
- Do not pull excessively on the protection wire or the EZ Retrieval Sheath to avoid filter bag tears, filter loop detachment or other protection wire damage.

PRECAUTIONS

- Do not autoclave or expose the FilterWire EZ™ Embolic Protection System to organic solvents.
- Confirm that the interventional devices to be used are compatible with the 0.014 in (0.36 mm) diameter protection wire before actual use.
- For SVG applications, prior to use, heparinize the patient to achieve an Activated Clotting Time (ACT) of > 300 seconds (> 200 seconds if using GP IIb/IIIa inhibitors) or in absence of ACT measurement, administer an appropriate weight based bolus of heparin (approximately 125 units/kg).
- For carotid procedures, administer heparin and monitor ACT to maintain > 275 seconds throughout the procedure.
- Venous access should be available in order to manage bradycardia and/or hypotension by either pharmaceutical intervention or placement of a temporary pacemaker, if needed.
- Prior to use, inspect the protection wire and sheaths for damage (e.g. bends, kinks, filter tears). Do not use the FilterWire EZ System if damaged. Guidewires are delicate instruments and should be handled carefully.
- Maintain the sterile contents of the pouch (including the coils and accessory tool kit) in the sterile field during the procedure. Do not store the accessory tool kit for future use.
- Confirm angiographically that the vessel diameter is appropriate for treatment with the FilterWire EZ System and ensure that the correct device size is selected. The reference vessel diameter at the filter loop deployment site must be at least 3.5 mm but no greater than 5.5 mm for the FilterWire EZ.
- To prevent overstretching and potential dissection of the artery, ensure the device selected matches the reference vessel diameter.
- Confirm angiographically that there is sufficient vessel length to place the filter to allow for and maintain adequate distance between the protection wire's catheter stop and the stent delivery system or other compatible interventional devices (see Table 3 and Figure 2).
- Minimum inner diameter of a guide catheter or guide sheath must be 0.066 in (1.68 mm).
- Minimum inner diameter bore of a rotating hemostasis valve must be 0.075 in (1.91 mm) to accommodate use of the peel-away introducer supplied with the FilterWire EZ Embolic Protection System.
- Use of a guide catheter with a rotating hemostasis valve is recommended when performing a procedure with the FilterWire EZ Embolic Protection System.
- When using a rotating hemostasis valve, ensure that the hemostasis valve is sufficiently opened prior to insertion of any device.
- A stent that is not well apposed to the vessel wall can lead to difficulty crossing with the retrieval sheath and may increase the risk of filter entanglement with the stent.
- Ensure the stent selected matches the diameter of the vessel being treated to maximize the potential for adequate stent apposition to the vessel wall.
- Ensure adequate guide catheter or guide sheath support and position during the interventional procedure. Guide catheter or guide sheath movement can cause the filter to be inadvertently pulled proximally in the vessel.
- During the procedure, always attempt to keep the filter loop, lesion site, and radiopaque tip of the guide catheter or guide sheath in the field of view, minimizing the likelihood of the guide catheter or guide sheath backing out or prolapsing into the aortic arch.
- Prolapse of the guide catheter or guide sheath can result in movement of an open filter through an untreated lesion, filter/stent entanglement, filter basket detachment, proximal movement of the stent, or protection wire damage.
- Always advance or retract the FilterWire EZ Embolic Protection System slowly under fluoroscopic imaging. If resistance is felt and/or observed, determine the cause and take any necessary remedial action.
- Advancing or torquing the FilterWire EZ System against resistance could result in vessel trauma.
- Always tighten the wire torquer appropriately to the protection wire. Over-tightening of the wire torquer may cause difficulty/inability to remove the wire torquer from the protection wire. Under-tightening of the wire torquer may lead to improper filter deployment.
- Always advance/retract the wire torquer carefully along the protection wire to prevent kinking.
- Keep the wire position stable by using the wire torquer while the filter is deployed.
- In vitro testing for compatibility of the FilterWire EZ System with coronary drug-eluting stents (DES) indicates that the DES coating should not be compromised upon advancing or retracting the FilterWire EZ System; however, care should still be taken not to disrupt the drug coating if a coronary DES has been recently placed.
- To address lesions in subsequent vessels, use a new FilterWire EZ Embolic Protection System for each vessel.
- Do not use the FilterWire EZ System if the yellow protective housing is not clipped to the coil and the filter and spring coil tip are exposed.
- Do not use the FilterWire EZ System if the retention clip is not clipped onto the coil with the delivery sheath secured to it.
- If a buddy wire is used to facilitate crossing the lesion, ensure that the buddy wire does not wrap around the protection wire.
- If a buddy wire is used, leaving the buddy wire in place during lesion treatment can lead to vessel trauma or accidental entrapment of the protection wire.
- To avoid device damage, do not pre-dilate the lesion with the sheathed FilterWire EZ System in the vicinity of the lesion.
- During deployment, failure to position the wire torquer against the hemostasis valve and secure it to the protection wire may lead to wire kinking and/or filter movement during deployment.
- If blood flow is slowed during the procedure, the protection wire may be removed and replaced, if necessary. However, at the operator's discretion, the procedure may be completed before removing the protection wire.

Boston Scientific (Master Brand DFUTemplate 8in x 8in Global, 90106041AS), DFU, MB, FilterWire EZ, Global, 50607386-01A

Black (K) ΔE ≤5.0

- If there is no blood flow, the protection wire should be removed and replaced (unless the procedure is complete).
- During FilterWire EZ™ System retrieval, do not rotate the EZ Retrieval Sheath more than 90 degrees in either direction; this can cause the protection wire to wrap around the EZ Retrieval Sheath shaft and prevent capture of the filter.
- Always use a FilterWire EZ System compatible retrieval sheath to remove the protection wire. Pulling an open filter through a stent should be a last resort for retrieval as it may lead to entanglement with the stent and possible filter loop detachment.
- Use caution when advancing or retracting the FilterWire EZ System through a deployed stent as this may cause stent/filter entanglement or stent dislocation.
- Use of a guide sheath with a fixed hemostasis valve requires use of a valve dilator tool. The use of a guide sheath with a fixed hemostasis valve may cause the filter bag to tear at the hemostasis valve upon removal.

ADVERSE EVENTS

Possible adverse events associated with FilterWire EZ Embolic Protection System use and application procedure include, but are not limited to, the following:

- Angina
- Bleeding complications
- Bradycardia or arrhythmias, including ventricular fibrillation or tachycardia
- Congestive heart failure
- Damage to or dislocation of the implanted stent(s)
- Death
- Detachment and/or implantation of a component of the system
- Drug reactions, allergic reaction to contrast media, medications, or device materials
- Embolization of air, tissue, thrombus or other embolic debris
- Emergent Surgery
- End organ ischemia/infarction
- Headache
- Hypotension/hypertension
- Infection (local or systemic)
- Myocardial infarction
- No-reflow resulting from reduced blood flow through the FilterWire EZ System filter
- Pain
- Puncture site complications (i.e., vessel occlusion, hemorrhage, hematoma, pseudoaneurysm or arteriovenous fistula)
- Renal insufficiency, kidney failure, hematuria
- Stroke/Cerebrovascular accident (CVA), transient ischemic attack (TIA) or seizure
- Vessel damage, dissection, occlusion, aneurysm, perforation, rupture, injury, thrombosis, or spasm

The FilterWire EZ Embolic Protection System has been tested for compatibility with drug eluting stents (DES). The coating of the DES is not compromised when used in conjunction with this embolic protection system.

HOW SUPPLIED

- Do not use if package is opened or damaged.
- Do not use if labeling is incomplete or illegible.

Handling and Storage

- Store in a cool, dry, dark place.
- Use the device prior to the "Use By" date noted on the box and pouch.

INSTRUCTIONS FOR USE

Preparation of the Protection Wire and EZ Delivery Sheath:

1. Open the pouch using standard sterile handling procedures and place the packaging coils and accessory tool kit into the sterile field.
2. Unclip and carefully remove the yellow protective housing from the filter.

Note: The protection wire is preloaded into the EZ Delivery Sheath.

Caution: Do not use the FilterWire EZ System if the yellow protective housing is not clipped to the coil and the filter and spring coil tip are exposed.

3. Remove the preloaded protection wire from the retaining clip; advance the EZ Delivery Sheath until the clear section is exposed; then, grasp the clear section of the sheath, and remove the preloaded protection wire from the packaging coil.

Caution: Do not use the FilterWire EZ System if the retention clip is not clipped onto the coil with the delivery sheath secured to it.

4. Retain the packaging coil on the sterile field to ensure availability of the EZ Retrieval Sheath when needed.
5. Attach the wire torquer to the protection wire near the exit port. This will facilitate retraction of the filter into the EZ Delivery Sheath during preparation.
6. Grasp the EZ Delivery Sheath towards the very distal end of the clear section and submerge the filter and sheath in heparinized saline.
7. While the filter and delivery sheath are submerged in the heparinized saline, sheath the protection wire by slowly retracting it into the EZ Delivery Sheath until the nose cone is partially retracted into the sheath.

Note: Do not over-retract the filter.

Note: Do not sheath the filter until immediately prior to use. Make certain that the filter is thoroughly saturated with saline prior to retracting the filter into the EZ Delivery Sheath.

Insertion and Positioning of the FilterWire EZ System:

Warning: Introduce and advance devices slowly to prevent air embolism or trauma to the vasculature.

1. Carefully insert the sheathed protection wire tip into the peel-away introducer. If desired, shape the spring coil tip of the protection wire. In order to prevent spring coil tip damage, ensure the tip is retracted inside the peel-away introducer.
2. Advance the FilterWire EZ System and peel-away introducer assembly into the hemostasis valve attached to a guide catheter or guide sheath.
3. Carefully remove the peel-away introducer.
4. Advance the FilterWire EZ System through the guide catheter or guide sheath.

Boston Scientific (Master Brand DFUTemplate 8in x 8in Global, 90106041AS), DFU, MB, FilterWire EZ, Global, 50607386-01A

Black (K) ΔE ≤5.0

Note: Carefully advance both the EZ Delivery Sheath and protection wire together to prevent inadvertent deployment of the filter.

Warning: Do not attempt to move the protection wire without observing the resultant tip response.

Warning: All distal wire tips have the potential to cause vessel injury. Confirm that the wire tip is free within the vessel.

5. Under fluoroscopic imaging, steer the FilterWire EZ™ System into the target vessel. Using a two-handed technique, torque the protection wire with one hand and advance the EZ Delivery Sheath with the other hand.

Note: Always advance the protection wire and EZ Delivery Sheath together to prevent inadvertent deployment.

Note: Do not allow the protection wire to wrap around the EZ Delivery Sheath shaft to ensure the filter can be deployed smoothly.

Note: Consider using a "buddy wire" or pre-dilate the lesion if it is tight and difficult to pass.

Caution: If a buddy wire is used, ensure that the buddy wire does not wrap around the protection wire.

Caution: If a buddy wire is used, leaving the buddy wire in place during lesion treatment can lead to vessel trauma or accidental entrapment of the protection wire.

Caution: To avoid device damage, do not pre-dilate the lesion with the sheathed FilterWire EZ System in the vicinity of the lesion.

6. Advance the FilterWire EZ System across the lesion until the apex of the filter loop can be deployed in the recommended minimum landing zone (see Table 3 and Figure 2). This distance should be maintained throughout the procedure.

Warning: Do not use excessive force when attempting to cross the lesion with the FilterWire EZ System.

Note: If extra landing zone length is available, placing the filter loop beyond the minimum distance is desirable.

Table 3. Minimum Landing Zone Requirements

	Minimum Landing Zone should be:
FilterWire EZ System (3.5 mm - 5.5 mm)	≥ 3.0 cm


Figure 2. Minimum Landing Zone

Note: If an interventional device with a tip-to-shoulder length greater than 10 mm is to be used, a longer vessel length and landing zone may be required to place the protection wire in order to avoid contact between the protection wire's catheter stop and the tip of the interventional device.

7. Once the protection wire is advanced past the lesion, slide the wire torquer along the protection wire and secure it against the hemostasis valve.
8. Deploy the filter by holding the protection wire in place with the wire torquer pressed against the hemostasis valve while simultaneously retracting the EZ Delivery Sheath.

Caution: Failure to position the wire torquer against the hemostasis valve and secure it to the protection wire may lead to wire kinking and/or filter movement during deployment.

9. Continue retracting the EZ Delivery Sheath, peeling the sheath away until the clear section of the EZ Delivery Sheath is encountered. Remove the wire torquer and carefully remove the remaining distal segment of the EZ Delivery Sheath from the wire while maintaining protection wire stability.
10. Inject contrast media to verify the protection wire is in the proper position and that there is adequate flow.

Note: If repositioning is required, refer to "Repositioning the Filter Using the EZ Retrieval Sheath" section.

During the Procedure:

1. Use the protection wire like a traditional guidewire to track catheters and stent delivery systems to the target treatment site. Carefully maintain the position of the filter, particularly during device exchanges.

Warning: Observe all protection wire movement in the vessels under fluoroscopic imaging.

Warning: Always keep the open filter loop distal to the deployed stent. Pulling the filter loop into the stent area may lead to entanglement with the stent and possible filter loop detachment.

Warning: Ensure that the protection wire is stabilized throughout the procedure. Failure to stabilize the protection wire could lead to inadvertent movement of the filter, resulting in protection wire entanglement and/or delay in the procedure.

Boston Scientific (Master Brand DFUTemplate 8in x 8in Global, 90106041AS), DFU, MB, FilterWire EZ, Global, 50607386-01A

Black (K) ΔE ≤5.0

Note: Inject contrast after any intervention or exchange. Visually verify that blood flow is not obstructed and that filter loop/vessel wall apposition is maintained. Reposition if necessary.

Caution: If blood flow is slowed during the procedure, the protection wire may be removed and replaced. However, at the operator's discretion, the procedure may be completed before removing the protection wire.

Caution: If there is no blood flow, the protection wire should be removed and replaced if further intervention is required.

2. Carefully remove all interventional devices, ensuring that the protection wire is not kinked and its position is maintained during removal.

Preparation of the EZ Retrieval Sheath and Filter Capture:

1. Remove the EZ Retrieval Sheath from its packaging coil by grasping it at the white proximal handle and pulling it carefully out of the packaging coil, taking care not to damage the EZ Retrieval Sheath.
2. Flush the EZ Retrieval Sheath with heparinized saline. Be careful not to damage the EZ Retrieval Sheath tip with the syringe.

Note: If the EZ Retrieval Sheath was used for repositioning, flush it again with heparinized saline and clean the outer surface before using it for filter retrieval.

Warning: Observe all protection wire movement in the vessels under fluoroscopic imaging.

3. While maintaining the protection wire position distal to the lesion, advance the EZ Retrieval Sheath over the protection wire past any deployed stent(s) until the tip of the EZ Retrieval Sheath reaches the protection wire's catheter stop.

Caution: Do not rotate the EZ Retrieval Sheath more than 90 degrees in either direction; this can cause the protection wire to wrap around the EZ Retrieval Sheath shaft and prevent capture of the filter.

Caution: Always use a FilterWire EZ™ System compatible retrieval sheath to remove the protection wire. Pulling an open filter through a stent should be a last resort for retrieval as it may lead to entanglement with the stent and possible filter loop detachment.

Note: If it is difficult to cross a deployed stent or tortuous vessel anatomy with the EZ Retrieval Sheath, the EZ Bent Tip Retrieval Sheath (available separately) may be used.

4. Slide the wire torquer along the protection wire and secure it against the hemostasis valve. Gently and slowly retract the protection wire and filter loop back into the EZ Retrieval Sheath until resistance is felt.

Warning: Do not pull excessively on the protection wire to avoid filter bag tears, filter loop detachment, or other protection wire damage.

Note: The distal edge of the collapsed filter loop should align with, or be proximal to, the EZ Retrieval Sheath marker band (see Figure 3).


Figure 3. Retrieved Filter Position inside the EZ Retrieval Sheath

Note: If the distal edge of the collapsed filter loop is not aligned with, or proximal to, the EZ Retrieval Sheath marker band, it is possible that the filter is overfilled. Carefully retract the protection wire and EZ Retrieval Sheath together as a system.

Note: If repositioning is required, use the EZ Retrieval Sheath. Refer to "Repositioning the Filter Using the EZ Retrieval Sheath" section for instructions.

Removal of the Protection Wire Using the EZ Retrieval Sheath:

1. Slowly and carefully retract the entire system as a unit until the tip of the EZ Retrieval Sheath is adjacent to the tip of the guide catheter or guide sheath.

Caution: Use caution when advancing or retracting the FilterWire EZ System through a deployed stent as this may cause filter/stent entanglement or stent dislocation.

Warning: Do not pull excessively on the protection wire or the EZ Retrieval Sheath to avoid filter membrane tears, filter loop detachment, or other protection wire damage.

Note: If any resistance is met during retraction of the sheathed protection wire and filter, slightly advance the sheathed protection wire and rotate the EZ Retrieval Sheath before continuing to retract.

2. Carefully retract the FilterWire EZ System across the guide catheter or guide sheath up to the hemostasis valve.

Note: If there is any resistance at the guide catheter or guide sheath, retract the guide catheter or guide sheath and FilterWire EZ System together.

3. Remove the FilterWire EZ System from the patient, ensuring that the hemostasis valve is fully opened prior to FilterWire EZ System removal. If a fixed hemostasis valve is used, place the valve dilator tool on the shaft of the retrieval sheath and slide the tool into the valve. Then, remove the FilterWire EZ System from the patient through the valve dilator tool. Remove the valve dilator tool.

Caution: Use of a guide sheath with a fixed hemostasis valve requires use of a valve dilator tool. The use of a guide sheath with a fixed hemostasis valve may cause the filter bag to tear at the hemostasis valve upon removal.

Repositioning the Filter Using the EZ Retrieval Sheath:

1. Prepare the EZ Retrieval Sheath and capture the filter according to the instructions outlined in the "Preparation of the EZ Retrieval Sheath and Filter Capture" section.

Note: Observe all protection wire movement in the vessels under fluoroscopic imaging.

Boston Scientific (Master Brand DFUTemplate 8in x 8in Global, 90106041AS), DFU, MB, FilterWire EZ, Global, 50607386-01A

Black (K) ΔE ≤5.0

2. Guide the EZ Retrieval Sheath and protection wire assembly to the desired location. Make sure to advance both the EZ Retrieval Sheath and the protection wire as a system to ensure that the protection wire does not prematurely deploy.
3. If needed, the spring coil tip of the protection wire can be steered by turning the wire torquer previously installed on the protection wire.

Note: Do not allow the protection wire to wrap around the EZ Retrieval Sheath shaft to ensure that the filter can be deployed smoothly.

Note: If it is difficult to advance or reposition the protection wire, this may be the result of a full filter. Care must be taken not to damage the filter.

4. Advance the FilterWire EZ™ System until the apex of the filter loop can be deployed in the recommended minimum landing zone (see Table 3 and Figure 2). This distance should be maintained throughout the procedure.

Note: Deploy the filter as distal as possible to the lesion.

5. Once the protection wire is in the desired location, slide the wire torquer along the protection wire and secure it against the hemostasis valve.
6. Deploy the filter by holding the protection wire in place with the wire torquer pressed against the hemostasis valve while simultaneously retracting the EZ Retrieval Sheath.
7. Ensure that the filter has been deployed in the desired location. Remove the wire torquer and the EZ Retrieval Sheath from the patient while holding the protection wire in place.
8. Inspect the EZ Retrieval Sheath for damage. If no damage is noted, set it aside for eventual retrieval and removal of the protection wire.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose**. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments**.

Boston Scientific (Master Brand DFUTemplate 8in x 8in Global, 90106041AS), DFU, MB, FilterWire EZ, Global, 50607386-01A

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© 2017 Boston Scientific Corporation or its affiliates. All rights reserved.