BONA博納 SHENZHEN BONA PHARMA TECHNOLOGY CO., LTD

Bona smart plant, No. 6 Gaoxinyuan Guansheng 1st Road, Luhu community, Guanhu Subdistrict,
Longhua District, Shenzhen city, Guangdong province, P.R. China

Tel: +8675526994208 Fax: +8675526527737 www.bona-cn.com

INSTRUCTIONS FOR USING

1. Item:

Vaginal applicator in individual pack

2. Description:

The item includes: applicator and sealed poly bag.

3. How to use: (See below picture)

Step 1 Open the sealed poly bag.

Step 2 Screw on the applicator on fitting bottle or container which filled liquid medicine.

Step 3 Gently insert the applicator into the vagina as for as it will go comfortably.

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BONA博納 SHENZHEN BONA PHARMA TECHNOLOGY CO., LTD
Bona smart plant, No. 6 Gaoxinyuan Guansheng 1st Road, Luhu community, Guanhu Subdistrict,
Longhua District, Shenzhen city, Guangdong province, P.R. China
Tel: +8675526994208 Fax: +8675526527737 www.bona-cn.com

SPECIFICATION OF VAGINAL APPLICATOR

1. Item Description

Vaginal applicator in individual pack is a dispenser for vaginal medicine or rinser, include only one component.

2. Structure Features

With holes at the top of applicator, to dispensing liquid deeply into vaginal. Produced and packed in class 100,000 clean room.

3. Appearance

As per drawing

4. Component and raw material list

Component name	Material name	Material Grade
Applicator	PP	T03 (Material specificate as per attached)

5. Packaging and stock indication

Individual packing in poly bag and sealed, then pack 700pcs in two bigger poly bags and seal, and put in cardboard box and seal.

Suggest to stock in the place which should be clean, dry, airy and shade, and the validity period is 2 years before you use them to pack your goods. Suggested temperature of stock place is between 0℃ to 40℃ and humidity is between 35% with 95%.

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