MOZEC™ PTA — Balloon Dilatation Catheter — Instructions for use - English

1. **Device Name:**

The device brand name is Mozec™ PTA, the generic name is Percutaneous Transluminal Angioplasty (PTA) Balloon Dilatation Catheter.

2. **Device Description:**

Mozec™ PTA OTW Balloon Dilatation Catheter consists of a Balloon (dilatation element) near the distal tip with radiopaque marker bands, which represent the approximate balloon working length at nominal pressure. The proximal end has a Y luer hub for balloon inflation.

The catheter has a soft tip at the distal end. Dual lumens shaft which permit movement of guide wire and balloon inflation. The distal portion of the shaft has hydrophilic coating applied on it, in order to facilitate catheter advancing through the vascular district and the vessel stenosis.

The compliance chart on the inner and outer label indicates how the balloon diameter increases with increasing pressure. The compliance data is based on in vitro testing of balloons at 37°C temperature.

- Guide wire compatibility: 0.035"
- 0.035" OTW PTA Balloon Dilatation Catheter available with two usable length viz. 80 cm and 135 cm.

The Mozec™ PTA Balloon Dilatation Catheter is available in different balloon sizes. Nominal balloon diameters and lengths are printed on the hub. Refer package label for information on catheter for usable length, Nominal and Rated burst pressure and Sheath compatibility.

3. **Indications:**

The Mozec™ PTA Balloon Dilatation Catheter is intended to dilate stenosis in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteiovenous dialysis fistula.This device is also indicated for stent post-dilatation in the peripheral vasculature.

4. **Contraindications:**

The Mozec™ PTA Balloon Dilatation Catheter is contraindicated for the use in coronary arteries and neuro vasculature.

5. **Warnings:**

- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death.
- Reuse, Reprocessing or re-sterilization may also create a risk of contamination of the device and/ or cause patient infection or cross-infection, including, but not limited to the transmission of infectious disease (s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Percutaneous Transluminal Angioplasty should only be performed at hospitals where emergency artery bypass graft surgery can be quicklyperformed in the event of a potentially injurious or life-threatening complication.
- Use the product before the "Use By" date specified on the product label.
- Introducer sheath used must have lumen sizes that are suitable to accommodate the introduction of Mozec™ PTA Balloon Dilatation Catheter.
- Do not use if the inner package is opened or damaged. Carefully remove the PTA catheter from the pouch to prevent damage and premature removal of the balloon cover.
- Use only the appropriate Balloon inflation media. Do not use air/any other gaseous media/oil based inflation media/alcohol/organic solvents to inflate the PTA catheter as this may cause uneven expansion/catheter leakage/loss of lubrication.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract catheter unless the balloon is fully deflated under vacuum. If resistance is felt during manipulation, determine the cause of the resistance before proceeding.
- Balloon pressure should not exceed the Rated Burst Pressure (See Compliance Chart on the packing labels). Use of a pressure monitoring device is recommended to prevent over pressurization.
- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.

6. **Precautions:**

- Prior to angioplasty, examine the PTA Balloon Dilatation catheter to verify functionality. Ensure that its size and shape are suitable for the specific procedure for which it is to be used.
- Only the physicians thoroughly trained and educated in the performance of PTA should use this device. Physicians should keep themselves informed and updated on recent publications about PTA techniques.
- The balloon catheter should be used with caution for procedure involving calcified lesions due to the abrasive nature of these lesions.
- Precautions to prevent or reduce clotting should be taken when any catheter is used.
- Flush or rinse all products entering in the vascular system sterile isotonic saline or similar solution via guidewire access port prior to use. Consider the use of systemic heparinization.
- Prior to reinsertion or withdrawal of the PTA catheter, wipe the guide wire with saline-soaked gauze to remove blood or other residues.
- During the procedure, appropriate anti-coagulant therapy must be provided to the patient as needed, Anti-coagulant therapy should be continued for a period of time to be determined by the physician after the procedure.
- The catheter must always be introduced, moved and withdrawn over a guidewire.
- Never try to move the guide wire when the balloon is inflated.
- The size of the inflated balloon should be selected not to exceed the diameter of the artery immediately distal or proximal to the stenosis.
- After use, eliminate the product according to safety requirements related to the products contaminated by blood.

7. **Adverse Effects:**

Potential adverse events, which may be associated with the use of a PTA Balloon Dilatation Catheter, include but are not limited to:

- Spasm reaction of vessel
- Artery wall injuries of different degrees, including perforation, dissection, rupture.
- Restenosis of the dilated vessel
- Haemorrhage or haematoma
- Arteriovenous fistula
- Hypotension, Arrhythmias
- Pseudo-aneurysm
- Vascular thrombosis, systemic embolization
- Device embolization
- Pain and tenderness
- Sepsis, Infection
- Endocarditis
- Pyrogenic reaction
- Additional interventional or surgical procedure
- Amputation
- Death

Complications related to concomitant medication, e.g.:

- Drug reaction
- Bleeding from anticoagulant/antiplatelet medication
- Allergic reactions to contrast medium

8. **How Supplied:**

Sterile : This device is sterilized with Ethylene Oxide, Non pyrogenic.

Contents: One (1) Mozec™ PTA Balloon Dilatation Catheter housed in a blister tray , One (1) instructions for use.

Storage : Store in a cool, dry and dark place in its original packaging.

9. **Operator's Manual:**

The following instructions provide guidance but do not obviate the necessity of training in the use of the device.

- The PTA procedure should be carried out according to standard PTA guidelines.
- Read Instructions for Use.

9.1 Preparation:

a. Open product box to remove the inner Tyvek pouch containing blister tray with sterile Mozec™ PTA Balloon Dilatation Catheter.

b. Remove protective sheath and stylet from the tip.

c. Prior to using Mozec™ PTA Balloon Dilatation Catheter, carefully inspect for bends, kinks, and other damages. Ensure that the radiopaque marker(s) are placed in the balloon.

d. Flush the guidewire lumen of the Catheter through guide wire port of the Y luer hub with sterile heparinised solution by gently applying pressure with thesyringe for at least 10 seconds. The saline should exit from the guidewire port. Remove the syringe while flushing the needle.

e. Attach a stopcock to the Catheter's inflation hub.

f. Attach an appropriate size of syringe filled with a minimum of 3cc of contrast medium in sterile heparinised saline to the stopcock.

g. Open the stopcock and induce negative pressure, pulling the syringe plunger as far back as possible without dislodging it from the syringe barrel. Aspirate for 15 seconds till the air is fully evacuated.

h. Gently release the negative pressure to fill the Catheter lumen with the solution.

i. Close the stopcock, remove the syringe and purge the air.

j. Repeat steps (g) through (j) if, necessary till there are no air bubbles left.

k. Fill an Angioplasty Inflation Device with a minimum of 3cc of contrast medium in sterile heparinised saline and attach it to the purged Catheter.

Note: Purge the Angioplasty Inflation Device with air prior to connecting to the Catheter.

l. Induce negative pressure and close the stopcock. Purge any air from the Angioplasty Inflation Device through the stopcock. Leave the Catheter under negative pressure until it is ready to be used.

9.2 Assembly and Insertion Procedure for PTA:

a. Position the guide wire under fluoroscopic guidance in accordance with standard percutaneous transluminal angioplasty(PTA) .

b. Insert the distal tip of the dilatation catheter on to the proximal end of the guide wire and advance until the guide wire exits the proximal end of the catheter.

c. Carefully insert the catheter through the Introducer Sheath .

d. Advance the catheter over the guide wire towards the lesion.

e. Position the balloon across the lesion using the balloon radiopaque markers as reference points under fluoroscopy.

f. If resistance is encountered during advance, do not force the passage. Determine the cause of resistance and take corrective actions.

Caution:

1. During the advance, ensure that negative pressure is maintained in the Inflation Device to keep the balloon fully deflated.
2. Do not advance or withdraw the PTA Catheter within the vasculature unless it is preceded by a guide wire.

9.3 Inflation:

a. Dilate the stenotic lesion following standard PTA techniques.

Caution: Use of oversized balloon may cause dissection. Do not exceed the Rated Burst Pressure mentioned on the product labels.

b. If needed, carry out repeat inflations. Maintain negative pressure between two inflations.

c. After each subsequent inflation, determine the distal blood flow.

d. If significant stenosis persists, the lesion may require successive inflations/further treatment.

e. Determine the results under fluoroscopy.

9.4 Catheter Withdrawal and Disassembly Procedure:

a. Apply negative pressure to the inflation device to achieve complete deflation of the Balloon. Confirm this under fluoroscopy.

b. Withdraw the deflated PTA Catheter and guide wire from the Introducer sheath.

c. Remove Sheath from the vasculature following an accepted technique for removal.

Note: If the balloon cannot be withdrawn through Introducer sheath, withdraw both as a single unit.

9.5 Disposal:

After usage, dispose off the hardware in accordance with accepted regulations for medical waste management.

Note: Updated guidelines for the procedure have to be consulted through new publications.

10. **Disclaimer of Warranty and Limitation of Remedy:**

There is no express or implied warranty, including without limitation any implied warranty of merchantability of fitness for a particular purpose, on the Meril Life Sciences Pvt. Ltd. product(s) described in this publication. Under no circumstances shall Meril Life Sciences Pvt. Ltd. be liable for any direct, indirect, incidental or consequential damages. Some jurisdictions do not allow the exclusion of or limitations on an implied warranty. For jurisdictions which do not allow the exclusion or limitations of incidental or consequential damages, some of the above exclusions may not apply. The patient may also have other rights, which vary from jurisdiction to jurisdiction.

Descriptions or specifications in Meril Life Sciences Pvt. Ltd. printed matter, including this publications, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Meril Life Sciences Pvt. Ltd. will not be responsible for any direct, indirect, incidental, or consequential damages resulting from reuse of the product.

Symbols used in labeling

Balloon outer diameter | Balloon length | Contains one unit
Do not use if box open or damaged | Max. guide wire diameter | Keep dry
Consult instruction for use | Reference number | Lot number
For single use only do not reuse | Manufacturer | Manufacturing date
STERILE EO — Sterilized using ethylene oxide | Do not resterilize | Use by
EC REP — Authorized representative in the European community | | Non-Pyrogenic

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