SHENZHEN XINJITAI SPRAYER AND PUMP CO., LTD.
CM2, TANGWEI INDUSTRIAL ZONE, INDUSTRIAL COMPANY, GONGMING STREET,
GUANGMING NEW DISTRICT, SHENZHEN, GUANGDONG, CHINA.
TEL: +86-755-26728939 FAX: +86-755-26728819

Instruction for Usage of Products

1. General Information

Manufacturer:
Name: Shenzhen Xinjitai Sprayer & Pump Co., Ltd.
Address: No CM2, Tangwei Industrial Zone, Industrial Company, Gongming Street, Guangming New District, Shenzhen, Guangdong, China. 518107.
Tel: 0086- 755- 26984062 Fax: 0086-755 - 26728819

Product: Fine Mist Spray Bottle including Nasal Spray Bottle, Throat Spray Bottle, Topical Spray Bottle

Application: pharmaceutical fine mist spray for nose, throat or mouth, skin, etc.

Declaration of conformance:
We declare that the materials used to produce this product are approved materials and comply with FDA 21 CFR 177.1520, FDA 21 CFR 177.1640, and FDA 21 CFR 175.300(c)(3). All products meet the intend use. All supporting documentation are retained under the premise of manufacturer.

2. Instruction of Usage

The fine mist spray bottle is for the packaging of pharmaceutical fine mist spray. Its expiry date is not more than 3 years.

The product is not sterile. If sterilized, the valid period of sterilization is 2 years.

The product is for single use;

Store in a place where it avoids high temperature, high humidity, and does not have caustic gases, and should have good dryness and ventilation;

Do not use if package or product has been damaged or contaminated.

Manufacturing environment: Filling must be operated under controllable environment to ensure number of bacteria, dust and other matters are under the allowed standard in customer's production.

Manufacturing operation: Filling procedure should follow the operation instruction of the filling machine used.

Using the Final Medical Liquid: To use a new bottle, press the actuator to spray to a tissue or cotton a few time until the dosage reach the normal range. Then press the actuator to spray the applied area with the dosage determined by the customer pharmaceutical company.

General Manager: Cao Jinhui
Signature & Date: Stamp:
Cao Jinhui
3.2.2020

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SHENXHEN XINJITAI SPRAYER AND PUMP CO., LTD
CM2, Tangwei Industrial Zone, Industrial Company, Gongming Street,
Guangming New District, Shenzhen, Guangdong, China. 518107.
Tel: +86-755-26728939 Fax: +86-755-26728819

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