Boston Scientific

NaviGuide™

Percutaneous Access Needle

Directions for Use	**2**
Instrucciones de uso	**6**
Mode d'emploi	**10**
Gebrauchsanweisung	**14**
Istruzioni per l'uso	**18**
Gebruiksaanwijzing	**22**
Instruções de Utilização	**26**

50440648-01

2019-10

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AP), DFU, MB, NaviGuide, Global, 50440648-01B

Black (K) ΔE ≤5.0

TABLE OF CONTENTS

WARNING **3**

DEVICE DESCRIPTION **3**

INTENDED USE **3**

INDICATIONS FOR USE **3**

CONTRAINDICATIONS **3**

ADVERSE EVENTS **3**

PRIOR TO USE **3**

OPERATIONAL INSTRUCTIONS **4**

WARNING: **4**

WARNING: **4**

HOW SUPPLIED **5**

HANDLING AND STORAGE **5**

WARNING **5**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AP), DFU, MB, NaviGuide, Global, 50440648-01B

Black (K) ΔE ≤5.0

NaviGuide™

Percutaneous Access Needle

℞ ONLY

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DEVICE DESCRIPTION

The NaviGuide™ Percutaneous Access Needle is comprised of a hollow outer cannula and conically grouped Tri-Facet trocar with a locking hub. Some configurations contain an adjustable slider with a radiopaque ring.

Contents: (1) 18 ga. (1.3 mm) Echogenic Needle with Trocar and slider (if applicable).

Image of device indicating slider


Figure 1.

INTENDED USE/INDICATIONS FOR USE

The Percutaneous Access Needle is intended to facilitate access into the renal pelvis, and to advance guidewires into the collecting system in an antegrade manner.

CONTRAINDICATIONS

None known.

ADVERSE EVENTS

The following are some of the complications which can result from the use of access needles:

- Renal perforation
- Ureteral perforation
- Hemorrhage
- Edema
- Tissue Damage
- Infection

PRIOR TO USE

The Percutaneous Access Needle is supplied sterile in a sealed pouch and intended for single use only. Prior to use, inspect both the sterile seal and device for any damage.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AP), DFU, MB, NaviGuide, Global, 50440648-01B

Black (K) ΔE ≤5.0

1. Remove the protective sheath from the needle with hub or slider assembly.

Caution: Visually inspect the device for loose, bent or broken parts, cracks, or other abnormalities. Inspect the needle for any kinks or other damage. If an abnormality is detected, DO NOT USE. Broken parts, cracks, or kinks will hinder the mechanical operation of the Percutaneous Access Needle. If the device fails to operate properly in any way or shows any signs of damage, DO NOT USE IT. Call Boston Scientific Customer Service and return the product. If there is breach to the package seal or if the device is damaged, DO NOT USE.

OPERATIONAL INSTRUCTIONS

1. The Percutaneous Access Needle is supplied in two configurations, with or without a slider, and may be used in the following manners under fluoroscopic or ultrasound guidance:
 A. If using a slider, adjust the slider to the desired location.
 B. Select the optimal puncture site on the patient. If using a slider, the radiopaque ring may assist targeting when performing the bulls-eye technique using fluoroscopy.

Note: When grasping the device with forceps, ensure that appropriate grip is achieved.

 C. Insert the fully-assembled Percutaneous Access Needle into the appropriate area of the pelvicalyceal system. Entry may be confirmed using imaging and by removing the inner stylet to either aspirate or visualize urine from the renal pelvis.

Precaution: Before removing the inner stylet, comfirm the stylet hub is fully unlocked by turning the hub a ¼ counterclockwise turn (as denoted by the visual cue on the top of the stylet hub) to ensure needle access is not inadvertently compromised.

Precaution: If urine cannot be confirmed, STOP, determine fluoroscopic position and reintroduce assembled access needle if needed.

 D. Guidewire advancement can now be attempted through the cannula lumen.

Precaution: If resistance is encountered during advancement or withdrawal, STOP. DO NOT continue without first determining cause of resistance and taking remedial action.

WARNING:

If the guidewire must be withdrawn while the cannula is inserted, remove both cannula and wire as a unit to prevent the cannula from damaging the guidewire.

WARNING:

Use of this device should be restricted to use by or under the supervision of physicians trained in urologic and percutaneous access procedures. Care should be exercised to prevent accidental advancement into other non-intended puncture areas.

The Percutaneous Access Needle is for single patient use and has no user serviceable parts.

- Do not attempt to repair or to alter the Percutaneous Access Needle.
- Do not resterilize the device.
- Any use of this device, other than those indicated in these instructions is not recommended.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AP), DFU, MB, NaviGuide, Global, 50440648-01B

Black (K) ΔE ≤5.0

HOW SUPPLIED

The NaviGuide™ Percutaneous Access Needle is supplied as a sterile, single-use, fully-assembled device. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

HANDLING AND STORAGE

Store in a cool, dry, dark place. Rotate inventory so that product is used prior to the expiration date on package label.

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DEVICE DESCRIPTION

The Percutaneous Access Set is a single-use, sterile, disposable set used for percutaneous nephrostomy.

This Directions For Use is supplied with products using one of the following contents listings:

Percutaneous Access Set Contents Listing:

(1) 21 ga (0.81 mm) Trocar Point Needle

(1) 19 ga (1.1 mm) Trocar Point Needle with 5F (1.7 mm) PTFE Sheath

(6) Fascial Dilators [6F (2.0 mm) -16F (5.3 mm) Even French sizes]

(1) Catheter Retention Device

(1) 0.038 in (0.97 mm) x 150 cm "J" Tip Stiff Bodied PTFE Coated Guidewire

(1) Connecting Tube with Removable Stopcock

(1) #11 Disposable Scalpel

Percutaneous Access Set with 8F Pigtail Catheter Contents Listing:

(1) 8F (2.7 mm) Percutaneous Pigtail Catheter

(1) 21 ga (0.81 mm) Trocar Point Needle

(1) 19 ga (1.1 mm) Trocar Point Needle with 5F (1.7 mm) PTFE Sheath

(3) Fascial Dilators [6F (2.0 mm) -10F (3.3 mm) Even Sizes]

(1) Catheter Retention Device

(1) 0.038 in (0.97 mm) x 150 cm "J" Tip Stiff Bodied PTFE Coated Guidewire

(1) Connecting Tube with Removable Stopcock

(1) #11 Disposable Scalpel

Percutaneous Access Set with 10F Pigtail Catheter Contents Listing:

(1) 10F (3.3 mm) Percutaneous Pigtail Catheter

(1) 21 ga (0.81 mm) Trocar Point Needle

(1) 19 ga (1.1 mm) Trocar Point Needle with 5F (1.7 mm) PTFE Sheath

(4) Fascial Dilators [6F (2.0 mm) -12F (4.0 mm) Even Sizes]

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 AS), DFU, MB, Percutaneous Access Set, Global, 50549799-01B

Black (K) ΔE ≤5.0

(1) Catheter Retention Device

(1) 0.038 in (0.97 mm) x 150 cm "J" Tip Stiff Bodied PTFE Coated Guidewire

(1) Connecting Tube with Removable Stopcock

(1) #11 Disposable Scalpel

Percutaneous Access Set with 14F Malecot Catheter Contents Listing:

(1) 14F (4.7 mm) Malecot Nephrostomy Catheter with Stylet

(1) 21 ga (0.81 mm) Trocar Point Needle

(1) 19 ga (1.1 mm) Trocar Point Needle with 5F (1.7 mm) PTFE Sheath

(6) Fascial Dilators [6F (2.0 mm) -16F (5.3 mm) Even Sizes]

(1) Catheter Retention Device

(1) 0.038 in (0.97 mm) x 150 cm "J" Tip Stiff Bodied PTFE Coated Guidewire

(1) Connecting Tube with Removable Stopcock

(1) #11 Disposable Scalpel

(1) Removable Stylet

INTENDED USE/INDICATIONS FOR USE

Percutaneous Access Set is intended to establish percutaneous nephrostomy access and is indicated for the following:

Diagnostic Indications:

- Antegrade pyelography
- Pressure/perfusion study (Whitaker test)

Therapeutic Indications:

- Nephrostomy catheter drainage
- Perfusion chemolysis of renal stones
- Post-percutaneous nephrolithotomy
- Post-percutaneous resection and coagulation of urothelial tumors

CONTRAINDICATIONS

The Percutaneous Access Set is not indicated for use in patients with blood clotting anomalies due to coagulopathies or pharmacological anticoagulations.

WARNINGS

None known.

PRECAUTIONS

1. Bending or kinking during, prior to, or after placement (including external fixation) could damage the integrity of the catheter.
2. If resistance is encountered during advancement or withdrawal of the catheter, STOP. Do not continue without first determining the cause of the resistance and taking remedial action.

Note: Where long-term use is indicated, it is recommended that indwelling time for catheter not exceed 90 days*. This catheter should be evaluated by the physician on or before 90 days postplacement.

*Biocompatibility data on file.

3. Periodic examinations are recommended to assess catheter patency and function.
4. Catheters are not intended to be permanent implant devices.
5. These recommendations are meant to serve only as a basic guide to the utilization of this device. The performance of percutaneous nephrostomy should not be undertaken without comprehensive knowledge of the indications, techniques and risks of the procedure.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 AS), DFU, MB, Percutaneous Access Set, Global, 50549799-01B

Black (K) ΔE ≤5.0

ADVERSE EVENTS

- Hemorrhage
- Hematoma
- Fistula
- Edema
- Extravasation
- Loss of renal function
- Perforation
- Infection
- Encrustation
- Pneumothorax
- Catheter occlusion
- Catheter dislodgement
- Colon injury
- Tissue damage
- Inflammation
- Discomfort / pain
- Hydronephrosis

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Handling and Storage

Store product at controlled room temperature. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that products are used prior to the sterilization expiration date on the package label.

OPERATIONAL INSTRUCTIONS

Prior to Use

Carefully read all instructions **prior** to using this device.


RECOMMENDED PROCEDURE

1. The intrarenal collecting system must be localized appropriately for definitive puncture. Localization puncture for contrast injection, using the 21 gauge (0.81 mm) needle, is an optional step.
2. Select the optimal puncture site and make a small skin incision with the scalpel supplied.
3. Insert the 19 gauge (1.1 mm) (brown hub) needle/directional sheath combination into the appropriate area of the pelvocalyceal system. Entry is confirmed by flow or aspiration. Maintain intraluminal position by fractional advancement of the sheath as the needle component is removed.
4. Use the inserter to pass the flexible "J" tip of the PTFE coated guidewire into the sheath hub. Advance the wire into the collecting system and down ureter, if required, (rotating the directional sheath to control the course of the wire).
5. If the guidewire is inserted directly through a metal cannula/needle and needs to be withdrawn, remove the guidewire and metal cannula/needle as a unit. Failure to do so may result in destruction and/or separation of the outer polymer jacket requiring retrieval. If a needle is used for initial placement, a plastic entry needle is recommended when using the guidewire.
6. Hold the guidewire in position during removal of the directional sheath and also during progressive insertion of the dilators necessary for the passage of the catheter. Placement of a semirigid sheath of appropriate size into the access tract is a helpful adjunct at this point.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 AS), DFU, MB, Percutaneous Access Set, Global, 50549799-01B

Black (K) ΔE ≤5.0

7. Pass the desired catheter over the guidewire according to the procedure appropriate for that catheter.
8. Secure the catheter to the skin with the retention device. The self adhesive back or sutures may be used to firmly fix this device to the skin. The catheter, in turn, is secured to the retention device using pull ties. (See illustration.) Dressings may also be applied. The connector tube is attached to the catheter and to the drainage bag.

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Khasidy, L.R., Smith A.D.: The Re-Entry Nephrostomy Catheter for Endourological Application. Journal of Urology Vol. 133, (Feb.) 1985.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 AS), DFU, MB, Percutaneous Access Set, Global, 50549799-01B

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Lot
Lote
Lot
Charge
Lotto
Partij
Lote

Recyclable Package
Envase reciclable
Emballage recyclable
Wiederverwertbare Verpackung
Confezione riciclabile
Recyclebare verpakking
Embalagem Reciclável

Use By
Fecha de caducidad
Date limite d'utilisation
Verwendbar bis
Usare entro
Uiterste gebruiksdatum
Validade

Australian Sponsor Address
Dirección del patrocinador australiano
Adresse du promoteur australien
Adresse des australischen Sponsors
Indirizzo sponsor australiano
Adres Australische sponsor
Endereço do Patrocinador Australiano

Argentina Local Contact
Contacto local en Argentina
Contact local en Argentine
Lokaler Kontakt Argentinien
Contatto locale per l'Argentina
Contactpersoon Argentinië
Contacto local na Argentina

For single use only. Do not reuse.
Para un solo uso. No reutilizar.
À usage unique. Ne pas réutiliser.
Für den einmaligen Gebrauch. Nicht wieder verwenden.
Esclusivamente monouso. Non riutilizzare.
Uitsluitend bestemd voor eenmalig gebruik. Niet opnieuw gebruiken.
Apenas para uma única utilização. Não reutilize.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 AS), DFU, MB, Percutaneous Access Set, Global, 50549799-01B

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Recyclable Package

CE 2797

© 2020 Boston Scientific Corporation or its affiliates.
All rights reserved.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 AS), DFU, MB, Percutaneous Access Set, Global, 50549799-01B

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- The Stone Retrieval Basket has no other user serviceable parts.
- Do not attempt to repair or to alter the device.
- Do not resterilize the device.
- Any use of this device, other than those indicated in these instructions is not recommended.

DEVICE DESCRIPTION

The Stone Retrieval Basket consists of a wire shaft (or cannula) housed in a flexible sheath with a multi-wire 360°, self-expanding basket at the distal end. A handle is attached proximally via a luer lock connector to open and close the basket. The Stone Retrieval Basket is designed to be used through the small working channel of an endoscope, such as a ureteroscope or nephroscope during the procedure.

This product contains no DEHP.

INTENDED USE/INDICATIONS FOR USE

The Stone Retrieval Basket is intended to be used during urological procedures to endoscopically grasp, manipulate and remove calculi and other foreign objects.

CONTRAINDICATIONS

None known.

WARNINGS

- Some objects may be too large to be removed endoscopically with a retrieval basket. To avoid object impaction, fluoroscopy or x-ray should be used to determine the size of the object; do not use the Stone Retrieval Basket if the object is too large to be removed endoscopically or held in basket.
- This device must not come in contact with any electrified instrument.
- This device should not be directly fired upon by any lithotrite. To do so may damage the device and could result in patient injury.
- Use of device without introducer may cause damage to the basket.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 92238512 Rev/Ver. A), DFU, MB, Stone Retrieval Basket, Global, 50940729-01A

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- Perforation
- Evulsion
- Edema
- Basket inversion
- Entrapment
- Hemorrhage
- Inability to disengage from irretrievable object
- Pain/Discomfort
- Infection
- Abrasion
- Stone Migration

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OPERATIONAL INSTRUCTIONS

Preparation

The Stone Retrieval Basket is packaged sterile in a sealed blister pouch. To protect the basket assembly during shipment and storage, the open basket is secured in a protective casing.

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1. Pull apart the top two edges of the blister pouch located above the product label to break the seal and present the product to the sterile field.
2. Remove and discard the protective casing.
3. Inspect the basket assembly for any visible damage.
4. Prior to use, ensure that the handle is securely attached to the wire cable by testing the device as outlined below:

Thumb Slide Style Handle

Place handle in hand with thumb tab up. Push thumb tab forward to close basket and pull back to open the basket by advancing it out of the sheath. Test the basket by opening and closing it several times.

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Precaution: It is recommended that the ureter be predilated before endoscopy and object retrieval to prevent stone impaction when extracting the object.

Insertion

(Please refer to contents section of package label)

The Stone Retrieval Basket devices are introduced into an endoscope in either of two ways; i.e., via an introducer to facilitate introduction into the scope or directly into the working channel port of the scope.

(For retrieval basket devices supplied with an Introducer)

1. First, insert Introducer into the nipple of the working channel port of the endoscope.

 (For all retrieval basket devices)

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 92238512 Rev/Ver. A), DFU, MB, Stone Retrieval Basket, Global, 50940729-01A

Black (K) ΔE ≤5.0 V2

2. Maintain the basket cage in the closed position, while introducing the (distal) working end of the device into the working channel of the endoscope.
3. Carefully advance basket device through the working channel. Using deliberate short strokes (2 to 3 cm [20 to 30 mm]), advance the sheath through the scope channel while watching for the sheath as it emerges from the distal scope channel.

Capture

1. Advance the closed basket alongside or beyond the object and open the basket by sliding the thumb tab.
2. Opening basket facilitates stone entrapment.
3. Slowly sweep the open basket across the face of the object to engage it into the wire cage.
4. Sliding thumb tab back and forth to open and close the basket cage facilitate engagement of object.

Removal

Once the object is captured within the basket, carefully maintain pressure on thumb tab in closed position to hold the object.

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Caution: Objects may be too large once captured in the basket to withdraw the basket fully into the sheath or the scope. **Do not exert excessive force to withdraw basketed object into sheath or into scope.** Carefully maintain basket position to hold the object in place.

Objects that are too large to recover through the sheath or through the scope channel will require the scope and basket to be removed simultaneously from the urinary tract.

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Nitinol Stone Retrieval Device

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

- The Stone Retrieval Device has no user serviceable parts.
- Do not attempt to repair or to alter the device.
- Do not resterilize the device.
- Any use of this device, other than those indicated in these instructions is not recommended.

DEVICE DESCRIPTION

The DAKOTA Nitinol Stone Retrieval Device with OpenSure handle consists of a wire shaft (or cannula) housed in a flexible sheath with a multi-wire, open-ended, self-expanding basket at the distal end. A handle is attached proximally to open and close the basket. Incorporated into the handle of the device is a feature that enables the distal end of the device to open to a larger diameter beyond the baseline basket size. It is controlled via the trigger of the handle. This feature is activated by depressing the trigger on the distal end of the handle and is intended to aid in the release of calculi or foreign objects from the basket. The Stone Retrieval Device is designed to be used through the small working channel of an endoscope, such as a ureteroscope or nephroscope during the endoscopic procedure.

This product contains no DEHP (Diethylhexyl-phthalate).

INTENDED USE/INDICATIONS FOR USE

The Stone Retrieval Device is intended to be used during urological procedures to endoscopically grasp, manipulate, release and remove calculi and other foreign objects from the upper urinary tract (ureter, renal pelvis and calyces).

CONTRAINDICATIONS

Use in patients with an untreated urinary tract infection.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 Rev/Ver. AP), DFU, MB, DAKOTA, Global, 90894127-01C

Black (K) ΔE ≤5.0

WARNINGS

- The use of this device must be restricted to use by or under the supervision of physicians trained in urological endoscopic procedures.
- Exercise care to prevent perforation or trauma of the linings of the urinary tract and associated tissues.
- Some objects may be too large to be removed by the retrieval basket without being fragmented into smaller pieces. Appropriate imaging should be done in advance of the procedure to anticipate whether fragmentation will be necessary.
- This device must not come in contact with any electrified instrument.
- This device should not be directly fired upon by any lithotrite (laser, EHL, etc.). Doing so may damage the device and could result in a patient injury.
- If a basket wire is broken during laser lithotripsy, discontinue laser use immediately. Good surgical procedures should be followed for basket removal.
- Use of device without introducer may cause damage to the basket.
- This device contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

ADVERSE EVENTS

Possible complications generally recognized with stone retrieval procedures and/or the use of stone retrieval basket devices include, but may not be limited to:

- Ureteral perforation
- Ureteral avulsion
- Ureteral edema
- Ureteral laceration
- Ureteral stricture
- Urothelial damage
- Pain/Discomfort
- Infection/Sepsis
- Stone impaction
- Stone migration
- Entrapment of stone and basket
- Hemorrhage
- Inability to release a grasped object

HOW SUPPLIED

Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Preparation

The Stone Retrieval Device is packaged sterile in a sealed blister pouch. To protect the basket assembly during shipment and storage, the basket is secured in a protective casing.

Prior to Use

Precaution: Before using, inspect the blister pouch for any breach of the package to ensure a sterile product and inspect product for any damage. If seal has been broken or product is damaged DO NOT USE. Immediately return package and product to your Boston Scientific representative.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 Rev/Ver. AP), DFU, MB, DAKOTA, Global, 90894127-01C

Black (K) ΔE ≤5.0

1. Pull apart the top two edges of the blister pouch located above the product label to break the seal and present the product to the sterile field .
2. Open the protective casing, remove and discard the packing material.
3. Remove the basket assembly and discard the protective casing.
4. Inspect the basket assembly for any visible damage.
5. Prior to use, ensure that the handle is securely attached to the wire cable by testing the device as outlined below.


Figure 1.

Operation of basket

Place handle in hand with thumb tab up. Push thumb tab forward to close basket and pull back to open the basket. Test the basket by opening and closing it several times.

Operation of OpenSure Handle feature:

Activate the OpenSure Handle by depressing the trigger in the direction of the arrows on the sides of the trigger while the slide control is in the open position.

Precaution: Kinks in the sheath will hinder the mechanical operation of the basket.

Precaution: If the ureter is too small to accommodate a ureteroscope and to safely remove the stone or stone fragments then prestenting should be considered.

Insertion

An introducer is supplied with the DAKOTA™ Stone Retrieval Device. Use of the introducer both facilitates insertion of the device and reduces the risk of basket damage.

1. Ensure that the working channel port is open and insert the Introducer into the working channel through the endoscopic cap or adaptor.
2. Maintain the basket cage in the closed position, while introducing the (distal) working end of the device into the working channel of the endoscope.

Caution: Ensure the basket cage is in the closed position during insertion and advancement through the working channel of the endoscope. Insertion or advancement of an open basket may result in basket inversion or device damage.

3. Withdraw the introducer after basket insertion.
4. Carefully advance the basket device through the working channel. Using deliberate short strokes (2 to 3 cm [20 to 30 mm]), advance the basket device through the scope channel until the tip of the device is seen to just emerge from the scope.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 Rev/Ver. AP), DFU, MB, DAKOTA, Global, 90894127-01C

Black (K) ΔE ≤5.0

Capture

1. Under direct visualization, advance the closed basket near the object. Then open the basket by sliding the thumb tab back.
2. Slowly advance the opened basket up to the object until the basket wires come in contact with it.
3. Close the basket around the object by sliding the thumb tab forward.

Note: The OpenSure feature is to better facilitate the release of captured objects.

Removal

Once the object is captured within the basket, carefully maintain pressure on thumb tab in closed position to hold the object in place.

If the object is too large to be removed through the sheath or scope channel then either reorient the object or release the object and break it into fragments that are small enough to be retrieved. In some circumstances, the scope, basket and object can be removed simultaneously ONLY if the ureter is large enough to safely accommodate the object.

Caution: If resistance is encountered during advancement or withdrawal of the device, DO NOT EXERT EXCESSIVE FORCE.

Caution: Objects may be too large once captured in the basket to withdraw the basket fully into the ureter, sheath or the scope. Do not exert excessive force to withdraw basketed object into the ureter, sheath or into scope. Excessive force can either damage the basket or result in the object and basket becoming impacted into the ureter, the sheath or the scope. Excessive force can also result in ureteral laceration or avulsion. Carefully maintain the thumb tab position to hold the object in place

Repositioning and Release

Objects that are too large to recover may be repositioned and released to facilitate lithotripsy.

1. Under direct visualization, reposition the basket and object to the desired location for release.
2. Slide the thumb tab back to open the basket and release the object.
3. If the object cannot be released from the open basket, depress the trigger of the handle control in the direction of the arrow to activate the OpenSure™ Handle feature.

Handle Disassembly

The handle is permanently attached to the sheath. If handle removal is desired, the sheath and drivewire must be cut: Hold the handle and proximal portion of the sheath, keeping the sheath taught.

1. To facilitate cutting of the sheath components, keep the sheath and drivewire taught
2. Use surgical scissors to cut through the sheath and drivewire
3. Dispose of the device in accordance with local laws and policy.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 Rev/Ver. AP), DFU, MB, DAKOTA, Global, 90894127-01C

Black (K) ΔE ≤5.0


Figure 2.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040 Rev/Ver. AP), DFU, MB, DAKOTA, Global, 90894127-01C

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DEVICE DESCRIPTION

Contents

The Amplatz Type Renal Sheath Set includes the following:

CONTENTS:

4 - Sheaths: 24F (8.1mm) – 30F (10.1mm) I.D., 17 cm Length, Radiopaque PTFE

The Amplatz Type Renal Sheath is made of radiopaque PTFE. Renal sheaths are used to establish percutaneous access for surgical instruments.

Carefully read all the instructions prior to using this device.

INTENDED USE/INDICATIONS FOR USE

The Amplatz Type Renal Sheath is recommended for the establishment of percutaneous access.

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ADVERSE EVENTS

The adverse events which may result from an access procedure include but are not limited to:

- Tissue trama
- Tissue perforation
- Acute bleeding
- Injury to the kidney

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Handling and Storage

Store product at a controlled room temperature. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that products are used prior to the expiration date on package label.

OPERATIONAL INSTRUCTIONS

Recommended Procedures

1. Establishment of a nephrostomy tract and renal dilatation should always be performed under fluoroscopic control.
2. Percutaneous access is established and a 0.038 in (0.97 mm) diameter guidewire passed into a secure location within the intra-renal system. If feasible, guidewire should be passed down ureter and into the bladder.
3. Dilate the tract to the desired size either with a nephrostomy dilatation balloon catheter or with sequential rigid dilatation catheters.
4. Once the tract is established, the appropriate sheath is passed over the dilatation device and positioned as desired.
5. Remove the dilator from the sheath, maintaining the sheath and guidewire in position.
6. Stone manipulation may now be performed through the sheath.

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DEVICE DESCRIPTION

The Clear Renal Sheath is a 30F (10.0 mm) I.D., 34F (11.5 mm) O.D., 17 cm length sheath with a guidewire notch at the proximal end and a radiopaque marker at the distal tip. Renal sheaths are used to establish percutaneous access for instrumentation.

INTENDED USE/INDICATIONS FOR USE

The Clear Renal Sheath is recommended for the establishment of percutaneous access.

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None known.

PRECAUTIONS

A thorough understanding of the technical principles, clinical applications, and risks associated with percutaneous nephrostomy is necessary before using this product.

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ADVERSE EVENTS

The adverse events which may result from an access procedure include:

- Tissue trauma
- Tissue perforation
- Acute bleeding
- Injury to the kidney

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OPERATIONAL INSTRUCTIONS

Recommended Procedures

1. Establishment of a nephrostomy tract and renal dilation should always be performed under fluoroscopic control.
2. Percutaneous access is established and a 0.038 in (0.97 mm) diameter guidewire passed into a secure location within the intra-renal system. If feasible, guidewire should be passed down ureter and into bladder.
3. Dilate the tract to the desired size either with a nephrostomy dilatation balloon catheter or with sequential rigid dilatation catheters.
4. Once the tract is established, the appropriate sheath is passed over the dilatation device and positioned as desired.
5. Remove the dilator from the sheath, maintaining the sheath and guidewire in position.
6. Guidewire may be secured in guidewire notch.
7. Stone manipulation may now be performed through the sheath.

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240AE), DFU, MB, Clear Renal Sheath, Global, 90959786-01B

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Boston Scientific

Navigator™

Ureteral Access Sheath Set with Hydrophilic Coating

Directions for Use	**3**
Instrucciones de uso	**6**
Mode d'emploi	**9**
Gebrauchsanweisung	**12**
Istruzioni per l'uso	**15**
Gebruiksaanwijzing	**18**
Instruções de Utilização	**21**

90964948-01

2019-12

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, Navigator, Global, 90964948-01B

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, Navigator, Global, 90964948-01B

Navigator™

Ureteral Access Sheath Set with Hydrophilic Coating

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DEVICE DESCRIPTION

The Navigator Ureteral Access Sheath Set with Hydrophilic Coating consists of an inner tapered dilator and an outer sheath. The Navigator Ureteral Access Sheath Set may be advanced over a 0.035 in. (0.89 mm) or 0.038 in. (0.97 mm) guidewire. The sheath is designed to create a conduit for urological procedure instruments. Refer to package label for sheath dimensions.

INTENDED USE/INDICATION FOR USE

The Navigator Ureteral Access Sheath is indicated for use in endoscopic procedures to facilitate the passage of endoscopes, urological instruments and for the injection of fluids into the urinary tract.

CONTRAINDICATIONS

- Patients who are contraindicated for retrograde urological procedures.
- Patients who are contraindicated for antegrade urologic procedures, including, but not limited to patients with blood clotting anomalies due to coagulopathies or pharmacological anticoagulations.
- Patients who have the presence of tight strictures which would limit use of the device.
- Patients who have the presence of large obstructing distal ureteral calculi.

PRECAUTIONS

The recommendations given are meant to serve only as a basic guide to the utilization of this access sheath set. The ureteral sheath set should not be used without comprehensive knowledge of the indications, techniques and risks of the procedure.

To minimize resistance during advancement, ensure the hydrophilic coating on the dilator and sheath is activated with saline or sterile water prior to placement.

DO NOT bend or kink the dilator or sheath prior to placement; to do so could damage the integrity of the device and result in patient injury.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, Navigator, Global, 90964948-01B

If resistance is encountered during advancement or withdrawal of the device, **STOP**. **DO NOT** continue without first determining the cause of the resistance and taking remedial action.

DO NOT advance the sheath without the dilator in place; to do so could cause injury to the patient.

ADVERSE EVENTS

Potential adverse events associated with the use of the transurethral access device include, but are not limited to:

- Mucosal irritation, inflammation and edema
- Uretheral strictures
- Acute bleeding or hemorrhage
- Urethral, bladder, or ureteral perforation
- Other injury to the urinary tract

HOW SUPPLIED

Store in a cool, dry, dark place. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that products are used prior to the expiration date on package label. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

DIRECTIONS FOR USE

1. Place a 0.035 in. (0.89 mm) or 0.038 in. (0.97 mm) guidewire into the desired location within the urinary tract.
2. Prior to use, place the dilator and sheath into a container of saline or sterile water to activate the hydrophilic coating.
3. Insert the dilator fully inside the sheath and secure by pushing the dilator hub until it snaps into the sheath hub.
4. Advance the dilator/sheath assembly over the guidewire to the desired location.
5. Once positioned, gently withdraw the dilator by grasping the luer and the tab, while maintaining sheath position. DO NOT advance the sheath without the dilator in place.
6. If desired, the sheath can be secured to the surgical drapes by attaching sutures through the holes located on the sheath hub body.
7. Introduce desired instrumentation through the sheath as needed.
8. To perform a retrograde pyelogram, reinsert the dilator into the sheath (see Step 3) and inject contrast through the luer fitting of the dilator. Remove dilator as desired per Step 5.
9. Upon completion of the access procedure, gently withdraw the device.
10. Discard device upon completion of procedure in accordance with hospital procedures.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, Navigator, Global, 90964948-01B

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DEVICE DESCRIPTION

The Helical Stone Retrieval Basket consists of a wire cable housed in a flexible sheath with a round-wire, helix-shaped basket attached to the distal end. A handle is attached proximally via a luer lock connector to open and close the basket. The Helical Stone Retrieval Basket is designed to be used through the small working channel of an endoscope, such as a ureteroscope or nephroscope during the procedure.

INDICATIONS FOR USE/INTENDED USE

The Helical Stone Retrieval Basket is used endoscopically to entrap and remove calculi and other foreign objects from the ureter.

CONTRAINDICATIONS

None known.

WARNINGS

Use of this device should be restricted to use by or under the supervision of physicians trained in urologic endoscopic procedures. Care should be exercised to prevent perforation or trauma of the linings and associated tissue of vessels or ducts.

The Helical Stone Retrieval Basket is for single patient use and has no other user serviceable parts.

- Some objects may be too large to be removed endoscopically with a retrieval basket. To avoid object impaction, fluoroscopy or x-ray should be used to determine the size of the object; do not use the Helical Stone Retrieval Basket if the object is too large to be removed endoscopically or held in the basket.
- Do not resterilize the device.
- Any use of this device, other than those indicated in these instructions is not recommended.

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- Perforation
- Entrapment of calculi too large to be removed
- Evulsion

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, Helical Stone Retrieval Basket, Global, 90966738-01B

- Hemorrhage
- Inability to disengage from irretrievable stone
- Edema

HOW SUPPLIED

The Helical Stone Retrieval Basket is packaged sterile in a sealed blister pouch. To protect the basket assembly during shipment and storage, the open basket is secured in a protective casing.

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OPERATIONAL INSTRUCTIONS

Prior To Use

Precaution: Before using, inspect the blister pouch for any breach of the package to ensure a sterile product and inspect product for any damage. If seal has been broken or product is damaged **DO NOT USE.** Immediately return package and product to Boston Scientific.

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 A. Thumb Slide Style Handle

 Place handle in hand with thumb slide up. Push thumb tab forward to close basket and pull back to open the basket by advancing it out of the sheath. Test the basket by opening and closing it several times.

 Precaution: Kinks in the sheath will hinder the mechanical operation of the basket.

 Precaution: It is recommended that the ureter be predilated before endoscopy and object retrieval to prevent object impaction when retracting the object.

WARNING

- To avoid object impaction, use fluoroscopy or x-ray to determine the size of the object. Do not use basket if object is too large to be removed endoscopically or held in the basket.
- This device must not come in contact with any electrified instrument.
- This device should not be directly fired upon by a pulsed dye laser. To do so may damage the device and could result in patient injury.

Insertion

1. Close the basket and maintain light pressure.
2. Insert the distal end of the sheath into the channel port of the endoscope. Using deliberate short strokes (2-3 cm), advance the sheath through the scope channel while watching for the sheath as it emerges from the distal scope channel.

Capture

1. Advance the closed basket beyond the object and open the basket.
2. Slowly sweep the open basket across the face of the object to capture the object in the wire cage. The Helical Stone Retrieval Basket can be rotated to aid basket positioning and to enhance the sweeping action.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, Helical Stone Retrieval Basket, Global, 90966738-01B

Removal

Once the object is in the confines of the basket, carefully maintain pressure in closed position to hold the object.

Caution: If resistance is encountered during advancement or withdrawal of the device, **DO NOT EXERT EXCESSIVE FORCE**.

Precaution: Objects may be too large once captured in the basket to withdraw the basket fully into the sheath. **Do not exert excessive force to withdraw basket into sheath or scope.** Carefully maintain basket position to hold the object in place.

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Handle Disassembly

A. Thumb Slide Style Handle

1. Unscrew the luer lock connector at the sheath/ handle junction.
2. Unscrew the flange cap at the proximal end of the handle.
3. Slide the handle off the mandrel.

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DEVICE DESCRIPTION

Renal dilators and sheaths are used to dilate nephrosotomy tracts to allow instrumentation.

Contents

12 Dilators, 8F through 30F (2.7 mm through 10.1 mm) (even sizes) of 35 cm length

4 Sheaths, 24F through 30F (8.1 mm through 10.1 mm) (even sizes)

1 Stylet, 8F (2.7 mm) O.D., 70 cm length

1 Sheath, 8F (2.7 mm) I.D., 10F (3.4 mm) O.D. radiopaque PTFE, of 30 cm length

INTENDED USE/INDICATIONS FOR USE

The Amplatz Type Renal Dilator/Sheath Set is recommended for dilatation of the nephrostomy tract and sheath placement.

CONTRAINDICATIONS

This product is contraindicated in the presence of conditions which create unacceptable risk during catheterization.

WARNINGS

None known.

PRECAUTIONS

A thorough understanding of the technical principles, clinical applications, and risks associated with nephrostomy dilatation is necessary before using this product.

The recommendations given are meant to serve only as a basic guide to the utilization of this set. The performance of percutaneous nephrostomy should not be undertaken without comprehensive knowledge of the indications, techniques and risks of the procedure.

ADVERSE EVENTS

The adverse events which may result from a dilatation procedure include:

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Recommended Procedures

Prior to Use

1. Establishment of a nephrostomy tract and renal dilatation should always be performed under fluoroscopic control.
2. Percutaneous access is established and a 0.038 in (0.97 mm) diameter guidewire passed into a secure location within the intra-renal system. If feasible, guidewire should be passed down ureter and into bladder.
3. Progressively dilate the tract over the 0.038 in (0.97 mm) diameter guidewire with the 8 through 12F (2.7 mm - 4.1 mm) fascial dilators while holding the guidewire in place.
4. Pass the 8F (2.7 mm) catheter/stylet and sheath over guidewire into the tract.
5. Remove the 8F (2.7 mm) catheter/stylet while maintaining position of the guidewire and sheath.
6. Place additional 0.038 in (0.97 mm) diameter guidewire(s) through sheath into a secure position location in the intra renal system.
7. Remove the sheath maintaining position of guidewires. The safety guidewire(s) is then laid aside and affixed to the drape.
8. Pass the 8F (2.7 mm) catheter/stylet over the working guidewire.
9. The tract may now be dilated progressively from 14 through 30F (4.7 mm - 10.1 mm) by passing the dilators over the 8F (2.7 mm) catheter/stylet.
10. When dilating the tract, the 8F (2.7 mm) catheter/stylet must be held in position and the dilators passed into the tract using a twisting motion.
11. The working sheaths on the 24 through 30F (8.1 mm - 10.1 mm) dilators may be passed with the dilators or over the dilator after placement of the dilator in the tract.
12. Once the tract has been dilated to the desired size, the 8F (2.7 mm) catheter/stylet is then removed, maintaining the guidewire, dilator and sheath in position.
13. Remove the dilator from the sheath, maintaining the sheath and guidewire in position.
14. Stone manipulation may now be performed through the sheath.

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Boston Scientific

Amplatz Type Graduated Renal Dilator Set

Directions for Use	2
Instrucciones de uso	4
Mode d'emploi	7
Gebrauchsanweisung	10
Istruzioni per l'uso	13
Gebruiksaanwijzing	16
Instruções de Utilização	19

90967963-01

2019-10

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, Amplatz Grad Renal Dilation Set, Global, 90967963-01B

Amplatz Type Graduated Renal Dilator Set

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DEVICE DESCRIPTION

Renal dilators are used to dilate nephrostomy tracts.

Contents

12 Dilators, 8F through 30F (2.7 mm through 10.1 mm) (Even Sizes) of 35 cm length

1 Stylet, 8F (2.7 mm) O.D., 70 cm length

1 Sheath, 8F (2.7 mm) I.D., 10F (3.4 mm) O.D. Radiopaque PTFE, of 30 cm length

INTENDED USE/INDICATIONS FOR USE

The Amplatz Type Graduated Renal Dilator Set is recommended for dilatation of the nephrostomy tract.

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OPERATIONAL INSTRUCTIONS

Recommended Procedures

1. Establishment of a nephrostomy tract and renal dilatation should always be performed under fluoroscopic control.
2. Percutaneous access is established and a 0.038 in (0.97 mm) diameter guidewire passed into a secure location within the intra-renal system. If feasible, guidewire should be passed down ureter and into bladder.
3. If desired, pass the 8F (2.7 mm) catheter/stylet and sheath over guidewire into the tract.
4. Remove the 8F (2.7 mm) catheter/stylet while maintaining position of the guidewire and sheath.
5. Place additional 0.038 in (0.97 mm) diameter guidewire(s) through sheath into a secure position location in the intra renal system.
6. Remove the sheath maintaining position of guidewires. The safety guidewire(s) is then laid aside and affixed to the drape.
7. Pass the 8F (2.7 mm) catheter/stylet over the working guidewire.
8. The tract may now be dilated progressively from 14 through 30F (4.7 through 10.1 mm) by passing the dilators over the 8F (2.7 mm) catheter/stylet.
9. When dilating the tract, the 8F (2.7 mm) catheter/stylet must be held in position and the dilators passed into the tract using a twisting motion.
10. The working sheaths on the 24 through 30F (8.0 through 10.1 mm) dilators may be passed with the dilators or over the dilator after placement of the dilator in the tract.
11. Once the tract has been dilated to the desired size, the 8F (2.7 mm) catheter/stylet is then removed, maintaining the guidewire, dilator and sheath in position.
12. Remove the dilator from the sheath, maintaining the sheath and guidewire in position.
13. Stone manipulation may now be performed through the sheath.

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Amplatz Type Graduated Renal Dilator Set

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
DEVICE DESCRIPTION

The Dual Lumen Ureteral Catheter is a catheter that will provide for placement of a safety guidewire and second procedural guidewire, if necessary, drainage of urine and contrast injection.

INTENDED USE/INDICATIONS FOR USE

Ureteral Catheters are indicated for use for catheterization of the urinary tract during interoperative procedures, including but not limited to the following applications:

- Ureteral occlusion for injection of contrast media
- Delivery of irrigation fluids to the ureter
- Drainage of fluids from the urinary tract
- Access, advancement or exchange of ancillary devices including guidewires.

CONTRAINDICATIONS

Ureteral catheters are in general contraindicated for use in patients unsuitable for catheterization.

PRECAUTIONS

1. Do not use this product without reading and understanding the complete instructions herein.
2. Avoid bending or kinking the catheter prior to placement. To do so could damage the integrity of the catheter and result in patient injury.
3. If resistance is encountered during advancement or withdrawal of the catheter, **STOP. DO NOT CONTINUE** without first determining the cause of the resistance and taking remedial action.

WARNINGS

None known.

ADVERSE EVENTS

The adverse events which may result from ureteral catheter placement include, but are not limited to:

- Extravasation
- Hemorrhage
- Sepsis
- Edema
- Loss of renal function
- Perforation of kidney, renal pelvis, ureter and/or bladder
- Peritonitis
- Urinary tract infection

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Handling and Storage

Store at controlled room temperature. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that products are used prior to the sterilization expiration date on package label.

OPERATIONAL INSTRUCTIONS

Prior to Use

1. Cystoscopically place an 0.038 in (0.97 mm) dia. guidewire an optimal distance up the ureter.
2. Fill the yellow injection lumen with the desired solution prior to introduction to help eliminate air bubbles.

RECOMMENDED TECHNIQUES

Contrast Media Injection

1. Attach the luer injection hub to the end of the catheter.
2. Advance the tip end of the catheter through the cystoscope and into the ureteral orifice.
3. Properly place catheter tip to occlude the ureter. Inject contrast media through luer injection hub.
4. Visualize the urinary tract radiographically.
5. Withdrawal catheter through cystoscope.

Guidewire Placement/Introduction and Placement of Second Guidewire

1. To access the site for guidewire placement, pass the tip end of the catheter through the cystoscope and into the ureteral orifice.
2. Confirm catheter position radiographically and then advance the appropriate guidewire through the catheter.
3. To place a second guidewire, advance the second guidewire through the second or open port, and advance through the catheter.
4. Confirm guidewire(s) placement, then slowly remove the catheter while maintaining the guidewire position.

Guidewire Exchange

1. Back feed the existing guidewire into the cystoscope and out of the working port. Take care not to lose the guidewire access.
2. Slowly pass the cystoscope over the guidewire into the bladder.
3. Introduce the tip end of the catheter over the guidewire, through the scope, and advance to desired location.
4. Confirm the catheter position radiographically.
5. Maintain the catheter position with the cystoscope and remove the pre-placed guidewire.
6. Insert the exchange into the catheter until in proper position.
7. Withdraw the catheter through the cystoscope.

DRAINAGE

Caution: Indwelling time not to exceed 24 hours.

1. Introduce the tip end of the catheter over the appropriate guidewire, through the scope, and advance to desired location.
2. Carefully remove the guidewire and cystoscope.
3. Confirm catheter position radiographically and secure externally to prevent dislodgment.
4. Attach appropriate drainage system to the other end.
5. To remove catheter, grasp firmly and gently pull.

REFERENCES

1. Lyon, W.S., University of Chicago, Illinois.
2. Lupu, A.N., Fuchs, G.J., Chaussy, C. G., Depts. of Surgery and Urology, University of California, Los Angeles: A New Approach to Ureteral Stone Manipulation for ESWL, Endourology Publications, Inc. (Volume 1, Number 1, March 1986).

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose**. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments**.

Boston Scientific (Master Brand DFU Template 5in x 5.25in Global, 90348240AE), DFU, MB, Dual Lumen Ureteral Catheter, Global, 90968455-01B

Boston Scientific (Master Brand DFU Template 5in x 5.25in Global, 90348240AE), DFU, MB, Dual Lumen Ureteral Catheter, Global, 90968455-01B

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DEVICE DESCRIPTION

Contents

The Angled Tip Ureteral Catheter includes the following components:

Includes 2 - Catheter Connectors

Includes 1 - 0.035 in (0.89 mm) Bentson Tip Guidewire

Ureteral Catheters are made from biocompatible materials. A variety of tip configurations and sizes are available to facilitate ureteral access. A separate luer hub is provided for access and collection. Refer to recommended guidewire size on package label if a guidewire is required.

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- Ureteral occlusion for injection of contrast medium
- Delivery of irrigation fluids to the ureter
- Drainage of fluids from the urinary tract
- Access, advancement or exchange of ancillary devices including guidewires

CONTRAINDICATIONS

Ureteral Catheters are in general contraindicated for use in patients unsuitable for catheterization.

WARNINGS

None known.

PRECAUTIONS

- Avoid bending or kinking the catheter prior to placement. To do so could damage the integrity of the catheter and result in patient injury.
- If resistance is encountered during advancement or withdrawal of the catheter, **STOP. DO NOT CONTINUE** without first determining the cause of the resistance and taking remedial action.

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- Extravasation
- Hemorrhage
- Sepsis
- Edema
- Loss of renal function
- Perforation of kidney, renal pelvis, ureter and bladder
- Peritonitis
- Urinary tract infection

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Handling and Storage

Store catheters at controlled room temperature. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that catheters are used prior to the sterilization expiration date on package label.

OPERATIONAL INSTRUCTIONS

Prior to Use

Carefully read all instructions **Prior** to using this device.

RECOMMENDED PROCEDURES

Catheter Selection

Selection of appropriate catheter size and tip configuration is dependent on patient anatomy, application and user preference. Refer to package label for sizing information.

Boston Scientific (Master Brand DFU Template 5in x 5.25in Global, 90348240AG), DFU, MB, Ureteral Catheters, Global, 90971494-01C

Black (K) ΔE ≤5.0

Note: Remove product inner diameter protector prior to use.

Contrast Media Injection

1. Attach the luer injection hub to the end of the catheter.
2. Advance the tip end of the catheter through cystoscope and into the ureteral orifice.
3. Properly place the catheter tip to occlude the ureter. Inject contrast media through luer injection hub.
4. Visualize the urinary tract radiographically.
5. Withdraw catheter through cystoscope.

Guidewire Placement

1. To access the site for guidewire placement, pass the tip end of the catheter through the cystoscope and into ureteral orifice.
2. Confirm catheter position radiographically and then advance the appropriate guidewire through the catheter.
3. Confirm guidewire placement, then slowly remove the catheter while maintaining guidewire position.

Guidewire Exchange

1. Back feed the existing guidewire into the cystoscope and out of the working port. Take care not to lose the guidewire access.
2. Slowly pass the cystoscope over the guidewire into the bladder.
3. Introduce the tip end of the catheter over the guidewire, through the scope, and advance to desired location.
4. Confirm the catheter position radiographically.
5. Maintain the catheter position within the cystoscope and remove the preplaced guidewire.
6. Insert the exchange guidewire into the catheter until in proper position.
7. Withdraw the catheter through the cystoscope.

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1. Introduce the tip end of the catheter over the appropriate guidewire, through the scope, and advance to desired location.
2. Carefully remove the guidewire and the cystoscope.
3. Confirm catheter position radiographically and secure externally to prevent dislodgment.
4. Attach one end of luer hub to proximal end of catheter and attach appropriate drainage system to the other end.
5. To remove catheter, grasp firmly and gently pull.

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Recyclable Package

CE 2797

© 2019 Boston Scientific Corporation or its affiliates.
All rights reserved.

Boston Scientific (Master Brand DFU Template 8.2677in x 11.6929in A4, 90105918AN), DFU, MB, Ureteral Catheters, ru, 90971494-15C

Black (K) ΔE ≤5.0

Boston Scientific

Single Action Pumping System

For Standard & Optional*
Continuous Flow Models

*Available for Urology Applications

Directions for Use	2
Instrucciones de uso	4
Mode d'emploi	6
Gebrauchsanweisung	8
Istruzioni per l'uso	10
Gebruiksaanwijzing	12
Instruções de Utilização	14

90986318-01

2019-12

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240AE), DFU, MB, Single Action Pumping System, Global, 90986318-01B

Boston Scientific (Master Brand DFU Template 5in x 5.25in Global, 90348240AE), DFU, MB, Single Action Pumping System, Global, 90986318-01B

Single Action Pumping System

For Standard & Optional* Continuous Flow Models

*Available for Urology Applications

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DEVICE DESCRIPTION

The Single Action Pumping System is designed for controlled irrigation through an Endoscope or accessory during urological and/or endoscopic procedures. One-way valves provide automatic refill of the syringe for access to irrigation fluid.

Contents

A single package contains one single use pumping system.

INTENDED USE/INDICATIONS FOR USE

The Single Action Pumping System is intended for use with urological and/or endoscopic surgical procedures requiring continuous irrigation.

CONTRAINDICATIONS

None known.

PRECAUTIONS

Do not use this product without reading and understanding the complete instructions enclosed herein.

Note: This product contains DEHP (Bis(2-ethylhexyl) phthalate). BSC has conducted a Chemical Safety Assessment of DEHP in this device based on currently known acceptable exposure limits for sensitive populations (children and pregnant/nursing women), and has determined that the potential exposure is below the known acceptable exposure limits for the device when used in accordance with the indications and instructions for use.

ADVERSE EVENTS

Any irrigation product may cause the following adverse tissue reactions.

- Hydro-nephrosis of the Urinary Tract
- Excess Distention of the Kidney

HOW SUPPLIED


The Single Action Pumping System is supplied sterile by ethylene oxide gas. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place. Rotate inventory so that devices are used prior to the expiration date on package label.

DIRECTIONS FOR USE

1. Attach vacuum syringe "A" to dual check valve Luer Lock "B".
2. Attach side inlet tube "C" to fluid supply "D" by removing protective cap before inserting Fluid Spike (Fluid Spike may be detached from side tube inlet if male luer connector is desired).
3. *(If using Optional Continuous Flow Model)* Close Roller Clamp "E".
4. Depress Plunger "F" and release. This will expel air and begin to draw in liquid. Repeat the process to bleed all air out of the system. Point the tip of the vacuum syringe upward to complete this process.
5. Attach forward outlet tube "G" to the device receiving fluid and depress the plunger to expel all air from the system before using.
6. To pump, depress the plunger "F" to deliver the desired amount of fluid and release. The vacuum syringe will automatically refill.
7. *(If using Optional Continuous Flow Model)* To irrigate with gravity fed continuous flow: open Roller Clamp "E" to allow fluid to flow. Before removing device, close Roller Clamp "E".

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Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240AE), DFU, MB, Single Action Pumping System, Global, 90986318-01B

Boston Scientific (Master Brand DFU Template 5in x 5.25in Global, 90348240AE), DFU, MB, Single Action Pumping System, Global, 90986318-01B

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DEVICE DESCRIPTION

The Escape Nitinol Stone Retrieval Basket consists of a nitinol wire housed in a flexible sheath with a nitinol multi-wire, self-expanding basket at the distal end. A handle is attached proximally to open and close the basket. The retrieval basket is designed to be used through a working channel of an endoscope, such as a nephroscope or ureteroscope, during urological procedures.

This package contains one (1) basket, one (1) Y-Adaptor to facilitate the use of the basket with a maximum 200 µm holmium laser fiber, and one (1) t-connector to provide an alternate irrigation port for endoscopes without a separate irrigation port.


Figure 1. (A) Y-Adaptor, (B) Y-Adaptor with t-connector, (C) Retrieval Device

INDICATIONS FOR USE/INTENDED USE

This basket is intended to be used endoscopically to grasp, manipulate and remove calculi from the urinary tract.

CONTRAINDICATIONS

- Pregnancy
- Renal malignancy
- Sepsis
- Anticoagulation or bleeding disorders
- Presence of very large stones that in the past have failed lithotripsy with basket retrieval

WARNINGS

- The use of this device should be restricted to use by or under the supervision of physicians trained in urological endoscopic procedures.
- Exercise care to prevent perforation or trauma of the linings of the urinary tract and associated tissues.
- Some stones may be too large to be removed endoscopically with the stone retrieval basket. Use fluoroscopy or x-ray to determine the size of the stone prior to instrumentation. Do not use the basket if the stone is too large to be removed endoscopically. To do this may result in stone impaction and patient injury.
- If resistance is encountered while attempting to withdraw the basket, DO NOT EXERT EXCESSIVE FORCE.
- Do not fire the laser directly upon the basket.

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040AL), DFU, MB, Escape, Global, 90998359-01B

Black (K) ΔE ≤5.0

- If a basket wire is broken during lasing, discontinue lasing immediately. Good surgical procedures should be followed for basket removal.
- This basket has no serviceable parts. Do not attempt to repair or alter the device. Do not resterilize the device. Any use of the device other than those indicated is not recommended.
- Do not use this basket without the pre-packaged Y-adapter. Doing so may result in inability to advance the laser fiber through the working channel of the scope.
- This device must not come in contact with any electrified instrument.

ADVERSE EVENTS

Potential complications generally recognized with the use of stone retrieval basket devices include, but are not limited to:

- Edema
- Entrapment
- Evulsion
- Fever
- Hemorrhage
- Pain
- Stone impaction
- Stone migration
- Urinary frequency
- Urinary tract infection
- Ureteral trauma (including stripping, abrasion or perforation)

There may be other potential adverse events that are unforeseen at this time.

HOW SUPPLIED

The Escape™ Nitinol Stone Retrieval Basket is protected by a clam shell tray which is placed inside a sealed sterile pouch. The Escape Nitinol Stone Retrieval Basket is sold in single unit configurations.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Prior To Use

Precaution: Before using, inspect the blister pouch for any breach of the package to ensure a sterile product and inspect product for any damage. If seal has been broken or product is damaged. ***DO NOT USE.*** Immediately return package and product to your Boston Scientific representative.

1. Pull apart the top two edges of the blister pouch located above the product label to break the seal and present the product to the sterile field.
2. Remove and discard the protective casing.
3. Inspect the device for any visible damage.
4. Prior to use, ensure the handle is securely attached to the wire cable and the basket is working properly by using the push thumb tab on the handle to open and close the basket several times. Withdraw the closed basket into the sheath prior to inserting the device into the scope.

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Precaution: It is recommended that the ureter be pre-dilated before endoscopy and stone retrieval to prevent impaction when retracting the stone.

Insertion

This basket is introduced into the endoscope via the Y-adapter per Figure 2.

1. The Y-adapter is preloaded onto the sheath. Close the basket and pull the retrieval device out of the adaptor.
2. Attach the Y-adapter into the working channel port of the scope. For the scopes without an irrigation port, detach the blue seal on the Y-adapter and attach the t-connector per Figure 1.
3. Maintain the basket cage in the closed position while introducing the distal working end of the device into the working channel of the endoscope.
4. Carefully advance the basket device through the working channel of the endoscope and look for the sheath as it emerges from the distal end of the scope channel.


Figure 2. Device Setup

Boston Scientific, (Master Brand DFU Template 3in x 9in Global, 90106040AL), DFU, MB, Escape, Global, 90098359-01B

Black (K) ΔE ≤5.0

Capture

1. Advance the closed basket alongside or beyond the stone and open the basket by sliding the thumb tab.
2. Slowly sweep the open basket across the face of the stone to engage the stone into the wire cage.
3. Sliding the thumb tab back and forth to open and close the basket will also facilitate engagement of the stone.
4. Once the stone is captured within the basket, carefully maintain pressure on the thumb tab to hold the stone.

Laser Use

Warning: Refer to the Directions for Use supplied by the manufacturer of the laser and the laser fiber.

Warning: Do not fire the laser if you cannot see the aiming beam. Do not fire the laser directly upon the basket.

1. With the stone captured, insert the 200 µm laser fiber into the open port on the Y-Adaptor (second port) per Figure 2.
2. Carefully advance the fiber through to the end of the working channel of the endoscope.
3. Using the notches on the handle as a reference point, slightly open the basket cage by sliding back the blue thumb tab approximately 1–2 notches to further separate the basket wires prior to lasing, as shown in Figure 3.
4. Position the laser fiber in desired location in order to fragment stone. Ensure that the stone is visible between the wires and that the fiber is not in direct contact with the basket wire.
5. Carefully fire the laser to fragment the stone into pieces that can be safely removed.


Figure 3. Handle Detail

Removal

1. Slowly and carefully, withdraw the basket with the stone fragment(s).
2. If resistance is encountered during advancement or withdrawal of the basket, **DO NOT EXERT EXCESSIVE FORCE.**
3. Some stones may be too large once captured in the basket to withdraw the basket through the scope. In this situation, remove the scope and basket simultaneously from the urinary tract while carefully maintaining basket position to hold the stone in place.

Handle Disassembly

If handle removal is desired:

1. Unscrew the luer lock connector at the sheath/handle junction.
2. Unscrew the flange cap at the proximal end of the handle.
3. Slide the handle off the mandrel.

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DEVICE DESCRIPTION

The OptiFlex Nitinol Stone Retrieval Basket consists of a shaft housed in a flexible sheath with a nitinol multi-wire and a self-expanding basket at the distal end. A handle is attached proximally to open, close and rotate the basket. The retrieval basket is designed to be used through a working channel of an endoscope, such as a nephroscope or ureteroscope, during urological procedures.

Contents

This package contains a Nitinol Stone Retrieval Basket and an Introducer, See Figure 1.


Figure 1. (A) Retrieval Basket; (B) Introducer

INTENDED USE/INDICATION FOR USE

The OptiFlex Nitinol Stone Retrieval Basket is intended to be used endoscopically to grasp, manipulate and remove calculi and other foreign objects from the urinary tract.

CONTRAINDICATIONS

The OptiFlex Nitinol Stone Retrieval Basket is contraindicated for:

- Renal malignancy
- Sepsis
- Anticoagulation or bleeding disorders

WARNINGS

- The use of this device should be restricted to use by or under the supervision of physicians trained in urological endoscopic procedures.
- Exercise care to prevent perforation or trauma of the linings of the urinary tract and associated tissues.
- Some stones may be too large to be removed endoscopically with the nitinol stone retrieval basket. Use

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AL), DFU, MB, OptiFlex, Global, 90998360-01B

Black (K) ΔE ≤5.0

fluoroscopy or x-ray to determine the size of the stone prior to instrumentation. Do not use the basket if the stone is too large for endoscopic removal. To do so may result in stone impaction and patient injury.

- If resistance is encountered while attempting to withdraw the basket, **DO NOT EXERT EXCESSIVE FORCE.**
- This device should not be directly fired upon by a lithotrite. To do so may damage the device and could result in patient injury.
- If a basket wire is broken during use, discontinue use immediately. Good surgical procedures should be followed for basket removal.
- This basket has no serviceable parts. Do not attempt to repair or alter the device. Do not resterilize the device. Any use of the device other than those indicated is not recommended.
- Do not use the device without the Introducer as this may result in damage to the basket.
- This device must not come in contact with any electrified instruments.

ADVERSE EVENTS

Adverse events generally recognized with the use of stone retrieval basket devices include, but are not limited to:

- Edema
- Stone impaction
- Entrapment
- Stone migration
- Ureteral Avulsion
- Urinary frequency
- Fever
- Urinary tract infection
- Hemorrhage
- Pain
- Ureteral Intussusception
- Ureteral Stricture
- Ureteral trauma (including stripping, abrasion or perforation)

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HOW SUPPLIED

The OptiFlex™ Nitinol Stone Retrieval Basket is protected by a clam shell tray which is placed inside a sealed sterile pouch. The OptiFlex Nitinol Stone Retrieval Basket is sold in single unit configurations.

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Handling and Storage

Store in cool, dry, dark place. Product expiration date is displayed on the product label.

OPERATIONAL INSTRUCTIONS

Prior To Use

Precaution: Before using, inspect the product pouch for any breach of the package to ensure a sterile product and inspect product for any damage. If the seal or product is damaged, DO NOT USE. Immediately return package and product to your Boston Scientific representative.

1. Pull apart the top two edges of the blister pouch located above the product label to break the seal and present the product to the sterile field.
2. Remove and discard the clam shell tray.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AL), DFU, MB, OptiFlex, Global, 90998360-01B

Black (K) ΔE ≤5.0

3. Inspect the device for any visible damage.
4. Prior to use, ensure the handle is securely attached to the wire cable and the basket is working properly by using the thumb wheel on the handle to open and close the basket several times. Withdraw the closed basket into the sheath prior to inserting the device into the endoscope.

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Precaution: It is recommended that the ureter be pre-dilated before endoscopy and stone retrieval to prevent impaction when extracting the stone.

Insertion

The OptiFlex™ Nitinol Stone Retrieval Basket device is inserted into an endoscope via an Introducer to facilitate introduction into the endoscope.

1. First, insert the Introducer into the nipple of the working channel port of the endoscope.
2. Maintain the basket cage in the closed position, while introducing the (distal) working end of the device into the working channel of the endoscope.
3. Carefully advance basket device through the working channel. Using deliberate short strokes (2 to 3 cm), advance the sheath through the endoscope channel while watching for the sheath as it emerges from the distal scope channel.


Figure 2. Thumb Wheel to Open/Close the Basket

Capture

1. Advance the closed basket alongside or beyond the stone and open the basket turning the thumb wheel toward the "OPEN" position. (See Figure 2)
2. Slowly sweep the open basket across the face of the stone to engage the stone into the wire cage. Turning the thumb wheel back and forth to open and close the basket will also facilitate engagement of the stone.
3. The basket cage may also be rotated by turning the basket rotation wheel at the end of the handle. (See Figure 3)
4. Once the stone is captured within the basket, carefully turn the thumb wheel toward the "CLOSE" position to hold the stone. (See Figure 2)


Figure 3. Basket Rotation Wheel to Rotate the Basket

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AL), DFU, MB, OptiFlex, Global, 90998360-01B

Black (K) ΔE ≤5.0

Removal

1. Slowly and carefully, withdraw the basket with the stone fragment(s).
2. If resistance is encountered during advancement or withdrawal of the basket, **DO NOT EXERT EXCESSIVE FORCE**.
3. Some stones may be too large once captured in the basket to withdraw the basket through the endoscope. In this situation, simultaneously remove the endoscope and basket from the urinary tract while carefully maintaining basket position to hold the stone in place.

Handle Disassembly

If handle removal is desired:

1. While holding the basket rotating wheel still, unscrew the end cap at the proximal end of handle, See Figure 4. (A pair of forceps may be used to unscrew the end cap).
2. Slide the handle off the basket shaft.


Figure 4. Handle Disassembly

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