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EN

Instructions for use

NOTE: people with reduced mental/physical capacities should be assisted while using the device.

Here are the steps to follow to ensure accurate use of the pen needle:

- Check the expiry date of the device. Do not use if expired
- Check the compatibility between the pen needle and the pen to be used for administrating the medication, considering the information on the packaging
- Always mount a new needle before each injection, and dispose of it afterwards. As a matter of fact, never store the pen with a needle inserted as this would leave an open path of communication between the cartridge and the outside. The pen needle is a sterile disposable device, therefore, it must be changed after each insertion, otherwise you risk:
 - inaccurate dosage. If the needle is left on the pen between injections, the composition of the medication may be altered;
 - air getting into the cartridge through the needle. The air bubbles reduce the normal flow of medication during the injection, with the consequent risk of a reduced dose of the drug;
 - the crystallization of the medication. If the needle is left mounted on the pen, the medication remaining inside the pen can transform into crystals, and partially or completely obstruct the needle, causing difficulty in or inability to carry out the injection;
 - contracting infections. Once used, the pen needle is no longer sterile
 - experiencing a more painful injection. During the injection, the lubricant put on the needle is completely removed as the tip of the needle undergoes a progressive deformation, losing its original shape. Often the tip of the needle becomes a "hook" shape, causing increased pain and possible tearing of the skin
- Before each injection, make sure the medication comes out properly. Always keep a spare pen needle as a replacement
- Use the technique recommended by your doctor for administration. Once it is inserted, press the injection button all the way down, being careful not to change the angle of the pen. If, after the needle is removed from the skin, you notice medication dripping, keep the needle inserted longer the next time, at the same dosage.
- After the injection, remove the needle from the pen and dispose of it in an appropriate container respecting applicable regulations and laws.

TO REDUCE PAIN use new, short, small calibre needles and always follow your physician's instructions

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Instructions for use

N.B. people with reduced mental/physical capacities should be assisted while using the device.

1. Check the expiry date of the device. Do not use if expired
2. Check the compatibility between the pen needle and the pen to be used for administering the medication, considering the information on the packaging. As a matter of fact, never store the pen with a mounted needle as this would leave an open path of communication between the cartridge and the outside
3. Always mount a new needle before each injection and dispose of it afterwards. The pen needle is a sterile disposable device; therefore, it must be changed after each injection, otherwise you risk:
 - inaccurate dosage. If the needle is left on the pen between injections, the composition of the medication may be altered.
 - air getting into the cartridge through the needle. The air bubbles reduce the normal flow of medication during the injection, with the consequent risk of a reduced dose of the drug.
 - the crystallization of the medication. If the needle is left mounted on the pen, the medication remaining inside the pen can transform into crystals, and partially or completely obstruct the needle, causing difficulty in or inability to carry out the injection.
 - contracting infections. Once used, the pen needle is no longer sterile.
 - experiencing a more painful injection. During the injection, the lubricant, put on the needle is completely removed and the tip of the needle undergoes a progressive deformation, losing its original shape. Often the tip of the needle becomes a "hook" shape, causing increased pain and possible tearing of the skin.
4. Before each injection, make sure the medication comes out properly. Always keep a spare pen needle as a replacement
5. Use the technique recommended by your doctor for administration. Once it is inserted, press the injection button all the way down, being careful not to change the angle of the pen. If, after the needle is removed from the skin, you notice medication dripping, keep the needle inserted longer the next time, at the same dosage
6. After the injection, remove the needle from the pen and dispose of it in an appropriate container respecting applicable regulations and laws

TO REDUCE PAIN use new, short, small caliber needles and always follow your physicians instructions.

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