Occlusafe

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DESCRIPTION

<INDICATIONS>
The Occlusafe is intended to be used for gaining temporary occlusion in the peripheral vasculature. The Occlusafe provides temporary vascular occlusion which is useful in selectively stopping or controlling blood flow.
<DESCRIPTION>
The Occlusafe is a coaxial balloon catheter with an inner lumen used to facilitate the passage of a guidewire and an outer lumen for inflation and deflation of the balloon. The Occlusafe is compatible with 0.014"/0.36 mm or smaller guidewires. The Occlusafe has a proximal stiff part and a distal flexible part coated with hydrophilic polymer which generates lubricity when wet. On the distal tip two radiopaque marker bands are incorporated.

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ENGLISH

Please read all instructions prior to use.

CAUTIONS FOR USE

1. Contraindicated for use in:
- Coronary vessels
- Central circulatory system
 'Central circulatory system' means the following vessels: pulmonary artery, ascending aorta, aortal arch, descending aorta to the aortal bifurcation, coronary artery, common carotid artery, external carotid artery, internal carotid artery, cerebral artery, brachiocephalic artery, coronary vein, pulmonary vein, vena cava superior, vena cava inferior.
- Embolectomy and angioplasty procedures
- Pediatrics and neonates


2. Potential Complications

Complications associated with the use of the Occlusafe are similar to the ones associated with standard percutaneous interventional procedures. Possible complications may include, but are not limited to:
- Local or systemic infection
- Local hematoma
- Intimal disruption
- Arterial dissection, perforation, rupture or injury
- Arterial thrombosis
- Distal embolization of blood clots and plaque
- Arterial spasm
- Arteriovenous fistula formation, Pseudoaneurysm or Bleeding Complication at access site
- Arrhythmias, including life-threatening ventricular fibrillation
- Hemorrhage
- Hypotension
- Death

3. Important safety instructions

WARNINGS
- **If a guiding catheter is used with a stopcock, do not operate this stopcock while the Occlusafe is in the guiding catheter. It may result in breakage or separation of the Occlusafe.**
- **Verify the size of the vessel under fluoroscopy. Always visualize under fluoroscopy to ensure patient safety while inflating and deflating the balloon .**
- **Do not advance the Occlusafe while the balloon is inflated. This may cause vessel rupture and intimal dissection.**
- **Do not exceed the maximum recommended inflation volume as balloon rupture may occur.**
- **Do not torque the Occlusafe if it is bent in order to avoid separation or break of the Occlusafe.**
- **Do not apply high pressure infusion via a guiding catheter working lumen on this product while it is inserted into a guiding catheter. This may compress the catheter portion of the Occlusafe, resulting in unexpected balloon inflation.**
- **Perform appropriate anticoagulant or antiplatelet therapy according to the patient's condition in order to avoid thrombotic embolization.**
- **While this product is inserted into the guiding catheter, continuously flush the guiding catheter lumen with heparinized saline solution. Otherwise, vascular embolization may occur.**
- **Do not advance the Occlusafe or guidewire against resistance. If resistance is felt, determine the cause under fluoroscopy.**
- **The Occlusafe has not been tested for compatibility or use with liquid embolics.**

CAUTIONS
- Inspect the catheter prior to use for any bends or kinks. Do not use a damaged catheter.
- Exercise care during handling in order to avoid possible damage to the catheter. Avoid acute bending or kinking of the catheter itself.
- Any advancement after introduction of the catheter into the vessel should be done under fluoroscopy. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Excessive tightening of a hemostatic valve on the catheter shaft may result in damage to the catheter or could delay balloon inflation and deflation.
- Verify the Occlusafe compatibility when using other ancillary devices commonly used in intravascular procedures. Physician must be familiar with percutaneous, intravascular techniques and possible complications associated with the procedure.
- Do not reshape the Occlusafe by any means. Attempting to reshape the catheter may cause damage, resulting in the release of catheter fragments into the vessel.
- Always keep observing the inflation state of the balloon during injection under fluoroscopy.
- Use the syringe provided with this product to inflate/deflate the balloon. Use of syringes other than the syringe provided with the Occlusafe may cause variation in the volume used for inflation, or excessive negative pressure, resulting in damage to the catheter.
- A mixed solution containing the same proportions of contrast medium (containing 300 mg/mL of iodine) and heparinized saline solution shall be used as an inflation solution (herein after referred to as "inflation solution"). Otherwise, the balloon may not inflate/deflate correctly.
- Always hold the hub as shown in the figure below when attaching or removing the syringe.

Hold the hub when attaching or removing the syringe.

- If any excessive resistance is felt during injection, discontinue the manipulation.
- Do not soak or clean the Occlusafe in agents containing organic solvents, such as alcohol for disinfection. Contact with these agents may lead to damage or breakage of the Occlusafe and loss of lubricity.

4. Precautions for application.
- The Occlusafe is sterilized and intended for single use only.
- The Occlusafe should be used only by physicians trained in percutaneous, intravascular techniques and well trained in the use of the device.
- The Occlusafe has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not re-sterilize. Do not re-process. Re-processing may compromise the sterility, biocompatibility and functional integrity of the device.
- Sterile and non-pyrogenic in an unopened and undamaged unit package. Do not use if the unit package or the Occlusafe have been damaged or soiled.
- The Occlusafe should be used immediately after opening the package and be disposed of safely and properly after use in accordance with hospital, administrative and/or local government policy.
- The entire operation should be carried out aseptically.

DIRECTIONS FOR USE

The following instructions provide technical direction but do not take the place of formal training in the use of the Occlusafe. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

1 PREPARATION FOR USE

1.1 Before removing the Occlusafe from the holder tube, fully hydrate the hydrophilic segment of the device by flushing heparinized saline solution through the holder tube using a syringe attached to the flushing luer.


1.2 To remove the Occlusafe from the holder tube, gently pull the hub out from the holder tube. Remove the Occlusafe by pulling it from the holder tube. If resistance is observed, repeat the flushing procedure until the Occlusafe is well hydrated and can be easily removed from the holder tube. Inspect the Occlusafe thoroughly to ensure it is not damaged. Do not allow the Occlusafe to dry prior to introduction into the guiding catheter. Do not reinsert a hydrated balloon catheter into its holder tube.

2 Balloon Preparation

2.1 Insert the enclosed stylet into the inner lumen from the distal tip of the Occlusafe to avoid any damages of the Occlusafe during the balloon preparation.

2.2 Attach the three-way stopcock and a 20mL syringe (available in the hospital) to the balloon inflation port.


2.3 To purge air from the balloon, apply negative pressure to the balloon lumen using the syringe and completely deflate the balloon (left lower figure), and close the three-way stopcock (right lower figure).

2.4 To purge air from the three-way stopcock, fill the 1.0mL syringe (provided with the Occlusafe) with "inflation solution", attach the syringe to the three-way stopcock, and fill the three-way stopcock with the inflation solution (figure below).

CAUTION


A mixed solution containing the same proportions of contrast medium (containing 300 mg/mL of iodine) and heparinized saline solution shall be used as an inflation solution. Otherwise, the balloon may not inflate/deflate correctly.

2.5 When the three-way stopcock is opened, inflation solution equal to the inflation prime volume will flow into the catheter by negative pressure of the balloon lumen (figure below). When flow ceases, additionally infuse 0.2 mL of the inflation solution slowly, using the syringe.

CAUTION

Do not infuse rapidly. This may cause damage to the balloon or the catheter.

2.6 Holding the distal end of the catheter downward, apply negative pressure to the balloon with the syringe, to completely deflate the balloon and to purge the remaining air from balloon lumen, then close the three-way stopcock.

2.7 Open the three-way stopcock and the balloon lumen will be filled again with inflation solution equal to the inflation priming volume, then slowly infuse the recommended balloon inflation volume, and check the balloon for any abnormality.

Note: Ensure that there is no air left inside the balloon. Residual air may compromise fluoroscopic visualization.

2.8 With the balloon completely deflated, close the three-way stopcock.

2.9 Fill the syringe with the recommended balloon volume inflation solution, and attach the syringe to the three-way stopcock.

2.10 Then remove the stylet from the distal tip of the Occlusafe, flush the Occlusafe guide wire lumen lumen with heparinized normal saline.

2.11 Immerse this product in a tray filled with heparinized normal saline, and ensure that the lubricity of the outer surface of this product increases.

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3 Direction for use

3.1 Attach a hemostatic valve (marketed product) to a compatible guiding catheter (marketed product), and flush the guiding catheter lumen with heparinized saline solution.

WARNING

While this product is inserted into the guiding catheter, continuously flush the guiding catheter lumen with heparinized saline solution. Otherwise, vascular embolization may occur.

3.2 After advancing the guiding catheter into the target site, appropriately release the hemostatic valve to permit insertion of a compatible guidewire (marketed product), and then carefully insert the Occlusafe. When using a hemostatic valve, use an inserter to insert a compatible guidewire, and then insert the Occlusafe over the guide wire.

Note: Deflate the balloon completely. Otherwise, the Occlusafe will not pass easily through the guiding catheter.

Note: Sufficiently wet the outer surface of the Occlusafe to increase its lubricity. Otherwise, resistance will readily occur during insertion of the Occlusafe.

Note: Use a guidewire concomitantly when inserting the Occlusafe . Otherwise, kinking of this product occurs easily.

3.3 Advance the guidewire under fluoroscopic guidance, and then insert the Occlusafe to the target site, over the guidewire.

3.4 Adjust the position of the balloon while confirming the position of the radiopaque marker band.

3.5 Inflate the balloon slowly under fluoroscopy.

WARNING

Verify the size of the vessel under fluoroscopy. Always inflate and deflate the balloon under fluoroscopy to ensure patient safety.

WARNING

Do not move the Occlusafe while the balloon is inflated. This may cause vessel rupture and intimal detachment.

WARNING

Do not exceed the maximum balloon recommended inflation volume as balloon rupture may occur.

3.6 On completion of the procedure, deflate the balloon completely and remove it slowly.


WARNING

If a guiding catheter is used with a stopcock, do not operate this stopcock while the Occlusafe is in the guiding catheter. It may result in breakage or separation of the Occlusafe.

Note: Ensure that the balloon is deflated completely. Otherwise, the balloon may be unable to be removed from the guiding catheter.

Reference Values

Indication on Label	Meaning	Reference Values		
∅	Balloon diameter	4.0 mm		
—	Balloon length	10 mm		
L	Catheter effective length	110cm	130cm	150cm
<Vol>	Balloon recommented inflation volume	0.10 mL		
Max.Volume	Balloon maximum inflation volume	0.20 mL		
IPV	Inflation prime volume	0.34mL	0.36 mL	0.38mL
	Minimum inner diameter	0.43 mm (0.017 inch)		
—	Catheter distal outer diameter	2.7 Fr (0.90mm)		
—	Catheter proximal outer diameter	2.8 Fr (0.94mm)		
Min GC	Minimum guiding catheter inner diameter	1.02mm (0.040")		
Max Guide Wire O.D.	Maximum guide wire outer diameter	0.36mm (0.014")		

STORAGE

Avoid exposure to water, direct sunlight, extreme temperature, or high humidity during storage.

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
Indication on Label	Meaning	Reference Values		
∅	Balloon diameter	4.0 mm		
—	Balloon length	10 mm		
L	Catheter effective length	110cm	130cm	150cm
<Vol>	Balloon recommented inflation volume	0.10 mL		
Max.Volume	Balloon maximum inflation volume	0.20 mL		
IPV	Inflation prime volume	0.34mL		0.38mL
	Minimum inner diameter	0.43 mm (0.017 inch)		
—	Catheter distal outer diameter	2.7 Fr (0.90mm)		
—	Catheter proximal outer diameter	2.8 Fr (0.94mm)		
Min GC	Minimum guiding catheter inner diameter	1.02mm (0.040")		
Max Guide Wire O.D.	Maximum guide wire outer diameter	0.36mm (0.014")		

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© TERUMO CLINICAL SUPPLY Co., Ltd. 2015-09
OCCLUSAFE is a trademark of TERUMO CLINICAL SUPPLY Co., Ltd..
TERUMO is a trademark of TERUMO CORPORATION.
TERUMO CLINICAL SUPPLY Co., Ltd.