LIVEWIRE™ TC
CATHETER EXTENSION CABLE

English - Instructions for use 2
Deutsch - Gebrauchsanleitung 5
Français - Mode d'emploi 7
Italiano - Istruzioni per l'uso. 9
Español - Instrucciones de uso 11
Polski - Instrukcja obsługi 13
Türkçe - Kullanım talimatları 15
中文 - 使用说明 17
Αγγλικά - Οδηγίες χρήσης 19
Русский - Инструкция по эксплуатации 21
Magyar - Használati utasítás 23
Português - Instruções de Utilização 25
Nederlands - Gebruiksaanwijzing 27
Svenska - Bruksanvisning 29
Norsk - Bruksanvisning 31
Suomi - käyttöohjeet 33
Dansk - Brugsanvisning 35
Česky - Návod k použití 37
Eesti keeles - Kasutusjuhend. 39
Română - Instrucţiuni de utilizare 41
한국어 – 사용 설명서 43
اللغة العربية - تعليمات الاستخدام 45
繁體中文 - 使用說明書 47
SYMBOL 49

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LIVEWIRE™ TC
CATHETER EXTENSION CABLE
Instructions for use

Read Instructions for Use prior to use of this device.

SEE INDIVIDUAL STERILE PACKAGE LABEL FOR CONTENTS.

CAUTION: FEDERAL (U.S.) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

DESCRIPTION AND APPLICATION

The St. Jude Medical (SJM™) Livewire™ TC Catheter Extension Cable is an ablation extension cable with conductors available for temperature monitoring. The SJM Livewire™ TC Catheter Extension Cable is designed to connect SJM Livewire™ TC Steerable Electrophysiology Ablation Catheters equipped with temperature monitoring to a compatible RF generator. The special indexing connector located at the distal end of the cable is intended for connection to SJM Livewire™ TC Ablation Catheters. The proximal end connector with special indexing is intended for connection to the generator. Ensure that the connection is secure between the cable and connecting equipment.

The compatible generator for your SJM Livewire™ TC Catheter Extension Cable can be found in the table below.

Cable	Reorder Number	Length (cm)	RF Generator
SJM Livewire™ TC Catheter Extension Cable	402505	300	STOCKERT EP-Shuttle™ Generator
	402543	150	IBI-1500T Cardiac Ablation Generator or Ampere™ RF Ablation Generator
	402544	300	

CAUTION

- **Do not** alter this device.
- Observe polarity.
- Inspect components prior to use.
- **STORE IN A COOL, DARK, DRY PLACE.**

CONTRAINDICATIONS

- There are no known contraindications for this device.

CABLE REUSE

- Properly clean and re-sterilize cables within twenty-four (24) hours after use.
- Do not submerge cables in an enzymatic detergent or water bath during transport. Avoid contact with contaminated equipment or instruments during cable transport.
- Cables are intended to be reused a maximum of ten (10) times. Discard cables that have reached this number of uses.
- Inspect cables prior to reuse. Do not use cables if there is visible damage. Discard damaged cables.

Note: *To ensure cables are properly cleaned and externally sterile, use the specified cleaning and sterilization methods/parameters. Deviations from the specified methods/parameters may result in reduced instrument life or inadequate cleaning or sterilization. If alternate methods/parameters are used, it is the customer's responsibility to ensure cables are effectively cleaned and sterilized. The sterilization and cleaning methods listed in this manual are only valid when using properly installed, maintained, calibrated and compliant cleaning and sterilization equipment.*

CABLE CLEANING METHODS

MANUAL METHOD

1. Soak the cable for five (5) minutes, fully submerged in ENZOL™ Enzymatic Detergent (or equivalent), diluted according to the manufacturer's recommendations with warm 31°–43° C (88°–110° F) tap water. Remove visible soil with a soft-bristled brush.
2. Rinse the cable with clean warm running tap water for a minimum of one (1) minute.
3. Dry with a clean, lint-free cloth.
4. After cleaning, visually inspect cables for any remaining soil or moisture. Repeat cleaning and/or drying if necessary.

 Note: *Cable must be thoroughly dried before sterilization.*

AUTOMATED METHOD

Preparation for Automated Cleaning:

Prior to putting cables in the washer/disinfector, perform the following steps to ensure that any organic matter on the cables is rehydrated.

1. Soak the cable for one (1) minute, fully submerged in ENZOL™ Enzymatic Detergent (or equivalent), diluted according to the manufacturer's recommendations with warm 31°–43° C (88°–110° F) tap water.
2. Remove visible soil with a soft-bristled brush.
3. Rinse the cable with clean, warm 31°–43° C (88°–110° F), running tap water for a minimum of one (1) minute.

Automated Cleaning:

1. Load cables into a tray in a manner that is consistent with the manufacturer's instructions. Do not place heavy objects on top of cables or place cables next to sharp objects during an automated cleaning cycle. Doing so may damage the cables.
2. Set the automated washer/disinfector RELIANCE™ 444 (or legally marketed equivalent) to the following settings prior to running the cleaning cycle.

 Note: *No thermal disinfection is recommended.*

Cycle	Time	Water Temperature (Tap or Purified)	Detergent & Concentration
Pre-Wash	1 min	Cold/Warm Approx.16°–43°C (61°–110° F)	N/A
Enzymatic Wash	1 min 30 sec	Warm/Hot Approx. 43°–60°C (110°–140° F)	PROLYSTICA™ 2X Concentrate Enzymatic Presoak and Cleaner (or equivalent), diluted per manufacturer's instructions
Wash	2 min 30 sec	Heated 60°–70°C (140°–158°F)	PROLYSTICA™ 2X Concentrate Neutral Detergent (or equivalent), diluted per manufacturer's instructions

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Cycle	Time	Water Temperature (Tap or Purified)	Detergent & Concentration
Final Rinse	1 min 30 sec	Heated 60˚–70˚C (140˚–158˚F)	N/A
Drying	1 hour	Air temperature 70˚–80˚C (158˚–176˚F)	N/A

3. After cleaning, visually inspect cables for any remaining soil or moisture. Repeat cleaning and/or drying if necessary.

STERILIZATION METHODS

EtO STERILIZATION

Place cable into a breathable package suitable for EtO sterilization (e.g. TYVEK™ pouch or equivalent). Load cables into chamber and run equipment in a manner that is consistent with the manufacturer's instructions.

St. Jude Medical has validated the following cycle parameters.

Cycle Parameters	55˚C Temperature Cycle
Pre-Conditioning Chamber	
Pre-Conditioning Time	30–45 min
Relative Humidity	40%–80%
Ethylene Oxide (EtO) Process	
Exposure Time	1 hour
Vacuum Set	1.3 psia
EtO Gas	100% EtO
EtO Concentration	725 mg/L–880 mg/L
Aeration	
Time	12 hours

STERRAD™ STERILIZATION SYSTEMS

Ensure cables are thoroughly dried before loading into the STERRAD™ System to avoid a potential cycle cancellation. Package the cable per facility practice in an appropriately sized single pouch (e.g. TYVEK™ pouch or equivalent), tray, or container labeled as appropriate for the STERRAD™ hydrogen peroxide gas plasma systems. Do not overlap cables. Load cables into STERRAD™ System and run equipment in a manner that is consistent with the manufacturer's instructions.

The following STERRAD™ Systems and cycles are recommended for reprocessing.

STERRAD™ System and Cycle within the United States	STERRAD™ System and Cycle outside United States
STERRAD™ 50	STERRAD™ 50
STERRAD™ 200	STERRAD™ 200 (Short cycle)
STERRAD™ 100S	STERRAD™ 100S (Short cycle)
STERRAD™ NX™ (Standard cycle)	STERRAD™ NX™ (Standard cycle)
STERRAD™ 100NX™ (Standard cycle)	STERRAD™ 100NX™ (Standard cycle)

NOTE: *The symbols section contains all the symbols that may be used on product labels. Product is labeled as required.*

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USA ONLY

LIMITED WARRANTY AND DISCLAIMER

St. Jude Medical (SJM) hereby warrants that if any SJM product fails to perform within normal tolerances for a patient due to a defect in materials or workmanship, SJM will provide, at no charge, a replacement SJM product for the patient's use. This limited warranty applies only if each of the following conditions are met.

1. The product was packaged and labeled by SJM.
2. The failed product must be returned to SJM and becomes the property of SJM.
3. The product has not been mishandled, reprocessed or altered in any way.
4. The product was used before the "USE BY" date marked on the packaging of the product.

No representation or warranty is made that a SJM product will not fail. SJM disclaims responsibility for any medical complications, including death, resulting from the use of its products. Except as expressly provided by this limited warranty, SJM IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE, OR MALFUNCTION OF ITS PRODUCTS, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE. Some states do not allow the exclusion or limitation of incidental or consequential damages however, so the above limitation or exclusion may not apply to you.

Except as expressly provided by the limited warranty, SJM MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE.

STATE OF CALIFORNIA (USA ONLY)

WARNING: This product and its packaging have been sterilized with ethylene oxide. The packaging may expose you to ethylene oxide, a chemical known to the State of California to cause cancer or birth defects or other reproductive harm.

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Australian Sponsor:
St. Jude Medical Australia Pty Limited
17 Orion Rd
Lane Cove NSW 2066
AUSTRALIA
+61 (02) 9936 1200

St. Jude Medical
14901 DeVeau Place
Minnetonka, MN
55345-2126 USA
+1 855 478 5833
+1 651 756 5833
sjm.com

EC REP

St. Jude Medical
Coordination Center BVBA
The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem
Belgium
+32 2 774 68 11

IFU: 100108513
ART: ARTMT100108513 C
2019-06

EP-Shuttle is a trademark of Stockert Instrumente GmbH.
ENZOL, STERRAD, STERRAD NX, and STERRAD 100NX are trademarks of Johnson & Johnson Corporation or a related company.
TYVEK is a trademark of E.I. Du Pont de Nemours and Company Corporation.
RELIANCE and PROLYSTICA are trademarks of Steris Inc.
Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL and the nine-squares symbol are trademarks and services marks of St. Jude Medical, Inc. and its related companies.
© 2014 St. Jude Medical. All rights reserved.

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