DePuy Synthes
MITEK SPORTS MEDICINE
PART OF THE Johnson & Johnson FAMILY OF COMPANIES

HEALIX ADVANCE™ KNOTLESS Suture Anchor

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P/N: 116046
Rev: B Issued 12/2019
© DePuy Synthes 2018-2019.
All rights reserved.

ENGLISH

HEALIX ADVANCE™ KNOTLESS Suture Anchor

℞ Only

CONTENTS

Contents are STERILE unless packaging is damaged or opened.

DESCRIPTION

The HEALIX ADVANCE™ KNOTLESS Anchors are one piece implants, consisting of either:

- PEEK (Polyetheretherketone)

 OR

- BIOCRYL RAPIDE™ (BR), an absorbable composite material (PLGA/TCP).

See package label for anchor material.

The HEALIX ADVANCE KNOTLESS Anchors are sterile implants designed to secure soft tissue to bone when used in conjunction with high-strength #2 suture or PERMATAPE™ 2.5mm suture.

MATERIALS

Anchor: BIOCRYL RAPIDE (BR) or PEEK (Polyetheretherketone)

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Threader Tab: Polycarbonate and stainless steel

Stay Suture: Orthocord® Suture is a synthetic, sterile, braided, composite suture composed of dyed (D&C Violet #2 or D&C Blue #6) absorbable polydioxanone (PDS®) and un-dyed non-absorbable polyethylene. The partially absorbable suture is coated with a copolymer of 90% caprolactone and 10% glycolide.

INDICATIONS

The HEALIX ADVANCE KNOTLESS Anchors are indicated for use in the following procedures for reattachment of soft tissue to bone:

	Indication	BR	PEEK
Shoulder	Rotator Cuff	X	X
	Biceps Tenodesis	X	X
	Deltoid Repair	X	X
Foot/Ankle	Achilles Tendon	X	X

CONTRAINDICATIONS

- Procedures other than those listed in the INDICATIONS section.
- Pathologic conditions of bone such as cystic changes or severe osteopenia that would impair its ability to securely fix the HEALIX ADVANCE KNOTLESS Anchors.
- Pathologic changes in the soft tissues being fixated to bone that would prevent their secure fixation by the HEALIX ADVANCE KNOTLESS Anchors.
- Comminuted bone surface that would militate against secure fixation of the HEALIX ADVANCE KNOTLESS Anchors.
- Physical conditions that would eliminate or tend to eliminate adequate implant support or retard healing, i.e., blood supply limitation, previous infection, etc.
- Conditions which tend to limit the patient's ability to restrict activities or follow directions during the healing period.
- The HEALIX ADVANCE KNOTLESS Anchors are not designed for and should never be used to attach artificial ligaments.

PRECAUTIONS

- The HEALIX ADVANCE KNOTLESS Anchors are supplied STERILE, and are intended for single use only. Do not re-sterilize. Do not use if sterile packaging appears to be damaged.
- The product should be stored under cool, dry conditions.
- Inspect all instruments for damage before use. Do not attempt to repair.

Property of Medical Device
www.smaller.tn | info@smaller.tn

- A surgeon should not begin clinical use of the HEALIX ADVANCE KNOTLESS Anchors without reviewing the Instructions for Use and practicing the procedure in a skills laboratory.
- Ensure that Anchor is inserted axially to the bone hole (limit off-axis insertion).
- Ensure area of anchor insertion is free of soft tissue.
- Improper instrument use, axial misalignment or levering with the anchor upon insertion may result in anchor fracture or reduced performance.

WARNINGS

- DePuy Mitek anchors are designed to lock into cortical or cancellous bone. Bone stock must be adequate to allow proper and secure anchor placement. Incomplete insertion or over-tensioning or poor bone quality may result in anchor pullout.
- Immediate range of motion should be avoided to allow biological bony/soft tissue healing.
- This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the thoracic or lumbar spine.
- Users should be familiar with surgical procedures and techniques involving absorbable and non-absorbable suture before employing ORTHOCORD suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.
- This product is for single-use only. It has not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use instruments can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.

ADVERSE EFFECTS

Adverse effects of absorbable implanted devices include mild inflammatory and foreign body reactions.

MR STATEMENT

The HEALIX ADVANCE KNOTLESS Implant(s) and Sutures are MR SAFE.

INSTRUCTIONS FOR USE

1. Remove the HEALIX ADVANCE KNOTLESS Anchor from the sterile package.
2. Pass the suture through soft tissue with the surgeon's preferred suture-passing technique. Keep the sutures at equivalent lengths after passing/tying.
3. Clean the surface of excess tissue around the point of insertion, do not de-corticate surface.
4. Create the hole as follows:
 a. Awl – Use a mallet to impact the distal tip of the instrument into bone until the laser line is flush with bone. Remove the awl and set aside.

 OR

 b. Drill – Drill into the bone until the distal shoulder on the drill is flush with the bone. Remove the drill and set aside.
5. In hard bone, it is necessary to tap the pre-drilled bone hole. Insert the distal tip of the appropriate tap into the pre-drilled hole until the threads are flush with the bone. Rotate the tap in a clockwise direction until the laser line is flush with the bone. Rotate the tap in a counter-clockwise direction to remove.
6. Load the desired suture(s) through the Wire Kite(s). The suture tails should extend approximately 1" through the Wire Kite(s).

 NOTE: For Achilles Tendon repair with BR Anchor load 4 to 6 #2 high-strength suture tails or 2 to 3 PERMATAPE 2.5mm Suture Tails.
7. Pull Threader Tab, which is attached to the Wire Kite until all sutures have been pulled through the Anchor, Driver Collar, and out of the slot in the Driver. Discard the Threader Tab.

 NOTE: To reload the Threader Tab, slide the Driver Collar to the distal end of the slot in the Driver and thread the Wire Kite through the Driver Collar, into the slot of the Driver. Feed the Wire Kite until it appears through the distal tip of the Anchor.

Permission granted to reproduce for personal use only.
www.arthrex.com | info@arthrex.com

8. For open procedures, slide the collar towards the handle. Turn collar counter-clockwise until it locks into notch.
9. While holding Sutures under slight tension, slide Anchor Assembly down the Sutures toward the insertion site.

10. While keeping desired tension on the suture strands, place the distal nose of the Anchor into the bone hole.
11. If more tension is desired:
 a. Individually tension sutures, while holding anchor nose in hole, until proper tension is achieved.
 b. Remove anchor nose from the bone hole and tension sutures. Re-insert the Anchor nose back in hole.
12. Once proper tension is achieved, release Sutures.

 NOTE: Maintaining the tension on the sutures exiting the driver while inserting the anchor will result in reduced fixation.
13. Begin inserting the Anchor by turning the Handle clockwise and applying slight downward pressure until the laser line is flush with the bone.
14. Open the slide cover on the Handle and remove suture card.
15. Pull the handle upwards to remove the shaft from the anchor.

STORAGE

Store below 25°C (77°F). Do not use after expiration date.

STERILIZATION

This product is provided sterile. Do not re-sterilize. A DePuy Mitek ANCHOR MUST NEVER BE REUSED.

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