Sensor Enabled™ Ablation Connection Cable

REF A-FASE-CBL4

Instructions for Use

Abbott

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LIMITED WARRANTY AND DISCLAIMER

Abbott Medical hereby warrants that if any Abbott Medical product fails to perform within normal tolerances for a patient due to a defect in materials or workmanship, Abbott Medical will provide, at no charge, a replacement Abbott Medical product for the patient's use. This limited warranty applies only if each of the following conditions are met.

1. The product was packaged and labeled by Abbott Medical.

2. The failed product must be returned to Abbott Medical and becomes the property of Abbott Medical.

3. The product has not been mishandled, reprocessed or altered in any way.

4. The product was used before the "USE BY" date marked on the packaging of the product.

No representation or warranty is made that an Abbott Medical product will not fail. Abbott Medical disclaims responsibility for any medical complications, including death, resulting from the use of its products. Except as expressly provided by this limited warranty, ABBOTT MEDICAL IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE, OR MALFUNCTION OF ITS PRODUCTS, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE. Some states do not allow the exclusion or limitation of incidental or consequential damages however, so the above limitation or exclusion may not apply to you.

Except as expressly provided by the limited warranty, ABBOTT MEDICAL MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE.

Proposition 65, a State of California voter initiative, requires the following notice:

⚠ WARNING: This product can expose you to chemicals including ethylene oxide, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

™ Indicates a trademark of the Abbott group of companies.
‡ Indicates a third party trademark, which is property of its respective owner.

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Instructions for Use

Caution

- United States law restricts this device to sale by or on order of a physician.
- Read directions prior to use.

Description

Electrophysiology cables are flexible, insulated cables constructed of metal and plastic materials.

Indications

The connector cable is intended for use with a compatible radiofrequency generator and sensor enabled catheter positioning and navigation applications.

Directions

1. Verify the connector cable is at room temperature prior to use.
2. Inspect the connector cable carefully for integrity and overall condition.
3. Use the connection with the green strain relief to connect to the FlexAbility ™ Ablation Catheter, Sensor Enabled™.
4. When connecting to an RF generator, refer to the RF generator instructions for use for details and diagrams on how to connect the cable to the RF generator. Use the connection with the yellow strain relief to connect to the RF generator.
5. Use the connection with the white strain relief to connect to the navigation and positioning systems.
 - For the MediGuide™ System, connect the white strain relief to the MediGuide™ CathConnect, Sensor Enabled™.
 - For the EnSite Precision™ System, connect the white strain relief to the EnSite Precision™ Link, Sensor Enabled™.
6. Ensure that the connection is secure between the cable and connecting equipment.

Warnings

- Ensure the connector cable is at room temperature before use to ensure accurate temperature measurement.
- No modification of this equipment is allowed.
- Connected equipment must be patient isolated and current leakage must not exceed 10 microamps.

Precautions

- The connector cable should only be used by or under the supervision of physicians trained in ablation procedures using compatible navigation and visualization technology systems.
- Personnel handling the connector cable should wear gloves.
- Standard grounding procedures should be followed if electrosurgical instruments are used.
- The connector cable has been evaluated at a maximum voltage of 240 volts.

Packaging and Shelf Life

The packaging is designed to provide for aseptic product transfer. It is recommended that the product remain in the unopened package until time of use. Contents are sterile if package is unopened and undamaged. The expiration date is marked on the outside of the package. The product should be stored in a cool, dry location. The instructions for use is recyclable. Dispose of used product and packaging following standard solid biohazard waste procedures.

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Storage and Transport

Keep each connector cable in a sealed pouch in a cool, dry location away from sunlight, including UV light, during storage and transport.

Cable Reuse

- Properly clean and re-sterilize cables within twenty-four (24) hours after use.
- Do not submerge cables in an enzymatic detergent or water bath during transport. Avoid contact with contaminated equipment or instruments during cable transport.
- Cables are intended to be reused a maximum of ten (10) times. Discard cables that have reached this number of uses.

 Note: It is the customer's responsibility to track the number of reuses of the cable. Customers should follow their established procedures for reuse tracking.
- Inspect cables prior to reuse. Do not use cables if there is visible damage (e.g., cuts, bulges, loose strain reliefs, or connectors, etc.). Discard damaged cables.

 Note: To ensure cables are properly cleaned and externally sterile, use the specified cleaning and sterilization methods/parameters. Deviations from the specified methods/parameters may result in reduced instrument life or inadequate cleaning or sterilization. If alternate methods/parameters are used, it is the customer's responsibility to ensure cables are effectively cleaned and sterilized. The sterilization and cleaning methods listed in this manual are only valid when using properly installed, maintained, calibrated and compliant cleaning and sterilization equipment. Customers are responsible for ensuring that the cleaning and sterilization methods recommended can be duplicated in their environment.

Cable Cleaning Methods

Manual Method

1. Soak the cable for five (5) minutes, fully submerged in diluted ENZOL‡ Enzymatic Detergent (or equivalent). Dilute the detergent according to the manufacturer's recommendations with warm 31°–43°C (88°–110°F) tap water. Remove visible soil with a soft-bristled brush.
2. Rinse the cable with clean warm running tap water for a minimum of one (1) minute.
3. Dry with a clean, lint-free cloth.
4. After cleaning, visually inspect cables for any remaining soil or moisture. Repeat cleaning and/or drying if necessary.

 Note: Cable must be thoroughly dried before sterilization.

Automated Method

Preparation for Automated Cleaning:

Prior to putting cables in the washer/disinfector, perform the following steps to ensure that any organic matter on the cables is rehydrated.

1. Soak the cable for one (1) minute, fully submerged in diluted ENZOL‡ Enzymatic Detergent (or equivalent). Dilute the detergent according to the manufacturer's recommendations with warm 31°–43°C (88°–110°F) tap water.
2. Remove visible soil with a soft-bristled brush.
3. Rinse the cable with clean, warm 31°–43°C (88°–110°F), running tap water for a minimum of one (1) minute.

Automated Cleaning:

1. Load cables into a tray in a manner that is consistent with the manufacturer's instructions. Do not place heavy objects on top of cables or place cables next to sharp objects during an automated cleaning cycle. Doing so may damage the cables.
2. Set the automated washer/disinfector RELIANCE‡ 444 (or legally marketed equivalent) to the following settings prior to running the cleaning cycle.

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Note: No thermal disinfection is recommended.

Cycle	Time	Water Temperature (Tap or Purified)	Detergent & Concentration
Pre-Wash	1 min	Cold/Warm Approx.16°–43°C (61°–110°F)	N/A
Enzymatic Wash	1 min 30 sec	Warm/Hot Approx. 43°–60°C (110°–140°F)	PROLYSTICA‡ 2X Concentrate Enzymatic Presoak and Cleaner (or equivalent), diluted per manufacturer's instructions
Wash	2 min 30 sec	Heated 60°–70°C (140°–158°F)	PROLYSTICA‡ 2X Concentrate Neutral Detergent (or equivalent), diluted per manufacturer's instructions
Final Rinse	1 min 30 sec	Heated 60°–70°C (140°–158°F)	N/A
Drying	1 hour	Air temperature 70°–80°C (158°–176°F)	N/A

3. After cleaning, visually inspect cables for any remaining soil or moisture. Repeat cleaning and/or drying if necessary.

Sterilization Methods

EO Sterilization

1. Ensure cables are thoroughly dried before loading into the EO sterilizer.
2. Coil cable loosely and place into an appropriately-sized, breathable package suitable for EO sterilization (e.g., TYVEK‡ pouch or equivalent). Do not place labels on the breathable side of the package as this may impede sterilant penetration.
3. Load cables into chamber and run equipment in a manner that is consistent with the EO sterilizer manufacturer's instructions.

Abbott Medical has validated the following cycle parameters.

	Cycle Parameters	55°C Temperature Cycle
Pre-Conditioning Chamber	Pre-Conditioning Time	30–45 min
	Relative Humidity	40%–80%
Ethylene Oxide (EtO) Process	Exposure Time	1 hour
	Vacuum Set	1.3 psia
	EtO Gas	100% EtO
	EtO Concentration	725 mg/L–880 mg/L
Aeration	Time	12 hours

STERRAD‡ Sterilization Systems

1. Ensure cables are thoroughly dried before loading into the STERRAD‡ System to avoid a potential cycle cancellation.
2. Coil cable loosely and package the cable per facility practice in an appropriately sized pouch (e.g., TYVEK‡ pouch or equivalent), tray, or container labeled as appropriate for STERRAD‡ hydrogen peroxide gas plasma systems. Do not overlap cables.
3. Load cables into STERRAD‡ System and run equipment in a manner that is consistent with the STERRAD‡ sterilizer manufacturer's instructions.

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The following STERRAD‡ Systems and cycles are recommended for reprocessing.

STERRAD‡ System and Cycle within the United States	STERRAD‡ System and Cycle outside United States
STERRAD‡ 50	STERRAD‡ 50
STERRAD‡ 200	STERRAD‡ 200 (Short cycle)
STERRAD‡ 100S	STERRAD‡ 100S (Short cycle)
STERRAD‡ NX‡ (Standard cycle)	STERRAD‡ NX‡ (Standard cycle)
STERRAD‡ 100NX‡ (Standard cycle)	STERRAD‡ 100NX‡ (Standard cycle)

Steam Autoclave

1. Ensure cables are thoroughly dried before loading into the autoclave.
2. Handle cables in the same manner as other reusable instruments that require packaging or wrapping during sterilization.
3. Coil cable loosely and place in an appropriately sized package that is labeled for steam sterilization use (e.g. sterilization wrap, TYVEK‡ pouch, or equivalent).
4. Load cables and run equipment in a manner that is consistent with the autoclave sterilizer manufacturer's instructions.

Abbott Medical has validated the following cycle parameters. Do not place labels on the breathable side of the package as this may impede sterilant penetration.

	Cycle Parameters	Low Temperature Cycle	High Temperature Cycle
Dynamic/Pre-Vac	Temperature	121°C (250°F)	132°C–135°C (270°F–275°F)
	Time	40 min	4 min
Gravity Displacement	Temperature	121°C (250°F)	132°C–135°C (270°F–275°F)
	Time	40 min	10 min

Note: Sterilization temperatures must not exceed 135°C (275°F).

Symbols

Note: The symbols section contains all the symbols that may be used on product labels. Product is labeled as required.

Symbol	Description
	Caution, consult accompanying documents
	Quantity
	Keep dry
	Keep away from sunlight
	Catalog number
	Batch Code

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Symbol	Description
UDI	Unique Device Identification
X	Number of electrodes
LENGTH	Length
	Dispose of hardware in accordance with local law
STERILE EO	Sterilized with Ethylene Oxide Gas.
Rx ONLY	CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician.
	Manufacturer
	Date of Manufacture
EC REP	Authorized European representative
	Consult instructions for use
medical.abbott/manuals	Follow instructions for use on this website
Diagnostic Catheter Cable	Diagnostic Catheter Cable
Ablation Connection Cable	Ablation Connection Cable
For use with 12-pole Electrophysiology Catheters	For use with 12-pole electrophysiology catheters
	Do not use if package is damaged
	Use by
CE 2797	Conformité Européenne (European Conformity). Affixed in accordance with European Council Directive 93/42/EEC (NB 2797) and 2011/65/EU. Hereby, Abbott Medical declares that this device is in compliance with the essential requirements and other relevant provisions of this directive.
	Temperature limitations
	Humidity limitation

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© 2019 Abbott. All Rights Reserved.