Boston Scientific

Instructions for Use 2

51105206-01

2020-07

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 92238512A), IFU, MB, STARTER, Global, 51105206-01A_pretrans

Black (K) ΔE ≤5.0

TABLE OF CONTENTS

WARNING 3

DEVICE DESCRIPTION 3

TIP STYLE 3

INTENDED USE 3

CONTRAINDICATIONS 3

WARNINGS 3

PRECAUTIONS 4

ADVERSE EVENTS 4

HOW SUPPLIED 4

Handling and Storage 4

PREPARATION TECHNIQUE 4

OPERATIONAL INSTRUCTIONS 4

REFERENCES 6

WARRANTY 7

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 92238512A), IFU, MB, STARTER, Global, 51105206-01A_pretrans

Black (K) ΔE ≤5.0

Starter™
Guidewire

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

Starter guidewires are coated stainless steel guidewires manufactured utilizing proprietary construction processes. Starter guidewires are packaged in a spiral hoop fitted with a luer hub for flushing prior to use. "J" tip guidewires are packaged with a guidewire introducer. The guidewire introducer aids in the insertion of the guidewire through the entry needle.

TIP STYLE

INTENDED USE

Starter guidewires are designed for percutaneous peripheral vessel entry using the Seldinger technique to facilitate subsequent introduction(s) of an intravascular device.

CONTRAINDICATIONS

Starter guidewires are not intended for use in the cerebrovasculature or coronary arteries.

WARNINGS

Vessel trauma may result from the improper use of this device. Follow the enclosed directions carefully. When the guidewire is in the body, it should be manipulated only under flouroscopy. Do not attempt to move the wire without observing the resultant tip response. Never advance the guidewire against resistance without first determining the reason for resistance under fluoroscopy. Excessive force against resistance may result in separation of the guidewire tip, damage to the catheter, or vessel damage.

For Starter movable core guidewires, do not advance the movable core while it is positioned in a curved vessel. Advancement and excessive force may cause movable core penetration of the guidewire coil resulting in possible vessel damage.

Avoid withdrawing a guidewire through metal needles; guidewires may shear or coating may scrape off against the needle bevel.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 92238512A), IFU, MB, STARTER, Global, 51105206-01A_pretrans

Black (K) ΔE ≤5.0

PRECAUTIONS

Use the device prior to the "Use By" date noted on the package.

Starter™ guidewires should be used only by physicians thoroughly trained in interventional procedures.

Prior to use, all equipment to be used for the procedure should be carefully examined to verify proper function and integrity.

The Starter guidewires are supplied STERILE and in the unopened and undamaged package. Starter guidewires are sterilized by ethylene oxide gas. Verify that the sterility of the device has not been compromised by assuring the package integrity has been maintained.

Remove the guidewire carefully from the carrier tube to reduce the possibility of damage to the distal tip.

Exercise care in handling of the guidewire during a procedure to reduce the possibility of accidental breakage, bending, kinking, or coil separation. Do not use a guidewire that has been damaged. Do not attempt to straighten a wire that has been kinked or bent. Do not advance a kinked guidewire into a balloon catheter or guidewire to reduce the potential of wire breakage.

ADVERSE EVENTS

Other potential adverse reactions which may result from the improper use of this device include, but are not limited to:

- Embolism (air, plaque, thrombus, device, or other)
- Hematoma at the puncture site
- Hemorrhage
- Infection
- Vascular thrombosis
- Vessel spasm
- Vessel trauma and perforation of the vessel

HOW SUPPLIED

This product is supplied sterile and intended for single use.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.


PREPARATION TECHNIQUE

Every Starter guidewire is provided in a dispenser package. It is recommended the dispenser be filled with heparinized saline solution to flush the guidewire.

Secure a filled syringe to the luer hub on the spiral hoop and flush several times. After flushing, remove the guidewire from the hoop and follow the instructions for use.

OPERATIONAL INSTRUCTIONS

1. PUNCTURE VESSEL.


New Cannula in Vessel

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 92238512A), IFU, MB, STARTER, Global, 51105206-01A_pretrans

Black (K) ΔE ≤5.0


2. Pass the guidewire through the needle and advance 5-10 cm of the guidewire into the punctured vessel. It may be necessary to manipulate the guidewire in order to successfully advance it. Usually a gentle, rotating motion is sufficient to traverse an obstruction. Avoid rough or overly vigorous manipulation of the guidewire to prevent damage to the guide or the vessel.

 If "J" wire is used, prepare "J" guidewire by sliding the guidewire introducer over "J" tip to straighten it. Insert the "J" guidewire into the needle hub and gently advance it.


Guidewire in Vessel

3. REMOVE.


Needle Removed from Guidewire


4. DILATE TISSUE AND VESSEL WITH THE DILATOR USING A SLIGHT ROTARY MOTION.

5. REMOVE DILATOR.


6. INTRODUCE CATHETER BY SLIDING IT OVER THE GUIDEWIRE.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 92238512A), IFU, MB, STARTER, Global, 51105206-01A_pretrans

Black (K) ΔE ≤5.0

7. REMOVE GUIDEWIRE FROM THE CATHETER.

Boston Scientific does not recommend a particular technique for the use of this device. The steps contained in the preceding directions discuss the Seldinger Technique of percutaneous entry and are for information purposes only. Each physician should evaluate their appropriateness according to individual patient conditions and his or her medical training and experience.

REFERENCES

This procedure should be performed only after a thorough review of technical references which emphasize precautions, contraindications, and risks in much greater detail (refer to the references).

1. McAufey BJ: Oesterle S: Simpson JB
 Advances in guidewire technology.
 AMERICAN JOURNAL OF CARDIOLOGY, 53 (12)
 94C-96C, 1984.
2. McDermott EA
 Coronary angioplasty guidewire technology.
 ZEITSCHRIFT FUR KARDIOLOGIE, 76 Suppl 6 29-32, 1987.
3. Moreyre AE: Khew K: Wilson AC: Patel R
 Effect of prolonged blood contact time on deposition
 of cellular and emorphous material on Teflon-coated
 guidewires: a scanning electron microscopy study.
 CATHETERIZATION AND CARDIOVASCULAR DIAGNOSIS,
 38 (4) 355-9, 1995.
4. Nakhjevan FK
 Use of angioplasty guidewire for technically
 difficult angiography.
 CATHETERIZATION AND CARDIOVASCULAR DIAGNOSIS,
 14 (3) 213, 1988.
5. Roid KD: Smith DC
 Guidewire with a J-shaped tip and a tapered movable core.
 RADIOLOGY, 149 (1) 317, 1983.
6. Seldinger, Sven Iva
 Acta Radiology, 39: 368-376, May 1953.
7. Violini R
 Guidewire technique (letter)
 CATHETERIZATION AND CARDIOVASCULAR DIAGNOSIS,
 16 (1) 80, 1989.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 92238512A), IFU, MB, STARTER, Global, 51105206-01A_pretrans

Black (K) ΔE ≤5.0

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 92238512A), IFU, MB, STARTER, Global, 51105206-01A_pretrans

Black (K) ΔE ≤5.0

Catalog Number
Número de catálogo
Numéro de catalogue
Bestell-Nr.
Numero di catalogo
Catalogusnummer
Referência

Consult instructions for use.
Consultar las instrucciones de uso.
Consulter le mode d'emploi.
Gebrauchsanweisung beachten.
Consultare le istruzioni per l'uso.
Raadpleeg instructies voor gebruik.
Consulte as Instruções de Utilização

Contents
Contenido
Contenu
Inhalt
Contenuto
Inhoud
Conteúdo

EU Authorized Representative
Representante autorizado en la UE
Représentant agréé UE
Autorisierter Vertreter in der EU
Rappresentante autorizzato per l'UE
Erkend vertegenwoordiger in EU
Representante Autorizado na U.E.

Legal Manufacturer
Fabricante legal
Fabricant légal
Berechtigter Hersteller
Fabbricante legale
Wettelijke fabrikant
Fabricante Legal

Lot
Lote
Lot
Charge
Lotto
Partij
Lote

Recyclable Package
Envase reciclable
Emballage recyclable
Wiederverwertbare Verpackung
Confezione riciclabile
Recyclebare verpakking
Embalagem Reciclável

Use By
Fecha de caducidad
Date limite d'utilisation
Verwendbar bis
Usare entro
Uiterste gebruiksdatum
Validade

Australian Sponsor Address
Dirección del patrocinador australiano
Adresse du promoteur australien
Adresse des australischen Sponsors
Indirizzo sponsor australiano
Adres Australische sponsor
Endereço do Patrocinador Australiano

Argentina Local Contact
Contacto local en Argentina
Contact local en Argentine
Lokaler Kontakt Argentinien
Contatto locale per l'Argentina
Contactpersoon Argentinië
Contacto local na Argentina

Brazil Local Contact
Contacto local en Brasil
Contact local au Brésil
Lokaler Kontakt Brasilien
Contatto locale per il Brasile
Contactpersoon Brazilië
Contacto local no Brasil

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 92238512A), IFU, MB, STARTER, Global, 51105206-01A_pretrans

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Black (K) ΔE ≤5.0

EC REP **EU Authorized Representative**

Lake Region Medical Limited
Butlersland, New Ross
Co Wexford
IRELAND

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Legal Manufacturer

Manufactured by:
Lake Region Medical
340 Lake Hazeltine Drive
Chaska, MN 55318
USA

Distributed by:
Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
USA
USA Customer Service 888-272-1001

Do not use if package is damaged.

Recyclable Package

CE 0050

© 2020 Boston Scientific Corporation or its affiliates.
All rights reserved.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 92238512A), IFU, MB, STARTER, Global, 51105206-01A_pretrans

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