SKW RETRACTOR SYSTEM

INSTRUCTIONS FOR USE
(ENGLISH)

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R_x Only

DEVICE DESCRIPTION/INTENDED USE

The SKW Retractor System is a single-use, disposable device that is used to retract tissue from the surgical site. It is designed for use with an inflatable penile prosthesis (penoscrotal approach), or an artificial urinary sphincter (transverse scrotal approach).

SKW RETRACTOR SYSTEM COMPONENTS

- 1 Retractor Ring
- 1 Penile Strap
- 2 Baby Deavers
- 1 Catheter Clip
- Set of Hooks:
 - 7 Yellow, 12-mm Blunt Stay Hooks (3 Extra)
 - 1 Green, 5-mm Sharp Stay Hook
 - 2 Spira Three-Finger Hooks

INDICATIONS FOR USE

The SKW Retractor System is indicated for use to provide incision retraction and to aid in visualization during surgical procedures.

CONTRAINDICATIONS


There are no known contraindications associated with the use of this device.

WARNINGS

- Do not use if sterile barrier is damaged.
- Do not use if the device has been opened or damaged.
- For single use only. Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Dispose of in accordance with all applicable Federal, State and local Medical Hazardous waste practices.

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INSTRUCTIONS FOR USE

1. Place the Retractor Ring with the big end up and the little end down. Position the penile strap in the first slot of the retractor ring above the hinges. Only four beads are visible over each side of retractor ring. Penile strap must be placed at the base of the penis.
2. Place a sharp hook in meatus at 12:00. Penis should be stretched and elevated by penile strap.
3. Make a full thickness high scrotal transverse incision. Move upper scrotal incision to penoscrotal junction. Hold scrotal incision on penis while placing hooks.
4. Blunt hooks are placed at 1:00, 3:00, 9:00, and 11:00. Spira three-finger hooks (or blunt hooks) are placed at 5:00 and 7:00.
5. Dissect overlying tissue until tunica and urethra are exposed. Pass Metzenbaum's Scissors outside of tunica to ischial tuberosity.
6. For use with the **artificial urinary sphincter**, place Baby Deaver blades in dissection plane and completely incise septum for full exposure of bulbar urethra.

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EXPLANATION OF SYMBOLS

QTY	QUANTITY
REF	REORDER NUMBER
LOT	BATCH CODE
	USE-BY DATE
	DATE OF MANUFACTURE
STERILE EO	STERILIZED USING ETHYLENE OXIDE
	DO NOT REUSE
	DO NOT RESTERILIZE
	CONSULT INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGE IS DAMAGED
EC REP	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY.
R_x Only	U.S. FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED PRACTITIONER OR LICENSED MEDICAL FACILITY.
	MANUFACTURER

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Trumbull, CT 06611 U.S.A.
Made in the U.S.A.
Boston Scientific Corporation