Boston Scientific

CrossBoss™

Catheter

Directions for Use	2
Instrucciones de uso	4
Mode d'emploi	6
Gebrauchsanweisung	8
Istruzioni per l'uso	11
Gebruiksaanwijzing	13
Instruções de Utilização	15

50949163-01

2019-09

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240AM), DFU, MB, CrossBoss, Global, 50949163-01A

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WARNING

Contents supplied STERILE using a Radiation process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

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DEVICE DESCRIPTION

The CrossBoss Catheter facilitates the placement and support of a guidewire into discrete regions of the coronary vasculature through its central guidewire lumen. The CrossBoss Catheter contains a rounded distal tip (1 mm diameter) mounted to a flexible and torqueable proximal shaft. The CrossBoss Catheter is packaged with a dedicated torque device positioned at its proximal end. The distal end of the catheter is hydrophilic coated. The CrossBoss Catheter is compatible with 6F guide catheters with minimum inner diameter of 0.070 in (1.7 mm), and may be used with guidewires ≤0.014 in (0.36 mm). All components of the CrossBoss Catheter are provided sterile and non-pyrogenic. For dimensional information refer to label.

Contents

Quantity	Material
1	CrossBoss Catheter

INTENDED USE / INDICATIONS FOR USE

The CrossBoss Catheter is intended for use with a guidewire to access discrete regions of the coronary vasculature.

When used as part of the system consisting of the CrossBoss Catheter, Stingray™ LP Catheter, and Stingray Guidewire, the CrossBoss Catheter is indicated for use to facilitate the intraluminal placement of conventional guidewires beyond stenotic coronary lesions (including chronic total occlusions [CTOs]) prior to PTCA or stent intervention.

CONTRAINDICATIONS

- Do not use with guidewire extension systems with a coupling profile larger than 0.014 in (0.36 mm) diameter (i.e. wave pattern coupling mechanism).

WARNINGS

- Only physicians thoroughly trained in interventional procedures should use the CrossBoss Catheter.
- To reduce the potential for vessel damage, the CrossBoss Catheter should only be used in vessels that are ≥1 mm in diameter.
- Always use the included torque device during catheter advancement and manipulation especially during device rotation. Failure to use the included torque device may result in catheter failure and may result in patient injury.

PRECAUTIONS

- Use prior to the Expiration Date.
- The CrossBoss Catheter should only be used in hospitals where emergency coronary bypass surgery can be immediately performed in the event of a potentially injurious or life threatening complication.
- Before insertion of the CrossBoss Catheter, administer appropriate anticoagulant and vasodilator therapy.
- The CrossBoss Catheter should be handled with care. Prior to use and during the procedure, inspect the packaging and catheter for bends, kinks, or other damage. Discontinue use if the catheter becomes damaged.
- Do not expose to organic solvents.
- The CrossBoss Catheter should only be manipulated under fluoroscopic observation.
- Do not use a syringe smaller than 5 cm^3 (5 cc) when flushing the catheter.
- Do not use an inflation or power assist device when flushing the catheter.

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- Acute myocardial infarction
- Artery spasm
- Drug reactions, allergic reaction to contrast media
- Embolism

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240AM), DFU, MB, CrossBoss, Global, 50949163-01A

Black (K) ΔE ≤5.0

- Hemorrhage or hematoma
- Infection
- Neurological deficit
- Stroke
- Vessel trauma requiring surgical repair or intervention

OPERATIONAL INSTRUCTIONS

Supplies

The following supplies are not provided with the device and need to be available and prepped prior to use of the CrossBoss™ Catheter:

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Product Preparation

1. Remove the CrossBoss Catheter from the sterile packaging.
2. Flush the guidewire lumen with saline.
3. Ensure all device surfaces are exposed to heparanized and sterile saline.
4. Gently tighten the black knob of the CrossBoss Catheter torque device onto the proximal shaft at a position near the proximal hub.

Product Operation

1. Access the vascular system and place the appropriate size introducer sheath per manufacturer's instructions.
2. Engage the vessel to be treated with an appropriate guiding catheter per manufacturer's instructions. Connect a hemostatic Y-adapter to the guiding catheter hub.
3. Assure the CrossBoss Catheter has been prepped per instructions.
4. Remove the appropriate guidewire from its package and inspect for damage.
5. Advance the guidewire to the target location using fluoroscopic guidance.
6. Back-load the CrossBoss Catheter onto the proximal portion of the guidewire.
7. Loosen the knob on the hemostatic Y-adapter.
8. Track the CrossBoss Catheter over the wire and into the proximal guiding catheter lumen.
9. Tighten the hemostatic Y-adapter to create a seal around the CrossBoss Catheter. Do not over-tighten.
10. Advance the CrossBoss Catheter over the guidewire using fluoroscopic guidance.
11. Once the CrossBoss Catheter tip exits the guiding catheter, loosen the black knob on the torque device, slide the torque device to the appropriate position along the proximal shaft and firmly tighten the black knob onto the proximal shaft.
12. While holding the torque device, advance while rotating the CrossBoss Catheter in the clockwise or counterclockwise direction. The CrossBoss Catheter may be advanced over the wire or in front of the wire to the target location.
13. In the event the device encounters significant torsional resistance, the torque device is designed to slip while emitting an audible and tactile "click." Upon torque device slip in one direction, reverse the direction of rotation and continue advancement. Continue clockwise and counterclockwise rotations as needed to advance to the target location.
14. In addition, the CrossBoss Catheter can also be withdrawn slightly and re-advanced in the event of torque device slip.
15. Once the CrossBoss Catheter has been successfully delivered to the target location, advance an appropriate guidewire through the central guidewire lumen of the CrossBoss Catheter.
16. Withdraw the CrossBoss Catheter over the wire and discard. Catheter rotation may also be advantageous during withdrawal.
17. Upon completion of the use of the CrossBoss Catheter, discard according to hospital standards and national legal requirements.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose**. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments**.

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240AM), DFU, MB, CrossBoss, Global, 50949163-01A

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Stingray™ LP
Catheter

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Stingray LP Catheter facilitates the placement, support and steering of a guidewire into discrete regions of the coronary vasculature through the central guidewire lumen or through one of two side-ports. These side-ports are on opposite sides of the balloon and are identified by radiopaque markers. The side-ports communicate with the central guidewire lumen and facilitate guidewire steering (at an angle to the central lumen) by allowing the guidewire to exit the Stingray LP Catheter. The Stingray LP Catheter contains a small balloon used for fluoroscopic orientation on the distal tip of a flexible shaft. The distal end of the catheter is hydrophilic coated. The Stingray LP Catheter is compatible with 6F guide catheters with minimum inner diameter of 0.070 in (1.7 mm), and may be used with guidewires ≤0.014 in (0.36 mm). All components of the Stingray LP Catheter are provided sterile and non-pyrogenic. For dimensional information refer to label.

User Information

An interventional cardiologist performing percutaneous coronary interventions.

Contents

Quantity	Material
1	Stingray LP Catheter

INTENDED USE / INDICATIONS FOR USE

The Stingray LP Catheter is indicated for directing, steering, controlling, and supporting a guidewire in order to access discrete regions of the coronary vasculature.

When used as part of the System consisting of the CrossBoss™ Catheter, Stingray LP Catheter and Stingray Guidewire, the Stingray LP Catheter is indicated for use to facilitate the intraluminal placement of conventional guidewires beyond stenotic coronary lesions (including chronic total occlusions [CTOs]) prior to PTCA or stent intervention.

CONTRAINDICATIONS

- Not intended for use in the cerebral vasculature.
- Unprotected left main coronary artery.
- Coronary artery spasm in the absence of a significant stenosis.

WARNINGS

- Only physicians thoroughly trained in interventional procedures should use the Stingray LP Catheter.
- To reduce the potential for vessel damage, the inflated balloon should only be used in vessels that are ≥2.5 mm in diameter as measured proximal and distal to the inflation site.
- Balloon pressure should not exceed the rated burst pressure. Use of a pressure monitoring device is recommended to prevent over pressurization.
- Use only the recommended balloon inflation medium.
- Do not advance or retract the Stingray LP Catheter unless the balloon is fully deflated under vacuum.
- Pay special attention when directing a guidewire out a side port. Never push, withdraw or torque components that meet resistance. Device kinking/breaking or vessel damage may occur. Fluoroscopy should always be used to aid device manipulation.

PRECAUTIONS

- Use prior to the Expiration Date.
- The Stingray LP Catheter should only be used in hospitals where emergency coronary bypass surgery can be immediately performed in the event of a potentially injurious or life threatening complication.
- Before insertion of the catheter, administer appropriate anticoagulant and vasodilator therapy.
- The Stingray LP Catheter should be handled with care. Prior to use and during the procedure, inspect the packaging and catheter for bends, kinks, or other damage. Discontinue use if the catheter becomes damaged.
- Do not expose to organic solvents.
- Always use wire support for advancement, manipulations, and withdrawal of the Stingray LP Catheter.
- The Stingray LP Catheter should only be manipulated under fluoroscopic observation.

Boston Scientific (Master Brand DFUTemplate 8in x 8in Global, 90106041AS), DFU, MB, OUS, Stingray LP, Global, 50949160-01A

Black (K) ΔE ≤5.0

ADVERSE EVENTS

Potential adverse events include, but are not limited to, the following:

- Acute myocardial infarction
- Artery spasm
- Cardiac tamponade
- Drug reactions, allergic reaction to contrast media
- Embolism
- Hemorrhage or hematoma
- Infection
- Recurrence of angina
- Stroke/neurological deficit
- Vessel trauma that may require further intervention or surgical repair

HOW SUPPLIED

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Supplies

The following supplies are not provided and need to be available and prepped prior to use of the Stingray™ LP Catheter:

- Guiding catheter
- Angiographic imaging supplies (i.e. radiopaque contrast, manifold, tubing, etc.)
- Sterile heparinized normal saline
- Guidewire
- Other supplies as needed to complete the established laboratory protocol

Product Preparation

1. Remove the Stingray LP Catheter from the sterile packaging.
2. Remove the balloon protector and stylet.
3. Flush the guidewire lumen with saline.
4. Ensure all device surfaces are exposed to heparinized and sterile saline.
5. Prepare balloon with either method 1 or 2 below:

 Prep method 1:

 a. Prepare a 1:1 mixture of contrast medium and sterile saline.
 b. Aspirate 3 cm^3 to 4 cm^3 (3 cc to 4 cc) of contrast mixture into a 20 cm^3 (20 cc) or larger syringe and purge air from the syringe barrel.
 c. Attach the syringe to a female luer lock on a 3-way stopcock and fill the stopcock lumens with contrast medium mixture.
 d. Attach the 3-way stopcock and syringe to the Stingray LP Catheter inflation port.
 e. Hold the syringe with the nozzle pointing downward and aspirate for 5 seconds.
 f. Close the Stingray LP Catheter inflation port with the 3-way stopcock, remove the syringe, and evacuate air from the barrel.
 g. Reconnect the syringe to the 3-way stopcock and purge the 3-way stopcock.
 h. Open the stopcock, aspirate until bubbles no longer appear, release the plunger, and disconnect the syringe from the 3-way stopcock.

 Prep method 2:

 a. Connect a new/unused stopcock to a new/unused 20 cm^3 (20 cc) syringe and close stopcock side-port.
 b. Connect Stingray LP Catheter inflation port to stopcock.
 c. Hold the syringe with the nozzle pointing downward and aspirate for 15 seconds.
 d. Close stopcock to Stingray LP Catheter.
 e. Remove the syringe and evacuate all air from the barrel.
 f. Reconnect the syringe to the stopcock and close the side-port.
 g. Hold the syringe with the nozzle pointing downward and aspirate for 15 seconds.
 h. Close stopcock to Stingray LP Catheter and remove the syringe.
 i. Fill 3 cm^3 (3 cc) syringe with 3 cm^3 (3 cc) of 100% contrast.
 j. Purge air from 3 cm^3 (3 cc) syringe.
 k. Connect 3 cm^3 (3 cc) syringe to stopcock and purge air from stopcock.
 l. Open stopcock between 3 cm^3 (3 cc) syringe and Stingray LP Catheter to allow contrast to draw into Stingray LP Catheter.
 m. Close stopcock to Stingray LP Catheter and remove 3 cm^3 (3 cc) syringe.
6. Prepare a balloon inflation device according to the manufacturer's instructions.

Product Operation

1. Access the vascular system and place the appropriate size introducer sheath per manufacturer's instructions.
2. Engage the vessel to be treated with an appropriate guiding catheter per manufacturer's instructions. Connect a hemostatic Y-adapter to the guiding catheter hub.
3. Assure the Stingray LP Catheter has been prepped per instructions.
4. Remove the appropriate guidewire from its package and inspect for damage.
5. Advance the guidewire to the target location using fluoroscopic guidance.
6. Back-load the Stingray LP Catheter onto the proximal portion of the guidewire.
7. Loosen the knob on the hemostatic Y-adapter.
8. Track the Stingray LP Catheter over the wire and into the proximal guiding catheter lumen.
9. Tighten the hemostatic Y-adapter to create a seal around the Stingray LP Catheter. Do not over-tighten.
10. Advance the Stingray LP Catheter over the guidewire using fluoroscopic guidance. Position the Stingray LP Catheter ports (using the radiopaque markers) at the target location.
11. Exchange the guidewire for a shaped guidewire (~45° bend at the distal 1.5 mm segment).
12. Connect the balloon inflation device to the 3-way stopcock and purge the stopcock with contrast medium mixture.
13. Inflate the balloon to 3 atm - 4 atm (304 kPa-405 kPa).
14. Adjust the fluoroscopic view until the balloon is shown at its minimum width.

Boston Scientific (Master Brand DFUTemplate 8in x 8in Global, 90106041AS), DFU, MB, OUS, Stingray LP, Global, 50949160-01A

Black (K) ΔE ≤5.0

15. Steer the shaped guidewire through the desired port by slowly rotating the shaped guidewire immediately proximal of the marker band until the shaped guidewire begins to exit the port.
16. Assess proper guidewire position and direction.
17. If the guidewire position or direction is incorrect, target the other side port by steering the shaped guidewire through the desired port by slowly rotating the shaped guidewire immediately proximal of the marker band until the shaped guidewire begins to exit the port.
18. Advance the shaped guidewire through the port as needed.
19. After completing the procedure, deflate the balloon, withdraw the Stingray™ LP Catheter and discard according to hospital standards and national legal requirements.

GUIDEWIRE TRAPPING TECHNIQUE

Simultaneous use of two balloon catheters in a guide catheter for the use of guidewire securement (Guidewire Trapping Technique): Bench testing has shown that one Stingray LP Catheter and one 3.0 mm x 20 mm (or smaller) MR balloon catheter with a max shaft OD of 0.036 in (0.91 mm) can be inserted simultaneously into a 7F [minimum 0.080 in (2.0 mm) ID] guide catheter. Care should be used when attempting to use two balloon catheters simultaneously in a guide catheter; this technique was not clinically evaluated for safety and effectiveness in a clinical trial. Balloon catheters with a diameter greater than those mentioned have not been tested for simultaneous use in a single guide catheter.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Boston Scientific (Master Brand DFUTemplate 8in x 8in Global, 90106041AS), DFU, MB, OUS, Stingray LP, Global, 50949160-01A

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All rights reserved.
© 2019 Boston Scientific Corporation or its affiliates. All rights reserved.