Pleurofix® No. 1/No. 2

GB	Set for puncture of the pleural cavity	3
DE	Besteck zur Pleurapunktion	3
BG	Комплект за пункция на плевралната кухина	4
CZ	Souprava k punkci pleurálního prostoru	5
DK	Sæt til punktur af lungehulen	5
EE	Komplekt pleuraõõne punktsiooniks	6
ES	Set de punción para la cavidad pleural	7
FI	Pleuraontelon punktiosetti	7
FR	Nécessaire pour ponction de la cavité pleurale	8
GR	Σετ παρακέντησης υπεζωκοτικής κοιλότητας	9
HR	Set za punkciju pleuralne šupljine	10
HU	Mellkaspunkciós készlet	10
IT	Set completo per drenaggio pleurico	11
KZ	Өкпеқап қуысын шанышқылауға арналған жинақ	12
LT	Pleuros ertmės punkcijos rinkinys	12
LV	Pleiras dobuma punkcijas komplekts	13
NL	Set voor punctie van de pleuraholte	14
NO	Punksjonssett for plevral kavitet	14
PL	Zestaw do punkcji jamy opłucnej	15
PT	Conjunto para punção da cavidade pleural	16
RO	Set de puncţie a cavităţii pleurale	17
RS	Set za punkciju pleuralne šupljine	17
RU	Набор для пункции плевральной полости	18
SE	Pleuratappningsset	19
SI	Komplet za punkcijo plevralnega prostora	19
SK	Súprava na punciu pleurálneho priestoru	20
TR	Plevral boşluk ponksiyon seti	21

B BRAUN

12414204_Pleurofix_A5_IFU_1117.indd 1 LLDorder 3752 - G 171305 01.12.17 14:13

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GB Instructions for use

Contents
Pleurofix® No.1, consisting of:
Thin-wall puncture needle
Omnifix®-Locksyringe 1,8 x 80 mm, 50 ml (usable upto 60 ml)
Discofix®-3 three-way stopcock, yellow
Secretion collecting bag, 2.0

Pleurofix®No.2 consisting of:
Thin-wall puncture needle
Omnifix®-Lock syringe 1,8 x 80 mm, 50 ml (usable upto 60 ml)
Anti-reflux valve
Secretion collecting bag, 2.0 l

Materials used
ABS, IR, NR, PA/PE, PC, PE, PP, PVC, steel

Indications
- Removal of secretion from the pleural cavity, pneumothorax
- Irrigation of the pleural cavity, instillation of medicaments

Contraindications
- Coagulation disturbances (absolute or relative contraindications)

Risks
Depending on how the set is used the known risks apply, especially
- Haemorrhage
- Injury to the visceral pleura and lung
- Damage to the lung by excessive suction (use only instruments with reliable suction limiter)
- Circulatory distress if excessive amounts of pleural fluid are aspirated
- Infection
- Pneumothorax if improperly used

Duration of use
Only for single removal of secretion. Do not leave in situ because of risk of injury.

Warning
Animal experiments have shown DEHP to be potentially toxic to reproduction. Given the current state of scientific knowledge, a risk (in case of long-term use) especially for premature male babies cannot be completely excluded. As a precaution, the application of medical products containing DEHP should be restricted to short-term use for pregnant women, breastfeeding mothers, infants, and children.
Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient. Do not use if packaging is damaged. Do not resterilise.

Application
Caution:
Puncture must be effected under sterile conditions.
- Connect cannula and syringe to the three-way stopcock or the extension tubing
- Close the cannula by turning the three-way stopcock accordingly
- Puncture, advance the cannula into the cavity of the pleura
- Aspirate secretion into the syringe and conduct secretion into the secretion collecting bag (turn three-way stopcock*)
- Remove subsequent secretion as well
- Withdraw cannula out of the pleural cavity

*) using Pleurofix No. 2, the anti-reflux valve takes on the function of the three-way stopcock. Therefore it doesn't have to be turned.

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B. Braun Melsungen AG