PTFE Peelable Introducer Kit

Instructions For Use (en)	1
Mode d'emploi (fr)	2
Gebrauchsanleitung (de)	3
Gebruiksaanwijzing (nl)	4
Instrucciones de uso (es)	5
Istruzioni per l'uso (it)	6
Instruções de Utilização (pt EU)	7
Brugsvejledning (da)	8
Bruksanvisning (sv)	9
Οδηγίες χρήσης (el)	10
Návod k použití (cs)	11
Instruções de uso (pt BR)	12
使用方法 (ja)	13
Kullanım Talimatları (tr)	14
Инструкции по применению (ru)	15
使用说明 (zh CN)	16
사용 지침 (ko)	17
วิธีใช้ (th)	18
使用說明 (zh TW)	19
Bruksinstruksjoner (no)	20
Használati utasítás (hu)	21
Instrukcja użytkowania (pl)	22
Návod na použitie (sk)	23

Greatbatch™ Medical

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PTFE Peelable Introducer Kit (en)

Warnings and Precautions

- Product is sterile in unopened, undamaged package.
- For single use only. DO NOT REUSE. DO NOT REPROCESS.
- Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.
- Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure.
- Federal (USA) law restricts this device to sale by or on the order of a physician.
- The introducer shall be used only with components contained in this kit and is intended for use with leads or catheters with the same or smaller diameter than the inner diameter stated on the label.
- Do not attempt to re-insert a partially or completely withdrawn needle.
- Discard of Medical Device: Safely discard the medical device in accordance with local laws and regulations.

Indications for Use

Peelable introducer is intended for use in the percutaneous insertion of pacing leads or catheters in the venous system.

Contraindications

Use of this introducer is contraindicated if the patient has a known or suspected obstruction in the subclavian vein. There is increased risk of pneumothorax for a patient who has severe chronic lung disease. Poor healing may result in a patient who has had irradiation to the anterior chest.


Potential Adverse Events

The potential complications related to the use of the introducer include, but are not limited to the following:


- Subclavian artery puncture
- Pneumothorax
- Hemothorax
- Air Embolism
- Hematoma formation
- Brachial plexus injury
- Subclavian vein thrombosis
- Wound infection
- Mediastinal widening

Instructions for Use

1. Gain percutaneous access with a guidewire introducer needle.
2. Advance a guidewire through the needle.
3. Withdraw the needle while leaving the guidewire in place.
4. Advance the peelable sheath and dilator assembly over the guidewire.

5. Remove the vessel dilator and guide wire, leaving the sheath as a conduit into the subclavian vein and superior vena cava. Immediately place a finger over the sheath to prevent excessive bleeding or possible air aspiration.
6. Advance a catheter or pacemaker lead through the sheath.
7. Once the catheter or pacemaker lead is properly positioned, peel the sheath off the catheter or pacemaker lead by pulling the handles so the sheath separates longitudinally.

Storage Conditions

Store device in ambient conditions.

Limited Warranty and Disclaimer

Greatbatch Medical hereby warrants that if the Greatbatch Medical product fails to perform within Greatbatch Medical's normal tolerance due to defects in materials or workmanship, Greatbatch Medical in its sole discretion will either provide, at no charge, a replacement to the failed product or a refund of the purchase price. This Limited Warranty applies only if:
1. Greatbatch Medical packaged and labeled the product;
2. The product was used before the applicable "Use By" or expiration date;
3. The failed product is returned to Greatbatch Medical;
4. The product was used solely for the purpose for which it was intended by Greatbatch Medical; and
5. The product has not been mishandled, reprocessed, altered or damaged.

Greatbatch Medical makes no representation of warranty (whether written, oral, express or implied, herein or otherwise) that a Greatbatch Medical product will not fail. Greatbatch Medical is not responsible for any medical complication, including death, resulting from the use of its products.

GREATBATCH MEDICAL SHALL NOT BE RESPONSIBLE FOR ANY DIRECT, INDIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON OR RELATING TO ANY DEFECT, FAILURE, OR MALFUNCTION OF ITS PRODUCTS, WHETHER BASED ON WARRANTY, CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE. FURTHERMORE, EXCEPT AS EXPRESSLY PROVIDED BY THE LIMITED WARRANTY ABOVE, GREATBATCH MEDICAL MAKES NO WARRANTY (WHETHER WRITTEN, ORAL, EXPRESS OR IMPLIED, HEREIN OR OTHERWISE) INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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© 2016 Greatbatch Medical. All rights reserved.
Medical Product Service GmbH