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
EN

Instructions for Use

Guide Wires for Central Circulatory Applications

Intended Use

Guide Wires are intended to be placed in the human body invasively or surgically-invasively through a natural or artificially produced body orifice or mini-mally-invasively through a port in the blood vessels of the central circulatory system, to be guided by controllable distal and proximal ends to a certain area of the human body to be used as a guide for catheters or other medical devices which are intended to be used in this area of the body.

Indications

Angioplasty (PTA, PTCA), Percutaneous Dilatationstracheotomy (PDT), Chronic Total Occlusion (CTO), Stent Implantation (peripheral, coronary, AAA), Carotid Entarterectomy (CEA), Carotid Artery Stenting (CAS); Transcatheter Aortic Valve Implantation (TAVI)

Contraindications

As the guide wire is only an accessory for the above mentioned indications, the contraindications are limited to the contraindications for these applications.

- Applications in the central nervous system
- Contraindications for guide wires > 0.014" (0.35 mm) can be Chronic Total Occlusions (CTOs).

Application

General

The application of the product is exclusively restricted to medical specialists and experienced users. The interventional combination device has to be prepared according to the instructions of the manufacturer of the device. The guide wire lumen of the combination device has to be flushed before insertion of the guide wire.

Guide Wires with Hydrophilic Coating

Prior to the withdrawal of the full-length hydrophilic guide wire from the dispenser, physiological saline solution has to be injected into the hub end of the dispenser in order to wet the complete hydrophilic surface of the guide wire.

After the injection of the saline solution, the guide wire has to be carefully withdrawn from the dispenser. If the guide wire cannot be easily withdrawn from the dispenser further saline solution has to be injected into the dispenser.

After withdrawal, the guide wire with hydrophilic coating must not be introduced into the dispenser again.

Warning Notices

If not used correctly guide wires can lead to arterial perforation. The damages caused by arterial perforation are very often extremely serious and can lead to death of the patient. Arterial perforation primarily occurs due to kinking of the guide wire. During the use of guide wires an insertion technique must be used that avoids kinking of the guide wire. An insertion guide has to be used. Furthermore, there have been cases of overinsertion of guide wires; this is generally considered a serious risk. Overinsertion of guide wires can lead to injury and/or perforation of the vessel being treated.

Guide wires contain ferromagnetic materials and must therefore not be used under MR conditions, as heating and movement of the guide wire can occur in magnetic fields. This can lead to severe complications.

ATTENTION: With the exception of the correspondingly marked products, guide wires must not be used under MR conditions!!

The guide wires marked as „MR conditional" (see symbol description) were tested for safety in applications under MR conditions and were classified as MR conditional. They can be used under conditions specified below under MRT conditions:.

- Static magnetic fields of 1.5 T and 3 T
- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force produced of 231 T2/m
- Theoretically expected maximum whole body average (WBA) specific absorption rate (SAR) of 2 W/kg (normal operation mode) or 4 W/kg (first level operation mode).

Under the conditions described above, after 15 minutes of continuous scanning, a maximum temperature rise is to be expected of less than

- 1,6 °C (2 W/kg, 1.5 Tesla) RF-dependent temperature increase with a background temperature rise of ≈ 1.4 °C (2 W/kg, 1.5 Tesla) and 3.2 °C (4 W/kg, 1.5 Tesla) RF-dependent temperature increase with a background temperature rise of ≈ 2.8 °C (4 W/kg, 1.5 Tesla)
- 1,2 °C (2 W/kg, 3 Tesla) RF-dependent temperature increase with a background temperature rise of ≈ 1,2 °C (2 W/kg, 3 Tesla) and 2.5 °C (4 W/kg, 3 Tesla) RF-dependent temperature rise with a background temperature rise of ≈ 2.5 °C (4 W/kg, 3 Tesla)

In non-clinical tests, the image artifacts generated by the guide wire expanded by about 9.62 mm with a gradient echo pulse sequence and a 3 Tesla system.

The guide wires, which are marked as MR conditional are non-conductive and may therefore also be used together with HF devices in the field of gastroenterology.

According to clinical data found in literature, a passing through strictures of greater than 1 cm in length may lead to damage or breakage of soft guide wires. This should be taken into consideration during surgical planning. In rare cases, repeated leading back and forth of the guide wire may lead to knot building at the distal tip of the wire. Therefore, the user should have appropriate knowledge in such techniques in order to manage serious complications related to such events.

If coated guide wires are used with other medical devices (e.g. catheter), special attention has to be paid to an appropriate combination to avoid damages to the integrity of the coating. Certain cases of entrapped or breakdown guide wires have been reported as well as possible techniques to retrieve the entrapped or broken parts. The user should have knowledge in such techniques in order to prevent serious consequences related to such events.

Like almost every product that is introduced into blood vessels, guide wires may lead to thrombosis. The user is responsible for taking adequate measures to prevent thrombosis during the application of the guide wire.

The guide wire should be manipulated slowly and carefully. The behavior and advancement of the guide wire tip should be controlled under fluoroscopy. Manipulation of the guide wire without control by fluoroscopy can result in injury to the patient.

Likewise, the guide wire must be held in position when changing or withdrawing the catheter in order to prevent injury to the patient.

A scale on the guide wire or scaled radiopaque markers are for orientation only.

The measuring results should be confirmed by other methods.

If any resistance is felt during use, the manipulation of the guide wire and/or the combination device used such as catheter must be stopped immediately and the cause for the resistance must be determined by use of various visual methods. Otherwise, overexpansion, kinking, breaks or partial detachment of the outer polymer jacket or the coating material of the guide wire can occur. Remove, if required, the guide wire and the combination devices used as a complete unit to avoid complications.

Never insert, advance and/or withdraw the guide wire through a metal cannula or needle or other sharp-edged devices. This may damage the guide wire and may lead to destruction and/or detachment, in particular of the outer polymer jacket or the coating, and therefore may require subsequent retrieval of the fragments.

Use extreme caution when using a device that emits energy (laser, pressure, ultrasound,...), and retract the guide wire into a position where it will not be damaged. Direct contact may cause damage to the wire and/or sever the wire.

Preventive Measures

- Guide wires are fragile devices which should be used very carefully. Before and, if possible, during its use the guide wire should be checked visually for damages, misalignment, kinks, abrasions of the outer polymer jacket or the coating and for other deformations. Damaged, misaligned, kinked and other deformed guide wires must not be used any longer and must be discarded due to potential risk of arterial perforation.
- The combination ability between the guide wire and interventional medical device must be verified before use.
- Avoid withdrawing the guide wire through a metal cannula. The sharp edges of the cannula can damage or shear the PTFE or polymer coating. The cannula should be replaced as soon as possible after the insertion of the guide wire into the vessel by a catheter, introducer sheath or vessel dilator.
- Due to variations in stiffness, inner diameter or shape of certain catheter tips, abrasions of the outer polymer jacket or the coating may occur when manipulating. If any resistance is felt during introduction of the catheter, such catheters should not be used.
- Guide wires with partial or full-length hydrophilic coating must not come into contact with alcohol, antiseptic substances or substances containing solvents.

Reprocessing and Sterilisation

Guide wires are delivered in sterile condition. Due to the design of the devices, after one use, devices **CANNOT BE EFFECTIVELY** cleaned, disinfected and sterilised after successful application, as body fluids are migrated into the inner part of the wire and cannot be removed afterwards. Used once, the guide wire must be discarded.

The devices only remain sterile if packaging has not been damaged or opened and the use-by date has not been exceeded. In case of damaged sterile packaging or exceeded use-by date, the devices must be discarded.

Disposal and Reprocessing

After use, this product may present a biological hazard. Handling and disposal must be carried out according to recognised medical procedures and in accordance with applicable legal requirements.

All EPflex guide wires are disposable products and must not be reprocessed and resterilised, as damage to the materials and the coatings cannot be ruled out.

Annotations

Symbols used on the packaging and labeling of these devices have the following meanings:

Date of manufacture	Do not use if package is damaged
Manufacturer	Do not reuse
Use by[1]	Caution, consult accompanying documents
LOT Batch code	Consult instructions for use
STERILE EO Sterilized by ethylene oxide	CE 0297 CE-sign and ID-number of the Notified Body DQS Medizintechnik GmbH
Keep dry	MR MR unsafe Do not use this device in MR fields
Keep away from sunlight	MR MR Conditional Used in combination with 1.5 Tesla and 3 Tesla devices under the conditions described above

[1] The symbol is accompanied by a date to indivate that the device must not be used after the end of the year, month, or day shown

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