Johnson&Johnson VISION

IFU in English and Russian for the IOL Model ZKB00

TECNIS® Multifocal 1-piece Intraocular Lens (IOL)

Z311276 Rev. 02
Revision Date: 06/2019
Sheet 1 of 2
EN

Rx Only DESCRIPTION:

The TECNIS® multifocal foldable acrylic 1-piece lens, Model ZKB00, is an ultraviolet light-absorbing posterior chamber intraocular lens (IOL). It is designed to be positioned in the lens capsule where the lens should replace the optical function of the natural crystalline lens. The TECNIS® multifocal foldable acrylic 1-piece lens incorporates a proprietary wavefront-designed aspheric optic with a squared posterior edge designed to provide a 360-degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects.

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1. Some visual effects associated with multifocal IOLs may be expected because of the superposition of focused and unfocused images. These may include a perception of halos or glare around lights under nighttime conditions. It is expected that, in a small percentage of patients, the observation of such phenomena will be annoying and may be perceived as a hindrance, particularly in low illumination conditions. On rare occasions these visual effects may be significant enough that the patient will request removal of the multifocal IOL.
2. Under low-contrast conditions, contrast sensitivity is reduced with a multifocal lens compared to a monofocal lens. Therefore, subjects with multifocal lenses should exercise caution when driving at night or in poor visibility conditions.
3. Children under the age of 2 years are not suitable candidates for intraocular lenses.
4. Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the lens may exacerbate an existing condition, may interfere with diagnosis or treatment of a condition or may pose an unreasonable risk to the patient's eyesight:
 a. Patients with recurrent severe anterior or posterior segment inflammation or uveitis.
 b. Patients in whom the intraocular lens may interfere with the ability to observe, diagnose or treat posterior segment diseases.
 c. Surgical difficulties at the time of cataract extraction and/or intraocular lens implantation that might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
 d. A compromised eye due to previous trauma or developmental defects in which appropriate support of the IOL is not possible.
 e. Circumstance that would result in damage to the endothelium during implantation.
 f. Suspected microbial infection.
 g. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support for the IOL.
 h. Congenital bilateral cataracts.
 i. Recurrent severe anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye.
 j. Previous history of, or a predisposition to, retinal detachment.
 k. Patients with only one eye with potentially good vision.
 l. Medically uncontrollable glaucoma.
 m. Corneal endothelial dystrophy.
 n. Proliferative diabetic retinopathy.
5. The TECNIS® Multifocal 1-piece IOL should be placed entirely in the capsular bag. Do not place the lens in the ciliary sulcus.
6. The splitting of the light into more than one focus may affect image quality and lead to some reduction of contrast sensitivity.
7. Well-informed patients with well-defined visual needs and preferences should be selected for TECNIS® Multifocal 1-piece lens implantation. The patients should be informed about the possibility that a decrease in contrast sensitivity and an increase of visual disturbances may affect their ability to drive a car under certain environmental conditions, such as driving at night or in poor visibility conditions.
8. Patients with a predicted postoperative astigmatism greater than 1.0 diopter may not be suitable candidates for multifocal IOL implantation since they may not fully benefit from a multifocal IOL in terms of potential spectacle independence.
9. Care should be taken to achieve IOL centration, as lens decentration may result in patients experiencing visual disturbances, particularly in patients with large pupils under mesopic conditions.
10. AMO IOLs are single-use medical devices and are labeled with instructions for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The reuse/resterilization/ reprocessing of AMO single-use medical devices may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/ or illness due to product contamination, transmission of infection, and lack of product sterility.

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Additional complications that may occur following implantation of a multifocal intraocular lens include: objectionable visual quality under certain lighting conditions (e.g., chromatic aberrations, halos, night glare, decreased contrast),

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hazy/filmy images, “ghost” imaging of near and/or distant objects, diplopia, and poor stereopsis. Some of these optical effects may be mitigated upon patients’ adaptation to the multifocal optic. However, some of the above listed complications may require secondary surgical intervention for intraocular lens exchange or explantation.

DETAILED DEVICE DESCRIPTION:

The TECNIS® Multifocal 1-piece lens, Model ZKB00, is a one-piece, foldable, posterior chamber lens made of soft acrylic material. The lens has a diffractive multifocal surface on the posterior side and a modified prolate (aspheric) surface on the anterior side of the lens. The optic is 6.0 mm in diameter and the lens has an overall diameter of 13.0 mm. The add power is +2.75 diopters, corresponding to +2.01 diopters in the spectacle plane.

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LENS POWER CALCULATIONS:

Accurate keratometry and biometry are essential to successful visual outcomes. Preoperative calculation of the required dioptric lens power for these posterior chamber intraocular lenses should be determined by the surgeon’s experience, preference, and intended lens placement. **Emmetropia should be targeted.** The appropriate A-constants listed on the outer label are presented as a guideline and is a starting point for implant power calculations. The physician should determine preoperatively the dioptric power of the lens to be implanted.

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- Holladay JT, Musgrove KH, Prager TC, Lewis JW, Chandler TY and Ruiz RS. A three-part system for refining intraocular lens power calculations. J Cataract Refract Surg. 14:17-24, 1988.
- Retzlaff JA, Sanders DR and Kraff MC. Development of the SRK/T intraocular lens implant power calculation formula. J Cataract Refract Surg. 16:333- 340, 1990; ERRATA, 16:528, 1990.
- Olsen T, Olesen H, Thim K and Corydon L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. J Cataract Refract Surg. 18:280-285, 1992.
- Hoffer KJ. The Hoffer Q formula: A comparison of theoretic and regression formulas. J Cataract Refract Surg. 19:700-712, 1993; ERRATA, 20:677, 1994.
- Holladay JT. Standardizing constants for ultrasonic biometry, keratometry and intraocular lens power calculations. J Cataract Refract Surg. 23:1356- 1370, 1997.

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• Norrby NES. Unfortunate discrepancies. Letter to the editor and reply by Holladay JT. J Cataract Refract Surg. 24:433-434, 1998.
• Norrby S, Lydahl E, Koranyi G, Taube M. Reduction of trend errors in power calculation by linear transformation of measured axial lengths. J Cataract Refract Surg. 2003; 29:100-105.
• http://ocusoft.de/ulib/c1.htm is in particular useful for Zeiss IOLMaster users.

PATIENT CARD:
An implant identification card, to be supplied to the patient, is included in the package. The patient should be instructed to keep the card as a permanent record of his/her implant and to show the card to any eye care practitioner he/ she may see in the future.

REPORTING:
Adverse events and/or potentially sight-threatening complications that may reasonably be regarded as lens-related and that were not previously expected in nature, severity or rate of occurrence must be reported to AMO. This information is being requested from all surgeons in order to document potential long-term effects of intraocular lens implantation, especially in younger patients.
Physicians are required to report these events in order to aid in identifying emerging or potential problems with posterior chamber lenses. These problems may be related to a specific lot of lenses or may be indicative of long-term problems associated with these lenses or with intraocular lenses in general.

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Note: *The cut-off wavelengths and the spectral transmittance curves represent the range of the transmittance of IOLs (5-34 diopter) made with this material. Spectral transmission measurements were taken in water at room temperature.*

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1. Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the lens may exacerbate an existing condition, may interfere with diagnosis or treatment of a condition or may pose an unreasonable risk to the patient's eyesight. These conditions are not specific to the design of the lens and are attributed to cataract surgery and IOL implantation in general:
 a) Patients with recurrent severe anterior or posterior segment inflammation or uveitis of unknown etiology, or any disease producing an inflammatory reaction in the eye.
 b) Patients in whom the intraocular lens may affect the ability to observe, diagnose or treat posterior segment diseases.
 c) Surgical difficulties at the time of cataract extraction, which may increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure or significant vitreous prolapse or loss).
 d) A compromised eye due to previous trauma or developmental defects in which appropriate support of the IOL is not possible.
 e) Circumstances that would result in damage to the endothelium during implantation.
 f) Suspected microbial infection.
 g) Patients in whom neither the posterior capsule nor the zonules are intact enough to provide support for the IOL.
 h) Congenital bilateral cataracts.
 i) Previous history of, or a predisposition to, retinal detachment.
 j) Patients with only one good eye with potentially good vision.
 k) Medically uncontrollable glaucoma.
 l) Corneal endothelial dystrophy.
 m) Proliferative diabetic retinopathy.
 n) Children under the age of 2 years are not suitable candidates for intraocular lenses.
2. The TECNIS Eyhance IOL, model ICB00, should be placed entirely in the capsular bag and should not be placed in the ciliary sulcus.
3. Johnson & Johnson Surgical Vision, Inc. IOLs are labeled with instructions for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The resterilization/reprocessing of the IOLs may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

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EN

Device Description

The **TECNIS Eyhance IOL (intraocular lens), model ICB00**, is an ultraviolet-light absorbing posterior chamber IOL which provides distance vision and a dysphotopsia profile comparable to a standard aspheric monofocal IOL. The Model ICB00 lens has a modified aspheric anterior surface that is designed to extend depth of focus, which improves intermediate vision as compared to a standard aspheric monofocal IOL. Accommodation will not be restored. The TECNIS Eyhance IOL, model ICB00, is designed to be positioned in the lens capsule to replace the optical function of the natural crystalline lens. The IOL contains a squared posterior optic edge that provides a 360-degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects.

Indications For Use

The TECNIS Eyhance IOL, model ICB00, is indicated for the visual correction of aphakia in adult patients in whom a cataractous lens has been removed by extracapsular cataract extraction. The lens extends the depth of focus, which improves vision for intermediate tasks, and provide similar distance vision as compared to a standard aspheric monofocal IOL. The lens is indicated for placement in the capsular bag only.

Precautions

1. Autorefractors may not provide optimal postoperative refraction of patients with TECNIS Eyhance IOL, model ICB00, lenses. Manual refraction with maximum plus technique is strongly recommended.
2. Recent contact lens usage may affect the patient's refraction; therefore, for patients that wear contact lenses, surgeons should establish corneal stability without contact lenses prior to determining IOL power.
3. Do not resterilize the lens. Most sterilizers are not equipped to sterilize the soft acrylic material without producing undesirable side effects.
4. Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.
5. Do not store the lens in direct sunlight or at a temperature greater than 113°F (45°C). Do not autoclave the intraocular lens.
6. Do not re-use the intraocular lens.
7. Prior to implanting, examine the lens package for proper lens model, dioptric power, and expiration date.
8. The lens is designed for optimum visual performance when emmetropia is targeted.
9. Please refer to the specific instructions for use provided with the insertion instrument or system for the amount of time the IOL can remain folded before the IOL must be discarded.
10. When the insertion system is used improperly, the haptics of the TECNIS Eyhance IOL, model ICB00, may become broken. Please refer to the specific instructions for use provided with the insertion instrument or system.

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1. Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the lens may exacerbate an existing condition, may interfere with diagnosis or treatment of a condition or may pose an unreasonable risk to the patient's eyesight. These conditions are not specific to the design of the lens and are attributed to cataract surgery and IOL implantation in general:
 a) Patients with recurrent severe anterior or posterior segment inflammation or uveitis of unknown etiology, or any disease producing an inflammatory reaction in the eye.
 b) Patients in whom the intraocular lens may affect the ability to observe, diagnose or treat posterior segment diseases.
 c) Surgical difficulties at the time of cataract extraction, which may increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure or significant vitreous prolapse or loss).
 d) A compromised eye due to previous trauma or developmental defects in which appropriate support of the IOL is not possible.
 e) Circumstances that would result in damage to the endothelium during implantation.
 f) Suspected microbial infection.
 g) Patients in whom neither the posterior capsule nor the zonules are intact enough to provide support for the IOL.
 h) Congenital bilateral cataracts.
 i) Previous tiistory of, or a predisposition to, retinal detachment.
 j) Patients with only one good eye with potentially good vision.
 k) Medically uncontrollable glaucoma.
 l) Corneal endothelial dystrophy.
 m) Proliferative diabetic retinopathy.
 n) Children under the age of 2 years are not suitable candidates for intraocular lenses.
2. The TECNIS Eyhance IOL, model ICB00, should be placed entirely in the capsular bag and should not be placed in the ciliary sulcus.
3. Johnson & Johnson Surgical Vision, Inc. IOLs are labeled with instructions for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The resterilization/reprocessing of the IOLs may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

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1. Endophthalmitis/intraocular infection
2. Hypopyon
3. Hyphema
4. IOL dislocation
5. Cystoid macular edema
6. Pupillary block
7. Retinal detachment/tear
8. Persistent corneal stromal edema
9. Persistent iritis
10. Persistent raised IOP (intraocular pressure) requiring treatment
11. Acute corneal decompensation
12. Secondary intraocular surgical intervention (including implant repositioning, removal, AC tap performed later then one week after cataract surgery, or other surgical procedure)
13. Any other adverse event that leads to permanent visual impairment or requires surgical or medical intervention to prevent permanent visual impairment

Detailed Device Description:

The TECNIS Eyhance IOL, model ICB00, is a one-piece, foldable, posterior chamber lens with an overall diameter of 13.0 mm and an optic diameter of 6.0 mm. It incorporates a spherical posterior surface and a modified aspheric anterior surface that is designed to provide distance vision and a dysphotopsia profile comparable to a standard aspheric monofocal IOL. Additionally, the lens extends the depth of focus, which provides improved intermediate vision compared to a standard aspheric monofocal IOL. The lens also compensates for corneal spherical aberration, similar to the TECNIS® monofocal IOL, model ZCB00.

Lens Optic

1. Material: Optically clear, soft foldable hydrophobic acrylic with a covalently bound UV absorber. Full transmission of the blue wavelength light for optimal scotopic sensitivity.
2. Power: +5.0 to +34.0 diopter powers in 0.5 diopter increments.
3. Optic Center Thickness: 0.72 mm (+20.0D).
4. Optic Edge Design: PROTEC 360 square posterior edge.
5. Index of Refraction: 1.47 at 35°C.
6. Light Transmittance: UV cut-off at 10% T for a +5.0 diopter lens (thinnest) and a +34.0 diopter lens (thickest) are shown in **Figure 1.**

Haptics

1. Material: Soft foldable hydrophobic acrylic with a covalently bound UV absorber.
2. One-piece lens.
3. Configuration: TRI-FIX design Modified C, integral with optic.
4. Haptic thickness: 0.46 mm.

Directions for Use

1. Prior to implanting, examine the lens package for IOL type, power, proper configuration and expiration date.
2. Open the peel pouches and remove the lens in a sterile environment. Verify the dioptric power of the lens.
3. Examine the lens thoroughly to ensure particles have not become attached to it, and examine the lens' optical surfaces for other defects.
4. The lens may be soaked or rinsed in sterile balanced salt solution until ready for implantation.
5. Handle the lens by the haptic portion. Do not grasp the optical area with forceps.
6. Transfer the lens, using a sterile technique, to an appropriate loading device.
7. The surgeon should target emmetropia, as this lens is designed for optimal visual performance when emmetropia is achieved.
8. Johnson & Johnson Surgical Vision, Inc. recommends using The UNFOLDER Platinum 1 Series implantation system or an equivalent insertion instrument or system to insert the TECNIS Eyhance IOL, model ICB00, lenses. Only insertion instruments that have been validated and approved for use with this lens should be used. Please refer to the directions for use with the insertion instrument or system for additional information on how to use the insertion instrument.

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Lens Power Calculations:

Accurate keratometry and biometry are essential to successful visual outcomes. Preoperative calculation of the required diopter power for this posterior chamber intraocular lens should be determined by the surgeon's experience, preference, and intended lens placement. The lens is designed for optimum visual performance when emmetropia is targeted. The A-constants listed on the outer label are presented as a guideline and serve as a starting point for implant power calculations. The physician should determine preoperatively the diopter power of the lens to be implanted.

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How Supplied

The TECNIS Eyhance IOL, model ICB00, is supplied sterile in a lens case within a double aseptic transfer peel pouch. The double aseptic transfer peel pouch is sterilized with ethylene oxide and should be opened only under sterile conditions. The pouch and product labels are enclosed in a shell pack. The external surfaces of the outer pouch are not sterile.

Expiration Date

The Use By date on the lens package is the sterility expiration date. This lens should not be implanted after the indicated sterility expiration date.

Return/Exchange Policy

Contact your local Johnson & Johnson Surgical Vision, Inc. representative for the return/exchange policy. Return the lens with proper identification and the reason for the return. Label the return as a biohazard. Do not attempt to resterilize the lens.

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Symbol/Explanation

SYMBOL	EXPLANATION
STERILE EO	Sterilized using ethylene oxide
	Do not re-use
	Use-by date (YYYY-MM-DD: Year-Month-Day)
	Consult instructions for use
	Manufacturer
EC REP	Authorized representative in the European Community
	Do not resterilize
45°C (113°F)	Upper limit of temperature
	Keep away from sunlight
	Date of manufacture (YYYY-MM-DD: Year-Month-Day)
	Do not use if package is damaged
REF	Catalogue Number
SN	Serial number

TECNIS, TECNIS Eyhance, PROTEC, TRI-FIX, and UNFOLDER are trademarks of Johnson & Johnson Surgical Vision, Inc.

Product of USA
AMO Puerto Rico Manufacturing, Inc.
Road 402, Km. 4.2
Añasco, Puerto Rico 00610

Johnson & Johnson Surgical Vision, Inc.,
1700 East St. Andrew Pl.,
Santa Ana, CA 92705 USA

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
© Johnson & Johnson Surgical Vision, Inc. 2020

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Note: ***The cut-off wavelengths and the spectral transmittance curves represent the range of the transmittance values of IOLs (5-34 diopter) made with this material.***

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TECNIS Synergy™ OptiBlue® IOL, Model ZFR00V

DEVICE DESCRIPTION:

The TECNIS Synergy™ OptiBlue® foldable acrylic 1-piece lens, Model ZFR00V, is an ultraviolet (UV) and violet light-filtering posterior chamber intraocular lens (IOL), which is intended to provide high-quality vision over a full range of distances, including far, intermediate and near to increase spectacle independence. The lens optic is biconvex and consists of a soft acrylic material with UVAM (a supplemental UV-light absorbing material) and proprietary violet light-filtering chromophore which reduces transmittance of violet light wavelengths.

The lens is designed to be positioned in the lens capsule to replace the optical function of the natural crystalline lens. The lens model ZFR00V incorporates a proprietary wavefront-designed aspheric optic with a squared posterior edge that provides a 360-degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects. The posterior optic of the lens has a proprietary diffractive surface designed to correct chromatic aberration for enhanced contrast sensitivity and provide high-quality vision over a full range of distances. The diffractive design of this lens leads to strong visual acuities across all distances and a continuous range of vision from intermediate through near, resulting in increased spectacle independence. Accommodation will not be restored.

INDICATIONS FOR USE:

The lens model ZFR00V is indicated for primary implantation for the visual correction of: [1] aphakia in adults in whom a cataractous lens has been removed, and [2] aphakia following refractive lensectomy in presbyopic adults who may benefit from useful near vision and reduced spectacle dependence across a range of distances. The lens is intended to be placed in the capsular bag.

PRECAUTIONS:

1. When performing refraction in patients implanted with the lens model ZFR00V, interpret results with caution when using autorefractors or wavefront aberrometers that utilize infrared light, or when performing a duochrome test. Confirmation of refraction with maximum plus manifest refraction technique is strongly recommended.
2. The ability to perform some eye treatments (e.g., retinal photocoagulation) may be affected by the IOL optical design.
3. Recent contact lens usage may affect the patient's refraction; therefore, in contact lens wearers, surgeons should establish corneal stability without contact lenses prior to determining IOL power.
4. Do not resterilize the lens. Most sterilizers are not equipped to sterilize the soft acrylic material without producing undesirable side effects.
5. Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.
6. Do not store the lens in direct sunlight or at a temperature greater than 45°C (113°F). Do not autoclave the intraocular lens.
7. Do not re-use the lens.
8. Prior to implanting, examine the lens package for proper lens model, dioptric power, and expiration date.
9. The surgeon should target emmetropia, as this lens is designed for optimum visual performance when emmetropia is achieved.
10. Care should be taken to achieve centration of the intraocular lens.
11. Please refer to the specific instructions for use provided with the insertion instrument or system for the amount of time the IOL can remain folded before the IOL must be discarded.
12. When the insertion system is used improperly, the haptics of the IOL may become broken. Please refer to the specific instructions for use provided with the insertion instrument or system.
13. Prior to surgery, the surgeon must inform prospective patients of the possible risks and benefits associated with the use of this device and provide a copy of the patient information brochure to the patient.
14. Safety and effectiveness in patients 21 years or younger have not been established in clinical studies.
15. A patient implanted with an intraocular lens should be monitored on a regular basis for long-term postoperative follow-up.
16. The intraocular pressure of implanted patients with glaucoma should be carefully monitored for postoperative changes.

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1. Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the lens may exacerbate an existing condition, may interfere with diagnosis or treatment of a condition or may pose an unreasonable risk to the patient's eyesight:
 a. Recurrent severe anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye.
 b. Patients in whom the intraocular lens may interfere with the ability to observe, diagnose or treat posterior segment diseases.
 c. Surgical difficulties at the time of cataract extraction and/or intraocular lens implantation that might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
 d. A compromised eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
 e. Circumstance that would result in damage to the endothelium during implantation.
 f. Suspected microbial infection.
 g. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support for the IOL.
 h. Congenital bilateral cataracts.
 i. Previous history of, or a predisposition to, retinal detachment.
 j. Patients with potentially good vision in only one eye.
 k. Medically uncontrollable glaucoma.
 l. Corneal endothelial dystrophy.
 m. Proliferative diabetic retinopathy.
 n. Children under the age of 2 years are not suitable candidates for intraocular lenses.
2. The lens model ZFR00V should be placed entirely in the capsular bag. Do not place the lens in the ciliary sulcus.
3. Well-informed patients with well-defined visual needs and preferences should be selected for lens model ZFR00V implantation. The patients should be informed of the possibility of visual effects (such as halo or glare) in nighttime or poor visibility conditions. Patients may perceive these visual effects as an annoyance or hindrance, which, on rare occasions, may be significant enough for the patient to request removal of the IOL.
4. Patients with a predicted postoperative astigmatism greater than 1.0 diopter may not be suitable candidates for lens model ZFR00V implantation since they may not fully benefit in terms of potential spectacle independence.
5. The lens model ZFR00V may affect image quality and lead to some reduction of contrast sensitivity compared to a monofocal lens. Therefore, patients should exercise caution when driving at night or in poor visibility conditions.
6. Johnson & Johnson Surgical Vision, Inc. IOLs are labeled with instructions for use and handling to minimize exposure to conditions that may compromise the product, patient, or the user. The resterilization/reprocessing of Johnson & Johnson Surgical Vision, Inc. IOLs may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

ADVERSE EVENTS:

General Adverse Events for IOLs

Potential adverse events during or following cataract surgery with implantation of the IOL may include but are not limited to:

1. Endophthalmitis/intraocular infection
2. Hypopyon
3. Hyphema
4. IOL dislocation
5. Persistent cystoid macular edema
6. Pupillary block
7. Retinal detachment/tear
8. Persistent corneal stromal edema
9. Persistent iritis
10. Persistent raised intraocular pressure (IOP) requiring treatment.
11. Secondary surgical intervention (including implant repositioning, removal, AC tap, or other surgical procedure).

Adverse events can lead to permanent visual impairment and may require secondary surgical intervention, including intraocular lens exchange or explantation.

DETAILED DEVICE DESCRIPTION:

The lens model ZFR00V is a one-piece foldable posterior chamber lens with an overall diameter of 13.0 mm and an optic diameter of 6.0 mm. The optic and haptics are made of OptiBlue® hydrophobic acrylic, a soft acrylic material with UVAM (a supplemental UV-light absorbing material) and proprietary violet light-filtering chromophore which reduces transmittance of violet light wavelengths. The lens incorporates a proprietary wavefront-designed aspheric optic design on the anterior optic surface that reduces spherical aberration to near zero, and a proprietary diffractive design on the posterior surface designed to compensate for the eye's chromatic aberrations, for enhanced contrast sensitivity. The proprietary diffractive design on the posterior surface leads to strong visual acuities across all distances and a continuous range of vision from intermediate through near, resulting in increased spectacle independence. Lens model ZFR00V is designed for pupil-independent lens performance in any lighting condition.

Lens Optic:

1. Material: Soft acrylic with UVAM and proprietary violet light-filtering chromophore
2. Power: +5.0 to +34.0 diopter powers in 0.5 diopter increments
3. Optic Center Thickness: 0.7 mm (+20.0D)
4. Optic Edge Design: PROTEC 360 square posterior edge
5. Index of Refraction: 1.471 at 35°C
6. Light Transmittance: UV cut-off at 10%T for a +5.0 diopter lens (thinnest), +20.0 diopter lens (mid) and a +34.0 diopter lens (thickest) are shown in **Figure 1**.

Haptics:

1. Material: Soft foldable acrylic with UVAM and proprietary violet light-filtering chromophore
2. One-piece lens
3. Configuration: TRI-FIX design, Modified C, integral with optic
4. Haptic Thickness: 0.46 mm

DIRECTIONS FOR USE:

1. Prior to implanting, examine the lens package for proper lens model, dioptric power, and expiration date.
2. Open the peel pouches and remove the lens in a sterile environment. Verify the dioptric power of the lens.
3. Examine the lens thoroughly to ensure particles have not become attached to it and examine the lens optical surfaces for other defects.
4. The lens may be soaked or rinsed in sterile balanced salt solution until ready for implantation.
5. Handle the lens by the haptic portion. Do not grasp the optical area with forceps.
6. Transfer the lens, using a sterile technique, to an appropriate loading device.
7. The physician should consider the following points:
 - The surgeon should target emmetropia, as this lens is designed for optimum visual performance when emmetropia is achieved.
 - Care should be taken to achieve centration of the intraocular lens.
8. Johnson & Johnson Surgical Vision, Inc. recommends using The UNFOLDER Platinum 1 Series implantation system or an equivalent insertion system or instrument to insert the lens model ZFR00V. Only insertion instruments that have been validated and approved for use with this lens should be used. Please refer to the directions for use with the insertion instrument or system for additional information.

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LENS POWER CALCULATIONS:

Accurate keratometry and biometry are essential to successful visual outcomes. Preoperative calculation of the required spherical equivalent lens power for these posterior chamber intraocular lenses should be determined by the surgeon's experience, preference, and intended lens placement. Emmetropia should be targeted. Accuracy of IOL power calculation is particularly important with the lens model ZFR00V, as spectacle independence is the goal of IOL implantation. The A-constants listed on the outer label are presented as a guideline and serve as a starting point for power calculations. The physician should preoperatively determine the spherical equivalent of the lens to be implanted.

Physicians requiring additional information on lens power calculations may contact the local Johnson & Johnson Surgical Vision, Inc. representative.

Lens power calculation methods are described in the following references:

- Haigis W. "The Haigis formula". In: Shammas HJ, ed, Intraocular Lens Power Calculations. Thorofare, NJ, Slack, 2004; 41-57.
- Hoffer K.J., "The Hoffer Q formula: a comparison of theoretic and regression formulas", J Cataract Refract Surg, **19**, 700-712 (1993). Erratum in: J Cataract Refract Surg 1994;20:677. Erratum in: J Cataract Refract Surg 2007;33:2-3
- Holladay J.T., Prager T.C., Chandler T.Y., Musgrove K.H., Lewis J.W. and Ruiz R.S., "A three-part system for refining intraocular lens power calculations", J Cataract Refract Surg. **14**,17-24 (1988).
- Retzlaff J.A, Sanders D.R. and Kraff M.C., "Development of the SRK/T intraocular lens implant power calculation formula", J Cataract Refract Surg. **16**, 333-340 (1990). Erratum in: J Cataract Refract Surg. 1990;16:528.
- Olsen T. "The Olsen formula". In: Shammas HJ, ed, Intraocular Lens Power Calculations. Thorofare, NJ, Slack, 2004; 27–40
- Canovas C., Artal P. "Customized eye models for determining optimized intraocular lenses power". Biomed. Opt. Express 2011;2:1649-1662

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REPORTING:

Adverse events and/or potentially sight-threatening complications that may reasonably be regarded as lens-related and that were not previously expected in nature, severity, or rate of occurrence must be reported to Johnson & Johnson Surgical Vision, Inc. This information is being requested from all surgeons in order to document potential long-term effects of IOL implantation, especially in younger patients.

Physicians are required to report these events to aid in identifying emerging or potential problems with posterior chamber lenses. These problems may be related to a specific lot of lenses or may be indicative of long-term problems associated with these lenses or with IOLs in general.

HOW SUPPLIED:

Each ZFR00V lens is supplied sterile, in a lens case within a double aseptic transfer peel pouch. The double aseptic transfer peel pouch is sterilized with ethylene oxide and should be opened only under sterile conditions. The pouch and product labels are enclosed in a shelf pack. The external surfaces of the outer pouch are not sterile.

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RETURN/EXCHANGE POLICY:

Contact the local Johnson & Johnson Surgical Vision, Inc. representative for the return/exchange policy. Return the lens with proper identification and the reason for the return. Label the return as a biohazard. Do not attempt to resterilize the lens.

Patient Information

Each patient should receive information regarding intraocular lenses, including information on different IOL options available to them, prior to the decision to implant an intraocular lens.

PRODUCT OF THE NETHERLANDS

Symbol/Explanation

SYMBOL	EXPLANATION
STERILE EO	Sterilized using ethylene oxide
	Do not reuse
	Use-by date (YYYY-MM-DD: Year-Month-Day)
	Consult instructions for use
	Manufacturer
EC REP	Authorized representative in the European Community
	Do not resterilize
45°C (113°F)	Upper limit of temperature
	Keep away from sunlight
	Date of manufacture (YYYY-MM-DD: Year-Month-Day)
	Do not use if package is damaged
REF	Catalogue Number
SN	Serial number

AMO Groningen BV, van Swietenlaan 5, 9728 NX Groningen, Netherlands


TECNIS, PROTEC, TRI-FIX, UNFOLDER and OptiBlue are trademarks of Johnson & Johnson Surgical Vision, Inc.

Johnson & Johnson
Surgical Vision, Inc.
1700 East St. Andrew Pl.,
Santa Ana, CA 92705 USA

EC REP
AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

CE
0344

Figure 1

OptiBlue® spectral transmittance curve is for spherical equivalent (SE) 5-, 20-, and 34-diopter lenses, compared to that of a 53-year-old phakic eye.

Sources: TR8026 (Percent Transmission Spectra in Aqueous for Low-, Mid-, and High-Diopter Power ZXR00V IOLs); Boettner EA and Wolter JR. Transmission of the Ocular Media. Investigative Ophthalmology 1962;1:776-783.

© Johnson & Johnson Surgical Vision, Inc. 2018

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EN

DEVICE DESCRIPTION:

The TECNIS® Symfony® Extended Range of Vision Intraocular Lens (IOL) (model ZXR00 and toric lens models ZXT100, ZXT150, ZXT225, ZXT300, ZXT375, ZXT450, ZXT525, and ZXT600) are ultraviolet-light absorbing posterior chamber IOLs which are intended to provide a continuous range of high-quality vision for far through intermediate distances in the non-toric constant condition for aphakia, the toric IOLs compensate for corneal astigmatism.

TECNIS® Symfony Extended Range of Vision IOLs are designed to be positioned in the lens capsule to replace the optical function of the natural crystalline lens. TECNIS® Symfony IOLs incorporate a proprietary wavefront-designed aspheric or toric optiless optic (ZXR00, ZXT100, ZXT150, ZXT225, ZXT300, ZXT375, ZXT450, ZXT525, ZXT600) with a 360-degree barrier. One end of the optic has a frosted design to reduce unwanted edge glare effects. The posterior optic of the lenses has a proprietary achromatic diffractive surface designed to correct chromatic aberration for enhanced contrast sensitivity and a unique echelette feature to extend the range of vision. TECNIS® Symfony IOLs are manufactured from a hydrophobic, foldable, acrylic material and are designed to minimize glistenings.

INDICATIONS FOR USE:

The TECNIS® Symfony IOL, Model ZXR00, is indicated for primary implantation for the visual correction of aphakia in adult patients with and without presbyopia in whom a cataractous lens has been removed by extracapsular cataract extraction, and optional following refractive lenectomy in presbyopic adults, who desire useful vision over a continuous range of distances (far through intermediate) and who may have low to moderate spectacle independence. The device is intended to be placed in the capsular bag.

The TECNIS® Symfony IOLs, Models ZXT100, ZXT150, ZXT225, ZXT300, ZXT375, ZXT450, ZXT525, and ZXT600, are indicated for primary implantation for the visual correction of aphakia and pre-existing corneal astigmatism in adult patients with and without presbyopia in whom a cataractous lens has been removed by extracapsular cataract extraction, and aphakia following refractive lenectomy in presbyopic adults, who desire useful vision over a continuous range of distances (far through intermediate) and who may have low to moderate spectacle independence. The devices are intended to be placed in the capsular bag.

PRECAUTIONS:

- When performing refractive in patients implanted with the TECNIS® Symfony IOL, interpret results with caution when using autorefractors or wavefront aberrometers that utilize infrared light, or when performing a duochrome test. Confirmation of refractive results with subjective refraction is recommended.
- The ability to perform future retinal laser photocoagulation may be affected by the TECNIS® Symfony IOL optical design.
- Recent cataract surgery may affect the patient's refraction; therefore, in contact lens wearers, surgeons should establish corneal stability without contact lenses before implanting the IOL.
- Do not resterilize the lens. Most sterilizers use an autoclave cycle that uses hot, moist aseptic material without producing undesirable side effects.
- Do not soak or store the intraocular lens in any fluid or solution other than sterile balanced salt solution or sterile normal saline.
- Do not use the lens if the temperature is less than a temperature greater than 113 °F (45 °C). Do not subject the lens to temperature shock.
- Prior to implanting, examine the lens package for proper lens model, dioptric power and expiration date.
- The surgeon should target emmetropia as this lens is designed for optimum visual performance in the emmetropic patient (plano target).
- Care should be taken to achieve centration of the intraocular lens.
- Please refer to the specific instructions for use provided with the inserter instrument or cartridge for the amount of time the IOL can remain loaded into the IOL must be discarded.
- If the TECNIS® Injector is used improperly, the haptics of the TECNIS® Symfony IOL may become damaged. Please refer to the specific instructions for use provided with the insertion instrument or system.

WARNINGS:

Physicians considering lens implantation under any of the following circumstances should weigh the potential benefits against the benefit ratio:

- Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the lens may exacerbate or complicate the condition, may interfere with diagnosis of treatment of a condition or may pose an unreasonable risk of surgical complications:
 - Patients with recurrent severe anterior or posterior segment inflammation or uveitis of unknown etiology, or any disease producing an inflammatory reaction in the eye.
 - Patients in whom the intraocular lens may affect the ability to observe, diagnose, or treat posterior segment disease.
 - Surgical difficulties at the time of cataract extraction, which may increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure or significant intraoperative miosis).
 - A compromised eye due to previous trauma or developmental defects in which the implantation of an IOL may be contraindicated.
 - Circumstances that would result in damage to the endothelium during implantation.
 - Suspected microbial infection.
 - Patients in whom the integrity of the posterior capsule nor the zonules are intact to provide support for the IOL.
 - Congenital bilateral cataracts.
 - Previous history of, or a predisposition to, retinal detachment.
 - Patients with only one eye with potentially good vision.
 - Clinically uncontrollable glaucoma.
 - Corneal guttata or dystrophy.
 - Proliferative diabetic retinopathy.
 - Children under the age of 2 years are not suitable candidates for intraocular lenses.
- The TECNIS® Symfony IOL should be placed entirely in the capsular bag and should not be placed in the ciliary sulcus.
- Patients should be informed of the possibility of visual disturbances.
- The TECNIS® Symfony IOL, Model ZXR00, has not been evaluated for patients with pre-existing corneal astigmatism and may cause visual disturbances.
- Rotation of the TECNIS® Symfony IOL away from the intended axis may result in their astigmatic correction. Misalignment greater than 30° may increase postoperative residual refractive cylinder. If necessary, lens repositioning should occur as early as possible prior to lens encapsulation.
- The TECNIS® Symfony IOL, Model ZXR00, is an optiless design. Residual viscoelastic may allow the lens to rotate, causing misalignment of the toric TECNIS® Symfony IOL with the intended axis of placement.
- ABO IOLs are supplied non-sterile and intended for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The manufacturer recommends that the IOL be used as soon as possible after it is removed from the package. The medical device is intended as a single-use device. Damage, illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

ADVERSE EVENTS:

General Undesirable Events for IOLs

Potential adverse events during or following cataract surgery with implantation of an IOL may include but are not limited to:

- Endophthalmitis/intraocular infection
- Vitreous loss
- IOL dislocation
- Pupillary block
- Retinal detachment
- Persistent corneal stromal edema
- Chronic uveitis/intraocular inflammation
- Secondary surgical intervention including implant repositioning, removal, AC tap, or other surgical procedures)
- Any other adverse event that leads to permanent visual impairment or requires surgical or medical intervention to prevent permanent visual impairment

DETAILED DESCRIPTION:

The TECNIS® Symfony IOLs (lens model ZXR00 and toric lens models ZXT100, ZXT150, ZXT225, ZXT300, ZXT375, ZXT450, ZXT525, and ZXT600) have a 6.0 mm posterior optic diameter (lenses with an overall diameter of 13.0 mm and an optic diameter of 6.0 mm). They incorporate a proprietary aspheric optic or toric aspheric optic design on an anterior optic surface that reduces spherical aberration to near zero, and a diffractive design on the posterior surface designed to compensate for the eye's natural spherical aberration. The lenses also have a proprietary achromatic diffractive surface and increase spectacle independence, compared to a standard monofocal IOL that does not have these design features. In addition, the toric IOLs compensate for corneal astigmatism while achieving the ANSI Standard for toric intraocular lens rotational stability requirement (+30% and having ±5° axis change between consecutive visits approximately three months apart).

Lens Optics:

- Material: Optically clear, soft foldable hydrophobic acrylic with a covalently bound UV absorber
- Optic: 6.0 to +34.0 diopter powers in 0.5 diopter increments
- Optic Power: toric lens models ZXT100, ZXT150, ZXT225, ZXT300, ZXT375, ZXT450, ZXT525, and ZXT600: 1.00 diopter, 1.50 diopters, 2.25 diopters, 3.00 diopters, 3.75 diopters, 4.50 diopters, 5.25 diopters and 6.00 diopters (as measured at the IOL plane)

Conversion table for cylinder powers:

	Cylinder Power (D)								
IOL Plane (Spectacle)	1.00	1.50	2.25	3.00	3.75	4.50	5.25	6.00	
Corneal Plane*	0.69	1.03	1.54	2.06	2.57	3.08	3.60	4.11	

*The corresponding cylinder values at the corneal plane have been calculated based on the average postoperative ACD of 3.15 mm.

- Optic Center Thickness: 0.7 mm (+20.0 D)
- Optic Edge Design: 360° barrier, square posterior edge
- Index of Refraction: 1.47 at 35 °C
- Light Transmission: UV cut-off at 10% T for a +3.0 diopter lens (thickest) and a 1.47 diopter lens (thickest) are shown in Figure 1.

Haptics:

- Material: Soft foldable acrylic with a covalently bound UV absorber.
- Two-piece lens
- Configuration: TR-Fit design Modified C, integral with optic
- Haptic thickness: 0.46 mm

DIRECTIONS FOR USE:

- Prior to implanting, examine the lens package for IOL type, power, proper configuration and expiration date.
- Do not use the lens if the package shows that the lens is a sterile environment. Verify the dioptric power of the lens.
- Examine the lens thoroughly to ensure particles have not become attached to it, and examine the lens' optical surfaces for other defects.
- The lens may be soaked or rinsed in sterile balanced salt solution until ready for implantation.
- Handle the lens by the haptic portion. Do not grasp the optical area with forceps.
- Transfer the lens, using a sterile technique, to an appropriate loading device.
- The procedure should not compromise the following points:
 - the surgeon should target emmetropia, as this lens is designed for optimum visual performance in the emmetropic patient (plano target).
 - Care should be taken to achieve centration of the intraocular lens.
- AMO recommends using the UNFOLDER Platinum 1 Series implantation system or an equivalent insertion instrument or system to insert the TECNIS® Symfony lenses. Only insertion instruments that have been validated and approved for use with the TECNIS® Symfony lens should be used. Please refer to the directions for use with the insertion instrument or system for additional information.

CAUTION:

Do not use the lens if the package has been damaged. The sterility of the lens may have been compromised.

SELECTION AND PLACEMENT OF THE TORIC TECNIS® Symfony IOL (MODELS ZXT100, ZXT150, ZXT225, ZXT300, ZXT375, ZXT450, ZXT525, AND ZXT600):

The astigmatism to be corrected should be determined from keratometry and biometry rather than refraction data since the presence of refractive astigmatism in the cataractous lens to be removed may influence results. The size and location of the incision(s) can affect the amount of astigmatism present postoperatively. Prior to selection IOL selection and placement, AMO provides a free calculation tool for the surgeon. Preoperative keratometry and biometry data, incision location, and the patient's estimated surgically induced corneal astigmatism are used to determine the appropriate toric TECNIS® Symfony IOL model, cylinder power, and axis of placement in the eye.

The TECNIS® Symfony toric IOL is indicated for primary implantation for the placement of the lens within the capsular bag. The anterior surface of the IOL is marked with indentations that are aligned with the haptic-optic junction that identify the flat meridian of the toric TECNIS® Symfony lens optic. These indentations form an alignment axis that is used to orient the lens in the eye. The lens should be implanted so that the lens TECNIS® Symfony IOL markings are aligned with the intended steep axis of the post-incision steep corneal meridian (intended axis of placement).

Prior to surgery, the operative eye should be marked in the following manner: With the patient sitting upright, precisely mark the twelve o'clock and/or the six o'clock position with an 11-mm-wide surgical skin marker, or a surgical pen, centered on the ophthalmic axis. Using these marks as reference points, an axis marker can be used to transfer the 0°–180° axis to the patient's eye. Ensure that the patient remains upright during the axis marking to avoid the risk of cyclotorsion to the ophthalmic axis of placement.

After the lens is inserted, precisely align the axis marking indentations on the toric TECNIS® Symfony IOL with the marked axis of lens placement. Carefully remove all viscoelastic from the capsular bag. This may be accomplished by manipulating the IOL optic with a Simcoe cannula or by irrigation/aspiration. Ensure that all residual viscoelastic from the lens' trifocal transition zone has been removed. If preferred, to ensure removal of viscoelastic from behind the lens implant. Special care should be taken to ensure proper positioning of the lens TECNIS® Symfony IOL optic edge after following viscoelastic removal. Residual viscoelastic may allow the lens to rotate, causing misalignment of the toric TECNIS® Symfony IOL with the intended axis of placement.

Misalignment of the axis of the lens with the intended axis of placement may compromise the astigmatic correction. Such misalignment can result from inaccurate keratometry or marking of the cornea, inaccurate placement of the toric TECNIS® Symfony IOL's axis mark, or rotation of the implant after placement in the eye. Misalignment of the AMO toric calculator to determine the optimum axis of placement. To minimize this effect, the surgeon should be careful to ensure that preoperative keratometry and biometry is accurate and that the IOL is properly oriented prior to the end of surgery.

POWER CALCULATIONS:

Accurate keratometry and biometry is essential to successful visual outcomes. Keratometry and biometry should be performed on each patient to determine the posterior chamber intraocular lenses should be determined by the surgeon's preferred method. Recommended A-constants are provided below. Optimized A-constants may be determined by the surgeon's preferred method. Use of the AMO toric calculator as spectacle independence is a goal of TECNIS® Symfony IOL implantation. The A-constants listed on the outer label are presented as a guideline and serve as a starting point for implant power calculations. The physician should determine preoperatively the appropriate lens power to target for the patient's vision correction.

For toric lens models: Use of the AMO unveiled toric calculator tool is recommended for determining the appropriate toric TECNIS® Symfony IOL model, optimal axis of IOL placement and cylinder power. The toric TECNIS® Symfony IOL is indexed with the IOL distance spherical equivalent power.

Please refer to the calculator instructional information on lens power calculations may contact the local AMO representative. Lens power calculation methods are described in the following publications:

- Hoffer K.J., "The Hoffer Q formula: A comparison of theoretic and regression Calculations, Thorofare, NJ, Slack, 2004: 41-57.
- Hoffer K.J., "The Hoffer Q formula: a comparison of theoretic and regression formulas," J Cataract Refract Surg, 19: 700-712 (1993). Erratum in: J Cataract Refract Surg 1994;20:339.
- Holladay J.T., Prager T.C., Chandler T.Y., Musgrove K.H., Lewis J.W. and Ruiz R.S., "A three-part system for refining intraocular lens power calculations", J Cataract Refract Surg, 14:17-24 (1988).
- Retzlaff J.A, Sanders D.R. and Kraff M.C., "Development of the SRK/T intraocular lens implant power calculation formula", J Cataract Refract Surg, 16: 333-340 (1990). Erratum in: J Cataract Refract Surg 1990;16:528.
- Olsen T., "The Olsen formula", In: Shammas HJ, ed. Intraocular Lens Power Calculations. Thorofare, NJ, Slack, 2004: 27-40.
- Canovas C, Arbelaez J, "Customized eye models for determining optimized intraocular lens power", Biomed. Opt. Express 2011;2:1649-1662

PATIENT CARD:

Along with the implantation card, to be supplied to the patient, is included in the package. The patient should be instructed to keep the card as a permanent record. It is important and to show the card to any eye care practitioner he/she may see in the future.

REPORT:

Any adverse events and/or potentially sight-threatening complications that may reasonably be regarded as lens-related and that were not previously reported (expected in nature, severity or rate of occurrence must be reported to AMO. This information is being requested to assist AMO in the ongoing evaluation of the safety and performance of its IOL products, especially its elongate patients.

AMO is also requesting that these events be laid in identifying emerging or potential problems with posterior chamber lenses. These patients may be related to a specific lot of lenses or may be indicative of long-term problems associated with these lenses or with IOLs in general.

HOW SUPPLIED:

The TECNIS® Symfony lenses are supplied sterile in a lens case within a double aseptic transfer peel pouch. The double aseptic transfer pouch is sterilized with ethylene oxide and should be opened only under sterile conditions. Each lens case and pouch label are enclosed in a shelf pack. The external surfaces of the outer pouch are not sterile.

EXPIRATION DATE:

The expiration date on the lens package is the sterility expiration date. This lens should not be implanted if the package is damaged or past the expiration date.

RETURN/EXCHANGE POLICY:

Each patient should receive information regarding intraocular lenses prior to the decision to undergo an intraocular lens.

SYMBOLS/EXPLANATION:

SYMBOL	EXPLANATION
STERILE EO	Sterilized Using Ethylene Oxide
	Do Not Reuse
	Use By (YYYY-MM-DD: Year-Month-Day)
	Consult Instructions for Use
	Manufacturer
EC REP	Authorized Representative in the European Community
	Do Not Resterilize
	Upper Limit of Temperature
	Keep Away From Sunlight
	Date of Manufacture
	Do not use if package is damaged

PRODUCT OF THE NETHERLANDS

Groningen BV, 9728 NX Groningen, Netherlands

Johnson & Johnson Surgical Vision, Inc.
1700 E. St. Andrew Place
Santa Ana, CA 92705 USA

0344

AMO Ireland
Block B
Liffey Valley Office Campus
Dublin 22, Co. Dublin, Ireland

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
TECNIS, TECNIS Symfony, PROTEC, TR-FIX and UNFOLDER are trademarks of Johnson & Johnson Surgical Vision, Inc.

© Johnson & Johnson Surgical Vision Inc. 2018

Figure 1
Light Transmission

LEGEND:
Spectral transmittance curve of a typical 5-diopter IOL (thinnest). UV cut-off at 10%T is 374 nm
Spectral transmittance curve of a typical 20-diopter IOL. UV cut-off at 10%T is 375 nm
Spectral transmittance curve of a typical 34-diopter IOL (thickest). UV cut-off at 10%T is 375 nm
Spectral transmittance curve of a 53-year-old phakic eye, from *Boettner, E.A., and Wolter J.R. Transmission of the Ocular Media. Investigative Ophthalmology. 1962;1:776-783.*

Note: The cut-off wavelengths and the spectral transmittance curves represent the range of the transmittance of the IOLs (s-340) made with this material. Spectral transmission measurements were taken in water at room temperature (ref: TR7475).

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Rx Only

Device Description

The TECNIS® 1-piece lens is an ultraviolet-light absorbing posterior chamber intraocular lens (IOL) which compensates for corneal spherical aberrations. Accommodation will not be restored. It is designed to be positioned in the lens capsule where the lens should replace the optical function of the natural crystalline lens. AMO 1-piece IOLs incorporate a proprietary wavefront-designed aspheric optic with a squared posterior optic edge designed to provide a 360 degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects. The image quality of the TECNIS® optic is illustrated in Figure 1.

Indications For Use

AMO TECNIS® 1-piece lenses are indicated for the visual correction of aphakia in adult patients in whom a cataractous lens has been removed by extracapsular cataract extraction. These devices are intended to be placed in the capsular bag.

Precautions

Do not resterilize the lens. Most sterilizers are not equipped to sterilize the soft acrylic material without producing undesirable side effects.

Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.

Do not store the lens in direct sunlight or at a temperature greater than 113°F (45°C). Do not autoclave the intraocular lens.

Do not reuse lens.

Please refer to the specific instructions for use provided with the insertion instrument or system for the amount of time the IOL can remain folded before the IOL must be discarded. When the insertion system is used improperly, the haptics of the TECNIS® 1-piece lens may become damaged. Please refer to the specific instructions for use provided with the insertion instrument or system.

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1. Patients with recurrent severe anterior or posterior segment inflammation or uveitis.
2. Patients in whom the intraocular lens may affect the ability to observe, diagnose or treat posterior segment diseases.
3. Surgical difficulties at the time of cataract extraction, which may increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure or significant vitreous prolapse or loss).
4. A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
5. Circumstances that would result in damage to the endothelium during implantation.
6. Suspected microbial infection.
7. Patients in whom neither the posterior capsule nor the zonules are intact enough to provide support for the IOL.
8. Children under the age of 2 years are not suitable candidates for intraocular lenses.
9. AMO IOLs are single-use medical devices and are labeled with instructions for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The reuse/resterilization/reprocessing of AMO single-use medical devices may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

The TECNIS® 1-piece IOL should be placed entirely in the capsular bag and should not be placed in the ciliary sulcus.

Adverse Events

As of August 10, 2007, the incidence of adverse events experienced during the clinical trial for Model AAB00 is similar to or less than those of the historic control population (FDA Grid for posterior chamber IOLs) as shown in Table 1.

Clinical Trial (SENSAR 1-Piece Lens, Model AAB00):

The parent lens of the TECNIS® ZCB00 Lens is the Model AAB00. The difference between lens model ZCB00 and lens model AAB00 is the modified prolate anterior optic surface of lens model ZCB00. The clinical results of the investigation of lens model AAB00 are pertinent to lens model ZCB00. A multi-center, unilateral, open-label, non-comparative clinical evaluation was performed for the SENSAR 1-piece lens, Model AAB00. The clinical trial was initiated on November 30, 2005. The purpose of the study was to evaluate the safety and effectiveness of lens model AAB00 in subjects undergoing cataract removal and intraocular lens implantation. Following routine cataract removal by extracapsular cataract extraction, all IOLs were implanted in the capsular bag with a continuous curvilinear capsulorhexis.

The results achieved by 117 patients followed for one year provide the basis for the data supporting the use of this lens design for visual correction of aphakia. In the total study population (123 patients), 56.9% of the patients were female and 43.1% were male; 93.5% were Caucasian, 4.1% were Black and 2.4% were Asian. The best corrected distance visual acuity results for the "best case" patients at 1 year (330-420 days) postoperatively are provided in Table 2. In addition the data compared to the FDA Grid values (historical control) are presented in Table 3.

Clinical Trial (TECNIS® Lens, Model Z9000):

In a controlled, multi-center, intra-individual clinical investigation of the Z9000 lens (modified prolate anterior surface) and an acrylic lens (spherical optic), ocular spherical aberrations were significantly less with the TECNIS® lens than with the acrylic lens. The simulated night driving results (functional vision) under several of the conditions tested and the visual acuity results were statistically significantly better in eyes implanted with the TECNIS® lens (TECNIS® eyes). The clinical significance of the reduction of ocular spherical aberration was to be demonstrated using contrast sensitivity testing. Functional significance was to be demonstrated using simulated night driving.

Spherical Aberration

The mean ocular spherical aberration of the TECNIS® eyes was not significantly different from zero. This was not true for eyes implanted with the lens with the spherical optic. The mean difference in ocular spherical aberration between the two eyes of patients was statistically significantly different from zero. Figure 2 provides the mean spherical aberration measurements of all eyes with evaluable wavefront measurements. As evidenced in the figure, the significant reduction in the spherical aberration in the TECNIS® eyes was independent of age.

Figure 3 present the wavefront measurements of the 22 patients for whom evaluable data were available for both eyes.

Visual Acuity

The monocular visual acuity results (90 ± 15 days postoperatively) of each patient in the Safety Population and in the subset of patients who underwent wavefront measurement and night driving simulation are presented in Table 4.

Contrast Sensitivity

The primary objective of the clinical investigation was to demonstrate the mesopic (6 cd/m²) intra-individual difference in the postoperative quality of vision using sine-wave contrast sensitivity testing between the TECNIS® lens (Z9000) and a lens with a spherical optic. In this clinical investigation, the contrast sensitivity results were not significantly different. The mesopic log contrast sensitivity results at all spatial frequencies tested for the Z9000 and the control lens are presented in Figure 4.

Simulated Night Driving

A subset of patients (29) randomly selected from all of the investigational sties underwent testing in a validated night driving simulator. Patients were tested monocularly in simulated city normal, city glare, rural normal and rural glare lighting conditions.

The night driving simulator consisted of an automobile cab/frame with a windshield, video scene and target projectors, glare sources, a display screen and a computer. The front cab included a windshield with a rear-view mirror, a non-functioning dashboard, a door-mounted side-view mirror, front seats with contoured headrests, seat belts and a steering wheel. The ambient lighting of the simulator was similar to average nighttime scenes.

The nighttime city driving scene was of a long, straight city street with a simulated traveling speed of 35 miles per hour with a variety of street lights, cars, store lights and signs creating a high degree of ambient lighting. The nighttime rural driving scene was of a long, straight country road with a traveling speed of 55 miles per hour and minimum ambient lighting. Each driving scene was about 30 seconds in duration.

For testing under glare conditions, the constant size glare source was a simulation of a real-life headlight disability glare from a following vehicle reflected in rear- and side-view mirrors adjusted to shine in the eyes of the patient. The amount of glare was set to produce a 10% loss in detection distance.

Patients were asked to detect and identify targets, including white-on-green information highway signs, black-on-yellow warning signs and pedestrian hazards. They were asked to respond when the sign or the hazard was first detected, and the detection distances were recorded. Patients were then asked to respond when the sign or hazard could first be identified, i.e., what did the sign say, what direction was the pedestrian walking, and the identification distances were recorded. The patient responses for each target set and visibility condition were averaged.

Figures 5 and 6 present the average difference between the detection and identification distances with testing of the Z9000 eye and the detection and identification distances with testing of the spherical optic lens of each subject (the mean of the intra-individual differences).

The Z9000 eyes performed functionally better than the control eyes in 21 of the 24 conditions tested. This means the Z9000 lens improves both detection and identification distances across the driving scenes (city and rural) and visibility conditions (with/without glare) compared to the control lens. Z9000 eyes performed statistically significantly better than the control eyes in 9 of the test conditions. The greatest advantage of the Z9000 lens is for increased detection and identification of the pedestrian hazard under rural visibility conditions with and without glare. Under these conditions, the increased visibility distance at 55 mile per hour provides for an average of about 0.5 seconds more time to perceive and react to the pedestrian hazard. A 0.5 second increase in perception and reaction time is functionally significant in increasing the time to take evasive action, time to stop or effect of impact.

These findings suggest there is likely to be a meaningful safety benefit to elderly drivers with TECNIS® lenses, and to the drivers and pedestrians with whom they share the road. The results of this performance/functional test demonstrate that the TECNIS® lens improved functional vision, which in turn may improve patient safety for other life situations under low visibility conditions.

Detailed Device Description

Lens Optic

1. Optic Material: Optically clear, soft foldable hydrophobic acrylic with a covalently bound UV absorber. Full transmission of the blue wavelength light for optimal scotopic sensitivity.
2. Power: +5.0 to +34.0 diopter powers in 0.5 diopter increments.
3. Optic Center Thickness: 0.722 mm (+20.0D)
4. Optic Edge Design: PROTEC 360 Square posterior edge
5. Index of Refraction: 1.47 at 35°C.
6. Light Transmittance: UV cut-off at 10% T for a +5.0 diopter lens (thinnest) and a +34.0 diopter lens (thickest) are shown in Figure 7.

Haptics

1. Material: Soft foldable hydrophobic acrylic with a covalently bound UV absorber.
2. One-piece lens.
3. Configuration: TRI-FIX design, Modified C, integral with optic
4. Haptic Thickness: 0.46 mm

Directions For Use

1. Prior to implanting, examine the lens package for IOL type, power, proper configuration and expiration date.
2. Open the peel pouch and remove the lens in a sterile environment. Verify the dioptric power of the lens.
3. Examine the lens thoroughly to ensure particles have not become attached to it, and examine the lens optical surfaces for other defects.
4. If desired, the lens may be soaked or rinsed in sterile balanced salt solution until ready for implantation.
5. AMO recommends using Duckworth & Kent 7722 folding forceps with the 7741 insertion forceps or an equivalent insertion instrument or system to insert the TECNIS® 1-Piece lens. Only insertion instruments that have been validated and approved for use with this lens should be used. Please refer to the directions for use with the insertion instrument or system for additional information.

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Lens Power Calculations

The physician should determine preoperatively the power of the lens to be implanted. Lens power calculation methods are described in the following references:

1. Hoffer, K.J. The Hoffer Q formula: a comparison of theoretic and regression formulas. Journal of Cataract and Refractive Surgery. 1993; 19:700-712; ERRATA; 1994; 20:677.
2. Holladay, J.T., Musgrove, K.H., Prager, T.C., Lewis, J.W., Chandler, T.Y., and Ruiz, R.S.. A three-part system for refining intraocular lens power calculations. Journal of Cataract and Refractive Surgery. 1988; 14:17-24.
3. Holladay, J.T. Standardizing constants for ultrasonic biometry, keratometry and intraocular lens power calculations. Journal of Cataract and Refractive Surgery. 1997; 23:1356-1370.
4. Norrby NES. Unfortunate discrepancies. Letter to the editor and reply by Holladay, J.T. Journal of Cataract and Refractive Surgery. 1998; 24:433-434.
5. Olsen, T., Olesen, H., Thim, K., and Corydon, L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. Journal of Cataract and Refractive Surgery. 1992; 18: 280-285.
6. Retzlaff, J.A., Sanders, D.R. and Kraff, M.C. Development of the SRK/T intraocular lens implant power calculation formula. Journal of Cataract and Refractive Surgery. 1990; 16:333-340; ERRATA, 1990; 16:528.

Physicians requiring additional information on lens power calculation may contact AMO.

Patient Registration Section (For US)

Each patient who receives a TECNIS® 1-piece IOL must be registered with AMO at the time of lens implantation.

Registration is accomplished by completing the Implant Registration Card that is enclosed in the lens package and mailing it to AMO. Patient registration is essential for AMO's long-term patient follow-up program and will assist AMO in responding to Adverse Reaction Reports and/or potentially sight-threatening complications.

An Implant Identification Card is supplied in the lens package. This card should be given to the patient with instructions to keep it as a permanent record of the implant and to show the card to any eye care practitioner seen in the future.

Reporting

All adverse events, regardless of severity and whether or not attributed to the implant, are to be reported to AMO at 1-877-266-4543 (U.S. only) or by contacting your local AMO representative. In the event of a life-threatening incident or serious adverse event, AMO must be notified immediately (no later than 48 hours upon detection) by phone and by faxing a completed adverse event form.

How Supplied

The TECNIS® 1-piece lenses are supplied sterile in a lens case within a double aseptic transfer peel pouch. The double aseptic transfer peel pouch is sterilized with ethylene oxide and should be opened only under sterile conditions. The pouch and product labels are enclosed in a shelf pack. The external surfaces of the outer pouch are not sterile.

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Return/Exchange Policy

Please contact your local AMO office regarding lens return or exchange.

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AMO, TRI-FIX and PROTEC, SENSAR and TECNIS are trademarks of Johnson & Johnson Surgical Vision, Inc.

Symbol/Explanation:

SYMBOL	ENGLISH
STERILE EO	STERILIZED USING ETHYLENE OXIDE
	DO NOT REUSE
	USE BY (YYYY-MM-DD: year-month-day)
	CONSULT INSTRUCTIONS FOR USE
	DO NOT RESTERILIZE
	MANUFACTURER
EC REP	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN UNION
45°C (113°F)	UPPER LIMIT OF TEMPERATURE
	DO NOT USE IF PACKAGE IS DAMAGED
	KEEP AWAY FROM SUNLIGHT
	DATE OF MANUFACTURE (YYYY-MM-DD: year-month-day)

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Table 1
Adverse Events Model AAB00
All Subjects (N = 123)

ADVERSE EVENTS	CUMULATIVE		PERSISTENT AT 1 YEAR		FDA GRID	
	N	%	N	%	CUM %	PER %
Persistent Corneal Edema	-	-	0	0.0	-	0.3
Cystoid Macular Edema (CME)	4	3.3*	1	0.9†	3.0	0.5
Endophthalmitis	0	0.0	-	-	0.1	-
Hyphema	0	0.0	-	-	2.2	-
Hypopyon	0	0.0	-	-	0.3	-
Persistent Iritis	-	-	0	0.0	-	0.3
Secondary Surgical Intervention - Pars Plana Vitrectomy with Membrane Peel	1	0.8	-	-	0.8	-
Lens Dislocation	0	0.0	-	-	0.1	-
Pupillary Block	0	0.0	-	-	0.1	-
Retinal Detachment	0	0.0	-	-	0.3	-
Persistent Raised IOP Requiring Treatment	-	-	0	0.0	-	0.4
Lens Exchange - Torn Haptic related to improper loading technique	1	0.8	-	-	-	-

* This rate is not statistically significantly higher than the FDA Grid cumulative rate for posterior chamber IOLs of 3.0% (p=0.5060).
† This rate is not statistically significantly higher than the FDA Grid rate for posterior chamber IOLs of 0.5% (p=0.4437).

TABLE 2
Best Corrected Distance Visual Acuity (Snellen Equivalent) at 1 Year
Best Case Subjects* (N = 110)

AGE GROUP	N	20/20 or BETTER		20/25 to 20/40		20/50 to 20/100		20/125 or WORSE	
		n	%	n	%	n	%	n	%
< 60	11	11	100.0	0	0.0	0	0.0	0	0.0
60-69	35	29	82.9	6	17.1	0	0.0	0	0.0
70-79	46	39	84.8	7	15.2	0	0.0	0	0.0
> 80	18	14	77.8	4	22.2	0	0.0	0	0.0
TOTAL†	**110**	**93**	**84.5**	**17**	**15.5**	**0**	**0.0**	**0**	**0.0**


* Excludes subjects with macular degeneration at any time during the study
† Includes three subjects who experienced a Nd:YAG posterior capsulotomy

TABLE 3
Best Corrected Distance Visual Acuity (Snellen Equivalent) at 1 Year
Best Case Subjects* (N = 110) vs. FDA Grid

AGE DECADE	TOTAL		VISUAL ACUITY 20/40 OR BETTER		FDA GRID
	N	%	N	%	%
< 60	11	10.0	11	100.0	98.5
60 - 69	35	31.8	35	100.0	96.5
70 - 79	46	41.8	46	100.0	97.5
> 80	18	16.4	18	100.0	94.8
TOTAL†	**110**	**100.0**	**110**	**100.0**	**96.7**

* Excludes subjects with macular degeneration at any time during the study
† Includes three subjects who experienced a Nd:YAG posterior capsulotomy

Figure 1

* The image quality of the TECNIS® 1-Piece lens was characterized by measuring modulation transfer function (MTF) in the Average Cornea Eye (ACE) model. The ACE model was designed using wavefront technology to reproduce the average amount of spherical aberration of the cornea. Measurements taken using the ACE eye model reflect the optical benefit of the TECNIS® aspheric design feature.
† MTF of the TECNIS® 1-Piece lens as measured using a model eye as specified in ISO 11979-2, Ophthalmic implants - Intraocular lenses, Part 2: Optical properties and test methods.

Figure 2



Figure 3

* ITT - intent-to-treat population: All randomized patients undergoing bilateral cataract surgery for whom postoperative data were available

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* One-sided paired t-test
** Difference between logMAR visual acuity in the eye receiving Z9000 and the eye receiving the control lens (spherical optic) for each patient.

Figure 4


Figure 5


Figure 6



Figure 7

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Note: *The cut-off wavelengths and the spectral transmittance curves represent the range of the transmittance values of IOLs (5-34 diopter) made with this material.*

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1. Hoffer, K.J. The Hoffer Q formula: a comparison of theoretic and regression formulas. Journal of Cataract and Refractive Surgery. 1993; 19:700-712; ERRATA; 1994; 20:677.
2. Holladay, J.T., Musgrove, K.H., Prager, T.C., Lewis, J.W., Chandler, T.Y., and Ruiz, R.S.. A three-part system for refining intraocular lens power calculations. Journal of Cataract and Refractive Surgery. 1988; 14;17-24.
3. Holladay, J.T. Standardizing constants for ultrasonic biometry, keratometry and intraocular lens power calculations. Journal of Cataract and Refractive Surgery. 1997; 23:1356-1370.
4. Norrby NES. Unfortunate discrepancies. Letter to the editor and reply by Holladay, J.T. Journal of Cataract and Refractive Surgery. 1998; 24:433-434.
5. Olsen, T., Olesen, H., Thim, K., and Corydon, L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. Journal of Cataract and Refractive Surgery. 1992; 18: 280-285.
6. Retzlaff, J.A., Sanders, D.R. and Kraff, M.C. Development of the SRK/T intraocular lens implant power calculation formula. Journal of Cataract and Refractive Surgery. 1990; 16:333-340; ERRATA, 1990; 16:528.

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DEVICE DESCRIPTION:

The TECNIS® Toric 1-piece lens is an ultraviolet-light absorbing posterior chamber intraocular lens (IOL) which compensates for corneal spherical aberrations and corneal astigmatism. It is designed to be positioned in the lens capsule where the lens should replace the optical function of the natural crystalline lens. The TECNIS® Toric 1-piece IOL incorporates a proprietary wavefront-designed toric aspheric optic with a squared posterior optic edge designed to provide a 360 degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects. The anteriorly located cylinder axis marks denote the meridian with the lowest power.

INDICATIONS FOR USE:

The TECNIS® Toric 1-piece lens is indicated for the visual correction of aphakia and pre-existing corneal astigmatism in adult patients with or without presbyopia, in whom a cataractous lens has been removed by extracapsular cataract extraction and who desire improved uncorrected distance vision, reduction of residual refractive cylinder and increased spectacle independence for distance vision. This device is intended to be placed in the capsular bag.

PRECAUTIONS:

1. Do not resterilize the lens. Most sterilizers are not equipped to sterilize the soft acrylic material without producing undesirable side effects.
2. Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.
3. Do not store the lens in direct sunlight or at a temperature greater than 113°F (45°C). Do not autoclave the intraocular lens.
4. Please refer to the specific instructions for use provided with the insertion instrument or system for the amount of time the IOL can remain folded before the IOL must be discarded.
5. When the insertion system is used improperly, the haptics of the TECNIS® Toric 1-piece lens may become broken. Please refer to the specific instructions for use provided with the insertion instrument or system.

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1. Patients with recurrent severe anterior or posterior segment inflammation or uveitis.
2. Patients in whom the intraocular lens may affect the ability to observe, diagnose or treat posterior segment diseases.
3. Surgical difficulties at the time of cataract extraction, which may increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure or significant vitreous prolapse or loss).
4. A compromised eye due to previous trauma or developmental defects in which appropriate support of the IOL is not possible.
5. Circumstances that would result in damage to the endothelium during implantation.
6. Suspected microbial infection.
7. Patients in whom neither the posterior capsule nor the zonules are intact enough to provide support for the IOL.
8. Children under the age of 2 years are not suitable candidates for intraocular lenses.
9. The TECNIS® Toric 1-piece IOL should be placed entirely in the capsular bag and should not be placed in the ciliary sulcus.
10. Rotation of TECNIS® Toric 1-piece IOLs away from their intended axis can reduce their astigmatic correction. Misalignment greater than 30° may increase postoperative refractive cylinder. If necessary, lens repositioning should occur as early as possible prior to lens encapsulation.
11. Carefully remove all viscoelastic from the capsular bag. Residual viscoelastic may allow the lens to rotate, causing misalignment of the TECNIS® Toric 1-piece IOL with the intended axis of placement.
12. AMO IOLs are single-use medical devices and are labeled with instructions for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The reuse/resterilization/reprocessing of AMO single-use medical devices may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

ADVERSE EVENTS:

Potential adverse events during or following surgery of the TECNIS® Toric 1-piece lens may include but are not limited to:

- Endophthalmitis
- Hypopyon
- Pupillary block
- Retinal detachment
- Acute corneal decompensation
- Secondary surgical intervention (including implant repositioning, removal, AC tap or other surgical procedure).

DETAILED DEVICE DESCRIPTION:

Lens Optic:

- Optic Material: Optically clear, soft foldable acrylic with a covalently bound UV absorber.
- Power: +5.0 to +34.0 diopter powers in 0.5 diopter increments.
- Cylinder Power: 1.00 diopter, 1.50 diopter, 2.25 diopter, 3.00 diopter, 3.75 diopter, 4.00 diopter, 4.50 diopter, 5.25 diopter, 6.00 diopter, 7.00 diopter and 8.00 diopter (as measured at the IOL plane).

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- Optic Center Thickness: 0.7 mm (+20.0D).
- Optic Edge Design: PROTEC 360 Square posterior edge.
- Index of Refraction: 1.47 at 35°C.
- Light Transmittance: UV cut-off at 10% T for a +5.0 diopter lens (thinnest) and a +34.0 diopter lens (thickest) are shown in Figure 1.

Haptics:

- Material: Soft foldable acrylic with a covalently bound UV absorber.
- One-piece lens.
- Configuration: TRI-FIX design Modified C, integral with optic.
- Haptic Thickness: 0.46 mm.

LENS POWER CALCULATIONS:

Accurate keratometry and biometry is essential to successful visual outcomes. Preoperative calculation of the required spherical equivalent lens power for these posterior chamber intraocular lenses should be determined by the surgeon's experience, preference, and intended lens placement. The A-constant listed on the outer label is presented as a guideline and is a starting point for implant power calculations. The physician should determine preoperatively the spherical equivalent and cylindrical power of the lens to be implanted.

TECNIS® Toric 1-piece lenses are labeled with the IOL spherical equivalent power. Lens power calculation methods are described in the following references:

- Hoffer, K.J. The Hoffer Q formula: a comparison of theoretic and regression formulas. Journal of Cataract and Refractive Surgery. 1993; 19:700-712; ERRATA; 1994; 20:677.
- Holladay, J.T., Musgrove, K.H., Prager, T.C., Lewis, J.W., Chandler, T.Y., and Ruiz, R.S.. A three-part system for refining intraocular lens power calculations. Journal of Cataract and Refractive Surgery. 1988; 14:17-24.
- Holladay, J.T. Standardizing constants for ultrasonic biometry, keratometry and intraocular lens power calculations. Journal of Cataract and Refractive Surgery. 1997; 23:1356-1370.
- Norrby NES. Unfortunate discrepancies. Letter to the editor and reply by Holladay, J.T. Journal of Cataract and Refractive Surgery. 1998; 24:433-434.
- Olsen, T., Olesen, H., Thim, K., and Corydon, L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. Journal of Cataract and Refractive Surgery. 1992; 18: 280-285.
- Retzlaff, J.A., Sanders, D.R. and Kraff, M.C. Development of the SRK/T intraocular lens implant power calculation formula. Journal of Cataract and Refractive Surgery. 1990; 16:333-340; ERRATA, 1990; 16:528.

DIRECTIONS FOR USE:

1. Prior to implanting, examine the lens package for IOL type, power, proper configuration and expiration date.
2. Open the peel pouches and remove the lens in a sterile environment. Verify the dioptric power and cylinder power of the lens.
3. Examine the lens thoroughly to ensure particles have not become attached to it, and examine the lens' optical surfaces for other defects.
4. If desired, the lens may be soaked or rinsed in sterile balanced salt solution until ready for implantation.
5. AMO recommends using the UNFOLDER Platinum 1-Series implantation system with the 1MTEC30 cartridge, the ONE SERIES Ultra implantation system with the 1VIPR30 cartridge, the UNFOLDER Emerald AR Series implantation system with the 1CART30 cartridge, or an equivalent insertion instrument or system to insert the TECNIS® Toric 1-piece lens. Only insertion instruments that have been validated or approved for use with this lens should be used. Please refer to the directions for use with the insertion instrument or system for additional information.

SELECTION AND PLACEMENT OF THE TECNIS® TORIC 1-PIECE IOL:

The astigmatism to be corrected should be determined from keratometry and biometry data rather than refractive data since the presence of lenticular astigmatism in the crystalline lens to be removed may influence results. The size and location of the surgical incision may affect the amount and axis of corneal astigmatism. In order to optimize IOL selection and axis placement, AMO provides a web-based tool (www.TecnisToricCalc.com) for the surgeon. Preoperative keratometry and biometry data, incision location, and the surgeon's estimated surgically induced corneal astigmatism are used to determine the appropriate TECNIS® Toric 1-piece IOL model, spherical equivalent lens power, and axis of placement in the eye.

For optimal results, the surgeon must ensure the correct placement and orientation of the lens within the capsular bag. The anterior surface of the IOL is marked with indentations (four at opposite sides) at the haptic/optic junction that identify the flat meridian of the TECNIS® Toric 1-piece optic. These indentations form an imaginary line representing the plus cylinder axis (note: IOL cylinder steep meridian is 90° away). The TECNIS® Toric 1-piece IOL cylinder axis marks should be aligned with the post-incision steep corneal meridian (intended axis of placement).

Prior to surgery the operative eye should be marked in the following manner:
With the patient sitting upright, precisely mark the twelve o'clock and/or the six o'clock position with a T marker, a surgical skin marker, or a marking pencil indicated for ophthalmic use. Using these marks as reference points, an axis marker can be used immediately prior to or during surgery to mark the axis of lens placement following the use of the web-based www.TecnisToricCalc.com to determine the optimal axis of placement.

After the lens is inserted, precisely align the axis marking indentations on the TECNIS® Toric 1-piece IOL with the marked axis of lens placement. Carefully remove all viscoelastic from the capsular bag. This may be accomplished by manipulating the IOL optic with the I/A tip and using standard irrigation/aspiration techniques to remove all viscoelastic from the eye. Bimanual techniques may be used, if preferred, to ensure removal of viscoelastic from behind the lens implant. Special care should be taken to ensure proper positioning of the TECNIS® Toric 1-piece IOL at the intended axis following viscoelastic removal. Residual viscoelastic may allow the lens to rotate, causing misalignment of the TECNIS® Toric 1-piece IOL with the intended axis of placement.

Misalignment of the axis of the lens with the intended axis of placement may compromise its astigmatic correction. Such misalignment can result from inaccurate keratometry or marking of the cornea, inaccurate placement of the TECNIS® Toric 1-piece IOL axis during surgery, an unanticipated surgically induced change in the cornea, or physical rotation of the TECNIS® Toric 1-piece IOL after implantation. In order to minimize this effect, the surgeon should be careful to ensure that preoperative keratometry and biometry is accurate and that the IOL is properly oriented prior to the end or surgery.

CAUTION: Do not use the lens if the package has been damaged. The sterility of the lens may have been compromised.

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HOW SUPPLIED:

The TECNIS® Toric 1-piece lenses are supplied sterile in a lens case within a double aseptic transfer peel pouch. The double aseptic transfer peel pouch is sterilized with ethylene oxide and should be opened only under sterile conditions. The pouch and product labels are enclosed in a shelf pack. The external surfaces of the outer pouch are not sterile.

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RETURN / EXCHANGE POLICY:

Contact the local AMO representative for the return policy. Return the lens with proper identification and the reason for the return. Label the return as a biohazard. Do not attempt to resterilize the lens.

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Symbol/Explanation:

SYMBOL	EXPLANATION
STERILE EO	Sterilized using Ethylene Oxide
	Do Not Reuse
	Use By (YYYY-MM: Year-Month)
	Caution: See Instructions for Use
	Manufacturer
EC REP	Authorized Representative in the European Union
	Do Not Resterilize
45°C (113°F)	Upper Limit of Temperature
	Keep Away from Sunlight
	Do Not Use if Package is Damaged

AMO, TECNIS, TRI-FIX, UNFOLDER, PROTEC and ONE SERIES are trademarks of Johnson & Johnson Surgical Vision, Inc.

Manufactured in the Netherlands:
AMO Groningen BV, 9728 NX Groningen, The Netherlands

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
Figure 1
Light Transmittance

Legend:
Curve 1: Spectral Transmittance curve of a typical 5 diopter IOL (thinnest). UV cut-off at 10% T is 375 nm
Curve 2: Spectral Transmittance curve of a typical 34 diopter IOL (thickest). UV cut-off at 10% T is 380 nm
Curve 3: Spectral Transmittance (T) Curve* Corresponding to 53 year-old Phakic Eye

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EN

Rx Only

DESCRIPTION:
The TECNIS® multifocal foldable acrylic 1-piece lens, Model ZKB00, is an ultraviolet light-absorbing posterior chamber intraocular lens (IOL). It is designed to be positioned in the lens capsule where the lens should replace the optical function of the natural crystalline lens. The TECNIS® multifocal foldable acrylic 1-piece lens incorporates a proprietary wavefront-designed aspheric optic with a squared posterior edge designed to provide a 360-degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects.

The lens incorporates a diffractive multifocal optic pattern designed to provide both near and distance vision and thereby reduce spectacle dependency. The light distribution between the distance and near focus is approximately 50/50. The labeled power of the lens is the distance power. The near power represents a +2.75 diopter add in actual lens power. However, accommodation will not be restored.

INDICATIONS FOR USE:
The TECNIS® multifocal foldable acrylic 1-piece intraocular lens is indicated for primary implantation for the visual correction of [1] aphakia in adults with and without presbyopia in whom a cataractous lens has been removed by phacoemulsification and [2] aphakia following refractive lensectomy in presbyopic adults who desire improved uncorrected near, intermediate, and distance vision and increased spectacle independence across a range of distances. The intraocular lens is intended to be placed in the capsular bag.

PRECAUTIONS:
1. Prior to surgery, the surgeon must inform prospective patients of the possible risks and benefits associated with the use of this device and provide a copy of the patient information brochure to the patient.
2. There were no patients 21 years old or younger included in the clinical study of parent lens TECNIS® Multifocal IOL, Model ZM900 (N=348); therefore there are insufficient clinical data to demonstrate safety and effectiveness in this age group.
3. The central one millimeter area of the TECNIS® Multifocal 1-piece IOL creates a far image focus in accordance with the labeled power of the IOL, so patients with abnormally small pupils (~1mm) should achieve, at a minimum, the prescribed distance vision under photopic conditions; however, because this multifocal design has not been tested in patients with abnormally small pupils, it is unclear whether such patients will derive any near vision benefit.
4. Autorefractors may not provide optimal postoperative refraction of patients with multifocal lenses. Manual refraction is strongly recommended.
5. Recent contact lens usage may affect the patient's refraction; therefore in contact lens wearers, surgeons should establish corneal stability without contact lenses prior to determining IOL power.
6. When performing wavefront measurements on a patient with a multifocal lens, two different wavefronts are produced. One wavefront will be in focus (either far or near) and the other wavefront will be out of focus. In this situation, incorrect interpretation of the wavefront measurements is possible.
7. The long-term effects of intraocular lens implantation have not been determined. Therefore the physician should continue to monitor implant patients postoperatively on a regular basis.
8. Secondary glaucoma has been reported occasionally in patients with controlled glaucoma who received lens implants. The intraocular pressure of implant patients with glaucoma should be carefully monitored postoperatively.
9. Do not resterilize this intraocular lens by any method.
10. Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.
11. Do not store the lens in direct sunlight or at a temperature greater than 45 °C (113 °F). Do not autoclave the intraocular lens.
12. Prior to implanting, examine the lens package for proper lens model, dioptric power, and expiration date.
13. The surgeon should target emmetropia as this lens is designed for optimum visual performance when emmetropia is achieved.
14. Care should be taken to achieve centration of the intraocular lens.
15. Please refer to the specific instructions for use provided with the insertion instrument or system for the amount of time the IOL can remain folded before the IOL must be discarded.
16. When the insertion system is used improperly, the haptics of the TECNIS® Multifocal 1-piece lens may become broken. Please refer to the specific instructions for use provided with the insertion instrument or system.

WARNINGS:
Physicians considering lens implantation under any of the following circumstances should weigh the potential risk/benefit ratio:
1. Some visual effects associated with multifocal IOLs may be expected because of the superposition of focused and unfocused images. These may include a perception of halos or glare around lights under nighttime conditions. It is expected that, in a small percentage of patients, the observation of such phenomena will be annoying and may be perceived as a hindrance, particularly in low illumination conditions. On rare occasions these visual effects may be significant enough that the patient will request removal of the multifocal IOL.
2. Under low-contrast conditions, contrast sensitivity is reduced with a multifocal lens compared to a monofocal lens. Therefore, subjects with multifocal lenses should exercise caution when driving at night or in poor visibility conditions.
3. Children under the age of 2 years are not suitable candidates for intraocular lenses.
4. Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the lens may exacerbate an existing condition, may interfere with diagnosis or treatment of a condition or may pose an unreasonable risk to the patient's eyesight:
 a. Patients with recurrent severe anterior or posterior segment inflammation or uveitis.
 b. Patients in whom the intraocular lens may interfere with the ability to observe, diagnose or treat posterior segment diseases.
 c. Surgical difficulties at the time of cataract extraction and/or intraocular lens implantation that might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
 d. A compromised eye due to previous trauma or developmental defects in which appropriate support of the IOL is not possible.
 e. Circumstance that would result in damage to the endothelium during implantation.
 f. Suspected microbial infection.
 g. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support for the IOL.
 h. Congenital bilateral cataracts.
 i. Recurrent severe anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye.
 j. Previous history of, or a predisposition to, retinal detachment.
 k. Patients with only one eye with potentially good vision.
 l. Medically uncontrollable glaucoma.
 m. Corneal endothelial dystrophy.
 n. Proliferative diabetic retinopathy.
5. The TECNIS® Multifocal 1-piece IOL should be placed entirely in the capsular bag. Do not place the lens in the ciliary sulcus.
6. The splitting of the light into more than one focus may affect image quality and lead to some reduction of contrast sensitivity.
7. Well-informed patients with well-defined visual needs and preferences should be selected for TECNIS® Multifocal 1-piece lens implantation. The patients should be informed about the possibility that a decrease in contrast sensitivity and an increase of visual disturbances may affect their ability to drive a car under certain environmental conditions, such as driving at night or in poor visibility conditions.
8. Patients with a predicted postoperative astigmatism greater than 1.0 diopter may not be suitable candidates for multifocal IOL implantation since they may not fully benefit from a multifocal IOL in terms of potential spectacle independence.
9. Care should be taken to achieve IOL centration, as lens decentration may result in patients experiencing visual disturbances, particularly in patients with large pupils under mesopic conditions.
10. AMO IOLs are single-use medical devices and are labeled with instructions for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The reuse/resterilization/reprocessing of AMO single-use medical devices may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

ADVERSE EVENTS:
Adverse events that occurred during the clinical trial of parent lens TECNIS® Multifocal IOL, Model ZM900 (N=348) include the following: macular edema, endophthalmitis and hypopyon (both followed by non lens- related surgical re-interventions), and surgical re-intervention. Only three subjects out of 348 experienced lens-related events (3/348; 0.9%) in which the IOLs were ultimately removed (two due to halos and one due to blurry/filmy vision). Only two secondary surgical re-intervention events for multifocal first eyes were lens-related (lens removal due to halos/glare and lens repositioning); the majority of secondary surgical re-intervention events were not lens-related (n=11, including iris prolapse/wound repair, lens exchange for refractive error and incorrect lens type, retinal repair, trabeculectomy and subsequent filtration bleb revisions, and treatment injections for cystoid macular edema). The only persistent adverse event was raised intraocular pressure (IOP) due to secondary glaucoma requiring treatment.

Adverse events that occurred during the clinical trial of parent lens SENSAR 1-piece IOL, Model AAB00, (N=123) include the following: macular edema, secondary surgical intervention (pars plana vitrectomy with membrane peel), and a torn haptic related to improper loading technique, resulting in lens exchange.

Additional complications that may occur following implantation of a multifocal intraocular lens include: objectionable visual quality under certain lighting conditions (e.g., chromatic aberrations, halos, night glare, decreased contrast), hazy/filmy images, "ghost" imaging of near and/or distant objects, diplopia, and poor stereopsis. Some of these optical effects may be mitigated upon patients' adaptation to the multifocal optic. However, some of the above listed complications may require secondary surgical intervention for intraocular lens exchange or explantation.

DETAILED DEVICE DESCRIPTION:
The TECNIS® Multifocal 1-piece lens, Model ZKB00, is a one-piece, foldable, posterior chamber lens made of soft acrylic material. The lens has a diffractive multifocal surface on the posterior side and a modified prolate (aspheric) surface on the anterior side of the lens. The optic is 6.0 mm in diameter and the lens has an overall diameter of 13.0 mm. The add power is +2.75 diopters, corresponding to +2.01 diopters in the spectacle plane.

Lens Optic:
1. Optic Material: Optically clear, soft foldable hydrophobic acrylic with a covalently bound UV absorber
2. Dioptric Power: +5.0 D to +34.0 D in 0.5 D increments
3. Optic Center Thickness: 0.7 mm (+20.0D)
4. Optic Edge Design: PROTEC 360 Square posterior edge
5. Index of Refraction: 1.47 at 35 °C
6. Light Transmittance: UV cut-off at 10%T for a +5.0 diopter lens (thinnest) and a +34.0 diopter lens (thickest) are shown in **Figure 1**.
7. Modulation Transfer Function (MTF) Through-Focus Response: MTF values are shown in **Figure 2**.

Haptics:
1. Material: Soft foldable acrylic with a covalently bound UV absorber
2. One-piece lens
3. Configuration: TRI-FIX design, Modified C, integral with optic
4. Haptic Thickness: 0.46 mm

DIRECTIONS FOR USE:
1. Prior to implanting, examine the lens package for proper lens model, dioptric power, and expiration date.
2. Open the peel pouches and remove the lens in a sterile environment. Verify the dioptric and near add power of the lens.
3. Examine the lens thoroughly to ensure particles have not become attached to it, and examine the lens optical surfaces for other defects.
4. If desired, the lens may be soaked in sterile balanced salt solution or sterile normal saline until ready for implantation.
5. Handle the lens by the haptic portion. Do not grasp the optical area with forceps.
6. Transfer the lens, using sterile technique, to an appropriate loading device.
7. The physician should consider the following points:
 - The surgeon should target emmetropia as this lens is designed for optimum visual performance when emmetropia is achieved.
 - Care should be taken to achieve centration of the intraocular lens.
8. AMO recommends using the UNFOLDER Platinum 1-Series Implantation system with 1MTEC30 cartridge, the ONE SERIES Ultra Implantation system with the 1VIPR30 cartridge, the UNFOLDER Emerald AR Series Implantation system with the 1CART30 cartridge, or an equivalent qualified insertion instrument or system to insert the TECNIS® Multifocal 1-piece lens. Only insertion instruments that have been validated and approved for use with this lens should be used. Please refer to the directions for use with the insertion instrument or system for additional information.

CAUTION:
Do not use the lens if the package has been damaged. The sterility of the lens may have been compromised.

LENS POWER CALCULATIONS:
Accurate keratometry and biometry are essential to successful visual outcomes. Preoperative calculation of the required dioptric lens power for these posterior chamber intraocular lenses should be determined by the surgeon's experience, preference, and intended lens placement. **Emmetropia should be targeted.** The appropriate A-constants listed on the outer label are presented as a guideline and is a starting point for implant power calculations. The physician should determine preoperatively the dioptric power of the lens to be implanted.

Physicians requiring additional information on lens power calculations may contact the local Johnson & Johnson Surgical Vision, Inc. representative. Lens power calculation methods are described in the following references:

- Holladay JT, Musgrove KH, Prager TC, Lewis JW, Chandler TY and Ruiz RS. A three-part system for refining intraocular lens power calculations. J Cataract Refract Surg. 14:17-24, 1988.
- Retzlaff JA, Sanders DR and Kraff MC. Development of the SRK/T intraocular lens implant power calculation formula. J Cataract Refract Surg. 16:333-340, 1990; ERRATA, 16:528, 1990.
- Olsen T, Olesen H, Thim K and Corydon L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. J Cataract Refract Surg. 18:280-285, 1992.
- Hoffer KJ. The Hoffer Q formula: A comparison of theoretic and regression formulas. J Cataract Refract Surg. 19:700-712, 1993; ERRATA, 20:677, 1994.
- Holladay JT. Standardizing constants for ultrasonic biometry, keratometry and intraocular lens power calculations. J Cataract Refract Surg. 23:1356-1370, 1997.
- Norrby NES. Unfortunate discrepancies. Letter to the editor and reply by Holladay JT. J Cataract Refract Surg. 24:433-434, 1998.
- Norrby S, Lydahl E, Koranyi G, Taube M. Reduction of trend errors in power calculation by linear transformation of measured axial lengths. J Cataract Refract Surg. 2003; 29:100-105.
- http://ocusoft.de/ulib/c1.htm is in particular useful for Zeiss IOLMaster users.

PATIENT CARD:
An implant identification card, to be supplied to the patient, is included in the package. The patient should be instructed to keep the card as a permanent record of his/her implant and to show the card to any eye care practitioner he/she may see in the future.

REPORTING:
Adverse events and/or potentially sight-threatening complications that may reasonably be regarded as lens-related and that were not previously expected in nature, severity or rate of occurrence must be reported to AMO. This information is being requested from all surgeons in order to document potential long-term effects of intraocular lens implantation, especially in younger patients.

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HOW SUPPLIED:
Each TECNIS® Multifocal 1-piece lens is supplied sterile, in a lens case within a double aseptic transfer peel pouch. The double aseptic transfer peel pouch is sterilized with ethylene oxide and should be opened only under sterile conditions. The pouch and product labels are enclosed in a shelf pack. The external surfaces of the outer pouch are not sterile.

EXPIRATION DATE:
The expiration date on the lens package is the sterility expiration date. The lens should not be implanted after the indicated sterility expiration date.

RETURN/EXCHANGE POLICY:
Contact the local AMO representative for the return lens policy. Return the lens with proper identification and the reason for the return. Label the return as a biohazard. Do not attempt to resterilize the lens.

PATIENT INFORMATION
Each patient should receive information regarding intraocular lenses prior to the decision to implant an intraocular lens.

SYMBOL/EXPLANATION

SYMBOL	EXPLANATION
STERILE EO	Sterilized Using Ethylene Oxide
	Do Not Reuse
	Use By (YYYY-MM-DD: year-month-day)
	Consult Instructions for Use
	Manufacturer
EC REP	Authorized Representative in the European Union
	Do Not Resterilize
45 °C (113 °F)	Upper Limit of Temperature
	Do Not Use if Package is Damaged
	Keep Away from Sunlight
	Date of Manufacture (YYYY-MM-DD: year-month-day)

Product of the Netherlands
AMO Groningen BV
Van Swietenlaan 5
9728 NX Groningen
The Netherlands

Johnson & Johnson Surgical Vision, Inc.
1700 E. St. Andrew Place
Santa Ana, CA 92705 USA

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The CE marked IOLs comply with the European Council Directive 93/42/EEC of June 14, 1993.

AMO, TECNIS, SENSAR, TRI-FIX, UNFOLDER, ONE SERIES and PROTEC are trademarks of Johnson & Johnson Surgical Vision, Inc.


MUKEESH SABARAD
FEB 13TH, 2020

Regulatory Affairs
Surgical Vision Inc.

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Figure 1
Light Transmittance curves


LEGEND:
Curve 1: Spectral Transmittance curve of a typical 5 diopter IOL (thinnest), UV cut-off at 10% T is 375nm.
Curve 2: Spectral Transmittance curve of a typical 34 diopter IOL (thickest), UV cut-off at 10% T is 380nm.
Curve 3: Spectral Transmittance (T) Curve* Corresponding to 53 year-old Phakic Eye

Note: ***The cut-off wavelengths and the spectral transmittance curves represent the range of the transmittance of IOLs (5-34 diopter) made with this material. Spectral transmission measurements were taken in water at room temperature.***

*** Boettner, E.A., and Wolter J.R. Transmission of the Ocular Media. Investigative Ophthalmology. 1962; 1:776-783.**

Figure 2
ACE Model – MTF Through-Focus for a 3mm Pupil, 20-Diopter Lenses

Note: The figure above describes the optical performance of the TECNIS® Multifocal 1-piece lens, Model ZKB00, in ACE* model eyes at 50 cycles/mm using green light. The graph depicts the focus as it is gradually shifted from that of a far object to increasingly nearer objects, with higher numbers typically indicating better performance. The ACE* model has the spherical aberration of an average natural cornea, which the TECNIS® lens is designed to compensate. The combination is aberration-free.

* Norrby S, Piers P, Campbell C, van der Mooren M. "Model eyes for evaluation of intraocular lenses." Appl Opt. 2007 Sep 10;46(26):6595-605.

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- Holladay JT, Musgrove KH, Prager TC, Lewis JW, Chandler TY and Ruiz RS. A three-part system for refining intraocular lens power calculations. J Cataract Refract Surg. 14:17-24, 1988.
- Retzlaff JA, Sanders DR and Kraff MC. Development of the SRK/T intraocular lens implant power calculation formula. J Cataract Refract Surg. 16:333-340, 1990; ERRATA, 16:528, 1990.
- Olsen T, Olesen H, Thim K and Corydon L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. J Cataract Refract Surg. 18:280-285, 1992.
- Hoffer KJ. The Hoffer Q formula: A comparison of theoretic and regression formulas. J Cataract Refract Surg. 19:700-712, 1993; ERRATA, 20:677, 1994.
- Holladay JT. Standardizing constants for ultrasonic biometry, keratometry and intraocular lens power calculations. J Cataract Refract Surg. 23:1356-1370, 1997.
- Norrby NES. Unfortunate discrepancies. Letter to the editor and reply by Holladay JT. J Cataract Refract Surg. 24:433-434, 1998.
- Norrby S, Lydahl E, Koranyi G, Taube M. Reduction of trend errors in power calculation by linear transformation of measured axial lengths. J Cataract Refract Surg. 2003; 29:100-105.
- http://ocusoft.de/ulib/c1.htm — ресурс, особенно полезный для пользователей системы Zeiss IOLMaster.

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Rx Only

DESCRIPTION:

The TECNIS® multifocal foldable acrylic 1-piece lens, Model ZLB00, is an ultraviolet light-absorbing posterior chamber intraocular lens (IOL). It is designed to be positioned in the lens capsule where the lens should replace the optical function of the natural crystalline lens. The TECNIS® multifocal foldable acrylic 1-piece lens incorporates a proprietary wavefront-designed aspheric optic with a squared posterior edge designed to provide a 360-degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects.

The lens incorporates a diffractive multifocal optic pattern designed to provide both near and distance vision and thereby reduce spectacle dependency. The light distribution between the distance and near focus is approximately 50/50. The labeled power of the lens is the distance power. The near power represents a +3.25 diopter add in actual lens power. However, accommodation will not be restored.

INDICATIONS FOR USE:

The TECNIS® multifocal foldable acrylic 1-piece intraocular lens is indicated for primary implantation for the visual correction of [1] aphakia in adults with and without presbyopia in whom a cataractous lens has been removed by phacoemulsification and [2] aphakia following refractive lensectomy in presbyopic adults who desire improved uncorrected near, intermediate, and distance vision and increased spectacle independence across a range of distances. The intraocular lens is intended to be placed in the capsular bag.

PRECAUTIONS:

1. Prior to surgery, the surgeon must inform prospective patients of the possible risks and benefits associated with the use of this device and provide a copy of the patient information brochure to the patient.
2. There were no patients 21 years old or younger included in the clinical study of parent lens TECNIS® Multifocal IOL, Model ZM900 (N=348); therefore there are insufficient clinical data to demonstrate safety and effectiveness in this age group.
3. The central one millimeter area of the TECNIS® Multifocal 1-piece IOL creates a far image focus in accordance with the labeled power of the IOL, so patients with abnormally small pupils (~1mm) should achieve, at a minimum, the prescribed distance vision under photopic conditions; however, because this multifocal design has not been tested in patients with abnormally small pupils, it is unclear whether such patients will derive any near vision benefit.
4. Autorefractors may not provide optimal postoperative refraction of patients with multifocal lenses. Manual refraction is strongly recommended.
5. Recent contact lens usage may affect the patient's refraction; therefore in contact lens wearers, surgeons should establish corneal stability without contact lenses prior to determining IOL power.
6. When performing wavefront measurements on a patient with a multifocal lens, two different wavefronts are produced. One wavefront will be in focus (either far or near) and the other wavefront will be out of focus. In this situation, incorrect interpretation of the wavefront measurements is possible.
7. The long-term effects of intraocular lens implantation have not been determined. Therefore the physician should continue to monitor implant patients postoperatively on a regular basis.
8. Secondary glaucoma has been reported occasionally in patients with controlled glaucoma who received lens implants. The intraocular pressure of implant patients with glaucoma should be carefully monitored postoperatively.
9. Do not resterilize this intraocular lens by any method.
10. Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.
11. Do not store the lens in direct sunlight or at a temperature greater than 45 °C (113 °F). Do not autoclave the intraocular lens.
12. Prior to implanting, examine the lens package for proper lens model, dioptric power, and expiration date.
13. The surgeon should target emmetropia as this lens is designed for optimum visual performance when emmetropia is achieved.
14. Care should be taken to achieve centration of the intraocular lens.
15. Please refer to the specific instructions for use provided with the insertion instrument or system for the amount of time the IOL can remain folded before the IOL must be discarded.
16. When the insertion system is used improperly, the haptics of the TECNIS® Multifocal 1-piece lens may become broken. Please refer to the specific instructions for use provided with the insertion instrument or system.

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1. Some visual effects associated with multifocal IOLs may be expected because of the superposition of focused and unfocused images. These may include a perception of halos or glare around lights under nighttime conditions. It is expected that, in a small percentage of patients, the observation of such phenomena will be annoying and may be perceived as a hindrance, particularly in low illumination conditions. On rare occasions these visual effects may be significant enough that the patient will request removal of the multifocal IOL.
2. Under low-contrast conditions, contrast sensitivity is reduced with a multifocal lens compared to a monofocal lens. Therefore, subjects with multifocal lenses should exercise caution when driving at night or in poor visibility conditions.
3. Children under the age of 2 years are not suitable candidates for intraocular lenses.
4. Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the lens may exacerbate an existing condition, may interfere with diagnosis or treatment of a condition or may pose an unreasonable risk to the patient's eyesight:
 a. Patients with recurrent severe anterior or posterior segment inflammation or uveitis.
 b. Patients in whom the intraocular lens may interfere with the ability to observe, diagnose or treat posterior segment diseases.
 c. Surgical difficulties at the time of cataract extraction and/or intraocular lens implantation that might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
 d. A compromised eye due to previous trauma or developmental defects in which appropriate support of the IOL is not possible.
 e. Circumstance that would result in damage to the endothelium during implantation.
 f. Suspected microbial infection.
 g. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support for the IOL.
 h. Congenital bilateral cataracts.
 i. Recurrent severe anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye.
 j. Previous history of, or a predisposition to, retinal detachment.
 k. Patients with only one eye with potentially good vision.
 l. Medically uncontrollable glaucoma.
 m. Corneal endothelial dystrophy.
 n. Proliferative diabetic retinopathy.
5. The TECNIS® Multifocal 1-piece IOL should be placed entirely in the capsular bag. Do not place the lens in the ciliary sulcus.
6. The splitting of the light into more than one focus may affect image quality and lead to some reduction of contrast sensitivity.
7. Well-informed patients with well-defined visual needs and preferences should be selected for TECNIS® Multifocal 1-piece lens implantation. The patients should be informed about the possibility that a decrease in contrast sensitivity and an increase of visual disturbances may affect their ability to drive a car under certain environmental conditions, such as driving at night or in poor visibility conditions.
8. Patients with a predicted postoperative astigmatism greater than 1.0 diopter may not be suitable candidates for multifocal IOL implantation since they may not fully benefit from a multifocal IOL in terms of potential spectacle independence.
9. Care should be taken to achieve IOL centration, as lens decentration may result in patients experiencing visual disturbances, particularly in patients with large pupils under mesopic conditions.
10. AMO IOLs are single-use medical devices and are labeled with instructions for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The reuse/resterilization/reprocessing of AMO single-use medical devices may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

ADVERSE EVENTS:

Adverse events that occurred during the clinical trial of parent lens TECNIS® Multifocal IOL, Model ZM900 (N=348) include the following: macular edema, endophthalmitis and hypopyon (both followed by non lens- related surgical re-interventions), and surgical re-intervention. Only three subjects out of 348 experienced lens-related events (3/348; 0.9%) in which the IOLs were ultimately removed (two due to halos and one due to blurry/filmy vision). Only two secondary surgical re-intervention events for multifocal first eyes were lens-related (lens removal due to halos/glare and lens repositioning); the majority of secondary surgical re-intervention events were not lens-related (n=11, including iris prolapse/wound repair, lens exchange for refractive error and incorrect lens type, retinal repair, trabeculectomy and subsequent filtration bleb revisions, and treatment injections for cystoid macular edema). The only persistent adverse event was raised intraocular pressure (IOP) due to secondary glaucoma requiring treatment.

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DETAILED DEVICE DESCRIPTION:

The TECNIS® Multifocal 1-piece lens, Model ZLB00, is a one-piece, foldable, posterior chamber lens made of soft acrylic material. The lens has a diffractive multifocal surface on the posterior side and a modified prolate (aspheric) surface on the anterior side of the lens. The optic is 6.0 mm in diameter and the lens has an overall diameter of 13.0 mm. The add power is +3.25 diopters, corresponding to +2.37 diopters in the spectacle plane.

Lens Optic:

1. Optic Material: Optically clear, soft foldable hydrophobic acrylic with a covalently bound UV absorber
2. Dioptric Power: +5.0 D to +34.0 D in 0.5 D increments
3. Optic Center Thickness: 0.7 mm (+20.0D)
4. Optic Edge Design: PROTEC 360 Square posterior edge
5. Index of Refraction: 1.47 at 35 °C
6. Light Transmittance: UV cut-off at 10%T for a +5.0 diopter lens (thinnest) and a +34.0 diopter lens (thickest) are shown in **Figure 1**.
7. Modulation Transfer Function (MTF) Through-Focus Response: MTF values are shown in **Figure 2**.

Haptics:

1. Material: Soft foldable acrylic with a covalently bound UV absorber
2. One-piece lens
3. Configuration: TRI-FIX design, Modified C, integral with optic
4. Haptic Thickness: 0.46 mm

DIRECTIONS FOR USE:

1. Prior to implanting, examine the lens package for proper lens model, dioptric power, and expiration date.
2. Open the peel pouches and remove the lens in a sterile environment. Verify the dioptric and near add power of the lens.
3. Examine the lens thoroughly to ensure particles have not become attached to it, and examine the lens optical surfaces for other defects.
4. If desired, the lens may be soaked in sterile balanced salt solution or sterile normal saline until ready for implantation.
5. Handle the lens by the haptic portion. Do not grasp the optical area with forceps.
6. Transfer the lens, using sterile technique, to an appropriate loading device.
7. The physician should consider the following points:
 - The surgeon should target emmetropia as this lens is designed for optimum visual performance when emmetropia is achieved.
 - Care should be taken to achieve centration of the intraocular lens.
8. AMO recommends using the UNFOLDER Platinum 1-Series Implantation system with 1MTEC30 cartridge, the ONE SERIES Ultra Implantation system with the 1VIPR30 cartridge, the UNFOLDER Emerald AR Series Implantation system with the 1CART30 cartridge, or an equivalent qualified insertion instrument or system to insert the TECNIS® Multifocal 1-piece lens. Only insertion instruments that have been validated and approved for use with this lens should be used. Please refer to the directions for use with the insertion instrument or system for additional information.

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LENS POWER CALCULATIONS:

Accurate keratometry and biometry are essential to successful visual outcomes. Preoperative calculation of the required dioptric lens power for these posterior chamber intraocular lenses should be determined by the surgeon's experience, preference, and intended lens placement. **Emmetropia should be targeted.** The appropriate A-constants listed on the outer label are presented as a guideline and is a starting point for implant power calculations. The physician should determine preoperatively the dioptric power of the lens to be implanted.

Physicians requiring additional information on lens power calculations may contact the local Johnson & Johnson Surgical Vision, Inc. (AMO) representative. Lens power calculation methods are described in the following references:

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REPORTING:

Adverse events and/or potentially sight-threatening complications that may reasonably be regarded as lens-related and that were not previously expected in nature, severity or rate of occurrence must be reported to AMO. This information is being requested from all surgeons in order to document potential long-term effects of intraocular lens implantation, especially in younger patients.

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HOW SUPPLIED:

Each TECNIS® Multifocal 1-piece lens is supplied sterile, in a lens case within a double aseptic transfer peel pouch. The double aseptic transfer peel pouch is sterilized with ethylene oxide and should be opened only under sterile conditions. The pouch and product labels are enclosed in a shelf pack. The external surfaces of the outer pouch are not sterile.

EXPIRATION DATE:

The expiration date on the lens package is the sterility expiration date. The lens should not be implanted after the indicated sterility expiration date.

RETURN/EXCHANGE POLICY:

Contact the local AMO representative for the return lens policy. Return the lens with proper identification and the reason for the return. Label the return as a biohazard. Do not attempt to resterilize the lens.

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SYMBOL/EXPLANATION

SYMBOL	EXPLANATION
STERILE EO	Sterilized Using Ethylene Oxide
	Do Not Reuse
	Use By (YYYY-MM-DD: year-month-day)
	Consult Instructions for Use
	Manufacturer
EC REP	Authorized Representative in the European Union
	Do Not Resterilize
45°C (113°F)	Upper Limit of Temperature
	Do Not Use if Package is Damaged
	Keep Away from Sunlight
	Date of Manufacture (YYYY-MM-DD: year-month-day)

Product of the Netherlands
AMO Groningen BV
Van Swietenlaan 5
9728 NX Groningen
The Netherlands

Johnson & Johnson
Surgical Vision, Inc.
1700 E. St. Andrew Place
Santa Ana, CA 92705 USA

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

Note: The figure above describes the optical performance of the TECNIS® Multifocal 1-piece lens, Model ZLB00, in ACE* model eyes at 50 cycles/mm using green light. The graph depicts the focus as it is gradually shifted from that of a far object to increasingly nearer objects, with higher numbers typically indicating better performance. The ACE* model has the spherical aberration of an average natural cornea, which the TECNIS® lens is designed to compensate. The combination is aberration-free.

* Norrby S, Piers P, Campbell C, van der Mooren M. "Model eyes for evaluation of intraocular lenses." Appl Opt. 2007 Sep 10;46(26):6595-605.

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The TECNIS® Multifocal 1-Piece IOL

DESCRIPTION:

The TECNIS® multifocal foldable acrylic 1-piece lens, Model ZMB00, is an ultraviolet light-absorbing posterior chamber intraocular lens (IOL). It is designed to be positioned posterior to the iris where the lens should replace the optical function of the natural crystalline lens. The TECNIS® multifocal foldable acrylic 1-piece lens incorporates a proprietary wavefront-designed aspheric optic with a squared posterior edge designed to provide a 360-degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects.

The lens also incorporates a diffractive multifocal optic pattern designed to provide both near and distance vision and thereby reduce spectacle dependency. The light distribution between the distance and near focus is approximately 50/50. The labeled power of the lens is the distance power. The near power represents a +4 diopter add in actual lens power. However, accommodation will not be restored.

INDICATIONS FOR USE:

TECNIS® multifocal foldable acrylic 1-piece intraocular lenses are indicated for primary implantation for the visual correction of: [1] aphakia in adults in whom a cataractous lens has been removed, and [2] aphakia following refractive lensectomy in presbyopic adults who may benefit from useful near vision and reduced spectacle dependence across a range of distances. The lens is intended to be placed in the capsular bag.

PRECAUTIONS:

1. Prior to surgery, the surgeon must inform prospective patients of the possible risks and benefits associated with the use of this device and provide a copy of the patient information brochure to the patient.
2. There were no patients 21 years old or younger included in the clinical study; therefore there are insufficient clinical data to demonstrate safety and effectiveness in this age group.
3. The central one millimeter area of the TECNIS® Multifocal 1-Piece IOL creates a far image focus in accordance with the labeled power of the IOL, so patients with abnormally small pupils (~1mm) should achieve, at a minimum, the prescribed distance vision under photopic conditions; however, because this multifocal design has not been tested in patients with abnormally small pupils, it is unclear whether such patients will derive any near vision benefit.
4. Autorefractors may not provide optimal postoperative refraction of patients with multifocal lenses. Manual refraction is strongly recommended.
5. Recent contact lens usage may affect the patient's refraction; therefore in contact lens wearers, surgeons should establish corneal stability without contact lenses prior to determining IOL power.
6. When performing wavefront measurements on a patient with a multifocal lens, two different wavefronts are produced. One wavefront will be in focus (either far or near) and the other wavefront will be out of focus. In this situation, incorrect interpretation of the wavefront measurements is possible.
7. The long-term effects of intraocular lens implantation have not been determined. Therefore the physician should continue to monitor implant patients postoperatively on a regular basis.
8. Secondary glaucoma has been reported occasionally in patients with controlled glaucoma who received lens implants. The intraocular pressure of implant patients with glaucoma should be carefully monitored postoperatively.
9. Do not resterilize this intraocular lens by any method.
10. Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.
11. Do not store the lens in direct sunlight or at a temperature greater than 45°C (113°F). Do not autoclave the intraocular lens.
12. Prior to implanting, examine the lens package for proper lens model, dioptric power, and expiration date.
13. The surgeon should target emmetropia as this lens is designed for optimum visual performance when emmetropia is achieved.
14. Care should be taken to achieve centration of the intraocular lens.
15. Please refer to the specific instructions for use provided with the insertion instrument or system for the amount of time the IOL can remain folded before the IOL must be discarded. When the insertion system is used improperly, the haptics of the TECNIS® Multifocal 1-Piece lens may become broken. Please refer to the specific instructions for use provided with the insertion instrument or system.

WARNINGS:

1. Some visual effects associated with multifocal IOLs may be expected because of the superposition of focused and unfocused images. These may include a perception of halos or glare around lights under nighttime conditions. It is expected that, in a small percentage of patients, the observation of such phenomena will be annoying and may be perceived as a hindrance, particularly in low illumination conditions. On rare occasions these visual effects may be significant enough that the patient will request removal of the multifocal IOL.
2. Under low-contrast conditions, contrast sensitivity is reduced with a multifocal lens compared to a monofocal lens. Therefore, subjects with multifocal lenses should exercise caution when driving at night or in poor visibility conditions.
3. Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the lens may exacerbate an existing condition, may interfere with diagnosis or treatment of a condition or may pose an unreasonable risk to the patient's eyesight:
 a. Patients in whom the intraocular lens may interfere with the ability to observe, diagnose or treat posterior segment diseases.
 b. Surgical difficulties at the time of cataract extraction and/or intraocular lens implantation that might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
 c. A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
 d. Circumstance that would result in damage to the endothelium during implantation.
 e. Suspected microbial infection.
 f. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support.
 g. Congenital bilateral cataracts.
 h. Recurrent severe anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye.
 i. Previous history of, or a predisposition to, retinal detachment.
 j. Patients with only one eye with potentially good vision.
 k. Medically uncontrollable glaucoma.
 l. Corneal endothelial dystrophy.
 m. Proliferative diabetic retinopathy.
4. The TECNIS® Multifocal 1-Piece IOL should be placed entirely in the capsular bag. Do not place the lens in the ciliary sulcus.
5. The splitting of the light into more than one focus may affect image quality and lead to some reduction of contrast sensitivity.
6. Well-informed patients with well-defined visual needs and preferences should be selected for TECNIS® Multifocal 1-Piece lens implantation. The patients should be informed about the possibility that a decrease in contrast sensitivity and an increase of visual disturbances may affect their ability to drive a car under certain environmental conditions, such as driving at night or in poor visibility conditions.
7. Patients with a predicted postoperative astigmatism greater than 1.0 diopter may not be suitable candidates for multifocal IOL implantation since they may not fully benefit from a multifocal IOL in terms of potential spectacle independence.
8. AMO IOLs are single-use medical devices that are labeled with instructions for use and handling to minimize exposure to conditions that may compromise the product, patient, or the user. The reuse/resterilization/reprocessing of AMO single-use medical devices may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

ADVERSE EVENTS:

Adverse events that occurred during the clinical trial of parent lens TECNIS® Multifocal IOL, Model ZM900 (N=348) include the following: macular edema, endophthalmitis and hypopyon (both followed by non-lens-related surgical re-interventions), and surgical re-intervention. Only three subjects out of 348 experienced lens-related events (3/348; 0.9%) in which the IOLs were ultimately removed (two due to halos and one due to blurry/filmy vision). Only two secondary surgical re-intervention events for multifocal first eyes were lens-related (lens removal due to halos/glare and lens repositioning); the majority of secondary surgical re-intervention events were not lens-related (n=11, including iris prolapse/wound repair, lens exchange for refractive error and incorrect lens type, retinal repair, trabeculectomy and subsequent filtration bleb revisions, and treatment injections for cystoid macular edema). The only persistent adverse event was raised intraocular pressure (IOP) due to secondary glaucoma requiring treatment.

Adverse events that occurred during the clinical trial of parent lens SENSAR® 1-Piece IOL, Model AAB00, (N=123) include the following: macular edema, secondary surgical intervention (pars plana vitrectomy with membrane peel), and a torn haptic related to improper loading technique, resulting in lens exchange.

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DETAILED DEVICE DESCRIPTION:

The TECNIS® Multifocal 1-piece lens, Model ZMB00, is a one-piece foldable posterior chamber lens. The optic and haptics are made of soft acrylic material. The lens has a diffractive multifocal surface on the posterior side of the lens and a modified prolate (aspheric) surface on the anterior side. The optic is 6.0 mm in diameter and the lens has an overall diameter of 13.0 mm. The add power is +4 diopters, corresponding to +3 diopters in the spectacle plane.

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Lens Optic:

1. Optic Material: Optically clear, soft foldable hydrophobic acrylic with a covalently bound UV absorber
2. Power: +5.0 to +34.0 diopter powers in 0.5 diopter increments
3. Optic Center Thickness: 0.722 mm (+20.0D)
4. Optic Edge Design: PROTEC 360 Square posterior edge
5. Index of Refraction: 1.47 at 35°C
6. Light Transmittance: UV cut-off at 10%T for a +5.0 diopter lens (thinnest) and a +34.0 diopter lens (thickest) are shown in **Figure 1**.
7. Modulation Transfer Function (MTF) Through-Focus Response: MTF values are shown in **Figure 2**.

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DIRECTIONS FOR USE:

1. Prior to implanting, examine the lens package for proper lens model, dioptric power, and expiration date.
2. Open the package and remove the lens in a sterile environment.
3. Examine the lens thoroughly to ensure particles have not become attached to it, and examine the lens optical surfaces for other defects.
4. The lens may be soaked in sterile balanced salt solution or sterile normal saline until ready for implantation.
5. Handle the lens by the haptic portion. Do not grasp the optical area with forceps.
6. Transfer the lens, using sterile technique, to an appropriate loading device.
7. The physician should consider the following points:
 - The surgeon should target emmetropia as this lens is designed for optimum visual performance when emmetropia is achieved.
 - Care should be taken to achieve centration of the intraocular lens.
8. AMO recommends using Duckworth & Kent 7722 folding forceps with the 7741 insertion forceps or an equivalent insertion instrument or system to insert the TECNIS® Multifocal 1-Piece lens. Only insertion instruments that have been validated and approved for use with this lens should be used. Please refer to the directions for use with the insertion instrument or system for additional information.

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LENS POWER CALCULATIONS:

The physician should determine preoperatively the power of the lens to be implanted. **Emmetropia should be targeted.** The estimated A-constant for this lens is provided on the lens box. Accuracy of IOL power calculation is particularly important with multifocal IOLs as spectacle independence is the goal of multifocal IOL implantation.

Physicians requiring additional information on lens power calculations may contact the local AMO representative. Lens power calculation methods are described in the following references:

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HOW SUPPLIED:

Each TECNIS® Multifocal 1-Piece lens is supplied sterile, in a lens case within a double aseptic transfer peel pouch. The double aseptic transfer peel pouch is sterilized with ethylene oxide and should be opened only under sterile conditions.

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Symbol/Explanation

SYMBOL	EXPLANATION
STERILE EO	Sterilized using Ethylene Oxide
	Do Not Reuse
	Use By (YYYY-MM-DD: year-month-day)
	Consult Instructions for Use
	Manufacturer
EC REP	Authorized Representative in the European Union
	Do Not Resterilize
	Do Not Use if Package is Damaged
	Keep Away from Sunlight
45°C (113°F)	Upper Limit of Temperature
	Date of Manufacture (YYYY-MM-DD: year-month-day)

Product of the Netherlands
AMO Groningen BV
Van Swietenlaan 5
9728 NX Groningen
The Netherlands

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
TECNIS, SENSAR, TRI-FIX and PROTEC are trademarks of Johnson & Johnson Surgical Vision, Inc.

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Figure 1
Light Transmittance



Legend:
Curve 1: Spectral Transmittance curve of a typical 5 diopter IOL (thinnest), UV cut-off at 10% T is 375 nm
Curve 2: Spectral Transmittance curve of a typical 34 diopter IOL (thickest), UV cut-off at 10% T is 380 nm
Curve 3: Spectral Transmittance (T) Curve* Corresponding to 53 year-old Phakic Eye

Note: The cut-off wavelengths and the spectral transmittance curves represent the range of the transmittance of IOLs (5-34 diopter) made with this material. Spectral transmission measurements were taken in water (aqueous) at room temperature.

*Boettner, E.A., and Wolter J.R. Transmission of the Ocular Media. Investigative Ophthalmology. 1962; 1:776-783.

Figure 2
MTF Through-Focus Response for Multiple Pupil Sizes, 20-Diopter Lens

Note: The figures above describe the optical performance of the TECNIS® Multifocal 1-piece lens, Model ZMB00, in both the ISO and ACE* model eyes at 50 cycles/mm as the focus is gradually shifted from that of a far object to increasingly nearer objects, with higher numbers typically indicating better performance. However, this may not be true for an aberration-correcting IOL such as the TECNIS® lens.

The ISO eye model has an aberration-free cornea, so what is measured is the aberration of the lens in the converging light from the cornea. As a result, MTF will vary with pupil sizes. Without aberration, MTF is independent of pupil size. The natural cornea is not aberration-free. The ACE* eye model has the spherical aberration of an average natural cornea, which the TECNIS® lens is designed to compensate. The combination is aberration-free.

*Norrby S, Piers P, Campbell C, van der Mooren M. "Model eyes for evaluation of intraocular lenses." Appl Opt. 2007 Sep 10;46(26):6595-605.

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- Holladay JT, Musgrove KH, Prager TC, Lewis JW, Chandler TY and Ruiz RS. A three-part system for refining intraocular lens power calculations. J Cataract Refract Surg. 14:17-24, 1988.
- Retzlaff JA, Sanders DR and Kraff MC. Development of the SRK/T intraocular lens implant power calculation formula. J Cataract Refract Surg. 16:333-340, 1990; ERRATA, 16:528, 1990.
- Olsen T, Olesen H, Thim K and Corydon L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. J Cataract Refract Surg. 18:280-285, 1992.
- Hoffer KJ. The Hoffer Q formula: A comparison of theoretic and regression formulas. J Cataract Refract Surg. 19:700-712, 1993; ERRATA, 20:677, 1994.
- Holladay JT. Standardizing constants for ultrasonic biometry, keratometry and intraocular lens power calculations. J Cataract Refract Surg. 23:1356-1370, 1997.
- Norrby NES. Unfortunate discrepancies. Letter to the editor and reply by Holladay JT. J Cataract Refract Surg. 24:433-434, 1998.
- Norrby S, Lydahl E, Koranyi G, Taube M. Reduction of trend errors in power calculation by linear transformation of measured axial lengths. J Cataract Refract Surg. 2003; 29:100-105.
- http://ocusoft.de/ulib/c1.htm — ресурс, особенно полезный для пользователей диагностической системы Zeiss IOLMaster.

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The TECNIS® Multifocal Toric 1-Piece IOL

Device Description

The TECNIS® Multifocal Toric 1-piece lens, Model ZMT, is an ultraviolet-light absorbing posterior chamber intraocular lens (IOL) which compensates for corneal spherical aberrations and corneal astigmatism and is intended to provide both near and distance vision and reduce spectacle dependence. It is designed to be positioned in the lens capsule where the lens should replace the optical function of the natural crystalline lens.

The TECNIS® Multifocal Toric 1-piece IOL incorporates a proprietary wavefront-designed toric aspheric optic with a squared posterior optic edge designed to provide a 360 degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects. The anterior optic compensates for corneal spherical aberrations and corneal astigmatism; the anteriorly located cylinder axis marks denote the meridian with the lowest power. On the posterior side the lens incorporates a diffractive multifocal surface designed to provide both near and distance vision and thereby reduce spectacle dependency. The light distribution between the distance and near focus is approximately 50/50. The labeled power of the lens is the distance spherical equivalent power. The near power represents a +4 diopter add in actual lens power. However, accommodation will not be restored.

Indications For Use

The TECNIS® Multifocal Toric 1-piece lens is indicated for the primary implantation for the visual correction of aphakia and pre-existing corneal astigmatism in [1] astigmatic adult patients with and without presbyopia in whom a cataractous lens has been removed by phacoemulsification and [2] aphakia following refractive lensectomy in astigmatic presbyopic adults, who desire improved uncorrected vision, reduction of residual refractive cylinder, useful near vision and reduced spectacle dependence across a range of distances. The intraocular lenses are intended to maintain rotational stability after implantation in the capsular bag.

Precautions

1. Prior to surgery, the surgeon must inform prospective patients of the possible risks and benefits associated with the use of this device and provide a copy of the patient information brochure to the patient.
2. There were no patients 21 years old or younger included in the clinical studies; therefore there are insufficient clinical data to demonstrate safety and effectiveness in this age group.
3. The central one millimeter area of the TECNIS® Multifocal Toric IOL creates a far image focus in accordance with the labeled power of the IOL, so patients with abnormally small pupils (~1mm) should achieve, at a minimum, the prescribed distance vision under photopic conditions; however, because this multifocal design has not been tested in patients with abnormally small pupils, it is unclear whether such patients will derive any near vision benefit.
4. Autorefractors may not provide optimal postoperative refraction of multifocal patients. Manual refraction is strongly recommended.
5. Recent contact lens usage may affect the patient's refraction; therefore in contact lens wearers, surgeons should establish corneal stability without contact lenses prior to determining IOL power.
6. When performing wavefront measurements on a patient with a multifocal toric lens, two different wavefronts are produced. One wavefront will be in focus (either far or near) and the other wavefront will be out of focus. In this situation, incorrect interpretation of the wavefront measurements is possible.
7. The long-term effects of intraocular lens implantation have not been determined. Therefore the physician should continue to monitor implant patients postoperatively on a regular basis.
8. Secondary glaucoma has been reported occasionally in patients with controlled glaucoma who received lens implants. The intraocular pressure of implant patients with glaucoma should be carefully monitored postoperatively.
9. Do not resterilize this intraocular lens by any method.
10. Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.
11. Do not store the lens in direct sunlight or at a temperature greater than 45°C (113°F). Do not autoclave the intraocular lens.
12. Prior to implanting, examine the lens package for proper lens model, dioptric power, and expiration date.
13. The surgeon should target emmetropia as this lens is designed for optimum visual performance when emmetropia is achieved.
14. Care should be taken to achieve centration of the intraocular lens.
15. Please refer to the specific instructions for use provided with the insertion instrument or system for the amount of time the IOL can remain folded before the IOL must be discarded.
16. When the insertion system is used improperly, the haptics of the TECNIS® Multifocal Toric 1-piece lens may become broken. Please refer to the specific instructions for use provided with the insertion instrument or system.

Warnings

Physicians considering lens implantation under any of the following circumstances should weigh the potential risk/benefit ratio:

1. Mild visual effects associated with multifocal IOLs may be expected because of the superposition of focused and unfocused images. These may include a perception of halos or glare around lights under nighttime conditions. It is expected that, in a small percentage of patients, the observation of such phenomena may be perceived as a hindrance, particularly in low illumination conditions. On rare occasions these visual effects may be significant enough that the patient will request removal of the TECNIS® Multifocal Toric 1-piece IOL.
2. Under low-contrast conditions, contrast sensitivity is reduced with a multifocal lens compared to a monofocal lens. Therefore, subjects with multifocal toric lenses should exercise caution when driving at night or in poor visibility conditions.
3. The splitting of the light into more than one focus may affect image quality and lead to some reduction of contrast sensitivity.
4. Well-informed patients with well-defined visual needs and preferences should be selected for TECNIS® Multifocal Toric lens implantation. The patients should be informed about the possibility that a decrease in contrast sensitivity and an increase of visual disturbances may affect their ability to drive a car under certain environmental conditions, such as driving at night or in poor visibility conditions.
5. Rotation of TECNIS® Multifocal Toric 1-piece IOLs away from their intended axis can reduce their astigmatic correction. Misalignment greater than 30° may increase postoperative refractive cylinder. If necessary, lens repositioning should occur as early as possible prior to lens encapsulation.
6. Carefully remove all viscoelastic from the capsular bag. Residual viscoelastic may allow the lens to rotate, causing misalignment of the TECNIS® Multifocal Toric 1-piece IOL with the intended axis of placement.
7. Children under the age of 2 years are not suitable candidates for intraocular lenses.
8. The TECNIS® Multifocal Toric 1-piece IOL should be placed entirely in the capsular bag and should not be placed in the ciliary sulcus.
9. Patients with any of the following conditions may not be suitable candidates for an intraocular lens because the lens may exacerbate an existing condition, may interfere with diagnosis or treatment of a condition or may pose an unreasonable risk to the patient's eyesight:
 a. Patients with recurrent severe anterior or posterior segment inflammation or uveitis.
 b. Patients in whom the intraocular lens may interfere with the ability to observe, diagnose or treat posterior segment diseases.
 c. Surgical difficulties at the time of cataract extraction and/or intraocular lens implantation which may increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
 d. A compromised eye due to previous trauma or developmental defects in which appropriate support of the IOL is not possible.
 e. Circumstance that would result in damage to the endothelium during implantation.
 f. Suspected microbial infection.
 g. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support for the IOL.
 h. Congenital bilateral cataracts.
 i. Recurrent severe anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye.
 j. Previous history of, or a predisposition to, retinal detachment.
 k. Patients with only one eye with potentially good vision.
 l. Medically uncontrollable glaucoma.
 m. Corneal endothelial dystrophy.
 n. Proliferative diabetic retinopathy.
10. AMO IOLs are single-use medical devices and are labeled with instructions for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The reuse/resterilization/reprocessing of AMO single-use medical devices may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

Adverse Events

Potential adverse events during or following surgery of the TECNIS® Multifocal Toric IOL 1-piece lens may include but are not limited to:
Endophthalmitis/Intraocular infection
Hyphema
Hypopyon
Pupillary block
Retinal detachment
Acute corneal decompensation
Secondary surgical intervention (including implant repositioning, removal, AC tap or other surgical procedure).
Any sight threatening complications

Adverse events that occurred during the clinical trial of parent lens TECNIS® Multifocal IOL, Model ZM900 (N=348) include the following: macular edema, endophthalmitis, and hypopyon (followed by non-lens related surgical re-interventions), and surgical re-intervention. Only three subjects out of 348 experienced lens-related events (3/348; 0.9%) in which the IOLs were ultimately removed (two due to halos and one due to blurry/filmy vision). Only two secondary surgical re-intervention events for the multifocal first eyes were lens-related (lens removal due to halos/glare and lens repositioning); the majority of secondary surgical re-intervention events were not lens-related (n=11, including iris prolapse/wound repair, lens exchange for refractive error and incorrect lens type, retinal repair, trabeculectomy and subsequent filtration bleb revisions and treatment injections for cystoid macular edema). The only persistent adverse event was raised intraocular pressure (IOP) due to secondary glaucoma requiring treatment.

Adverse events that occurred during the clinical trial of parent lens SENSAR 1 Piece IOL, Model AAB00, (N=123) include the following: macular edema, secondary surgical intervention (pars plana vitrectomy with membrane peel), and a torn haptic related to improper loading technique, resulting in lens exchange.

Additional complications that may occur following implantation of a multifocal intraocular lens include: objectionable visual quality under certain lighting conditions (e.g., chromatic aberrations, halos, night glare, decreased contrast), hazy/filmy images, "ghost" imaging of near and/or distant objects, diplopia and poor stereopsis. Some of these optical effects may be mitigated upon patients adaptation to the multifocal optic. However, some of the above listed complications may require secondary surgical intervention for intraocular lens exchange or explantation.

Detailed Device Description

The TECNIS® Multifocal Toric 1-piece IOL Model ZMT, is a one-piece, foldable, posterior chamber lens made of soft acrylic material. It incorporates a proprietary toric modified prolate (aspheric) design with cylinder axis markers on the anterior optic surface and a diffractive multifocal design on the posterior optic surface. The optic is 6.0 mm in diameter and the lens has an overall diameter of 13.0 mm. The add power is +4 diopters, corresponding to +3 diopters in the spectacle plane.

Lens Optic

- Optic Material: Optically clear, soft, foldable hydrophobic acrylic with a covalently bound UV absorber.
- Diopter power: 5.0 D to 34.0 D in 0.5 D increments
- Cylinder Power: 1.50 diopter, 2.25 diopter, 3.00 diopter, and 4.00 diopter (as measured at the IOL plane)

Conversion table for Cylinder powers

	Cylinder Powers (D)			
IOL Plane (labeled)	1.50	2.25	3.00	4.00
Corneal Plane	1.03	1.54	2.06	2.74

- Optic Center Thickness: 0.7 mm (+20.0D)
- Optic Edge Design: PROTEC 360 Square posterior edge
- Index of refraction: 1.47 at 35°C
- Light Transmittance: UV cut-off at 10% T for a +5.0 diopter lens (thinnest) and a +34.0 diopter lens (thickest) are shown in **Figure 1**.
- Modulation Transfer Function (MTF) Through-Focus Response: MTF values are shown in **Figure 2.**

Haptics

- Material: Soft foldable acrylic with a covalently bound UV absorber
- One-piece lens.
- Configuration: TRI-FIX design Modified C, integral with optic
- Thickness: 0.46 mm

© Johnson & Johnson Surgical Vision, Inc. 2019

Lens Power Calculations

Accurate keratometry and biometry are essential to successful visual outcomes. Preoperative calculation of the required spherical equivalent lens power for these posterior chamber intraocular lenses should be determined by the surgeon's experience, preference, and intended lens placement. Emmetropia should be targeted. The appropriate A-constant listed on the outer label is presented as a guideline and is a starting point for implant power calculations. The physician should determine preoperatively the spherical equivalent and cylindrical power of the lens to be implanted.

Use of the AMO provided Toric calculator is recommended for determining the appropriate Toric IOL spherical equivalent (SE) power, optimal axis of Toric IOL placement and Toric IOL cylinder (Cyl) power. The TECNIS® Multifocal Toric 1-piece lenses are labeled with the IOL distance spherical equivalent power.

Physicians requiring additional information on lens power calculations may contact the local AMO representative. Lens power calculation methods are described in following references:

- Hoffer KJ. The Hoffer Q formula: A comparison of theoretic and regression formulas. J Cataract Refract Surg. 19:700-712, 1993; ERRATA 20:677, 1994.
- Holladay JT. Standardizing constants for ultrasonic biometry, keratometry and intraocular lens power calculations. J Cataract Refract Surg. 23-1356-1370, 1997.
- Holladay JT, Musgrove KH, Prager TC, Lewis JW, Chandler TY and Ruiz RS. A three-part system for refining intraocular lens power calculations. J Cataract Refract Surg. 19:17-24 1988.
- Norrby NES. Unfortunate discrepancies. Letter to the editor and reply by Holladay JT. J Cataract Refract Surg. 24:433-434, 1998.
- Norrby S, Lydahl E, Koranyi G, Taube M. Reduction of trend errors in power calculation by linear transformation of measured axial lengths. J Cataract Refract Surg 2003; 29:100-105
- Olsen T, Olesen H, Thim K and Corydon L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. J. Cataract Refract Surg. 18:280-285, 1992.
- Retzlaff JA, Sanders DR and Kraff MC. Development of the SRK/T intraocular lens implant power calculation formula. J. Cataract Refract Surg. 16:333-340, 1990; ERRATA, 16:528, 1990.
- http://www.augenklinik.uni-wuerzburg.de/ulib/index.htm is in particular useful for Zeiss IOLMaster users.

Directions For Use

1. Prior to implanting, examine the lens package for IOL type, power, cylinder, proper configuration and expiration date.
2. Open the peel pouches and remove the lens in a sterile environment. Verify the dioptric power and cylinder power of the lens.
3. Examine the lens thoroughly to ensure particles have not become attached to it, and examine the lens' optical surfaces for other defects.
4. If desired, the lens may be soaked or rinsed in sterile balanced salt solution until ready for implantation.
5. Handle the lens by the haptic portion. Do not grasp the optical area with forceps.
6. Transfer the lens, using sterile technique, to an appropriate loading device.
7. The physician should consider the following points:
 - The surgeon should target emmetropia as this lens is designed for optimum visual performance when emmetropia is achieved.
 - Care should be taken to achieve centration of the intraocular lens.
8. AMO recommends using the UNFOLDER Emerald AR Series Implantation system with 1CART30 Cartridge, the ONE SERIES Ultra Implantation system with the 1VIPR30 cartridge or the UNFOLDER Platinum 1-Series Implantation system with the 1MTEC30 cartridge or an equivalent insertion instrument or system to insert the TECNIS® Multifocal Toric 1-piece lens. Only insertion instruments that have been validated and approved for use with this lens should be used. Please refer to the directions for use with the insertion instrument or system for additional information.

Caution: Do not use the lens if the package has been damaged. The sterility of the lens may have been compromised.

Selection and Placement of the TECNIS® Multifocal Toric 1-piece IOL

The astigmatism to be corrected should be determined from keratometry and biometry data rather than refractive data since the presence of lenticular astigmatism in the crystalline lens to be removed may influence results. The size and location of the surgical incision may affect the amount and axis of corneal astigmatism. In order to optimize IOL selection and axis placement, AMO provides a web-based tool (www.TecnisToricCalc.com) for the surgeon. Pre-operative keratometry and biometry data, incision location, and the surgeon's estimated surgically induced corneal astigmatism are used to determine the appropriate TECNIS® Multifocal Toric 1-piece IOL model, spherical equivalent lens power, and axis of placement in the eye.

For optimal results, the surgeon must ensure the correct placement and orientation of the lens within the capsular bag. The anterior surface of the IOL is marked with indentations (four at opposite sides) at the haptic/optic junction that identify the flat meridian of the TECNIS® Multifocal Toric 1-piece optic. These indentations form an imaginary line representing the plus cylinder axis (note: IOL cylinder steep meridian is 90° away). TECNIS® Multifocal Toric 1-piece IOL cylinder axis marks should be aligned with the post-incision steep corneal meridian (intended axis of placement).

Prior to surgery the operative eye should be market in the following manner:
With the patient sitting upright, precisely mark the twelve o'clock and/or the six o'clock position with a T marker, a surgical skin marker, or a marking pencil indicated for ophthalmic use. Using these marks as reference points, an axis marker can be used immediately prior to or during surgery to mark the axis of lens placement following the use of the web-based www.TecnisToricCalc.com to determine the optimal axis of placement.

After the lens is inserted, precisely align the axis marking indentations on the TECNIS® Multifocal Toric 1-piece IOL with the marked axis of lens placement. Carefully remove all viscoelastic from the capsular bag. This may be accomplished by manipulating the IOL optic with the I/A tip and using standard irrigation/aspiration techniques to remove all viscoelastic from the eye. Bimanual techniques may be used, if preferred, to ensure removal of viscoelastic from behind the lens implant. Special care should be taken to ensure proper positioning of the TECNIS® Multifocal Toric 1-piece IOL at the intended axis following viscoelastic removal. Residual viscoelastic may allow the lens to rotate, causing misalignment of the TECNIS® Multifocal Toric 1-piece IOL with the intended axis of placement.

Misalignment of the axis of the lens with the intended axis of placement may compromise its astigmatic correction. Such misalignment can result from inaccurate keratometry or marking of the cornea, inaccurate placement of the TECNIS® Multifocal Toric 1-piece IOL axis during surgery, an unanticipated surgically induced change in the cornea, or physical rotation of the TECNIS® Multifocal Toric 1-piece IOL. In order to minimize this effect, the surgeon should be careful to ensure that preoperative keratometry and biometry is accurate and that the IOL is properly oriented prior to the end or surgery.

Patient Card:

An implant identification card, to be supplied to the patient, is included in the package. The patient should be instructed to keep the card as a permanent record of his/her implant and to show the card to any eye care practitioner he/she may see in the future.

Reporting:

Adverse events and/or potentially sight-threatening complications that may reasonably be regarded as lens related and that were not previously expected in nature, severity or rate of occurrence must be reported to AMO. This information is being requested from all surgeons in order to document potential long-term effects of intraocular lens implantation, especially in younger patients.

Physicians are required to report these events in order to aid in identifying emerging or potential problems with posterior chamber lenses. These problems may be related to a specific lot of lenses or may be indicative of long-term problems associated with these lenses or with intraocular lenses in general.

How Supplied:

The TECNIS® Multifocal Toric 1-piece lenses are supplied sterile in a lens case within a double aseptic transfer peel pouch. The double aseptic transfer peel pouch is sterilized with ethylene oxide and should be opened only under sterile conditions. The pouch and product labels are enclosed in a shelf pack. The external surfaces of the outer pouch are not sterile.

Expiration Date

The expiration date on the lens package is the sterility expiration date. This lens should not be implanted after the indicated sterility expiration date.

Return/Exchange Policy

Contact the local AMO representative for the return lens policy. Return the lens with the proper identification and the reason for the return. Label the return as a biohazard. Do not attempt to resterilize the lens.

Patient Information

Each patient should receive information regarding intraocular lenses prior to the decision to implant an intraocular lens.

Manufactured in the Netherlands:

AMO Groningen BV, 9728 NX Groningen, The Netherlands

Manufactured for:

Johnson & Johnson Surgical Vision, Inc., 1700 E. St. Andrew Pl., Santa Ana, CA 92705 USA

Johnson & Johnson Surgical Vision, Inc.
1700 E. St. Andrew Place
Santa Ana, CA 92705 USA

CE 0344

EC REP AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland


Symbol/Explanation

SYMBOL	EXPLANATION
STERILE EO	Sterilized using Ethylene Oxide
	Do Not Reuse
	Use By (YYYY-MM-DD: year-month-day)
	Caution; See Instructions for Use
EC REP	Manufacturer
	Authorized Representative in the European Union
	Do Not Resterilize
45°C 113°F	Upper Limit Of Temperature
	Keep Away From Sunlight
	Do Not Use if Package is Damaged
	Date of Manufacture (YYYY-MM-DD: year-month-day)

AMO, TECNIS, TRI-FIX, SENSAR, PROTEC, UNFOLDER and ONE SERIES are trademarks of Johnson & Johnson Surgical Vision, Inc.

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Figure 1
Light Transmittance


Legend:
Curve 1: Spectral transmittance curve of a typical 5 diopter IOL (thinnest), UV cut-off at 10% T is 375.4 nm.
Curve 2: Spectral transmittance curve of a typical 34 diopter IOL (thickest), UV-cut-off at 10% T is 380.4 nm.
Curve 3: Spectral transmittance curve corresponding to 53-year-old Phakic eye*.

Note: The cut-off wavelengths and the spectral transmittance curves represent the range of the transmittance values of IOLs (5-34 diopter) made with this material.

* Boettner, E.A., and Wolter J.R. Transmission of the Ocular Media. Investigative Ophthalmology. 1962; 1:776-783.

Figure 2
ACE Model- MTF Through-Focus for Multiple Pupil Sizes, 20-Diopter, 3.0 Cylinder Power Lens

Note: The figure above describes the optical performance of the TECNIS® Multifocal Toric 1 piece lens - Model ZMT300 in ACE model eyes at 50 cycles/mm using white light. The graph depicts the focus as it is gradually shifted from that of a far object to increasingly nearer objects, with higher numbers typically indicating better performance. The ACE eye model has the spherical aberration of an average natural cornea, which the TECNIS® lens is designed to compensate. The combination is aberration free.

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- Hoffer KJ. The Hoffer Q formula: A comparison of theoretic and regression formulas. J Cataract Refract Surg. 19:700-712, 1993; ERRATA 20:677, 1994.
- Holladay JT. Standardizing constants for ultrasonic biometry, keratometry and intraocular lens power calculations. J Cataract Refract Surg. 23-1356-1370, 1997.
- Holladay JT, Musgrove KH, Prager TC, Lewis JW, Chandler TY and Ruiz RS. A three-part system for refining intraocular lens power calculations. J Cataract Refract Surg. 19:17-24 1988.
- Norrby NES. Unfortunate discrepancies. Letter to the editor and reply by Holladay JT. J Cataract Refract Surg. 24:433-434, 1998.
- Norrby S, Lydahl E, Koranyi G, Taube M. Reduction of trend errors in power calculation by linear transformation of measured axial lengths. J Cataract Refract Surg 2003; 29:100-105.
- Olsen T, Olesen H, Thim K and Corydon L. Prediction of pseudophakic anterior chamber depth with the newer IOL calculation formulas. J. Cataract Refract Surg. 18:280-285, 1992.
- Retzlaff JA, Sanders DR and Kraff MC. Development of the SRK/T intraocular lens implant power calculation formula. J. Cataract Refract Surg. 16:333-340, 1990; ERRATA, 16:528, 1990.
- http://www.augenklinik.uni-wuerzburg.de/ulib/index.htm — ресурс, особенно полезный для пользователей диагностической системы Zeiss IOLMaster.

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