The Unfolder Platinum 1 Series Implantation System
with 1MTEC30 Cartridge and AMO Acrylic 1-Piece IOLs

Z311145 Rev. C
Revision Date: 09/2019

EN

RX Only

Description
The UNFOLDER Platinum 1 Series Implantation System is designed for use with AMO Acrylic 1-Piece intraocular lenses. The UNFOLDER Platinum 1 Series Implantation System consists of the following:

- The UNFOLDER Platinum 1 Series handpiece, Model DK7796: an autoclaveable, reusable titanium handpiece which is provided non-sterile (Fig. 1).
- The UNFOLDER Platinum 1 Series cartridge, Model 1MTEC30: a single-use, non-pyrogenic cartridge which is provided sterile (Fig. 2).

Indication For Use
The Model DK7796 handpiece is used in combination with the Model 1MTEC30 cartridge to fold and assist in inserting AMO Acrylic 1-Piece intraocular lenses, ONLY into the capsular bag.

Contraindications
Do not use the handpiece if the rod tip appears nicked or damaged in any way.

Instructions For Use
Read the entire instructions before using this device.

Handpiece Preparation
1. Clean and sterilize the handpiece according to the method specified under **"Instructions for Reprocessing."** Sterility assurance is the responsibility of the user.
2. Confirm that the handpiece is engraved with the designation "The UNFOLDER Platinum 1 Series."
3. Prior to each use visually inspect the rod tip to establish that there are no material deposits or damage.
4. Ensure the handpiece plunger is fully retracted for cartridge loading.

Cartridge Preparation & Lens Loading
1. Inspect the cartridge tray before each use. Ensure that the sterilized, outer sealed tray contains the inner tray which holds one cartridge. Do not use the cartridge if the tray is damaged or the seal is broken. Remove the inner tray from the sealed tray.
2. Remove the cartridge from the inner tray and fill the entire inside of the cartridge with viscoelastic (Fig 3: Filling Viscoelastic).
3. Using forceps grasp the lens by the optic edge only anterior side up. Use the IOL diagram on the cartridge for correct IOL orientation (Fig. 2:a. IOL Diagram). Do not grasp the trailing haptic.
4. Hold the cartridge with the IOL diagram facing up. Engage the lead haptic with the canopy to initiate tucking of the haptic over the lens optic body, (Fig. 4:a. area to engage lead haptic). Continue advancing the lens into the cartridge so that the lead haptic sweeps over the optic body in one motion (Fig. 5). Stop when half the optic is inside the cartridge. Ensure that the lead haptic is fully tucked over the optic body before proceeding (Fig. 6).
5. Using forceps grasp the trailing haptic and tuck it over the optic body (Fig. 7). Advance the lens until it is positioned under the lens diagram on the cartridge so that the lens can be held in the folded position until it is ready to be delivered (Fig. 2:a. lens diagram, Fig. 8: a. Leading and trailing haptics tucked over optic). Once the lens is under the diagram, ensure that the lens and haptics remain in place after retracting the forceps.

Lens Delivery
1. Insert the cartridge into the handpiece with the IOL diagram and canopy facing up and the cartridge tip bevel down (Fig. 2:a. IOL Diagram, b. Canopy, c. Cartridge tip). Insert the tip first so that the cartridge slides into the cartridge slots in the handpiece (Fig. 9:a. Cartridge slot) and then push down firmly on the back end of the cartridge to securely snap it into the handpiece (Fig. 10). Both sides of the cartridge should snap into place. If the rod tip makes contact with the cartridge upon entry the cartridge has not been inserted correctly. Retract the rod and ensure that the cartridge is fully snapped into the handpiece before moving forward.
2. Insert the cartridge bevel down into the incision and slowly advance the plunger into the cartridge.
3. Rotate the knob to advance the plunger to push the lens forward. Ensure that the rod tip engages the lens optic body and remains behind the optic body through the entire delivery.
4. Slowly advance the lens, ensuring that the haptics remain tucked over the optic, until the lens is fully released from the cartridge tube.
5. Slowly reverse screw the knob to withdraw the plunger.

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Instructions For Reprocessing
First Use
1. Handpiece is supplied "non-sterile" and must be cleaned and sterilized before first use.

Dis-assembly
1. Remove and dispose of cartridge. Disassemble the handpiece before reprocessing. **Do not autoclave or resterilize the cartridge by any method.**
2. The handpiece is comprised of two components, the main body and the plunger (Fig. 11:a. Main Body, b. Plunger). To separate the two, (1) push the plunger from the knurled ring into the main body and (2) rotate the plunger in a counter-clockwise direction (Fig. 12).

Cleaning
1. It is recommended that the instrument be cleaned as soon as possible after use.
2. It is recommended that a non-ionic detergent suitable for cleaning metallic surgical instruments such as Lancerzyme detergent is used.
3. Fully immerse the two components of the disassembled handpiece into a detergent solution not exceeding 30°C. Vigorously wash and scrub the two components with an appropriate brush for at least one minute, applying the detergent solution to all surfaces. It is recommended that the instrument be cleaned as soon after use as possible; however where blood, tissue, saline or viscoelastic has been left to dry it is recommended that the two components are left to soak for 30 minutes in the detergent solution.
4. After manual cleaning, rinse the components in clean water for a minimum of 3 times.
5. Place the two separated components in a suitable container and subject to an automatic dishwasher cycle, consisting of 3 key elements, a 2-minute pre-wash, a 3-minute detergent wash at 93°C and a drying stage sufficient to ensure complete absence of moisture on the two components.
6. Avoid the use of acid base products when cleaning the instrument and always follow the guidelines set by the detergent and mechanical cleaner manufacturer.
7. If possible, sterilizing trays should be used to hold the instrument in a secure position during cleaning and storage to prevent damage. Store handpiece in a dry area.

Sterilization
Moist heat (steam) in autoclavable bags is the preferred method of sterilization. Please ensure system manufacturer's guidelines are followed at all times.

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Please reference the ANSI/AAMI ST79: 2006 Guideline "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities," Section 8.6.1, "Sterilization Parameters for Wrapped or Containerized Items."

Another possible method for steam sterilization is the following:

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Re-Assembly
1. To re-assemble the handpiece gently push the plunger into the main body until the pin (Fig. 13:a) is flush with the back of the main body.
2. Rotate plunger until the pin engages the slot (Fig. 13:b). Continue rotating plunger in clockwise direction until pin locks into position.

General Care
- Do not allow blood, tissue, saline or viscoelastic to dry on the handpiece.
- Never use saline or balanced salt solution for rinsing the handpiece.

Reprocessing validation
All cleaning and sterilization processes require validation at the point of use. Due to the different reprocessing options, Duckworth & Kent (D&K) cannot validate all possible process variations. The user is responsible for qualifying any method that deviates from D&K's recommended method of cleaning and sterilizing The UNFOLDER Platinum 1 Series handpiece.

Limitations on reprocessing
Repeated processing has minimal effect on the instruments. End of life is normally determined by wear and damage due to use.

Precautions
1. The use of viscoelastics is required when loading the IOL into the cartridge. For optimal performance use the AMO HEALON® family of viscoelastics. Do not use balanced salt solution.
2. The combination of low operating room temperatures and high IOL diopter powers may require a slower delivery.
3. Do not use if any component of The UNFOLDER Platinum 1 Series Implantation System has been dropped or inadvertently struck while outside of the shipping case.
4. Do not store the cartridges at temperatures under 41°F (5°C) or over 86°F (30°C).

Warnings
1. The UNFOLDER Platinum 1 Series Implantation System should be used ONLY with AMO Acrylic 1-Piece IOLs.
2. Do not use if the cartridge tip is cracked or split prior to implantation.
3. Never release the plunger until the optic body has been completely released from the cartridge tube.
4. The lens and cartridge should be discarded if the lens has been folded within the cartridge for more than 5 minutes.
5. If the IOL is not properly placed in the cartridge, the IOL may be damaged and/or implanted upside down.
6. Do not attempt to modify or alter this device or any of the components, as this can significantly affect the function and/or structural integrity of the design.
7. Use of methyl cellulose viscoelastics is not recommended as they have not been validated with the Model 1MTEC30 Cartridge and AMO Acrylic 1-Piece IOLs.
8. Do not implant lens if rod tip becomes jammed in the cartridge.
9. AMO single-use medical devices are labeled with instructions for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The reuse/resterilization/reprocessing of AMO single-use medical devices may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.

Symbols used on Sterile Packaging

SYMBOL	ENGLISH
	SEE INSTRUCTIONS FOR USE
	DO NOT REUSE
STERILE EO	STERILIZED USING ETHYLENE OXIDE
REF	MODEL NUMBER
LOT	LOT/BATCH NUMBER
	USE BY (YYYY-MM-DD YEAR-MONTH-DAY)
5°C (41°F) 30°C (86°F)	STORE AT 5-30°C (41-86°F)
	DO NOT RESTERILIZE
	MANUFACTURER
EC REP	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN UNION
	DO NOT USE IF PACKAGE IS DAMAGED
	DATE OF MANUFACTURE (YYYY-MM-DD YEAR-MONTH-DAY)

UNFOLDER, AMO, and HEALON are trademarks of Johnson & Johnson Surgical Vision, Inc.

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The DK7786/DK7791 Implantation System
with ONE SERIES Ultra Cartridge and AMO Acrylic 1-Piece IOLs

IOL Data
Description
The DK7786/DK7791 implantation System is designed for use with AMO Acrylic 1-Piece intraocular lenses.
The DK7786/DK7791 implantation System consists of the following components:
- A single-use, sterile, preloaded ONE SERIES Ultra Cartridge in which the IOL/viscoelastic is autoclaved, movable, and provided non-sterile.
- A reusable non-sterile ONE SERIES Ultra Cartridge which is assembled by the user.

Sodium Hyaluronate (HA), used in the ONE SERIES Ultra Cartridge coating, is produced by microbiological fermentation method.

Instruction For Use
Either the DK7786 or DK7791 handpiece is used in combination with the ONE SERIES Ultra Cartridge to fold and assist in inserting the AMO Acrylic 1-Piece intraocular lenses (IOLs) into the capsular bag.

Contraindications
Do not use the handpiece if the nozzle appears nicked or damaged in any way.

Instruction For Use
Read the entire instructions before using this device.
Preparation
Clean and sterilize the handpiece according to the method specified under "Instructions for Reprocessing".
Sterility assurance is the responsibility of the user.
Prepare the handpiece for use as follows (designation DK7786 or DK7791):
- Prior to each use visually inspect for any cracks that may lead to oxidation and flash them on for material deposits or damage.
- Ensure the handpiece plunger is fully retracted for cartridge loading.

Cartridge Preparation and Loading
1. Inspect the cartridge tip before each use. Ensure that the electrolyte outer sealed tray contains the inner tray which holds one cartridge. Open the outer tray and lift the inner tray out. Open the inner tray and lift the cartridge out by the cartridge body.
2. Remove the cartridge from the inner tray and fill the entire inside of the cartridge with viscoelastic (Fig. 4 Ring Segment).
3. Using bumpy grasp the lens by the optic edge only anterior side up. Use the IOL diagram on the cartridge for correct IOL orientation. Load the IOL into the cartridge (Fig. 5 Lens Loading).
4. Fill the cartridge with the IOL diagram facing up. Engage the load nozzle with the cartridge. It rotates locking of the nozzle over the lens optic body (Fig. 6), a lens is integra ted haptic). Continue advancing the lens into the cartridge so that the load haptic is positioned in the cartridge tip. Ensure that the lens is loaded into the cartridge (Fig. 7 Lens Load Tip Haptic) and subsequently advance the lens into the cartridge tip.
5. Using bumpy grasp the trailing haptic and tuck it over the optic body (Fig. 8). Advance the lens until it is positioned under the cartridge tip (Fig. 9). Ensure that the lens is loaded into the cartridge tip (Fig. 10 Lens Load Tip Haptic) and that the lens diagram (Fig. 9 a. Loading and Folding Haptics tucked over optic). Once the lens is under the diagram, ensure that the lens and haptics remain in place after retracting the forceps.

Lens Delivery
1. Insert the cartridge into the handpiece with the IOL diagram and cranjo facing up and the cartridge tip bevel down (Fig. 2 Cartridge Tip). Push the cartridge into the handpiece until it is fully seated and the cartridge tip is in the eye (Fig. 3 Cartridge Tip). Fig 10 a. Cartridge seal and flush guide shown from the back end of the cartridge to accurately assess the IOL position (Fig. 11 Lens Load Tip Haptic). Fig 11. Both sides of the cartridge should snap into place if the tool is fully seated. Fig 12. Cartridge seal and flush guide shown from the back end of the cartridge to accurately assess the IOL position (Fig. 13 Lens Load Tip Haptic). Fig 13. Both sides of the cartridge should snap into place if the tool is fully seated.
2. Insert the cartridge tip bevel down into the incision and slowly advance the plunger into the cartridge. For the DK7786 handpiece, to avoid damage to the IOL, do not advance the plunger beyond the end of the cartridge body (Fig. 14).
3. The DK7786 handpiece is designed to be used with a single hand. Use the thumb of the hand that holds the handpiece to engage the plunger. Slowly advance the lens, ensuring that the haptics remain tucked over the optic, until the lens is fully released from the cartridge (Fig. 15).
4. The DK7786 allows the delivery to commence slowly, always release the plunger so that the rod tip fully retracts before removing the cartridge from the eye. For DK7791, slowly reverse screw the knob to withdraw the plunger.

Caution: Do not use the handpiece if the nozzle appears nicked or damaged in any way. If the IOL is not delivered, rotate handpiece body with cartridge still inside to extend the rod tip and remove the cartridge tip to the right to re-orient the haptic fix.

Instructions For Reprocessing
Disposal
- Handpiece is supplied "non-sterile" and must be cleaned and sterilized before first use.

Cleaning
Remove and dispose of cartridge. Clean handpiece immediately after surgery, the metal body and the plunger (Fig. 1). a. Main Body, b. Plunger). To separate the two, (1) push the plunger from the hexagonal flange into the main body and (2) rotate the plunger one a counter clockwise direction to release the plunger.
The DK7786 handpiece is designed to be used with a single hand. Use the thumb of the hand that holds the handpiece to engage the plunger (Fig. 15).

Cleaning
It is recommended that the instrument be cleaned as soon as possible after use.
It is recommended that a soft cloth moistened with cleaning solution, ultrasonic cleaning solution, surgical instruments such as Lancetronic detergent is used.
Ensure that the components of the disassembled handpiece onto a detergent solution not exceeding 30°C. Vigorously shake the components in the detergent solution with an appropriate brush for at least one minute, applying the detergent to all surfaces. It is recommended that the instrument be cleaned as soon after use as possible, however where blood, tissue, or other debris has dried on the surface, it may be necessary to soak the components until the debris is soft enough to be removed with a brush.
After manual cleaning, rinse the components in clean water for a minimum of 3 times.
Place the two components into a suitable container and subject to an automatic dishwasher cycle, consisting of 3 key elements, a 2 minute dry wash, a 2 minute detergent wash at 55°C and a drying stage sufficient to ensure complete absence of moisture from the components.
Avoid the use of acid based detergents.
Inspect the components for cleanliness and ensure that the instrument and plunger follow the assembly order and the mechanical steam manufacturer.

Sterilization
Steam sterilization in autoclavable bags is the preferred method of sterilization. Please ensure system manufacturer's guidelines are followed at all times.

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UNFOLDER® Platinum Push Implantation System with 1VTEC30 Cartridge and TECNIS® and SENSAR Acrylic 1-Piece IOLs

Z311354 Rev. A
Revision Date: 10/2018

EN

Description

The UNFOLDER® Platinum Push Implantation System is designed for use with TECNIS® and SENSAR Acrylic 1-Piece intraocular lenses.

The UNFOLDER® Platinum Push Implantation System consists of the following:

- The UNFOLDER® Platinum Push handpiece, Model DK7798, an autoclavable, reusable titanium handpiece which is provided non-sterile (Fig. 1).
- The UNFOLDER® Platinum Push cartridge, Model 1VTEC30: a single-use, non-pyrogenic cartridge which is provided sterile (Fig. 2).

This syringe style system is for use with TECNIS® and SENSAR Acrylic 1-Piece intraocular lenses up to 28.0 Diopter.

Indication For Use

The Model DK7798 handpiece is used in combination with the Model 1VTEC30 cartridge to fold and assist in inserting TECNIS® and SENSAR Acrylic 1-Piece intraocular lenses, ONLY into the capsular bag.

Contraindications

Do not use the handpiece if the rod tip appears nicked or damaged in any way

Instructions For Use

Read the entire instructions before using this device.

Handpiece Preparation

1. Clean and sterilize the handpiece according to the method specified under **"Instructions for Reprocessing."** Sterility assurance is the responsibility of the user.
2. Prior to each use visually inspect the rod tip and O-ring to establish that they are clean and undamaged
3. Replace the O-ring, if damaged from the packet of replacement O-rings supplied with this device. The packet of replacement O-rings should be retained for future use. Contact the local sales representative for additional replacement O-rings quoting reference number DK0010.

- To remove the O-ring place the plunger assembly in the middle of a dry gloved hand (Fig. 14). Grip the plunger assembly firmly and pull the whole assembly out of the closed hand in the direction shown (Fig. 15).
- To install, place the O-ring on the plunger, grip with thumb and index finger and pull toward the handle until the O-ring seats in the first groove.

4. Ensure the handpiece plunger is fully retracted for cartridge loading.

Cartridge Preparation & Lens Loading

1. Inspect the cartridge tray before each use. Ensure that the sterilized, outer sealed tray contains the inner tray which holds one cartridge. Do not use the cartridge if the tray is damaged or the seal is broken. Remove the inner tray from the sealed tray.
2. Remove the cartridge from the inner tray and fill the entire inside of the cartridge with viscoelastic (Fig 3: Filling Viscoelastic).
3. Using forceps grasp the lens by the optic edge only anterior side up. Use the IOL diagram on the cartridge for correct IOL orientation (Fig. 2:a. IOL Diagram). Do not grasp the trailing haptic.
4. Hold the cartridge with the IOL diagram facing up. Engage the lead haptic with the canopy to initiate tucking of the haptic over the lens optic body, (Fig. 4:a. area to engage lead haptic). Continue advancing the lens into the cartridge so that the lead haptic sweeps over the optic body in one motion (Fig. 5). Stop when half the optic is inside the cartridge. Ensure that the lead haptic is fully tucked over the optic body before proceeding (Fig. 6).
5. Using forceps grasp the trailing haptic and tuck it over the optic body (Fig. 7). Advance the lens until it is positioned under the lens diagram on the cartridge so that the lens can be held in the folded position until it is ready to be delivered (Fig. 2: a. lens diagram, Fig. 8: a. Leading and trailing haptics tucked over optic). Once the lens is under the diagram, ensure that the lens and haptics remain in place after retracting the forceps.

Lens Delivery

1. Insert the cartridge into the handpiece with the IOL diagram and canopy facing up and the cartridge tip bevel down (Fig. 2:a. IOL Diagram, b. Canopy, c. Cartridge tip). Insert the tip first so that the cartridge slides into the cartridge slots in the handpiece (Fig. 9:a. Cartridge slot) and then push down firmly on the back end of the cartridge to securely snap it into the handpiece (Fig. 10). Both sides of the cartridge should snap into place. If the rod tip makes contact with the cartridge upon entry the cartridge has not been inserted correctly. Retract the rod and ensure that the cartridge is fully snapped into the handpiece before moving forward.
2. Insert the cartridge bevel down into the incision and slowly advance the plunger into the cartridge.
3. Advance the plunger to push the lens forward. Ensure that the rod tip engages the lens optic body and remains behind the optic body through the entire delivery.
4. Slowly advance the lens for a minimum of 12 seconds to deliver the IOL (no faster than a rate of 2mm/sec), ensuring that the haptics remain tucked over the optic, until the lens is fully released from the cartridge tube.
5. Slowly withdraw the plunger.

Note: If the lead haptic is not folded within the lens or pointing to the left as the IOL is delivered, rotate handpiece body with cartridge until haptic is pointed to the left or sweep the cartridge tip to the right to re-orient the haptic tip.

Instructions For Reprocessing

First Use

1. Handpiece is supplied "non-sterile" and must be cleaned and sterilized before first use.

Dis-assembly

1. Remove and dispose of cartridge. Disassemble the handpiece before reprocessing. **Do not autoclave or resterilize the cartridge by any method.**
2. The handpiece is comprised of two components, the main body and the plunger (Fig. 11:a. Main Body, b. Plunger). To separate the two, (1) push the plunger from the knurled ring into the main body and (2) rotate the plunger in a counter-clockwise direction (Fig. 12).

Cleaning (Model DK7798)

1. It is recommended that the instrument be cleaned as soon as possible after use.
2. It is recommended that a non-ionic detergent suitable for cleaning metallic surgical instruments such as Lancerzyme detergent is used.
3. Fully immerse the two components of the disassembled handpiece into a detergent solution not exceeding 30°C. Vigorously wash and scrub the two components with an appropriate brush for at least one minute, applying the detergent solution to all surfaces. It is recommended that the instrument be cleaned as soon after use as possible; however where blood, tissue, saline or viscoelastic has been left to dry it is recommended that the two components are left to soak for 30 minutes in the detergent solution.
4. After manual cleaning, rinse the components in clean water for a minimum of 3 times.
5. Place the two separated components in a suitable container and subject to an automatic dishwasher cycle, consisting of 3 key elements, a 2-minute pre-wash, a 3-minute detergent wash at 93°C and a drying stage sufficient to ensure complete absence of moisture on the two components.
6. Avoid the use of acid base products when cleaning the instrument and always follow the guidelines set by the detergent and mechanical cleaner manufacturer.
7. If possible, sterilizing trays should be used to hold the instrument in a secure position during cleaning and storage to prevent damage. Store handpiece in a dry area.

Sterilization (Model DK7798)

Moist heat (steam) in autoclavable bags is the preferred method of sterilization. Please ensure system manufacturer's guidelines are followed at all times.

Possible methods for steam sterilization are the following:

- Gravity displacement cycle suitable for wrapped items at 270°F-275°F (132°C for 15 minutes or 135°C for 10 minutes), or
- Prevacuum cycle consisting of wrapped items at 270°F-275°F (132°C-135°C) for 4 minutes.

Please reference the ANSI/AAMI ST79: 2006 Guideline "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities," Section 8.6.1, "Sterilization Parameters for Wrapped or Containerized Items." Another possible method for steam sterilization is the following:

- Cycle for wrapped goods at a sterilization temperature of 134°C with a maximum allowable temperature of 137°C for a minimum holding time of 3 min.

Please reference HTM 2010 Part 4, Section 5.0, "Operation of porous load sterilizers."

A sterilization cycle of 126°C with a holding time of 26 minutes has been validated.

A sterilization cycle of 134°C with a holding time of 3 minutes has also been validated.

Re-Assembly

1. To re-assemble the handpiece gently push the plunger into the main body until the pin (Fig. 13:a) is flush with the back of the main body.
2. Rotate plunger until the pin engages the slot (Fig. 13:b). Continue rotating plunger in clockwise direction until pin locks into position.

General Care

- Do not allow blood, tissue, saline or viscoelastic to dry on the handpiece.
- Never use saline or balanced salt solution for rinsing the handpiece.

Reprocessing validation (Model DK7798)

All cleaning and sterilization processes require validation at the point of use. Due to the different reprocessing options, Duckworth & Kent (D&K) cannot validate all possible process variations. The user is responsible for qualifying any method that deviates from D&K's recommended method of cleaning and sterilizing the UNFOLDER® Platinum Push handpiece.

Limitations on reprocessing

Repeated processing has minimal effect on the instruments. End of life is normally determined by wear and damage due to use.

Precautions

1. The use of viscoelastics is required when loading the IOL into the cartridge. For optimal performance use the HEALON® family of viscoelastics. Do not use balanced salt solution.
2. The combination of low operating room temperatures and high IOL diopter powers may require a slower delivery.
3. Do not use if any component of the UNFOLDER® Platinum Push Implantation System has been dropped or inadvertently struck while outside of the shipping case.
4. Do not store the cartridges at temperatures under 41°F (5°C) or over 86°F (30°C).
5. The model DK7798 handpiece should be used ONLY with TECNIS® and SENSAR Acrylic 1-Piece IOLs with diopter sizes of 28.0D or less at insertion speeds no greater than 2mm/second (minimum 12 seconds to deliver the IOL).
6. Evaluate O-ring for wear, replace if damaged and then resterilize the device.

Warnings

1. The UNFOLDER® Platinum Push Implantation System should be used ONLY with TECNIS® and SENSAR Acrylic 1-Piece IOLs.
2. Do not use if the cartridge tip is cracked or split prior to implantation.
3. The lens and cartridge should be discarded if the lens has been folded within the cartridge for more than 5 minutes.
4. If the IOL is not properly placed in the cartridge, the IOL may be damaged and/or implanted upside down.
5. Do not attempt to modify or alter this device or any of the components, as this can significantly affect the function and/or structural integrity of the design.
6. Use of methyl cellulose viscoelastics is not recommended as they have not been validated with the Model 1VTEC30 Cartridge and TECNIS® and SENSAR Acrylic 1-Piece IOLs.
7. Do not implant lens if rod tip becomes jammed in the cartridge.
8. Johnson & Johnson Surgical Vision, Inc. single-use medical devices are labeled with instructions for use and handling to minimize exposure to conditions which may compromise the product, patient, or the user. The reuse/resterilization/reprocessing of Johnson & Johnson Surgical Vision, Inc. single-use medical devices may result in physical damage to the medical device, failure of the medical device to perform as intended, and patient illness or injury due to infection, inflammation, and/or illness due to product contamination, transmission of infection, and lack of product sterility.
9. Advancing the plunger in less than the recommended 12 seconds may increase the potential for lens acceleration upon delivery.

Symbols used on Sterile Packaging

SYMBOL	ENGLISH
	Consult Instruction for use
	Do not re-use
STERILE EO	Sterilized using ethylene oxide
REF	Catalogue number
LOT	Batch code
	Use by date (YYYY-MM-DD) year-month-day
5°C (41°F) 30°C (86°F)	Temperature limit store at 5-30°C (41-86°F)
	Do not resterilize
	Manufacturer
EC REP	Authorized Representative in the European Community
	Do not use if package is damaged
	Date of manufacture (YYYY-MM-DD) year-month-day

UNFOLDER, TECNIS, SENSAR, and HEALON are trademarks of Johnson & Johnson Surgical Vision, Inc.

Lancerzyme is a trademark of Lancer UK.

PRODUCT OF THE USA:
AMO Puerto Rico Manufacturing, Inc.
Road 402, Km. 4.2
Añasco, Puerto Rico 00610 USA

CE 0344

Johnson & Johnson Surgical Vision, Inc.
1700 E. St. Andrew Place
Santa Ana, CA 92705 USA
1-877-266-4543

EC REP AMO Ireland
Block B
Liffey Valley Office Campus
Quarryvale, Co. Dublin, Ireland

Johnson & Johnson Regulatory Affairs Surgical Vision

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The UNFOLDER Emerald Series Implantation System

Description
The UNFOLDER Emerald Series Implantation System consists of the following components:
- The UNFOLDER Emerald Series Cartridge, Model EMERALDC
- The UNFOLDER Emerald Series Implant, Model EMERALDIMP (Fig. 2a), EMERALDIMP-L (Fig. 2b), or EMERALDIMP-S (Fig. 2c)

The UNFOLDER Emerald Series Cartridge is provided sterile and non-pyrogenic and is available (Fig. 2) in the following configurations:
The UNFOLDER Emerald Series Implant is provided sterile and is available (Fig. 2):
The UNFOLDER Emerald Series Implantation System is designed for use ONLY with ARC-ALIGN Navigation Software.

Prescription Statement
Rx Only – In the U.S., prescription only device.

Contraindications, Warnings
Contraindications
The UNFOLDER Emerald Series Implantation System is used to find and assist in inserting the UNFOLDER Emerald Series Implant. There are no known contraindications.

Warnings
The UNFOLDER Emerald Series Implantation System is provided non-sterile and must be cleaned and sterilized before use. Failure to properly clean and sterilize the device before each use could result in infection.
- Thoroughly rinse and flush the handpiece exposure with demineralized or distilled water to remove all visible debris. Failure to remove all debris could result in contamination and/or sterilization water. Follow the manufacturer's instructions for the ultrasonic cleaning bath.
- After ultrasonic cleaning, thoroughly flush and rinse the instrument under running water and then dry it with a soft, lint-free cloth or sterilization wrap. The device should be dried with a hot air blower or lint-free cloth before returning to storage.
- Immediately after each use, thoroughly flush and rinse the instrument under running water and then dry it with a soft, lint-free cloth or sterilization wrap.
- Clean the handpiece with a mild soap solution and gently scrub with a soft bristle brush paying particular attention to the inside and the rim of the instrument.
- Thoroughly rinse and flush the handpiece exposure with demineralized or distilled water to remove all visible debris. Failure to remove all debris could result in contamination and/or sterilization water. Follow the manufacturer's instructions for the ultrasonic cleaning bath.
- After ultrasonic cleaning, thoroughly flush and rinse the instrument under running water and then dry it with a soft, lint-free cloth or sterilization wrap. The device should be dried with a hot air blower or lint-free cloth before returning to storage.
- Immediately after each use, thoroughly flush and rinse the instrument under running water and then dry it with a soft, lint-free cloth or sterilization wrap.
- Clean the handpiece with a mild soap solution and gently scrub with a soft bristle brush paying particular attention to the inside and the rim of the instrument.
- Thoroughly flush the instrument under running water and then do a final rinse using demineralized or distilled water. Dry the device with a hot air blower or lint-free cloth before returning to storage.

AUTOMATED CLEANING
- An initial rinse of the numerous of reprocessing procedures is determined by the function / wear of the device. In case of damage the device should be reprocessed immediately. Do not let the device dry out.
- Clean the device immediately after use. Don't use a flushing detergent or hot water (>40°C) as this can cause protein fixation and make cleaning more difficult. Do not use cleaning agents containing aldehyde as this can lead to permanent fixation of residual blood and damage or contamination to the environment.
- Do not use abrasive cleaning agents for manual cleaning as a protective layer for the instrument can be damaged. Follow the instructions of the disinfectant and start the program.
 - 4 minutes pre-washing with cold water (<40°C)
 - 5 minutes cleaning with alkali cleaner (55°C)
 - 2 minutes neutralizing with neutralizer (<40°C)
 - 3 minutes intermediate rinsing with warm water (<40°C)
- Dry the instrument with a validated automatic instrument and surface disinfection machine must be followed.

MANUAL CLEANING / Disinfection
- Manual Cleaning and Disinfection in a washer/disinfector, consider validated instruments in regard to the effectiveness of the cleaning and disinfection process.
- Do not allow blood, tissue, or saline to dry on the instrument.
- Do not use cleaning agents containing fixatives such as bleach or glutaraldehyde.
- Clean one instrument or multiple ones for free range.
- Do not use metal brushes or abrasive powders to clean the instrument.

3. Silicon sterilize the UNFOLDER Emerald Series Cartridge in the sterilization box/pouch at least once and dry (thread end stopper) according to one of the following methods:
- Steam sterilization at least once suitable for wrapped items at 132°C–135°C (270°F–275°F) for 15 minutes.

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