CE2460

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CE2460

Disposable Scalpel with Plastic Handle -IFU

1.Component

Disposable Scalpel with Plastic Handle is surgical blade firmly worked in Plastic Handle of medical application, applied to incision soft tissues on human body.

2.Configuration

Disposable Scalpel with Plastic Handle is made of stainless steel and carbon steel material applied to medical with sharpness and tenacity through punch, heat treatment ,grind and burnish process .etc. Disposable Scalpel with Plastic Handle are assembled on the surgical Plastic Handles for the surgeons to cut the soft tissues. After packaged, the products shall be sterilized with Gamma ray for disposable clinical use.

3 . Capability :

a) Blade sharpness : Incision force ≤0.8N

b) Hardness :Carbon steel ≥800HV10 Stainless steel ≥700HV10

c) Sterilization level :SAL=10^{-6}

4. Specification

9,10,11,12,12d,13,14,15,15c,16,18,19,20,21,22,23,24,24d,25,26,27,34,36,36d
Note：In addition ,our company also developed many special blade for different demands in clinic.

5. Plate .

Interpretation of the symbols used for packing:

	Do not re-use	STERILE R	Sterilization using irradiation
	Use by		Consult instructions for use
	Do not use if package damaged		Do not resterilize
EC REP	Authorized Representative of European		Manufacturer
	Date of manufacture	LOT	Batch code
	Keep away from sunlight		Keep dry
	Fragile. Handle with care	CE2460	CE mark and Identity NO. Of Notified Body.

6. Notice :

a) Users should be familiar with surgical procedures and techniques. Please choose proper specifications surgical Scalpel with Plastic Handle according to the demand of surgery.

b) Please tear open according to the indicated directions of the words on the single and middle packing, then take out the rolled scalpel.

c) It have been sterilized with single package not be re-sterilized. Reuse of the Sterile scalpel with plastic handle cause the following situation during surgery, can't be applied normally、the incision time extends, and for the patient, more risks after surgery, like fever, infection, etc.

d) Do not use if inner package damaged, or there'll be risk of fever, infection, etc.

Tel：+86-517-83806532 Fax：+86-517-83806511 PC：223003
A an, Jiangsu,China. Web：http://www.hy-medical.com

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e) The castoff after operation would be dispose in order to prohibited hurt .

f) Surgical Scalpel with Plastic Handle after used should be discarded into separated 'sharps' container to prevent edge damage.

g) Dispose of contaminated devices and packaging materials utilizing standard hospital procedures and universal precautions for biological hazardous waste.

h) Do not use if the blade has rust .

i) The shelf life is 5 years .

7. Storage

a) Disposable Scalpel with Plastic Handle should be prevented pressure , sunshine and rain .

b) This product should be stored in clean and well ventilated places with suggested relative humidity not higher than 80% and storage temperature between 5℃and 25 ℃; avoiding direct sunlight or heat sources and free of erosive gases.

c) In transportation process, Should prevent break and humid .

8. Effective date

This edition of the User's Manual takes effect on 2017.11.01 ED.:D/2

: Huaiyin Medical Instruments Co., Ltd.
Address：No.8 West Ming Yuan Road, 223003，Huaian, China
Tel：0517-83806499
Fax：0517-83806511
ZIP：223003
E-mail:hy-medical@hy-medical.com
website：http://hy-medical.com

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CE2460

Sterile Surgical Blade -Introduction for Use

1.Component

The sterile surgical blade needs to be used together with plastic handle to cut tissues in surgeries. It can also be used in human anatomy teaching.

2.Configuration

The surgical blade used in our products is made of stainless steel or carbon steel through many processes such as blanking , heat treatment, grind and polish, etc. The blade is sharp. The product adopts the process of sealed package of polyester aluminum complex material .The product is sterilized by gamma irradiation.

3. Performance :

a) Blade sharpness : Cutting force ≤0.8N
b) Hardness :Carbon steel ≥800HV10; Stainless steel ≥700HV10
c) Sterilization level : SAL=10^{-6}

4. Specification

9,10,11,12,12d,13,14,15,15c,16,18,19,20,21,22,23,24,24d,25,26,27,34,36,36d
Note：In addition, our company also developed many special blade for different demands in clinic.

5. Illustration .

Interpretation of the packing symbols:

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6. Notice :

a) Users should be familiar with surgical procedures and techniques.Please choose proper surgical blade according to the demand of surgery.
b) Please tear open the package according to the indicated directions of the words on the single packing and middle packing, then take out the rolled blade.
c) Do not re-sterilize! And do not re-or-cross use! Reuse will increase the risks for the patient after surgery, like fever, infection, etc.
d) Do not use if inner package damaged, or there' ll be risk of fever, infection .
e) The blade shall not be used when the pack is damaged.
f) Surgical blade after used should be discarded into separated 'sharps' container to prevent edge damage.
g) Dispose of contaminated devices and packaging materials utilizing standard hospital procedures and universal precautions for biological hazardous waste,to prevent environmental pollution.
h) Do not use if the blade has rust .
I) The shelf life is 5 years .

7. Storage

a) Surgical blade should be prevented pressure , sunshine and rain .

Tel：+86-517-83806532 Fax：+86-517-83806511 PC：223003
ian, Jiangsu,China. Web：http://www.hy-medical.com

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b) This product should be stored in clean and well ventilated places with suggested relative humidity not higher than 80% and storage temperature between 5℃ and 25 ℃; avoiding direct sunlight or heat sources and free of erosive gases.
c) During transportation, damage and humidity should be avoided .

8. Effective date

This edition of the User's Manual takes effect on 2019.01.20 ED.: D/3

: Huaiyin Medical Instruments Co.,Ltd.
Address: No.8 West Ming Yuan Road, 223003, Huaian,China
Tel: 0517-83806499 Fax: 0517-83806511 ZIP: 223003
E-mail:hy-medical@hy-medical.com website: http://hy-medical.com

EC REP Company Name : Wellkang Ltd
Company Address: Suite B 29 Harley Street
LONDON W1G9QR, England, United Kingdom
Tel: +44 (20) 30869438,32876300
Fox: +44 (20) 76811874

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