
ROTEM® *delta*

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©2002, 2005, 2007, 2008, 2009, 2010, 2012, 2013, 2014, 2016 Tem Innovations GmbH. All rights reserved.

No part of this manual may be reproduced in any way nor edited, copied or distributed electronically (by printing, photocopying or any other method) without written consent of the publisher.

ROTEM® is a trademark of Tem Innovations GmbH in Europe, in the USA and in other countries.

Other product names marked with ® are trademarks of their trademark owners.

Revision History

Manual Version	Software Version	Date	Changes
1.5.3.01.EN	1.5.3	08-2008	first version
1.6.0.01.EN	1.6.0	04-2009	actualization with new SW
1.6.0.02.EN	1.6.0	11-2010	actualization (including new company name and logo) and update with "important note EN_091013"
2.2.0.01.EN	2.2.0	05-2012	actualization with new software version
2.6.0.01.EN	2.6.0	06-2014	actualization with new software version
2.6.0.02.EN	2.6.0 or higher	11-2014	actualization with new optimal range for the service amplitude
2.8.0.01.EN	2.8.0 or higher	08-2016	actualization with new software version

Manual Version: 2.8.0.01.EN
Software Version: 2.8.0
Release Date: 08/2016
Reference Number: REF 200137-EN

This User Manual is published by the manufacturer.
Tem Innovations GmbH
Martin-Kollar-Str. 13-15
D-81829 Munich/Germany

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This manual reflects the technical specifications at the time of printing. Major amendments will be included in a new issue of the Operating Manual. Any data sheets included with this product represent the latest information update.

The right to make amendments of the technical development without explicit reference is reserved.

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3.1 Introduction 3-1

3.2 The ROTEM® measuring principle 3-1

3.3 Reaction curve and parameters in thromboelastometry 3-3

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5 Measurement 5-1

5.1 Introduction 5-2

5.1.1 Usage of the pipette 5-2
5.1.2 Sample handling 5-3

5.2 Measurement preparation 5-4

5.2.1 Switching on the ROTEM® system 5-4
5.2.2 Reagent preparation 5-5
5.2.3 Refill cups and tips 5-6
5.2.4 Measuring cell preparation 5-6
5.2.5 Sample preparation 5-7
5.2.6 Adjust measurement temperature at the ROTEM® 5-7
5.2.7 Select Channel 5-8
5.2.8 Select test 5-8
5.2.9 Enter patient data 5-8

5.3 Performing measurements 5-10

5.3.1 Start measurement 5-11
5.3.2 Graphics of measurement results 5-12
5.3.3 Warning messages during measurement 5-13
5.3.4 Interpretation of measurement results 5-14

5.4 Finish measurement 5-14

5.4.1 Stop measurement 5-14
5.4.2 Print measurement results 5-14
5.4.3 Save and clear measurement 5-14
5.4.4 Discard measurement 5-15
5.4.5 Remove measuring cells 5-15
5.4.6 Log off the system 5-16
5.4.7 Daily maintenance 5-16
5.4.8 Switch off 5-16

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5.5 Quality control 5-17
5.5.1 Control intervals 5-17
5.5.2 Continuous self-monitoring when powered on 5-18
5.5.3 Preparation of quality control material 5-18
5.5.4 Performance of the quality control 5-19
5.5.5 Check quality control results 5-20
5.5.6 Suggested actions for QC out of range 5-21
5.6 Managing measurement results in the database 5-22
5.6.1 Find a patient / QC result 5-22
5.6.2 Print patient / QC results 5-24
5.6.3 Export patient or QC results to USB stick 5-24
5.6.4 Backup 5-25
5.7 Troubleshooting 5-26
5.7.1 Warning screen errors 5-26
5.7.2 Errors during pipetting 5-27
5.7.3 Incorrect TEMograms 5-27
6 System software 6-1
6.1 Operating elements 6-1
6.2 Measurement module 6-2
6.2.1 Preparation 6-2
6.2.2 Measurement 6-3
6.2.3 Multi-TEM 6-4
6.2.4 Screenshots 6-4
6.2.5 Overlays 6-4
6.3 Database 6-4
6.3.1 Database user interface 6-4
6.3.2 Manual selection of patient / QC results 6-6
6.3.3 Selecting measurements (patients or QC) using the search form 6-6
6.3.4 Making measurements (patients or QC) available for processing 6-7
6.3.5 Printing 6-7
6.3.6 Exporting 6-7
6.3.7 Evaluating patient or QC results 6-10
6.4 Setup 6-14
6.4.1 User Setup 6-15
6.4.2 Online setup 6-16
6.4.3 Remote access 6-18
6.4.4 Printer setup 6-20
6.4.5 Test setup 6-21
6.4.6 Parameter setup 6-24

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6.4.7 Comment setup 6-24
6.4.8 QC setup 6-25
6.5 Service 6-28
6.5.1 Information 6-28
6.5.2 Settings 6-28
6.5.3 Log book 6-30
6.6 Liquitrans 6-31
6.7 File system 6-32
6.7.1 Admin Tools 6-32
6.7.2 User 6-37
6.8 USB stick 6-38
6.8.1 Update 6-38
6.8.2 Importing data from the USB stick into Excel® 6-39
7 Initiation and installation 7-1
7.1 Initiation 7-1
7.1.1 Open package 7-1
7.1.2 Touch screen preparation 7-1
7.1.3 Site installation requirements 7-2
7.2 Installation 7-3
8 Maintenance, cleaning, disinfection 8-1
8.1 Introduction 8-1
8.2 Maintenance plan 8-2
8.2.1 Daily maintenance 8-3
8.2.2 Weekly maintenance 8-3
8.2.3 Quarterly maintenance 8-3
8.2.4 Semi-annual maintenance 8-7
8.2.5 Yearly maintenance 8-9
8.2.6 Every 3 years 8-10
8.3 Maintenance of the ROTEM® electronic eLine pipette 8-10
8.3.1 Replacement of the pipette filter 8-10
8.3.2 Cleaning of the pipette surface 8-11
8.3.3 Cleaning and disinfection of individual parts 8-11
8.3.4 Check and calibrate ROTEM® eLine pipette 8-14
8.3.5 Functional check of the ROTEM® eLine pipette 8-14
8.3.6 New calibration 8-14
8.3.7 Troubleshooting 8-15
9 Placing out of operation, transport and disposal 9-1

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9.1 Placing out of operation 9-1
9.2 Transport 9-1
9.3 Disposal 9-1
10 Appendix 10-1
11 Index 11-1

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1 Introduction

This operating manual applies to the ROTEM® *delta* system with the respective ROTEM® software. It is part of the product and describes the intended use and safe exposure during the whole product life cycle.

In this document the terms ROTEM® and ROTEM® system always refer to ROTEM® *delta*.

1.1 Subject to change and validity

The information in this manual concerns the approved technical specifications at the time of printing. Significant amendments will be considered in a new edition of the manual. Any data sheets accompanying this manual contain the most current information regarding the product. The document and version number of this manual are mentioned in thefooter.

1.2 Completeness

This operationd manual is only complete with further applicable documents for the respective application.

1.3 Depository

This manual as well as further applicable documents for the respective application have to be kept within reach and always accessible near by the ROTEM® system.

1.4 Further applicable documents

A range of system reagents and consumables are available for the ROTEM® *delta* Thromboelastometry system. Always refer to the instructions for use (IFU = *Instruction for use*) for the reagents and control kits regarding the most recent information concerning specific sample requirements, execution of tests, limits of the method and of interpretation of the results.

Furthermore please refer to the universally valid, legal or other binding regulations of european and national law, as well as to the effective regulations for prevention of accidents and environmental protection applicable in your country.

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1.5 Typographical conventions

This section shows the applied abbrevations, symbols and textual accentuations.

1.5.1 Used abbrevations

Table 1-1 shows used abbrevations in this document.

Abbreviation	Signification
A(x)	Amplitude after (x) minutes
ACF	Actual Clot Firmness
AIDS	Acquired Immune Deficiency Syndrome
AR(x)	ARea under the curve at time (x)
AUC	Area under the Curve (from the first derivative)
CFR	Clot Formation Rate
CFT	Clot Formation Time
CLR	Clot Lysis Rate (Lyse Rate)
CT	Clotting Time
DB	Database
G	Elastic shear modulus
HBsAg	Hepatitis B Surface Antigen
HCV	Hepatitis C Virus
HIV	Humane Immunodeficiency Virus
ID	Identification (e.g.: patient ID, sample ID, etc.)
IFU	Instruction For Use
INR	International NormaliSed Ratio
IVD	In Vitro Diagnostic
LAN	Local Area Network
LI(x)	Lysis Index (after (x) minutes)
LIS	Laboratory Information System
LOT	Lysis Onset Time
LT	Lysis Time
MA	Maximum Amplitude
maxVel	Maximum Velocity
MCE	Maximum Clot Elasticity
MCF	Maximum Clot Firmness
ML	Maximum Lysis
N	Normal (for QC or blood samples)
P	Pathological (for QC or blood samples)
POC	Point Of Care
QC	Quality Control
REF	Reference Number
RO	Read Only
RW	Read and Write
SN	Serial Number
TEM	TEMogram
t-MaxVel	Time to Maximum Velocity
TPI	Thrombodynamic Potential Index
VNC	Virtual Network Computing
∠	Alpha angle

Table 1-1: Used abbrevations

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1.5.2 Text conventions

Table 1-2 shows applied accentuations in this document.

Indentification	Applied for	Example
<>	Keyboard entries	<Enter>
Small capitals	Surface texts	DATABASE
Bold	Buttons	**Exit**
Italic	Input field	*Patient ID*

Table 1-2: Text conventions

1.5.3 Warning concept

Warning messages in this document are constructed as follows:

⚠ SIGNAL WORD

Danger

Consequences

➥ Precautions

According pictogram is shown if it is a matter of biohazard (Table 1-3).

Symbol	Applied in	Relevance
	Warning message	Refers to danger. The level of danger is labelled by the signal word.
	Warning message	Refers to biohazard.

Table 1-3: Applied symbols in the warning concept

1.5.4 Strategic concept

Operation procedure instructions in this document are constructed as follows:

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1.6 Help function

A help file is available in the MEASUREMENT MODULE and in the modules DATABASE, SETUP, SERVICE and LIQUITRANS. It explains, additional to this document, the handling of the equipment, handling of the reagents and in particular the interpretation of the results. A troubleshooting guide is also available in this menu.

1.7 Target group

This document is adressed to qualified hospital staff, working in the field of haemostasis.

- Physicians (mostly anaesthesists and laboratory physicians)
- Laboratory technicians
- Perfusionists and others

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2 Safety

2.1 General safety instructions

This manual contains important information for the safe and appropriate installation, operation, transport, storage, and maintenance of the ROTEM®system.

- Please read the operating manual carefully to acquaint yourself with the equipment.
- Special attention should be paid to the information provided in this chapter.
- Please note the instructions mentioned in Appendix A regarding the applied symbols for the equipment and reagents.
- Before performing measurements and quality controls, please read the instructions of the reagents and controls.

2.2 Intended use

ROTEM® Thromboelastometry systems are designed for *in vitro* diagnostics in the POC (Point of Care) field and in hospital laboratories. The ROTEM® system is intended to be used to provide a quantitative and qualitative indication of the coagulation state of a blood sample. Different parameters of coagulations are measured, analysed, monitored and interpreted (constitution of clot) for this purpose.

Additionally, the ROTEM® system records the kinetic changes in a sample of citrated whole blood as the sample clot retracts and/or lyses (breaks apart).

Results from the ROTEM® system should not be the sole basis for a patient diagnosis; ROTEM® results should be considered along with a clinical assessment of the patient's condition and - if necessary - other laboratory tests. The equipment is not evaluated to measure veterinary or other non-human samples.

Non-intended usage means all other usage of the ROTEM® *delta* than described in this manual.

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2.3 Ambient conditions

For the proper operating of the ROTEM® *delta* please note the requirements within the work environment as follows:

Operating temperature	+15 °C to +30 °C
Storage temperature	+0 °C to +50 °C
Relative humidity	20% to 85%
Pollution degree	2
Air pressure (Operating, Warehousing, Transport)	700 hPa to 1060 hPa
Operable to height above sea level	2000 m
Distance of examined, treated or monitored patients	1.5 m

Table 2-1: Ambient conditions

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2.4 Safety standards

The in-vitro diagnostic equipment of the ROTEM® category has been manufactured in compliance with EC Directive 98/79/EC of October 27, 1998 (including Corrigendum to Directive 98/79/EC) regarding in-vitro diagnostic equipment (IVD guideline).

The equipment has been manufactured in accordance with the following standards:

IEC 61010-1:2010 Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1: General requirements (EN 61010-1:2010)

IEC 61010-2-081:2001 + A1:2003 Safety requirements for electrical equipment for measurement, control and laboratory use. Part 2-081: Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes

IEC 61010-2-101:2002 Safety requirements for electrical equipment for measurement, control and laboratory use. Part 2-101: Particular requirements for *in vitro* diagnostic (IVD) medical equipment (EN 61010-2-101:2002)

IEC 60601-1-2:2007 (Modified) Medical electrical equipment. Part1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests excluding chapter 5: Markings, Inscriptions and Accompanying Documents (EN 60601-1-2:2007)

IEC 61326-1:2012 Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements (EN 61326-1:2013)

IEC 61326-2-6:2005 Electrical equipment for measurement, control and laboratory use - EMC requirements. Part 2-6: Particular requirements - *In vitro* diagnostic (IVD) medical equipment (EN 61326-2-6:2006)

The equipment uses RF energy only for its internal function. Therefore, its RF-emissions are very low and are not likely to cause any interference in nearby electronic equipment.

This equipment is suitable for use in all establishments other than domestic.

The equipment fulfils this safety requirements when it is leaving the manufacturing site[1].

[1] See Declaration of conformity (Appendix D)

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2.5 Personal protection

2.5.1 Biohazards

Always treat test samples as a potential biohazard and ensure proper handling and disposal of the samples.

- Always use personal protective equipment, such as laboratory coat, gloves and goggles when handling samples. Adhere to local statutory provisions for proper disposal (general recommendation, available in *Laboratory Biosafety Manual World Health Organization, 1984*).

2.5.2 Electrical hazard

⚠ WARNING

Danger of life

Misuse of electrical equipment can cause electrocution or death by burning.

- Always adhere to basic safety precautions for electrical equipment.
- Don’t touch the patient and the device at the same time.
- Any adjustments, service or repairs that have to be carried out while the equipment is energised, may only be performed by qualified service personnel for ROTEM® systems.
- Do not leave the equipment unattended while it is plugged in.
- The equipment is supplied with an external grounded power cable. Connect the equipment only to an grounded power socket. Utilised adaptors have to be grounded, equally.
- Keep equipment dry. If equipment gets wet, unplug it before touching.
- Utilise only original accessories (Appendix G).
- Do not switch the equipment on if it is damaged. Contact qualified service personnel.
- If the power cable (instrument plug) is lost or damaged, it must be replaced by an original power cable (Appendix G).

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2.6 Residual risks

Adhere to following instructions for the safe usage of the ROTEM® *delta*:

- Do not use the equipment if it is not working properly or damaged.
- Do not leave the equipment unattended if it is used near children or disabled persons.
- Use the equipment only for purposes described in this operating manual (Chapter 2.2).
- Liquids may damage the internal components of the equipment. Do not spill liquids over or into the ROTEM® system. Do not immerse the equipment in water or other liquids.
- Only qualified service personnel may access internal components of the ROTEM® systems.
- If safe operation is not guaranteed, usage of the ROTEM® system has to be discontinued and further usage has to be prevented.
- Do not let the equipment or its power cable come into contact with surfaces that are too warm to touch.
- Do not cover air openings nor place equipment on a soft surface which might block them. Keep air openings free from lint, hair, dust etc.
- Do not place anything on top of the equipment, except the accessories recommended by the manufacturer.
- Unless specified in the operating manual, do not stuff or put anything into the openings, tubes or couplings of the equipment.
- Do not use the ROTEM® *delta* near flammable or explosive atmospheres.
- The equipment is designed for indoor use only.
- Only software explicitly supplied by Tem Innovations GmbH may be installed on the ROTEM® *delta* system.
- Do not drop the equipment or expose to shocks.
- Do not touch or stop the rotating axis with your hand.

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3 Basics

3.1 Introduction

The ROTEM® system is a further development of classical thromboelastography developed by Hartert.

The ROTEM® *delta* Thromboelastometry systems are designed for *in vitro* diagnostics in the POC *(Point of Care)* field or in hospital laboratories. The system is intended to provide a quantitative and qualitative indication of the coagulation state of a blood sample. The system records the kinetic changes in a sample of citrated whole blood during clot formation as well as when the sample clot retracts and/or lyses (breaks apart). Different parameters of the clotting are measured, analysed, monitored, interpreted, and charted for this purpose. The graphical presentation reflects the various physiological results, which describe the interaction between components like coagulation factors and inhibitors, fibrinogen, thrombocytes, and the fibrinolysis system. Additionally, the effect of certain drugs influencing haemostasis, in particular anticoagulants, can be detected.

3.2 The ROTEM® measuring principle

The patented ROTEM® technology is based on a fixed cylindrical cup and a permanently oscillating vertical axis (Figure 3-1).

The axis is supported by a high precision ball bearing and oscillates to the left and to the right through an angle of 4.75 °. The rotation of the axis is driven by a motor that is connected to the axis via an elastic spring.

For the measurement, a disposable plastic pin with 6 mm diameter is placed firmly on the axis and the blood sample is filled into a disposable 8 mm diameter cup and is then uplifted onto the measurement channel. Hence, the plastic pin is immersed into the blood sample.

The rotation is detected optically via a mirror plate at the upper end of the axis, a diode as light source and a light sensitive sensor (CCD Chip). If no clotting takes place, the movement is not obstructed. When a clot is formed and attaches itself between pin and cup surfaces, the movement is obstructed.

The result is a balance between the spring tension and the tension of the clot. As the clot becomes firmer, the rotational amplitude of the axis is reduced.

The results of the measurement are interpreted with special software (Chapter 6).

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
Figure 3-1: Principle of thromboelastometry with ROTEM® *delta*

1 Axis (~4.75 °)
2 Spring
3 Light source/diode
4 Mirror
5 Detection device (electric camera)
6 Sensor pin
7 Cup filled with blood
8 Fibrin fibres and thrombocyte aggregate
9 Heated cup holder
10 Ball bearings
11 Data processing

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3.3 Reaction curve and parameters in thromboelastometry

Like in classical thrombolastography, the ROTEM®system generates a reaction curve and calculates different numeric (kinetic and firmness) parameters in a mathematical analysis of that curve.

3.3.1 Reaction curve

The ROTEM® software utilises a highly developed curve smoothing algorithm and noise filter which prevents potential errors due to mechanical or electronic noise. The parameters are determined in real time during the tests, calculated, and are represented graphically in TEMograms.

The rotation amplitude of the pin is converted into a graphical amplitude, thus the following definition applies for the ROTEM® system:

Graphical amplitude in mm	Rotation of the pin	Clotting
0	free rotation of the pin	no clotting
100	no rotation of the pin	maximum possible firmness of the clot

Table 3-1: Definition of the amplitude

Figure 3-2 shows the most important phases of clotting with the according parameters in a TEMogram.


Figure 3-2: Example of a ROTEM® *delta* reaction curve (TEMogram)

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3.3.2 The most important routine parameters

Generally, the measurement results are interpreted with routine parameters. Figure 3-3 shows the most important routine parameters. The y-axis shows the amplitude in mm, the x-axis shows the time in minutes.


Figure 3-3: The most important routine parameters

The routine parameters are described in the following sections, in Table 3-9 a compact overview is given.

Clotting Time, CT, [s]

Definition	The CT is the time from the beginning of the test by adding the clotting activator until the time when an amplitude of 2 mm is achieved.
Description	The CT describes how fast the formation of fibrin starts. This parameter is analogical to the clotting time in a classical clotting test in a laboratory. However, they are not identical, as more fibrin is created and has to be stabilised in order to achieve a certain firmness of the clot that is sufficient to connect the two moving parts of the measuring cell.
Main influencing factors	Clotting factors, anticoagulants (sensitivity is dependent on the test)
Clinical application	The CT parameter facilitates the decision to substitute clotting factors (e. g. through fresh frozen plasma, factor concentrates, activated factor concentrates or anticoagulant inhibitors (e. g. Protamin)).

Table 3-2: Routine parameter – CT

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Clot Formation Time, CFT, [s]

Definition	The CFT is the time between a 2 mm amplitude and a 20 mm amplitude of the clotting signal.
Description	The CFT describes the next phase of the clotting: the kinetics of the formation of a stable clot through both activated thrombocytes and fibrin.
Main influencing factors	The amount of thrombocytes and their contribution to the clot firmness. Fibrinogen level and its ability to polymerise.
Clinical application	The CFT parameter facilitates the decision to substitute with thrombocyte concentrate or fibrinogen (as cryoprecipitate, fresh frozen plasma, fibrinogen concentrate) or both. A shortened CFT is indicative for hypercoagulation (as well as a high MCF parameter and alpha angle). Samples with very low clot formation, the clot formation time may not be achieved and therefore not be indicated.

Table 3-3: Routine parameter – CFT

Alpha angle (α, [°])

Definition	The alpha angle is defined as the angle between the middle axis and the tangent to the clotting curve through the 2 mm amplitude point. It describes the kinetic of clotting.
Description	The diagnostic information of this parameter is similar to CFT.
Clinical application	A reduced alpha angle indicates a hypocoagulable state (refer to CFT).

Table 3-4: Routine parameter – Alpha angle

Maximum Clot Firmness, MCF, [mm]

Definition	The MCF is the measure for the firmness of the clot and therefore the clot quality. It is the maximum amplitude that is reached before the clot is dissolved by fibrinolysis and the clot firmness falls again.
Main influencing factors	Thrombocytes, fibrinogen (concentration and the ability to polymerise),F XIII, the status of fibrinolysis.
Clinical application	A low MCF indicates a low clot firmness and therefore a potential bleeding risk. The MCF value is used to facilitate the decision for substitution therapy with thrombocyte concentrate or fibrinogen (concentrate, cyroprecipitate or fresh frozen plasma, depending on availability). Hyperfibrinolysis should necessarily be ruled out before treating with a fibrinogen source, because hyperfibrinolysis may lead to an unstable clot. A high MCF value may indicate a hypercoagulable state.

Table 3-5: Routine parameter – MCF

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A(x)-values ([mm])

Definition	A(x)-values represent the clot firmness. An A(x)-value is the amplitude after a certain time x after CT (e.g. A10 after 10 min).
Main influencing factors	Thrombocytes, fibrinogen (concentration, ability to polymerise), F XIII.
Clinical application	see MCF.

Table 3-6: Routine parameter – A(x)-values

Lysis index at 30 min (LI30) and related parameters ([%])

Definition	The LI30 value represents the fibrinolysis 30 min after CT. It is the relation of the amplitude to the maximum clot firmness (% remaining clot firmness). The LI45 and LI60 parameters describe the according remaining clot firmness 45 and 60 min after CT.
Clinical application	Due to the high concentration of fibrinolysis inhibitors, almost no fibrinolysis can be observed in samples from healthy persons. An abnormal LI30 value mostly indicates hyperfibrinolysis. Therefore the LI30 parameter may facilitate the decision for or against a therapy with anti-fibrinolytic drugs. In certain cases, hyperfibrinolysis may develop relatively late. In such cases, also LI 45 or LI60 may be used for a decision.

Table 3-7: Routine parameter – LI

Maximum lysis (ML, [%])

Definition	The parameter of maximum lysis (ML) describes the degree of fibrinolysis relative to maximum clot firmness (MCF) achieved during the measurement. (% clot firmness lost).
Description	A ML of 5% means that, at the period of observation, the MCF has decreased by 5%. As the maximum lysis is not calculated at a fixed time point, but is defined as % lysis at the end of the measurement, the total runtime and the time after maximum clot formation should always be considered.
Clinical application	see LI(x)

Table 3-8: Routine parameter – ML

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3.3.3 Research parameters

An important extension of classical thrombolastography has been introduced by Sørensen et al. in 2003[1]. This group has successfully applied a new method of data analysis to the ROTEM® system, primarily focusing on the rate of coagulation via derivative curves (Figure 3-4). The research parameters described by Sørensen *et al*., MaxVel (MAXV), t,MaxVel (MAXV-t) and AUC, have been implemented into the ROTEM® *delta* software and can be used for research applications.

Note: Research parameters are not validated for the routine application in the patient diagnostics (with the exception of MCE).


Figure 3-4: TEMogram and 1st derivative curve

[1] **Literature:** Sørensen B., Johansen P., Christiansen K., Woelke M., Ingerslev J.: Whole blood coagulation thrombelastographic profiles employing minimal tissue factor activation. J Thromb Haemost 2003;1:551-8

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The lower part of Table 3-9 gives an overview of the research parameters. A summary of the routine parameters described in Chapter 3.3.2 is shown in the grey shaded area.

Acronym	Parameter	Definition	Unit
ROUTINE PARAMETERS			
Coagulation activation and clot polymerisation parameters			
CT	Coagulation Time (synonym r)	The time from test start until an amplitude of 2 mm is reached	s
CFT	Clot Formation Time (synonym k)	The time between 2 mm amplitude and 20 mm amplitude	s
α	α-Angle	Angle between the baseline and a tangent to the clotting curve through the 2 mm point	degree (°)
Clot firmness parameters			
MCF	Maximum Clot Firmness (synonym MA)	The maximum amplitude reached during the test	mm
A(x)	Amplitude (firmness) at time x	Clot firmness (in mm amplitude) at the respective time point after CT	mm
A5	Firmness at time 5 minutes		
A10	Firmness at time 10 minutes		
A15	Firmness at time 15 minutes		
A20	Firmness at time 20 minutes		
A25	Firmness at time 25 minutes		
A30	Firmness at time 30 minutes		
Clot lysis parameters			
ML	Maximum Lysis	Maximum lysis detected during the run time, described as the difference between MCF and the lowest amplitude, in % of MCF (% lost firmness)	%
LI(x)	Lysis Index at time x minutes	Ratio of the amplitude and MCF at a given time point after CT (% firmness remaining) Calculated: A/MCF*100	%
RESEARCH PARAMETERS			
Research parameters for derivative curves (according to Sørensen et al. 2003)			
maxV	Maximum velocity	The maximum of the 1st derivative of the curve	mm/min
maxV-t	Time to maximum velocity	Time from reaction start until the maximum of the 1st derivative of the curve is reached	s
AUC	Area under 1st derivative curve	Area under the 1st derivative curve from start of the derivative curve until MCF is reached	mm x 100
Other research parameters			
CFR	Clot Formation Rate	Angle between the baseline and the tangent at the maximum slope	degree (°)
ACF	Actual Clot Firmness or Last Clot Firmness	Clot firmness at the actual time point of observation	mm
MCF-t	MCF-Time	The time from CT until MCF is reached	s
MCE	Maximum Clot Elasticity MCE=100*MCF/(100-MCF)	MCE is a parameter calculated from MCF. The spreading (solution) of this parameter at high amplitudes is better as compared to MCF	
AR(x)	Area under the curve at time x	The area under the curve from CT to the respective time point (minutes)	mm^2
LOT	Lysis Onset Time	The time span from CT to the start of significant lysis. Significant lysis is defined as a decrease of the amplitude of 15% as compared to MCF	s
CLR	Clot Lysis Rate	The strongest lysis, described by the angle between the baseline and the tangent to the declining firmness curve at the minimum of the 1st derivative	degree (°)
LT	Lysis Time	The time from CT until the clot firmness is decreased to 10% of the MCF	s
G	G = 5000MCF/(100 – MCF) Shear Elastic Modulus Strength	G is a calculated parameter. It increases exponentially as compared to the amplitude. It allows a more sensitive resolution at high amplitudes as compared to MCF	
TPI	Thrombodynamic Potential Index TPI = EMX/K EMX = (100*MCF)/(100-MCF)	According to Raby 1975, the Thrombodynamic Potential Index describes the patient's global coagulation	

Routine parameters are displayed in grey

Table 3-9: Routine and research parameter

3.3.4 Reference ranges

Reaction curves and the kinetic parameters of each sample of a patient can be compared with the results of healthy reference persons.

The reference values are described in a multicenter study for all ROTEM® system tests[1]. Results for each of the test reagents are indicated in the respective instructions for use. However, pre-analytical factors (sampling, storage period and physical stress such as shaking of the tubes) may influence the results. These reference values should be regarded as a guideline and should be checked in each hospital on the basis of an appropriate reference population.

[1] Literatur: Lang T., Bauters A., Braun S.L., Potzsch B., von Pape K.W., Kolde H.J., Lakner M.: Multi-centre investigation on reference ranges for ROTEM® thromboelastometry. Blood Coagul Fibrinolysis 2005;16:301–310.

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3.4 ROTEM® differential analysis

The ROTEM® analysis is a diagnostic approach which aims to detect all relevant alterations in haemostasis which may have an impact on coagulation rapidly. It uses thromboelastometry of citrated blood in combination with a variety of specific reagents.

In individual cases, the curve shape of the TEMogram shows characteristics that make an immediate diagnosis possible (e. g. typical spindle shape in the case of hyperfibrinolysis). Also a hypercoagulative state can be characterised by a typical curve shape with a rapid start of clotting, an explosive clot formation and a rapidly achieved maximum firmness.

In the case of complex haemostasis disorders, only a comparison between reference values and the results obtained through the application of specific reagents, can result in a better understanding of haemostasis disorders. This differential diagnostic procedure results in targeted therapy.

The ROTEM® software HELP menu contains a clear tutorial for differential diagnosis (Chapter 1.6).

3.4.1 The citrated blood sample

In classical thromboelastography, non-anticoagulated whole blood or citrated blood (after recalcification) was used. However, non anticoagulated blood is unstable and may be used only for a limited period after sampling. The measurement results may vary rapidly and there is a risk of producing artefacts. Citrated whole blood has the advantage that it is more stable than non-anticoagulated whole blood, especially as platelet activity can change rapidly over time. That is why the ROTEM® tests are performed with citrated whole blood (3.2 or 3.8% citrate, same tube as for coagulation tests in the laboratory). The pipetting programme of the ROTEM® system is designed for citrated whole blood.

Sample stability

The stability of the sample varies according to the tests. Generally, within 4 hours no influence on the measurement results of activated tests is noticed, with samples from healthy people[1]. The instructions for use of each test contain detailed information.

Note: Normally, the sample stability is reduced for non activated thromboelastometry (30 min), hyperfibrinolytic samples and samples from patients under heparin therapy.

Influence of the sample temperature

The temperature of a fresh blood sample adopts the room temperature. At a room temperature of 23 °C, the temperature decreases within 5 min to 30 °C and within 20 min to 25 °C.

Experiments have shown that the impact of temperature on the measurement results of the most ROTEM® parameters is quite small as long as the ambient temperature is >20 °C. However, in order to get reproducible and precise results, especially for CT and CFT, the sample should preferably always have the same temperature.

[1] **Literature:** Bohner J., von Pape K.-W.: Method Dependent Reference Values and Pre-analytical Influences in Rotation Thrombelastography. Ann Haematol 2003;82(S1):13.

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For this purpose, the ROTEM® *delta* disposes of a temperature controlled work area (incl. a sample preheating station). Samples of patients that are not analysed immediately after sampling, have to be preheated 5–10 min (and not longer than 30 min) in the sample preheating station before analysing.

Also, control samples (ROTROL N or ROTROL P) must be kept on the temperature controlled work area for 5 min before measurement starts.

Table 3-10 shows different values for the CT and the CFT for EXTEM when not adhering to this instruction. ROTROL N was transferred from a vial at the indicated temperature to a cup in a cup holder at 37 °C. The measurement with the EXTEM test was started immediately.

Temperature	CT	CFT
20 °C	47 s	42 s
25 °C	43 s	34 s
30 °C	37 s	29 s
37 °C	33 s	28 s

Table 3-10: Measurement results for the CT and the CFT at different sample temperatures.

3.4.2 ROTEM® tests

Introduction

Acute bleeding during or after surgery requires rapid differentiation between surgical induced bleeding and haemostasis disorders. The combination of the ROTEM® *delta* with additional diagnostic methods, considering the given limitations (platelet aggregation is not measured by ROTEM® *delta* for example, but needs to be measured with ROTEM® *platelet*), facilitates further differential treatment strategies.


Figure 3-5: Differentiation of surgical bleeding and haemostasis disorder with the ROTEM® analysis

Contrary to the classical thromboelastography, executed without reagents, the ROTEM® whole blood haemostasis analysis with the specific system reagents enables an extensive diagnostics as basis for therapeutical decisions [1].

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The ROTEM® analysis extends the resulting diagnostic power by a number of additional tests and parameters, which

- shorten the reaction time relevantly,
- increase the precision,
- inhibit certain factors (e.g. heparin),
- allow differentiation between fibrin and platelet contributions to the overall clot.


Figure 3-6: The principle of the ROTEM® analysis

The ROTEM® tests

You will find detailed information regarding the specific requirements for the samples, limitations, current test characteristics, limits of the method and interpretation of results, in the instructions for use of the ROTEM® system reagents.

Following name conventions are applied:

- The tests are referred to as INTEM(S),EXTEM(S), FIBTEM (S), HEPTEM(S) and APTEM(S).
- The associated reagents are referred to as in-tem®(S), ex-tem®(S) and r ex-tem® , fib-tem®(S), hep-tem®(S), ap-tem®(S) and (t) ap-tem®.

There are two product lines for ROTEM® reagents:

- For frequent users: stable liquid reagents with approx. 5 – 10 tests per vial.
 Those reagents are pipetted using the INTEM, EXTEM etc. pipetting sequences.
- For POC users: Single Use Reagent.
 Those reagents are pipetted using the INTEM S, EXTEM S etc. pipetting sequences.

Note: Always use the appropriate pipetting sequence for the type of reagents. The sequences have to activated in the SETUP menu, so that they are available in the system (Chapter 6.4.5).

[1] **Literature**: Calatzis A., Haas S. et al.: Thromboelastographic coagulation monitoring during cardiovascular surgery with the ROTEG coagulation analyzer. In: Pfarr, R., ed. Management of Bleeding in Cardiovascular Surgery, 2000, Hanley Belfus Inc, Philadelphia.

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The following sections describe the general principles of the available tests in ROTEM® *delta*.

INTEM (S)		
ROTEM® Reagents	*Activator*	*Recalcifier*
Liquid reagents	in-tem®	star-tem®
Single Use Reagent	in-tem® S	
Principle:	mild intrinsic coagulation activation	
Application area:	global analysis of coagulation in combination with HEPTEM (see below)	
Sensitive to:		*Most affected ROTEM® parameters:*
Factor deficiency (intrinsic and common pathways)		CT
Anticoagulant effects (e. g. heparin, direct thrombin inhibitors)		CT
Thrombocyte contribution to clot firmness		MCF, A(x), (CFT, α)
Fibrin polymerisation (and fibrinogen concentration)		MCF, A(x), (CFT, α)
Hyperfibrinolysis		ML, LI(x)
(F-XIII-deficiency)		MCF, A(x), ML, LI(x)
Limitations:	• Not very sensitive to mild coagulation factor deficiencies • Insensitive to defects of primary haemostasis (platelet aggregation)	




Table 3-11: INTEM (S)

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HEPTEM (S)		
ROTEM® Reagents	*Activator*	*Heparinase and Recalcifier*
Liquid reagents	in-tem®	hep-tem®
Single Use Reagent	hep-tem® S	
Principle:	Mild intrinsic coagulation activation in the presence of a heparin degrading enzyme (heparinase I).	
Application area:	• global analysis of coagulation after the elimination of the influence of heparin • In comparison to INTEM: Qualitative assay for the presence of heparin	
Sensitive to:		*Most affected ROTEM® parameters:*
Factor deficiency (intrinsic and common pathways)		CT
Anticoagulants (e.g. direct thrombin inhibitors, not heparin!)		CT
Platelet contribution to clot firmness		MCF, A(x), (CFT, α)
Fibrin polymerisation (and fibrinogen concentration)		MCF, A(x), (CFT, α)
Hyperfibrinolysis		ML, LI(x)
(F-XIII-deficiency)		MCF, A(x), ML, LI(x)
Note:	The combination of INTEM and HEPTEM confirms the existence of heparin in the sample. It also provides information on alterations in the coagulation potential of the sample under elimination of the heparin therapy.	
Limitations:	• Not very sensitive to mild coagulation factor deficiencies • Insensitive to defects of primary haemostasis (platelet aggregation)	

The prolonged CT in the INTEM test is significantly shortened in the HEPTEM test (to normal value). Also CFT, A(x) and MCF parameters return to normal values. The sample therefore contains heparin.

Table 3-12: HEPTEM (S)

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EXTEM (S, L)		
ROTEM®Reagents	*Activator*	*Recalcifier*
Liquid reagents	r ex-tem®, ex-tem® L	star-tem®
Single Use Reagent	ex-tem® S	
Principle:	mild extrinsic coagulation activation	
Application area:	• global analysis of coagulation largely heparin insensitive (does not apply to EXTEM S!) • in conjunction with FIBTEM and APTEM (see below)	
Sensitive to:		*Most affected ROTEM® parameters:*
Factor deficiency (extrinsic and common pathways)		CT
Platelet contribution to clot firmness		MCF, A(x), (CFT, α)
Fibrin polymerisation (and fibrinogen concentration)		MCF, A(x), (CFT, α)
Hyperfibrinolysis		ML, LI(x)
(F-XIII-deficiency)		MCF, A(x), ML, LI(x)
Limitations:	• Not very sensitive to mild coagulation factor deficiencies. • Insensitive to defects of primary haemostasis (platelet aggregation) • May still be normal when INR *(International Normalized Ratio)* is elevated. • May show pathological values caused by very high heparin levels.	



Table 3-13: EXTEM (S, L)

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FIBTEM (S)		
ROTEM®Reagents	*Activator*	*Thrombocyte inhibitor and Recalcifier*
Liquid reagents	r ex-tem® , ex-tem® L	fib-tem®
Single Use Reagent	fib-tem® S	
Principle:	Mild extrinsic coagulation activation in the presence of a thrombocyte inhibitor (cytochalasin D). The TEMogram represents the fibrin part of the clot only.	
Application area:	• Detection of fibrinogen deficiency and fibrin polymerisation disorders • Compared to EXTEM: indirect evaluation of the thrombocyte component of the coagulation • Largely heparin insensitive (does not apply to FIBTEM S!)	
Sensitive to:		Mostly affected ROTEM® parameters:
Fibrin polymerisation (and fibrinogen concentration)		MCF, A(x), (CFT, α)
Hyperfibrinolysis		ML, LI(x)
(F-XIII-deficiency)		MCF, A(x), ML, LI(x)
Limitations:	• May be influenced by very high heparin levels. • FIBTEM represents only the fibrin component of the clot.	
Note:	In contrast to the measurement of only platelet count or fibrinogen concentration, the ROTEM® system captures the ability of fibrin and the cellular components to form a stable clot. The FIBTEM clot firmness generally correlates with the fibrinogen concentration. But a fibrin polymerisation disorder will give a pathological result in spite of a reasonably high fibrinogen concentration in the sample.[1]	



Table 3-14: FIBTEM (S)

[1] **Literature**: Innerhofer P., Fries D., Margreiter J., Klingler A., Kuhbacher G., Wachter B., Oswald E., Salner E., Frischhut B., Schobersberger W.:The effects of perioperatively administered colloids and crystalloids on primary platelet-mediated hemostasis and clot formation. Anesth Analg 2002;95(4):858–65.

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APTEM (S)		
ROTEM® Reagents	*Activator*	*Fibrinolysis inhibitor and Recalcifier*
Liquid reagents	r ex-tem®, ex-tem® L	t ap-tem® or ap-tem®
Single Use Reagent	ap-tem® S	
Principle:	Mild extrinsic coagulation activation in the presence of a fibrinolysis inhibitor (tranexamic acid or aprotinin).	
Application area:	• Confirmation of hyperfibrinolysis in comparison to EXTEM • Assessment of coagulation after hypothetical patient treatment with anti-fibrinolytic drugs • Largely heparin insensitive (does not apply to APTEM S!)	
Sensitive to:		*Most affected ROTEM® parameters:*
Factor deficiency (extrinsic and common pathways)		CT
Platelet contribution to clot firmness		MCF, A(x), (CFT, α)
Fibrin polymerisation (and fibrinogen concentration)		MCF, A(x), (CFT, α)
(F-XIII-deficiency)		MCF, A(x), ML, LI(x)
Note:	If in APTEM, with the occurrence of the typical pattern of hyperfibrinolysis (spindle shaped, total lysis of the clot firmness) in EXTEM (as in INTEM and FIBTEM), hyperfibrinolysis is not present, then hyperfibrinolysis is confirmed. The APTEM result represents *in vitro* the coagulation situation of the patient after a possible therapy with antifibrinolytics. Up to a certain extent , APTEM can identify *in vitro*, if antifibrinolytic therapy will restore the coagulation system of the patient or if further therapy is necessary.	
Limitations:	See EXTEM	


Table 3-15: APTEM (S)

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NATEM		
ROTEM® Reagents	*Activator*	*Recalcifier*
Liquid reagents	none	star-tem®
Principle:	Contact coagulation activation	
Application area:	global analysis of coagulation in combination with NA HEPTEM (see below)	
Sensitive to:		*Most affected ROTEM® parameters:*
Factor deficiency		CT
Anticoagulant effects (e. g. heparin, direct thrombin inhibitors)		CT
Thrombocyte contribution to clot firmness		MCF, A(x), (CFT, α)
Fibrin polymerisation (and fibrinogen concentration)		MCF, A(x), (CFT, α)
Hyperfibrinolysis		ML, LI(x)
(F-XIII-deficiency)		MCF, A(x), ML, LI(x)
Limitations:	• Insensitive to defects of primary haemostasis (platelet aggregation)	



Table 3-16: NATEM

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NA HEPTEM		
ROTEM® Reagents	*Activator*	*Heparinase and Recalcifier*
Liquid reagents	none	hep-tem®
Principle:	Contact coagulation activation in the presence of a heparin degrading enzyme (heparinase I).	
Application area:	• global analysis of coagulation after the elimination of the influence of heparin (e.g. to detect tissue factor expression on circulating cells in sepsis) • In comparison to NATEM: Qualitative assay for the presence of heparin	
Sensitive to:		*Most affected ROTEM® parameters:*
Tissue factor expression on circulating cells in sepsis		CT
Factor deficiency (intrinsic and common pathways)		CT
Anticoagulants (e.g. direct thrombin inhibitors, not heparin!)		CT
Platelet contribution to clot firmness		MCF, A(x), (CFT, α)
Fibrin polymerisation (and fibrinogen concentration)		MCF, A(x), (CFT, α)
Hyperfibrinolysis		ML, LI(x)
(F-XIII-deficiency)		MCF, A(x), ML, LI(x)
Note:	The combination of NATEM and NA HEPTEM confirms the existence of heparin (LMWH or UFH) in the sample.	
Limitations:	• Insensitive to defects of primary haemostasis (platelet aggregation)	

The prolonged CT in the NATEM test is significantly shortened in the NA HEPTEM test (to normal value). CFT, A(x) and MCF parameters also improved. The sample therefore contains heparin.

Table 3-17: NA HEPTEM

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⚠ NOTE

Test selection mistakes

Test pipetting programmes RUO1 and RUO2 are meant to be used only for research.

- Use RUO1 and RUO2 for research purposes only.
- Configure your system without RUO1 and RUO2 if you do not need them for research use.

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4 System description

4.1 Scope of supply

The ROTEM® *delta* equipment is an integrated complete system containing the ROTEM® *delta* device and the respective system software. Table 4-1 lists components in the scope of supply.

The package ROTEM® *delta* (REF 200100-XX*) contains	**Amount**	**REF**
Accessories box	1	400472
Barcode scanner 3000/4000 (mounted)	1	400566
Barcode scanner 4000 holder (mounted)	1	208705
Barcode scanner 3000/4000 cable	1	400568
Clean Ring Rod	1	400596
Colour printer, USB (optional)	1	400222
Connecting cable eLine	1	400580
Connection cable for printer, USB (optional)	1	400223
Cup holder pro	4	400016
Disc for cup ejector pin	1	400586
eLine electronic pipette	1	700030
eLine pipette tips for the system pipette	1x96	400041
Keyboard	1	400475-EN
Keyboard protection foil	1	400476
MC Rod pro	1	400045
Set MC Rod Holder *delta*	1	200175
Measuring cells (Cup&Pin pro)	1x20	200011
Operating manual ROTEM® *delta*	1	200137-EN
Parameter list	1	200086-EN
Pipette cable guide	1	500220
Pipette and Clean Ring Rod holder	1	200206
Plastic clip for cup holder plus	8 (4 on holders)	400482-4
Reagent rack *delta*	1	400484
ROTEM® *delta*	1	200090
ROTEM® *delta* power cable	1	700048-XX*
ROTEM® *delta* Secure Viewer CD	1	200178
ROTEM® *delta* software	1	200140

Table 4-1: Scope of supply ROTEM® *delta*

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The package ROTEM® *delta* (REF 200100-XX*) contains	Amount	REF
ROTEM® system USB Stick	1	400221
Short user manual Single Use Reagents	1	200176-EN
Short user manual ROTEM® *delta*	1	200124-EN
Touch screen, glas	1	408707
Touch screen, glas, connecting cable USB	1	200185
Touch screen, glas, connecting cable VGA	1	400584
Touch screen, glas, connecting cable power	1	407107

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4.2 Hardware

4.2.1 Components

Figure 4-1 shows components of the ROTEM® *delta*, Table 4-2 shows their description.


Figure 4-1: Components of the ROTEM® *delta*

No.	Component	Description
1	Temperature controlled work area with ejector	The work area serves as a preheating station for cup holders and additionally contains a sample preheating station. The temperature in the work area is maintained at the preset measurement temperature (generally 37 °C ± 1.0 °C, physiological temperature). To facilitate removal of the blood-filled measurement cells (cups with pins), a metal pin as measuring cell ejector (with surrounding disc) is integrated into the work area.
2	Accessories box	Cups&pins within their package are stored in the drawer. The box of pipette tips can be placed on top.
3	Reagent rack	Holds the reagents during measurement.

Table 4-2: Components of the ROTEM® *delta*

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No.	Component	Description
4	Sockets	Sockets on the rear side of the device for the connection of touch screen cables, keyboard cable etc.
5	Measurement station with cup holders	Every ROTEM® system has four independent channels. The channels are tempered to 37 ± 1.0 °C (default, physiologic temperature). The cup holders are of massive, stainless material and contain two high power magnets for fixing.
6	Type label	The type label provides information on model, serial number, date of manufacturing, required patient distance, power data, IVD designation, CE mark, SGS mark and manufacturer. The symbols are described in Appendix A.
7	Touch screen with barcode reader	The software is controlled by the touch screen. The barcode reader scans Single Use Reagent, ROTROL and sample barcodes.
8	System pipette	Using the electronic system pipette, semi-automatic and software controlled pipetting and mixing of liquids is possible.
9	Keyboard with integrated track ball and protection foil	The keyboard with the integrated track ball is connected to one of the four USB slots at the rear side of the instrument. The keyboard is placed in front of the touch screen.
10	Axis	The pins are placed on the central, partly slotted axis.
11	Pipette and Clean Ring Rod holder with pipette cable guide	Holds the pipette, the pipette cable guide and the Clean Ring Rod.
12	On/off switch	The blue button on the right side of the device turns the ROTEM® *delta* on (the button glows blue, software is starting).
13	Pipette socket	Connection of the system pipette.
14	Guiding rods and rails	Guides the cup holder when placing it on the axis and ensures the proper positioning of the cup under the axis.
15	Main switch (optional)	This switch starts up the ROTEM® *delta* or disconnects it from the grid. Switch I – equipment is switched on Switch O – equipment is switched off
16	Power cable	Utilise only the original power cable. It is also available as an accessory (Appendix G).
17	two 6 pin plugs	1 for the connection with the touch screen monitor 1 for the connection with the ROTEM® *platelet*
18	MC Rod Holder	Holds the MC Rod pro

Table 4-2: Components of the ROTEM® *delta*

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4.2.2 Accessories

Accessories	Description
Operating manual	This manual contains comprehensive information on handling of the equipment.
Short user manual	This brief instruction is laminated to ensure high durability under laboratory conditions. Please store it near the instrument.
Pipette tips	Each ROTEM® system is supplied with a box with 96 single use disposable pipette tips (Tip Tray) suitable for the system pipette. Replacement tip trays are available as accessories (Appendix G).
Pipette filters	Each ROTEM® pipette is supplied with 5 replacement filters suitable for the system pipette. Additional replacement filters are available as accessories (Appendix G).
Measuring cells	The measuring cells (Cup&Pin pro) consist of one cup and one pin which are made of special plastic material for use with the ROTEM® *delta*. They are supplied as single-use disposable products (Appendix G).
Clean Ring Rod	Tool to remove the Clean Ring around the axis (Chapter 8.2.3).
MC Rod pro	With the MC Rod pro, the cup is pushed into the cup holder and fixed.
Integrated pin remover (plastic clip cup holder)	The cup holders are equipped with an integrated pin remover.
Reagent rack	ROTEM® *delta* is equipped with a solid reagent rack. There you can place the reagents during the tests. At the ROTEM® *delta*, the position for the reagent rack is given. There is no space for blood samples on the reagent rack. They should be stored on the sample preheating station ofthe system.

Table 4-3: Accessories. Please refer to Appendix G.

Printer (optional)

In case that a printer is connected to the system, it is possible to print numeric results and TEMograms. The installation of the printer is described in Chapter 6.4.4.

Barcode reader

The barcode scanner (REF 400 566) is mounted on the touch screen and connected on the rear of the device via one of the USB ports. It allows you to record the *patient* and/or the *sample ID* precisely, if the blood sample vials have an adhesive label with a barcode on it. Likewise, barcodes on Single User Reagents and ROTROL boxes / control sheets / vials can be scanned.

Note: always hold the vial, the ROTROL box / control sheet with the barcode to be scanned at the level of the arrow on the touch screen. This is the best position for scanning.

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4.3 Software

The user management is preset to show the log in screen when switching on the device. After the log in process the measurement module is shown. After quitting the measurement module, the main screen is shown.

The ROTEM® *delta* software supports you through intuitive user guidance during the performance, the interpretation and the management of measurements. The central screen of the software is the main screen. From there you can open the eight modules using the touch screen (Figure 4-2). Table 4-4 shows an overview of the modules. Chapter 6 gives a detailed description of the modules.

Figure 4-2: Main screen of the system software

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Module	Description	Details
Measurement module	In the MEASUREMENT MODULE, you can select, start, and interpret tests.	Chapter 5
Database	In the module DATABASE, you can manage the sample measurements and QC results and further evaluate them.	Chapter 6.3
USB stick	Data in the ROTEM® system is usually transferred via a USB stick. The USB stick is used for backup/restore and the export of TEMograms and numerical data.	Chapter 6.8
File system	In the module FILE SYSTEM and in the sub modules ADMIN TOOLS and USER, you can manage some sytem settings, the authorized users and the exported data, subject to the condition that you have appropriate user rights.	Chapter 6.7
Liquitrans	The LIQUITRANS module guides you step by step through the pipetting of up to 2000 µl liquid with the automatic pipette.	Chapter 6.6
Service	In the module SERVICE, you influence the recording of raw data, subject to the condition that you have appropriate user rights. The module also gives an overview of the technical interventions on the system.	Chapter 6.5
Setup	In the module SETUP, you change system settings, subject to the condition that you have appropriate user rights.	Chapter 6.4

Table 4-4: Modules of the system software

4.3.1 Technical data

Measurements and weight	Length in mm	545
	Width in mm	335 without pipette holder
	Height in mm	205 without touch screen/pipette
		approx. 590 including touch screen
	Weight in kg	approx. 15 without touch screen
		approx. 22 including touch screen
Power supply instrument	Input	115/230 V_{AC}; 2.3/1.15 A; 60/50 Hz
	Power input	Max. 265 VA
Power supply touch screen (powered by mains supply)	Input	115/230 V_{AC}; 0.5/1 A; 60/50 Hz
	Power input	Max. 50 VA
Power supply touch screen (powered by ROTEM®)	Input	24 V_{DC}; 0.7 A
	Power input	Max. 16 VA

Table 4-5: Technical data

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Installation	Installation category	I
Data input	USB	See respective operating manual
Data output	USB, LAN	See respective operating manual
	optional printer	See operating manual for the printer
Ambient conditions	Operating temperature	+15 °C to +30 °C
	Storage temperature	+0 °C to +50 °C
	Relative humidity	20% to 85%
	Operable to height above sea level	2000 m
Temperature control	Measuring station	+30 °C ± 1.0 °C to + 40 °C ± 1.0 °C Reference measurement temperature +37 °C ± 1.0 °C
Specifications	Measuring principle	Thromboelastometry
	Angular deflection detection	Optical system with CCD sensor
	Heater	Independent heaters with independent temperature sensors: - for each measuring position (altogether 4) - for the temperature controlled work area
	Number of channels	Four independent channels
Connections	Power supply	Original ROTEM® power cable (country specific)
	Touch screen	Three original ROTEM® connection cables for the equipment
	Pipette connection	Original ROTEM® system pipette with helix cable and RJ12/Western connector (6-pin)
Data backup	USB stick 1GB	Allows to export data to other media.
Printer (optional)	supported printer	Printers with USB port are supported. Please ensure that required drivers are available in the system (Chapter 6.7.1). The printer must follow IEC 60950 standards to comply with DIN EN 61010.

Table 4-5: Technical data

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⚠ WARNING

Danger of infection

Inappropriate handling of contaminated samples may cause infections.

- Always wear personal protective equipment such as laboratory coat, gloves and goggles when handling samples.
- Always follow guidelines for infective materials[1].

[1] Laboratory Biosafety Manual World Health Organization, 1984

⚠ CAUTION

Measurement errors

Inappropriate storage of the materials may lead to measurement errors.

- Always store measurement cells, pipette tips, reagents and samples according to the environmental conditions as recommended on the package.

⚠ NOTE

Damage of the touch screen

Inappropriate handling of the touch screen may lead to damage.

- Do not press the keys of the keypad with sharp or pointed objects (e.g. pencils).

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5.1 Introduction

The ROTEM® *delta* system software supports you during the measurements through the step-by-step menu. You will find additional instructions on how to interpret measurement results in the HELP menu, which can be opened by the **Help** button on the screen.

This chapter defines how the different ROTEM® tests are conducted. However, it is essential to read the other chapters of this manual thoroughly.

5.1.1 Usage of the pipette

Figure 5-1 shows the system pipette, Table 5-1 lists its buttons.

Note: The buttons are only active in the pipette mode of the measurement module.

Figure 5-1: System pipette

Button	Function
SELECT	none
ENTER	none
➢ \| ➢	none
⏻	tip eject
CONFIRM	this central blue button triggers the next pipetting step
RESET	small reset button on the back side of the pipette

Table 5-1: Buttons of the system pipette

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Note the following during pipette handling:

- Hold the pipette vertically during aspiration of reagents.
- Immerge the tip only a few millimetres into the liquids.
- To ensure that no air is absorbed during pipetting of blood, immerge tip approximately 10 millimetres.
- Only press briefly the confirm button.
- Use the integrated pipette holder to put down the filled pipette in order to avoid contamination of the mechanical part of the pipette.

5.1.2 Sample handling

Blood sample collection

Take the blood sample carefully according to the given references:

NCCLS Document H3-A4 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture Fourth Edition; Approved Standard (1998)

NCCLS Document H21-A2. Collection, transport, and processing of blood specimens for coagulation testing and performance of coagulation assays, 3rd ed. Approved Guideline 1998

Transport of samples

⚠ CAUTION

Measurement errors

Inappropriate transport of samples may lead to measurement errors.

- Do not shake samples.
- Do not roll the sample container.
- Prior to use, pneumatic tube transportation systems need to be validated.
- Do not expose the blood samples to cold temperatures.
- Never store samples on ice!

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5.2 Measurement preparation

Please note the following during measurement preparation:

- Prepare tests according to test specific instructions.
- Gain and process samples according to the instructions.

5.2.1 Switching on the ROTEM® system

✓ There are no cup holders present on any of the four measurement channels.

1. If required, activate the ROTEM® system with the optional main switch on the back of the device (switch setting „I“).
2. Switch on the printer (optional).
3. Push the blue on/off button on the right hand side of the instrument.
 ↳The ROTEM® system software starts automatically and checks the system.
 ↳The control unit is initialised.
 ↳The connection between the ROTEM® system, the connected pipette and the optional printer is established.
 ↳The log in screen or the measurement module screen is loaded (depending on system configuration).
4. Log in to the system.
 - Touch the screen if the screen saver is active.
 - Select user.
 - Enter password.

 ↳Measurement module screen is shown (Figure 5-2).
 ↳Heating up of the device to operating temperature. This procedure may take several minutes. The colour of the 5 temperature indicators down to the right of the screen (1) changes from dark blue (very cold) over light blue (less cold) to white (target temperature).
5. The system displays QC and maintenance reminders according to the presetQC (Chapter 6.4.8) and maintenance schemes (Chapter 8.2).
6. Perform QC measurements (Chapter 5.5.4) and maintenance as prompted by the reminder.

Figure 5-2: Details of the measurement screen

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The colour of the status line (2) informs about the status of the measurement for each channel.

Colour	Status
Grey	Channel ready to use
Blue	Channel not ready, initialisation may be active
Green	Measurement is running
Yellow	Warning message during measurement
Red	Error, measurement is stopped

Table 5-2: Colours of the status line

5.2.2 Reagent preparation

hep-tem® and **ex-tem® L** (liquid)

Before performing a HEPTEM or EXTEM L test, dissolve the lyophilised reagent (one bottle contains reagent for several tests).

1. Withdraw one small Lyo bottle (hep-tem® or ex-tem® L) and a big diluent bottle (hep-tem® or ex-tem® L Dil) from the package.
2. Open both bottles.
3. **Quit** measurement module and confirm with **Yes** to exit.
4. Select LIQUITRANS module.
5. For hep-tem® enter 200 *µ litre* in the entry field; for ex-tem® L enter 230 *µ litre* in the entry field.
6. Confirm with **Start**.
7. Take a new pipette tip.
8. Immerge the pipette into the diluent (Dil).
9. Push the blue start button of the pipette.
 ↳200 or 230 µl liquid is aspirated.
10. Hold the pipette tip over the target vial (Lyo).
11. Push the blue confirm button of the pipette.
 ↳200 or 230 µl liquid is dispensed again.
12. Choose **Done**.
13. Exit LIQUITRANS module with **Quit**.
14. Close both bottles with their caps.
15. Place the diluent in its package into the refrigerator.
16. Gently swirl the dissolved Lyo bottle.
17. Write the date of dissolving onto the given field on the new bottle.
 ↳The reagent is ready for use.
 ↳The hep-tem® reagent may be used for a period of 30 days; the ex-tem® L reagent may be used for a period of 7 days.
18. Place it in the reagent rack.
19. Place the reagent rack to its position in the front of the ROTEM® *delta*.
20. Before first use of the hep-tem® or ex-tem® L reagent, let it rest for 10 minutes.

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Note: It is also possible to reconstitute these reagents directly from the pipetting sequence in the MEASUREMENT menu.

All other reagents

Prepare all other reagents (except hep-tem® or ex-tem® L) as follows:

1. Take reagent rack with reagents already in use out of the refrigerator.
2. Withdraw additional fresh reagents for the work day from the packages inthe refrigerator.
3. Write the date into the given field on any new liquid reagent bottle.
4. Refill the reagent rack.
5. Place the reagent rack to its position in the front of the ROTEM® *delta*.
6. Wait until the reagent reaches room temperature.

Note: During the workday, reagents may be kept at room temperature but vials must be closed tightly immediately after each use.

5.2.3 Refill cups and tips

Refill tips and measurement cells in the accessories box if necessary.

5.2.4 Measuring cell preparation

⚠ CAUTION

Measurement errors

Inappropriate handling of the measurement cells may lead to measurement errors and error messages.

- Guide the axis into the upper opening of the pin.
- Place the pin firmly all the way up the axis.
- Push the cup using the MC Rod down until the cup is flush with the surface of the cup holder. Otherwise the cup is positioned too high and the pin will touch the bottom of the cup during measurement.
- Always touch measuring cells at the exterior of the cup, never touch the pin.

The steps for the preparation of the measurement cell are summarised in the submenu **Preparation** of the measurement module screen. Please follow these instructions (Figure 5-3).

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
Figure 5-3: Instructions for the preparation of measuring cells

1. Take the cup with the pin in it from the storage box.
2. Push the pin in the cup onto the axis chosen for the measurement.
 Note: The status line under the channel is grey when the channel is ready. In case the axis has been moved heavily when attaching the pin, the channel becomes inactive and the status line turns blue during the time of initialisation.
3. Place the cup with its opening facing upwards into the appropriate preheated cup holder.
4. Place the cup holder onto the temperature controlled work area again.
5. Push and fix the cup in the cup holder using the MC Rod.

The measuring cell is prepared.

5.2.5 Sample preparation

⚠ CAUTION

Measurement errors

The temperature of the sample may influence the measurement results and/or reproducibility.

- Measure the blood sample directly after sampling.
- If this is not possible, preheat the blood sample for 5–10 min before measurement in the sample preheating station of the ROTEM® *delta*.
- Do not leave the sample more than 30 min in the sample preheating station.

5.2.6 Adjust measurement temperature at the ROTEM®

Normally the samples are measured at a temperature of 37 °C. Under certain circumstances, it may be necessary to adjust the measurement temperature in the module SETUP, e.g. in case of an undercooled patient.

Note: This requires the user profile "Administrator" (admin).

1. Leave measurement module with the **Quit** button.
2. Confirm with **Yes**.

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3. Select SETUP module.
4. Enter required temperature in the USER SETUP module (Chapter 6.4.1).
5. Exit module with **Quit**.
6. Confirm with **Yes**.
7. Select MEASUREMENT MODULE.
 ↳The status line now displays the new temperature (highlighted in red).
 ↳The 5 temperature indicators at the bottom right of the screen shortly turn to red or blue until the new temperature is reached.

Note: After measurement, set temperature to 37° C again.

5.2.7 Select Channel

The measurement screens for the four channels appear identical. The lower part of the screen constantly displays all channels. With channel 1 to 4 left to right. The active selected channel is highlighted in blue. The upper half of the screen displays the patient demographic data fields as soon as a channel is selected. While a test is pipetted, the pipetting instructions are displayed. During measurement, a larger TEMogram and its numerical results are shown.

1. Each of the four channels can be selected to start a measurement by touching the screen on the corresponding channel.
2. To obtain a view of the 4 channels press **Multi-TEM**.

5.2.8 Select test

For each channel a test is preset in the menu SETUP. These presets may be overwritten and another test may be chosen.

Note: This choice cannot be changed during a running measurement.

1. Touch the test name on top of the selected channel.
 ↳A list with all possible tests pops up.
2. Select required test by touching it.
 ↳ The new test name above the channel is highlighted in the color code of the ROTEM® reagent bottle.
 ↳The colour coded reagent(s) necessary for this test are displayed under the *comment* field of the patient data.

If necessary, you may adjust the order of the tests in the SETUP menu. Frequently applied tests are listed at the top of the list (Chapter 6.4.5).

5.2.9 Enter patient data

Patient data is assigned to each channel and therefore each measurement. The patient data contains patient ID, patient name (first name, last name), sample ID, comment, date of birth and gender. In the module SETUP it is defined which fields are available and if these are mandatory or just available (Chapter 6.4.1). It is also possible to select a normal range profile specific for the patient (e.g. pediatric).

1. Touch one of the four channels.
 ↳In the upper part of the screen entry fields for patient data are displayed (Figure 5-4).

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
Figure 5-4: Entry fields for patient data

2. Touch the respective entry field.
3. Enter patient data.

 Note: Alternatively patient data may be read in, using a barcode reader.

 Note: Mandatory fields are displayed in blue.
4. Choose a normal range profile for the patient.

 ↳ The normal values from this profile are used for this patient.

 ↳The profile overlay TEMogram is formed with these normal values.

Conventions for entry

Consider the following conventions when entering patient data to sort the measurements in the database logically and to avoid failures.

- Enter names with first and last name. Datasets can be sorted consistently by names.
- Do not use special characters like /, \, #, * ,, „, ?, < . Errors may occur when exporting data containing these characters.
- Do not use spaces within the patient ID. Errors may occur, especially during patient overlays.

Simplification when entering data

The system supports you during the entry of patient data.

- If a patient ID contained already in the database is re-entered or read in, the fields containing *Patient name*, *Date of birth*, and *Gender* are filled in automatically with the according deposited patient data when leaving the field with the <**Tab**> key. The entries in the fields, containing automatically filled out data, cannot be changed.

 Note: Patient data is transmitted to the database when the according measurement is saved/cleared (Chapter 5.4.3).
- Patient data can be copied from one channel to another by **Copy** and **Paste** or **Copy to all**.
- All patient data of a channel can be deleted with the <**Esc**> key or with the **Reset data** button.
- When **Single Patient Mode** is selected in SETUP, the patient specific entry of the channel is automatically copied in all other channels.

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⚠ CAUTION

Measurement errors

Inappropriate handling of the pipette may lead to measurement errors.

➥ Follow instructions for the handling and cleaning of the pipette.

⚠ CAUTION

Measurement errors

User errors like spilling the reagent or the usage of a wrong reagent may lead to measurement errors.

➥ Cancel the measurement in this case.

➥ Read the instructions for use of the reagents.

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5.3.1 Start measurement

1. For liquid reagents: mix reagents gently and thoroughly and open vials.
2. Mix the blood sample by repeated slow tilting of the sample.
3. Start the pipetting sequence for the selected channel by using the **Start** button.
4. If required, scan the (single use) reagent vial barcode upon prompt.

 Note: always hold the vial near the arrow on the screen for scanning.

 Note: alternatively, enter manually the reagent's expiry date, lot number and TEM Code (the 4 numbers just above the barcode on the vial) with **Manual input** and press **Check.**

 ↳The instructions (incl. pictogram) for the next step to be performed are displayed on the screen (Figure 5-5).
5. Open the single use reagent vial after scanning.

Figure 5-5: Pipetting sequence

6. Follow the menu navigation of the pipetting programme and confirm each step by pressing the blue confirm button of the pipette.

 Note: If reconstitution of the hep-tem® or ex-tem® L reagent is required, press **Reconstitute.**

 Note: If the pipetting sequence requires a pipette tip change, the blue confirm button stays inactive until the pipette tip is ejected by using the eject button.
7. The pipetting steps can be repeated in the event of pipetting failure.

 Note: Always hold the pipette tip above the waste container before you select **Step back** (**<**). The pipette will automatically empty the liquids in the pipette tip when performing this step back.

 - Press **Step back** (**<**) button (arrow to the left).

 ↳The pipetting program blows out the liquids in the pipette tip.

 ↳The previous pipetting step is shown.
 - Follow the pipetting sequence again.

 ↳After the blood has been added to the reagents, the reaction starts.
8. Place the cup holder onto the measuring position using the guiding rods.

 Note: The cup holder is held in measuring position with magnets. In case that the cup holder is not positioned in time, the measurement is aborted.
9. If required, press <enter> or touch the screen to confirm that the cup holder is in place.

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5.3.2 Graphics of measurement results

During measurement, the large TEMogram is shown at the left upper side of the screen. In the upper right part of the screen the current measurement results of the test parameters are shown (Chapter 3). Some of the coagulation parameters are shown as “preliminary values” in the display while a clot is developing. These values are based on a limited amount of data and are replaced by the final values as soon as the appropriate amount of data points are available. Preliminary results are marked with an asterisk (*).


Figure 5-6: Measurement module with profile overlay

With **Profile overlay,** a TEMogram within the normal range of the current patient profile is displayed over the measured TEMogram.

By using the **Patient overlay** button, the 4 newest TEMograms of the same test and the same patient from the database are displayed. With each overlayed TEMogram the starting time of the according measurement is displayed in colours. By using the **Patient data** button, you can switch between the TEM view and the patient data (Figure 5-6).

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By using the **Multi-TEM** button an overview of all four channels is displayed (Figure 5-7).

For measurements run > 60 minutes, it is possible to view the graph from 60-180 minutes by double-touching the graph.

Figure 5-7: 4 TEMograms in the measurement module

5.3.3 Warning messages during measurement

Usually the **Warning** button is hidden. It only appears if a system error occurs. The status line directly under the according channel turns yellow, the **Warning** button flashes red.

1. Touch **Warning** button.
 ↳The channel with the error is displayed in red.
2. Note the caption down for eventual requests of the service provider (e. g. CCD invalid).
3. Read the descriptions and instructions under the caption.

 Note: The messages are shown with a short coded error message and a longer description of the cause and troubleshooting of the problem. The messages in the warning windows are self explanatory. See Table 5-5 for additional help for troubleshooting.
4. Check if TEMogram and pameter results are plausible and decide if the measurement shall be continued (**Ignore**) or cancelled (**Stop**).

 Note: In case of any doubt, stop and repeat measurement.
5. Please contact your service provider if the error occurs again during the next measurement.

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5.3.4 Interpretation of measurement results

The system software helps you to interpret the measurement results.

1. Select **Help** button.
2. Select **Interpretation Cases** and/or **Expert Recommendations**.

5.4 Finish measurement

A measurement is finished in several steps described in the sections from Chapter 5.4.1 to Chapter 5.4.8.

5.4.1 Stop measurement

A measurement can be stopped automatically or manually.

- The measurement is stopped automatically when the maximum measurement time determined in the module SETUP is reached.
- **Stop channel** stops the measurement manually.

5.4.2 Print measurement results

✓ Channel is not yet cleared/discarded.

1. Select **Print** button.
2. Select number of TEMograms (measurements) per page.
 - 1 TEM (1 TEMogram over 60 min and additionally the same TEMogram over 8 hours)
 - 4 TEM (4 different TEMograms of 60 min each)

 ↳If normal values (or QC target ranges) have been entered, they are printed as well.

 ↳Overlays are printed when selected on the screen and when pre-set in SETUP - PRINTER SETUP.

 ↳The normal range profile for the patient is printed.

Note: If a warning message was displayed during a run, the respective printout will be marked with “check plausibility”. Check the data for irregularities and completeness and consider to discard the measurement results.

Note: Patient or control results that are out of range are flagged with an arrow up or down. Patient measurements run after missing yearly maintenance are flagged. Patient measurements run after missing or out of range QC are equally flagged.

5.4.3 Save and clear measurement

When a channel is saved and cleared, the data is transmitted to the database and only remains available there for printing (Chapter 6.3.5) or exporting. Check especially the patient ID and the sample ID before transmitting the measurement to the database. For safety reasons, data cannot be changed after a measurement has been finished and transmitted to the database.

1. Select **Save/clear channel**.

 ↳The measurement is removed from the screen.

 ↳The measurement is saved in the database.

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5.4.4 Discard measurement

An erroneous measurement may be discarded. For safety reasons, it will still be present in the database after discarding. It will however only be highlighted in light grey. In case of a connection with the LIS or HIS, it will not be transferred. Moreover, it is not used for patient overlays.

1. Select **Discard channel**.
2. Enter comment.
3. Press **Accept.**
 ↳The measurement is removed from the screen.
 ↳The measurement is saved in the database but greyed out on the database display.

5.4.5 Remove measuring cells

⚠ CAUTION

Wrong calculations of the parameters

Removing the cup holder before measurement is finished may lead to wrong calculations, especially of the lysis parameters.

➥ Remove cup holder only after the measurement has finished.

The cup holder is equipped with an integrated blue pin remover. The cup and the pin are removed together with the cup holder from the measuring position.

1. Hold the cup holder with one hand and push, with the thumb of the other hand, the blue pin remover up towards the device (Figure 5-8).
 ↳The pin is pushed into the cup.
2. Remove cup holder.
 ↳Cup holder, pin and cup are removed together.


Figure 5-8: Remove cup holder

Remove cup and pin

1. Insert metal ejector pin at the right edge of the work area into the bottom of thecup holder.
2. Push the cup holder down the ejector pin.

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↳Cup is being removed (Figure 5-9).

3. Dispose of cup and pin according to the effective regulations for biohazardous material.


Figure 5-9: Remove cup from the cup holder

5.4.6 Log off the system

1. Leave measurement module with the **Quit** button.
2. Confirm with **Yes**.
3. Select **Logout**.
 ↳The log in screen is shown.

5.4.7 Daily maintenance

A detailed description of the necessary maintenance actions can be found in Chapter 8. If yearly maintenance is not performed, all subsequent measurement results will be flagged on printouts and in the database (Chapter 8.2).

The system automatically reminds the user to perform maintenance. The user is responsible for the correct following of the maintenance plan.

Following maintenance must be performed daily on the ROTEM® *delta*.

1. Clean and disinfect cup holder.
2. Clean and disinfect the lower surface of the measuring block (around the axis).
3. Check the pipette filter for blood stains and replace if required (Figure 8-10).

5.4.8 Switch off

Normally the device remains ready for use.

Note: Finish running software before switching the device off, to avoid damage of the database.

Note: Switching off the system other than described below (e.g. pressing the blue on/off button long) can lead to damages to the hard disk or wrong time settings.

1. Select **Quit** button in the measurement module.
 ↳The main menu is shown.

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2. Press blue on/off button on the right side of the device.
 ↳The software is being shut down.

In addition, if the device has to be transported, the main power switch next to the power cable has to be switched off. The monitor also needs to be switched off and should be secured before transportation.

5.5 Quality control

The proper function of the measuring mechanism and electronics are controlled via the equipment's auto-monitoring. Channel status information and warnings are displayed automatically as soon as a self-test failure is detected.

Errors	Display in the status line	Consequence
Channels are not ready for use before measurement	blue	**Start** button is not active. Measurement cannot be started.
Channel fails during measurement	yellow	Warning is displayed. Measurement is marked with "check plausibility" on the printout and with an error code in the database. Continue or stop measurement.

Table 5-3: Checking the measurement channels

Like any other *in vitro* diagnostic system, the ROTEM® system requires quality control by performing tests with standardised quality control (QC) materials. The ROTEM® standardised system controls ROTROL N and ROTROL P are based upon human plasma and will show reaction curves at two different levels.

The ROTEM® quality controls ROTROL N and ROTROL P are run as performance tests to identify potential problems with

- reagents
- consumables (tips, cups)
- user technique
- pipetting volumes

5.5.1 Control intervals

The frequency of the quality controls with ROTROL N and ROTROL P has to be chosen according to local guidelines. Thanks to the continuous control of the function of the measuring mechanics and the measuring electronics directly on the equipment, a weekly quality control is recommended beeing sufficient. The system reminds the user automatically when to perform QC according the preset QC-frequency (Chapter 6.4.8). The user is reponsible for the correct following of the QC plan.

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5.5.2 Continuous self-monitoring when powered on

The values included in the table below are representative of the continuous monitoring system that monitors the performance of the axis movement, sensors and temperatures. The user accesses the SERVICE menu (SERVICE > SETTINGS) and evaluates the following values per channel (see also Chapter 8.2.3):

Value	Definition	Target range
Amplitude	Measurement confirming the movement of the axis and proper functioning of the supporting ball bearings	4400 to 5100
Center	Measurement indicating the location of the amplitude on the CCD chip	-832 to 708
Variance	Measurement providing information about the mechanical status of each channel	<12,000
Temperature	Measurement of the instrument temperature settings as well as the pre-warming station (temperature controlled work area)	35°C to 38.5°C

Table 5-4: Acceptable range for the internal quality control values

The internal monitoring validates the proper functioning of the axis on each channel and the performance of the instrument optics depicted in the SERVICE menu channel graphs:

1. Black curve depicts the actual light beam from the CCD.
2. Red curve depicts the calculated light beam from the CCD.
 Note: The values are documented in the menu SERVICE-LOG BOOK-DATALOGGER

5.5.3 Preparation of quality control material

⚠ WARNING

Danger of infection

The raw material for each batch of ROTROL was tested using appropriate methods and was found non-reactive with HBsAg and antibody-negative for HIV 1, HIV 2, HCV and *Treponema pallidum*.

There is no known test method that is able to provide information with 100% certainty that products obtained from human blood are free from hepatitis, AIDS or other infectious diseases. These products should be regarded as potentially infectious, just as all materials of human origin.

- Always wear personal protective equipment such as laboratory coat, gloves and goggles when handling samples.
- Always follow guidelines for infective materials.

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The control materials are easy to prepare and no accurate pipetting is required.

✓ The control material has reached room temperature.

1. Reconstitute ROTROL N (or P) by transferring the content of a vial with diluent into a vial of lyophilisate (Figure 5-10).

 Note: The vial should not be emptied completely with a pipette.
2. Tightly close the vial with rubber stopper and screw cap.
3. Gently swirl the lyophilisate to dissolve it.

 Note: Ensure the lyophilisate is dissolved completely.
4. Let plasma reconstitute in the closed vial for 15 minutes (ROTROL N) or 30 minutes (ROTROL P).
5. Carefully mix it.
6. Check measurement temperature in the status line (37 °C).
7. Adjust temperature under **Setup - User Setup,** if it is not set at the appropriate temperature.
8. Prewarm for 5 minutes on the temperature controlled work area.
9. Mix again carefully by gently inverting and swirling.

The control is ready for use.

Figure 5-10: Reconstitution of ROTROL N (or P)

5.5.4 Performance of the quality control

⚠ CAUTION

Error during quality control

Systematic errors may possibly not be recognised by inexperienced users.

➥ Quality control has to be performed regularly by every ROTEM® user.

The control is usually performed as a QC test combination on all 4 channels e.g. 2 channels INTEM (S) and 2 EXTEM (S). INTEM (S) and EXTEM (S) are the most commonly performed tests and contain the base activators for FIBTEM (S), APTEM (S) and HEPTEM (S).

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Note: The required QC test combination appears automatically on the channels and can be configured in SETUP - QC SETUP - REMINDER (Chapter 6.4.8).

The only difference between the pipetting sequences for the routine tests and the ones for the control tests, is that instead of blood samples, controls are pipetted. The control measurements are then saved in the special QC database (Chapter 5.6.1).

Control material	Range of control	QC Test	Equivalent patient test
ROTROL N	level 1	QCinN, QCinN S	INTEM, INTEM S
ROTROL N	level 1	QCexN, QCexN S	EXTEM, EXTEM S
ROTROL P	level 2	QCinP, QCinP S	INTEM, INTEM S
ROTROL P	level 2	QCexP, QCexP S	EXTEM, EXTEM S

When the QC reminder is displayed on the screen:

1. Check the proposed combination of QC tests.
 - Confirm the proposed tests with **OK**
 - Or select another QC test combination with **Other QC combination**.

When the QC reminder is not displayed:

1. Choose **ROTROL (P/N)** in the test bar of any channel.
 ↳The QC mode is opened.
 - Confirm the proposed QC test combination with **OK**
 - Or select another QC test combination with **Other QC combination**.
2. Select a channel.
3. Touch the field *Lot number*.
4. Scan the ROTROL vial (or select the ROTROL lot from the drop down list).
5. If prompted (by new lot), scan the ROTROL box or the control sheet.
 ↳The target values for this new lot are automatically saved.
6. Press **Start.**
7. Follow pipetting sequence.
8. Let the QC measurement run until parameter A20 has been reached.

Note: with every new ROTROL lot, the corresponding target ranges displayed in SETUP-QC SETUP-ROTROL have to be compared to the ones on the control sheet.

5.5.5 Check quality control results

The ROTEM® QC materials have to produce repeatable results in the range of the ROTROL control sheet during the whole quality control.

The target values for a specific ROTROL lot are automatically saved in the system after scanning the corresponding ROTROL box or control sheet (either in the MEASUREMENT menu or in SETUP - QC SETUP). Alternatively, the target values of the according tests for each new ROTROL lot can be entered manually in SETUP - QC SETUP (e. g. for QCinN).

↳QC results deviating from the target values are highlighted in red on screen.

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↳All patient results obtained after any out of range (or missing) QC result are flagged in the database and on printouts. In the database they are also highlighted yellow.

Note: Flags will appear the day after the QC reminder due date.

If no target values are available in the system, the QC results have to be compared with the control sheet in the package insert.

Complying with QC reminder rules

- The required QC tests should be run as soon as the QC reminder is displayed on the screen (e.g. every Monday morning at 8 am).
- Otherwise QC tests should be run at the latest before midnight that same day (e.g. Monday).

 ↳In case the controls aren't run before midnight, all following patient measurement results will be flagged with "QC missing" in the database and on the printout.
- QC is passed when within range QC results are obtained for all tests from one of the combinations activated in SETUP - QC SETUP - REMINDER (Chapter 6.4.8)

 ↳If any QC result of a QC test combination is out of range, all following patient measurement results will be flagged with "QC out of range" in the database and on the printout.
- Flags (QC missing or QC out of range) will not be posted anymore as soon as a successful QC combination is run.

Note: The current QC status (OK, out of range, missing, blocked) can be viewed in SERVICE - INFORMATION.

5.5.6 Suggested actions for QC out of range

If a QC result from the required QC test combination is outside the expected range, examine the potential reasons.

1. **Discard** the out of range QC measurement in the MEASUREMENT module.

 ↳Only the missing QC test needs to be re-run.

 Note: The missing QC test should be run within 4 hours. After that, a complete QC combination needs to be run in order to pass QC.
2. Rerun the missing QC test (on any available channel).
3. In case the new QC test is within the target range, patient samples may be tested again and valid results can be expected.[1]

 ↳patient measurements will not be flagged anymore with "QC out of range"
4. In case the new QC test is again outside the target range, contact your service provider for further troubleshooting or follow the instructions from Figure 5-11.

[1] **Literature:** NCCLS. Statistical quality control for quantitative measurements: principles and definitions; Approved Guideline – Second Edition. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA, 1999.

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
Figure 5-11: Possible causes for out of range QC measurements

If the out of range result is confirmed in the original channel but the result of the second channel is error free, probably a channel specific problem or a calibration problem exists. Do not use this channel for further measurements and contact your technical service provider.

5.6 Managing measurement results in the database

The following sections describe the most important steps to find, print and export measurement and QC results in the database. Chapter 6.3 contains detailed information on these topics.

Patient results should be interpreted as soon as they become available in the MEASUREMENT menu.

The database should only be used as a documentation tool, not for result interpretation.

5.6.1 Find a patient / QC result

1. **Quit** measurement module.
2. Confirm with **Yes**.
3. Open DATABASE module.
4. Select either the patient database or the QC database with the **Patients/QC** button.

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5. Select requested patient or QC result either by entering the search criteria into the search form or by touching its line in the list on the upper left side (Figure 5-12).


Figure 5-12: Select measurement

6. To search the patient/QC results using the search form, select **Start search**.
 ↳The located measurements are highlighted in the list.
 ↳The according TEMograms are overlayed in the TEMogram window.
7. Add located patient/QC results with the **Add** button to the selection box on the lower left.
8. To search further patient/QC results, clear search criteria by **clear search** and search and add further measurements.
9. To remove single patient/QC results, highlight the measurements in the selection window and select **Remove**.
10. To clear all patient/QC results from the selection window, select **Reset**.

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5.6.2 Print patient / QC results

Normally measurement results are printed immediately after the measurement. Additionally there is the possibility to print data directly from the database.

1. Select and add patient / QC results.
2. Select **Print**.
3. Determine number of TEMograms per page.
 - Printout 1 (1 measurement result: the TEMogram of 60 min and the same one of 8 hours).
 - Printout 4 (4 measurement results: 4 different TEMograms of 60 min each).

 ↳Flags and errors are printed together with TEMogram and numerical results.

5.6.3 Export patient or QC results to USB stick

Special characters (*, /, &, etc.) are not accepted in file names when copied from the internal hard disk to an USB stick. Please do not use any such characters when renaming a file during the export function.

1. Select and add patient / QC results.
2. Select **Accept** button.

 ↳Evaluation screen is opened (Figure 5-13)


Figure 5-13: Evaluation screen

3. Touch and select a maximum of 4 TEMograms for the export.

 ↳The 4 TEMograms are displayed overlayed on the right side.
4. Select **Export 4 TEM**.

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↳The TEMograms are saved in a picture on the hard drive in the file system.

5. Enter a title if required.
6. Confirm with **OK**.
7. Optionally, change file name of the exported file.
8. Confirm with **OK**.
9. Repeat action until all required TEMograms are exported.
10. **Quit** DATABASE module.
11. Confirm with **Yes**.
12. Open FILE SYSTEM.
13. Select USER.
14. Select GRAPHEXPORT.

 ↳A preview of the stored and exported TEMograms is displayed (max 4 per screen).
15. Insert a USB stick to the USB port at the right side of the device.
16. Open USB STICK module.
17. Drag&drop the exported pictures with your finger or the trackball to the left side into the USB stick window.
18. **Copy** or **Move** the pictures.
19. Close the file system (right) and USB stick (left) with **Quit**.

 ↳The pictures are now saved onto the USB stick.
20. Remove USB stick.

The TEMograms are now available for further processing on a Windows computer.

5.6.4 Backup

In case the remote backup feature has not been implemented (Chapter 6.7.1), take a backup of the database and the settings at least once a week.

1. **Quit** measurement module.
2. Confirm with **Yes**.
3. Open FILE SYSTEM.
4. Select ADMIN_TOOLS.
5. Select DB BACKUP.
6. Insert USB stick to the USB port at the right side of the device.
7. Enter START BACKUP.

 ↳Message appears: Backup in progress, please wait.

 ↳Message appears: Backup has been successfully copied to USB stick!
8. Confirm with **OK**.
9. Remove USB stick.

Chapter 6.7.1 describes how to restore the backup data in case of data loss.

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5.7 Troubleshootin

In case of regular errors, contact your serviceprovider.

5.7.1 Warning screen

Probable cause of the error	Warning message and code
Misaligned optical detection system	CCD_DARK 1003
Stray light inside the device	CCD_BRIGHT 1013
Failure in optical detection (LED, mirror, CCD)	CCD_NOISE 1023
Temporary malfunction of axis	CCD_INVALID 1033
Dried sample still in place in holder	POS_TIMEOUT 2013
Obstacle in the optical path	CARRIER_SPIKE 3003
Interrupted data stream (computer/board)	CARRIER_LOST 3013
Severe device failure (hardware clock, communication firmware)	CARRIER_INVALID 3023
Significant mechanical shock or vibration	CARRIER_TIMEOUT 3033
Drive failure (motor, voltage supply)	MOTION_VANISH 4003
Cuvette or pin not correctly placed	MOTION_DRIFT 4013
Worn ball bearing on the axis	MOTION_BAD 4023
Axis jammed	MOTION_TIMEOUT 4033
Mechanical problem (excenter loose, slider grinding)	AMPL_JUMP_HIGH 5003
Firmness cannot be determined due to mechanical failure	AMPL_JUMP_LOW 5013
Sample was removed during measurement	AMPL_BAD_IDLE_LOW 5023
Clot is torn-off from the cuvette surface	AMPL_BAD_IDLE_HIGH 5033
Axis heavily disturbed during placing of pin	AMPL_BAD_SAMPLE 5043
Failure when smoothing the TEMogram curve	TGM_SPLINE 6003
Clotting already started before positioning cuvette holder	TGM_BASELINE 6013
Running measurement influenced by drying of sample	TGM_NOSE 6023
Signal lost during pipetting sequence	TGM_DRYING 6033
Yearly maintenance missing	TGM_DRY_PARA 6043
Quarterly maintenance missing	SYNC_LOST 10093
Weekly maintenance missing	YEARLY M OVERDUE 10083
QC out of range	QUARTERLY M OVERDUE 10073
QC missing	WEEKLY M OVERDUE 10063
Measurement discarded	QC OUT OF RANGE 10053
Temperature low during measurement	QC MISSING 10043
Temperature high during measurement	MEAS_DISCARDED 10033
	TEMP_LOW 10023
	TEMP_HIGH 10013

Table 5-5: Error codes

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5.7.2 Errors during pipetting

Errors	Status bar	Channel	Repair
Start button is not active	blue	initialising	Wait 30 sec until status line turns grey.
	stays blue	defect	Call technical service provider.
Message window: "Error of pipetting programme occurred, pipette will be checked. Please cancel test. Please push **Reset** at the back side of the pipette."			Push the small blue button at the back side of the pipette. Pipette is initialising shortly (ejecting movement). Repeat pipetting sequence.
The pipetting programme is not available after starting.			Abort sequence (x on screen). Wait until status line turns grey. Repeat **Start**.

Table 5-6: Errors during pipetting

5.7.3 Incorrect TEMograms

Errors	Repair
The TEMogram is interrupted and displayed again after a few minutes.	No warning signal. Check the temperature in the SERVICE module. If the temperature differs more than 3°C from the adjusted target temperature (mostly 37°C), no measurement data is recorded. Check the service values in the SERVICE module (Chapter 8.2.3). Contact the technical service provider.
Other abnormal TEMograms	A warning message is displayed for the affected channel. Follow the instructions of the warning menu.

Table 5-7: Incorrect TEMograms

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6 System software

The ROTEM® *delta* software provides you with intuitive user guidance when performing, evaluation and managing measurements. General operating elements and the modules are described elaborately in this chapter. The main screen is described in Chapter 4.3.

6.1 Operating elements

On-screen keyboard

When data needs to be entered in a menu of the software (e.g. patient data), a touch on the field to be filled out will display the on-screen keyboard.

Track ball

The user interface is optimised for touch screen operation in connection with a trackball (comparable to a mouse) on the top right side of the keyboard and the corresponding buttons on the top left side of the keyboard.

Keyboard

Additionally, some keyboard shortcuts facilitate the navigation in the program (Table 6-1).

Shortcut	Function
<**Tab**> or <**Enter**>	Move forwards between sub menus, buttons and entry fields
<**Shift**> + <**Tab**>	Move backwards between sub menus, buttons and entry fields
<**left arrow**>, <**right arrow**>	Navigation within an entry field with a date

Table 6-1: Shortcuts for the navigation in the programme

On the screen, list boxes to select elements like the list box COMMENT, where you can select an appropriate comment, are displayed.

Table 6-2 gives an overview about keyboard shortcuts that facilitate the work with list boxes.

Shortcut	Function
<**Shift**> + <**Space**>	Open highlighted list box and show entries
<**Enter**>	Select highlighted element
<**arrow down**>, <**arrow up**>	Scroll in list boxes
<**Esc**>	Close list box

Table 6-2: Shortcuts for the navigation in list boxes

Buttons

Buttons can be controlled by touching the screen, with the trackball and with its corresponding buttons and with <**Enter**> or <**Space**>. With the <**Tab**> key you navigate to the required button.

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6.2 Measurement module

The measurement module consists of three sub areas: preparation, measurement and graphical representation of the measurement results. The following sections describe the user interface of these three sub areas.

This chapter solely deals with the user interface. Chapter 5 gives a detailed description of the measurement.

6.2.1 Preparation

When opening the measurement module the start screen is opened (Figure 6-1) and the channel check is displayed. Touching the button **Preparation** returns from other windows of the measurement module back to the preparation window.


Figure 6-1: Measurement module – preparation window

The preparation window shows in the upper area, the steps to carry out before starting a measurement. The lower area shows the four channels.

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6.2.2 Measurement

Touching one of the channel fields opens the measurement window with its six areas (Figure 6-2), (Table 6-3):


Figure 6-2: Measurement module – measurement

	Area	Content
1	Action area	Entry fields for patient data (upper half of the screen) and as well as the buttons **Reset**, **Copy, Copy to all** and **Paste**. While running a measurement: current TEMogram is displayed.
2	Main menu	Buttons **Preparation** (or **ROTEM® *platelet*** in case a ROTEM® *platelet* is connected), **Multi-TEM**, **Screenshot**, **Profile overlay**, **Patient overlay**, **Help** and **Quit**. The button **Warning** is only visible when errors occur.
3	Submenu	Functional buttons depending on state of measurement: During start of test: Pipetting instructions (upper half of screen). While running a measurement: Numeric parameter values. The order of the parameters is determined in the menu SETUP - PARAMETER SETUP (Chapter 6.5.2). With the arrow keys on the right, it is possible to scroll in the list.
4	Status line	Date and time, software version, user, target temperature as well as temperature indicators for the work area and the channels. The colour of the test indicators is blue (too cold), white (ready for use) or red (too warm).
5	TEMogram display	Channel 1 to 4 from left to right.The selected and therefore active channel is highlighted in blue, the appropriate parameter values are shown under the particular channel.
6	Test line	Colour coded selected test. Selection of other tests by touching the line.

Table 6-3: Measurement module – sub areas of the measurement window

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6.2.3 Multi-TEM

The button **Multi-TEM** opens an enlarged view of all four channels. This view is especially suitable for remote transmission of the results (please refer to the ROTEM® Secure Viewer manual).

Touching any TEMogram once on the screen displays the measurement window again.

Touching any TEMogram twice on the screen displays measurements run up to 180 minutes.

6.2.4 Screenshots

The measurement and service modules provide the possibility to capture the current screen content as a graphic file in .jpg format. These screenshots can be viewed and processed in every graphical software that can read .jpg files.

- Select **Screenshot**.
 ↳The file is saved in the folder FILESYSTEM - USER - SCREENSHOTS under the name Screenshot_yyyymmtt**T**hhmmss.jpg (e. g. Screenshot_20070413T174652.jpg). These files are automatically deleted 60 days after creation. Chapter 6.7.2 contains further information regarding the FILESYSTEM.

6.2.5 Overlays

With **Profile overlay,** a TEMogram within the normal range of the current patient profile is displayed over the measured TEMogram.

By using the **Patient overlay** button, the 4 newest TEMograms of the same test and from the same patient ID in the database (of the last 48 hours) are displayed. With each overlayed TEMogram, the starting time of the according measurement is displayed in colour.

6.3 Database

All measurements (patient and QC) are saved in the database. They are available for later retrieval and for export to other programs. Measurements may be selected manually or using a search form (Chapter 6.3.2), (Chapter 6.3.3).

Note: To avoid the loss or unwanted editing of data, it is not possible to delete, change or add data sets to the database. Measurements that have been previously discarded (instead of saved/cleared) in the measurement menu are displayed grey but still present.

Note: Measurements with flags are highlighted in yellow.

6.3.1 Database user interface

The database module contains 2 databases:

1. the patient database
2. the QC database

The required database can be selected by toggling the button **Patients/QC**

Note: In case a ROTEM® *platelet* is connected to the ROTEM® *delta*, the ROTEM® *platelet* database can be selected with **Platelet** (see ROTEM® *platelet* operating manual).

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The database module user interface consists of five areas (Figure 6-3), (Table 6-4):


Figure 6-3: Database module

	Area	Content
1	Database buttons	Selection of QC or patient database, show or hide discarded measurements (if connected, ROTEM® *platelet* database)
2	Database window	List of performed patient or QC results
3	TEMogram window	According to the selection in the database window, either one single or overlayed TEMograms. The time pattern can be switched with the right trackball key from a 10 to a 60 minute interval.
4	Search form	Entry fields for the filtering of measurement results Search buttons
5	Selection screen	List of selected measurements

Table 6-4: Database module

Additionally there are further buttons for the handling of the database module available at the right margin. These are described in the particular sections.

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6.3.2 Manual selection of patient / QC results

Sorting data sets

By default, the data sets are sorted according to the start time of the measurements, beginning with the newest one.

With a click to the column caption, a sorting criterion is determined. With another click on the same column, the order is inverted.

The view selection can be adjusted with the horizontal and vertical scroll bar.

Hiding/showing discarded measurements

With the button **Hide discarded** it is possible to hide all the discarded measurements from the database view. With the button **Show all** they are shown again.

Highlighting data sets

A data set is highlighted by touching its line. Sliding down will highlight the data sets below. Pressing the <**Ctrl**> key allows to select multiple, non successive data sets.

↳The lines of the data sets are highlighted.
↳The selected measurements are shown overlayed in the TEMogram window.
Proceed as described in Chapter 6.3.4 for further processing of the selected results.

6.3.3 Selecting measurements (patients or QC) using the search form

The database module provides the possibility to filter the measurements and to select them for further applications, using the search form.

It is possible to enter wildcard characters "?" for single characters and "*" for multiple characters in the search form of the database. For example searching for "M??er" retrieves the results for Meier, Meyer, Maier etc. Searching for "*er" retrieves the results for all patient names ending with "er".

1. In the patient database: Specify at least one of the search criteria *patient ID, sample ID, patient name, comment, lot, data range or test* in the SIMPLE SEARCH FORM.
2. In the QC database: Specify at least one of the search criteria *QC lot, expiry date, comment, data range or test* in the SIMPLE SEARCH FORM.
3. Select **Start search**.
 ↳The database is scanned for the required entries.
 ↳The results complying with the given criteria are highlighted in the database window.
 ↳The results are shown in the TEMogram window.

 Note: If the search does not provide the requested results, delete the sorting criteria you have entered using the button **Clear search**.

 Note: If necessary, selected QC measurements can be discarded with **Discard QC.**
 ↳The discarded QC measurements are displayed in grey.
4. Proceed as described in Chapter 6.3.4 for processing of the results.

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6.3.4 Making measurements (patients or QC) available for processing

1. Select **Add** button.
 ↳The highlighted measurement results are copied to the selection screen and are available for printing and exporting (Chapter 6.3.5), (Chapter 6.3.6).

Table 6-5 shows an overview of the buttons used in the selection screen.

Button	Function
Remove	Remove the highlighted search results from the selection screen
Reset	Remove all search results from the selection screen

Table 6-5: Buttons for selection screen

6.3.5 Printing

The settings for the printer are determined in the menu SETUP - PRINTER SETUP.

Patient measurements and QC results can be printed from the respective databases.

✓ A printer is connected to the system.
✓ The measurements to be printed are copied to the selection screen.

1. Select **Print** button.
 ↳The SELECT PRINT STYLE window is opened.
2. Determine number of TEMograms to be printed per page.
3. Confirm with **OK**.

6.3.6 Exporting

A timeout can occur when transferring extensive data sets via network to an external information system. Do not export more than thousand measurements at once.

✓ The measurements to be exported have been copied to the selection screen.

1. Select **Export** button.
 ↳A configuration window is shown (Figure 6-4).
2. Determine export destination:
 - Online (Section ***Online export)***:
 This export function allows measurement results to be batched and transmitted to the laboratory information system (LIS) if no real-time automatic transmission is possible.
 - Hard drive (Section ***Hard disk)***:
 Measurement results can be exported to the FILESYSTEM on the hard drive. From there the results can be transferred to a USB stick and made available for evaluation in other software like Excel®. This also acts as a data backup feature (Chapter 6.7.1), (Chapter 6.8.2).

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
Figure 6-4: Export options

Online export

Appendix E gives a detailed description of the data transmission protocol and possible error sources.

✓ Export destination ONLINE is selected.

1. Confirm with **OK**.
 ↳The settings are saved.
 ↳The window is closed.
2. Select **Export**.
 ↳The measurement results are sent to the interface for transmission to the LIS.

Hard disk

✓ Export destination HARD DISK is selected.

1. Mark the option EXPORT WITH TEMOGRAM, when not only the numeric measurement results but also the corresponding TEMograms should be exported.
 ↳The two options for selecting the graphical format are active (Figure 6-5).

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
Figure 6-5: Export options for TEMograms

2. Select requested graphical format:
 - .bmp for maximum resolution (450 kB/file)
 - .jpg if multiple TEMograms should be exported (70 kB/file)

 ↳A window for file name determination is opened.
3. Enter a new name for the export file if necessary.
4. Confirm with **OK**.

 ↳The measurement results are exported.

 ↳A pop up confirms the successful export.

Name conventions during the export of the measurement results to the file system

The numeric measurement results are saved in a .txt file, the TEMograms are saved in separate single files, in a folder with the same name as the .txt file.

By default the system assigns names for the .txt files according to following conventions:
TEM-Export date (yyyy-mm-tt)_Export time (hh-mm-ss).txt
(e.g. TEM-2008-06-18_13-44-53.txt).

The file name may be changed but must have the extension ".txt" for the import into Excel®.

The file names for the single TEMograms are composed of:
Patient name (15 characters)_Name of the test_Date (yymmtt)time (hhmmss).jpg or
Patient name (15 characters)_Name of the test_Date (yymmtt)time (hhmmss).bmp.
e. g. Miller-John_INTEM_080618134453.jpg.

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Exported parameter of the numeric measurement results

With the export of numeric measurement results, 43 parameters are exported (Table 6-6). If TEMograms are exported additionally, the 44th parameter is the name of the folder, containing those TEMograms.

No.	Content	No.	Content	No.	Content	No.	Content
1	Patient ID	12	CT	23	LI30	34	AR25
2	Sample ID	13	A5	24	LI45	35	AR30
3	Patient name	14	A10	25	LI60	36	MCE
4	Sex	15	A15	26	MI	37	ACF
5	Birth date	16	A20	27	CFR	38	G
6	Comment	17	A25	28	LOT	39	TPI
7	Test name	18	A30	29	CLR	40	MAXV
8	Start time	19	CFT	30	AR5	41	MAXV-t
9	Runtime	20	MCF	31	AR10	42	AUC
10	Channel	21	MCF-t	32	AR15	43	LT
11	Error code	22	alpha	33	AR20		

Table 6-6: Exported parameters

Note: The exported .txt file is available for measurement result evaluation only. In case of multiple error codes during the measurement, only 1 error code is reported. Flags, normal range profiles and lot numbers (single use reagents and QC) are not reported. Error codes, flags and lots can be viewed in the database module.

6.3.7 Evaluating patient or QC results

The database module offers the possibility to view and compare the stored measurement results in the patient or the QC database.

✓ The required measurement results are copied into the selection screen.

1. Select **Accept**.
 ↳The measurement results are prepared for evaluation.

During the data transfer, you are informed about the progress in a status pop up. The transmission can be stopped with **Cancel**.

↳The measurement results can be interpreted in a new screen (Figure 6-6).

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
Figure 6-6: Measurement result evaluation screen

User interface of the evaluation screen

The evaluation screen consists of three areas:

- Area 1: Thumbnails of the TEMograms
- Area 2: TEMograms in larger view (1h or 8h view possible with a right click)
- Area 3: numeric measurement results (colour coded)

By using the **Back** button, the evaluation screen is closed and the first database screen is shown.

Selecting TEMograms

In area 1 the TEMograms for evaluation are selected:

- Touch the required TEMograms in area 1
 ↳The background changes from grey to white.
 ↳Explicit colours are assigned to the TEMograms.

If more then 4 TEMograms are selected, all other TEMograms are shown in black by default. With the **Colour** button at the lower part of the screen, further colours can be defined and assigned (Section *Define and assign colours)*.

↳ In area 2 the TEMograms are shown in these colours.
↳The numeric measurement results are shown in area 3 and are highlighted with colour bars in these colours.

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Define and assign colours

Additionally to the automatic assignment of colours, there are multiple possibilities to assign colours to TEMograms:

1. Select TEMogram to which the colour is to be assigned.
2. Select **Color** button.
 ↳The colour selection window is opened.
3. Select one of the basic colours

 or

 Add a user-defined colour
4. Confirm with **OK**.
 ↳The colour is assigned to the TEMogram.

Evaluating numeric measurement results

According to the settings in the menu SETUP - PARAMETER SETUP - LARGE TEM, the numeric measurement results are shown in area 3 (Chapter 6.4.6).

The heading shows the patient name and the executed test from the measurement .

If normal values have been entered for the parameters, all values outside the reference range are highlighted with a red background.

If multiple TEMograms are selected, it is possible to scroll between the different measurements with the arrow keys at the bottom of the screen.

Evaluating TEMograms

For research purposes, it may be necessary to overlay not only multiple TEMograms but also (with the appropriate buttons) their first or second derivative.

Table 6-7 shows possible combinations for the overlay of curves.

Number of selected TEMograms	TEMogram	1. Derivative	2. Derivative
1	X	X	X
> 1		X	

Table 6-7: Combinations for overlay

The first and second derivative are scaled proportionally to the scale of the firmness curves (y-axis). Both scale factors are shown at the bottom left of the TEMograms.

Printing TEMograms and overlays

The evaluation screen offers the possibility to print TEMograms and overlays.

For the purpose of clarity, the number of printable curves is limited to four curves.

✓ At least one but not more than four curves are selected (TEMograms and/or derivative).

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or

- Select **Print 4 TEM** and enter title.

Figure 6-7 shows an example from a printout of 4 overlay TEMograms.


Figure 6-7: Printout of 4 overlay TEMograms

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Export TEMograms and overlay

Selecting **Export overlay** or **Export 4 TEM** copies the curves as .jpg files into the FILESYSTEM. From the file system they can be transmitted to a USB stick.

6.4 Setup

The ROTEM® system has preset factory settings and is ready for testing after installation. Some ROTEM® *delta* settings can be changed according to internal standards.

The authorisation to change settings require appropriate user rights. If this user rights are not given, the settings can only be viewed.

1. Select **Setup** button in the main menu.
 ↳ The ROTEM® SETUP module is opened (Figure 6-8).
2. Customise settings. The setting possibilities are described in the designated sections.
3. Save changes using the **Quit** button.
4. Confirm to save the changes.


Figure 6-8: Setup - User setup

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6.4.1 User Setup

The following entry and list fields are displayed in the upper area of the screen:

Field	Content
Device name	Name for identification on printouts when multiple devices are utilised
Location, institute	Information for the correct inscriptions on the printouts
Date, time	System date and time, save with **Set clock**
Temperature	Temperature of the sample preheating station and of the temperature controlled work area set to 37 °C. **Note:** The temperature set here has to be at least 5 °C above ambient temperature, otherwise the system cannot maintain the correct working temperature.
Max. measurement time	Runtime of complete patient sample measurement. After this time period, the measurement is stopped (normally < 2 hours, for research purposes only: in a range of up to 8 hours). **Note:** Independently from this setting, QC measurements are stopped after 60 min.

Table 6-8: General settings

The lower part of the screen allows additional configurations:

- In the area PIPETTE you can select between manual and automatic pipette. Normally the automatic pipette is used, the manual pipette is used for research purposes only.
- In the area ADDITIONAL SETTINGS I, you determine which patient data can be entered (AVAILABLE) and which has to be entered (MANDATORY).

 Note: At least one field has to be selected for the identification of the measurement in the database (preferably the *patient ID* or the *patient name*). Not selected fields are not available in the measurement module.
- In the area ADDITIONAL SETTINGS II, *Single Patient Mode, Results barcode, On-screen keyboard default*, *Identification name (for patient query), User Lockout* can be selected.
 - ↳with *Single Patient Mode* selected, patient data entered in the first selected channel, are automatically copied in all other channels.
 - ↳with *single patient mode* selected, the *Results barcode* can be activated and its update interval determined. The barcode includes the current parameter results from all channels as well as the ROTEM® serial number, model and patient ID.
 - ↳with *On-screen keyboard default* selected, the on-screen keyboard is always displayed when touching an entry field.
 - ↳with *Identification name (for patient query)* selected, LIS connected systems display a unique patient identification entry field in the MEASUREMENT menu. The name of this field can be entered in the SETUP menu field. It is always used in combination with the option *Patient data are automatically retrieved from HIS/LIS when patient ID is scanned* (Chapter 6.4.2) and *Single Patient Mode*.
 - ↳with *User Lockout* selected, only trained and authorized personnel can log in and the user is locked out when his training has expired.

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6.4.2 Online setup

ONLINE SETUP contains the settings for the connection between the ROTEM® *delta* and the laboratory information system (LIS).

The window consists of four areas (Figure 6-9):


Figure 6-9: Setup - Online setup

Transfer of measurement results (area 1)

In area 1 you can determine:

- if the measurement results should be transmitted automatically to the LIS after clearing the measurement (Chapter 5.4.3).
- whether the patient data is automatically retrieved from the LIS when the patient ID is scanned (only possible in combination with the ROTEM® Connectivity Kit and the ROTEM XMLv2 protocol).
- to use the NET OVERLAY feature which allows patient overlays with patient measurements run on a different ROTEM® system (only possible in combination with the ROTEM® Connectivity Kit or ROTEM® *live*).

Note: the overlay will be formed with data stored in the ROTEM® Connectivity Kit, not in the database.

- to allow a patient ID query from the ROTEM® Secure Viewer, showing the required patient on the right ROTEM® system (only possible in combination with *Single Patient Mode* and the ROTEM® Connectivity Kit).
- to make XMLv2 protocol selectable in area 4.
- if the serial port is used instead of the usual TCP/IP port.

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Note: The serial interface is only available on ROTEM® *delta* series 2000 and 3000.

Select transmission protocol (area 4)

In area 4 you select the protocol for data transmission. Currently following communication protocols are implemented (Appendix E):

- DOS compatible: exact copy of the protocol in the version of the ROTEG® 05 Software Version 2.98
- ROTEMXMLv1: With this protocol all demographic data and all numeric measurement values are transmitted.
- ROTEMXMLv1+jpeg: With this protocol additional to the demographic data and the numeric values, the TEMograms are transmitted.
- ROTEMXMLv2: With this protocol demographic data, numeric values and TEMograms are transmitted. This protocol is compatible with the ROTEM® Connectivity Kit (which includes ASTM low level, ASTM high level and HL7 protocols). Further more, several devices can be connected together, using only 1 LIS connection. For more information about this kit, please refer to its user manual.

Determining transmission parameters for serial interface (area 3)

In area 3 you can determine the transmission parameters for serial interfaces. Standard settings:

Parameter	Default
SERIAL INTERFACE	COM 1
CHARACTER SET	ANSI
BAUD RATE	115200
PARITY	none
DATABITS	8
STOPBITS	1

Table 6-9: Transmission parameters

Determining connection parameters for the ROTEM® Connectivity Kit or ROTEM® *connect* (area 2)

In area 2 you may enter the IP address used for the ROTEM® Concentrator (part of the Connectivity Kit, please refer to its manual) or ROTEM® *connect* (please refer to its manual).

Note: Do not change the port standard setting.

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6.4.3 Remote access

The system software of the ROTEM® *delta* offers the possibility to view the current measurements via the ROTEM® *live* or ROTEM® Secure Viewer to another computer within the same Local Area Network (LAN).

Contact your network administrator before activating ROTEM® Secure Viewer or ROTEM® *live* for the first time. Please refer to the corresponding manual.

In the menu REMOTE ACCESS settings for the connection can be found (Figure 6-11).

ROTEM® *live*

1. If the server is running, select **Stop server** button to change the settings.
2. Select ROTEM® *live*.
3. If required, adjust the the ROTEM® *live* update interval (default 10 s / adjustable from 5 - 10 s).


Figure 6-10: Setup – Remote access - ROTEM® *live*

Note: In case a measurement from a patient is saved/cleared on the ROTEM® *delta* while that same patient is being viewed on an external computer with ROTEM® *live*, the measurement that was just cleared will not be available for patient overlays in ROTEM® *live*. In order to become available, the patient search in ROTEM® *live* has to be repeated.

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ROTEM® Secure Viewer (discontinued)

1. If the server is running, select **Stop server** button.
2. Select ROTEM® Secure Viewer.
3. The ROTEM® Secure Viewer can be started with the displayed default settings. If only one computer shall be authorized for the remote view, enter the IP address of this computer into the entry field *IP range* and select the option SINGLE ACCESS.
4. If the connection shall be established automatically when starting the ROTEM® *delta*, check the relevant option.

Note: The ROTEM® Secure Viewer password can be changed in FILESYSTEM - ADMIN_TOOLS - DELTA SETTINGS. It can only contain lower case characters and numbers. It may not contain any upper case nor special characters.


Figure 6-11: Setup – Remote access - ROTEM® Secure Viewer

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6.4.4 Printer setup

There is the option to export the measurement results and interpretations to a printer or a file. Under the menu item PRINTER SETUP the according setting window opens (Figure 6-12).

Figure 6-12: Setup – Printer setup

Printing

✓ A printer is connected to the system.

1. Highlight the option CHOOSE PRINTER.
2. Select printer from the list.
3. Select printouts in COLOUR (recommended) or in GREY SCALE.

Including normal values (or target ranges) on the printout

✓ In the settings for the test (or QC test), normal values (or target ranges) are provided.

1. Highlight the option INCLUDE NORMAL VALUES.
 ↳The normal values (or target ranges) of the tests (or QC tests) are displayed together with the measurement values.
 ↳Measurement values outside the ranges are marked with an arrow up or down.

Including overlays on the printout

✓ Patient overlays have been selected during measurement.

1. Highlight the option INCLUDE OVERLAYS.
 ↳The Patient or Profile overlays of the measurements are displayed on the printout, when overlays are selected in the MEASUREMENT MODULE.

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6.4.5 Test setup

In the menu **Test setup** the lists with available tests can be configured and normal values can be entered.

Determining available tests

In the measurement module, the desired tests for each measurement are determined. In the menu TEST SETUP it can be determined which tests are available in this list and in which order they appear in the list.

This setting affects only the measurement module, not the database.

1. Select **delta Test selection**.
 ↳The window to select the tests opens (Figure 6-13).
2. Mark the desired test in the list TESTPOOL.
 ↳The selected test will be highlighted.
3. Move test with **<Arrow right>** to the list TEST CHOSEN.
4. Repeat steps 2 and 3 until the list TEST CHOSEN contains the desired entries.
 ↳These tests can be selected in the measurement module.
5. To move an entry within the list, highlight this entry and move it with the **<Arrow up>** or **<Arrow down>** to the required position.
 Note: The first entry in the list will appear in the measurement module on channel 1, the second on channel 2 etc. If less than four tests are selected, the tests are repeated in the set order in the following channels.

With the **Clear all** button, the list of the selected tests is cleared. With **<Arrow left>**, highlighted tests are removed from the list.

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
Figure 6-13: Setup – Test setup

Defining normal values

For the evaluation of measurement results, comparison with the reference ranges of the tests may be helpful. By default, normal values of the ROTEM® parameters are preset here. It is possible to enter normal values for different patient profiles (e.g. pregnant women, pediatric, ...). These normal values will be taken in account to form the profile overlay in the MEASUREMENT menu.

1. Select **delta Normal values**
 ↳The window to select the tests is opened (Figure 6-14)
2. Select or create the required patient profile (e.g. "pediatric") in the list CURRENT PROFILE:
3. Select the required test in the list NORMAL VALUES FOR:
4. Enter the reference values for each parameter in the columns LOWER LIMIT and UPPER LIMIT
5. Confirm each entry with <**Enter**> (next line) or <**Tab**> (next column)

By using the **Clear** button all normal values are deleted for the selected test.

6. Select **Enable overlay** if the entered normal values should be used for the profile overlay in the MEASUREMENT menu
 Note: for this feature, values for all parameters between brackets must be filled out. All values and the curve must be checked for correctness before confirmation
 ↳The patient profile can now be selected in the MEASUREMENT menu

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↳The corresponding profile overlay can be displayed over the running measurement
↳The profile name and the corresponding profile ranges are printed on the printout
↳The corresponding profile ranges are sent automatically to the LIS
↳The profile name is added to the export file

Creating a new patient normal range profile

1. Press **New**
2. Enter the new profile name (max. 10 characters)
3. Confirm with **OK**
4. Select the new name in the drop down list
5. Enter the new normal values for each required test as described above
6. Check the correctness of the entered values
7. Check the correctness of the curve
8. Select **Enable overlay** if the entered normal values should be used for the standard overlay in the MEASUREMENT menu
9. If required, preselect the current profile as default to be used in the MEASUREMENT MENU.

It is possible to create a new profile with a copy of existing normal values with **Copy.**
After selecting an existing profile in the drop down list, it is possible to **Rename** or **Delete** it.


Figure 6-14: Setup – Test setup - delta Normal values

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Note: If the instructions for use of the reagents contain updated reference ranges, a new patient normal range profile needs to be created and the corresponding new values have to be entered.

6.4.6 Parameter setup

When displaying the TEMograms, the appropriate parameters are shown in the order determined in the menu PARAMETER SETUP.

1. Select **Parameter setup**
 ↳The window for the determination of the order of the parameters for the SMALL TEM, the LARGE TEM, the MULTI-TEM and the PRINTER TEM is opened (Figure 6-15).
2. With **<Arrow up>** and **<Arrow down>**, highlight parameters and move as required.

Note: The maximum number of parameters shown at the screen or the printout is stated in brackets.


Figure 6-15: Setup – Parameter setup

6.4.7 Comment setup

Under the menu item COMMENT SETUP, additional information for reutilisation, like clinical studies and diagnoses, may be entered. These text items become then available in the measurement module in the list box COMMENT.

1. Select **Comment setup**.
 ↳ The window for the entry of comments opens.
2. Enter text (max. 50 characters).

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3. Confirm with <**Enter**>.
 ↳The text is adopted to the list.
4. With the <**Arrow up**> and <**Arrow down**>, highlight comments and move as required.

By using the **Delete** button, the highlighted comments are removed from the list.

6.4.8 QC setup

In QC SETUP - ROTROL (Figure 6-16), control material can be:

- scanned (control name, lot number, expiration date and target ranges are transferred to the system)
- added (the scanned control is added to the control list and becomes available in the MEASUREMENT menu)
- deactivated (when the control should not be available in the MEASUREMENT menu)
- activated again (in order to become available in the MEASUREMENT menu)


Figure 6-16: Setup - QC setup - ROTROL

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The list in the upper part of the screen displays all controls, their different lots and their actual status (activated or deactivated). Viewing just one control type is possible by selecting the concerned control type in the upper *Control* field. The list in the lower part of the screen displays the history of all actions (when, by which user etc) taken for each lot. Pressing the required column header, allows to sort the data in both lists. The lower list can be printed out with **Print Lot History.**

Entering a new control

- ✓ The control has not expired
- ✓ The control has not been entered yet

1. Press **SCAN**

2. Scan the control box or control sheet barcode by holding the barcode under the scanner and near the arrow on the screen

 ↳control type, lots, expiry date and TEM Code are displayed on top of the screen

3. Press **Add Control**

 ↳the control lot is now added to the upper list

 ↳in its STATUS column, it is shown as "ACTIVATED"

 ↳in the lower history list ACTION column, it is shown as "added"

 ↳the target ranges for the control lot are saved and can be viewed in SETUP, QC SETUP, TARGET RANGES

 ↳the lot becomes available for the user in the MEASUREMENT screen

 Note: scanned but expired controls will not be added to the list.

 Note: with every new ROTROL lot, the corresponding target ranges displayed in SETUP-QC SETUP-ROTROL have to be compared to the ones on the control sheet.

Entering a new control manually

In case the barcode scanner is not available, control data can be entered the following way:

1. Select control type
2. Type lot numbers, expiration date and TEM code (4 numbers on the right of the large barcode on the box or on the control sheet) in the corresponding fields

 Note: Each ROTROL batch has 3 lot numbers: 1 for the complete kit, 1 for the lyophilisate and 1 for the diluent.

3. Press **Add Control**

 ↳the new control now appears in the upper list

4. Select the new control

5. Press **Target ranges**

6. Select a test in the drop down menu
7. Enter the lower and upper limit for each test result parameter (to be found on the control sheet in the ROTROL box)
8. Repeat steps 6 and 7 for the other tests

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QC reminder

With the QC reminder, the user is reminded automatically to run quality control in the MEASUREMENT menu.

1. Select the QC reminder interval (daily or weekly on a particular day)
2. Select the time the user wishes to be reminded
3. Select if the system should be blocked for patient measurements after 3 consecutive QC runs with out of range results.
 ↳the system will only allow to run patient measurements again after a successful QC run with in range results

 Note: This feature can only be enabled or disabled by service engineers with the appropriate user rights.
4. Check the QC tests for the different channels in the proposed combinations.
5. If required, select other QC tests.

 Note: in case <4 QC tests should be run, select NO TEST when appropriate

 Note: per active combination at least one test needs to be selected
6. If required, deactivate any unwanted combinations

 ↳the QC reminder automatically reminds the user in MEASUREMENT when the next QC run is due.
 ↳the QC reminder will propose each combination alternately
 ↳the user always has the option to swich to the other combination.

In order to comply with local quality control regulations, it is strongly advised to leave the QC reminder activated!

Only service engineers with the appropriate user rights can disable the QC reminder.

↳the user will not be prompted automatically to run controls
↳patient results run after out of range QC will not be flagged
↳the user is responsible to run the controls timely and to check all QC results

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6.5 Service

The authorisation to change settings requires appropriate user rights. If this user rights are not given, the settings can only be viewed.

6.5.1 Information

The menu item INFORMATION contains a list of the installed software components and features installed on the system.

The menu also indicates the QC status:

1. OK: all QC tests have been successfully performed
2. Out of range (1x, 2x or 3x): at least one parameter from the last (1, 2 or 3)QC combination(s) is out of range
3. Missing: QC has not been performed yet
4. Blocked: the system has been blocked after 3 consecutive QC runs with out of range results

Figure 6-17: Service – Information

6.5.2 Settings

In the menu SETTINGS, the graphical illustration of the recorded data streams allows an evaluation of the raw data generation. If the curve is displayed in red (superimposed on black), the data stream is recorded correctly. If the red and black lines are not superimposed exactly, it is a sign of a technical defect.

Note: the technical values are balanced and ready for evaluation about 3 minutes after entering the screen.

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Calibration values (gains) and temperature correction values (offsets) of the channels are factory set. Phase shift, amplitude, center, variance and period are parameters providing information on the mechanical and optical functions of the system. They are helpful for experts to judge if ball bearings or other components of the system are in need of exchange or repair (Chapter 8.2.3).

Screenshots can be made by using the **Screenshot** button (Chapter 6.2.4).

The button **Print Screenshot** captures the screen on a printout that can be forwarded to the ROTEM® technical service for further evaluation.

The buttons **ADM 1, ADM 2, ADM 3 and ADM 4** show technical data concerning each channel module.

Figure 6-18: Service – Settings

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6.5.3 Log book

The log book contains 4 options:

- The user can view the **History** of service interventions (e.g. repairs, SW installation, maintenance...) (Figure 6-19).
- In **Datalogger,** the user can view trends of the service values (variance, center, amplitude and period). A report of this statistical information can be sent to the service representative (**Export**).
- The service engineer can enter his service actions under **Report.**
- The **Reminder** informs the user when the next yearly preventive maintenance is due. The date can be set by the service representative or is set automatically after adding a *yearly maintenance* in **Report**.

 ↳The user will be informed through a reminder in the MEASUREMENT menu one month before the yearly preventive maintenance is due.

 ↳All measurements run when the yearly maintenance is 30 days overdue will be flagged and coloured yellow in the database and flagged on printouts.


Figure 6-19: Service - Log book - History

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6.6 Liquitrans

The LIQUITRANS function helps to reconstitute reagents (see according instructions for use). It allows the pipetting of any volume between 10 and 2000 µl with the automatic pipette.

1. Select **Liquitrans**.
 ↳The window for the entry of the pipetting volume opens (Figure 6-20).
2. Enter pipetting volume in µ litres.
3. Confirm with **Start**.
4. Pipette the volume according to the instructions on the screen.
 Note: Volumes up to 300 µ litres may be pipetted in one step, larger volumes are divided into more steps.
5. Complete action by pressing **Done**.

Figure 6-20: Liquitrans – pipetting volume


Figure 6-21: Liquitrans – process

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6.7 File system

All system, patient and measurement data as well as administration tools are stored on the hard drive in a FILESYSTEM, which is divided into ADMIN TOOLS and USER.

The authorisation to work with the admin tools requires according user rights.

6.7.1 Admin Tools

The ROTEM® *delta* offers 11 tools for the administration of the system.

- DB-Backup
- DB-Restore
- Delta Settings
- Feature Activation
- Keyboard Layout
- Printer Setup
- Remote Backup
- Set your Password
- Touchscreen
- Show Logs
- User Management

DB-Backup

A database backup should be performed at regular intervals.

Please follow the local rules for backup frequency. The two last sequential backups should always be stored for safety reasons.

Automatic bakups can be performed using the **Remote Backup** tool.

1. Select **DB-Backup**.
2. Insert USB stick.
3. Enter *START BACKUP*.
4. Confirm with <**Enter**> or **OK**.
 ↳A backup of the database and all settings of the device is made.
 ↳A window informs you about the successful backup.
 ↳The backup is saved under the name "Backup_Serialnumber.tgz" on the USB stick.
5. Finish action with <**Enter**> or **OK**.
6. Remove USB stick.

Note: Backups of identical names are overwritten. Transmit the backup files from the USB stick to a computer and rename them, e. g. in "Backup Date_Serialnumber.tgz".

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DB-Restore

The function DB-RESTORE permits restoration of a saved database.

Restoration of the database (DB Restore) will overwrite the existing data and is always associated with a loss of the data generated between the last backup and the restoration. Therefore it should only be applied if absolutely necessary.

Databases can only be restored if the backup file has its original name "Backup_Serialnumber.tgz". If necessary rename the file, remove the USB stick and reconnect it.

- ✓ The backup file to be restored belongs to the same ROTEM® *delta* system.
- ✓ If the backup file is not from the same system (e.g. replacement system during repair period), the file has been taken from a ROTEM® *delta* with the same SW version.

1. Touch **DB-Restore**.
2. Insert USB stick.
3. Enter *START RESTORE*.
4. Confirm with <**Enter**> or **OK**.
5. Select one of the following ROTEM® restore options:
 - Restore all

 Note: this option also includes data from FILESYSTEM-SHOW LOGS and SERVICE-LOG BOOK-DATALOGGER
 - Restore all except device settings
 - Restore medical data only
 - Restore user management only
 - Restore medical data of backup [1.5.X/1.6.X]
6. Confirm with <**Enter**> or **OK**.
 ↳The database is restored.
 ↳A window informs you about the successful saving.
7. Finish action with <**Enter**> or **OK** .
8. Remove USB stick.
9. Reboot system.

Delta Settings

Under this menu item you can change some instrument settings of the ROTEM® *delta* (according user rights required).

Note: All changes require a reboot of the system to become effective.

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Note: The ROTEM® remote view password can only contain lower case characters and numbers. It may not contain any upper case nor special characters

The instrument settings may only be changed by experienced users.

Feature Activation

This menu will allow to activate/deactivate barcode scannning and cup holder recognition on the system. A service engineer or an administrator will be able to access this protected menu with a special activation key.

Keyboard Layout

This menu allows the selection of different types of (country specific) keyboards.

Printer Setup: Install new printer driver

For printers provided by Tem Innovations GmbH, the required printer driver will be found under the drivers installed by default (Chapter 6.4.4).

If a required driver is not listed, please request information from the manufacturer of the printer for an alternative driver, search the database under http://www.linuxprinting.org/. Always contact your local service provider before purchasing a new printer.

Proceed as follows to install a new printer driver.

1. Touch **Printer Setup**.
 ↳ The KDE control module starts.
2. Click **Add** button.
3. Select menu item ADD PRINTER/CLASS.
 ↳The installation wizard starts.
4. Click **Next**.
5. Select interface according to printer connection.
6. Click **Next**.
 ✓ A local printer is selected.
7. Select USB PRINTER or PARALLEL PORT .
 Note: Some printers have both interfaces. If a USB printer is installed, only the USB connection can be used (and vice versa).
8. Click **Next**.
9. Select manufacturer and type of the printer.
10. If required for your printer, select POSTSCRIPT PRINTER or RAW PRINTER.
11. Click **Next**.
 ↳The installation of the printer driver is finished.
12. Click **Settings** and make required changes.
13. Confirm entries with **OK** and the following windows with **Next**.
14. Assign name of the printer in the window GENERAL INFORMATION.
15. Click **Next**.

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16. Check entries in CONFIRMATION window.
 Note: If further amendments are required, click **Back**.
17. **Finish** action.
18. Define printer as standard printer:
 - Click with the right trackball button on the name of the printer.
 - Select SET AS USER DEFAULT .
 - Click with the right trackball button on the name of the printer.
 - Select SET AS LOCAL DEFAULT .
19. Confirm with **OK**.

Remote Backup

This feature allows automatic backup transfer to a remote computer via the network. A backup is automatically transferred when shutting down the system and whenever the screen saver is activated for 4 hours.

Note: During shut down, it is possible to cancel the automatic backup.

Note: When the screen saver is active for 4 hours, the remote backup is only performed if new (cleared or discarded) measurements are available in the database.

To setup this feature, configuration by IT personnel is required:

✓ The ROTEM® *delta* is connected to the local area network.

1. Create a password protected network folder in your local network.
2. Give access rights to the users who may access this folder.
3. Configure the REMOTE BACKUP settings.
 - Enter the *user name* of the network folder
 - Enter the *password* to access the network folder
 - Enter the *path* to find the network folder in the following format: // IP address / folder name

 Note: use a "slash", not a "backslash"
 - Enter the *number* of backup files (per ROTEM® *delta*) you wish to be saved in the network folder

 ↳the oldest backup file will be overwritten, once the maximum number is reached
 - Select **Save and exit.**
 - Confirm with **OK.**
4. Reboot the system a first time.
 ↳The configuration changes are saved.
5. Reboot the system a second time.
 ↳The remote backup is performed.
6. Access the network folder on the remote computer.
 ↳The backup file (e.g. backup-SN-date-time.tgz) is now present in the folder.

If the Remote Backup feature is not required, the feature needs to be deactivated.

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To deactivate the feature, configuration by the administrator is required:

1. Select the REMOTE BACKUP menu.
2. Remove the entry from the field *path*.
3. Select **Save and exit.**
4. Confirm with **OK.**
5. Reboot the system.
 ↳any configuration change is saved.

Touch screen

If the entry fields of the touch screen cannot be operated smoothly, the screen has to be recalibrated.

1. Follow the instructions on the screen.

Show logs

The log in list gives information about the people who used the equipment, any configuration changes and backup or restore procedures at a given time.

With the button **Filter**, it is possible to sort the logged data.

With **Export to USB-Stick,** the data are copied directly onto a USB stick.

Set your Password

Each user may change his password in this menu.

1. Type the current password and confirm with **OK**.
2. Type the new password and confirm with **OK**.
3. Type the new password again and confirm with **OK**.
 ↳A status message confirms the successful change.

In case a password has been forgotten, please contact the administrator. Only the administrator can reassign new passwords.

User management

Under USER MANAGEMENT, new users and their passwords can be set up. Passwords of existing users can be reset. Table 6-10 gives an overview of the hierarchy of user rights and the management of passwords.

In **Lockout** the administrator can set the training status for each user (date of last training, required training interval etc):

↳If the user's training is expiring soon, a message will be displayed upon log in.
↳If the user's training has expired, the user will be locked out after 30 days.

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With the button **Options** it can be defined if the password rules should apply (**ON**) or not (**OFF**). The password rule requires a password of minimum 8 characters in which at least three of the following four requirements have to be fulfilled:

- Uppercase
- Lower case
- Numbers
- Special Characters ($, #, *, etc.)

Furthermore the encryption and the password protection of all exported data can be selected. Either the use of the default password may be selected or a new password may be entered.

	(super)user	administrator	service1
Measurement	x	x	x
Database	x	x	x
Liquitrans	x	x	x
Setup (read only)	x	x	x
Setup (read & write)		x	x
Setup (system blocking after QC OR)			x
Service (read only)	x	x	x
Service (gains)		x	x
Service (read & write)			x
USB	x	x	x
FileSystem-User	x	x	x
FileSystem-Backup	x	x	x
FileSystem-Restore		x	x
FileSystem-Printer		x	x
FS-Deltasettings: everything except SN		x	x
FS-Deltasettings: serial number			x
FileSystem-Touchscreen calibration	x	x	x
FileSystem-Set your password	x	x	x
FileSystem-User Management		x	x

Table 6-10: Hierarchy of user rights

Note: Only administrators can add new users, set new passwords and remove users.

6.7.2 User

The area **User** contains the subfolders **DataExport**, **GraphExport** and **Screenshots**.

DataExport

All exported data is saved in this folder (Chapter 6.3.6). The complete file name is shown when you position the mouse pointer on the data symbol.

GraphExport

All evaluations, created in and exported from the evaluation screen of the database, are saved in this folder (Chapter 6.3.7).

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Screenshot

All screenshots are saved in this folder (Chapter 6.2.4).

Service

All exports from SERVICE-LOG BOOK are saved in this folder.

6.8 USB stick

The menu USB STICK offers the possibility to transfer data from the file system to a USB stick.

1. Connect USB stick to one of the USB ports from the device.
2. Select USB STICK.
 ↳ The USB conqueror is opened.
3. Select FILE SYSTEM - USER.
4. With the left trackball button, select the files or folders to be transmitted from the file system, hold them and drag&drop them with the trackball onto the USB stick.
 ↳A window opens and the user can copy or move the file.
5. Select one of the options.
 ↳ The file is transmitted.

In case the files you intend to export do not fit on one USB stick, a message appears to insert a new USB stick. If the copying has to be aborted, a list will be shown with all files that have not been transferred. Before starting the next data transfer, you will be asked whether to transfer the data not yet transferred or to delete them. But please note that after 60 days, data not yet transferred will be cleared automatically with no alert. They will however still be present in the database.

6.8.1 Update

The USB stick is also used to update the system with new software. Smaller updates are usually copied to a USB stick.

✓ A software update has been copied to the USB stick.

1. Connect USB stick to one of the USB interfaces.
2. Select USB STICK.
 ↳ The USB conqueror is opened.
3. Double touch the update file.
 ↳Update starts.
4. Please follow any operation instructions appearing on the screen.
5. Reboot the system after the update has been completed.

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6.8.2 Importing data from the USB stick into Excel®[1]

Note: Make backups on a regular basis. Not all system data can be reconstructed from exported data.

See Chapter 6.3.6 for more information on the data export.

Using the text conversion wizard is the most popular way to import files into an external programme like Excel®.

1. Open external programme on your computer, e.g. Excel® 2000.
2. Open menu **Data** - **Import external data** - **Import data**.
3. Highlight requested .txt file on the USB stick.
4. Click **Open** button.
 ↳The text conversion wizard is opened.
5. Adjust the settings according to Figure 6-22 to Figure 6-24 and confirm with **Next** or **Finish**.


Figure 6-22: Step 1 for the import of a .txt file to Excel®

[1] Excel® is a registered trademark of Microsoft Corporation 1985-2003

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
Figure 6-23: Step 2 for the import of a .txt file to Excel®


Figure 6-24: Step 3 for the import of a .txt file to Excel®

The exported files are now available for further processing in Excel®.

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7 Initiation and installation

7.1 Initiation

Before using the ROTEM® *delta* system for the first time, make sure you read and understand all information in the Operating Manual. Particular attention should be paid to Chapter 2, where important information is provided regarding risks, safety standards, operating conditions, and handling of reagents, samples and controls.

7.1.1 Open package

If your instrument is supplied and installed by a distributor authorised by Tem Innovations GmbH, you can skip this chapter and continue with Chapter 7.1.3. The ROTEM® *delta* system is supplied in two packages. Optional accessories are supplied in additional packages.

If damage is found to have occurred due to improper handling during transport, please inform your supplier or the shipper immediately.

Before installing the system, the package has to be opened and checked if the content is complete.

1. Open the upper side of the package in an environment suitable for operation.

 Note: The upper side is marked with the following symbol: 

2. To gain access to the equipment, remove the cardboard with the accessories and the foam parts.
3. Take the foam packaging by its cut-outs and lift the ROTEM® *delta* with the foam packaging out of the cardboard.
4. Take the equipment out of the plastic packaging with the drying agent.
5. Place the equipment on an even and clean surface.
6. Take power cable, operating manual, keyboard, keyboard foil, USB stick, USB cable, pipette cable, cup holder, MC Rod holder, pipette holder, MC Rod, Clean Ring Rod, tips, measurement cells (cup & pin) etc. out of the packaging.
7. Take the pipette (inc. quality certificate and accessories) out of the white packaging.
8. Check if the supply is complete using the checklist in Chapter 4.1.

7.1.2 Touch screen preparation

The touch screen and its cables are shipped in their original manufacturer packaging and have to be prepared for the installation to the system.

1. Take the touch screen and its foot out of the cardboard.
2. Assemble the parts according to the enclosed instructions.

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7.1.3 Site installation requirements

For the proper usage of the ROTEM® *delta*, the site has to fulfil the following requirements:

- Ensure that the on/off switch is always reachable and that the power cable can always be removed to disconnect the power socket in case of emergency.
- Place the ROTEM® *delta* on an even, clean and robust surface. The operation of the equipment on a robust trolley is possible but strong shocks during measurement should be avoided.
- Ensure that the place where the equipment is set up is capable of supporting the weight of the system (22 kg) and that at least two grounded power supplies are available (ROTEM® *delta*,optional printer).
- Avoid placing the equipment near open windows, radiators or other instruments that are producing heat or strong vibrations.
- Ensure sufficient ventilation of the ROTEM® *delta*, leaving a minimum distance of 20 cm between other instruments and walls.
- Ensure that the ventilation openings at the rear side of the equipment are not blocked

Additionally please note the instruction for ambient conditions listed in Chapter 2.3.

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⚠ CAUTION

Measurement errors

After the installation the functionality of the equipment is not guaranteed. This may cause measurement errors.

➥ Perform a control measurement after each installation of the equipment (Chapter 5.5).

⚠ CAUTION

Damage of equipment or loss of data

Inappropriate installation and operation of the ROTEM® may lead to damage of equipment or loss of data.

➥ Never connect pipette to the sockets of the PC on the rear side of the device which is identical in construction.

➥ Never remove the printer cable with the system switched on.

➥ Never press the main switch while software is running.

ROTEM® *delta* installation:

1. Connect the system pipette to the pipette connector on the right side of the ROTEM® *delta*:
 - Push down the plastic bolt at the connector.
 - Push the connector into the socket.
2. Place and screw pipette and Clean Ring Rod holder to right side of the device (Figure 7-1).

Figure 7-1: Adjusting of pipette and Clean Ring Rod holder

3. Push the pipette cable guide in the intended little borehole.
4. Place the pipette and the Clean Ring Rod in their holder.
5. Install, similarly to the pipette holder, the MC Rod holder on the left side of the device.
6. Place the MC Rod in its holder.

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7. Place touch screen on the ROTEM® *delta* and connect screen cable on the rear side of the device to the VGA socket.
8. Connect touch screen power cable with the power socket on the device.
9. Connect USB interface at the rear side of the touch screen with the USB cable with one of the USB ports ROTEM® *delta*.
10. Connect the barcode reader to one of the USB ports with the provided cable.
11. Connect the keyboard to one of the USB ports and place it on top of the device in front of the touch screen.
12. Place the keyboard foil on the keyboard.
13. Connect the country specific power cable to the rear side of the ROTEM® *delta* and connect the ROTEM® *delta* to the power socket.
 ↳110V or 220V is selected automatically.
14. Connect the optional printer to the USB port of the ROTEM® *delta* and then connect the printer to the power socket (See printer manual for details).
15. Place the reagent rack on its position on the front of the device.
16. Place the cup holders on the temperature controlled work area.
17. Switch on the ROTEM® *delta*, the (optional) printer and the touch screen.
18. If necessary calibrate the touch screen (Chapter 6.7.1).
19. If required configure the software (Chapter 6.4).

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8 Maintenance, cleaning, disinfection

8.1 Introduction

Maintenance must be performed in order to guarantee correct operation of the system!

The system reminds you automatically when to perform maintenance

The ROTEM® system requires only minor maintenance, if utilized according to the regulations:

- Cleaning of surfaces.
- Cleaning of the cup holder.
- Replacing the pipette filter.
- Cleaning of the pipette if necessary.
- Cleaning or replacing the pin remover.

Tem Innovations GmbH recommends to follow the maintenance schedule as outlined below (Chapter 8.2) to ensure serviceability.

⚠ WARNING

Danger of infection

Inappropriate handling of contaminated samples may lead to infections.

- Always wear gloves when cleaning the equipment.
- Always follow guidelines for infectious materials.
- Especially carefully disinfect surfaces that have been in contact with potential contaminated material.
- Follow local regulations for the cleaning of such surfaces.

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8.2 Maintenance plan

The maintenance plan contains daily, weekly, quarterly, semi-annual and yearly maintenance work scheme. The user is prompted to perform maintenance in the MEASUREMENT menu.

⚠ CAUTION

Damage of equipment

Inappropriate maintenance may lead to equipment damages through invading liquids.

- Switch off the device before performing maintenance.
- Ensure that no liquids are invading the equipment nor the keyboard.
- Use only damp, lint-free cloths for cleaning.
- Protect axis and ball bearings from cleansing liquids or disinfectants.

Table 8-1 shows various cleansers and disinfectants, tested and recommended for the different surfaces by Tem Innovations GmbH.

Before other disinfectants are used, it is recommended to contact Tem Innovations GmbH first or carry out appropriate compatibility tests.

Manufacturer	Name	Countries	Composition (active constituents)
Schülke&Mayr	Mikrozid-AF-Tücher	A, CH, D, I,	35% Propane-1-ol 25% Ethanol
Schülke&Mayr	Mikrozid AF Liquid	A, CH, CY, D, DK, EST, GR, I, IRL, L, LT, LV, PL, RUS, TR	35% Propane-1-ol 25% Ethanol
Schülke&Mayr	Terralin Liquid	D	35% Propane-1-ol 25% Ethanol
Schülke&Mayr	SEPTINOL SA	F	35% Propane-1-ol 25% Ethanol
Bode Chemie Hamburg	Bacillol®plus	All European countries including the Baltic states	40 g Propane-1-ol 20 g Propane-1-ol 0.1 g Glutaraldehyde

Table 8-1: Recommended cleansers and disinfectants

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Manufacturer	Name	Countries	Composition (active constituents)
Bode Chemie Hamburg	Bacillol®AF	All European countries including the Baltic states	45 g Propane-1-ol 25 g Propane-1-ol 4.7 g Ethanol
Bode Chemie Hamburg	Sterillium Tücher	All European countries including the Baltic states	45 g Propane-1-ol 30 g Propane-1-ol 0.2 g Mecetroniumetilsulphate
Sullivan-Schein Dental	Cavicide	USA	10-20% Isopropanol 1-5% Ethyleneglycolmonobutylether
Clorox Company DO NOT USE FOR TOUCH SCREEN	Clorox Regular Bleach	USA	5-10% Sodium hypochloride <1% Sodium hydroxide Use a 10% Clorox dilution or Clorox cloths for the disinfection of surfaces

Table 8-1: Recommended cleansers and disinfectants

8.2.1 Daily maintenance

Table 8-2 shows the maintenance to be performed daily

Device to maintain	Action
Outer surface of the ROTEM® *delta*	Clean and disinfect with lint-free cloth.
Cup holder	Cleaning and disinfection with lint-free cloth soaked with disinfectant. Cleaning in a dishwasher or autoclave affects the durability of the surface structure.
Pipette	Clean with a damp, lint-free cloth (Chapter 8.3.2).
Filter	Control and replace the filter if contaminated with blood (Chapter 8.3.1).

Table 8-2: Daily maintenance work

8.2.2 Weekly maintenance

In addition to the daily maintenance described in Table 8-2 , the pipette filter should be replaced at least once per week (Chapter 8.3).

8.2.3 Quarterly maintenance

In addition to the daily and weekly maintenance, the device temperature and CCD chip values have to be evaluated quarterly.

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Control of the device temperature

1. Switch on the ROTEM® system
2. Remove pins and cups from the axes.
3. Wait until the device reaches operating temperature (approx. 5 min).
4. Open menu SERVICE > SETTINGS.
5. Read and make a note of the temperature of the four channels. The temperature should be between 36.0 °C and 38.0 °C.

Inform the technical service personnel in any case if the values are out of the specification.

Control CCD chip values

The CCD chip values are a measure for the quality of the analytical system. Changes in this area refer to problems.

1. Open menu SERVICE > SETTINGS.
2. Remove pins and cups from the axes.
3. Wait until the curves appear and the values stabilize (approx. 3 min).
4. Read and take a note of the values for AMPLITUDE, CENTRE and VARIANCE.

If a printer is connected to the device, you can print these values with the button **Print Screenshot.** It is recommended to print out this report regularly and to file it. Further documentation and trend graphs are available in the menu SERVICE, LOG BOOK, DATALOGGER.

The values should be within the ranges below:

Parameter	Lower limit	Upper limit
Amplitude	4400	5100
Centre	-832	708
Variance	1	12000

Table 8-3: Acceptable range for the CCD chip values

In case the values for one particular channel are not within specifications, remove and clean the Clean Ring (Figure 8-1) around the channel axis as described below.

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
Figure 8-1: Clean Ring around channel axis

1. Power off the system
2. Place the Clean Ring Rod around the axis (Figure 8-2)

Figure 8-2: Fitting the Clean Ring Rod onto the axis

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3. Fit the two little pins of the Clean Ring Rod into the 2 little holes of the Clean Ring from the concerned channel (Figure 8-3)

Figure 8-3: Fitting the Clean Ring Rod's pins into the Clean Ring

4. Turn the Rod, unscrew and remove the Clean Ring (Figure 8-4)


Figure 8-4: Removing the Clean Ring with the Rod

5. Let the Clean ring soak in soapy water (Figure 8-5)

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Figure 8-5: Soaking the Clean Ring

6. Clean off any remaining dirt with a moist, cotton swab
7. Rinse with water
8. Let the Clean Ring dry
9. Clean the Clean Ring's position around the axis (in the measuring block) with moist cotton swab

 Note: do not spill any moisture inside the measuring
10. Let the Clean Ring's position dry
11. Replace the Clean Ring into its position around the channel axis with the Clean Ring Rod
12. Switch the system on
13. Check the CCD chip values again, as described above
14. Inform the technical service personnel if any values are still out of the specification

8.2.4 Semi-annual maintenance

Additionally to the daily, weekly and quarterly maintenance, the integrated pin remover on the cup holders must be replaced semi-anually or every 250 measurements, as the springs of the pin remover are subject to wear.

Replacement pin removers are available as accessories (Appendix G).

⚠ CAUTION

Measurement errors

Damage of the polished surface may lead to measurement errors.

- To remove the pin remover, always use the MC Rod.
- Do not use any metal tools on the polished surface of the cup holder.

Replacement of pin remover:

1. Push the flat part of the MC Rod into the opening at the lower side of the cup holder (Figure 8-6).

Figure 8-6: Remove pin remover

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Measurement errors

Damage of the sealing tongue may lead to measurement errors due to evaporation.

- Do not press or scratch the sealing tongue.
- Do not push the sealing tongue in place by brute force.
- Always lift the sealing tongue over the slot.

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2. Insert the sealing tongue deep into the central opening of the cup holder (Figure 8-7).


Figure 8-7: Insert new pin remover: Insert sealing tongue

3. Press the springs left and right at the pin remover carefully into their guides (Figure 8-8).

Figure 8-8: Insert new pin remover: Press the springs into their guides.

4. Press the pin remover until it clicks and the springs slide into their final position (Figure 8-9).


Figure 8-9: Insert new pin remover: Snap of the pin remover.

8.2.5 Yearly maintenance

The system reminds the user when to organize the yearly maintenance 30 days before the due date.

Maintenance of the ROTEM® *delta* system must be performed on a yearly basis by authorized service personnel. Please contact your local service provider for yearly maintenance or ask for a maintenance contract.

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If the service representative does not document in **SERVICE-LOG BOOK- REPORT** that yearly preventive maintenance was performed, all measurements performed 30 days after the due date will be flagged and highlighted in yellow in the database (yearly maintenance overdue). In the MEASUREMENT menu, *Service!* lights up in red in the lower left corner of the screen.

8.2.6 Every 3 years

The system internal hard disk should be exchanged every 3 years by authorized service personnel.

8.3 Maintenance of the ROTEM® electronic eLine pipette

Please refer to the maintenance plan for the periodic maintenance of the pipette (Chapter 8.2).

After the maintenance of the pipette, the pipetting accuracy has to be checked (Chapter 8.3.5).

If the accuracy check isn't possible, at least two control measurements have to be performed with ROTROL N (in-tem® or QCinN; ex-tem® or QCexN). The measurement results have to lie within the target range of the parameters for ROTROL N.

8.3.1 Replacement of the pipette filter

- Remove the filter with the provided filter tweezers (Figure 8-10).
- Insert new filter with the provided filter tweezers.

Figure 8-10: Remove the pipette filter

Pipette filters are available as accessories (Appendix G).

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8.3.2 Cleaning of the pipette surface

Clean the pipette surface with a lint-free cloth and one of the disinfectants listed in Table 8-1. Alternatively also Ethanol (70 %), Isopropanol (60 %) or mild chemical cleansers may be used.

The cleansing liquid must not invade the gaps of the pipette handle! The liquid should not humidify the filter.

8.3.3 Cleaning and disinfection of individual parts

Dismantle pipette

1 Piston
2 Spring
3 Pipette head with cable
4 Tip cone
5 Tip ejector
6 Tip holder

Figure 8-11: Individual parts of the pipette

Dismantle pipette for cleaning and disinfection completely (Figure 8-11):

1. Disconnect ROTEM® *delta* from the power outlet.
2. Remove pipette cable:
 - Press the locking (round buttons on both sides).
 - Remove cover.

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- Remove plug.
- Mount cover again.

3. Unscrew and remove tip ejector.
4. Unscrew the tip cone.
5. Remove tip cone, tip holder and spring.
6. Remove pipette filter
 ↳Piston is exposed.
7. Unscrew and remove piston.

The pipette is prepared for cleaning and disinfection.

Cleaning

1. Clean individual parts with an appropriate disinfectant or mild chemical cleanser. Clean the interior of the tip ejector and the tip holder with a cotton cleansing tissue.
2. Rinse thoroughly with distilled water.
3. Let parts dry.
4. Grease the piston very thinly with the grease provided.

Disinfection

1. Place individual parts (except the pipette head) for at least 30 minutes into a beaker with appropriate disinfectant.
2. Rinse thoroughly with distilled water.
3. Let parts dry.
4. Grease the piston very thinly with the grease provided.

Figure 8-12: Disinfection of individual parts

Autoclaving

The individual parts of the eLine pipette (except the pipette head) are autoclaved separately or as a unit (121 °C, 1 bar over pressure for 20 min).

Adequate disinfectants are available from the company Biohit (e.g. Biocontrol Biohit REF 724004, 5 litres). For more details about suppliers in the specific countries please refer to http://www.biohit.fi.

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Before autoclaving, rinse the interior of the tip ejector and clean it with cotton swabs.

It is also possible to clean the individual parts and grease the piston before autoclaving.

Only use very little grease.

Utilise only the grease provided with the pipette.

Reassemble pipette

⚠ CAUTION

Damage

Inappropriate assembling may lead to damage of the pipette.

- Before reassembling, remove any lint or particles on the surface of the piston.
- When reassembling the pipette, do not overtighten the parts.

Reassemble pipette:

- ✓ All parts are dry.
- ✓ The piston is greased very thinly.

1. Screw in the piston (1).
2. Place the spring (2) on the piston.
3. Push the tip holder (6) into position.
4. Screw tip cone (4) back to its position.
 ↳Tip holder is fixed.
5. Insert new pipette filter using the filter tweezers.
6. Fix tip ejector (5).
7. Open cover from pipette head (3).
8. Connect pipette cable.
9. Close cover.
10. Connect pipette cable to the socket on the side of the device.
11. Start ROTEM® system.
12. Open LIQUITRANS menu.
13. Check pipette function by pipetting different small and large volumes.
14. Check function of tip ejector.

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8.3.4 Check and calibrate ROTEM® eLine pipette

The ROTEM® system may only be used with a calibrated eLine pipette. The calibrating of the electric pipette has been factory checked and certified at 22 °C using distilled water (grade 3), according to ISO 8655.

The specifications of the pipette are only guaranteed using the original Tem Innovations GmbH 320 µl pipette tips.

8.3.5 Functional check of the ROTEM® eLine pipette

Check the ROTEM® eLine pipette, e. g. every 6-12 month for accuracy and precision, regularly, according the internal regulations (complying with ISO8655).

Keep in consideration that the requirements for accuracy of the tests, usage frequency, number of users, characteristics of the liquids to be pipetted and the maximum frequency of acceptable errors (ISO 8655-1).

If a pipette has been dropped or stressed mechanically and inexplicably doubtful results have been measured afterwards, the pipette performance should be checked.

The pipette may be sent to an authorised service station for checking. Please contact your local ROTEM® supplier regarding the checking of the pipette.

Tem Innovations GmbH offers a description (according to ISO8655), on which basis you can check your pipette. If the pipette does not fulfil the specifications of accuracy and precision, it has to be calibrated.

8.3.6 New calibration

The pipette should be sent to an authorised service station for new calibration.

Please contact your local ROTEM® supplier regarding the calibration of the pipette.

Tem Innovations GmbH offers two kinds of calibration certificates:

- Quality Control Certificate
- Accredited calibration service (e. g. FINAS)

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8.3.7 Troubleshooting

Table 8-4 lists possible errors when using the pipette and their repair.

Symptom	Possible cause	Solution
Droplets left inside the tip	Non conform tip	Use original tips
Leakage or pipetted volume too small	Irregular wetting of the pipette tip	Attach new tip
	Tip attached loosely	Attach firmly
	Unsuitable tip	Use original tips
	Foreign particles between tip and tip holder	Clean the tip holder and attach new tip
	Pipette contaminated	Clean and grease the piston and the tip holder
	Insufficient amount of grease on a piston	Grease according to maintenance instructions
Pipetting volume out of specifications	Pipette damaged	Send pipette for maintenance to the Tem Innovations service
Pipette is blocked, aspired volume too low.	Liquids (blood) got into the tip holder and dried up.	Clean and grease the piston and the tip holder
Tip ejector is damaged and works irregularly.	Tip ejector or tip holder is damaged or contaminated.	Remove and clean tip ejector and tip holder.

Table 8-4: Troubleshooting

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9 Placing out of operation, transport and disposal

9.1 Placing out of operation

Disinfect the ROTEM® *delta* when placing it out of operation (Chapter 8.2) and protect it from dust until reutilisation.

9.2 Transport

⚠ CAUTION

Damage of device

Inappropriate transport may cause damage of the device.

- Do not lift the device by the pipette holder, axes or guiding rods.
- Lift the device only by holding the lower side.
- Transport the ROTEM® *delta* only in the original card board transport boxes.
- Transport touch screen separately.

9.3 Disposal

The life cycle of the ROTEM® *delta* system is five years when maintained regularly.

The following directive: *Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE)* applies to electrical and electronic equipment. Equipment that can be disposed according to this directive is marked with the WEEE symbol and may be disposed of according to local directions.

Figure 9-1:WEEE disposal symbol

For equipment which came into contact with blood or other biological contaminants prior to end-of-life, seperate directives are applied. This equipment is marked at the front side and has to be disposed of according to the directives effective for bio hazardous equipment.

Figure 9-2:Symbol for biohazard

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Further information on the disposal of peripheral equipment is listed in the according operating manual.

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10 Appendix

Appendix A: Marks, labels and symbols

Symbol	Description
CE	Mark according to EU guideline
2	For single use only
IVD	*In vitro* diagnostic agent
2016-08-01	Date of manufacturing e.g. here: August 1st 2016
2010-10	Expiry date yyyy-mm
~	AC voltage
LOT xxx	Lot number, batch number (e.g. xxx)
SN xxx	Serial number (e.g. xxxx)
REF xxx	Order number, article number (e.g. xxxxxx)
	1. Note operating instructions 2. Note the information in the paragraphs and chapters of the user manual marked with this symbol
Σ n	Kit sufficient for (n) tests
2 … 8°C 2°C 8°C	Symbols for storage indicated at a temperature range (e.g. +2 to +8 °C)
	Caution, fragile Fragile packaged goods
	Protect from moisture
	This side up during transport and upon opening
	Do not touch pin with hands
115V/T3,15AH 230V/T1,60AH	Fuse (technical specifications with 115V or 230V).
	Biohazard

Table 10-1: Markings, symbols and labels

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
Symbol	Description
	Minimum Distance of 1,5m / 4.9 feet between device and patient is obligatory
	Devices marked like this must not be recycled as unsorted municipal waste. Instead, these devices must be collected separately in accordance the local recycling regulations.

Table 10-1: Markings, symbols and labels

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Appendix B: Warranty

Tem Innovations GmbH declares to the original purchaser that the instrument sold by Tem Innovations GmbH or an authorized distributor of Tem Innovations GmbH shall be free of defects in material and manufacture.

Unless otherwise expressly agreed, the warranty period for new goods as well as for used goods shall be limited to the warranty period requested by the local legal authorities of the according country. In case the customer is a contractor the warranty period for parts subject to wear and tear shall be six months.

In case the customer is a contractor he shall be obliged to inspect received goods with respect to defects, guaranteed qualities, deficiencies and identity and to inform us within one week at the latest in writing. Hidden defects must be reported immediately upon their discovery. In case of exceeding those terms all warranty claims shall expire.

In case the customer is not a contractor he shall be obliged to inform Tem Innovations GmbH within a period of 2 months after discovery of the lack of conformity of the goods with the contract about apparent defects in writing; in case of omission of this information the warranty claims shall expire. This shall not apply in case of malice by seller. The customer shall carry the burden of proof for the time the defects were discovered.

In case the customer is a contractor, Tem Innovations GmbH shall in case of justified complaints to our discretion deliver a replacement or repair the goods or supply deficiencies. The customer shall grant an appropriate period of time for the repair of the goods or the compensation delivery. To a reasonable extent, repeated approaches for amendment shall be admissible. In case the customer should deny amendment, this shall void our warranty responsibility.

In case the customer is not a contractor he shall be entitled in the first instance to choose whether the supplementary performance shall be made as an amendment or replacement. However, Tem Innovations GmbH shall be entitled to deny the chosen supplementary performance if it would only be possible at unreasonable costs, provided the alternative way of supplementary performance causes no substantial disadvantages to the customer.

In order to enable supplementary performance, the goods must be returned to Tem Innovations GmbH's principal office. In case the notice of defect is justified, Tem Innovations GmbH shall deliver the goods free of charge after the supplementary performance. Otherwise the customer shall bear costs of freight.

In case our supplementary performance or replacement delivery does not lead to success within an appropriate period of time, the customer shall be entitled to demand decrease of the payment or rescission of the agreement. However, in case of an insignificant breach of the contract, especially in case of insignificant defects, the customer shall not be entitled to rescission.

Warranty shall not cover normal wear and tear, damages due to careless handling, excessive strain, use of non-original consumables and accessories, inappropriate equipment and chemical, electro-chemical or electric influences not to be assumed under the agreement.

In case of slight neglect of duty, Tem Innovations GmbH liability shall be limited to the average predictable, typical to contract and immediate damage for that kind of goods. This shall also apply in case of slight neglect of duty by our statutory representatives and vicarious agents. In case the customer is a contractor we shall only be held liable in case of slight neglect of duty to the extent violation of essential contractual obligations is concerned.

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The preceding limitations of liability shall not apply with respect to claims by customers for product liability. Furthermore, those limitations of liability shall not be applied with respect to damages to body or health or to loss of life of the customer to the extent imputable to us.

Wearing parts such as e.g.

- ball bearings
- complete axis

are subject to a reduced warranty of 6 months.

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Appendix C: Software licence

The ROTEM® system contains proprietary software modules of Tem Innovations GmbH (measurement module, database module, liquitrans module, service module, setup module, hardware-software interface, calculation module). Furthermore it contains free software/open source software. Only the free software/open source software is licensed under the free software licences. The source code of all modules licensed under free software licenses and the respective license texts can be found on the CD-ROM that is delivered along with the ROTEM®.

As long as the Software is only used as embedded software of the ROTEM®, there is no requirement to comply with the additional conditions as set forth by the free software licences. In order to obtain additional rights, the user may, at his option, agree to these license agreements with the licensors of the respective software modules.

Compliance with the license conditions of the free software licenses entitles the user to use the free software as provided by the relevant license. Any modification of the Free Software embedded in the ROTEM® system has the effect that any warranty for the ROTEM® system is lost and that the user becomes liable for any damage arising in connection with the use of the ROTEM® system. The ROTEM® system is not designed for any kind of modifications by the user, including but not limited to modifications of the embedded software. All public authorisations will automatically terminate with any intervention in the integrity of the ROTEM® system.

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Appendix D: Declaration of conformity

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Appendix E: Description of the interface

Tem Innovations GmbH reserves the right for technical amendments, if they are necessary for further development and technical advances. Please ask for the current version of the interface description before programming the access to the laboratory information system under following address:

E-Mail: support@tem-international.de

The ROTEM® program can transmit a part of the medical parameters of the ROTEM® measurements to a host computer (LIS) according to the protocol described below.

Communication path

Communication takes place via the RS-232 interface of the ROTEM® system or can also take place through the LAN port (TCP/IP) connecting to the ROTEM® Concentrator from the ROTEM® Connectivity kit.

Operating mode, only for the RS-232 interface

Interface parameter:	1200/9600/19200/57600 or 115200 baud, 8 data bits, 1 stop bit, no parity
Flow control:	Software (the RTS/CTS lines are not supported)

Settings in the ROTEM® program

Data transmission has to be activated and configured in the **Setup** - **Online** menu. For details of configuration see Chapter 6.4.2.

Activate the data transfer

Data transfer is exclusively activated on the ROTEM®. Data can be transmitted in real time or batch mode. For "real time" mode, please select the option "measurement results are sent automatically when measurement is finished". For "batch mode", do not select this option (Figure 6-9).

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Label	Relevance
<or>	logical “or”: The one *or* the other sign is sent.

Table 10-3: Logical combination

Communication protocol

Action	ROTEM® is sending	Host PC is sending
Enquiry from ROTEM®	<ENQ>	
Acknowledgement or no acknowledgement from host		<ACK> <or> <NACK>

Table 10-4: Communication protocol (1)

Depending on what the host PC is answering, the protocol proceeds.

- The host PC is denying with <NACK>
 - The protocol is interrupted.
 - A message informs you about denial of the measurement by the host PC .
 Note: you can still transmit the measurement from the database.
 - The transmission mode remains active. The ROTEM® tries to transmit the next completed measurement.
- The host PC affirms with <ACK>
 - Protocol continues as described in Table 10-5.

Action	ROTEM®is sending	Host PC is sending
Data from ROTEM®	<STX> <CR> <LF> <RS> <data record> <GS> <check sum> <CR> <LF> <ETX>	
Status of the check sum by host		<ACK> <or> <NACK>

Table 10-5: Communication protocol (2)

Depending on the confirmation of the check sum by the host PC, the protocol proceeds:

- The host PC does not confirm the check sum and sends <NACK>.
 - The ROTEM® repeats this process twice and waits for the confirmation of the host PC.
 - A third <NACK> response by the host PC leads to a final termination of the protocol and to a check sum error detection.
- The host PC confirms the check sum with <ACK>.
 - The protocol is concluded correctly.

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Structure of the data record ROTEMXMLv2

The data record is structured as follows:

```
<?xml version="1.0" encoding="UTF-8"?>
<Rotem>

 <Type>MEDREC</Type>
 <Version>v2.0</Version>

 <Device>
  <Name>Device Name</Name>
  <Location>Location</Location>
  <Department>Institute</Department>
  <Type>DELTA</Type>
  <SN>2115</SN>
  <IN>---</IN>
  <SVersion>FrontEndD 1.6.0</SVersion>
  <Lang>russian</Lang>
 </Device>

 <Patient>
  <PID>patientid</PID>
  <Surname>lastname</Surname>
  <Firstname>firstname</Firstname>
  <Birthdate>19700101</Birthdate>
  <Sex>M</Sex>
  <Measurement>
   <SID>sampleid</SID>
   <Starter>admin</Starter>
   <StartDate>20090416</StartDate>
   <StartTime>115026</StartTime>
   <RunTime>000001</RunTime>
   <Status>database</Status>
   <Channel>1</Channel>
   <Test>
    <Name>extem</Name>
    <Color>ff0000</Color>
    <Pipseq>---</Pipseq>
    <LOT>---</LOT>
   </Test>
   <ChannelErrors>
    <Count>1</Count>
    <Error>
     <Code>0</Code>
     <Message>OK.</Message>
     <Runtime>---</Runtime>
    </Error>
   </ChannelErrors>
   <MedicalComment>
    <Count>1</Count>
    <Comment>
     <Text>hyperfibrinolysis</Text>
     <Runtime>---</Runtime>
     <Writer>---</Writer>
    </Comment>
   </MedicalComment>
   <MedicalParameters>
    <Count>32</Count>
    <Paramètre>
     <Name>CT</Name>
     <Value>59</Value>
     <Status>final_ok</Status>
     <RefMin>38</RefMin>
     <RefMax>79</RefMax>
     <Unit>s</Unit>
```

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This protocol is fully compatible with the ROTEM® Connectivity Kit.

The Concentrator of this kit is capable to translate this protocol into ASTM and HL7. Further more, several devices can be connected together, using only 1 LIS connection.

For more information, please refer to the user manual of the ROTEM® Connectivity Kit.

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Appendix F: Passwords

Table 10-6 lists the ROTEM® standard passwords.

If the user management has been deactivated, some functions are protected with the following passwords:

Function	Password
Setup	admin
Service	admin (only for gain values, other settings can only be changed by service personnel).
DB backup	pentagamma
DB restore	pentagamma
Printer setup	admin
Delta settings	pentagamma
Access to ROTEM® screen from a remote access PC via ROTEM® Secure Viewer.	rotemremote

Table 10-6: Standard passwords

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Appendix G: Accessories

Following components and accessories can be ordered with their according catalogue numbers (Table 10-7).

Article	Description	REF No. Tem Innovations
Accessories box	Additional rack for tips and measurement cells	400472
Barcode scanner 3000/4000 (mounted)	Barcode scanner to read patient or sample ID and to scan single use reagents	400566
Barcode scanner 4000 holder (mounted)	The barcode scanner is mounted directly on the touch screen with the holder	208705
Barcode scanner 3000/4000 cable	Connection cable for the barcode scanner	400568
Clean Ring Rod	Tool to remove the Clean Ring around the axis	400596
Colour printer	Colour printer (USB)	400222
Connection cable eLine	Connection cable for electronic eLine pipette	400580
Connection cable eLine for left handed people	Connection cable for electronic eLine pipette for left handed people	200163
Connection cable for printer	Printer cable 2m USB	400223
Cup & Pin pro	Measuring cells for the ROTEM® *delta* (cup and pin) (10x20 pieces)	200011
Cup holder pro	Metal receptacle for the cup	400016
Disc for cup ejector pin	Noise absorbing disc around the cup ejector pin on the temperature controlled work area	400586
eLine pipette filter	Replacement filters (50 pieces)	100048
eLine pipette	Automatic, electronic eLine pipette	700030
eLine pipette tip tray	Pipette tips box (10x96 pieces)	400041
eLine pipette tip tray refill	Pipette tips refill (10x96 pieces)	400040
Keyboard	Keyboard with track ball for the ROTEM® *delta* system (specified language XX*)	400475-XX*
Keyboard protection foil	Protection foil for keyboard	400476
Kit MC Rod holder *delta*	Holder for the MC Rod	200175
MC Rod pro	Cup positioner	400045
Operating manual ROTEM® *delta*	Operating manual (specified language XX*)	200137-XX*
Parameter list	List with all parameters and their definitions	200086-EN

Table 10-7: Available accessories

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Article	Description	REF No. Tem Inno-vations
Pipette cable guide	Cable guide for eLine pipette cable	500220
Pipette and Clean Ring Rod holder	Replacement holder for pipette and Clean Ring Rod	200206
Plastic clip cup holder plus	Cup remover integrated in cup holder	400482
Reagent rack *delta*	Reagent rack for the ROTEM® *delta* system	400484
ROTEM® *delta*	Thromboelastometry equipment	200090
ROTEM® *delta* power cable	Power cable for the ROTEM® *delta*	700048-XX*
ROTEM® *delta* Software	Software for ROTEM® *delta*	200140
ROTEM® *delta* Secure Viewer CD	Software and manual for the ROTEM® *delta* Secure Viewer	200178
ROTEM® system USB Stick	ROTEM® system USB stick	400221
ROTEM® Connectivity Kit	LIS/HIS connectivity Kit	400509
ROTEM® Database Commander	View, print and export of the database on a personal computer	200172
Rubber band sample holder	Rubber band that keeps the sample in place	400646
Short user manual ROTEM® *delta*	Brief instruction (specified language XX*)	200124-XX*
Short user manual Single Use Reagents	Brief instruction (specified language XX*)	200176-XX*
Touch screen, glas	Touch screen for the ROTEM® *delta* system	408707
Touch screen, glas, connecting cable, power	Connecting cable between touch screen and ROTEM® *delta*, power	407107
Touch screen, glas, connecting cable, USB	Connecting cable between touch screen and ROTEM® *delta*, USB	200185
Touch screen, glas, connecting cable, VGA	Connecting cable between touch screen and ROTEM® *delta*, VGA	400584

* Country or language specific (e.g. -EN for English)

Table 10-7: Available accessories

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Appendix H: Documentation sheets

Here you find documentation sheets for maintenance and quality control. A detailed maintenance plan can be found in Chapter 8.2. In Chapter 5.5 it is described how to perform quality control.

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MAINTENANCE ROTEM® SN: _____

Daily Maintenance

Clean and disinfect cup holders after each use if necessary.
Check the lower side of the measurement block and clean and disinfect if required: For cleaning, a disinfectant cloth or a damp lint-free cloth with a disinfectant mentioned in this operator manual should preferably be used. Do not use water to avoid corrosion.
Check the pipette filter for stains and replace if required

Weekly Maintenance

Clean the ROTEM® system.
Change the pipette filter. If the ROTEM® system is used infrequently only, the filter may be exchanged whenever necessary (contaminated by blood).

DATE										
ROTROL N run performed										
ROTROL P run performed										
Measuring block disinfected										
Cup Holders disinfected										
ROTEM® surface disinfected										
Pipette filter changed										
SIGNATURE										

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MAINTENANCE ROTEM® SN: _____

Quaterly Maintenance

DATE					DIVERS
Check of Temperature					
Check of Optical Signals					
Check of Axes					
SIGNATURE					

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Quality control sheet

ROTEM® Quality Control Sheet: ROTROL N

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ROTEM® Quality Control Sheet: ROTROL P

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Maximal lysis 3-6

Maximum clot firmness 3-5

Measurement
- Finish 5-14
- Prepare 5-4
- Start 5-10
- Stop 5-14

Measurement module 6-2

Measurement results
- Export 5-24
 - Conventions 6-9
 - Parameters 6-10
- Find 5-22
- Graphics 5-12
- Interpretation 5-14
- Overlay 5-13
- Print 5-14, 6-7
- Print out 5-24

Measuring module
- 4-TEM 6-4
- Measurement 6-3
- Preparation 6-2
- Screenshot 6-4

N

Normal values 6-22

O

Operating instructions 1-3

P

Parameter 6-24

Patient data
- Enter 5-8
 - Conventions 5-9
 - Simplifications 5-9

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Error 8-15
Handling 5-2
Individual parts 8-11
Reassemble 8-13

Placing out of operation 9-1

Preparation of reagents 5-5
hep-TEM® 5-5
Other reagents 5-6

Prepare measuring cells. 5-6

Print
Measurements 5-24
Measurment results 6-7

Print report 6-29

Printer Setup 6-20, 6-34

Q

QC-Test
QCexN 5-20
QCexP 5-20
QCinN 5-20
QCinP 5-20

Quality control 5-17

R

Reaction curve 3-3

Reagents
ap-TEM® 3-16
ex-TEM® 3-14, 3-17
fib-TEM® 3-15
hep-TEM® 3-13, 3-18
in-TEM® 3-12, 3-17
Product line 3-11

Reference ranges 3-8

Remote access 6-17

Remove measuring cells 5-15

Replace pin remover 8-7

Research parameters 3-7

ROTEM®
Differential analysis 3-9
Measuring principle 3-1
Reagents 3-11
Tests 3-11

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Sampling 5-3
Transport 5-3
Sample preheating station 3-10, 4-3
Sample stability 3-9
Sample temperature 3-9
Scope of supply 4-1
Screenshots 6-4
Search form 6-6
Select test 5-10, 6-21
Serial interface 6-8
Single Shot Reagent 3-11
Software licence 10-5
Sorting data sets 6-6
Standard overlay 5-12, 6-4
Status line 5-5
Stop channel 5-14
Stop server 6-19
Switch off 5-16
Switch on 5-4
Symbols 10-1
System controls
ROTROL N 5-17
ROTROL P 5-17

T

Target group 1-4
TEMogram 3-3
Define colours 6-12
Print 6-12
Temperature 6-15
Temperature controlled work area 3-10
Temperature indicators 5-4
Tests
APTEM 3-16
EXTEM 3-14
FIBTEM 3-15
HEPTEM 3-13, 3-18
INTEM 3-12, 3-17
Thrombin inhibitor 3-12, 3-13, 3-17, 3-18
Thromboelastometry 3-2
Touch screen calibration 6-35, 6-36
Touch screen preparation 7-1
Trackball 6-1
Transmission parameters 6-17
Transmission protocol 6-17
Transmit measurement results in real time 6-16
Transport 9-1
Troubleshooting 5-26
Type label 4-4

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