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INSTRUCTIONS FOR USE
Dallop® NM ULTRALIGHT

Urological mesh for surgical treatment of urinary incontinence

CHARACTERISTICS

Non-resorbable urological surgical mesh **Dallop® NM ULTRALIGHT**, made of knitted polypropylene monofilament, of medical grade lightweight. The mesh has a diamond-like shape with elongated cells. The mesh has a marker for identification. The ends of the mesh are equipped with needles. Depending on the conditions of the surgical procedure, the needles can be bent. The product is intended for implantation in the human body. The product is sterile. The product is intended for implantation in urethral slings or mesh slings. The product is intended for implantation in urethral slings or mesh slings. The product is intended for implantation in urethral slings or mesh slings.

INDICATIONS

Urological mesh **Dallop® NM ULTRALIGHT** is indicated for the surgical treatment of female stress urinary incontinence (SUI) caused by urethral hypermobility and/or intrinsic sphincter deficiency (ISD), including:

- Primary or recurrent SUI
- SUI associated with mild pelvic organ prolapse (POP)
- SUI in patients undergoing concomitant pelvic surgery
- SUI in patients with a history of failed conservative treatment

The mesh is used in the following surgical procedures:

- Retropubic tension-free vaginal tape (TVT)
- Transobturator tape (TOT) or transobturator tape-outside-in (TOT-O)
- Mini-sling procedures

CONTRAINDICATIONS

Urological mesh **Dallop® NM ULTRALIGHT** is contraindicated in patients with:

- Active urinary tract infection
- Pregnancy
- Known hypersensitivity to polypropylene
- Connective tissue disorders
- Severe immunosuppression
- Uncontrolled diabetes mellitus
- Bleeding disorders or anticoagulant therapy
- Pelvic malignancy
- Previous pelvic radiation therapy
- Inability to undergo surgical procedures
- Known allergy to any component of the mesh

WARNINGS AND PRECAUTIONS

- The mesh should be implanted only by surgeons trained in pelvic floor reconstructive surgery and familiar with the specific risks and complications associated with mesh implantation.
- Careful patient selection is crucial. Preoperative evaluation should include assessment of urinary symptoms, pelvic examination, and urodynamic testing when indicated.
- Meticulous surgical technique is essential to minimize the risk of complications such as bladder perforation, urethral injury, vaginal erosion, and mesh exposure.
- Avoid excessive tension when placing the mesh, as this can lead to voiding dysfunction, pain, or mesh erosion.
- Ensure proper positioning of the mesh under the urethra, avoiding placement too close to the bladder neck or urethral meatus.
- Postoperative monitoring for signs of infection, hematoma, or urinary retention is important.
- Patients should be advised to report any unusual symptoms, such as persistent pain, fever, or changes in urinary function.
- Long-term follow-up is recommended to assess the durability of the repair and detect late complications.
- The mesh is not intended for use in patients with active infections, including urinary tract infections.
- Use with caution in patients with a history of bleeding disorders or those receiving anticoagulant therapy.
- Avoid use in patients with severe immunosuppression or uncontrolled diabetes mellitus.
- Do not use in patients with a history of pelvic malignancy or previous pelvic radiation therapy.
- Ensure the patient is able to undergo surgical procedures and has no contraindications to surgery.
- Do not use in patients with known allergy to any component of the mesh.

PRECAUTIONS

- The mesh should be implanted only by surgeons trained in pelvic floor reconstructive surgery and familiar with the specific risks and complications associated with mesh implantation.
- Careful patient selection is crucial. Preoperative evaluation should include assessment of urinary symptoms, pelvic examination, and urodynamic testing when indicated.
- Meticulous surgical technique is essential to minimize the risk of complications such as bladder perforation, urethral injury, vaginal erosion, and mesh exposure.
- Avoid excessive tension when placing the mesh, as this can lead to voiding dysfunction, pain, or mesh erosion.
- Ensure proper positioning of the mesh under the urethra, avoiding placement too close to the bladder neck or urethral meatus.
- Postoperative monitoring for signs of infection, hematoma, or urinary retention is important.
- Patients should be advised to report any unusual symptoms, such as persistent pain, fever, or changes in urinary function.
- Long-term follow-up is recommended to assess the durability of the repair and detect late complications.
- The mesh is not intended for use in patients with active infections, including urinary tract infections.
- Use with caution in patients with a history of bleeding disorders or those receiving anticoagulant therapy.
- Avoid use in patients with severe immunosuppression or uncontrolled diabetes mellitus.
- Do not use in patients with a history of pelvic malignancy or previous pelvic radiation therapy.
- Ensure the patient is able to undergo surgical procedures and has no contraindications to surgery.
- Do not use in patients with known allergy to any component of the mesh.

USE

1. Before starting the procedure, ensure you have all the necessary instruments and equipment, including the urological mesh, needles, syringes, and other materials.
2. Prepare the patient for the procedure, including anesthesia and positioning.
3. Carefully open the mesh package, avoiding damage to the sterile barrier.
4. Using the needles, carefully guide the urological mesh through the tissues, following the recommended technique and anatomical landmarks.
5. Ensure the mesh is correctly positioned and tensioned before finishing the procedure.
6. Close the incisions carefully, using appropriate suture material and technique.
7. After the procedure, monitor the patient for complications, such as bleeding, infection, or pain.
8. Provide the patient with instructions for post-operative care and activity restrictions during recovery.
9. Contact a healthcare professional if the patient has any concerns or complications.

PRECAUTIONS

1. The urological mesh **Dallop® NM ULTRALIGHT** is intended for use by one patient only and should never be reused. Reuse can lead to infection, inflammation, tissue damage, and other complications.
2. To ensure the sterility of the product, the package must be intact and undamaged. Do not use the product if the package is damaged or if there are signs that it has been opened previously.
3. Before using the product, carefully inspect it for any damage or defects. Do not use the product if it appears damaged or defective.
4. The product should be used in accordance with the manufacturer's instructions and medical protocols. Do not use the product for purposes not intended by the manufacturer.
5. Medical personnel should be properly trained in the use of the product and aware of the potential risks and complications associated with its use.
6. The patient should be fully informed about the potential risks and complications associated with the use of the product before giving consent to the procedure.
7. The product should not be used in patients who have an allergy to the materials from which it is made.
8. The product should not be used in patients who have a urinary tract infection, until the infection has been cured.
9. The product should not be used in patients who have a bleeding disorder, until the disorder has been cured.
10. The product should not be used in patients who have a kidney disorder, until the disorder has been cured.
11. The product should not be used in patients who have a liver disorder, until the disorder has been cured.
12. The product should not be used in patients who have a heart disorder, until the disorder has been cured.
13. The product should not be used in patients who have a lung disorder, until the disorder has been cured.
14. The product should not be used in patients who have a gastrointestinal disorder, until the disorder has been cured.
15. The product should not be used in patients who have a urinary tract disorder, until the disorder has been cured.
16. The product should not be used in patients who have a genital disorder, until the disorder has been cured.
17. The product should not be used in patients who have a nervous system disorder, until the disorder has been cured.
18. The product should not be used in patients who have a musculoskeletal disorder, until the disorder has been cured.
19. The product should not be used in patients who have a skin disorder, until the disorder has been cured.
20. The product should not be used in patients who have an immune system disorder, until the disorder has been cured.
21. The product should not be used in patients who have an endocrine system disorder, until the disorder has been cured.
22. The product should not be used in patients who have a metabolic disorder, until the disorder has been cured.
23. The product should not be used in patients who have a digestive disorder, until the disorder has been cured.
24. The product should not be used in patients who have a respiratory disorder, until the disorder has been cured.
25. The product should not be used in patients who have a cardiovascular disorder, until the disorder has been cured.
26. The product should not be used in patients who have a lymphatic system disorder, until the disorder has been cured.

PROCESS STERILIZATION

Sterilization is carried out using ethylene oxide (EO), within a validated process.

The product remains sterile, in the specified shelf life, only in its original packaging, under appropriate storage conditions.

STORAGE CONDITIONS

Dallop® NM ULTRALIGHT medical products should be stored:

- in cleaned rooms, at humidity 25-65%.
- at temperature 15°C-30°C.
- in conditions excluding exposure to mechanical damage or chemical contamination.
- in the original packaging (carton and blister packaging). Any mechanical damage to the packaging can expose the product to sterility loss.

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