	Do not reuse		Expiry date		Keep away from rain
	Do not use if the package is damaged		Manufacturer		Fragile, please handle it with great care
CE	CE Marking		Consult instructions for use	EC REP	Authorized representative in the European community

Tianjin Huahong Technology Co., Ltd.
A01, Plant B, No.278, Hangkong Road, Tianjin Pilot Free, Trade Zone (Air Port Industrial Park), 300308 Tianjin, China
Tel: +86 022-58775333 +86 4008-587-753
Fax: +86 022-58775353
E-mail: huahongsales@163.com

EC REP **Shanghai International Holding Corp. GmbH (Europe)**
Eiffestrasse 80, 20537 Hamburg, Germany
Tel:49-40-2513175

File No.: HH-IFU-01 Ver.: 2020.06.01

Instructions for use

CE 0197

Sterile lancets (Twist-off lancets)

[Structure] The Twist-off Lancets consist of needle, main body and protective cap.

[Models]

Model	Total Length (L_1)	Main body Length(L_2)	Exposed Needle Length(L_3)	Main body Width(W)
I	32±1.5mm	22.9±0.3mm	3.2±0.3mm	6.4±0.3mm
II	32 ±1.5mm	21.9±0.3mm	3.2±0.3mm	5.0±0.3mm
III	34.8±0.3mm	27.2±0.3mm	3.2±0.3mm	6.4±0.3mm
IV	29.5±1.5mm	20.0±0.3mm	3.2±0.3mm	6.4±0.3mm

[Needle Specifications]

Gauge	16G	18G	19G	20G	21G	22G	23G	24G	25G	26G
Diameter	1.50	1.20	1.07	0.91	0.82	0.72	0.64	0.57	0.51	0.46
Gauge	27G	28G	29G	30G	31G	32G	33G	34G	36G	38G
Diameter	0.41	0.36	0.34	0.31	0.26	0.24	0.21	0.19	0.17	0.15

[Performance] The device is sterilized by gamma irradiation.
[Intended use] The twist-off lancet(use with lancing device) is a single use, sterile medical device, which is designed for use of micro blood sampling puncture to obtain capillary blood samples from the fingertip.
[Contraindications] None.
[Instructions]
1. Unscrew the lancing device cover from the body of the lancing device;
2. Insert a sterile lancet into the lancing device;
3. Hold the lancet firmly in the lancing device and twist the cap of the lancet until it loosens, then pull the cap off the lancet;
4. Carefully screw the cover back onto the lancing device;
5. Adjust the puncture depth by setting the lancing device;
6. Press the lancing device against the side of the finger to be lanced with the cover resting on the finger. Push the release button to prick your fingertip;
7. Apply drop of blood to test strip;
8. Recap lancet, discard in appropriate container;
9. Place the lancing device cover back on the lancing device.

⚠	The combination products must be compatible with each other. When using Huahong's sterile lancets with other manufacturers' lancing devices, the operator must be sure that these products are compatible in the field of application.
⚠	We recommend that you shall consult a health care professional before use.

[Storage] The products should be stored in a well-ventilated room and kept away from light and moisture and corrosive gases.
[Shelf life] 5 years.
[Cautions]
1. Read the instructions carefully before use;
2. Using the lancets within the expiry date;
3. Do not use if protective cap is detached or broken;
4. Do not store lancet in Lancing device, discard Lancet after use, do not re-use;
5. For your safety, never use a lancet that has been used by someone else. Lancets are not intended for use on multiple individuals.
6. In the event of any serious adverse event, please report to the manufacturer and the competent authority.

[Symbol descriptions]

LOT	Lot number	⚠	Caution		Keep away from sunshine
STERILE R	Radiation sterilization		Manufacturing date		Do not re-sterilize

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