ENGLISH

Terumo Cardiovascular System - CDI™ Products
CDI™ Shunt sensor with heparin treatments, Model CDI510H

THE CDI™ BLOOD PARAMETER MONITORING SYSTEM OPERATOR'S MANUAL CONTAINS COMPLETE INSTRUCTIONS FOR ASSEMBLY AND USE OF THE SENSORS WITH THE CDI BLOOD PARAMETER MONITORING SYSTEM. IT MUST BE READ IN ITS ENTIRETY BEFORE ATTEMPTING TO ASSEMBLE AND USE THE SYSTEM.

Indications
The disposable CDI™ Shunt Sensor with heparin treatment is intended for use with the CDI™ Monitor during cardiopulmonary bypass procedures when continuous blood gas, pH, and potassium monitoring is desired.

Contraindications
The CDI510H is not intended for use in situations in which there is no flow of blood through the extracorporeal circuit. A minimum blood flow of 35 ml/min is recommended for optimal measurement perfor- mance of the shunt sensor. Restoration of adequate blood flow through the pump line after an interruption will restore optimal performance of the system.

The CDI Shunt Sensor with heparin treatment, Model CDI510H (Figure I)

The Shunt Sensor contains porcine derived heparin.

The shunt sensors are sterile, heparin-coated, non-toxic, non-pyrogen- ic, and for single use only. This device is sterilized by γ-irradiation.

The heparin treatment provides a non-leaching heparinized surface. It is applied by a process which covalently binds heparin to the blood contacting surfaces of this device. The heparin treatment improves the thromboresistance of the non-biological materials. Surfaces with heparin treatment have demonstrated stability in the presence of blood.

WARNING: FOR SINGLE USE ONLY. DO NOT REUSE. DO NOT RESTERILIZE. DO NOT REPROCESS. REPROCESSING MAY COMPROMISE THE STERILITY, BIOCOMPATIBILITY AND FUNCTIONAL INTEGRITY OF THE DEVICE.

CAUTION: DO NOT USE THE SHUNT SENSOR AFTER THE DATE PRINTED ON THE PACKAGE LABEL. USE BEYOND THE DATE MAY RESULT IN INACCURATE PERFORMANCE.

CAUTION: DO NOT USE THE SHUNT SENSOR IF THE FOIL POUCH HAS BEEN DAMAGED OR DEFLATED. DAMAGE TO THE POUCH MAY LEAD TO THE LOSS OF PROTECTIVE GAS. DEFLA- TION OF THE POUCH WILL INDICATE THE LOSS OF STORAGE GAS. LOSS OF STORAGE GAS MAY RESULT IN INACCURATE PERFORMANCE.

CAUTION: CHECK THE FREEZE INDICATOR ON THE BOX/ CARTON BEFORE USING THE ENCLOSED SHUNT SENSORS. STAINED BACKGROUND PAPER INDICATES EXPOSURE TO FREEZING TEMPERATURES HAS OCCURRED AND SENSORS SHOULD NOT BE USED.

CAUTION: THIS PRODUCT CONTAINS GERMALL II IN THE CALI- BRATION FLUID. A POTENTIAL FORMALDEHYDE EXPOSURE MAY BE FORMALDEHYDE. EXPOSURE MAY CAUSE ADVERSE REACTIONS IN PATIENTS WITH FORMALDEHYDE SENSITIVITY.

WARNING: DO NOT FREEZE SHUNT SENSOR. STORE ONLY AT SHOWN TEMPERATURE RANGE. FREEZING OR INCORRECT STORAGE MAY RESULT IN INACCURATE PERFORMANCE.

WARNING: DO NOT REUSE SENSORS. USED SENSORS ARE CONTAMINATED AND CANNOT BE STERILIZED. RESTERILIZATION WILL DAMAGE THE MICROSENSORS AND MAY RESULT IN INACCURATE PERFORMANCE.

WARNING: DO NOT USE ON PATIENTS WITH HEPARIN SENSITIVITY. DEVICES WITH HEPARIN TREATED SURFACES MAY CAUSE AN ADVERSE REACTION IN PATIENTS WITH HEPARIN SENSITIVITY.

CAUTION: FAILURE TO FOLLOW THE INSTRUCTIONS CAN CAUSE THE MONITORING SYSTEM TO DISPLAY INACCURATE VALUES. THE ACCURACY OF THE RESULTS IS DEPENDENT UPON THE FOLLOWING:
- **READING AND UNDERSTANDING THE INSTRUCTIONS FOR USE**
- **PROPER SET-UP AND CALIBRATION OF THE SYSTEM**
- **USE OF ALL AVAILABLE SYSTEM FEATURES**
- **PERIODIC COMPARISON TO A LABORATORY REFERENCE SAMPLE**

WHEN A DISPLAYED VALUE IS SIGNIFICANTLY DIFFERENT FROM EXPECTATIONS BASED ON THE CLINICAL SITUATION, VERIFY ITS ACCURACY BY INDEPENDENT MEANS BEFORE INITIATING TREATMENT.

WARNING: WHEN ADMINISTERING INTRAVASCULAR DYES AND NOVEL PHARMACOLOGICAL AGENTS, OR WHEN DYSHEMO- GLOBINS OR ELEVATED BILIRUBIN LEVELS ARE PRESENT, IN- DEPENDENT EXTERNAL BLOOD GAS AND BLOOD CHEMISTRY ANALYSIS ARE REQUIRED FOR ACCURATE DETERMINATION OF ALL MEASURED PARAMETERS NEEDED TO GUIDE THERAPEU- TIC DECISIONS. ALL THEY CAN POTENTIALLY CAUSE INACCURA- CIES IN DISPLAYED VALUE. REFER TO THE OPERATOR'S MANUAL FOR THE CDI BLOOD PARAMETER MONITORING SYSTEM FOR FURTHER SPECIFICATION OF AGENTS KNOWN TO INTERFERE WITH THE SYSTEM, THE VALUES AFFECTED AND ADDITIONAL INSTRUCTIONS.

WARNING: THE SHUNT SENSORS ARE STERILE, HEPARIN- COATED, NON-TOXIC, NON-PYROGENIC DISPOSABLE DEVICES, AND FOR USE IN CARDIOPULMONARY BYPASS PROCEDURES FOR UP TO 6 HOURS.

WARNING: MAINTAIN ADEQUATE LEVELS OF ANTICOAGULATION DURING EXTRACORPOREAL CIRCULATION BY MONITORING ACTIVATED CLOTTING TIME (ACT) OR OTHER APPROPRIATE MEASUREMENT. USE OF A HEPARIN TREATED DEVICE DOES NOT SUBSTITUTE FOR ADEQUATE ANTICOAGULATION LEVELS.

Calibration (Figure II)

Connect the CDI™ Model 540 Calibrator cable to the monitor. Using the system map key, select the CALIBRATE mode on the CDI monitor. Calibration of the shunt sensor is achieved by performing the following steps:

1) Remove the shunt sensor from the foil pouch and perform a visual inspection for damage to the sensor. Note the 3 digit K+ calibration value printed on the label.
2) Using sterile technique, remove the small top luer cap (white) from the shunt sensor, discard calibration fluid just prior to use, and flush sensor lumen prior to one end of your shunt line circuit.
3) Using sterile technique, remove the sensor filter/sparger assembly (clear) from the end of the shunt sensor and attach the other end of the shunt line circuit.
4) Prime and deblubble purge/shunt line, inspecting the line and sensor for bubbles. Clamp shunt line primed line for leaks; do not use the sensor if any leaks are detected.

The CDI Blood Parameter Monitoring System is now ready for use.

After selecting the OPERATE mode, the monitor will begin reporting blood parameter values.

Shunt Sensor Replacement

Shunt Sensor Replacement During Cardiopulmonary Bypass:
Terumo does not recommend removing the shunt sensor during cardiopulmonary bypass. In the event that the CDI510H shunt sensor(s) or the CDI™ monitor fails, the case should be completed using routine laboratory blood gas analysis.

Disposal

Dispose of this device according to hospital procedure for contaminated materials.

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