WaveLight®

Alcon®

WaveLight® FS200

Patient Interface

Only use with WaveLight® FS200 laser system

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Tel.: +49 1805 / 52 62 62
Fax: +49 9131 / 61 86 2221

Manufactured for:

Alcon Laboratories, Inc.
6201 South Freeway,
Fort Worth, Texas
76134 USA

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
Page

English 5
German 6
French 7
Spanish 8
Portuguese 9
Italian 10
Greek 11
Turkish 12
Dutch 13
Danish 14
Swedish 15
Norwegian 16
Finnish 17
Russian 18
Ukraine 19
Hungarian 20
Czech 21
Polish 22
Chinese 23

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Indication for Use

The WaveLight® FS200 Patient Interface is indicated to be used only with the WaveLight® FS200 laser system, consistent with the cleared Indications for Use for this ophthalmic surgical laser.

Instructions for Use

The WaveLight® FS200 Patient Interface is a sterile disposable product intended for single use only.

Never recycle or reuse the applanation cone and the suction ring /vacuum tube, for this can damage the components and thus may lead to unforeseeable cutting results.

The mechanical characteristics are changed when using the tubing and the applanation cone several times, which can lead to abnormal pressure (insufficient leak tightness).

In addition, sterility is not assured when reusing the cone, the tubing and the suction ring which can lead to contamination of the patients eye.

Before use, please check the expiration date marked on the package.

Prior to using the WaveLight® FS200 Patient Interface, please make sure that the package has not been damaged.

For further detailed information, please refer to the WaveLight® FS200 User Manual.

The plasticizer DEHP (di(2-ethylhexyl) phthalate) may be released from this product. For sensitive populations such as male babies, women during pregnancy or lactation, adolescent males, etc., health concerns about DEHP risk should be taken into consideration when medical professionals evaluate factors involving clinical treatment before using this product.

Disposal


Material and waste possibly being contaminated with pathogenic agents or substances must be collected in tear proof, humidity resistant and tight containers directly at their site of occurrence and subsequently - without refilling or sorting them - transported in appropriate, safely closed containers to the central collection point under the surveillance of the hygienist in charge.

Package Content of the WaveLight® FS200 Disposable Set

20 Patient Interfaces, each containing one application cone, one suction ring with vacuum tubing system

Conditions of Return

Returning the WaveLight® FS200 Patient Interface is excluded.

How to use the WaveLight® FS200 Patient Interface

Register the WaveLight® FS200 Patient Interface by means of the dataMATRIX* code at the WaveLight® FS200 laser system.

Open the sterile packaging and place the WaveLight® FS200 Patient Interface onto a sterile base.

Lock in place the male connector of the vacuum tubing system at the vacuum unit connection.

Prior to using the WaveLight® FS200 Patient Interface, please make sure that there is no bend or damage visible at any point of the complete vacuum tube (see figure B on page 3). If the structure of the tubes is incorrect, do not use the WaveLight® FS200 Patient Interface and please send it back completely to the address shown on page 2 of this instruction leaflet.

For further detailed information, please refer to the WaveLight® FS200 User Manual.

LOT	Lot number		Single use only	STERILE EO	EO sterilized
	Use by YYYY-MM-DD		Consult insturctions for use		Do not use in case of damaged enclosure
SN	Serial No.				

*dataMATRIX is a trademark of Meridios, Ltd.

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© 2012 Alcon, Inc.