LIGASORB

Protein A column for immunoadsorption

de	Gebrauchsanweisung	REF F00001564
en	Package Insert	
cs	Příbalová informace	
da	Indlægsseddel	
es	Prospecto	
[illegible]	[illegible]oste	
fr	Notice d'utilisation	
hr bih	Upute za uporabu	
hu	Használati utasítás	
[illegible]	[illegible]	
nl	Bijsluiter	
pl	Ulotka dołączona do opakowania	
[illegible]	[illegible]	
sv	Bipacksedel	
[illegible]	[illegible]	

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en

Package insert – please read carefully

LIGASORB® – Protein A column for immunoadsorption

Note!
Before using the LIGASORB column, please read the ADAsorb (medicap clinic GmbH) instructions for use, the respective package inserts from the devices approved for the primary plasma separation in combination with ADAsorb, as well as the package inserts for the respective disposables, solutions and sets.
The information provided in the Ligasorb package insert does not replace the information provided in the aforementioned package inserts.
Physicians should perform a risk-benefit analysis for each individual case before treating the patient.

Indication:
The LIGASORB immunoadsorber is specified for the quantitative elimination of immunoglobulins from plasma in patients whose disease is associated with specific antibodies and/or immune complexes. LIGASORB is specified for intervention in the acute exacerbation of autoimmune diseases. Treatment with LIGASORB must be prescribed and monitored for each individual patient by a qualified physician.

1. Intended use of the LIGASORB immunoadsorber

LIGASORB is a sterile medical product. It contains an adsorbent made of recombinant protein A which is immobilised covalently to a porous Polymethacrylate matrix. Protein A binds human IgG molecules (particularly IgG1, IgG2, IgG4) and immune complexes. Other immune globulins (IgG3, IgM and IgA) are partially bound to these ligands.

The housing consists of polycarbonate and contains the adsorbent between 2 polyester sieves. The connections to the adsorber are integrated into the housing cover. The internal housing volume is about 162 ml.

The binding capacity of the adsorber is > 2.0 g per loading cycle, under saturation conditions. The system is designed to treat 1-1.5x plasma volume. One adsorber is used for each treatment which is regenerated multiple times during treatment. The adsorber is designed for single use only.

LIGASORB may only be used

- if packaging is undamaged,
- if both housing and adsorber are undamaged,
- if all connections are undamaged,
- if protective caps on the connectors are not loose or have fallen off,
- if the expiry date has not been exceeded.

Note:
LIGASORB must be stored after receipt at +2 °C to +8 °C. The adsorber may not be frozen.

This product must be taken from the styrofoam box before proceeding to all steps that follow, at least 30 minutes before preparation in order to bring the adsorber to room temperature before use. Also, the solutions must be brought to the treatment room before use so that they can equilibrate to room temperature. The adsorber's connectors may only be opened when used.

The adsorber is filled with a non-physiological solution at delivery, which contains ascorbic acid. The solution must be completely rinsed out in one rinsing step before treatment is commenced. This process is automatically controlled by ADAsorb.

2. Application

2.1 Anticoagulation

Suitable anticoagulation is strictly required for extracorporeal treatment with LIGASORB. The anticoagulation must be prescribed and monitored by the responsible physician depending upon coagulation status and the history of the patient. It is recommended to check blood coagulation at regular intervals during treatment. For patients with no risk of bleeding, a heparin bolus of 70 to 80 IU/kg IV as well as a maintenance dose of 15 to 20 IU/kg/hour is administered. The dose must always be adjusted to the status of the patient and his/her risk of bleeding.

The recommendations of the primary plasma separation device must be considered for anticoagulation. If the anticoagulation is done with heparin, a supplemental administration of ACD-A should be given in a ratio of at least 1:20 in plasma before loading the LIGASORB column. The ACD-A serves to prevent coagulation as well as to suppress possible activation of the complement system.

Attention!
Administration of ACD-A requires close monitoring of blood ionised calcium levels.

2.2 Use of LIGASORB with the extracorporeal system

LIGASORB is used in combination with ADAsorb (to control apheresis) with an approved primary plasma separation device as well as appropriate accessories and solutions. ADAsorb

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controls preparation, loading and regeneration of the adsorber. The adsorber is regenerated after loading the adsorber with plasma, i.e. the plasma is discharged to neutral buffer; afterwards, the adsorber is eluted with an elution solution, and then neutralised with neutral buffer before a new plasma loading cycle is started. After the last load cycle at the end of treatment, the plasma is simply discharged. Since this is a single-column system, the plasma is sent back to the patient over the ADAsorb SB02 tubing system during regeneration of the adsorber. When the treatment is finished, follow the ADAsorb instructions for plasma return to the patient.

All machines are equipped and prepared independently according to the instructions in the respective package inserts. They are connected only before being connected to the patient. When the treatment is completed, the relevant instructions of the primary plasma separation process are to be followed to return blood and plasma to the patient.

Use of the COM.TEC/P1R set

Attention!

During filling, the clamp for the P1R set bypass line must be closed, including the tubing of the plasma lines or plasma return lines which to the ingoing or outgoing ADAsorb lines, including the Heidelberg storage with NaCl,.

The P1R set bypass line clamp must be closed during treatment.

Use of ADAsorb/SB02/LIGASORB

Attention!

3000 IU of heparin must be administered into 1 litre NaCl solution bag to prepare LIGASORB

The adsorber must be shaken for 10 seconds before use in order to be as homogeneous as possible. Afterwards, the adsorber is to be put into a provided adsorber stand.


Connecting all disposables and machines

Attention!

Observe aseptic working conditions.

Before connecting the plasma inlet line to the primary plasma separation device or the Heidelberg extension to the drip chamber of the SB02 tubing system, one must ensure that the clamp is closed at its vent.

When connecting the machine, ensure that the particle filter is installed in the patient return, as described in the ADAsorb package insert.

One must ensure that the adsorber is correctly connected (see drawing).

3. Side effects

3.1 Non-specific side effects due to extracorporeal treatment

Every extracorporeal treatment procedure is connected to stresses specific to each individual patient which can appear as non-specific side effects. The physician must inform the patient of the possibility of individual performance limitations for a period of time after treatment.

Side effects can occur due to interactions between blood or plasma and the artificial surfaces of the extracorporeal circuit. These reactions can be accelerated and/or worsened by specific factors for the respective patient's disease.

Side effects which can generally occur in an extracorporeal circulating system include:

Immune system:

- Hypersensitivity and allergic reactions with the following symptoms: Headache, nausea, vomiting, chest and stomach pains, diarrhoea, lower back pain, hypertension, erythema, flushes, eyelid oedema, tachycardia, fever, chills, increased sweating, muscle spasm, ear noise, pruritus, malaise, coughing.
 - Anaphylactic reaction with obstruction of the respiratory tract, bronchospasm and shock.
- Activation of the bradykinin system (e.g. flushes, back and abdominal pain, hypotension, globus sensation)

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Blood / coagulation:

- Anaemia
- Thrombocytopenia
- Coagulation disruption, bleeding tendencies
- Temporary change in the concentration of electrolytes and other plasma components

A side effect of venous puncture can be the formation of a haematoma, thrombosis or inflammation of the area around the vessel. Due to the concomitant administration of anticoagulants, there can be a temporarily increased generalised haematoma formation as well as bleeding.

If anticoagulation is done with continuous administration of ACD-A, a temporary systemic reduction of ionised blood calcium is possible; this can lead to an increase in neuromuscular excitability (with, for example, paraesthesia, tonic muscle cramps, tetany, and paralysis symptoms). For patients who already have hypocalcaemia prior to treatment, calcium administration to normalise the calcium level before or continuously during treatment should beconsidered.

Depending upon the degree of the clinical signs of hypocalcaemia, action must be taken to counteract further derailment of the calcium balance. At first, blood flow must be reduced. If this is insufficient, the mixing ratio of ACD-A to blood may be lowered at the discretion of the physician. If needed, calcium can be substituted through the administration of, for example, calcium gluconate 10 % IV according to directions and with a physician's prescription.

Biocompatibility of the materials is a particular focus during the development of extracorporeal systems. However, an aphaeresis procedure can lead to a more or less pronounced influence on the cellular components as well as activation of clotting and complement system due to contact of the blood and plasma with a "foreign" surface. This activation can lead to an anaphylactic reaction.

3.2 Specific side effects

- Regular monitoring of fibrinogen levels is recommended due to the reduction of fibrinogen under treatment and the concomitant risk of impaired blood clotting. This must be especially considered in the case of patients with low starting values.
- The therapeutically proposed elimination of immunoglobulins may cause increased susceptibility to infection due to the transient weakening of the humoral immune system.
- Reduction of fibrinogen and elimination of immunoglobulins can result in hypoproteinaemia.
- Temporary leucocytosis can be experienced under treatment.
- Temporary phosphate increase can be experienced under treatment.

4. Drug interactions

There are no known interactions between LIGASORB and other drugs or treatments. Drugs may be taken out of the blood or diluted with this procedure; therefore, the time when this is used should be adjusted, to the extent that it is medically reasonable, to the apheresis treatment (e.g. administer drugs after apheresis treatment). The type of accompanying treatment, dosage and time of the drug administration must be determined by the treating physician.

It is basically possible that bradykinin synthesis can be activated in an apheresis procedure. Bradykinins are quickly decomposed by ACE (Angiotensin Converting Enzyme). Treatment with ACE inhibitors can lead to excessive bradykinin and anaphylactic reactions.

5. Contraindications

Use of LIGASORB is contraindicated for

- patients who are unable to withstand the stress of an extra-corporeal treatment procedure due to their age, their physical developments or their clinical constitution,
- patients who cannot have suitable anticoagulation treatment, such as due to known hypersensitivity to heparin or ACD-A,
- patients who suffer from severe cardiovascular disease, so that extracorporeal treatment is not possible,
- patients who suffer from an acute, systematic infection,
- patients with pronounced allergic diathesis, especially with known hypersensitivity to the materials used in the extra-corporeal circuit (e.g. adsorber material, ascorbic acid)
- patients with haemorrhagic diathesis in whom extracorporeal apheresis procedures and anticoagulation performed have a high bleeding hazard

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6. Precautions

Responsible staff must continuously monitor patient parameters during treatment, especially the fluid and calcium balance, anti-coagulation and blood and plasma flow rates. A loss of platelets or platelet aggregation as well as loss of erythrocytes can occur with the apheresis procedure.

For first use of LIGASORB in patients under ACE inhibitor treatment, one must watch out for a possible anaphylactic reaction. In patients who are treated with ACE inhibitors and show reactions, they must stop taking it in time before an apheresis treatment. One must therefore consider the half-life of the substances used. The longer the half-life, the earlier the taking of these drugs must be stopped.

Treatment with LIGASORB should, as for all extracorporeal procedures, be used with care for patients who cannot tolerate therapeutic apheresis and have already shown an over-sensitivity reaction in connection with therapeutic apheresis or other extracorporeal treatments.

One must consider the present permanent treatment in patients with oral anticoagulants (e.g. ASA or phenprocoumon) for anticoagulation.

A particularly careful benefit-risk evaluation is advised before the application of apheresis in the following conditions:

- Patients who tend to hypocalcaemia due to a concomitant disease.
- Patients with clinically-relevant hypotension for whom a further drop in blood pressure could be damaging.
- Patients with restricted kidney function and cardiac performance who can react sensitively to fluid and electrolyte shifts.
- Patients with decompensated liver failure who have reduced protein synthesis.
- Patients with assured or suspected intracranial disease for whom a fluid imbalance or change in pressure should be avoided.
- Patients with unexplained bleeding tendency.

Note: There is no experience with immunadsorption in pregnant and nursing mothers or children. Treatment should only be carried out under strict indications in these cases.

Furthermore, the following must be observed:

- Before treatment, the particle filter must be installed in the patient return line as indicated in the ADAsorb directions for use. In order to avoid endangering the patient, if there is suspicion of or a confirmed blockage of the particle filter during treatment, do not continue the treatment under any circumstance without a filter.
- The products are for one-time use. Further use could pose a health risk to patients and users. Cleaning materials, disinfectants and preservatives can damage the material. Safety and performance capabilities can no longer be assured, and the manufacturer assumes no liability.
- One must always work in clean, aseptic conditions.
- Blood and plasma can be contaminated with pathogens of communicable diseases. All components of the treatment system must be regarded as potentially infectious after contact with patient blood or plasma. All components must be disposed afterwards in accordance with respective regional disposal legal provisions.

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