Medtronic

Duran AnCore™

Annuloplasty Ring and Band
Анулопластични пръстен и лента
Anuloplastický prstenec a proužek
Annuloplastikring og -bånd
Anuloplastie-Ring und -Band
Δακτύλιος και ανοιχτός δακτύλιος δακτυλιοπλαστικής
Anillo y banda para anuloplastia
Annuloplastika rõngas ja võru
Anuloplastiarengas ja -nauha
Anneau et bande d'annuloplastie
Prsten i traka za anuloplastiku
Anuloplasztikai gyűrű és pánt
Anello ed emianello per anuloplastica
Anuloplastikos žiedas ir lankas
Anuloplastijas gredzens un aploce
Прстен и обрач за анулопластика
Annuloplastiekring en -band
Annuloplastikkring og -bånd
Pierścień zamknięty i pierścień otwarty do annuloplastyki
Anel e banda de anuloplastia
Anel e banda de anuloplastia
Inel și bandă de anuloplastie
Кольцо и полукольцо для аннулопластики
Anuloplastický prstenec a pásik
Obroč in trak za anuloplastiko
Prsten i traka za anuloplastiku
Anuloplastikring och band
Anüloplasti Halkası ve Bandı
Кільце і півкільце для анулопластики
620R, 620B

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Istruzioni per l'uso • Naudojimo instrukcijos • Lietošanas pamācība • Упатства за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instruções de utilização • Instrucţiuni de utilizare • инструкция по применению • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з експлуатації

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Annuloplasty Ring and Band

1. Device Description

The annuloplasty ring and band are available in the following models:

- CG Future™ annuloplasty ring, Model 638R
- CG Future™ annuloplasty band, Model 638B

The ring and band each consist of an MP35N wire with a silicone overmold and a polyester fabric cover.

Both the ring and band must be implanted in the mitral position. The ring and band are marked at 3 points by colored sutures; 2 markers correspond to the trigones of the mitral valve and 1 marker identifies the midpoint of the device. The device size is identified by the inside diameter of the ring or band at its widest point. The MP35N wire allows radiographic visualization of the devices. The ring also includes a silicone marker, impregnated with barium sulfate, for radiographic visualization.

2. Indications for use

The CG Future annuloplasty ring and the CG Future annuloplasty band are indicated for the reconstruction or remodeling of pathological mitral valves. Valvular insufficiency or stenosis may be corrected by appropriate repair and annular remodeling.

3. Contraindications

- Heavily calcified valves
- Valvular retraction with severely reduced mobility
- Active bacterial endocarditis

4. Warnings and Precautions

4.1. Warnings

- This device was designed for single use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Only surgeons who have received adequate training to determine whether incompetent, stenotic, or diseased heart valves are capable of being repaired or replaced should use this device.
- Only surgeons who have received appropriate training in valve repair, including ring and band implant and sizing techniques, should use this device.
- Do not use rings and bands that have been damaged or contaminated from patient contact.
- Correct annuloplasty ring and band sizing is an important element of a successful valve repair. Using a ring or band that is too small can result in dehiscence or fracture, or valve stenosis. Using a ring or band that is too large can result in fracture, or valve regurgitation.
- Ensure that sutures are not placed in atrial tissue. Such placement may result in impairment of the cardiac conduction system.
- Secure the ring or band into fibrous trigone tissue to maintain permanent attachment.
- Securely tie all suture knots. Loose knots and long suture tails may be a source for hemolysis, thrombosis, or thromboembolism.
- Do not cut the ring or band. Loose threads can be a source of hemolysis, thrombosis, or thromboembolism.
- Use intraoperative or postoperative echocardiography to evaluate the effectiveness of the valve repair. Minimizing regurgitation and preventing systolic anterior motion (SAM) are important elements of an effective repair.
- Surgeons who use annuloplasty rings and bands should be current on all anticoagulation regimens.
- Carefully monitor the patient's anticoagulant status when using postoperative anticoagulant therapy.
- Patients with intra-atrial thrombi or a giant left atrium may benefit from long-term anticoagulation therapy.
- Patients in atrial fibrillation may remain on anticoagulation therapy until sinus rhythm is established.
- Do not alter or deform the ring or band to conform to annular anatomy. Altering or deforming the ring or band could lead to fracture of the ring or band and possible mitral regurgitation.
- Do not squeeze the ring or band with sharp instruments. Squeezing the ring or band may damage the surface of the stiffener, and could result in fracture of the ring or band and possible mitral regurgitation.
- The ring and band are indicated for use in the mitral position only. Use in the tricuspid position may result in fracture or dehiscence of the ring or band, conduction system damage, or tricuspid regurgitation.

4.2. Precautions

- Do not use cutting edge needles. Cutting edge needles may damage the ring or band and potentially lead to dehiscence, fracture, or possible mitral regurgitation.

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5. Potential adverse events

While infrequent, certain complications have been reported when using annuloplasty rings and bands. Reported complications include:

- Uncorrected or recurrent regurgitation
- Stenosis
- Ring or band dehiscence
- Hemolysis (even with mild regurgitation)
- Low cardiac output
- Heart block
- Systolic anterior motion (SAM) and left ventricular outflow tract obstruction (LVOTO)
- Damage to coronary arteries
- Damage to coronary sinus
- Endocarditis
- Thrombosis
- Thromboembolism
- Anticoagulant-related hemorrhage
- Ring or band fracture
- Leaflet perforation

Consider the potential for these complications when selecting the most beneficial surgical procedure for each patient.

Use the sizing and implantation methods described in these Instructions for Use to avoid or minimize occurrence of these adverse events. Ensure that the physician has appropriate training and experience in valve repair.

6. Instructions for use

6.1. Sizing

Select the appropriate ring or band size to help restore proper function of the valve. Use the CG Future™ sizer set, Model 7638, and the Medtronic annuloplasty handle, Model 7615, to select the proper ring or band size for repair of the mitral valve. The reusable sizers and handle are provided nonsterile and must be thoroughly cleaned and sterilized prior to use. Refer to the appropriate instructions for use for further information.

Align and insert the handle into the snap cavity of the sizer (Figure 2). The thin section of the handle is malleable to align the sizers with the valve annulus.

Warning: Do not use other manufacturers' annuloplasty sizers or sizers from other Medtronic annuloplasty products to size the CG Future annuloplasty ring or band. Other sizers may not indicate the appropriate CG Future ring or band size.

To determine the proper ring or band size, both the distance between the annular trigones and the area of the anterior leaflet must be measured.

1. Lower the sizer onto the valve annulus and align the sizer notches with the annular trigones.
2. Gently extend the anterior (A) leaflet and cover its surface with the selected sizer.

Note: The sizer that has a notch spacing most nearly matching the intertrigonal distance and a surface area most nearly matching that of the anterior leaflet corresponds to the appropriate size of the ring or band (Figure 3).

6.2. Handling and preparation instructions

1. Open the box. Remove the product literature and Patient Registration Form.
2. Remove the double-aseptic transfer pouch containing the ring or band mounted to a holder.
3. Inspect the outer pouch. If the outer pouch is damaged, the exterior surface of the inner pouch may not be sterile.
4. Open the outer transfer pouch, and while still holding the bottom of the outer pouch, pass the inner pouch into the sterile field.
5. Inspect the inner pouch. The ring and band assemblies are sterile as long as the inner pouch has not been compromised. If the inner pouch is damaged, do not implant the ring or band.
6. Open the inner pouch only in the sterile field.

6.3. Device implantation

The ring and band assembly may be used with or without the Medtronic annuloplasty handle, Model 7615, provided separately.

1. Thoroughly clean and sterilize the handle prior to each use.
2. With the printed side of the holder facing up, align and engage the sterile handle into the snap fit cavity of the holder (Figure 4). The thin section of the handle is malleable to align the device with the annulus.

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3. Remove the serial number identification tag (Figure 5). Record the serial number in the patient's record. Verify that the serial number matches the serial number on the Patient Registration Form.

 Warning: The serial tag must be removed from the ring or band for proper function. Do not cut or tear the fabric of the ring or band during removal of the serial tag.

Mitral annulus suture placement

1. Place sutures in each trigone, approximately 4 mm in width.
2. Place additional interrupted sutures, approximately 4 mm in width, in the posterior portion of the mitral annulus following the line of cusp insertion (Figure 6). Place a total of 10 to 14 sutures for the ring or band.

Warning: Do not place sutures in the circumflex coronary artery.

Note: Pledgets may be used for trigone sutures and midposterior annulus sutures to reduce the possibility of dehiscence to the ring or band.

Band suture placement

1. Orient the band assembly on the mitral (M) annulus on the holder (Figure 7).
2. Pass sutures through the band, approximately 2 mm in width, entering at the bottom of the band and exiting at the periphery of the band (Figure 8).

Note: In the trigonal tissue, pass the suture through the eyelets of the band stiffener to ensure a secure attachment (Figure 9). Pass only 1 suture through the trigonal marker on the band.

Ring suture placement

1. Orient the ring assembly on the mitral (M) annulus on the holder (Figure 7).
2. Pass sutures through the ring, approximately 2 mm in width, entering at the bottom of the ring and exiting at the periphery of the ring (Figure 8).

Note: Place sutures completely around the device, including the anterior mitral annulus region. In the trigonal tissue, pass the suture through the eyelets of the ring stiffener to ensure a secure attachment (Figure 9). Pass only 1 suture through the trigonal marker on the ring.

Holder and handle removal

1. Push the ring or band down onto the valve annulus while pulling back on the sutures.
2. Use a scalpel to cut the retention sutures on the holder in the areas indicated (Figure 10).
3. Remove the holder and handle from the valve annulus (Figure 11). Dispose of the holder.

Warning: Do not cut the fabric of the ring or band while cutting the holder retention sutures.

Knots and testing

1. Tie all knots around the device securely (Figure 12).
2. Trim all excess sutures.
3. Test valvular competency after removing the holder.

Warning: The holder must be removed from the ring or band at the end of the procedure for proper function. Never leave the holder attached to the ring or band.

Warning: Suture knots must be securely tied. Loose knots and long suture tails could cause hemolysis, thrombosis, or thromboembolism.

7. Individualization of treatment

Consider postoperative anticoagulation therapy for at least 6 weeks to allow for healing and incorporation of the annuloplasty ring or band by host tissue, regardless of cardiac rhythm.

8. Patient counseling information

Patients with annuloplasty rings or bands are at increased risk for endocarditis due to bacteremia (for example, when undergoing dental procedures). Advise patients about prophylactic antibiotic therapy.

Encourage patients to carry the Implanted Device Identification Card, provided by Medtronic, with them at all times.

9. HOW SUPPLIED

9.1. Packaging

A single annuloplasty ring or band assembly (device and holder) is packaged in sterile, double-aseptic transfer pouches (Figure 1). The packaging system is designed to ease placement of the assembly into the sterile field. The assemblies are

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sterile if the pouches are undamaged and unopened. The outer surfaces of the outer pouch are nonsterile and must not be placed in the sterile field.

9.2. Available sizes

The CG Future annuloplasty ring, Model 638R, is available in the following sizes for the mitral position: 24 mm, 26 mm, 28 mm, 30 mm, 32 mm, 34 mm, 36 mm, and 38 mm.

The CG Future annuloplasty band, Model 638B, is available in the following sizes for the mitral position: 26 mm, 28 mm, 30 mm, 32 mm, 34 mm, 36 mm, and 38 mm.

9.3. Storage

Store the product in its original packaging, including the outer shelf carton, in a clean, cool, and dry area to protect the product and minimize the potential for contamination.

The sterility and nonpyrogenicity of the CG Future annuloplasty ring and band are validated to remain unaffected until the Use-by date identified on the shelf carton, provided the pouches are not opened or damaged.

10. Patient Information

10.1. Registration Information

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the European Union.

A patient registration form is included in each device package. After implantation, please complete all requested information. The serial number may be found on the package and on the identification tag attached to the device. Return the original form to the Medtronic address indicated on the form and provide the temporary identification card to the patient prior to discharge.

An Implanted Device Identification Card is provided to the patient. The card contains the name and telephone number of the patient's physician, as well as information that medical personnel would require in the event of an emergency.

11. Postoperative Information

11.1. MRI safety information 

Nonclinical testing has demonstrated the CG Future annuloplasty ring and band are MR Conditional. A patient with this device can be safely scanned immediately after implantation in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T or 3 T
- Maximum spatial field gradient of 2500 gauss/cm (25 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 4 W/Kg (First Level Operating Mode)

Under the scan conditions defined above, the CG Future ring and band is expected to produce a maximum temperature rise of less than 2°C after 15 minutes of continuous scanning.

In nonclinical testing, the image artifact caused by the device extends approximately 20 mm from the CG Future ring and band when imaged with a gradient echo pulse sequence and a 3 T MRI system.

12. Accessories

Use the CG Future sizer set, Model 7638, to determine the appropriate ring or band size.

If using a handle, only use the Medtronic annuloplasty handle, Model 7615, to interface with the holder and sizers.

Caution: Do not use other manufacturers' annuloplasty sizers, or sizers for another Medtronic annuloplasty product, to size the CG Future annuloplasty ring or band.

Caution: Do not use the accessories until they have been thoroughly cleaned and sterilized. Refer to the appropriate instructions for use for further instructions.

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ALTHOUGH THE CG FUTURE ANNULOPLASTY RING, MODEL 638R, AND CG FUTURE ANNULOPLASTY BAND, MODEL 638B, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED

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UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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ALTHOUGH THE CG FUTURE ANNULOPLASTY RING, MODEL 638R, AND CG FUTURE ANNULOPLASTY BAND, MODEL 638B, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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Table of Contents

1 DEVICE DESCRIPTION 11
2 INDICATIONS 11
3 CONTRAINDICATIONS 11
4 WARNINGS 11
5 PRECAUTIONS 12
6 POTENTIAL ADVERSE EVENTS 12
7 INDIVIDUALIZATION OF TREATMENT 12
8 PATIENT COUNSELING INFORMATION 12
8.1 Patient Dental Care 12
8.2 MRI Conditions 12

9 HOW SUPPLIED 13
9.1 Packaging 13
9.2 Storage 13

10 INSTRUCTIONS FOR USE 13
10.1 Sizing 13
10.1.1 Contour 3D Ring Sizer Set, Single-Use, Model 7690SD 14
10.1.2 Contour 3D Ring Sizer Set, Reusable, Model 7690S 14
10.2 Device Handling and Preparation Instructions 14
10.3 Device Implantation 14
10.4 Accessories 15
10.5 Sterilization 15

11 REGISTRATION INFORMATION 16
12 DISCLAIMER OF WARRANTIES 17

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1 DEVICE DESCRIPTION

The Contour 3D Annuloplasty Ring, Model 690R, (hereafter referred to as Contour 3D Ring) consists of a titanium core encapsulated within silicone and covered with polyester fabric. The ring must be implanted in the tricuspid position. The ring is marked at 3 points by colored sutures. The 3 markers correspond to the commissures of the tricuspid valve. The device size is identified by the inside diameter of the ring at its widest point. The titanium core enables radiographic visualization of the device.

2 INDICATIONS

The Contour 3D Ring is indicated for the reconstruction and/or remodeling of pathological tricuspid valves.

3 CONTRAINDICATIONS

The Contour 3D Ring is contraindicated for patients presenting with any of the following conditions:

- heavily calcified valves
- valvular retraction with severely reduced mobility
- active bacterial endocarditis

4 WARNINGS

- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Only surgeons who have received adequate training to determine whether incompetent, stenotic, or diseased heart valves are capable of being repaired or replaced should use this device.
- Only surgeons who have received appropriate training in valve repair, including ring implant and sizing techniques, should use this device.
- Correct annuloplasty ring sizing is an important element of a successful valve repair. Undersizing the ring can result in valve stenosis, ring dehiscence, and/or ring fracture. Oversizing the ring can result in valve regurgitation and/or ring fracture.
- Do not place sutures in the coronary sinus, right coronary artery, AV node or other conduction tissue, as this may result in impairment of the cardiac conduction system.
- Suture knots must be securely tied. Loose knots and long suture tails may be a source for hemolysis, thrombosis, and/or thromboembolism.
- Do not cut the ring. Resultant loose threads may be a source of hemolysis, thrombosis, and/or thromboembolism. In addition, resultant jagged edges of the stiffener may lead to tissue trauma.
- Do not alter or deform the ring to conform to annular anatomy as this could lead to ring fracture, possible valve regurgitation, or stenosis.
- Do not squeeze the ring with sharp instruments as this may damage the surface of the stiffener, which may result in ring fracture and possible tricuspid regurgitation.
- Intraoperative and/or postoperative echocardiography should be used to evaluate the effectiveness of the valve repair. Minimizing regurgitation is an important element of an effective repair.
- Surgeons who use annuloplasty rings should be current on all anticoagulation regimens.
- When postoperative anticoagulant therapy is used, the patient's anticoagulation status should be carefully monitored.
- Patients with intra-atrial thrombi or a giant left atrium may benefit from long-term anticoagulant therapy.

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- The surgeon may desire that patients in atrial fibrillation remain on anticoagulation therapy until sinus rhythm is established.
- The ring is indicated for use in the tricuspid position only. Use in the mitral position may result in ring fracture, conduction system damage, mitral regurgitation, and/or ring dehiscence.

5 PRECAUTIONS

Do not use cutting edge needles, as they may damage the annuloplasty device potentially leading to ring dehiscence, ring fracture, and possible tricuspid regurgitation.

6 POTENTIAL ADVERSE EVENTS

While infrequent, certain complications have been reported when using annuloplasty rings. These include the following:

- uncorrected or recurrent regurgitation
- stenosis
- ring dehiscence
- hemolysis (even with mild regurgitation)
- low cardiac output
- heart block
- damage to coronary arteries
- endocarditis
- thrombosis
- thromboembolism
- anticoagulant-related hemorrhage
- ring fracture
- leaflet perforation

The potential for these complications should be considered when selecting the most beneficial surgical procedure for each patient.

To avoid or minimize occurrence of these adverse events, the annuloplasty repair, including sizing and implantation, should be conducted in accordance with the methods prescribed in these Instructions for Use by surgeons with appropriate training and experience in valve repair.

7 INDIVIDUALIZATION OF TREATMENT

To allow for healing and incorporation of the annuloplasty ring by host tissue, regardless of cardiac rhythm, postoperative anticoagulation therapy should be considered for at least 6 weeks following surgery.

8 PATIENT COUNSELING INFORMATION

8.1 Patient Dental Care

Patients with annuloplasty rings who undergo dental or other potentially bacteremic procedures must be considered for prophylactic antibiotic therapy.

8.2 MRI Conditions

Nonclinical testing has demonstrated the Contour 3D Ring is Magnetic Resonance (MR) Conditional. It can be scanned safely according to the conditions presented in Table 1. Patients with the Contour 3D Ring may safely undergo MRI for Normal or First Level Controlled Mode of the MR System as defined in the IEC Standard 60601-2-33.

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MR image quality may be compromised if the area of interest is in the exact area or relatively close to the position of the Contour 3D Ring. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant. See Table 1 for maximum distortion distances for spin echo and gradient echo sequences. These values are based upon nonclinical testing at 3T according to ASTM F2119-01.

Table 1. *Nonclinical MR Test Conditions and Results*

	1.5 Tesla System	**3.0 Tesla System**
Spatial Gradient Field (Units: Gauss/cm)	2200	2200
Maximum whole-body-averaged specific absorption rate (SAR) (Units: W/kg)	4	4
Scanning Time (Units: minutes)	15	15
Temperature Rise[a] (Units: °C)	< 2	< 2
Spin Echo Sequence Maximum Distortion Distance (Units: mm)	≤ 10	≤ 10
Gradient Echo Sequence Maximum Distortion Distance (Units: mm)	≤ 15	≤ 15

[a] In nonclinical testing according to ASTM F2182-02a when the local SAR is normalized to 1 W/kg for 15 minutes of MR scanning in a GE Signa®[1] 1.5T whole body transmit coil, and a GE Signa® HDxt 3T whole body transmit coil.

9 HOW SUPPLIED

9.1 Packaging

The Contour 3D Ring is available in the following sizes for the tricuspid position: 26, 28, 30, 32, 34, and 36 mm. The package contains a single annuloplasty ring assembly consisting of the device and holder (Figure 1) packaged in sterile, double-aseptic transfer trays. The packaging system is designed to ease placement of the device into the sterile field. The ring assembly is sterile if the trays are undamaged and unopened. The outer surfaces of the outer tray are NONSTERILE and must not be placed in the sterile field.

9.2 Storage

Store the product in its original packaging, including the outer shelf carton, in a clean, cool, and dry area to protect the product and minimize the potential for contamination. The sterility and nonpyrogenicity of the Contour 3D Ring are validated to remain unaffected until the Use By date identified on the shelf carton, provided the trays are not opened or damaged.

10 INSTRUCTIONS FOR USE

10.1 Sizing

Proper ring size selection is an important part of valvular annuloplasty to help restore proper valve function.

Warning: Do not use other manufacturers' annuloplasty sizers or sizers from other Medtronic annuloplasty products to size the Contour 3D Ring. Other annuloplasty sizers may not indicate the appropriate Contour 3D Ring size.

Two Contour 3D Ring Sizer Sets are available, and either set may be used to select the proper Contour 3D Ring size for repair of the tricuspid valve. One set is provided sterile and is for single use only. The other set is reusable and must be cleaned and sterilized prior to use.

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[1] Signa® is a registered trademark of General Electric Company.

10.1.1 Contour 3D Ring Sizer Set, Single-Use, Model 7690SD

The Contour 3D Ring Sizer Set, Model 7690SD, is for single use only. The sizers are provided STERILE and must not be resterilized. Single-use sizers that have been damaged or contaminated should not be used. Single-use sizers must be properly disposed of after use.

Use a sterile grasping surgical instrument to align the sizer to the valve annulus.

Warning: This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

10.1.2 Contour 3D Ring Sizer Set, Reusable, Model 7690S

The Contour 3D Ring Sizer Set, Model 7690S, is provided NONSTERILE. The sizers are reusable and are used with the Annuloplasty Handle, Model 7686. Both the reusable sizers and the handle must be thoroughly cleaned and sterilized by autoclave (steam) prior to each use.

Align and insert the handle into the snap cavity of the sizer. The thin section of the handle is malleable, allowing the surgeon to align the sizers to the valve annulus.

Warning: The reusable sizers and handle are provided NONSTERILE and must be thoroughly cleaned and sterilized prior to use.

Tricuspid Valve

To determine the proper Contour 3D Ring size, both the distance between the commissures of the septal leaflet (Figure 2) and the area of the anterior leaflet must be measured. First, lower the sizer onto the valve annulus and align the sizer notches with the commissures of the septal leaflet (Figure 3). Second, gently extend the anterior leaflet and cover its surface with the selected sizer (Figure 4). The sizer that has a notch spacing most nearly matching the intercommissural distance of the septal leaflet and a surface area most nearly matching that of the anterior leaflet corresponds to the size of the ring that should be selected (Figure 5).

10.2 Device Handling and Preparation Instructions

- Open the box and remove the product literature and Patient Registration Form.
- Remove the double-aseptic transfer tray containing the ring mounted to the holder.
- Inspect the trays ensuring they have not been opened or damaged. The ring assembly is sterile as long as the inner tray has not been compromised. If the inner tray is damaged, do not implant the ring.
- If the outer tray is damaged, the exterior surface of the inner tray may not be sterile.
- Open the outer transfer tray, and while still holding the bottom of the outer tray, pass the inner tray into the sterile field.
- The inner tray should be opened only in the sterile field.

10.3 Device Implantation

The Contour 3D Ring assembly may be used with or without the Annuloplasty Handle, Model 7686 (provided separately).

Warning: The handle is provided NONSTERILE and must be thoroughly cleaned and sterilized prior to use.

To use the handle, align and engage the sterile handle into the snap fit cavity of the holder (Figure 6). The thin section of the handle is malleable, allowing the surgeon to align the device with the annulus.

Note: The handle must be inserted into the side of the holder with the snap fit cavity.

Remove the serial number identification tag and record the serial number in the patient's record. Verify that the serial number matches the serial number on the Patient Registration Form.

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Warning: The serial tag must be removed from the ring for proper function. Do not cut or tear the ring fabric during removal of the serial tag.

Tricuspid Suture Placement

Caution: The Contour 3D Ring is designed for implantation only with interrupted suture techniques.

The first annular suture should be placed at the midpoint of the septal leaflet. Place additional interrupted sutures, approximately 4 mm in width, continuing along the tricuspid annulus ending at the anteroseptal commissure (Figure 7).

Care should be taken to avoid placing a suture through area of the bundle of His.

Warning: Do not place sutures in the coronary sinus, right coronary artery, AV node, bundle of His, or other conduction tissue.

Note: Pledgets may be used with annular sutures to reduce the possibility of ring dehiscence.

Ring Suture Placement

Orient the ring assembly to the tricuspid annulus.

Pass sutures through the ring, approximately 2 to 4 mm in width, entering at the bottom of the ring and exiting the periphery of the ring (Figure 8).

Holder/Handle Removal

Use the holder/handle to push the ring down onto the valve annulus while pulling back on the sutures (Figure 9 and Figure 10). Using a scalpel, cut the retention sutures on the holder in the areas indicated (Figure 11 and Figure 12). Remove the disposable ring holder and handle from the valve annulus (Figure 13). Dispose of the ring holder.

Warning: Do not cut the ring fabric while cutting the holder retention sutures.

The Contour 3D Ring holder must be removed from the ring at the end of the procedure for proper function. Under no circumstances should the holder remain attached to the ring.

Knots, Testing

Tie all knots around the ring securely (Figure 14). Trim all excess sutures and test valvular competency after removing the holder (Figure 15).

Warning: Suture knots must be securely tied. Loose knots and long suture tails may be a source of hemolysis, thrombosis, and/or thromboembolism.

10.4 Accessories

Refer to the Instructions for Use provided with each accessory for detailed information on its use and sterilization.

Warning: Do not use other manufacturers' annuloplasty sizers or sizers from other Medtronic annuloplasty products to size the Contour 3D Ring. Other annuloplasty sizers may not indicate the appropriate Contour 3D Ring size.

Contour 3D Ring Sizer Set, Single-Use, Model 7690SD

The Contour 3D Ring Sizer Set, Model 7690SD, is for single use only. The sizers are provided STERILE and must not be resterilized. Single-use sizers that have been damaged or contaminated should not be used. Single-use sizers must be properly disposed of after use.

Contour 3D Ring Sizer Set, Reusable, Model 7690S

The Contour 3D Ring Sizer Set, Model 7690S, is provided NONSTERILE. The sizers are reusable, however, they must be thoroughly cleaned and sterilized by autoclave (steam) prior to each use.

Annuloplasty Handle, Model 7686

Use only the Annuloplasty Handle, Model 7686 to interface with the ring holder and Contour 3D sizers. The handle is provided NONSTERILE and must be thoroughly cleaned and sterilized prior to use.

10.5 Sterilization

The ring is provided sterile (steam) on the holder and must not be resterilized. Rings that have been contaminated from patient contact should not be used.

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11 REGISTRATION INFORMATION

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the EU.

A Patient Registration Form is included in each device package. After implantation, please complete all requested information. The serial number may be found on the package and on the identification tag attached to the device. Return the original form to the Medtronic address indicated on the form and provide the temporary identification card to the patient prior to discharge.

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12 DISCLAIMER OF WARRANTIES

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO UNITED STATES CUSTOMERS ONLY:

DISCLAIMER OF WARRANTY

ALTHOUGH THE CONTOUR 3D ANNULOPLASTY RING, MODEL 690R, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portions of this DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO CUSTOMERS OUTSIDE THE UNITED STATES:

DISCLAIMER OF WARRANTY

ALTHOUGH THE CONTOUR 3D ANNULOPLASTY RING, MODEL 690R, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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Trademarks may be registered and are the property of their respective owners.
All Rights Reserved
Contour 3D™ is a trademark of Medtronic, Inc.
Medtronic B.V.
Medtronic Australasia Pty Ltd
Medtronic of Canada Ltd