STORZ

KARL STORZ — ENDOSKOPE

Fiberskope
Serien 110xx ..., 111xx ..., 112xx ..., 113xx..., 13304C*)

Fiberscopes
Series 110xx ..., 111xx ..., 112xx ..., 113xx..., 13304C*)

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Content

14. Cleaning, disinfection, sterilization.....31
14.1 Reprocessing Overview 1, German.....34
14.2 Reprocessing Overview 1, English.....35
14.3 Reprocessing Overview 1, Spanish.....36
14.4 Reprocessing Overview 2, German.....37
14.5 Reprocessing Overview 2, English.....38
14.6 Reprocessing Overview 2, Spanish.....39
14.7 Reprocessing Overview 3, German.....40
14.8 Reprocessing Overview 3, English.....41
14.9 Reprocessing Overview 3, Spanish.....42
15. Preliminary cleaning.....43
15.1 Disassembling suction valve insert.....45
16. Manual reprocessing.....46
16.1 Manual cleaning.....46
16.2 Manual disinfection.....48
17. Machine reprocessing.....50
18. Inspection after reprocessing.....58
18.1 Assembling suction valve insert.....59
19. Gas sterilization/ plasma sterilization.....60
20. Chemical sterilization.....63
21. Storage.....64
22. Accessories.....65
23. Overview: Cleaning brushes.....68
24. Normatives.....70
24.1 Servicing and repair.....70
24.2 Limitation of liability.....70
24.3 Warranty.....70
24.4 Directive compliance.....70
25. Repair and replacement program.....71
25.1 Important information.....71
Subsidiaries.....72

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It has been determined that Broncho-Fiberscopes and Intubation fiberscopes (see table next page) incorporate a lumen material that is only approved for use as an external surface material and not included in the STERRAD® sterilization systems FDA-cleared claims for use in flexible endoscope lumens. This could result in lower than expected levels of hydrogen peroxide within the lumen during the exposure phase of the sterilization cycle and compromise sterilization efficacy.

Based on the risk factors discussed above, the probability of a patient infection occurring from a cleaned and disinfected bronchoscope or intubation scope processed in the STERRAD® sterilization systems is defined as "remote". Since bronchoscopes and intubation scopes are semi-critical (B) devices and need not necessarily sterilization, KARL STORZ has validated many sterilization methods.

Available validated sterilization methods include:

- STERRAD® 100S with Booster,
- STERIS® System 1,
- STERIS® System 1E,
- STERIS® AMSCO® V-PRO® maX and
- ethylene oxide (EO) gas.

The instructions „Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments", Version 1.0.0, No. 96216003DER, describe in detail the workflows and process-relevant parameters of the individual validated procedures. The procedure must be chosen taking into account applicable national requirements and in consultation with the device manufacturers.

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Important Information

The methodologies for manual and automated cleaning and disinfection are already described in details in this instruction manual. When selecting the disinfectant the range of action should be taken into account with regard to microorganisms. Use only chemicals approved by KARL STORZ for processing. For a complete list, please visit www.karlstorz.com. When preparing and using solution, strictly observe the chemical manufacturer's instructions on concentration, temperature, exposure time, and shelf life.

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Thank you for your expession of confidence in the KARL STORZ brand name. Like all of our other products, these products are the result of years of experience and great care in manufacture. You and your organization have decided in favor of modern, high-quality items of equipment from KARL STORZ.

This manual is intended to serve as an aid in the proper use, cleaning , desinfection and sterilization of fiberscopes manufactured by KARL STORZ. The fiberscopes are sensitive precision instruments which require particular care during reprocessing. To prevent damage during reprocessing, support is provided here for disassembly, cleaning, sterilization and assembly of the KARL STORZ fiberscopes. For this reason, read these instructions carefully and keep them as a ready reference.

This instruction manual describes three fiberscope versions, distinguished from one another by the imaging system. If a section refers exclusively to one of these versions, it is thus appropriately labeled.

PDD OPT	Telescope	
PDD	Telescope with special filters for PDD/AF	
DCI®	DCI® fiberscopes (Direct Coupled Interface)	

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope for supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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2. Controls

Connectors, accessories and their uses

495 F	Light cable adaptor for Wolf
495 G	Light cable adaptor for KARL STORZ, Olympus, Winter & Ibe
11025 E	Pressure compensation cap (red)
11301 CA	Lip valve, disposable version, pack of 20
13242 XA	Adaptor for leakage test in cleaning and disinfection machines (green)
29100	Plug for instrument channel, pack of 10, reusable
11301 CB	Suction valve insert, reusable version
5917800	Seal ring (suction valve insert)
6885391	Cap (suction valve insert)
5923110	Spring (suction valve insert)
6927791	Housing (suction valve insert)
8567290	Plunger (suction valve insert)
11301 CE	Suction valve insert, disposable version, pack of 20
11301 D1	LED battery-powered light source with fine screw thread
11301 D3	LED battery-powered source with quick-release thread
6927691	Three-way double stopcock adaptor with LUER-Lock
6011590	LUER cap
600007	Barbed tube, thin, with LUER-Lock

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3. Intended Use

Fiberscopes are used in diagnostic and surgical endoscopy. They must only be used by appropriately trained and qualified technical staff.

Indications: Examination, diagnosis and/or therapy in conjunction with endoscopic accessories in various medical fields.

Contraindications: No contraindications are known which directly relate to the product. The fiberscope must not be used if, in the opinion of a qualified physician, such an application would endanger the patient, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

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3.1 Safety precautions when using the fiberscopes

Fiberscopes must be used according to the medical rules and procedures of endoscopy recognized for this method.

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4. User qualification

KARL STORZ fiberscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the endoscopic technique.

The information given in these instructions only serves to instruct in the correct handling, cleaning, disinfection and sterilization of the fiberscopes. It is not intended as an introduction to the endoscopic technique.

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5. Warnings and cautions

Please read this manual and follow its instructions carefully. The words, **Warning, Caution** and **Note** convey special meanings. Wherever they are used in this manual, the text should be carefully reviewed to ensure the safe and effective operation of this product. To make the words **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

WARNING: *A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.*

CAUTION: *A Caution indicates that particular service procedures or safety precautions must be followed to avoid possible damage to the product.*

NOTE: *A Note indicates special information about operating the product, or clarifies important information.*

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WARNING: *Before using the fiberscope, read the following safety instructions carefully, to avoid putting your patients, personnel or yourself at risk.*

WARNING: *Inspect the fiberscope prior to each surgical procedure to ensure that it is functioning correctly and that it has been correctly cleaned, disinfected and (where applicable) sterilized.*

WARNING: *This equipment may not be used during the discharge of a cardiac defibrillator. Remove all endoscopic equipment from patient contact and from the surgical site.*

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Warnings and cautions

WARNING: *Verify the sound condition of the fiberscope as well as of accessories used in combination before every application. Damaged fiberscopes or damaged accessories must not be used.*

WARNING: *Combinations of medical devices are only then assured to be safe if*
- *they are identified as such in the respective instruction manual or*
- *the intended purpose and the interface specifications of the devices used in combination permit this.*

Pay careful attention to the instructions and interface specifications of medical products used in combination.

WARNING: *Lengthy direct contact between the endoscope tip and surrounding tissue may cause tissue damage due to the potentially high temperatures on the tip of the endoscope and must thus be avoided.*

WARNING: *The fiberscope is delivered unsterilized, and must therefore be cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.*

WARNING: *Remove the red pressure compensation cap from the vent port*
- *for every endoscopic intervention*
- *before immersion in liquid/chemical sterilization.*

WARNING: *Before the fiberscope is removed from the patient, the locking mechanism ⑮ must be released and the distal tip of the fiberscope put in the neutral, non-angled position. Otherwise the instrument may be damaged or the patient injured.*

WARNING: *If a malfunction should occur during application on a patient, stop the application immediately. Place the distal tip of the fiberscope in the neutral, non-angled position and remove the fiberscope slowly and carefully from the patient.*

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Warnings and cautions

CAUTION: *Do not steam sterilize (autoclave) the fiberscopes.*

CAUTION: *The maximum permissible temperature for cleaning, disinfection, sterilization, transport and storage is 65 °C.*

CAUTION: *Never immerse fiberscopes in normal saline solution, since even only brief contact can lead to pitting and corrosion.*

CAUTION: *Material damage: if the outer casing of the fiberscope is cleaned with alcohol, this will damage the material.*

CAUTION: *Fiberscopes must not be cleaned in an ultrasound bath under any circunstances.*

CAUTION: *Place the red pressure compensation cap on the vent port before*
- *gas sterilization/plasma sterilization*
- *shipping*

CAUTION: *Only use KARL STORZ accessories (guide probes, cleaning brushes etc.). Accessories from other manufacturers may damage the fiberscopes.*

CAUTION: *Only use guide probes with flexible ends. Other types of guide probe may cause considerable damage to fiberscopes.*

CAUTION: *Make certain to read the ‘Cleaning, sterilization and care of KARL STORZ instruments’ instructions. Methods of cleaning, disinfection and sterilization are explained in detail there.*

NOTE: *Any damage to the fiberscope resulting from incorrect operation is not covered by the guarantee.*

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NOTE: *If the instrument is in the working channel, its additional rigidity may limit angling.*

NOTE: *Separate instruction manuals are available for flexible videoscopes; they are therefore not described in this brochure.*

NOTE: *Federal law restricts this device to sale by or on the order of a physician (USA only).*

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6. Unpacking

Check for missing items and evidence of shipping damage. File any complaints to KARL STORZ or the supplier immediately.

6.1 Basic equipment

1	Fiberscope	Art. no.
1	Suction valve insert, disposable version, pack of 20 inserts	11301 CE
1	Leakage tester	13242 XL
1	Pressure compensation cap	11025 E
1	Plug for instrument channel, pack of 10	29100
2	Adaptor for cleaning the working channel in cleaning machines (green)	11301 CD
1	Cleaning brush for 2.8 x 50 and for 2.8 x 70	27651 A*)
	for 3.7 x 65	27651 AL*)
	for 5.0 x 65	27651 B
	for 5.6 x 65	27651 C
	for 6.4 x 65	27651 C
1	Biopsy forceps (only bronchofiberscopes)	11001 KL
1	Grasping forceps (only bronchofiberscopes)	11002 KS
1	Three-way double stopcock adaptor with LUER-Lock for O_2 connection (only intubation fiberscopes)	6927691
1	LIPP tube holder (only intubation fiberscopes)	11301 CF
1	Mouth piece (only gastro-, intubation and bronchofiberscopes)	13272

*) For the fiberscopes 11278 A/AI/AU/AIU use the cleaning brush 27651 A (27651 AL is not suitable for these fiberscopes).

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6.1.1 Accessories

1 Suction valve insert, reusable version 11301 CB
1 Battery light source 11301 D
1 Adaptor for leakage test in cleaning and disinfection machines (green) 13242 XA
4 Cleaning adaptors for cleaning macines 11301 GG2/FF2/HH/KK2

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7.1 Fitting the suction valve

Place the suction valve insert ② on the valve nozzle ③ and turn until the projections on the side of the insert lie over the nozzle recesses.

Push gently until the suction valve insert clicks into the valve nozzle.

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7.2 Fitting the tube holder (only intubation fiberscopes)

Push the tube holder onto the fiberscope from distal and secure it by turning the two parts in opposite directions.

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7.3 Fitting the lip valve (only bronchofiberscopes)

Place the disposable lip valve onto the instrument channel ⑧ and tighten gently by turning it in a clockwise direction (plastic part).

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8. Handling the fiberscopes

The optical section of KARL STORZ fiberscopes consists of a fiberoptic bundle in a flexible and watertight casing, a light guide, as well as objective and eyepiece lenses. Like all optical systems, these parts are very delicate. For this reason, careful handling of the fiberscopes is essential for a long service life and a consistent quality. Therefore be aware of the following points:

- Never kink or tightly coil the fiberscopes and light guides, so that no optical fibers break off. It is therefore best to suspend the fiberscopes for storage or keep in large containers.
- Do not jerk the up/down control ④ or move forcefully.
- Do not insert the instruments or cleaning brushes forcefully through the channels. This may cause damage to the casing, especially if the tip of the instrument is deflected.
- Carefully protect the objective lens ⑨ at the distal end from damage. Do not hit. Never remove soiling with sharp-edged objects.
- Never store or disinfect the fiberscope together with sharp-edged objects to prevent damage to the exterior.
- Never place heavy objects on top of the fiberscope.

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NOTE: *Every fiberscope supplied by KARL STORZ includes a photo showing the condition of the fiber bundle upon delivery. A copy of this photo is kept on file in our quality control office. The wedge-shaped marking (at 12 o'clock) serves as a guide for the user.*

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9. Using the fiberscopes

The following instructions are recommended for the proper operation of the fiberscope and accessories. This is not a reference for endoscopic techniques.

Consult the appropriate medical literature for information on endoscopic techniques.

You can find accessories for your fiberscope (catheters, flexible scissors and forceps, guide wires, stone baskets, lithotriptor probes, laser fibers, electrodes etc.) in the appropriate KARL STORZ catalog.

WARNING: *Always adjust the light source to the minimum illumination intensity necessary to achieve optimum visualization of the operating field. The higher the intensity setting of the light source, the greater the heat energy that will be generated at the tip of the fiberscope.*

WARNING: *Never lay the end of the fiber optic cable or the distal tip of the fiberscope on or under a sterile cover cloth. The light intensity may lead to burns on the patient and/or sterile cover cloth.*

- When connecting and operating light sources and suction/irrigation equipment, be certain to read the appropriate instruction manuals for the equipment. The irrigation equipment recommended by KARL STORZ is set to a maximum pressure of 54 kpa (400 mmHg/0,54 bar/7.8 psi).
- Endoscopic examinations must be carried out under aseptic conditions.
- When introducing the fiberscope into the patient, the distal tip of the fiberscope must be in a straight, nondeflected position (control lock-in can be felt and heard).
 The fiberscope must be handled according to the recognized medical rules and procedures for endoscopic methods.

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Using the fiberscopes

- If no instrument is inserted, close the instrument channel proximally with a plug (Art. no. 29100).
- If an instrument is to be inserted, position the disposable lip valve (Art. no. 11301 CA) to seal off the instrument channel proximally.
- Instruments especially designed for this may be inserted through the instrument channel of the fiberscopes. Make certain the accessories have a suitable length and diameter. Unsuitable accessories may damage the fiberscopes. For questions on accessories, please contact KARL STORZ directly or an authorized dealer.

ⓘ ***NOTE:*** *When the lip valve is first used, the instrument is subjected to increased resistance. This resistance will decrease the more often the lip valve is used.*

- When inserting an instrument into the instrument channel, the fiberscope must be in a straight, nondeflected position.

❗ ***CAUTION:*** *When inserting instruments, increasing resistance is felt with increasing deflection. Forced introduction of instruments can lead to damage to the fiberscope or the instrument channel.*

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Using the fiberscopes

CAUTION: *If undue resistance is encountered when inserting an instrument, do not attempt to force the instrument. Withdraw the instrument and apply a water soluble lubricant to the tip of the instrument and attempt to insert again.*
Prior to removal of the instrument from the patient, the distal tip must be in a straight position.

- When using forceps, scissors or other accessories that can be opened or closed, ensure that the instrument is completely closed prior to insertion into the instrument channel. Do not attempt to open the jaws until the distal tip is completely free of the channel. Opening the jaws inside the instrument channel can lead to damage to the fiberscope.

WARNING: *When using HF electrodes or laser probes, make sure that the active parts are clear of the instrument channel and in full view of the operating surgeon before power is applied. Activation within the instrument channel may damage the fiberscope and/or injure the patient.*

WARNING: *Make certain to read the instructions for the laser unit and/or HF generator.*

CAUTION: *If accessories with sharp tips, such as EHL probes or laser probes, are inserted into the instrument channel of a deflected fiberscope, this can damage the fiberscope.*

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Using the fiberscopes

CAUTION: *When using a stone basket, the stone must not be retracted through the instrument channel but instead the entire fiberscope must be removed. Once the stone is removed, the stone basket may then be taken out of the instrument channel.*

CAUTION: *The fiberscopes must be cleaned and then reprocessed immediately after the intervention in order to prevent incrustations. Dried-on incrustations may jeopardize successful decontamination and disinfection.*

- A leakage test must be carried out after every intervention, to prevent damage during reprocessing.

CAUTION: *Use only KARL STORZ accesories (guide probes, cleaning brushes etc.). Accessories from other manufacturers may damage the fiberscopes.*

NOTE: *If the instrument is in the working channel, its additional rigidity may limit maneuverability.*

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10. Preparation of the fiberscope

Before using the fiberscope on a patient, exterior damage and operational problems must be eliminated.

Check the outer sheath of the fiberscope. Be especially aware of visible damage to the casing, particularly in the area of the fiberscope tip.

Look through the eyepiece and make certain that you have a clear image. The grid (Moiré effect) should be easily recognizable. If the image is not clear but cloudy, clean the eyepiece and lens with a cotton applicator soaked in alcohol, maintaining hygienic conditions at all times.

DCI® fiberscope image quality can be controlled on screen. If necessary, clean image waveguides, light cables and eyepieces as described above.

If it is still cloudy, it is possible that the objective lens is scratched or that there is moisture inside the telescope. The fiberscope must not be used in this case, but instead sent in for repairs.

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CAUTION: *Before every application on the patient, check the instrument for leaks.*

Check the fiberscope for leaks before every application on a patient.

Make certain to read the specifications on page 43 of this instruction manual.

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10.1 Connect these supply lines

- If necessary, remove the adaptor 495 G or 495 F from the light inlet piece and connect the fiberoptic light cable or battery-powered light source to the light inlet piece ⑦.

NOTE: *Not necessary for battery-powered light sources with coarse threads (11301 D3).*

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Preparation of the fiberscope

- suction to suction outlet ⑥
- or O_2 supply (under anesthesia) to instrument channel ⑧ via the appropriate adaptor

CAUTION: *The instruction manuals of the supply units used must be observed.*

Switch on the respective supply units.

Insert the DCI® fiberscope into the DCI® camera head until it locks into place.

Insert the light cable (Art. no. 495 DV) into the DCI® camera head as far as it will go.

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Connect the light cable and the video signal cable up to the appropriate supply units.

10.2 Controlling the distal tip

The up/down movement of the distal tip on the fiberscope can be adjusted at the control ④: move smoothly and slowly.

Some fiberscopes have:
- contrapositive deflection direction (urology);
- a locking mechanism ⑮ for securing the fiberscope tip (gynecology).

Models with locking mechanism ⑮:

To secure the deflection postition of the distal tip, turn the locking mechanism ⑮ slowly to distal until resistance can be felt.

CAUTION: *Forcing the locking mechanism ⑮ further in the same direction may damage the locking mechanism ⑮.*

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Preparation of the fiberscope

WARNING: *Before the fiberscope is removed from the patient, the locking mechanism ⑮ must be released and the distal tip of the fiberscope put in the neutral, non-deflected position. Otherwise the instrument may be damaged or the patient injured.*

To release the locking mechanism ⑮ from the secure position, turn it slowly to proximal.

10.3 Suction

Press the suction valve ② to start suction (when the suction pump has been connected and switched on).

CAUTION: *There should be no instrument in the suction channel during suction. The instrument channel inlet should be closed so that no secondary air is drawn in.*

10.4 Insufflation

Under anesthesia, connect the O_2 hose to the instrument channel:

- using a barbed tube (Art. no. 600007) on the instrument channel or
- if you are using the three-way double stopcock adapter with LUER-Lock (Art. no. 6927691), attach the hose to the tube holder on the angled part of the adapter. This ensures that the straight-line passage of the adapter remains free for instrument insertion.

10.5 Focusing

The image is focused with the focusing ring ①.

NOTE: *When connecting photo or video systems recommended by KARL STORZ, the fiberscope focusing ring must be at the zero position (mark).*

DCI® fiberscopes are focused using the adjustment wheel of the DCI® camera (see page 28).

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10.6 Instrument channel

Flexible instruments can be inserted through the instrument channel ⑧ as required.

If the instrument channel is not being used: replace the sealing cap.

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10.7 PDD version

Fiberscopes intended for use in Photodynamic Diagnosis (PDD) or autofluorescence (AF) are illuminated with a special diagnosis light (D-Light C). Adjust the PDD filter changer ⑫ for special filters designed for this light and possibly for marker substances:

ALA	for a special tumor marker based on aminolevulinic acid
W	for white light
AF	autofluorescence mode

ⓘ ***NOTE:*** *If the filter changer is set to ‘ALA’ or ‘AF’, a yellowish image is produced when viewing through the integrated special filter.*

ⓘ ***NOTE:*** *No filter changer is available on PDD fiberscope art. no. 11272 CI. The filter (ALA) is permanently installed.*

For optimal illumination, use the fluid light cable 495 FP (2.5 m) or 495 FO (1.8 m) between the light source and PDD fiberscope.

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
11. DCI® fiberscope with camera

11.1 Assembly

Slide the DCI® into the snap-in DCI® coupling mechanism of the DCI® camera until the coupling snaps in automatically.

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Connecting the special light cable 495 DC or 459 DV to the DCI® camera head:

Insert the light cable plug into the DCI® camera head until it locks into place.

While doing this, hold the light cable by the handle piece. Never pull on the cable!

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
NOTE: *The entrance and exit faces of the light cable and endoscope must be kept clean if optimal light delivery is to be maintained. Clean the surfaces with a cotton applicator soaked in cleaning solution, then disinfect it with a cotton applicator soaked in 70 % isopropyl alcohol and dry.*

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DCI® fiberscopes are focused using the adjustment wheel on the DCI® camera.

The camera control unit is controlled using the two programmable control keys (silver, blue) on the DCI® camera head.

NOTE: *Further controls are described in the instruction manual of the camera control unit used.*

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11.2 Disassembly

Hold the DCI® camera and the DCI® fiberscope still. Push the release button on the DCI® camera: the 'snap-in' DCI® coupling mechanism is released.

Pull the light cable out of the DCI® camera head.

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12. DIN eyepiece on DCI® fiberscope

The DCI® eyepiece adaptor 20260031 ensures the connection of DCI® telescopes to camera heads with DIN eyepiece mounts.

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12.1 Assembly

Slide the DCI® fiberscope into the 'snap-in' DCI® coupling mechanism of the DCI® eyepiece adaptor until the coupling snaps in automatically.

Connect a light cable or a battery-powered light source.

WARNING: *The instruction manuals of the supply units used must be observed.*

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12.2 Disassembly

Remove the light cable or battery-powered light source.

Hold the DCI® eyepiece adaptor and the DCI® fiberscope. Press the release button on the DCI® eyepiece adaptor. The 'snap-in' DCI® coupling mechanism is released.

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13. DCI® camera on DIN eyepiece

13.1 Assembly

Put the focusing ring (if available) of the fiberscope in the zero position (mark).

Slide the DCI® adaptor 20260030 into the 'snap-in' DCI® coupling mechanism of the DCI® camera until the coupling snaps in automatically.

Connect the standard light cable at the light inlet piece ⑦ of the fiberscope.

NOTE: *The special light cable 495 DC or 495 DV cannot be used and may have to be removed from the DCI® camera.*

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Opening the instrument coupling: turn the outer ring against the spring force, insert the fiberscope DIN eyepiece and let go of the outer ring. The instrument coupling closes automatically as a result of the spring force.

CAUTION: *No further tightening of the outer ring to position the eyepiece is required. This may, among other things, cause damage to the instrument coupler!*

13.2 Disassembly

Hold the DCI® camera and the fiberscope still. Open the instrument coupler for the DCI® adaptor and remove the fiberscope.

Hold the DCI® camera and the DCI® adaptor still. Push the release button on the DCI® camera: the 'snap-in' DCI® coupling mechanism is released. Remove the DCI® adaptor.

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14. Cleaning, disinfection, sterilization

Cleaning, disinfection, care and sterilization

CAUTION: *The manual ‘Cleaning, sterilization and care of KARL STORZ instruments’ must be followed. There, the cleaning, disinfection and sterilization procedures are explained in detail.*

CAUTION: *For cleaning and disinfection, use a product from the KARL STORZ approved chemicals list. The current list of approved products can be found on our website (www.karlstorz.com).*

CAUTION: *When preparing and using the solutions, follow the chemical manufacturers’ instructions paying close attention to proper concentration and exposure times. Prolonged immersion and incorrect concentrations may result in damage. Bear in mind the microbiological range of action of the chemicals used.*

WARNING: *The stated cleaning, disinfection and, where necessary, sterilization methods alone do not guarantee proper reprocessing or sterility. This can only be achieved if the staff follows recognized and validated reprocessing methods.*

WARNING: *For all cleaning and disinfection of contaminated endoscopes and instruments, oberserve the guidelines of the Employers’ Liability Insurance Association and equivalent organizations.*

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Cleaning, disinfection, sterilization

WARNING: *Remember that sterilization can only be performed successfully if*
- *all organic and chemical residue on the instruments has been eliminated. Sterilization is only effective on clean, disinfected surfaces!*
- *work proceeds according to recognized and validated methods.*
- *any deviation from recommended procedures must be validated by the user.*
- *the recommended sterilization parameters are maintained. These are only valid with sterilization equipment that is properly maintained and calibrated.*

CAUTION: *Never clean fiberscopes with a metal brush.*

CAUTION: *Constantly changing between different reprocessing methods causes particular wear on the materials and must therefore be avoided. A once-only change in methods however, e.g., when acquiring new cleaning equipment, is perfectly safe.*

CAUTION: *Fiberscopes must not be cleaned in an ultrasonic bath.*

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CAUTION: *Never immerse fiberscopes in normal saline solution, since even brief contact can lead to pitting and corrosion.*

CAUTION: *Final rinsing during manual reprocessing of fiberscopes should be carried out using microbiologically pure/sterile water.*

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Reprocessing Overview 1

Product	Cleaning			Disinfection	Sterilization					Remarks
	manual	Cleaning machine 60 °C	Cleaning machine 93 °C	Chemical	Steam 134 °C, fractionated	Gas EO	Gas FO	Plasma STERRAD	Chemical STERIS	
110xx ...,111xx ...,112xx ...,113xx ... Fiberscopes	✓	✓	-	✓	-	✓	✓	✓	✓	
495 F/G; Light cable daptor	✓	✓	✓	✓	✓	✓	✓	✓	✓	
11025 E; Pressure compensation cap	-	-	-	-	-	✓	✓	✓	-	
11301 CD; Adaptor for cleaning working channel in cleaning machines	✓	✓	✓	✓	-	✓	✓	✓	✓	
11301 CF; LIPP tube holder	[1]✓)	[1]✓)	[1]✓)	✓	✓	✓	✓	✓	✓	[1]) dismantle for cleaning
11301 D1, 11301 D3; Battery light source (accessory)	[5]✓)	-	-	[5]✓)	-	✓	✓	✓	✓	[5]) Battery compartment with closed cover; liquid <15 cm deep; no pressure impact: rinse only with unpressurized liquids; do not dry with compressed air
11301 FF2; Adaptor for leakage test in ETD2/ETD3 cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	

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Reprocessing Overview 2

Product	Cleaning			Disinfection	Sterilization					Remarks
	manual	Cleaning machine 60 °C	Cleaning machine 93 °C	Chemical	Steam 134 °C, fractionated	Gas EO	Gas FO	Plasma STERRAD	Chemical STERIS	
11301 GG2; Adaptor for cleaning irrigation and working channel in ETD2/ETD3 cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	
11301 HH; Adaptor for leakage test in BHT cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	
11301 KK2; Adaptor for cleaning working channel in ETD3 cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	
13242 XA; Adaptor for leakage test in cleaning und disinfection machines	[2]✓)	[2]✓)	[2]✓)	✓	✓	✓	✓	[3]✓)	✓	[2]) only PH neutral cleaners [3]) change in color possible
13242 XL; Leakage tester	-	-	-	-	-	-	-	-	-	
13272; Mouth piece	✓	✓	✓	✓	✓	✓	✓	✓	✓	
29100; Plug for instrument channel	✓	✓	✓	✓	✓	✓	✓	✓	-	
600007; Barbed tube, thin, with LUER-Lock	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6011590; LUER cap	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6927691; Three-way double stopcock adaptor with LUER-Lock	✓	✓	✓	✓	✓	✓	✓	✓	✓	Disassemble stopcocks before cleaning and sterilizing

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Reprocessing Overview 3

Product	Cleaning			Disinfection	Sterilization					Remarks
	manual	Cleaning machine 60 °C	Cleaning machine 93 °C	Chemical	Steam 134 °C, fractionated	Gas EO	Gas FO	Plasma STERRAD	Chemical STERIS	
11301 CB; Suction valve insert, reusable version	[4]✓)	[4]✓)	[4]✓)	✓	✓	✓	✓	✓	✓	[4]) dismantle for cleaning
5917800; Seal ring (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
5923110; Spring (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6885391; Cap (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6927791; Housing (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
8567290; Plunger (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
20260030; DCI® adaptor	✓	✓	-	✓	-	✓	✓	✓	✓	
20260031; DCI® eyepiece adaptor	✓	✓	-	✓	-	✓	✓	✓	✓	

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15. Preliminary cleaning

Remove any heavy soiling (secretions, etc.) from the outer casing of the fiberscope immediately after use. Use a disposable cloth moistened with cleaning solution (no solutions predominantly containing alcohol). Submerge the distal end in a container with cleaning solution and activate the suction to clean the working channel.

CAUTION: *Before immersing the instrument in cleaning and disinfectant solution, check that it is leakproof, as otherwise penetration of liquid may cause damage.*

CAUTION: *For miniature endoscopes with a removable eyepiece, the eyepiece must either be sealed with a plug or firmly attached to the endoscope before sterilization.*

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Connect the leakage tester (art. no. 13242 XL) to the valve for pressure compensation ⑤ and build up pressure until the pointer is in the blue marked range. Watch the pointer for approx. 30 seconds: the pressure must not drop during this time.Press the red button on the pressure valve and release the pressure. Remove the test unit from the fiberscope.

NOTE: *During inflation, make certain that the red pressure relief button is completely closed.*

If there is a loss in pressure, the fiberscope has a leak and must neither be used nor immersed in liquid. In this event, send the fiberscope in for repairs.

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Preliminary cleaning

WARNING: Risk of infection: *In order to prevent infection, it is strictly forbidden to ship contaminated medical devices!*
All medical devices must be decontaminated on-site*, to avoid contact and aerogenous infections among staff. To eliminate the risk of infection, before the fiberscope is dispatched it should be* wiped-down desinfected and sterilized in the NT procedure with the pressure compensation cap in place*. We reserve the right to return contaminated instruments/equipment to the sender.*

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- Remove all instruments from the instrument channel ⑧.
- Remove adapter and other removable parts from the light inlet piece ⑦ and instrument channel ⑧.
- Remove suction valve insert ② from the valve nozzle ③: press the projections on the side together and pull out the suction valve insert.
- Dispose of disposable parts removed. If necessary, disassemble other parts further and reprocess separately.

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15.1 Disassembling suction valve insert

ⓘ ***NOTE:*** *Only the multi-use suction valve insert (silver-coloured cap, metal plunger) is suitable for reprocessing.*

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
16. Manual reprocessing

NOTE: *Before cleaning, remove both adaptors from the light inlet piece ⑦ to permit free access to the fiber surface.*

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16.1 Manual cleaning

Submerge the instrument in a cleaning solution approved by KARL STORZ and clean it using a cleaning agent (brush or cleaning gun, Art. no. 27660). Avoid splashing with the cleaning agent.

CAUTION: *Pass the brush through the channel in only one direction and pull back only after the brush head has completely left the distal end. The brush must not be pulled back and forth within the channel, since this can lead to damage.*

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Remove the brush from the channel carefully. Then thoroughly rinse the channels repeatedly using cleaning and disinfection agents. If required, add a cleaning agent booster.

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Detach the adaptor, clean the stopcocks and grease (see the ‘Cleaning, sterilization and care of KARL STORZ instruments’ instructions).

The adaptor 11301 CD can be used for manual rinsing of the suction/irrigation channel.

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Rinse to remove the cleaning solution.

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Clean the eyepiece and objective lens with a soft cloth, sponge or cotton swab, soaked with 70 % alcohol.

CAUTION: *Cleaning with the special cleaning paste should only be carried out if the image as viewed through the endoscope is visibly cloudy or if there are visible deposits on the optical surfaces, it should certainly not be used routinely during reprocessing.*

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Stubborn soiling and deposits may be removed with cleaning paste (art. no. 27661). To do this, a small amount of cleaning paste is placed on a clean, moistened cleaning swab and distributed over the glass surface. Through polishing, even stubborn deposits can be removed.

NOTE: *The instrument must be reprocessed again following treatment with the cleaning paste!*

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Finally, dry the surfaces of the fibers and the optical end faces with a cotton-tipped applicator soaked in 70 % alcohol.

NOTE: *The fiber surface in the light inlet piece ⑦ must be dried off especially thoroughly with alcohol. Residue from disinfectants and detergents in the light inlet piece could bake on if the fiber optic light cable is connected and severely impair light transmission.*

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16.2 Manual disinfection

All parts of the KARL STORZ fiberscope can be disinfected. For this purpose, materials approved by KARL STORZ and recommended for use with endoscopes should be used. For a complete list, please visit www.karlstorz.com.

For disinfection, we recommend a plastic basin (Art. no. 27647 V), in order to avoid scratches on the optical end faces in particular. A plastic basin also prevents electrolytic corrosion, which can occur if different metals are present in a solution.

For disinfection, submerge the fiberscope completely. Ensure that the lumina of the channels and sheaths are completely filled with the solution. There must be no more air (air bubbles) in the channels. Make sure that all the stopcocks, the valve insert and the valve cap are open or have been disassembled.
Immerse the valve insert, adaptors etc. separately.

Rinse out the instrument channel ⑧ with a syringe filled with disinfectant solution. Follow the chemical manufacturers' instructions for concentration and exposure times.

For each rinse, use new microbiologically pure/sterile water. It is absolutely necessary to rinse the fiberscope with microbiologically pure/sterile water in order to prevent potential toxic impact from any remaining chemical residue, which could cause the instrument to become conductive. This is also required to optimize the results of sterilization and avoid damage to the surface.

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17. Machine reprocessing

After appropriate preliminary cleaning

After the appropriate preliminary cleaning (see p. 43), all KARL STORZ fiberscopes may be connected to devices for semi-automatic or automatic reprocessing. The choice of reprocessing programs must be agreed with the machine manufacturer. There are usually special programs available for fiberscopes which are suitable for cleaning and which are gentle on materials. For machine reprocessing, make sure that the fiberscopes are not heated above 65 °C.

CAUTION: *Check the fiberscope for leaks before every reprocessing in semi- or fully automatic equipment.*

Leakage tester connection on washer/ disinfector:

- ETD3/ETD2 washers:
 Adaptor for leakage test (Art. no. 11301 FF2, see pp. 52/55)
- BHT washer/disinfector:
 Adaptor for leakage test (Art. no. 11301 HH)
- Belimed washer/disinfector:
 Adaptor for leakage test (Art. no. 11301 BM)
- On the fiberscope using the green adaptor included in delivery (Art. no 13 242 XA);
 On the machine: consult the manufacturer of the washer/disinfector.

If the reprocessing machine cannot be connected to the leakage testing unit for any reason, check for leakage manually as described on page 43 or check the adaptor connection using the 13242 XA.

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For machine reprocessing, attach the green cleaning adapter (Art. no. 11301 CD) for the reprocessing of suction and irrigation channels in cleaning machines to a reprocessing tube and position on the valve connections. Use the cleaning adapter (Art. no. 11301 GG2, see page 8) for reprocessing in ETD2/ETD3 machines. In addition, use the special connecting adapter (Art. no. 11301 KK2, see page 8) for cleaning in ETD3 machines. For fiberscopes with a working channel without a suction valve, connect the tube of the washer and disinfector up to the working channel ⑧. To avoid pressure losses in the washer and disinfector, lock unused connections with the LUER-Lock (Art. no. 6011590).

WARNING: *Please observe the instruction manual of the washer and disinfector.*

Lay the fiberscope carefully in the tray of the machine. Connect the rinsing tubes to the channels. Proceed with the endoscope cleaning and disinfection program according to the regulations of the machine manufacturer.

CAUTION: *After the program has finished, the fiberscopes (including all adapters except 11301 CD) are highly susceptible to corrosion due to the ambient conditions in the reprocessing machine.*

Remove the fiberscopes from the machine as soon as possible after the program has finished. Then dry the optical surfaces with a cotton-tipped applicator soaked in 70 % alcohol.

NOTE: *If, according to the high risk classification and/or the Infection Protection Act, sterilization of flexible fiberscopes is not required, the instruments should be stored safely in a hanging position free from dust (distal tip pointing down), so that any remaining moisture in the channels can drain out.*

NOTE: *For flexible fiberscopes without channels, a hanging storage position is not essential.*

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Connection for machine reprocessing in ETD3 machines

Dismantling the adaptor 11301 GG2:
Press the button corresponding to each tube and pull the tube off simultaneously.

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
1. Slide the free tube end ① of the adaptor (Art. no. 11301 GG2) onto the olive connector of the green adaptor (Art. no. 11301 CD).

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ⓘ ***NOTE:*** *For fiberscopes with a working channel without a suction valve, pull the tube not required, e.g. ① out of the metal block.*

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2. For machine reprocessing, attach the green adaptor (Art. no. 11301 CD) to the valve connections of the fiberscope.

3. Attach the tube with LUER-lock to the fiberscope instrument channel ⑧ and secure it.

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4. Insert the other tube ends in the metal part of the adaptor (Art. no. 11301 GG2).

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5. Connect the adaptor (Art. no. 11301 FF2) to the ‘valve for equalizing pressure and connecting leakage tester’.

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6. Connect the leakage tester to the ETD3 using the adaptor 11301 FF2.

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
7. Connect the cleaning adaptor for cleaning the instrument channel in EDT3 cleaning machines (Art. no. 11301 GG2) to the connector of the supply tube. Then connect up the special connecting adapter for ETD-3 machines (Art.-Nr. 11301 KK2).

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8. Lay the instrument carefully in the tray of the ETD3 machine.

9. Connect the rinsing tubes to the ETD3 machine.

WARNING: *To avoid larger pressure losses in the EDT3 machine, connections which are not in use should be locked with LUER locks.*

10. Screw the adaptor 11301 KK2 onto the designated connections of the ETD3 machine, according to the figure.

11. Use the cleaning and disinfection program for fiberscopes according to the ETD3 manufacturer’s instructions.

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
After cleaning has been completed, remove the green adaptor (11301 CD) for cleaning the working channel by pressing together the projections at the side.

ⓘ ***NOTE:*** *The green adaptor for cleaning the working channel may be left on the tube for later use.*

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Remove the connection from adapter 11301 FF2.

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
❗ ***CAUTION:*** *The fiberscope could be damaged by liquid which may enter the device via an adaptor for leakage testing (11301 FF2) left on the fiberscope after machine reprocessing. Therefore, you should always dismantle the adaptor (11301 FF2) immediately after cleaning is completed!*

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To disassemble adaptor 11301 GG2:
Press respective button and pull off tube simultaneously.

ⓘ ***NOTE:*** *If, according to the high risk classification and/or the Infection Protection Act, sterilization of flexible fiberscopes is not required, the instruments should be stored safely in a hanging position free from dust (distal tip pointing down), so that any remaining moisture in the channels can drain out.*

ⓘ ***NOTE:*** *For flexible fiberscopes without channels, a hanging storage position is not essential.*

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Connection for machine reprocessing in ETD2 machines

Push free end of the adaptor tube (Art. no. 11301 GG2) onto the barbed tube (Art. no. 11301 CD).

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Place the adaptor 11301 CD on the valve nozzle ③ of the fiberscope.

Insert the other end of the tube into the metal part of adaptor 11301 GG2.

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Place the LUER-Lock on the working channel ⑧ of the fiberscope.

Insert the other end of the tube into the metal part of adaptor 11301 GG2.

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Insert the supply tube with the fitting for ETD2 cleaning machines into the metal part of adaptor 11301 GG2.

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Machine reprocessing

Attach adaptor 11301 FF2 to the valve ⑤ for pressure compensation and the leakage tester connection of the fiberscope.

Connect leakage tester connection of the ETD2 cleaning machine.

Lay the endoscope carefully in the machine basket. Connect the irrigation tubing to the channels. Proceed with the endoscope cleaning program according to the regulations of the machine manufacturer.

CAUTION: *After the program run, the fiberscopes are greatly at risk of corrosion in the reprocessing machine.*

The fiberscopes (including all adapters except 11301 CD) must be removed from the machine as soon as possible once the program has finished.

Then if necessary dry off the optical surfaces with a cotton swab soaked with 70 % alcohol.

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18. Inspection after reprocessing

Before assembly it should be checked that all parts and surfaces are undamaged, free from any organic or chemical residue as well as completely clean and dry.

Replace defective parts or send in for repair. Before being sent in for repair, all parts must be cleaned, wipe-down disinfected and sterilized in the NT procedure if necessary.

If residue is still present after reprocessing, reprocess the parts again, subject them to special treatment, discard or repair. Remove stubborn deposits on glass surfaces with a special cleaning paste (Art. no. 27661).

The part must then be submitted to a complete reprocessing cycle.

Check plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their new condition, they must be replaced. This applies in particular to seals and rubber caps.

Test the instruments for proper operation: check the telescope for a clear view, the instrument for steerability and all canals for patency.

WARNING: *Damage to the tube: The use of ULTRASTOP as an anti-fog solution on the optically distal end faces of the flexible endoscopes can damage the tube connector.*

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19. Gas sterilization/plasma sterilization

Sterilization using formaldehyde or ethylene oxide and plasma sterilization are gentle and therefore recommended methods of sterilization.

NOTE: *Observe national laws and regulations for the gas sterilization method.*

NOTE: *Fiberscopes which are to be used for interventions in areas of the body which lack microbial colonies, must be sterilized in sterile packaging and stored in a closed cupboard free from contamination and dust.*

CAUTION: *If the instrument is to be sterilized, the entire system must be opened by fitting the red cap on the pressure compensation part to ensure pressure compensation.*

WARNING: *Before commencing gas sterilization, the instrument must be completely dry or else sterile conditions will not be achieved.*

- Place the instrument securely on a sterilization sieve.

Ethylene oxide gas sterilization has been validated under the following parameters:

Gas mixture	EO: HCFC*
Temperature (°F, °C)	130 °F ± 5 °F; 54 °C ± 2 °C
Relative air humidity	60 ± 20%
Pressure (p_e, kpa)	58 … 70 (0.56 … 0.7 bar/8 … 10 psi)
Exposure time	120 mins
EO concentration	600 ±30 mg/L

* Gas mixture consisting of 10 % ethylene oxide and 90 % chlorotetrafluoroethane (HCFC-124; percent by weight)

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	at 20 °C	at 42 °C	at 55 °C
Telescopes	24 hrs	18 hrs	12 hrs
Rubber items	7 days	4 days	3 days
Fiberscopes/ Plastics	5 days	2 days	30 hrs

(minimum time in simple transparent packaging and an application time on the patient of more than 30 minutes.)

Maximum acceptable levels (per 21 CFR 812.100; Federal Register, Vol. 43, No. 122, June 23 1978) of residues following EO sterilization are as follows:

Ethylene oxide: 25 ppm* / 250 ppm**
Ethylene chlorohydrine: 25 ppm* / 250 ppm**
Ethylene glycol: 250 ppm* / 5000 ppm**

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The instruments can be used immediately after sterilization with formaldehyde or plasma sterilization. Extra airing time is not necessary.

Remember if you use plasma sterilization that there are limits regarding the size of the lumen and material in regard to instruments which may be sterilized using STERRAD®.

Sterilization using STERRAD NX is also possible.

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Gas sterilization/plasma sterilization

Statements as to the efficiency of sterilization for the respective type of instrument (depending on the length and diameter of the lumina) may only be made by the manufacturer.

WARNING: *Sterility is only assured if the manufacturer's instructions for the STERRAD® procedure are followed.*

WARNING: *All deviations from the sterilization parameters recommended for the STERRAD® system must be validated by the user.*

CAUTION: *Attention must be paid to the information on reprocessing, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. The methods for cleaning, disinfection and sterilization are explained there in detail.*

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20. Chemical sterilization

Chemical sterilization (STERIS®) is a special method which is particularly suited for thermolabile goods such as fiberscopes.
Except for telescopes with a HAMOU® eyepiece drive mechanism, all KARL STORZ telescopes are suitable for STERIS® sterilization. Color changes to the surface of aluminum alloys could occur during sterlization, which do not however affect the function of the instrument.
It must be pointed out in this context that anodized aluminum parts are not always recognizable as such.

WARNING: *Sterility is only guaranteed if the manufacturer's instructions for the STERIS System have been followed. This also applies to use of the sterilization containers.*

WARNING: *Chemical sterilization is not possible on greased or oiled surfaces.*

WARNING: *Parts with a lumen must be reprocessed together with the STERIS Quick Connect kit. All parts must be thoroughly cleaned before sterilization using the STERIS PROCESS™, as also applies to the other sterilization methods.*

CAUTION: *Attention must be paid to the information on cleaning, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. The methods for cleaning, disinfection and sterilization are explained there in detail.*

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21. Storage

CAUTION: *Do not expose fiberscopes to direct sunlight or excessive heat.*

Sterilized fiberscopes should be packaged and stored in sterile cloths or suitable containers in accordance with EN 868, and must be absolutely dry before being put into storage.

Storage must take place in facilities with normal temperature and air humidity levels.

You can find special holders for this in the ‘Cleaning, sterilization and care of KARL STORZ instruments’ instructions.

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22. Accessories

Flexible forceps

Carefully wipe the outer surfaces of the flexible forceps with a disposable cloth soaked in disinfectant. Clean the forceps jaws with a brush.

Next, clean in an ultrasound bath. Afterwards, rinse with microbiologically pure/sterile water.

Manual/machine cleaning/disinfection

Both processes can run parallel to decontaminate flexible forceps.

Exception: If the flexible forceps are thermoresistant, they can be cleaned and disinfected in the 93 °C program.

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For handles which can be dismantled:

Dismantle handle as follows:

- Undo the thumbscrew on the handle and remove the rear part of the handle.
- Undo the thumbscrew on the connecting piece and remove the front part of the handle to the rear.

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Also oil the joints of the jaws by holding the instrument so that the oil can flow into the instrument channel and wet the pulling wire. Then reassemble in reverse order.

The sterilization of flexible forceps is compulsory!

ⓘ ***NOTE:*** *Only use special oils whose lubricating effect remains even in autoclaves (art. no. 27656).*

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Papillotome/Polypectomy snarn

Dismantle the device completely for reprocessing.

- Unscrew thumb handle.
- Undo contact socket.
- Unscrew grip piece.
- Remove inner guide sleeve.
- Pull papillotomy or polypectomy loop out of the plastic probe at the front.

Then reassemble the polypectomy snare in reverse order.

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Both processes can run parallel to decontaminate papillotome/polypectomy snares.

Exception: If the papillotome/polypectomy snares are thermoresistant, they can be cleaned and disinfected in the 93 °C program.

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Wipe the exterior of the instruments with a cloth soaked in disinfectant. Clean the forceps jaws with a disinfected brush.

Then clean in an ultrasound bath.

Afterwards, rinse with microbiologically pure/ sterile water.

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
Assembly of papillotomes:

- Screw inner guide sleeve into the grip piece
- Fit contact socket.
- Insert loop into plastic probe as follows: insert long end of loop into the rear opening of the probe and push to the rear, keeping the probe stretched.
- Insert shorter end of loop into the front opening of the probe and slide both parts towards the rear keeping them parallel.

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In polypectomy snares, the loop end is already shaped. Insert the proximal loop end into the plastic probe from the distal end and push towards rear.

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24. Normatives

24.1 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

24.2 Limitation of liability

As the supplier of this instrument, KARL STORZ shall be liable for the safety, operational reliability and performance of this instrument subject to the following conditions:

assembly, system expansion, readjustment, modification, or repair have been performed by persons authorized by KARL STORZ and that the instrument is used in accordance with the instruction manual.

24.3 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see “Subsidiaries” section), even during the warranty period.

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25. Repair and replacement program

KARL STORZ maintains a repair and replacement warehouse which is normally adequate to ensure prompt replacement of damaged telescopes and instruments. Under the repair and replacement plan, you receive an identical as-new instrument and are only charged the repair costs for the defective instrument. For telescopes, you receive a warranty of 1 year, and for instruments 6 months.

For fiberscopes and equipment, individual repair is necessary. To bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Straße 34
D-78532 Tuttlingen/Germany

In other countries please contact your local KARL STORZ branch or authorized dealer.

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25.1 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havanna, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsätragränd 14
12739 Skärholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoskopy Kasachstan
6, Saryarka str. BC „Arman", office314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd .
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang Intenational Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: , Fax: +63 2 31745-00
E-Mail: phillippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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KARL STORZ — ENDOSKOPE

GEBRAUCHSANWEISUNG
SILVER SCOPE® Video-Koloskope,
Modelle 1392x

INSTRUCTION MANUAL
SILVER SCOPE® Video colonoscopes,
Models 1392x

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Important information for users of KARL STORZ instruments

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of modern, high-quality items of equipment from KARL STORZ.

This instruction manual is designed to serve as an aid in the proper use, cleaning, disinfection and, if relevant, sterilization of videoscopes manufactured by KARL STORZ. The videoscopes are sensitive precision instruments which require particularly careful reprocessing. In order to avoid damaging the instruments during reprocessing, help is provided here for the disassembly, cleaning, disinfection, gas sterilization and assembly of videoscopes from KARL STORZ. Therefore, read the instructions carefully and keep them to hand for future reference.

KARL STORZ is continuously working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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Content

1 Safety instructions....11
2 Intended use11
3 User qualification....11
4 Warnings and cautions12
5 Safety precautions when using the videoscopes13
6 Handling....18
6.1 Compatible video processor units/light sources18
6.2 Connecting the videoscope to a light source illustrated with the XENON 100 SCB19
6.3 Operating the videoscope....20
6.3.1 Preparing the light source....20
6.3.2 Preparing the processor unit....20
6.3.3 Connecting the light conductor, air/water connections and the suction channel.....20
6.3.4 Deflection of the distal tip21
6.3.5 Suction21
6.3.6 Insufflation and light21
6.3.7 Water irrigation for lens cleaning22
6.4 Preparation before use22
6.4.1 Testing the various functions....22
6.4.2 Handling the endoscope correctly25
6.4.3 Additional recommendations for use26
6.4.4 Operating the endoscope buttons27
7 Reprocessing28
7.1 General information....28
7.2 General cleaning....30
7.2.1 Pre-cleaning31
7.2.2 Leakage test....32
7.3 Channel system of a KARL STORZ SILVER SCOPE® video colonoscope....35
7.4 Channel system of a KARL STORZ SILVER SCOPE® video colonoscope with additional rinsing channel (water-jet channel)....36

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
Contents

7.5 Cleaning/disinfection of the videoscope: Manual cleaning 37
7.6 Manual disinfection 41
7.7 Machine reprocessing of the videoscope 43
7.7.1 Cleaning adaptor set 13991 CK 43
7.8 Preliminary cleaning directly after the examination (see page 31) 49
7.9 Cleaning in the reprocessing place 49
7.10 Reprocessing in the washer/disinfector 50
7.11 Cleaning of accessories 50
7.12 Reprocessing the water bottle 13992 BS 51
7.13 Checking the reprocessing procedure 51
7.14 Sterilization 52
7.14.1 Gas sterilization 53
7.15 Storage 55
8 Storage/ transport conditions 55
9 Maintenance 56
10 Important information 56
11 Servicing and repair 57
12 Limitation of liability 57
13 Warranty 57
14 Directive compliance 57
15 Appendix: Troubleshooting 58

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Symbols employed

Symbol	Meaning
	Consult instructions for use
	Serial number
	Manufacturer
	Fragile, handle with care
	Keep dry
	Temperature limit
	Do not reuse
	Applied part type BF, defibrillation proof, IEC 60601-1
	Potential equalization connector
	DANGER: Risk of explosion if used in the presence of flammable anesthetics.
	Protective earth (ground)
	Alternating current
	Device is subject to the requirements of the WEEE Directive, 2012/19/EU.

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Symbols employed

Symbol	Meaning
	CAUTION: Do not open the device! Risk of electric shock. Refer servicing to qualified service personnel. Do not store liquids on or above the device. Keep the device out of reach of patients.
	On
	Off
	Video input
	Video output
	Color video camera
	Color video monitor
	White balance
	Output
	Input

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
Symbols employed

Symbol	Meaning
	Setup Wizard
	Signifies that this feature/option is available in the Options Menu.
	Signifies that this feature/option can be controlled via the User Menu.
	Signifies that this feature/option can be controlled via the buttons on the videoscope.
RxOnly	Federal (USA) law restricts this device to sale by or on the order of a physician.

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




Images
SILVER SCOPE® Series
Colonoscopes

① Up/down deflection wheels
② Up/down deflection brake
③ Programmable remote control buttons
④ Suction valve
⑤ Air/water valve
⑥ Working channel with rubber lip valve/biopsy valve
⑦ Right/left deflection wheels
⑧ Right/left control brake
⑨ Additional irrigation channel with LUER lock connection (only 13925 xKS, 13920 PKS/ NKS, 13926 PKS/NKS, 13927 PKS/NK, 13925 HPKS/NKS, 13924 HPKS/NKS
⑩ Nozzles for lens cleaning/irrigation/insufflation
⑪ Light cable
⑫ Lens
⑬ Distal water jet channel (only for 13925 PKS/NKS, 13920 PKS/NKS, 13926 PKS/NKS, 13927 PKS/NKS, 13925 HPKS/NKS)
⑭ Supply tube
⑮ Connection for leakage test
⑯ Connection for water bottle
⑰ Light cable connection
⑱ Port for video processor/electronic connector
⑲ Cap with connection for leakage tester
⑳ Suction/aspiration port
㉑ Air pin
㉒ Lanyard
㉓ Working channel

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1 Safety instructions

Please read these safety instructions carefully before using the device. Before using the instrument on the patient for the first time, it is imperative that you be acquainted with how the instrument operates and is controlled.

2 Intended use

The KARL STORZ video colonoscopes are optical instruments for medical, diagnostic and endoscopic treatment in the lower gastrointestinal tract of humans.

They are introduced through the anus, via the rectum and the entire colon into the lowest sections of the small intestine.

For use with accessories specified by KARL STORZ for this purpose.

Contraindications:

No contraindications relating directly to the product are currently known. Videoscopes must not be used if, in the opinion of the responsible physician, the patient's health is endangered through their use, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

3 User qualification

KARL STORZ videoscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the endoscopic technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the videoscopes.

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4 Warnings and cautions

Please read this manual carefully and follow the instructions exactly to avoid putting your patients, personnel or yourself at risk. The words **Warning, Caution** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be endangered. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or safety precautions must be followed to avoid any damage to the instrument.

NOTE: A Note indicates special information about operating the instrument, or clarifies important information.

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WARNING: *Combinations of medical devices are only assured to be safe if*

- *they are identified as such in the respective instruction manuals or*
- *the intended use and the interface specifications of the devices used in combination permit this.*

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WARNING: *The pressure compensation/EO cap 11025 XE must be removed from the connection piece during:*

- *endoscopic interventions*
- *manual/machine reprocessing*

WARNING: *Do not allow any live components or instruments to come into contact with the videoscope.*

WARNING: *TThe videoscope must not be used during the discharge of a defibrillator. Remove the videoscope from the operation area before discharge occurs.*

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Safety instructions

WARNING: *High energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the fiberoptic light cable and at the tip of the endoscope.*

Prolonged direct contact of light-emitting surfaces on tissue can lead to burns; this must therefore be avoided at all costs.

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WARNING: *Always unplug the videoscope before carrying out any maintenance and cleaning procedures on it.*

WARNING: *Risk of infection during shipment! Contaminated videoscopes may cause serious or life-threatening illnesses during transportation or on delivery to the recipient.*
Please cover the carrying case with the protective sheet 13990 SFN before sending contaminated instruments back to KARL STORZ.

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Safety instructions

WARNING: *If the endoscope should malfunction during an intervention on the patient, stop the intervention immediately. Place the distal tip of the videoscope in a neutral, nondeflected position and remove the videoscope slowly and carefully from the patient.*

WARNING: *The use of HF probes, HF electrodes, HF instruments or other energy-emitting probes is not permitted if explosive gases (for example, explosive/flammable anesthetics) are present.*

WARNING: *Patient leakage currents can accumulate if endoscopes and electrically powered endoscope accessories are used together.*

WARNING: *Before turning on the high voltage/electrical supply, ensure that there is no contact with metal parts of the endoscope or any other electricity-conducting accessory, and that the suction unit has no contact with conductive liquids.*

WARNING: *HF probes, HF electrodes, laser probes and other energy-emitting probes should only be used under constant, visual, endoscopic control.*
Activation should only occur when the tip of the energy-emitting probe is outside the working channel and within the endoscope's field of vision.

CAUTION: *Only use accessories (guide probes, cleaning brushes, etc.) that are specifically suited for the intended use with videoscopes and which can be employed in line with generally accepted technical standards and safety guidelines.*

CAUTION: *Only use accessories with flexible ends. Other types of accessories may cause damage to videoscopes.*

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Safety instructions

CAUTION: *Place the pressure compensation/EO cap 11025XE on the connection piece prior to:*

- *gas sterilization*
- *shipping*

CAUTION: *Do not expose the videoscope to direct sunlight or excessive heat, to avoid damage to the device.*

CAUTION: *Videoscopes must never be steam sterilized (autoclave).*

CAUTION: *The maximum permissible temperature for decontamination, sterilization, transport and storage is 65°C.*

CAUTION: *The manual ‘Reprocessing and care’ must be followed. The cleaning, disinfection, and gas sterilization procedures are explained therein in detail.*

CAUTION: *Do not hit the distal tip of the videoscope against hard objects, and never use a sharp object to remove dirt.*

CAUTION: *Never place heavy objects on the videoscope.*

CAUTION: *Do not use any accessories with sharp edges or points in or on the endoscope.*

NOTE: *The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts. In this case, the reduction of the patient leakage current should be effected using only CF application parts as far as this is possible.*

NOTE: *The endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against high frequency currents should be ensured through the use of correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.*

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Safety instructions

NOTE: *Any damage to the videoscope resulting from incorrect operation is not covered by the guarantee.*

NOTE: *Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided.*

NOTE: *Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.*

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6 Handling

6.1 Compatible video processor units/light sources

NOTE: *The KARL STORZ videoscope should not be used in combination with video processor units and light sources from other manufacturers. The use of non-compatible units can lead to patient injury or unit damage.*
It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other units at any time.

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
6.2 Connecting the videoscope to a light source illustrated with the XENON 100 SCB

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Connection of supply connector with light source and water bottle (model 13925 PKS/NKS,13925 HPKS/NKS,13927 PKS/NKS, 13926 PKS/NKS, 13920 PKS/NKS) is equipped with an additional connection/ irrigation channel (water jet channel) on the control piece.

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6.3 Operating the videoscope

6.3.1 Preparing the light source

Please read the full instructions in the instruction manual of the respective light source.

1. Install the light source on a firm, level surface.
2. Plug the power cord into the specified power socket.

6.3.2 Preparing the processor unit

Refer to the relevant instruction manual for the processor unit used.

6.3.3 Connecting the light conductor, air/water connections and the suction channel

1. Insert the supply connector into the light source until it audibly clicks into place (Fig. 1).
2. Connect the videoscope adaptor cable (VEA) to the endoscope. When doing so, take note of the white markings on the endoscope and the VEA cable (Fig. 2).
3. Attach the connector from the water bottle tube to the endoscope supply plug (irrigation bottle connector) (Fig. 3). The water bottle holder automatically positions itself correctly and cannot be adjusted manually.
 IMPORTANT: To fill the water bottle with sterile water, remove the bottle from the light source and unscrew the lid. Only fill the water bottle 2/3 full with sterile water (1/3 is required to cushion air pressure).
4. Connect the suction tube (patient tube) from the suction pump or central suction system to the supply plug connector (Fig. 4).

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6.3.4 Deflection of the distal tip

Rotating the large deflection wheel, labeled ‘U/D’, (up/down) deflects the distal tip. Rotating clockwise or forwards lowers the tip. Rotating counterclockwise or backwards raises the tip.

Right and left movements are controlled using the smaller deflection wheel (labeled ‘L/R’). Rotating clockwise or forwards = right/rotating counterclockwise or backwards = left (Fig. 5)

WARNING: *Do not rotate the deflection wheels forcibly; this could damage the instrument or traumatize the patient! Before attempting to introduce or remove the instrument, release the deflection wheel brakes!*

6.3.5 Suction

Switch on the suction pump. To perform suction, press the suction valve (with the red ring) on the endoscope (Fig. 6).

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6.3.6 Insufflation and light

Start the light source, switch on the insufflation pump.

For insufflation, cover the opening on the blue valve with a finger (do not depress the valve, since this starts water irrigation) (Fig. 7).

NOTE: *To avoid patient discomfort, use insufflation sparingly. Excessive insufflation may cause cramps, spasms and injuries!*

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6.3.7 Water irrigation for lens cleaning

Pressing on the air/water valve starts the water flow and/or lens cleaning. (Fig. 8).

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6.4 Preparation before use

WARNING: *Before being used for the first time, the instrument should undergo a complete reprocessing cycle, including preliminary cleaning.*

6.4.1 Testing the various functions


1. Prepare light source as described in 2.1.
2. Connect supply connector, water bottle and suction tube according to 2.3.
3. Inspect the supply tube for damage (such as breaks, cracks, twisted or squashed sections...) and check the supply connector for loose connections.
4. Examine the surface of the flexible insertion shaft visually and by touch for cracks, cuts, bulges, foreign bodies or other damage.

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5. Testing deflection – upward/downward control

The large deflection wheel controls up and down movements. To test up and down movements, release the up and down brake. Slowly rotate the large deflection wheel clockwise or forwards; the distal end should move downwards. Rotating the deflection wheel counterclockwise or backwards causes upward movement.

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Handling

6. Testing deflection – Right/left control
7. The smaller/outer deflection wheel controls right/left movement. To test right/left movement, release the right/left brake. Slowly rotate the small deflection wheel clockwise or forwards; ***the tip of the endoscope*** should move to the right. Rotating the deflection wheel counterclockwise or backwards causes movement to the left.
8. Testing deflection – brake function
 The deflection wheels should turn freely when the brakes are released.
9. Release the brakes and test their action (see 5/6). Activate the brakes by moving the brake lever counterclockwise – ensure that the angled position of the distal endpiece is fixed in the desired position.

Do not twist or bend the deflection piece by hand.

WARNING: *If the deflection and/or brake mechanisms do not function properly, do not use the device. Contact the KARL STORZ service department.*

10. Check the air supply. With the insufflation pump switched on, hold the distal end of the instrument under sterile water. Cover the air control opening in the air valve (blue ring) with a finger; bubbles should emerge from the air nozzle at the distal end.

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11. Check the function of the water supply: Depress the air/water valve (blue ring) fully. Water should flow out of the water nozzle (distal) and rinse the lens clean. IMPORTANT: When you first actuate the valve, it may take a couple of seconds for water to start flowing.

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12. Test the suction function. With the device correctly installed and the suction pump or central suction system switched on, immerse the distal end in a container of sterile water once again. Depress the suction valve (red ring) fully. Water should be suctioned through the suction channel of the instrument via the connected suction pump and into the collecting bottle of the suction system.

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13. Test the biopsy channel by introducing biopsy forceps.

CAUTION: *It should be possible to guide the biopsy forceps smoothly and without problems through the entire channel.*

CAUTION: *Full (>180°) deflection can lead to limitation in movement, depending on the accessories. If necessary, adjust the deflection slightly.*

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CAUTION: *Accessories such as biopsy forceps that are kinked or bent will not function properly, are difficult to guide through the channel and may cause severe damage to it. They must therefore be replaced and not used further! 9. Test the biopsy channel by introducing biopsy forceps.*

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Handling the endoscope correctly

Hold the control handle of the endoscope in the left hand in such a way that the thumb, index and middle fingers can easily control the large deflection wheel. In this position, the suction valve can be depressed with either the index or the ring finger. The air/water valve can also be covered with these fingers. The remote control button (e.g., Stills) is also within reach.

The small, outer deflection wheel can be controlled with the fingers of the left hand, or alternatively with the thumb and fingers of the right hand.

6.4.3 Additional recommendations for use

Observe the lumen the whole time. If you become disoriented, withdraw the instrument gently until vision is restored.

ⓘ ***NOTE:*** *If pollution causes a blurry picture during the procedure, this can normally be rectified by a combination of water and insufflation. Excessive water or air can be removed by applying simultaneous pressure on the suction valve.*

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6.4.4 Operating the endoscope buttons

Depending on the video processor used, the buttons can have two functions, which are initiated by means of short and long presses of the button. More detailed information on button assignment can be found in the manual for the camera control unit (as listed in section 6.1).

NOTE: *In case of 'red-out', a blurred red picture on the screen, caused by the distal endoscope tip resting directly on the mucosa, withdraw the instrument gently until vision is restored.*

WARNING: *Pushing the instrument forwards despite a 'red-out' effect may result in serious injury to the patient.*

NOTE: *If the angling of the distal tip can no longer be properly controlled, the examination must be ended immediately. Release the brakes and place the tip in the neutral position. Position the deflection wheels so that 'U' and 'L' are between the valves (see page 13/17); then carefully withdraw the device, whilst continuing to monitor the procedure.*

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7 Reprocessing

7.1 General information

WARNING: *The stated cleaning, disinfection and sterilization methods alone do not guarantee successful sterilization of the endoscopes. This can only be achieved if personnel follow recognized and validated reprocessing methods.*

WARNING: *The local occupational safety regulations must be followed and applied for all cleaning and disinfection work on contaminated endoscopes and instruments.*

WARNING. *Remember that decontamination is only then successful if*

- *the instruments have been freed of all organic material and chemical residues.*
- *only recognized and validated methods are used.*
- *any deviation from recommended procedures must be validated and tested by the user.*
- *the recommended sterilization parameters are maintained. These are only valid with cleaning, disinfection and sterilization equipment that is properly maintained and calibrated.*

CAUTION: *Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.*

CAUTION: *Never clean videoscopes with metal brushes.*

NOTE: *The manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’ (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.*

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CAUTION: *Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided. However, a one-off change of method, e.g., when acquiring new cleaning equipment, should be carried out in consultation with the manufacturers of the endoscope, the chemicals and the washer/disinfector.*
This change in method is then usually completely harmless.

CAUTION: *Videoscopes must not be cleaned in an ultrasonic cleaning bath.*

CAUTION: *Videoscopes must never be steam sterilized (autoclaved).*

CAUTION: *Never place videoscopes into physiological saline solution, since even brief contact may lead to pitting and corrosion.*

CAUTION: *During manual reprocessing, videoscopes may only be rinsed with microbiologically clean/sterile water.*

CAUTION: *Decontamination in washers and disinfectors is only possible using valid procedures.*

CAUTION: *Even parts that do not come into direct contact with the patient (e.g., video connecting cables) should be thoroughly disinfected following an intervention.*

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7.2 General cleaning

Cleaning should only be undertaken by authorized trained personnel. All safety regulations should be observed (gloves, protective glasses, mouth and nose protection, clothing, German Technical Rules for Biological Agents (TRBA 250)

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The national guidelines and recommendations for the reprocessing of instruments (for example RKI, ESGE/ESGENA, ASGE, LANGS, SAGES etc.) should be observed and the effectiveness verified.

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CAUTION: *Do not autoclave the endoscope or accessories (with the exception of accessories that are specifically designated as autoclavable, e.g. forceps and Albarran module).*

NOTE: *The endoscope must be cleaned immediately after every use. If organic material is allowed to remain and dry in the instrument, complete or partial blockage of the channel may result. Thorough cleaning is a prerequisite for effective reprocessing. This will avoid coagulation residue/blockages and the resulting repairs and instrument failures.*

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7.2.1 Pre-cleaning

1. The preliminary cleaning of the endoscope is performed at the place of operation. The main cleaning and disinfection is carried out in a separate reprocessing room.
2. Immediately after use, remove any heavy soiling (secretions, etc.) and organic material from the outer surfaces of the insertion tube with a disposable, lint-free cloth moistened with active cleaning solution. Then submerge the distal end in a container with active cleaning solution and activate the suction to rinse the working channel.
3. Actuate the suction valve and suction at least 100 ml of cleaning solution through the suction channel. (Fig. 2).
4. Use the cleaning valve 13991 SRV (Fig. 1) in place of the air/water valve (Fig. 2). Do not use this valve during examination.
5. Check that the water bottle is filled (max 2/3). Start the air supply (insufflation pump). Water and air are now passed through both the air channel and the water channel, depending on the position or actuation of the cleaning valve. (The process should run for at least 10 seconds).
 Air channel = Position 1
 Water channel = Position 2
6. It is important that both valve settings ('Actuated'/'Not actuated') are used multiple times for several seconds. This is the only way to ensure that both channels are rinsed through.
7. Switch off video system.
8. Pull out the videoscope connector from the videoscope adaptor (VEA) at the supply connector of the endoscope when the video processor unit is switched off. Put on the sealing cap for the video cable connector and lock it with a 1/4 turn.
9. Disconnect the endoscope from the video processor unit and the light source; place the cleaning valve in a cleaning solution.
10. Then transport in a protected case to the site of reprocessing.

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7.2.2 Leakage test

Preparing for the leakage test – dry

ⓘ ***NOTE:*** *Place the cap on the videoscope adaptor (VEA) connector and screw it on tightly. Ensure that the cap is secured properly.*

1. The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector to the supply connector and lock by turning clockwise 90°.

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2. Inflate the leakage tester until the maximum test pressure of 160 – 180 mmHg (blue area) is reached. The deflection rubber on the distal deflection part of the endoscope expands slightly. Move the deflection part with the deflection wheels on the control piece in all four directions as far as possible and check the manometer to see whether the pressure remains constant or falls (dry test).
3. Wait at least 1 minute.

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Performing the leakage test – in liquid

1. Immerse the device in water so that it is completely covered. It must be ensured that the leakage tester is not immersed but rather kept dry. Perform the leakage test as described above and observe whether bubbles emerge from the device. If after an observation period of at least 3 minutes no more bubbles can be seen emerging from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof and can be used.

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If this test is also negative and no bubbles can be seen rising in the liquid, the leakage tester can be removed. To do this, remove the supply connector from the liquid and press the vent button on the manometer to vent the endoscope.

NOTE: *Never release the pressure in the liquid. To do this, remove the connector from the liquid and press the vent button. The test pressure must fall to 0. Then remove the leakage tester.*

If a pressure loss occurs, this means that the videoscope leaks and may no longer be used or immersed or reprocessed by machine.

CAUTION: *In this event, fill out the accompanying repairs form (see Appendix), roughly clean and remove any organic soiling/residue from the endoscope, mark it as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the enclosed transport protection sheet 13990 SFN.*

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WARNING: *Risk of infection during shipment! Contaminated instruments can cause aerogenous or contact infections during transportation or at the recipient's premises.*
Shipping without the plastic film packaging must be avoided.

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
7.3 Kanalsystem eines KARL STORZ SILVER SCOPE® Video-Koloskops

7.3 Channel system of a KARL STORZ SILVER SCOPE® video colonoscope

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
7.4 Channel system of a KARL STORZ SILVER SCOPE® video colonoscope with additional rinsing channel (water-jet channel)

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7.5 Cleaning/disinfection of the videoscope: Manual cleaning

Please adhere to the national guidelines, e.g., ESGENA guidelines, DEGEA/RKI recommendations, ASGE, LANGS, SAGES etc. for cleaning, disinfection and sterilization, as well as to the instructions of the chemical manufacturer concerning concentration and exposure times.

1. Remove the rubber lip valve/biopsy valve, suction valve and air/water valve and place in an active cleaning solution and clean using a soft disinfected brush (13991 VB) or disposable brush.
2. Wipe the housing, insertion shaft, and supply tube carefully with a lint-free cloth soaked in cleaner (for suitable cleaning agents or disinfectants approved by KARL STORZ, see www.karlstorz.com). Avoid shocks and pulling movements so as not to damage the insertion shaft. Then rinse the endoscope with microbiologically clean/sterile water.
3. Immerse the endoscope in the cleaning solution and rinse the channels using a syringe and the irrigation tubes supplied (13991 SS, 13991 YS, 13991 SWJ).
4. The channels must be filled without any bubbles, to guarantee a complete moistening of the surfaces.

CAUTION: *Only immerse instrument parts once they have successfully passed a leakage test!*

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



5. When submerged, use a suitable cleaning brush (preferably disposable) (see accessories, page 9/10) to brush below the surface of the liquid in one direction starting from the supply plug (supply plug --> control body, control body --> distal end, biopsy channel --> distal end) (see Fig. 13).
 Repeat this procedure at least 3 times.

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
6. If you are using reusable brushes, clean the brush head at the distal end before retracting it. Do not brush backward and forward – brush in one direction only.

CAUTION: *Cleaning with brushes must be carried out with the instrument immersed.*

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
7. Thoroughly clean out the entire suction channel (shaft part and supply part) and the biopsy channel a number of times with a brush. Observe the exposure time for the cleaning solution!

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
8. Connect the Y-irrigation tube (13991 YS, 13991 SS and 13991 SWJ) and pump fresh cleaning fluid through the channels using a syringe. Afterwards, rinse with microbiologically clean/sterile water.

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9. Carefully clean all particles of dirt from the distal end. It may be necessary to clean the air/water valve with a brush (13991 VB) with soft bristles, in order to remove hardened contaminants. Never remove debris with a needle or other sharp objects, because these can cause damage.

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
Suitable brushes

Art. no.	Designation	Used for
110917-01	Cleaning brush dia. 0.5 mm	Water-jet channel
110910-01	Cleaning brush dia. 5 mm	Working channel from 2.8 to 4.2 mm
13991 VB	Valve brush	Valve seat

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
10. Remove the device from the cleaning solution and rinse with microbiologically clean/sterile water.
11. Carefully clean all valves with the cleaning brush (13991 VB) and then rinse.
12. The valves should also be cleaned in an ultrasound bath.

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Once the cleaning work has been completed, clean the cleaning instruments, for example the brushes, of residue and disinfect/sterilize them. Never reuse damaged or kinked brushes.

CAUTION: *Do not autoclave the endoscope or its accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps, cleaning brushes, Albarran module, etc.).*

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
7.6 Manual disinfection

1. After careful rinsing, immerse the endoscope in the disinfecting solution.

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CAUTION: *When preparing and using the solutions, follow the chemical manufacturers' instructions paying close attention to proper concentration and exposure times. Prolonged immersion can lead to damage. Bear in mind the microbiologic range of action of the chemicals used.*

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2. The endoscope must be completely immersed in a disinfectant solution. Using the irrigation tubes (13991 YS, 13991 SS, 13991 SWJ) and a syringe with LUER connector, fill the lumina in particular to ensure bubble-free wetting. (see illustration).

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3. At the end of the exposure time, the channels and surfaces of the video coloscope must be rinsed with microbiologically clean/sterile water.

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4. Finally, dry off all surfaces, joints, openings, channels and lumina, preferably with sterile compressed air at a max. of 0.5 bar/7.25 psi, in accordance with the national regulations.

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CAUTION: *The pressure must not exceed 50 kPa (0.5 bar/7.25 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.*

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If at all possible, store the videoscope in a dry cupboard or hang it up in a storage case for safe storage.

NOTE: *The valves of the videoscope must not be inserted; they must be stored separately.*

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7.7 Machine reprocessing of the videoscope

All flexible KARL STORZ endoscopes can be connected to approved washer/disinfectors for semi- or fully automatic cleaning and disinfection. Observe the national laws and regulations for cleaning and disinfection procedures. Leakage testing, manual pre-cleaning and (where necessary) rinsing with microbiologically clean water must be carried out beforehand. The manufacturer of the cleaning machine must be consulted when selecting the reprocessing programs. Usually, special programs for flexible endoscopes are available which are suited for decontamination and compatible with the materials. For machine reprocessing, make sure that the flexible endoscopes are not heated above 65°C.

CAUTION: *Check the videoscope for leaks before every disinfection in semi- or fully automatic equipment.*

CAUTION: *Ensure that only KARL STORZ approved chemicals are used for procedures involving peracetic acid.*

Lay the endoscope carefully in the tray or basket of the machine. Connect the irrigation tubes supplied with the cleaning adaptor set 13991 CK or the irrigation tubes of the manufacturer of the washer/disinfector to the endoscope.

7.7.1 Cleaning adaptor set 13991 CK

The cleaning adaptor set 13991 CK is used for validated reprocessing of KARL STORZ videoscopes from series 13xxx NKS/PKS (acc. to ISO 15883 part 4). The cleaning adaptor set is not intended for any other purpose as this may damage it.

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WARNING: *This set is not sterile when delivered, and must therefore be cleaned, disinfected, and (where applicable) sterilized prior to initial use and each subsequent reuse.*

WARNING: *Risk of infection from improper reprocessing: Reprocessing may only be carried out by experienced and trained personnel.*

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CAUTION: *Only use chemicals which have been tested for material compatibility by KARL STORZ. If you use other chemicals, get written confirmation from the chemical manufacturers that the the material is compatible with KARL STORZ endoscopes.*

CAUTION: *If the reprocessing chemicals are changed, it is essential to consult the manufacturer of the chemicals and of the washer/disinfector.*

CAUTION: *Please refer to our operating instructions for details on cleaning preparation under the section: Cleaning, disinfection and sterilization.*

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The following must be checked before use:

- Visual checks of the connection tubes (patency) and the O-rings (wear)
- Check the connections are correctly and securely in place
- Check for damage


In the event of noticeable problems or complaints, please contact our service department or send in the cleaning adaptor set for technical examination.

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
Assembly cleaning adaptor set

Plug the connector ❷ of the connection tube 13991 CL into the vent port ⑮ of the videoscope at the supply connector and secure by turning clockwise.

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
Plug the cleaning adaptor 13991 CR ❺ as far as it can go into the connector ㉒ of the air pin of the videoscope.

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
Then plug connector ❻ as far as it will go into the ‘Connection for water bottle’ ⑰ of the videoscope.

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Plug connector ⑨ of the cleaning adaptor 13991CS as far as it will go into the suction channel ㉑ at the supply connector of the videoscope.

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If necessary, remove the suction and irrigation/insufflation valves (④, ⑤) from the videoscope.

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
Insert the channel separator 13991 CD ⓮ with the valve pin ⓯ into the free opening of the irrigation/insufflation valve ④ as far as it will go so that the lock ⓰ lies on top of the free opening of the suction valve ④.

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
Press the slide mechanism of the channel separator 13991 CD ⓮ downward and slide it forward until the channel separator 13991 CD ⓮ audibly locks into place. Check that the channel separator is securely in place.
Connect the adaptor 13991 SWJ to the jet channel.

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ⓘ ***NOTE:*** *If necessary, remove the cap of the biopsy/working channel ⑥ from the videoscope.*

Connect the connection adaptor 13991 SS with the connector ⓬ to the biopsy/working channel ⑥ of the videoscope.


⚠ ***WARNING:*** *Danger of infection due to inadequate rinsing caused by kinked tubes (decreases cleaning performance): avoid kinking of the irrigation tubes.*

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Connect all open LUER connectors of the attached connection adaptors to the respective connections on the washer/disinfector.

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CAUTION: *When the program has finished, the videoscopes are at a risk of corrosion due to the conditions in the reprocessing machine.*

CAUTION: *Immediately after the end of the program cycle, the videoscopes are hot as they are heated up during chemical thermal reprocessing. Consequently, the instrument's material properties are different and, as such, the prepared endoscope should be allowed to cool down for at least 20 minutes at room temperature before being used again.*

CAUTION: *Ideally, remove the videoscopes from the machine immediately after the program has finished. Then dry off the optical surfaces, if necessary, with a cotton swab soaked with 70% alcohol.*

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After reprocessing, store the endoscopes freely suspended, dust-free and dry. Do not insert the valves – store them separately.

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Trocknung

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7.8 Preliminary cleaning directly after the examination (see page 31)

- The outer casing should be wiped down from top to bottom with a lint-free disposable cloth that is soaked with active cleaning solution.
- Suction ample amounts of cleaning solution through the distal end
- Apply the air-water supply for a few seconds every cycle
- Connect the cleaning adaptor
- Separate the flexible endoscope from the light source

CAUTION: *Close the video plug with the sealing cap.*

7.9 Cleaning in the reprocessing place

- Carry out leakage test
- Immerse the flexible endoscope
- Clean the channels with special cleaning brushes
- Clean the complete air/water system thoroughly (in the case of duodenoscopes, the Albarran unit must also be disassembled)
- Rinse the channels

Disinfection

- Immerse the endoscope and rinse the channels

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7.10 Reprocessing in the washer/ disinfector

- The machine must fulfill all the cleaning, disinfection and material compatibility requirements (validated process).
- Make sure that the environmental conditions with regard to temperature and pressure are not exceeded (+60°C (tolerance +/- 5% and 1 bar for the channels, 0.3 bar leakage test).
- Recommendation: Before placing the endoscope in the machine, clean all of the channels (working channels, water channels, Albarran channels, air/water channels) regardless of whether or not they have been used.
- We recommend sticking to one type of cleaning agent.

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7.11 Cleaning of accessories

1. Thorough cleaning is required before disinfection or sterilization. Carefully clean the biopsy forceps and brushes whilst they are immersed in a cleaning solution using a soft, disinfected brush. Take special care to remove residues such as blood, excrement and slime from the biopsy forceps and needles – an ultrasound bath is recommended.
2. Rinse in microbiologically clean/sterile water and dry (using sterile compressed air).
3. Disinfect accessories (work procedures as with cleaning).
4. Always test instruments after cleaning/ disinfection for malfunction and any damage.
5. Before sterilization, carefully oil moveable parts. For this, only use suitable steam-resistant, physiologically safe, residue-free oils (e.g.,27656 B). After oiling, open and close the forceps several times. For details on reprocessing, please read the instruction leaflet of the respective accessory.

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7.13 Checking the reprocessing procedure

Before assembly, it should be checked that all parts and surfaces are undamaged, clean and dry.

Replace defective parts or send them in for repair. Before sending them in for repair, all parts must be cleaned and disinfected/sterilized or clearly labeled as unsterilized.

If residues are still present after cleaning, reprocess the parts again, give them special treatment, discard or repair.

Check plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

Test the instruments for proper operation (deflection of the distal end, brake lever). Check the lens for a clear view, the instruments for appropriate deflection and all channels for patency (see section 6.4.1).

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7.14.1 Gas sterilization

Sterilization using ethylene oxide is a gentle and therefore recommended method of sterilization.

NOTE: Observe national laws and regulations for the gas sterilization method.

CAUTION: If the instrument is to be sterilized, the entire system must be opened by fitting the EO cap on the pressure compensation part to ensure pressure compensation.

WARNING: Before commencing gas sterilization, the instrument must be thoroughly cleaned and completely dry or else sterile conditions will not be achieved.

- Place the instrument securely on a sterilization tray.

	20:80 EO/ CO	10:90 EO/ HCFC
Temperature:	55°C	55°C
Relative humidity:	50%	50%
Vacuum:	533 mmHg effective	533 mmHg effective
Outlet pressure:	69 kPa 0.70 Kg/cm, 10 PSI	97 kPa 0.98 Kg/cm, 14 PSI
EO concentration:	450 mg/l	600 mg/l
Conditioning:	1 hour	1 hour
Gas exposure time:	5 hours	5 hours
Aeration:	12 h at 55°C	12 h at 55°C

An aeration period of 72 hours is required after EO sterilization. An aeration chamber (desorption chamber) can be used to reduce the aeration period to 12 hours. During sterilization, the temperature must not exceed 55°C.

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For gas sterilization with ethylene oxide, the following airing times must be observed for instruments due to gas absorption of the materials:

	at 20°C	at 42°C	at 55°C
Telescopes	24 hrs	18 hrs	12 hrs
Rubber items	7 days	4 days	3 days
Videoscopes, plastics	5 days	2 days	30 hours

(minimum time in simple transparent packaging and an application time on the patient of more than 30 minutes).

Maximum acceptable levels (per ISO 10993-7:2008) for residues following EO sterilization are as follows:

Ethylene oxide: 25 ppm* / 250 ppm**
Ethylene chlorohydrine: 25 ppm* / 250 ppm**
Ethylene glycol: 250 ppm* / 5000 ppm**

** for instruments in contact with skin and mucous membrane,
* for instruments in contact with blood.

WARNING: *The efficiency of the sterilization should be checked to see that it complies with the German Medical Devices Act (MPG).*

CAUTION: *Attention must be paid to the notes on cleaning, sterilization and storage in the ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’ instructions. The cleaning, disinfection, and sterilization procedures are explained therein in detail.*

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10 Important information

To avoid the spread of infectious illnesses among hospital personnel or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and, if necessary, disinfected/ sterilized before they are sent for repair. We reserve the right to return contaminated instruments, units or endoscopes to the sender.

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11 Servicing and repair

Defective endoscopes must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts

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12 Limitation of liability

As a supplier of this endoscope, KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:
All assembly operations, system expansions, readjustments, modifications, or repairs to the same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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13 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see ‘Subsidiaries’ section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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14 Directive compliance

In accordance with the Medical Device Directive (MDD): This medical device belongs to Class IIa

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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REPAIR DISPATCH NOTE

Sender Information:

Hospital / Practise | Postcode / City | Date

Dept. / Contact | Tel./Fax

Product / Modell | Part-No. | Your Ref.no.


Videoendoscope is gas sterilized ❑ yes ❑ no

Videoendoscope is disinfected ❑ yes ❑ no

Notes

Damage report

Hand piece
❑ Function buttons leaking
❑ Chassis leaking

Steering wheels
❑ loose
❑ leaking

Supply tube
❑ cracked / damaged

Working channel
❑ blocked
❑ leaking


Insertion sheath
❑ cracked
❑ worn out

Distal tip
❑ Lens scraped
❑ Insufflation disturbed
❑ Image transmission error

Deflection
❑ leaking
❑ blocked / damaged

Supply plug
❑ electrical contacts damaged / corrossive
❑ leaking

Light transmission
❑ damaged
❑ Illumination power not adequate

REPAIR DISPATCH NOTE

From customer enclosed parts:

Customer	Article / Part	KARL STORZ
❑	Valve / air-water, suction valve	❑
❑	Cleaning brush	❑
❑	Protection wrap	❑
❑	Biopsy forceps	❑
❑	Cleaning tubes	❑
❑	Instrument oil	❑
❑	ETO cap	❑
❑	Bite protector	❑
❑	Leakage tester	❑
❑	Seal cap for working channel	❑
❑	Video connection cable	❑

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havanna, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada, Del. Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsätragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kęstučio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

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KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: phillippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd .
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

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Contents

1 Important information 1
2 Pictures of the instruments 2
3 Controls, displays, connectors and their uses 4
4 Symbols employed 5
5 Intended use 6
5.1 Intended use 6
5.2 Indications 7
5.3 Contraindications 7
5.4 Patient population 7
5.5 User qualification 7
6 Safety instructions 8
6.1 Safety precautions when using flexible intubation videoscopes 8
7 Initial operation 14
7.1 Unpacking 14
7.2 Basic equipment for endoscopes with a working channel 14
7.2.1 Accessories (not included in the set) 15
7.3 Basic equipment for endoscopes without a working channel 15
7.3.1 Accessories (not included in the set) 15
7.4 Operating conditions 16
7.5 External system components 16
7.6 General information on use 17
7.7 Assembly before use 18
7.7.1 Fitting the sealing cap 18
7.7.2 Fitting the suction valve 18
7.7.3 Fitting the tube holder 18
7.8 Checking before use 19
7.8.1 Checking the suction (does not apply to 11301ABX) 19
7.8.2 Checking the insufflation (does not apply to 11301ABX) 20
7.8.3 Inspecting the surfaces 21
7.8.4 Inspecting the deflection 22
7.8.5 Inspecting the image quality 23

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Contents

7.9 Connecting and switching on..... 24
7.9.1 Connection to the C-MAC® monitor 24
7.9.2 Switching the flexible intubation videoscope On/Off 26
7.9.3 Functions of the C-MAC® monitor 8403ZX 27
7.9.4 Connection to the C-HUB® II..... 27
8 Handling 29
8.1 Holding the device correctly 29
8.2 Handling the sheath..... 29
8.3 Handling the deflection mechanism 30
8.4 Handling the working channel (does not apply to 11301ABX) 32
8.5 Handling the suction valve..... 33
8.6 Function buttons of the endoscope 34
9 Reprocessing 35
9.1 General information..... 36
9.2 Precleaning directly after use 38
9.3 Leakage test 38
9.3.1 Performing the leakage test – dry..... 38
9.3.2 Performing the leakage test – in liquid..... 39
9.4 Manual cleaning and disinfection 40
9.4.1 Manual cleaning 40
9.4.2 Manual disinfection..... 41
9.4.3 Manual High Level Disinfection (HLD)..... 41
9.4.4 Drying..... 41
9.5 Machine cleaning and disinfection 42
9.5.1 Manual preparation..... 42
9.5.2 Connecting..... 42
9.5.3 Machine cleaning/chemical-thermal disinfection..... 42
9.6 Assembly, inspection and care 43
9.7 Packaging 44

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11 Servicing and repair ... 60
12 Maintenance ... 60
13 Limitation of liability ... 60
14 Warranty ... 60
15 Directive compliance ... 61

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This instruction manual is intended to serve as an aid in the proper use, cleaning, disinfection and, if necessary, sterilization of the intubation videoscope. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

⚠ WARNING:

Risk of infection:

KARL STORZ instruments are not sterile when delivered, and must be cleaned, disinfected and, if necessary, sterilized prior to initial use and each subsequent reuse.

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3 Controls, displays, connectors and their uses

① Connecting cable
② Suction valve
③ Connection for leakage tester or pressure compensation cap (pressure compensation valve)
④ Connection for suction valve/irrigation adaptor 11301CB1/11301CE1/11301CD1/ 091011-01 (only for endoscopes with working channel)
⑤ Working channel with lip valve 29100 (only for endoscopes with working channel)
⑥ Tube holder
⑦ Sheath
⑧ Distal tip
⑨ Deflection lever
⑩ Video recording Start/Stop
⑪ White balance
⑫ Single image capture
⑬ Monitor 8403ZX
⑭ Second connection for videoscope
⑮ Primary connection for videoscope
⑯ Mains connection
⑰ Reset button
⑱ Holder for SD memory card
⑲ Status bar
⑳ Pilot lamp
㉑ On/Off switch button
㉒ Charging pilot lamp (orange)
㉓ Softkeys
㉔ Function symbols
㉕ HDMI port
㉖ USB port
㉗ Leakage tester 13242XL
㉘ Pressure compensation cap 11025E

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5 Intended use

5.1 Intended use

Intubation endoscopes, flexible, CMOS with working channel

Flexible intubation endoscopes with working channel (CMOS) are designed for both endotracheal and nasal intubation. Flexible intubation endoscopes with working channel (CMOS) are designed for endoscopic visualization and image transmission as well as for creating an access to introduce additional instruments for diagnosis and treatment of the upper and lower airways. Flexible intubation endoscopes with working channel (CMOS) are intended for temporary use in interventions performed through natural orifices.

Intubation endoscopes, flexible, CMOS without working channel

Flexible intubation endoscopes without working channel (CMOS) are designed for both endotracheal and nasal intubation. Flexible intubation endoscopes without working channel (CMOS) are designed for endoscopic visualization and image transmission for diagnosis of the upper and lower airways. Flexible intubation endoscopes without working channel (CMOS) are intended for temporary use in interventions performed through natural orifices.

NOTE: *The flexible intubation videoscope 11301BNX can also be used with the C-MAC® monitor 8402ZX. All other endoscopes in the FIVE series are only compatible with the C-MAC monitor 8403ZX.*

Flexible intubation videoscopes are applied parts, which must only be used with KARL STORZ supply units which are intended for this purpose.

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5.2 Indications

The use of intubation endoscopes is indicated if, in the opinion of the attending physician, a fiberoptic endotracheal intubation or inspection of the oropharyngeal, supra- or subglottic region is indicated or PCT is performed.

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- Endotracheal intubation with a difficult respiratory tract
- Endoscopic inspection and monitoring of the mouth, nose and throat region.
- Monitoring of percutaneous tracheostomy
- Anesthesiological bronchoscopy

5.3 Contraindications

The use of the intubation endoscope is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Intubation endoscopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

No other contraindications are known.

5.4 Patient population

Men and women as well as children and adolescents. Use is not limited to a certain patient profile (sex, age, weight, etc.).

NOTE: *The attending physician decides which patients to use the product on.*

5.5 User qualification

KARL STORZ flexible intubation videoscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the flexible laryngoscopy technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the flexible intubation videoscope. It is not intended as an introduction to the techniques of flexible laryngoscopy.

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6 Safety instructions

Please read this manual carefully and follow the instructions exactly. The words **Warning, Caution** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this instrument. To further emphasize the words **Warning, Caution** and **Note**, they are preceded by a pictogram.

⚠ WARNING:

‘Warning’ indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

⚠ ATTENTION:

‘Attention’ indicates that particular service procedures or safety precautions must be followed to avoid any damage to the product.

NOTE:

‘Notes’ indicate special information about operating the instrument, or clarify important information.

6.1 Safety precautions when using flexible intubation videoscopes

Flexible intubation videoscopes must only be used according to the medical rules and procedures of endoscopy recognized for this method.

NOTE: *The flexible intubation videoscopes do not offer any electric insulation against high frequency (HF) currents. The use of HF instrument sets can result in image faults and damage to the endoscope.*

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Safety instructions

⚠ ATTENTION:

Damage to product: Whenever the flexible intubation videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition.

➢ *Damaged flexible intubation videoscopes or damaged accessories must not be used.*

⚠ WARNING:

Combinations of medical devices are only assured to be safe if

- *they are identified as such in the respective instruction manuals or*
- *the intended use and interface specifications of the devices used in combination permit this.*

The instruction manuals and interface specifications of the medical devices used in combination must be observed precisely.

⚠ WARNING:

Risk of infection:

➢ *The flexible intubation videoscopes are not sterile when delivered, and must be cleaned and disinfected/sterilized prior to initial use and each subsequent reuse.*

⚠ WARNING:

The pressure compensation cap ㉘ must be removed from the pressure compensation valve ③ during:

- *endoscopic interventions*
- *immersion in liquid, e.g., disinfectant solution and for machine reprocessing*

⚠ WARNING:

Risk from electromagnetic incompatibility: Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Observe the stated EMC instructions during installation and operation.

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Safety instructions

⚠ ATTENTION:

Damage to product: The pressure compensation cap ㉘ must be placed on the pressure compensation valve ③ during:

- *gas sterilization (EO, FO, hydrogen peroxide)*
- *transport and shipping*

⚠ WARNING:

Risk of burns and glare: The LED light emitted through the flexible intubation videoscope may cause glare or high temperatures, and thereby burns, at the distal light outlet.

➢ *Never look directly into the light emitted by the flexible intubation videoscope.*

⚠ WARNING:

Risk of explosion: If flammable anesthetics are used in the immediate proximity of the flexible intubation videoscope.

➢ *The flexible intubation videoscope must not be used in the presence of flammable or combustible media.*

⚠ ATTENTION:

Only use accessories that are specifically suited for use with videoscopes in accordance with the intended use and which can be employed in line with generally accepted technical standards and safety guidelines.

⚠ ATTENTION:

Damage to product: The flexible intubation videoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for flexible intubation videoscopes that are repaired by unauthorized personnel.

NOTE: *Follow the 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions. The cleaning, disinfection and sterilization procedures are explained therein in detail.*

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NOTE: *Damage to the flexible intubation videoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.*

⚠ WARNING:

Always unplug the flexible intubation videoscope before carrying out any maintenance and cleaning procedures on it.

⚠ WARNING:

Risk of injury: The flexible intubation videoscope must not be used during the discharge of a cardiac defibrillator.

➢ *Remove the flexible intubation videoscope from the operation area before discharge takes place.*

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⚠ WARNING:

Risk of infection: In order to prevent infection and protect personnel, it is strictly forbidden to ship contaminated medical devices.

➢ *To avoid contact and aerogen infections, medical devices must be decontaminated on-site and wrapped in protective film.*

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⚠ WARNING:

Danger of electric shock:

➢ *Have servicing on these devices carried out by authorized personnel only. Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims. Exceptions are those that are described in the instruction manual.*

⚠ WARNING:

Only KARL STORZ devices may be used in combination with the flexible intubation videoscope. Disregarding this warning can endanger the patient and user.

⚠ WARNING:

Risk of injury: If the flexible intubation videoscope malfunctions during use on the patient, discontinue use immediately.

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⚠ ATTENTION:

Damage to product:

➢ *Only use accessories (guide probes, etc.) from KARL STORZ. Accessories from other manufacturers can damage the flexible intubation videoscope.*

⚠ ATTENTION:

Damage to product:

➢ *Flexible intubation videoscopes must never be steam sterilized (autoclaved).*

⚠ ATTENTION:

Damage to product:

➢ *The maximum permissible temperature for decontamination, sterilization, transport and storage is 65°C.*

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⚠ ATTENTION:

Damage to product:

➢ *Do not allow any live components or instruments to come into contact with the flexible intubation videoscope.*

⚠ ATTENTION:

Damage to product:

➢ *Also do not place the flexible intubation videoscope on other devices that become very hot.*

➢ *Do not expose the device to direct sunlight or excessive heat, to avoid damage to the device.*

⚠ ATTENTION:

Damage to product:

➢ *Flexible intubation videoscopes must never be exposed to X-rays, to avoid damage.*

⚠ ATTENTION:

Damage to product:

➢ *Do not hit the distal tip of the flexible intubation videoscope against hard objects.*

➢ *Do not remove soiling with sharp-edged objects under any circumstances.*

⚠ ATTENTION:

Damage to product:

➢ *Never lay heavy objects on the flexible intubation videoscope.*

⚠ ATTENTION:

Damage to product:

➢ *Never clamp or fold the sheath of the flexible intubation videoscope. This can result in damage to the videoscope.*

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⚠ WARNING:

Risk of infection:

➢ *The flexible intubation videoscope is not sterile when delivered, and must therefore be properly cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.*

Check for missing items and evidence of shipping damage.

File any complaints immediately with KARL STORZ or the supplier.

7.2 Basic equipment for endoscopes with a working channel

The 11302BDXK, 11301BNXK, 11303BNXK and 11304BCXK sets consist of:

1x Flexible intubation videoscope
11302BDX 4.0 x 65 cm or
11301BNX 5.5 x 65 cm or
11303BNX 5.5 x 65 cm or
11304BCX 6.5 x 65 cm
1x Case
1x Instruction manual
1x Leakage tester 13242XL
1x Pressure compensation cap 11025E
1x Suction valve insert for single use, pack of 20, 11301CE1
1x Cap for LUER lock irrigation port, black, pack of 10, autoclavable, 29100
2x Adaptors for cleaning the working channel in cleaning machines (black), 11301CD1
50x Cleaning brush
1x Tube holder for video endoscopes 11301CFX
1x Bronchoscope insertion tube, size 4, 10309
1x Bronchoscope insertion tube, size 2, 10310
1x Protective film 13990 SFN
1x Protective cap 8403YZ

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7.2.1 Accessories (not included in the set)

1x Suction valve insert, reusable 11301CB1
1x Adaptor for leakage test in washers and disinfectors (green), 13242XA
1x ‘Mainz Adaptor’ mask attachment, for endoscopic inspection with simultaneous respiration, sterile packaging, pack of 5, 11008C
1x Lip valve, non-sterile, for single use, pack of 20, 11301 CA
1x Biopsy Forceps with Oval Jaws, coated, working length 120 cm, for single use, sterile, package of 10, 110284-10

7.3 Basic equipment for endoscopes without a working channel

The 11301ABXK set consists of:
1x 11301ABXK flexible intubation videoscope 3.0 x 51.5 cm 11301ABX
1x Case
1x Instruction manual
1x Leakage tester 13242XL
1x Pressure compensation cap 11025E
1x Tube holder for video endoscopes 11301CFX
1x Bronchoscope insertion tube, size 4, 10309
1x Bronchoscope insertion tube, size 2, 10310
1x Protective film 13990 SFN
1x Protective cap 8403YZ

7.3.1 Accessories (not included in the set)

1x Adaptor for leakage test in washers and disinfectors (green), 13242XA

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7.6 General information on use

The flexible intubation videoscope must be used according to the medical rules and procedures recognized for laryngoscopy.

NOTE: *Any damage to the flexible intubation videoscope resulting from incorrect operation is not covered by the manufacturer's warranty.*

⚠ WARNING:

Risk of infection: The flexible intubation videoscope is not sterile when delivered, and must therefore be properly reprocessed prior to initial use and each subsequent reuse.

⚠ ATTENTION:

Damage to product:

➢ *In order to avoid damage to the exterior, do not store or disinfect the flexible intubation videoscope together with sharp objects.*

The following instructions outline the correct operation of the flexible intubation videoscope and the corresponding accessories. This is not an introduction to the techniques of laryngoscopy. Instructions on endoscopic technique can be found in the relevant medical literature.

The flexible intubation videoscope is used for endotracheal intubation and inspection of the oropharynx.

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7.7 Assembly before use

7.7.1 Fitting the sealing cap

Place the sealing cap (Art. no. 29100) on the proximal end of the working channel.

7.7.2 Fitting the suction valve

NOTE: *The suction valve only needs to be fitted if the flexible videoscope has a working channel.*

1. Place the suction valve (Art. no. 11301CE1) onto the valve connecting piece in such a way that the tube connection is positioned above the function buttons and the unit's side nubs engage with the recesses on the connecting piece.

2. Using gentle pressure, rotate the suction valve in a clockwise or counter-clockwise direction until it noticeably engages.

For additional information on handling the suction valve, see the instruction manual for the suction valve 97000282.

7.7.3 Fitting the tube holder

Push the tube holder ⑥ 11301CFX from distal onto the cone of the flexible intubation videoscope until the tube holder is securely in place.

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7.8 Checking before use

NOTE: *Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding precleaning.*

⚠ ATTENTION:

Damage to product: Damaged flexible intubation videoscopes or damaged accessories must not be used.

➢ *Whenever the flexible intubation videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition.*

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7.8.1 Checking the suction (does not apply to 11301ABX)

Suction can be triggered by pressing the suction valve ② (when the suction pump is connected and switched on).

NOTE: *There should be no instruments in the suction channel whilst suctioning is being performed. The instrument channel inlet should be sealed (Art. no. 29100) to prevent additional air from being drawn in.*

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7.8.2 Checking the insufflation (does not apply to 11301ABX)

Connect the O_2 tube to the instrument channel:

- Use either the olive connector (Art. no. 600007) on the instrument channel or,
- if you are using the three-way double stopcock with LUER lock (Art. no. 6927691), connect the tube to the tube holder on the angled part of the adaptor. This keeps the straight opening of the adaptor free for instrumentation.

⚠ WARNING:

Observe the instruction manuals of the supply units used.

➢ *Switch on the respective supply units.*

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⚠ ATTENTION:

Damage to product:

➢ *The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.*

⚠ ATTENTION:

Damage to product:

➢ *Do not squeeze the distal tip of the flexible intubation videoscope too tightly.*

Every time the flexible intubation videoscope is used, you must first check that its surface, and above all the surface of the sheath and distal tip, is in perfect condition and that there are no sharp edges. Intubation videoscopes with sharp edges must not be used.

They may cause serious injury to the patient or user.

Also check the connecting cable for damage (such as fractures, cracks, twisted or crushed areas, etc.) and the connecting plug for loose connections.

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7.8.5 Inspecting the image quality

⚠ WARNING:

Risk of injury: If the image fails or there is interference during use, this can result in injury to the patient.

- ➢ *Put the distal tip into the neutral position and remove it from the patient's body gently and carefully.*
- ➢ *If the intubation flexible videoscope continues to be used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.*

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7.9 Connecting and switching on

The flexible intubation videoscopes (FIVE) can be operated with the following monitors:

- C-MAC® monitor 8403ZX (11301BNX also with C-MAC® monitor 8402ZX)
- also with additional monitors via the camera control unit C-HUB® II 20290301

7.9.1 Connection to the C-MAC® monitor

The flexible intubation videoscopes 11301ABX, 11302BDX, 11303BNX and 11304BCX can be used together with the C-MAC® monitor 8403ZX. The instruction manual for the monitor must be observed (see manual 96076008DER).

NOTE: *Older monitor models (8402ZX/8401ZX) are not compatible.*

The flexible intubation videoscope 11301BNX can be connected to both the C-MAC® monitor 8402ZX and the C-MAC® monitor 8403ZX.

⚠ WARNING:

Risk of injury: It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other devices at any time.

⚠ WARNING:

Risk of injury: Optical radiation poses a risk of injury to eyes.

➢ *Never look directly into the light output at the distal tip.*

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Initial operation

Plug the connecting cable ① of the flexible intubation videoscope into one of the connection sockets ⑭ or ⑮ on the monitor.

The connecting cable ① can be connected or disconnected during operation. If necessary, connect the cord of the power supply up to the connection socket ⑯ on the C-MAC® monitor and plug it into the power supply. The power cord can be fixed to the serial port using the union nut.

The battery is charged automatically as soon as the C-MAC® monitor is connected to the power supply. This is indicated by the lighting up of the orange LED ㉒.

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7.9.2 Switching the flexible intubation videoscope On/Off

Switching On/Off: Pressing the on/off button ㉑ on the C-MAC® monitor switches on the flexible intubation videoscope with the integrated LED light source. Pressing the button a second time switches the unit off.

⚠ WARNING:

Danger of glare:

➢ *Do not look directly into the light beam.*

⚠ WARNING:

Risk of injury: When in operation, the integrated LED light source heats up the flexible intubation videoscope.

➢ *Switch off the flexible intubation videoscope before putting it down, and do not place it in the vicinity of the patient (as per IEC 60601-1-1).*

***NOTE:** The C-MAC® monitor can be switched on or off by pressing the On/Off button ㉑. The flexible intubation video endoscope itself does not have an On/Off button.*

Pressing and holding the On/Off button switches off the monitor together with the connected endoscope. The green indicator light ⑳ goes out.

If the On/Off button is actuated only briefly during operation, the monitor will switch to standby mode. This is indicated by the flashing indicator light on the monitor.

***NOTE:** When the monitor is switched off, the integrated light source in the flexible intubation videoscope also switches off.*

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7.9.3 Functions of the C-MAC® monitor 8403ZX

A more detailed description of the operation and functions of the monitor 8403ZX can be found in the corresponding instruction manual supplied with every monitor (see manual 96076008DER).

If you lose the instruction manual, you can also download it from our homepage at www.karlstorz.com.

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7.9.4 Connection to the C-HUB® II

The C-HUB® II (Art. no. 20290320) is a plug & play device for processing and transmitting image and control signals to external devices in medical video endoscopy. The flexible intubation videoscope can be used in combination with the C-HUB® II and an external monitor.

To do so, insert the connector of the flexible intubation videoscope into the serial port of the C-HUB® II.

A more detailed description of the operation and functions of the C-HUB® II can be found in the corresponding manual (manual 96206529DER).

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Initial operation

NOTE: *The flexible intubation videoscopes 11301ABX, 11302BDX, 11303BNX and 11304 BCX can only be connected to the C-HUB® II 20290301. Older models are not compatible. The flexible videoscope 11301BNX can be used with both the C-HUB® (Art. no. 20290101) and the C-HUB® II (Art. no. 20290320).*

NOTE: *Observe the instruction manual for the C-HUB® II (Art. no. 96206529DER).*

1. Connect the power supply unit to the C-HUB® II.
2. Connect the video connecting cable from the flexible intubation videoscope to the video input on the front of the C-HUB® II.
3. Connect the monitor to the S-Video or HDMI output on the rear of the C-HUB® II.

Focus

The image focus is adjusted automatically. Manual adjustment is not necessary.

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8 Handling

8.1 Holding the device correctly

Take the handle of the flexible intubation videoscope in your right or left hand and wrap your fingers around it to secure, as shown in the Figure. The thumb is located on the deflection lever ⑨. With this kind of control grip, the flexible intubation videoscope can be operated with one hand. (Fig. 1 Intubation videoscope with working channel; Fig. 2 Intubation videoscope without working channel).

8.2 Handling the sheath

⚠ ATTENTION:

Damage to product:

➢ *The flexible intubation videoscope must not be carried by the sheath. Do not pull, clamp or twist the sheath of the flexible intubation videoscope.*

⚠ ATTENTION:

Damage to product:

➢ *The sheath of the flexible intubation videoscope must never be clamped or coiled up too tightly. It is best to store the flexible intubation videoscope in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39406AS).*

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8.3 Handling the deflection mechanism

The deflection lever ⑨ can be used to angle the distal tip of the flexible intubation videoscope. To do so, the deflection lever needs to be moved upwards (tip angled downwards) or downwards (tip angled upwards) using the thumb. When the distal tip is moved back into the neutral position, there will be a click. This ensures that the user is aware that the distal tip is back in the neutral position.

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⚠ WARNING:

Risk of injury: When inserting the flexible intubation videoscope into the patient, the distal tip of the flexible intubation videoscope must be in the straight, non-deflected position.

⚠ WARNING:

Risk of injury: The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device!

- ➢ *Actuate the deflection lever gently and carefully as far as the limit stop.*
- ➢ *Move the deflection lever smoothly and slowly!*

⚠ WARNING:

Risk of injury:

- ➢ *Before the endotracheal tube is pushed via the sheath of the flexible intubation videoscope, the distal tip must be in the neutral position as the endoscope could otherwise be damaged.*

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⚠ WARNING:

Risk of injury:

➢ *Before removing the flexible intubation videoscope from the patient, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.*

⚠ WARNING:

Risk of injury:

➢ *If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible intubation videoscope must not be used.*

⚠ WARNING:

Risk of injury:

➢ *If the flexible intubation videoscope malfunctions during use on the patient, discontinue use immediately.*

➢ *Place the distal tip of the flexible intubation videoscope in the neutral, non-deflected position and remove the flexible intubation videoscope slowly and carefully from the patient.*

⚠ ATTENTION:

Damage to product: The resistance that can be felt when deflecting the distal tip using the deflection lever increases as the angle of deflection increases.

➢ *Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.*

⚠ ATTENTION:

Damage to product:

➢ *The lens at the distal tip of the flexible intubation videoscope may become detached, e.g., due to mechanical forces. For this reason, it is essential to check that the lens is undamaged and present before and after every decontamination.*

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8.4 Handling the working channel (does not apply to 11301ABX)

Do not force instruments or cleaning brushes through the channels. This can lead to damage to the working channel, especially if the tip of the sheath is deflected.

NOTE: *Check instruments, cleaning brushes and other accessories for compatibility before inserting these in the working channel.*

⚠ ATTENTION:

Damage to product: Unsuitable or damaged instruments may damage the working channel or cause the endoscope to malfunction.

➢ *Check instruments, cleaning brushes and other accessories for sharp edges and corners before inserting them in the working channel.*

When connecting and operating suction/irrigation devices, pay attention to the instruction manuals for the particular devices. The irrigation devices recommended by KARL STORZ are set to a maximum pressure of 54 kpa (400 mmHg/0.54 bar/7.8 psi).

- If no instrument is inserted, seal the proximal end of the working channel with a plug (Art. no. 29100).
- If an instrument is to be inserted, attach the lip valve (Art. no. 11301CA) for single use to seal the proximal end of the instrument channel.
- The working channel of the flexible intubation videoscope means that specially developed instruments can be introduced. Make sure that the accessories are of the appropriate length and diameter. Unsuitable accessories can damage fiberscopes and flexible intubation videoscopes. If you have any questions concerning accessories, please contact KARL STORZ directly or an authorized specialist dealer.

NOTE: *When using the lip valve, the instrument is subjected to increased resistance.*

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- The flexible intubation videoscope must be in a straight, undeflected position when inserting accessories into the working channel.

⚠ ATTENTION:

Damage to product: The resistance that can be felt when inserting the instruments increases as the angle of deflection increases.

➢ *If the instruments are inserted with force, this can damage the flexible intubation videoscope or the instrument channel.*

Forceps, scissors and other accessories whose jaw can be opened or closed must always be inserted into the instrument channel with closed jaws. The jaws may only be opened when the distal tip of the instrument is completely clear of the channel. Opening the jaws in the instrument channel can cause damage to the flexible endoscope.

⚠ ATTENTION:

Damage to product: Accessories from other manufacturers may damage the flexible endoscope.

➢ *Only use accessories (guide probes, etc.) from KARL STORZ.*

8.5 Handling the suction valve

When the suction valve is assembled, a suction pump must be connected to the suction valve via a tube. Suction is activated upon pressing the valve button. If the suction valve is released, the suction is interrupted.

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8.6 Function buttons of the endoscope

The endoscope is equipped with three buttons. These are used to activate the functions displayed as symbols on the endoscope.

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Video recording ⑩:
Starts or stops the recording of a video stream on the SD memory card inserted in the receptacle ⑱ of the C-MAC® monitor.

The video recording symbol appears in the top left-hand corner of the live image during video recording.

White balance ⑪:
Enables the user to carry out a manual white balance. To do this, point the tip of the videoscope at a white surface and press the button for 1 second.

NOTE: *With the flexible intubation videoscopes, the performance of the white balance procedure is visualized by a progress bar in the menu bar. The white balance can also be performed via the C-MAC® monitor menu.*

Single image capture ⑫:
Saves the live image on the screen to the SD memory card of the C-MAC® monitor.

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9 Reprocessing

NOTE: *Remove the connector from the C-MAC® monitor or the C-HUB® II. Remove the tube holder, suction valves, lip valves and other accessories from the endoscope, and reprocess them separately.*

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- ● This reprocessing step is part of the validated process.
- O This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility; however, it must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

NOTE: *Empty fields are not part of the validated process.*

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⚠ WARNING:

Risk of infection:

➢ *Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.*

⚠ WARNING:

Risk of infection from non-sterile delivery.

➢ *Reprocess the device before each use. Check the medical device for visible soiling before use.*

⚠ WARNING:

Damage to product:

➢ *Before reusing the device, replace any defective parts or send the flexible endoscope for repair.*

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⚠ ATTENTION:

Damage to product:

➢ *Never clean flexible endoscopes with metal brushes.*

⚠ ATTENTION:

Damage to product: Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.

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⚠ ATTENTION:

Damage to product:

➢ *Flexible endoscopes must never be steam sterilized (autoclaved).*

⚠ ATTENTION:

Damage to product:

➢ *To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.*

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NOTE: *We recommend using a plastic basin for disinfection to avoid damaging the endoscope. Furthermore, use of a plastic basin prevents electrolytic corrosion which can occur if there are different metals in a solution.*

NOTE: *The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.*

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9.3.2 Performing the leakage test – in liquid

NOTE: *Damage to product: To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.*

NOTE: *Any lumina must be covered with no air bubbles by filling in a targeted manner.*

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Remove the endoscope from the liquid and press the vent button on the manometer. The test pressure must drop to 0 before the leakage tester can be removed.

NOTE: *Never release the pressure in the liquid.*

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⚠ ATTENTION:

Damage to product: If a pressure loss occurs, this means that the endoscope leaks and may no longer be used, immersed in liquid or reprocessed by machine.

➢ *In this case, fill out the accompanying repairs form and send the endoscope, which has been cleaned by wiping to remove organic soiling/residue, for repair or, alternatively, contact KARL STORZ. The endoscope must not be labeled as disinfected. Line the storage/carrying case with the protective film 13990 SFN.*

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⚠ ATTENTION:

Damage to product:

➢ *Never remove dirt with a needle or other sharp objects, since these can cause severe damage.*

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9.4.2 Manual disinfection

The endoscope must be completely immersed in a disinfectant solution. The working channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time (as specified by the manufacturer), the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

9.4.3 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time (as specified by the manufacturer), the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

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9.4.4 Drying

Finally, the surfaces and channels/lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

⚠ ATTENTION:

Damage to product:

➢ *The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.*

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9.5 Machine cleaning and disinfection

The following method for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’:

9.5.1 Manual preparation

The endoscope must be completely immersed in a cleaning solution. The working channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channel is brushed until the brush head is visibly clean. Never brush back and forth within the channel.

Finally, it must be rinsed with cold water to ensure neutralization.

9.5.2 Connecting

In order to ensure effective machine cleaning and disinfection, the instrument must be connected to the suction and working channel of the endoscope in the washer and disinfector.

⚠ ATTENTION:

Damage to product:

➢ *During cleaning and disinfection the pressure compensation cap must not be in place on the pressure compensation valve.*

9.5.3 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C and pressure of 100 kPa (1 bar/14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer’s instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: *If necessary, the instrument must be dried off manually afterwards.*

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- Damaged or corroded medical devices must be withdrawn from use.

NOTE: *Inspect plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.*

- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surfaces (e.g., connection) and channels/lumina must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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➢ *Afterwards, a functional check must be carried out.*

NOTE: *Check for a clear view and steerability.*

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⚠ WARNING:

Risk of infection:

➢ *Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.*

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⚠ ATTENTION:

Damage to product:

➢ *The pressure compensation cap must be in place on the pressure compensation valve for sterilization and transport to ensure pressure compensation.*

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The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’. The procedure must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

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9.8.1 Hydrogen peroxide (H_2O_2) sterilization – ASP STERRAD®

⚠ ATTENTION:

Damage to product: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

NOTE: *For detailed information, please consult the user handbook of the respective device.*

NOTE: *Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.*

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- *STERRAD®* 100S Long Cycle (without channel)
- *STERRAD®* 100S Long Cycle with Boosters (with channel)
- *STERRAD®* NX® Advanced Cycle
- *STERRAD®* 100NX® Flex Cycle
- *STERRAD®* 100NX® DUO Cycle

9.8.2 Hydrogen peroxide (H_2O_2) sterilization – STERIS® V-PRO™

NOTE: *For detailed information, please consult the user handbook of the respective device.*

NOTE: *The 'STERIS Low Temperature Reprocessing Device Compatibility Matrix' can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.*

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⚠ WARNING:

Risk of infection:

➢ *Improperly reprocessed medical devices may pose a risk of infection for patients, users and third parties and may cause malfunction of the medical device. Observe the manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’ and the accompanying documentation.*

⚠ WARNING:

Risk of infection:

➢ *These medical devices are not supplied in sterile condition. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.*

⚠ WARNING:

Risk of infection:

➢ *When carrying out any work on contaminated medical devices, the guidelines of the Employers’ Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.*

⚠ ATTENTION:

Damage to product:

➢ *When preparing and using solutions, follow the chemical manufacturer’s specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.*

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⚠ ATTENTION:

Damage to product:

➢ *Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.*

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Brushes: 27652, 27650 B

9.11.3 Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer’s specifications must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth.

Preparation for cleaning and disinfection

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9.11.4 Manual precleaning

Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible, and/or must be opened.

Brushing the surfaces

Any visible contamination or heavy soiling must be removed by cleaning the surfaces under cold running water with the aid of a brush or sponge.

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9.11.5 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual, ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’ (Art. no. 96216003):

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9.11.6 Connecting

In order to ensure effective machine cleaning, disinfection and rinsing, if necessary, the instrument must be connected up to the washer and disinfector.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Disassembled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

NOTE: *During care procedures, use items from the catalog ‘Care, Sterilization and Storage Techniques’.*

9.11.8 Packaging systems

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

9.11.9 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’ (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

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Steam sterilization in the fractionated prevacuum procedure

NOTE: *Sterilize greased components in a disassembled state.*

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C – 137°C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

NOTE: *This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.*

Or:

- 134°C – 137°C with a minimum exposure time of 3 minutes and maximum of 18 minutes.

Hydrogen peroxide (H_2O_2) sterilization – ASP STERRAD®

⚠ WARNING:

Risk of infection:

➢ *Medical devices must be sterilized in dismantled condition.*

NOTE: *Refer to the '*STERRAD® *Sterility Guide' to ensure that the medical device concerned can be sterilized using the different* STERRAD® *devices.*

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The following STERRAD® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- STERRAD® 50, 100S, 200 without booster
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

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10 Technical data

10.1 Flexible intubation videoscope 11304BCX

Direction of view:	0°
Angle of view:	100°
Deflection:	up 180°, down 140°
Working length:	65 cm
Outer diameter distal tip:	6.5 mm
Working channel:	3.2 x 2.8 mm
Illumination:	Integrated LED light source
Electrical power consumption:	1 W
Sensor:	CMOS technology
Weight	425 g
Drip water protection:	IPX8
Classification:	IIa

Applied part: Provides Type BF protection when used in conjunction with a KARL STORZ camera control unit with a BF symbol on the connection port.

Operating/storage conditions

Operating conditions:	15°C ... +40°C 15–80% rel. humidity
Storage conditions:	-0°C ... +55°C 15–90% rel. humidity
Atmospheric pressure:	700–1080 hPa

⚠ ATTENTION:

The specifications in relation to atmospheric pressure can only be guaranteed if the pressure compensation cap is fitted to the pressure compensation valve during transportation.

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15 Directive compliance

In accordance with the Medical Devices Directive (MDD):

11302BDX, 11301BNX, 11303BNX and 11304BCX are Class II a medical devices.

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16 Electromagnetic Compatibility

The applied part described in this instruction manual has been tested as a system with the following basic devices. For relevant electromagnetic compatibility (EMC) information, refer to ‘Electromagnetic Compatibility (EMC) information’ in the instruction manuals for these basic devices:

Basic device	Instruction manual
C-MAC® monitor 8403ZX	96076008DER
C-HUB® II 20290320	96206529DER

The EMC warning statements, precautions, notes and Emission/Immunity limits specified in the instruction manuals for the basic devices also apply to the applied part described in this instruction manual.

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18 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the device as electronic scrap. Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your authorized dealer about the collection point in your area. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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29100 Sealing cap,
black sealing cap, for LUER lock irrigation port, black, pack of 10 autoclavable

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
11301CD1 Irrigation adaptor,
for machine cleaning, reusable, for use with flexible intubation videoscopes 11301BNX, 11302BDX, 11303BNX and 11304BCX

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
11301CE1 Suction valve,
for single use, pack of 20, for use with flexible intubation videoscopes 11301BNX, 11302BDX, 11303BNX and 11304BCX

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091011-20 Suction Valve,
sterile, for single use, package of 20, for use with flexible Intubation Video Endoscopes 11301BNX, 11303BNX,11302BDX, 11304BCX and flexible Video Cystoscopes 11272VH/VHU

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10309/10310 Bronchoscope insertion tube,
size 4, with integrated mouthpiece, for single use, sterile, insertion length 85 mm, made from EVA, pack of 10, for use with flexible intubation endoscopes

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
11025E Pressure compensation cap,
for venting during gas and plasma sterilization

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13242XL Leakage tester,
with blower and pressure gauge

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8403YZ Protection Cap,
for the C-MAC® System Interface

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100010-10 Endoscopic Seal,
for working channels for 3–7 Fr. instruments, for single use, sterile, package of 10

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11008C ‘Mainz Adaptor’ mask attachment,
for endoscopic inspection with simultaneous respiration, sterile packaging, pack of 5

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
110274-10 Biopsy Forceps,
with Oval Jaws, working length 120 cm, sterile, for single use, package of 10, for use with flexible endoscopes with working channels as of diameter 2.0 mm

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110282-10 Biopsy Forceps,
with Oval Jaws, coated, working length 180 cm, sterile, for single use, package of 10, for use with flexible endoscopes with working channels as of diameter 2.8 mm

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110284-10 Biopsy Forceps,
with Oval Jaws, coated, working length 120 cm, sterile, for single use, package of 10, for use with flexible endoscopes with working channels as of diameter 2.0 mm

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110345-10 Grasping Forceps Type Pelican,
with large opening angle, diameter 2.4 mm, working length 170 cm, sterile, for single use, package of 10, for use with flexible endoscopes with working channels as of diameter 2.8 mm

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110352-10 Grasping Forceps Type Alligator,
with ribbing, diameter 1.8 mm, working length 120 cm, sterile, for single use, package of 10, for use with flexible endoscopes with working channels as of diameter 2.0 mm

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110354-10 Grasping Forceps Type Alligator,
with ribbing, diameter 2.4 mm, working length 160 cm, sterile, for single use, package of 10, for use with flexible endoscopes with working channels as of diameter 2.8 mm

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INSTRUCTION MANUAL
SILVER SCOPE® Video gastroscopes,
Models 1382xxx

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, these products are the result of years of experience and great care in manufacturing. You and your organization have decided in favor of modern, high-quality items of equipment from KARL STORZ.

This manual is intended to serve as an aid in the proper use, cleaning, disinfection and sterilization of videoscopes manufactured by KARL STORZ. The videoscopes are sensitive precision instruments, which require particular care during preparation. To prevent damage during preparation, guidance is provided here for the disassembly, cleaning, sterilization and assembly of KARL STORZ videoscopes. For this reason, read these instructions carefully and keep them as a ready reference.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply, design, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, illustrations and descriptions in this manual.

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Contents

1 Symbols employed3
2 Images Gastroscope6
3 Images Accessories8
4 Safety instructions....10
5 Intended use10
6 User qualification....11
7 Warnings and cautions11
8 Safety precautions when using the videoscopes12
9 Handling....17
9.1 Compatible video processor units/light sources....17
9.2 Connecting to a light source, for example the KARL STORZ XENON 100 SCB....18
9.3 Operating the videoscope....19
9.3.1 Preparation of the light source....19
9.3.2 Preparation of the processor unit19
9.3.3 Connecting the light conductor, air/water connections and the suction channel....19
9.3.4 Deflection of the distal tip20
9.3.5 Suction20
9.3.6 Insufflation and light21
9.3.7 Water irrigation for lens cleaning21
9.4 Preparation before use21
9.4.1 Testing the various functions....22
9.4.2 Handling the endoscope correctly25
9.4.3 Additional recommendations for use26
9.4.4 Operating the endoscope buttons27
10 Reprocessing28
10.1 General information....28
10.2 General cleaning....30
10.2.1 Preparation for cleaning....31
10.2.2 Leakage test....32
10.3 Channel system of a KARL STORZ SILVER SCOPE® video gastroscope....35
10.4 Channel system of a KARL STORZ SILVER SCOPE® Interventional video gastroscope 13826 PKS/NKS....36

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Contents

10.5 Cleaning/disinfection of the videoscope: manual cleaning....37
10.6 Manual disinfection....42
10.7 Machine preparation of the videoscope....44
10.8 Preliminary cleaning directly after the examination (see page 31)....50
10.9 Cleaning in the preparation place....50
10.10 Preparation in the washer/disinfector (AER)....51
10.11 Cleaning of accessories....51
10.12 Preparation of water bottle 13992 BS....52
10.13 Checking preparation procedure....52
10.14 Sterilization....53
10.14.1 Ethylene oxide sterilization (EO)....54
10.15 Limits of reprocessing....54
10.16 Storage....54
11 Storage/transport conditions:....55
12 Disposal....55
13 Maintenance....56
14 Important information....56
15 Servicing and repair....57
16 Limitation of liability....57
17 Warranty....58
18 Directive compliance....58
19 Appendix: Troubleshootin....59

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1 Symbols employed

Symbol	Meaning
	Consult instructions for use
	Serial number
	Manufacturer
	Fragile, handle with care
	Keep dry
	Temperature limit
	Keep away from sunlight
	Do not reuse
	Applied part type BF, defibrillation proof, IEC 60601-1
	Potential equalization connector
	DANGER: Risk of explosion if used in the presence of flammable anesthetics.
	Protective earth (ground)
	Alternating current

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Symbols employed

Symbol	Meaning
	Device is subject to the requirements of the WEEE Directive, 2012/19/EU.
	CAUTION: To reduce the risk of electrical shock, do not remove cover. Refer servicing to qualified service personnel. Do not store liquids on or above the unit. Keep out of reach of patients.
	On
	Off
	Video input
	Video output
	Color video camera
	Color video monitor
	White balance
	Output

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Symbols employed

Symbol	Meaning
	Input
	Setup Wizard
	Signifies that this feature/option is available in the Options Menu
	Signifies that this feature/option can be controlled via the User Menu
	Signifies that this feature/option can be controlled via the buttons on the videoscope
RxOnly	Federal (USA) law restricts this device to sale by or on the order of a physician.

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



Images
SILVER SCOPE® Series
Gastroscopes

① Up/down deflection wheels
② Up/down deflection brake
③ Programmable remote control buttons
④ Suction valve
⑤ Air/water valve
⑥ Rubber lip valve/biopsy valve
⑦ Right/left deflection wheels
⑧ Right/left deflection brake
⑨ Nozzle for lens cleaning/irrigation/insufflation
⑩ Light cable
⑪ Lens
⑫ Supply tube
⑬ Connection for leakage test
⑭ Connection for water bottle
⑮ Light cable connection
⑯ Port for video processor/electronic connector
⑰ Cap with connection for leakage tester
⑱ Suction/aspiration port
⑲ Air pin
⑳ Lanyard
㉑ Working channel
㉒ Additional irrigation channel with LUER lock connection (only 13826 PKS/NKS)
㉓ Distal water-jet channel (only 13826 PKS/ NKS)

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4 Safety instructions

Please read these safety instructions carefully before using the instrument.
Before using the instrument on the patient for the first time, it is imperative that you be acquainted with how the instrument operates and is controlled.

5 Intended use

Gastroscopes are used for medical, diagnostic and endoscopic treatment in the gastrointestinal tract.

Gastroscopes are designed for transient use in endoscopic interventions performed through natural orifices.

Indications

The use of gastrocopes is indicated if, in the opinion of the attending physician, a gastroscopy is indicated.

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- Inexplicable abdominal pain
- Diagnostics
- Acid reflux disease
- Tumor removal, observation after tumor removal
- Treatment of polyps

Contraindications::

The use of gastroscopes is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated, or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Gastroscopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

Videoscopes are applied parts, which should only be used with the designated supply units approved by KARL STORZ (camera control units, light sources).

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6 User qualification

KARL STORZ videoscopes may only be used by persons with an appropriate medical qualification and who are familiar with the endoscopic technique. The information given in these instructions only serves to instruct in the correct handling and preparation of the videoscopes.

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7 Warnings and cautions

Please read this manual and follow its instructions carefully, to avoid putting your patients, personnel or yourself at risk. The words, **Warning, Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the instrument. To make the words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or safety precautions must be followed to avoid possible damage to the instrument.

NOTE: *A Note indicates special information about operating the instrument, or clarifies important information.*

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8 Safety precautions when using the videoscopes

Videoscopes must be used according to the medical rules and procedures of endoscopy recognized for this method.

WARNING: *Check the videoscope prior to each procedure to ensure that it functions correctly and has been properly cleaned and disinfected.*

WARNING: *Verify the sound condition of the videoscope, as well as of accessories used in combination, before every application. Damaged videoscopes or damaged accessories must not be used.*

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WARNING: *Remove the EO/pressure compensation cap 11025 XE from the connection piece during:*

- *endoscopic interventions*
- *manual/machine preparation.*

WARNING: *Do not allow any electrically active components or instruments to come into contact with the videoscope.*

WARNING: *The videoscope must not be used during the discharge of a defibrillator. Remove the videoscope from the operation area before discharge occurs.*

WARNING: *Danger of explosion if flammable gases are used in the immediate proximity of the videoscope. The videoscope must not be used in the presence of flammable or combustible materials.*

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Safety instructions

WARNING: *High-energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the fiberoptic light cable and at the distal tip of the endoscope.*

Prolonged direct contact of light-emitting surfaces on tissue can lead to burns; this must therefore be avoided.

Temperatures remain high even after the proximal connector is removed from the light source (risk of burns to the user from the light cable).

WARNING: *Before removal from the patient, the distal tip of the videoscope must be returned to the neutral, nondeflected position (L and U on top of each other and positioned between the valves) to prevent damage to the instrument and/or injury to the patient.*

WARNING: *The videoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for videoscopes that are repaired by unauthorized personnel.*

WARNING: *Damage to the videoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.*

WARNING: *Always unplug the videoscope before carrying out any maintenance and cleaning procedures.*

WARNING: *Risk of infection due to shipment! Contaminated videoscopes can cause serious or fatal diseases during transportation or at the recipient's premises. Please cover the carrying case with the protective sheet 13990 SFN before sending contaminated instruments back to KARL STORZ.*

WARNING: *If a malfunction of the videoscope should occur during an intervention on the patient, stop the intervention immediately. Place the distal tip of the videoscope in a neutral, nondeflected position and remove the videoscope slowly and carefully from the patient.*

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Safety instructions

CAUTION: *Place the EO/pressure compensation cap 11025XE on the connection piece prior to:*

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CAUTION: *Do not steam sterilize (autoclave) the videoscopes.*

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CAUTION: *The ‘Preparation and Care’ manual must be followed. The procedures for cleaning, disinfection and gas sterilization are explained therein in detail.*

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CAUTION: *Never lay heavy objects on the videoscope.*

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NOTE: *The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts.*
In this case, the reduction of the patient leakage current should be effected using only CF application parts whenever this is possible.

NOTE: *The endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against HF currents should be ensured through the use of the correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.*

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9.3.4 Deflection of the distal tip

Rotating the large deflection wheel, labeled ‘U/D’, (up/down) deflects the distal tip. Rotating clockwise (or forwards) lowers the tip. Rotating counterclockwise (or backwards) raises the tip.

Right and left movements are controlled using the smaller deflection wheel (labeled ‘L/R’, left/right). Rotating clockwise (or forwards) = right / rotating counterclockwise (or backwards) = left (Fig. 5).

WARNING: *Do not rotate the wheels forcibly, this could damage the instrument or cause injury to the patient! Before attempting to introduce or remove the instrument, release the deflection wheel brakes!*

9.3.5 Suction

Switch on the suction pump. To perform suction, press the suction valve (marked red) on the endoscope (Fig. 6a).

Suction can also be used when the biopsy forceps are placed in the biopsy channel and the rubber lip valve is mounted. During the procedure, the rubber lip valve/biopsy valve seals off the forceps or the accessories.

With the interventional gastroscope (13826 PKS/NKS), the suction channel can also be determined after the suction valve with the two lugs has been positioned in the direction of the distal end.

Should the pressed-down valve be aligned such that the mark points towards the distal end, suction will be performed via both working channels by pressing the suction valve (Fig 6b).

Should the pressed-down valve be turned anticlockwise, suction will be performed via the left working channel (dia. 2.8 mm) by pressing the suction valve (Fig 6c).

Should the pressed-down valve be turned clockwise, suction will be performed via the right working channel (dia. 3.8 mm) by pressing the suction valve (Fig 6d).

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9.3.6 Insufflation and light

Start the light source, switch on the insufflation pump.
For insufflation, cover the opening on the blue valve with a finger (do not depress the valve, since this starts water irrigation) (Fig.7).

ⓘ ***NOTE:*** *To avoid patient discomfort, use insufflation sparingly. Excessive insufflation may cause cramps, spasms and injuries!*

9.3.7 Water irrigation for lens cleaning

Pressing on the air/water valve starts the water flow and/or lens cleaning (Fig. 8).

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9.4 Preparation before use

⚠ ***WARNING:*** *Before being used for the first time, the device must undergo a complete preparation cycle.*

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9.4.1 Testing the various functions

1. Prepare light source as described in 9.3.1.
2. Connect supply connector, water bottle and suction tube according to 9.3.3.
3. Check the supply cable for damage (e.g., breaks, cracks, twisted or squashed areas...) and the supply connectors for loose connections.
4. Examine the surface of the flexible insertion shaft visually and by touch for cracks, cuts, bulges, foreign bodies or other damage.

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The large deflection wheel controls up and down movement. To test up and down movement, release the up and down brake. Slowly rotate the large deflection wheel clockwise (or forwards), the distal end should move downwards. Rotating the deflection wheel counterclockwise (or backwards) causes upward movement.

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Handling

6. Testing deflection – right/left control
7. The smaller/outer deflection wheel controls right/left movement. To test right/left movement, release the right/left brake. Slowly rotate the small deflection wheel clockwise or forwards, the distal tip should move to the right. Rotating the deflection wheel counterclockwise (or backwards) causes movement to the left.
8. Testing deflection – brake function
 The deflection wheels should turn freely when the brakes are released.
9. Release the brakes and test their action (see 5a/b). Activate the brakes by moving the brake lever counterclockwise – ensure that the deflected position of the distal endpiece is fixed in the desired position.

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10. Check the air supply. Switch on the insufflation pump and hold the distal end of the instrument under sterile water. Cover the air control opening in the air valve (blue ring) with a finger; bubbles should emerge from the air nozzle at the distal end.

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12. Test the suction function. With the unit correctly installed and the suction pump or central suction system switched on, immerse the distal end in a container of sterile water once again. Depress the suction valve (red ring) fully. Water should be suctioned through the suction channel of the instrument via the connected suction pump and into the collecting bottle of the suction equipment.

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9.4.3 Additional recommendations for use

1. Use a lubricant from KARL STORZ (see list on our website www.karlstorz.com under 'Hygiene') for the distal end of the insertion tube.
2. Observe the lumen the whole time. If you become disoriented, withdraw the instrument gently until vision is restored.

NOTE: *Do not put lubricant on the lens to avoid clouding!*

NOTE: *Do not use oil or petroleum, silicone, Vaseline or alcohol based media, to avoid damage to the endoscope material.*

NOTE: *If debris causes the lens to cloud during the procedure, this can normally be rectified by a combination of water and insufflation. Excessive water or air can be removed by applying simultaneous pressure on the suction valve.*

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9.4.4 Operating the endoscope buttons

Depending on the video processor used, the buttons can have two functions, which are initiated by means of short and long presses of the button. For more detailed information on button assignment, please consult the manual for the camera control unit used. A list of compatible camera control units can be found in section 9.1.

NOTE: *In case of ‘red-out’, a blurred red picture on the screen, caused by the distal endoscope tip resting directly on the mucosa, withdraw the instrument gently until vision is restored.*

WARNING: *Pushing the instrument forwards despite a ‘red-out’ effect may result in serious injury to the patient.*

NOTE: *If the angling of the distal tip can no longer be properly controlled the examination must be ended immediately. Release the brakes and place the tip in the neutral position. Position the deflection wheels so that ‘U’ and ‘L’ are between the valves (see page 13/15), then carefully withdraw the instrument, whilst continuing to monitor the procedure.*

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10 Reprocessing

10.1 General information

WARNING: *The stated cleaning, disinfection and sterilization methods alone do not guarantee successful sterilization of the endoscopes. This can only be achieved if personnel follow recognized and validated preparation methods.*

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WARNING: *Remember that decontamination is only then successful if*

- *the instruments have been freed of all organic material and chemical residues.*
- *work proceeds according to recognized and validated methods.*
- *any deviation from recommended procedures are validated and tested by the user.*
- *the recommended sterilization parameters are maintained. These are only valid with cleaning, disinfection and sterilization equipment that is properly maintained and calibrated.*

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CAUTION: *The manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’ must be followed. The cleaning, disinfection and sterilization procedures are explained therein in detail.*

CAUTION: *When preparing and using the cleaning and disinfection solutions, follow the manufacturer’s instructions paying close attention to proper concentration and exposure times. Disregarding the instructions can lead to damage. Bear in mind the microbiological range of action of the chemicals used.*

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***CAUTION:** Never clean videoscopes with metal brushes.*

***CAUTION:** Constantly changing between different preparation methods causes particular wear on the materials and must therefore be avoided. However, a once-only change in methods, e.g., when acquiring new cleaning equipment, should be carried out in consultation with the manufacturers of the endoscope, the chemicals and the washer/disinfector.*

This change in method is then usually completely harmless.

***CAUTION:** Do not clean videoscopes in an ultrasound bath.*

***CAUTION:** Do not steam sterilize (autoclave) videoscopes.*

***CAUTION:** Never immerse videoscopes in normal saline solution, since even brief contact can lead to pitting and corrosion.*

***CAUTION:** During manual preparation, videoscopes may only be rinsed with microbiologically clean/sterile water.*

***CAUTION:** Decontamination in washers and disinfectors is only possible using valid procedures.*

***CAUTION:** Even parts that do not come into direct contact with the patient (e.g., video connecting cables) should be thoroughly disinfected following an intervention.*

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10.2 General cleaning

Cleaning should only be performed by authorized, trained personnel. All safety regulations should be observed (gloves, protective glasses, mouth and nose protection, clothing, German Technical Rules for Biological Agents (TRBA 250).

Thorough cleaning of the endoscope and all its accessories after every use ensures problem-free working and the best results for disinfection or sterilization procedures.

The national guidelines and recommendations for the preparation of instruments (for example RKI, ESGE/ESGENA, ASGE, LANGS, SAGES etc.) should be observed and the effectiveness verified.

CAUTION: *Before immersion in cleaning solution, check the instrument for leaks. Otherwise fluid penetrating the videoscope could cause damage.*

CAUTION: *Do not autoclave the endoscope or accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps).*

NOTE: *The endoscope must be cleaned immediately after every use. If organic material is allowed to remain and dry in the instrument, complete or partial blockage of the channel may result. Thorough cleaning is a prerequisite for effective preparation. This will avoid coagulation residue / blockages and the resulting repairs and instrument failures.*

NOTE: *Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.*

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10.2.1 Preparation for cleaning

1. The preliminary cleaning of the endoscope is performed at the place of operation. The main cleaning and disinfection is carried out in a separate preparation room.
2. Immediately after use, remove any heavy soiling (secretions, etc.) and organic material from the outer surfaces of the insertion tube with a disposable cloth moistened with cleaning solution. Then submerge the distal end in a container with cleaning solution and activate the suction to rinse the working channel.
3. Actuate the suction valve and suction at least 100 ml of cleaning solution through the suction channel. (Fig. 14). Suction up at least 200 ml with the interventional gastroscope (13826 PKS/NKS). The suction valve must then be operated in all three positions as described in section 9.3.5 (Left --> left working channel dia. 2.8 mm, right --> right working channel dia. 3.8 mm and middle --> through both channels). Make sure that all the channels are free from blockages.
4. Use the cleaning valve 13991 SRV (Fig. 13) in place of the air/water valve. (Do not use this valve during examination!).
5. Check that the water bottle is filled (max 2/3). Start the air supply (insufflation pump). Water and air are now passed through both the air channel and the water channel, depending on the position of the cleaning valve. (The process should run for at least 10 seconds).
 Air channel = Position 1
 Water channel = Position 2
6. It is important that both valve settings ('Actuated'/'Not actuated') are used multiple times for several seconds. This is the only way to ensure that both channels are rinsed through.
7. Turn off video system.

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

8. Pull out the video connector from the videoscope adaptor (VEA) at the supply connector of the endoscope when the video processor unit is switched off. Put on the sealing cap for the video cable connector and lock it with a 1/4 turn.
9. Disconnect the endoscope from the video processor unit and light source. Place the cleaning valve in a cleaning solution.
10. Then transport in a protected case to the site of reprocessing.

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10.2.2 Leakage test

Preparing for the leakage test – dry

ⓘ ***NOTE:*** *Place the cap on the videoscope adaptor (VEA) connector and screw it on tightly. Ensure that the cap is secured properly (Fig. 16).*

1. The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector to the supply connector and lock by turning clockwise 90° (Fig. 17).

ⓘ ***NOTE:*** *The rubber seal in the connection adaptor of the leakage tester 13242 LX should be treated occasionally with instrument oil.*

Before every immersion in liquid, test the endoscope for leaks. Remove all valves.

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1. Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal deflection part of the endoscope expands slightly. Move the deflection part with the deflection wheels on the control piece in all four directions as far as possible and check the manometer to see whether the pressure remains constant or falls (dry test).

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1. Immerse the device in water so that it is completely covered. It must be ensured that the leakage tester is not immersed but rather kept dry. Perform the leakage test as described above, note whether bubbles emerge from the device. If after an observation period of at least 3 minutes no more bubbles can be seen emerging from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof and can be used.

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If this test is also negative and no bubbles can be seen rising in the liquid, the leakage tester can be removed.
To do this, remove the supply connector from the liquid and press the vent button on the manometer to vent the endoscope.

NOTE: *Never release the pressure in the liquid. For this, remove the connector from the liquid and press the vent button. The test pressure must fall to 0. Then remove the leakage tester.*

If a pressure loss occurs, this means that the videoscope leaks and may no longer be used, immersed or reprocessed by machine.

CAUTION: *In this case, fill out the accompanying repairs form (see Appendix), roughly clean and remove any organic soiling/ residue from the endoscope, mark it as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the enclosed transport protection sheet 13990 SFN.*

WARNING: *Risk of infection during shipment! Contaminated instruments can cause aerogenous or contact infections during transportation or at the recipient's premises. Shipping without the plastic film packaging must be avoided.*

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10.3 Kanalsystem eines KARL STORZ SILVER SCOPE® Video-Gastroskops


(13820 PKS/NKS, 13821 PKS/NKS, 13823 PKS/NKS)

10.3 Channel system of a KARL STORZ SILVER SCOPE® video gastroscope

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
10.4 Channel system of a KARL STORZ SILVER SCOPE® Interventional video gastroscope 13826 PKS/NKS

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
10.5 Cleaning/disinfection of the videoscope: manual cleaning

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1. Remove the rubber lip valve/biopsy valve, suction valve and air/water valve and place in an cleaning solution and clean using a soft disinfected brush (13991 VB) or disposable brush.
2. Wipe the housing, insertion shaft, and supply tube carefully with a lint-free cloth soaked in an cleaner (for suitable cleaning agents or disinfectants approved by KARL STORZ, see www.karlstorz.com). Avoid shocks and pulling movements so as not to damage the insertion shaft. Then rinse with microbiologically pure/ sterile water.
3. Immerse the endoscope in the cleaning solution and rinse the channels using a syringe and the supplied irrigation tubes -13820 PKS/NKS, 13821 PKS/NKS, 13823 PKS/NKS (Fig. 23a, 23b): 13991 SS, 13991 YS; -13826 PKS/NKS (Fig. 24a, 24b): 13991 DKSA, 13991 SWJ, 13991 JASA).
4. The channels must be filled without any bubbles, to guarantee a complete covering of the surfaces.

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


CAUTION: *Only immerse instrument parts once they have successfully passed a leakage test.*

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

5. When submerged, use a suitable brush (preferably disposable) (see accessories, page 9/10) to brush below the surface of the liquid in one direction starting from the supply plug (supply plug --> control body, control body --> distal tip, biopsy channel --> distal end) (Fig. 25a: 13820 PKS/NKS, 13821 PKS/NKS, 13823 PKS/NKS; Fig. 25b and 25c: 13826 PKS/NKS)). Repeat these procedures at least 3x.

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
6. Clean the brush head (if you are using reusable brushes) at the distal end before retracting it. Do not brush backward and forward, rather brush in one direction only.

CAUTION: *Cleaning with brushes must be carried out with the instrument immersed!*

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
7. Thoroughly clean out the entire suction channel (shaft part and supply part) and the biopsy channel a number of times with a brush. Observe the exposure time for the cleaning solution.

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
8. Carefully clean all particles of dirt from the distal tip. It may be necessary to clean the air/water valve with a brush (13991 VB) with soft bristles, in order to remove hardened contaminants. Never remove debris with a needle or other sharp objects, because these can cause damage.

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9. Connect the irrigation tubes:
- 13820 PKS/NKS, 13821 PKS/NKS, 13823 PKS/NKS: 13991 SS, 13991 YS;
- 13826 PKS/NKS: 13991 DKSA, 13991 SWJ, 13991 JASA) and pump fresh cleaning fluid through the channels using a syringe. Afterwards, rinse with microbiologically clean/ sterile water.

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

Suitable brushes


Art. no.	Designation	Used for
110915-01	Cleaning brush dia. 3 mm	Working channel 2 mm
110910-01	Cleaning brush dia. 5 mm	Working channel from 2.8 to 4.2 mm
13991 VB	Valve brush	Valve seat
110917-01	Cleaning brush dia. 1.5 mm	Water-jet channel 13826 PKS/NKS

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10. Remove the device from the cleaning solution and rinse with microbiologically pure/sterile water.
11. Carefully clean all valves with the disinfected cleaning brush (13991 VB) and then rinse.
12. The valves should also be cleaned in an ultrasound bath.
13. Flow control for the interventional gastroscope 13826 PKS/NKS:

- Remove the irrigation tubes from the device
- Close the working channel inlets on the handle using a biopsy valve/rubber lip valve ((Fig. 30a) 13991 DA).
- Connect the irrigation adaptor 13991 CS to the device's supply sheath (Fig. 30b)
- Connect the suction cleaning valve ((Fig. 30c) 13991 ARV) in the valve seat of the suction valve on the handle in two subsequent positions so that it engages perceptibly in place each time.

CAUTION: *The valve 13991ARV must be set in the direction of the right working channel. The marker on the valve then points towards the right working channel (Fig. 33d). Then valve must then be set in the direction of the left working channel. The marker on the valve points towards the left working channel (Fig. 33e). At least 50 ml of microbiologically clean / sterile water must be pumped through the irrigation tube in both valve settings.*

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It must be checked that the working channels are not blocked and that the amount of liquid pumped through emerges again at the distal tip.

Once the cleaning work has been completed, clean the cleaning instruments, for example the brushes and irrigation tubes and disinfect/sterilize them. Never reuse damaged or kinked brushes.

CAUTION: *Do not autoclave the endoscope or its accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps).*

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
10.6 Manual disinfection

1. After careful rinsing, fully immerse the endoscope in the disinfecting solution.

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2. Using the irrigation tubes and a syringe with LUER connector, rinse out the channels a number of times with disinfecting solution.

 Between rinses, allow the specified time (chemical manufacturers' instructions) for the disinfecting solution to take effect in the channel system.

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3. Then rinse the channels and surfaces of the video gastroscope with microbiologically pure/sterile water.

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4. After the final step in the cleaning process has been completed, dry the entire endoscope. Thoroughly dry off all external surfaces of the endoscope, preferably with sterile compressed air at a max. of 0.5 bar/7.25 psi, thereby removing any residues of cleaning and disinfection solution.

CAUTION: The pressure must not exceed 50 kPa (0.5 bar/7.25 psi). The cleaning gun might therefore have to be operated with a pressure reducer.

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10.7 Machine preparation of the videoscope

All flexible KARL STORZ endoscopes can be connected to approved washer/disinfectors for semi- or fully automatic cleaning and disinfection. Observe the national laws and regulations for cleaning and disinfection procedures. Leakage testing, manual cleaning and (where necessary) rinsing with microbiologically pure water must be carried out beforehand. The manufacturer of the cleaning machine must be consulted when selecting the preparation programs. Usually, special programs for flexible endoscopes are available which are suited for decontamination and compatible with the materials. For machine preparation, make sure that the flexible endoscopes are not heated above 65°C.

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CAUTION: *Ensure that only KARL STORZ approved agents are used for procedures involving peracetic acid.*

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
Cleaning adaptor set 13991 CK


The cleaning adaptor set 13991 CK is used for validated preparation of KARL STORZ videoscopes from series 13xxx NKS/PKS (acc. to ISO 15883 part 4). The cleaning adaptor set is not intended for any other purpose as this may damage it.

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WARNING: *This product is not sterile when delivered and must therefore be cleaned, disinfected and if necessary, sterilized prior to initial use and each subsequent reuse.*

WARNING: *Risk of infection by improper preparation: Preparation may only be carried out by experienced and trained personnel.*

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CAUTION: *Only use chemicals which have been tested for material compatability by KARL STORZ. If you use other chemicals, get written confirmation from the chemical manufacturers that the material is compatible with KARL STORZ endoscopes.*

CAUTION: *If the preparation chemicals are changed, it is essential that you consult the manufacturer of the chemicals and of the washer/ disinfector.*

CAUTION: *Please refer to our operating instructions for details on the cleaning preparation under the section: Cleaning, disinfection and sterilization.*

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- visual checks of the connection tubes (patency) and the O-ring (wear)
- check the connections are correctly and securely in place
- check for damage

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
Assembly cleaning adaptor set

Plug the connector ② of the connection tube 13991 CL into the vent port ⑮ of the videoscope at the supply connector and secure by turning clockwise.

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Plug the cleaning adaptor 13991 CR ⑤ as far as it can go into the connector ㉒ of the air pin of the videoscope.

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Then plug connector ⑥ as far as it will go into the ‘Connection for water bottle’ ⑰ of the videoscope.

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Plug connector ⑨ of the cleaning adaptor 13991 CS as far as it will go into the suction channel ㉑ at the supply connector of the videoscope.

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If necessary, remove the suction and irrigation/ insufflation valves (④, ⑤) from the videoscope.

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
Insert the channel separator 13991 CD ⓮ with the valve pin ⓯ into the free opening of the irrigation/ insufflation valve ④ as far as it will go so that the lock ⓰ lies on top of the free opening of the suction valve ④.

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
Press the slide mechanism of the channel separator 13991 CD ⓮ downward and slide it forward until the channel separator 13991 CD ⓮ audibly locks into place. Check that the channel separator is securely in place.

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ⓘ ***NOTE:*** *If necessary, remove the cap of the biopsy/working channel ⑥ of the videoscope.*

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
CAUTION: *When the program has finished, the videoscopes are at a risk of corrosion due to the conditions in the preparation machine.*

CAUTION: *Immediately after the end of the program cycle, the videoscopes are still hot, as they have been heated up during chemical thermal preparation. Consequently the instrument's material properties are different and, as such, the prepared endoscope should be allowed to cool down for at least 20 minutes at room temperature before being used again.*

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Store the endoscopes freely suspended, dust-free and dry. Store valves separately

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10.8 Preliminary cleaning directly after the examination (see page 31)

- The outer surfaces should be wiped down from top to bottom with a lint-free disposable cloth that is soaked with cleaning solution.
- Ample amounts of cleaning solution should be suctioned through the distal end.
- Apply the air-water supply for a few seconds every cycle
- Connect the cleaning adaptor
- Disconnect the flexible endoscope from the light source

CAUTION: *Close the video connector with the cap.*

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10.9 Cleaning in the preparation place

- Carry out leakage test
- Immerse the flexible endoscope
- Clean the channels with special cleaning brushes
- Clean the complete air/water system thoroughly
- Rinse the channels

Disinfection
- Immerse the endoscope and rinse channels

Drying
with sterile compressed air

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10.10 Preparation in the washer/ disinfector (AER)

- The machine must fulfill all the cleaning, disinfection and material compatibility requirements (validated process).
- Make sure that the environmental conditions with regard to temperature and pressure are not exceeded (+60°C (tolerance +/- 5%) and 1 bar for the channels, 0.3 bar leakage test).
- Recommendation: Before placing the endoscope in the machine, clean all of the channels (working channels, water channels, air/water channels) regardless of whether or not they have been used.
- We recommend sticking to one type of cleaning agent.

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10.11 Cleaning of accessories

1. Thorough cleaning is required before disinfection or sterilization. Carefully clean the biopsy forceps and brushes whilst they are immersed in a cleaning solution using a soft, disinfected brush. Take special care to remove residues such as blood, excrement and mucous from the biopsy forceps and needles – an ultrasound bath is recommended.
2. Rinse in microbiologically pure/sterile water and dry (using sterile compressed air).
3. Disinfect accessories (work procedures as with cleaning).
4. Always test instruments after cleaning/disinfection for proper functioning and possible damage.
5. Before sterilization, carefully oil movable parts. Use only suitable, steam-resistant, physiologically safe, residue-free oil (e.g., 27656 B). After oiling, open and close the forceps several times. For details on preparation, please read the instruction leaflet of the respective accessory.

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10.13 Checking preparation procedure

Before assembly it should be checked that all parts and surfaces are undamaged, clean and dry.

Replace defective parts or send them in for repair. Before sending them in for repair, all contaminated parts must be cleaned/disinfected and sterilized or clearly labeled as unsterilized.

If residues are still present after cleaning, prepare the parts again, subject them to special treatment, discard or repair.

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Test the instruments for proper operation (deflection of the distal end, brake lever). Check the lens for a clear view, the instrument for appropriate deflection and all channels for patency (see page 9).

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10.14 Sterilization

Sterilization of instruments to be used on patients is not necessary when coming into contact with the skin and mucous membrane, as well as when introduced through natural, unsterile orifices.

When used on patients via sterile tissue, a transluminal access or piercing of the mucosa, sterilization of the instruments is necessary.

WARNING: *Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.*

WARNING: *Risk of infection: Sterilize the medical device in a disassembled state.*

CAUTION: *Flexible endoscopes must never be steam sterilized (autoclaved).*

CAUTION: *The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.*

NOTE: *All deviations from the recommended sterilization parameters must be validated by the user.*

NOTE: *When using reprocessing equipment, please observe the equipment manufacturers' instructions.*

NOTE: *Sterilization is not possible on surfaces which have been lubricated and oiled.*

CAUTION: *The STERRAD process is not suitable for sterilization of videoscopes. This process damages the videoscope!*

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

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12 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary, or your authorized dealer about the collection point in your area.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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13 Maintenance

Preventive maintenance is not essential. Regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

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14 Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and, if necessary, disinfected/sterilized before they are sent for repair. We reserve the right to return contaminated instruments, units or endoscopes to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of units, endoscopes or instruments which have been repaired by unauthorized third parties.

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15 Servicing and repair

Defective endoscopes must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts.

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16 Limitation of liability

As a supplier of this endoscope, KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:
All assembly, operations, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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Important note for users of KARL STORZ devices and instruments

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This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and, where required, sterilization of the CMOS Video-Rhino-Laryngoscope. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

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Rx only *Federal (USA) law restricts this device to sale by or on the order of a physician.*

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Contents

1 Important note for users of KARL STORZ devices and instruments I
2 Pictures of the instruments 4
3 Symbol description 6
4 Safety information 8
5 Intended use 8
5.1 Indications 8
5.2 Contraindications 8
6 User qualification 9
7 Kompatibilität 9
8 Warnings and cautions 10
8.1 Safety precautions when using the CMOS Video-Rhino-Laryngoscopes 10
9 Initial operation 14
9.1 Unpacking the device 14
9.2 Basic equipment 14
9.3 Operating conditions 14
9.4 Electromagnetic compatibility 14
9.5 External system components 15
10 Handling 16
10.1 Operating 16
10.2 Holding correctly 16
10.3 Handling the shaft 17
10.4 Handling the deflection mechanism 18
10.5 Connection to the monitor 8403 ZX 20
10.6 Use of the HDMI output 20
10.7 Functions of the C-MAC® Monitor 8403 ZX 21
10.7.1 Focus 21
10.8 Function buttons of the endoscope 22
10.9 Monitor function keys 22
10.10 Preparation prior to use 23
10.10.1 Inspecting the surfaces 23
10.10.2 Checking the deflection 24
10.11 Switching the CMOS Video-Rhino-Laryngoscope on/off 25

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Contents

10.12 Inspecting the image quality 26
10.13 Switching off the C-MAC® Monitor 26
10.14 System reset 27
10.15 Monitor bag 8403 YD 27
11 Reprocessing 28
11.1 Overview of reprocessing procedures 28
11.2 Accessories 30
11.3 Preliminary cleaning directly after use .. 30
11.4 Leakage test 31
11.4.1 Performing the leakage test – dry 31
11.4.2 Performing the leakage test – in liquid 31
11.5 Disassembly 32
11.6 Manual reprocessing 32
11.6.1 Manual cleaning 32
11.6.2 Manual disinfection 33
11.6.3 Manual High Level Disinfection (HLD) .. 33
11.6.4 Drying 33
11.7 Automated reprocessing 34
11.7.1 Manual preparation 34
11.7.2 Machine cleaning/chemical-thermal disinfection 34
11.8 Assembly, inspection and care 34
11.9 Packaging 35
11.10 Sterilization 36
11.10.1 Hydrogen peroxide (H_2O_2) sterilization - ASP STERRAD® 36
11.10.2 Hydrogen peroxide (H_2O_2) sterilization – STERIS® AMSCO® V-PRO™ 37
11.10.3 Ethylene oxide sterilization (EO) 37
11.10.4 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1® 37
11.10.5 Low-temperature steam and formaldehyde process (LTSF) 38
11.11 Limits of reprocessing 38

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
Contents

12 Storage 38
13 System description 40
14 Technical description 41
15 Maintenance 43
15.1 Servicing and repair 43
15.2 Maintenance 43
15.3 Limitation of liability 43
15.4 Warranty 44
15.5 Standard compliance 44
15.6 Directive compliance 45
15.7 Important information 45
15.8 Disposal 45
16 Appendix 46
16.1 Troubleshooting 46
16.2 Accessories 48
16.3 Electromagnetic Compatibility (EMC) Information 49

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
ⓘ ***NOTE:*** *The illustrations in this instruction manual show model 11101 CM. However, all descriptions similarly apply to model 11102 CM.*

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Controls, displays, connectors and their uses

① Connecting cord
② Connection for leakage tester or pressure compensation cap (vent port)
③ Distal tip
④ Shaft
⑤ Handle
⑥ Deflection lever
⑦ Leakage tester 13242 XL
⑧ Pressure compensation cap 11025 E
⑨ Monitor
⑩ Connection 2 for video endoscope
⑪ Connection 1 for video endoscope
⑫ Mains connection
⑬ Reset button
⑭ Holder for SD memory card
⑮ Status bar
⑯ Pilot lamp
⑰ On/Off switch
⑱ Pilot lamp charging (orange)
⑲ Softkeys
⑳ Function symbols
㉑ USB port
㉒ HDMI port
㉓ Video recording Start/Stop
㉔ White balance
㉕ Still image recording
㉖ C-HUB®

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Please read these safety instructions carefully before using the device. Before using the instrument on the patient it is imperative that you be acquainted with how the instrument operates and is controlled.

5 Intended use

Video-Rhino-Laryngoscopes are used for endoscopic imaging and image transmission, thus enabling an endoscopic diagnosis of body cavities or hollow organs in ENT medicine.

Video-Rhino-Laryngoscopes are intended for temporary use in both surgical invasive interventions and interventions performed through natural orifices.

5.1 Indications

The use of the Video-Rhino-Laryngoscope is indicated if, in the view of the attending physician, an assessment of the nasal, paranasal sinuses, nasopharynx, laryngeal and esophageal structures is indicated with the flexible endoscope.

5.2 Contraindications

No contraindications relating directly to the product are currently known.

The use of the Video-Rhino-Laryngoscope is contraindicated if, in the opinion of the attending physician, the endoscopic surgical method as such is contraindicated or if the patient is not able to undergo surgery or anethesia due to his or her general condition.

The Video-Rhino-Laryngoscope must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

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6 User qualification

The KARL STORZ CMOS Video-Rhino-Laryngoscope may only be used by persons with an appropriate medical qualification and who are acquainted with the rhino-laryngoscopic technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the CMOS Video-Rhino-Laryngoscope. It is not intended as an introduction to the technique of rhino-laryngoscopy.

7 Compatibility

CAUTION! *Non-compatible accessories can lead to malfunctioning or damage to the CMOS video rhino-larynogscope.*

Only use the compatible accessories designated by KARL STORZ.

The CMOS video rhino-laryngoscope is exclusively intended for use with the following KARL STORZ devices.

- **C-MAC® monitor 8403 ZX** (art. no. of the instructions for use: 96076008D)
- **C-HUB® II 20 2903 20** (art. no. of the instructions for use: 96206529D)

The instructions for use of the device used in combination must be followed. The instructions for use can be requested from KARL STORZ or downloaded at www.karlstorz.com.

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8 Warnings and cautions

Please read this manual carefully and follow the instructions exactly. The words **Warning, Caution,** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

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CAUTION: *A Caution indicates that particular service procedures or safety precautions must be followed to avoid any damage to the product.*

NOTE: *A note indicates special information about operating the product, or clarifies important information.*

8.1 Safety precautions when using the CMOS Video-Rhino-Laryngoscopes

The CMOS Video Rhino-Laryngoscope must be used according to the recognized medical rules and procedures of endoscopy.

NOTE: *The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts.*

In this case, the reduction of the patient leakage current should be effected using only CF application parts as far as this is possible.

NOTE: *The flexible endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against high frequency currents should be ensured through the use of correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.*

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Application

WARNING: *Every time the CMOS Video-Rhino-Laryngoscope is used, you must first check that it, and any accessories used in combination with it, is in perfect condition. Damaged CMOS Video-Rhino-Laryngoscopes or damaged accessories must not be used.*

WARNING: *Combinations of medical devices are only then assured to be safe if*

- *they are identified as such in the respective instruction manual or*
- *the intended purpose and the interface specifications of the devices used in combination permit this.*

WARNING: *Danger of explosion if flammable gases are used in the immediate proximity of the CMOS Video-Rhino-Laryngoscope. The CMOS Video-Rhino-Laryngoscope must not be used in the presence of flammable or combustible media.*

WARNING: *Only the KARL STORZ devices listed in chapter 7 may be used. Accessories from other manufacturers can damage CMOS video rhino laryngoscopes.*

WARNING: *The CMOS Video-Rhino-Laryngoscope must not be used during the discharge of a cardiac defibrillator. Remove the CMOS Video-Rhino-Laryngoscope from the operation area before discharge takes place.*

WARNING: *When in operation, the integrated LED light source heats up the CMOS Video-Rhino-Laryngoscope. Switch off the Video-Rhino-Laryngoscope before putting it down and do not place it in the proximity of the patient (as per IEC 60601-1-1).*

WARNING: *Danger of glare and burns! Never look directly into the light of the CMOS Video-Rhino-Laryngoscope. The LED light radiated through the CMOS Video-Rhino-Laryngoscope may cause glare or high temperatures at the distal light outlet and therefore lead to burns.*

WARNING: *If the CMOS Video-Rhino-Laryngoscope should malfunction during application on a patient, stop the application immediately. Place the distal tip of the CMOS Video-Rhino-Laryngoscope in the neutral, nondeflected position and remove*

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the CMOS Video-Rhino-Laryngoscope slowly and carefully from the patient.

***CAUTION:** To avoid damaging the monitor, do not remove the SD card or connecting cord while a video is recording. Only change the SD card when the monitor is switched off.*

***CAUTION:** Do not allow any live components or instruments to come into contact with the CMOS Video-Rhino-Laryngoscope.*

***WARNING:** Damage to the CMOS Video-Rhino-Laryngoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.*

Reprocessing

***WARNING:** The CMOS Video-Rhino-Laryngoscope is delivered unsterilized and must be cleaned and disinfected/sterilized before initial and subsequent use.*

***CAUTION:** Do not steam sterilize (autoclave) the CMOS Video-Rhino-Laryngoscopes under any conditions.*

***CAUTION:** The maximum permissible temperature for decontamination, sterilization, transport and storage is 65 °C.*

***CAUTION:** Do not hit the distal tip of the CMOS Video-Rhino-Laryngoscope against hard objects, and never use a sharp object to remove dirt.*

***NOTE:** Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided.*

***WARNING:** The instructions 'Cleaning, Sterilization and Care of KARL STORZ Instruments' must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail.*

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Safety information

Storage

CAUTION: *Do not place the CMOS Video-Rhino-Laryngoscope on devices that become very hot. Do not expose it to direct sunlight or excessive heat, to avoid damage to the device.*

CAUTION: *Never lay heavy objects on the CMOS Video-Rhino-Laryngoscope.*

Maintenance

WARNING: *The CMOS Video-Rhino-Laryngoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for CMOS Video-Rhino-Laryngoscopes that are repaired by unauthorized personnel.*

WARNING: *Danger of electric shock! Refer servicing for this device to authorized personnel only. Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims. Exceptions are those that are described in the instruction manual.*

WARNING: *In order to prevent infection and protect personnel, it is strictly forbidden to ship contaminated medical devices. To avoid contact and aerogen infections, medical devices must be decontaminated on-site.*

WARNING: *The* ***pressure compensation cap*** ⑧ *must be* ***removed*** *from the vent port* ② *for:*

- *endoscopic interventions*
- *immersing in liquid, e.g. disinfectant solution*

WARNING: *The* ***pressure compensation cap*** ⑧ *must be* ***in place*** *on the vent port* ② *for:*

- *gas and plasma sterilization*
- *shipping*

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9.5 External system components

Charger

Using the medically approved charger (Art. No. ET27-30-0004515), the Li-Ion battery of the C-MAC® Monitor can be recharged extremely quickly.

Moreover, this charger enables the C-MAC® Monitor to be operated with a discharged Li-Ion battery.

The charger has a green LED display ㉘, which indicates that the output voltage is correct.

Country specific configuration

The charger features interchangeable primary adapters, which allow the system to be adapted to the country-specific power supply outlets.

a Primary adaptor for EUROPE:
ET27-30-0004 516
b Primary adaptor for UK:
ET27-30-0004 517
c Primary adaptor for US/Japan:
ET27-30-0004 518
d Primary adaptor for Australia:
ET27-30-0004 519
e Primary adaptor for USA:
ET27-30-0004 520
f Adaptor cable-ET27-30-0004 370

Press the button on the power supply unit to change the adaptor (fig. 1). This releases the adaptor which can then be lifted out.

When inserting a new adaptor, this must be placed in the guide (fig. 2) and pushed toward the cable until there is an audible click (fig. 2, see arrow).

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10 Handling

NOTE: *Damage to the CMOS Video-Rhino-Laryngoscope resulting from incorrect operation is not covered by the manufacturer's warranty.*

WARNING: *The CMOS Video-Rhino-Laryngoscope is delivered unsterilized, so it has to be reprocessed prior to initial use and before each subsequent reuse.*

CAUTION: *Do not store or disinfect the flexible CMOS Video-Rhino-Laryngoscope together with sharp objects in order to avoid causing damage to the outer casing.*

NOTE: *The C-MAC® monitor has passed a fall and crash test successfully.*

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10.1 Operating

The following instructions outline the correct operation of the CMOS Video-Rhino-Laryngoscope and the corresponding accessories. This is not an introduction to the techniques of rhino-/laryngoscopy. Instructions on the endoscopic technique can be found in the relevant medical literature.

The CMOS Video-Rhino-Laryngoscope is intended for visualisation and diagnostic purposes.

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10.2 Holding correctly

Grip the handle of the CMOS Video-Rhino-Laryngoscope in your left or right hand, as shown in the picture. With this kind of control grip, the CMOS Video-Rhino-Laryngoscope can be operated with one hand.

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10.3 Handling the shaft

CAUTION: *Do not carry the CMOS Video-Rhino-Laryngoscope by the shaft. Do not pull, clamp or twist the shaft of the CMOS Video-Rhino-Laryngoscope.*

CAUTION: *The shaft of the Video-Rhino-Laryngoscope must never be kinked or coiled up too tightly. It is best to store the video rhinolaryngoscope in a suspended position or to keep it in a tray specifically designed for sterilizing and storage (39405 AS).*

WARNING: *Do not apply lubricant to the lens, to avoid clouding the image and thus posing a risk to the patient!*

Use a suitable water-soluble lubricant for the distal section of the shaft (For materials recommended and approved insertion by KARL STORZ, please see www.karlstorz.com).

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the shaft of the CMOS Video-Rhino-Laryngoscope out of the patient's body.

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10.4 Handling the deflection mechanism

NOTE: *The neutral position of the distal tip is straight, without any deflection whatsoever.*

WARNING: *If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the CMOS Video-Rhino-Laryngoscope must not be used.*

WARNING: *The deflection lever must never be operated with force or abruptly to avoid perforation and injury to the patient or damaging the device irreparably.*

WARNING: *The distal tip of the CMOS Video-Rhino-Laryngoscope must be straight, not angled, when the videoscope is inserted into the patient*

WARNING: *Before removing the CMOS Video-Rhino-Laryngoscope from the patient, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.*

WARNING: *If a functional defect occurs with the CMOS Video-Rhino-Laryngoscope during use on the patient, discontinue use immediately. Put the distal tip of the CMOS Video-Rhino-Laryngoscope into the neutral, non-deflected position and remove the CMOS Video-Rhino-Laryngoscope slowly and carefully from the patient.*

CAUTION: *The resistance that can be felt when deflecting the tip using the deflection lever increases as the angle of deflection increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.*

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Insert the connecting cord ① of the flexible CMOS Video-Rhino-Laryngoscope into either connection socket ⑩ or ⑪ on the monitor.

The connecting cord ① can be connected or disconnected during operation.

Insert the cable from the power supply unit into the connection socket ⑫ of the C-MAC® Monitor and connect the power supply unit with the mains socket. The power cord can be fixed to the connection socket using the union nut.

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ⓘ ***NOTE:*** *The connected monitor must support a resolution of 1024 × 768 or 1280 × 720.*

ⓘ ***NOTE:*** *To connect the C-MAC® Monitor to a KARL STORZ monitor, the HDMI-DVI adaptor with article number **20**9190 10 is required.*

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⚠ ***WARNING:*** *The KARL STORZ flexible CMOS Video-Rhino-Laryngoscope can not be used in combination with camera control units and light sources. The use of non-compatible units can lead to patient injury or unit damage.*

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10.8 Function buttons of the endoscope

The endoscope is equipped with 3 keys, which are used to activate the functions displayed as symbols on the endoscope.

The following functions are available:

- Video recording ㉓: Starts or stops the recording of a video stream on the SD memory card inserted in the receptacle.
- White balance ㉔: Enables the user to carry out a manual white balance. Point the tip of the videoscope at a white surface and press the button for 1 sec.
- Single image capture ㉕: Saves the live image on the screen on the SD memory card.

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10.9 Monitor function keys

The monitor is equipped with 4 softkeys. These can be used to activate the functions displayed as symbols on the screen.

NOTE: *The* ***video recording*** *and* ***image storage*** *symbols are only visible if an SD memory card is inserted. The* ***change the camera unit*** *symbol is only visible if two video units are connected.*

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10.10 Preparation prior to use


NOTE: *Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding preliminary cleaning.*

WARNING: *Every time the CMOS Video-Rhino-Laryngoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged CMOS Video-Rhino-Laryngoscopes or damaged accessories must not be used.*

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10.10.1 Inspecting the surfaces

WARNING: *The surface of the flexible shaft must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.*

CAUTION: *Do not squeeze the distal tip of the CMOS Video-Rhino-Laryngoscope too tightly.*

Every time the CMOS Video-Rhino-Laryngoscope is used, you must first check that its surface, and above all the surface of the shaft and distal tip, are in perfect condition and that there are no sharp edges. CMOS Video-Rhino-Laryngoscopes with sharp edges must not be used. Otherwise
serious injuries may be caused to the patient or the user.

Inspect the connecting cord for damage (such as breaks, cracks, twisted or squashed sections...) and check the supply connector for loose connections.

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10.10.2 Checking the deflection

WARNING: *Actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated with force or abruptly, to avoid injuring the patient or damaging the device! Move the deflection lever smoothly and slowly!*

WARNING: *Before use, always check that the deflection lever and the distal tip connected to it function correctly. Defective CMOS Video-Rhino-Laryngoscopes must not be used.*

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10.11 Switching the CMOS Video-Rhino-Laryngoscope on/off

Switching on/off: The CMOS Video-Rhino-Laryngoscope with the integrated LED light source is switched on by pressing the power button ⑰ on the monitor. Pressing the key a second time switches the unit off.

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NOTE: *The function keys ⑰ on the monitor are used for switching the unit on/off. The CMOS Video-Rhino-Laryngoscope itself does not have any function keys.*

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Video recording: When the control key ‘Video recording’ on the C-MAC® Monitor or the flexible CMOS Video-Rhino-Laryngoscope is pressed, the flexible CMOS Video-Rhino-Laryngoscope records a video stream on the SD memory card. Pressing the key a second time stops the recording.

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Single image capture: When the control key ‘Single image capture’ on the C-MAC® Monitor or the flexible CMOS Video-Rhino-Laryngoscope is pressed, the current image on the monitor is saved on the SD memory card.

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10.12 Inspecting the image quality

WARNING: *If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the CMOS Video-Rhino-Laryngoscope from the patient's body gently and carefully. If the CMOS Video-Rhino-Laryngoscope is used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.*

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NOTE: *When the C-MAC® Monitor is switched off, the integrated light source in the flexible CMOS Video-Rhino-Laryngoscope also switches off.*

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10.14 System reset

If a system error occurs (unit no longer reacts, although connections and power supply are in order), a system reset has to be implemented. This is the equivalent of rebooting a computer. For this, a thin screwdriver or similar is inserted into the round opening next to the SD card holder ⑭ and carefully pushed downward in order to activate the reset button ⑬.

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10.15 Monitor bag 8403 YD

The optionally available protective bag 8403 YD facilitates the mobile use of the monitor.

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11 Reprocessing

11.1 Overview of reprocessing procedures

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- ● This reprocessing step is part of the validated procedure
- ❍ This reprocessing step may be carried out additionally in order to aid the cleaning. The sterilization procedure is approved in respect of material compatibility; however, it must be validated by the operator on site in respect of effectiveness.
- **X** On the grounds of potential damage, this reprocessing step must NEVER be performed

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Reprocessing

WARNING: *Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.*

WARNING: *Risk of infection from non-sterile delivery. Reprocess the device before each use. Check the medical device for visible soiling before use.*

WARNING: *Before reusing the device, replace any defective parts or send the flexible endoscope for repair.*

ATTENTION: *The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.*

CAUTION: *Never clean flexible endoscopes with metal brushes.*

CAUTION: *Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.*

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CAUTION: *To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.*

NOTE: *When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.*

NOTE: *Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.*

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ⓘ ***NOTE:*** *Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be found at www.karlstorz.com.*

ⓘ ***NOTE:*** *National laws and regulations must be observed.*

ⓘ ***NOTE:*** *We recommend using a plastic basin for disinfection to avoid damaging the endoscope. Furthermore, use of a plastic basin prevents electrolytic corrosion which can occur if there are different metals in a solution.*

ⓘ ***NOTE:*** *The manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’ can be downloaded or requested by visiting www.karlstorz.com.*

11.2 Accessories

Accessories required for carrying out reprocessing:

Accessories: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/ adaptor
Clean surface	27652
Ventilation connector	13242XA

11.3 Preliminary cleaning directly after use

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the endoscope immediately after use.

Wipe-clean the exterior surfaces of the endoscope with a cloth moistened with disinfectant.

ⓘ ***NOTE:*** *To avoid damaging the shaft, do not use pushing and pulling movements.*

ⓘ ***NOTE:*** *Do not use disinfectants for preliminary cleaning, these cause fixation of proteins and organic residues*

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11.4 Leakage test

The endoscope should be tested to ensure it is leakproof after use and before reprocessing.

11.4.1 Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector securely to the supply connector vent port and lock by turning the bayonet catch 25° clockwise.

ⓘ ***NOTE:*** *The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.*

Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer manometer. The test pressure can drop by up to 10 mmHg until the endoscope is completely filled with air.

After stabilizing the pressure on manometer wait at least 2 minutes to see whether the pressure remains constant or drops.

11.4.2 Performing the leakage test – in liquid

ⓘ ***NOTE:*** *To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.*

ⓘ ***NOTE:*** *Any lumina must be covered with no air bubbles by filling in a targeted manner.*

Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in ‘Leakage test – dry’ and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer

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does not show a fall in pressure, the device is leakproof.

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CAUTION: *In this case, fill out the accompanying repairs form (see Appendix or Art. no. 96096004DER), mark the endoscope, which has been cleaned by wiping to remove organic soiling/ residue, as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the transport protection sheet 13990 SFN.*

11.5 Disassembly

The valves must be removed in their individual parts prior to cleaning and disinfection.

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The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

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Horizontal storage is possible for sterilized endoscopes. KARL STORZ recommends the use of a specially designed tray for sterilization and storage.

CAUTION: *Do not expose flexible endoscopes to direct sunlight or excessive heat.*

CAUTION: *The distal tip of the flexible endoscope must be in the neutral position during storage.*

CAUTION: *Do not store the flexible endoscope in the case.*

CAUTION: *The shaft of the flexible endoscope must never be kinked or coiled up too tightly.*

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13 System description

The CMOS Video-Rhino-Laryngoscope from KARL STORZ has a 300 mm pellethane working shaft, 3.70 mm (11101 CM) or 2.9 mm (11102 CM) in diameter. The PEEK (polyetheretherketone) distal tip is 3.70 mm or 2.9 mm in diameter and houses the CMOS video module. The distal tip of the flexible CMOS Video-Rhino-Laryngoscope can be deflected 140° up or 140° down. The handle houses the deflection lever, the function keys and the LED light source.

The CMOS video rhino laryngoscope from KARL STORZ is connected to the recommended monitor 8403 ZX using the integrated 2 m supply cable. Images and videos can be saved using a SD memory card and the control keys on the endoscope or on the monitor, which are assigned, for example, the functions of video recording and image capture.

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C-MAC® Monitor 8403 ZX

Length:	230 mm
Width:	155 mm
Depth:	54 mm
Weight:	1250 g
Drip-water protection:	IP54

CAUTION: *Drip-water protection can only be ensured if the connection sockets are sealed using the designated rubber seals.*

Classification:	IEC Class II
Input voltage:	12 V
Current consumption:	1.5 A max.
Battery:	Lithium-Ion, 3.6 V 7800 mAh
Resolution:	1280 x 800
Screen width:	150 mm
Screen height:	94 mm
Memory card:	SD card
Video out signal:	The connected monitor must support the 720p mode.

Operating/storage conditions

Operating conditions:	0 °C to +40 °C 30 % to 70 % rel. humidity
Storage conditions:	-10 °C to +60 °C 5 % to 95 % rel. humidity
Atmospheric pressure:	700 to 1080 hPa

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Technical description

SD card

Throughput:	10 MB/s
Compatibility:	Compatible with the Sandisk Ultra II HD card and SDHC cards class 10

Color settings

Basic Setting

Brightness	50 %
Contrast	35 %
Saturation	31 %
Hue	50 %

Endoscope 11101 CM/11102 CM

Direction of view:	0°
Angle of view:	85° (11101 CM) 100° (11102 CM)
Deflection:	up 140°, down 140°
Working length:	30 cm
Outside diameter:	3,7 mm (11101 CM) 2,9 mm (11102 CM)
Illumination:	integrated LED light source
Electrical power consumption:	1W
Sensor:	CMOS technology
Weight	245 g (11101 CM) 325 g (11102 CM
Drip water protection:	IPX8
Classification:	IIa
Applied part:	Type BF

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15.3 Limitation of liability

KARL STORZ SE & Co. KG shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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15.6 Directive compliance

According to Medical Device Directive (MDD):

11101 CM /11102 CM:
This medical product belongs to Class IIa
8403 ZX: This medical product belongs to Class I

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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15.7 Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender.

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KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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15.8 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its service life, dispose of the device as electronic waste. Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information about your local collection point. Within the scope of application of this directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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Electromagnetic Compatibility user information

The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING: *Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.*

WARNING: *The use of an accessory or cable with the flexible CMOS Video-Rhino-Laryngoscope11101 CM/11102 CM together with the Monitor 8403 ZX other than those specified in this manual may result in increased emissions or decreased immunity of the flexible CMOS Video-Rhino-Laryngoscope11101 CM/11102 CM together with the Monitor 8403 ZX. When using an accessory or cable with the CMOS Video-Rhino-Laryngoscope11101 CM/11102 CM together with the Monitor 8403 ZX other than those specified in this manual, it becomes the responsibility of the user of the flexibleCMOS Video-Rhino-Laryngoscope11101 CM/11102 CM together with the Monitor 8403 ZX to determine compliance with IEC 60601-1-2 when using this item.*

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Appendix
Electromagnetic Compatibility (EMC) Information

WARNING: *Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the flexible CMOS Video-Rhino-Laryngoscope11101 CM/11102 CM together with the Monitor 8403 ZX including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.*

NOTE: *The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.*

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Accessories and cables which have been shown to comply with EN/IEC 60601-1-2:				
Cable type	**Shielded**	**Length [m]**	**Ferrite**	**Used for**
Charger ET27-30-0004515	No	2.0	No	Charging the internal battery

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Table 1
Guidelines and manufacturer's declaration – electromagnetic immunity

The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is intended for use in an electromagnetic environment specified below. The user of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines lines ± 1 kV for input and output lines 100 kH repetition	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	Voltage Dip: Dip to 0% for 1 cycles @ 0°phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100% for 250/300 cycles	Voltage Dip: Dip to 0% for 1 cycles @ 0°phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be me asured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3V/m 80 MHz to 2.7 GH	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	

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Table 2
Test levels for Proximity fields from RF wireless communications equipment

Test Frequency MHz	Band MHz	Service	Modulation	Immunity Test Level V/m	Compliance level V/m
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	27	27
450	430 – 470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz Sine Wave	28	28
710	704 – 787	LTE Band 13 & 17	Pulse modulation 217 Hz	9	9
745					
780					
810	800 – 960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	28	28
870					
930					
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28	28
1845					
1970					
2450	2400 – 2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	28	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
5500					
5785					

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Table 3
Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting

The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
			Portable and mobile RF communications equipment should be used no closer to any part of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF disturbances acc. to IEC 61000-4-6	3 V_{rms} 150 kHz to 80 MHz	3 V_{rms}	$d = 1.2\sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey[a] should be less than the $d = 1.2\sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz compliance level in each frequency range[b]. Interference may occur in the vicinity of equipment marked with the following symbol:
			

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.
Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.

[a] Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is used exceeds the above compliance levels, the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
[b] Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 **Guidelines and manufacturer's declaration – electromagnetic emissions**		
The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is intended for use in the electromagnetic environment specified below. The customer or user of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX uses RF energy only for its internal function. The customer or user of the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class B	The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is suitable for use in all establishments including domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class B	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5
Recommended separation distances between portable and mobile RF telecommunication devices and the flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX

The flexible CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CMOS Video-Rhino-Laryngoscope 11101 CM/11102 CM together with the Monitor 8403 ZX as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	**Separation distance d [m] as a function of the transmitter frequency**		
	150 kHz to 80 MHz $d = 1{,}2\sqrt{P}$	**80 MHz to 800 MHz** $d = 1{,}2\sqrt{P}$	**800 MHz to 2.5 GHz** $d = 2{,}3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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INSTRUCTION MANUAL
Flexible Cystoscope HD-VIEW 11272VH/VHU

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Content

1 Important information 1
2 Symbols employed 3
3 Intended use 6
3.1 Indications 6
3.2 Contraindications 7
4 User qualification 8
5 Warnings and cautions 8
6 Safety instructions 9
6.1 Unpacking 18
6.2 Basic equipment 19
6.3 Operating conditions 20
7 External system components 20
8 Preparation before use 21
9 Handling 21
9.1 Operation 21
9.2 Holding the device correctly 21
9.3 Operating the suction valve 22
9.4 Handling the sheath 22
9.5 Handling the deflection mechanism 23
9.6 Handling the working channel 24
9.7 Handling the accessories (forceps, cleaning brushes, HF electrodes etc.) 26
9.8 Inspecting the surfaces 28
9.9 Inspecting the deflection 29
9.10 Connecting to IMAGE1 S™ 29
9.11 Focus 29
9.12 Function buttons of the endoscope 30
9.13 Switching the flexible cystoscope HD VIEW ON/OFF 31
10 Reprocessing 32
10.1 General warnings 32
10.2 Accessories 33
10.3 Precleaning directly after use 34
10.4 Leakage test 35
10.5 Disassembly 37

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Content

10.6 Manual disinfection 40
10.6.1 Manual High Level Disinfection (HLD) 40
10.6.2 Drying 41
10.7 Machine cleaning and disinfection 41
10.7.1 Manual preparation 41
10.7.2 Connecting 42
10.7.3 Machine cleaning/chemical-thermal disinfection 42
10.8 Assembly, inspection, care and packaging 43
10.8.1 Packaging 44
10.9 Sterilization 44
10.10 Limits of reprocessing 46
10.11 Storage 46
11 Reprocessing of LUER port with stopcocks 48
11.1 Accessories 49
11.2 Preparation for cleaning and disinfection 49
11.3 Manual precleaning 49
11.4 Machine cleaning and disinfection 50
11.5 Connecting 51
11.6 Assembly, inspection and care 51
11.7 Packaging systems 51
11.8 Sterilization 52
11.9 Limits of reprocessing 52
12 Technical data 53
13 Servicing and repair 55
14 Maintenance 55
15 Limitation of liability 55
16 Warranty 55
17 Directive compliance 56
18 Important information 56
19 Disposal 56
20 Electromagnetic Compatibility user information 57

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This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and, if appropriate, sterilization of the flexible cystoscope HD VIEW. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

It is recommended that the suitability of the products for the intended procedure be checked prior to use.

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2 Symbols employed

	Manufacturer
	Consult instructions for use
REF	Catalogue number
QTY	Number of products in the product packaging
SN	Serial number
	Non-sterile
	Fragile, handle with care
	Keep dry
	Temperature limit
RxOnly	Federal (USA) law restricts this device to sale by or on the order of a physician.

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

Controls

① Connector
② Suction valve
③ Connection for leakage tester or pressure compensation cap (vent port)
④ Handle
⑤ Irrigation and working channel
⑥ Deflection lever
⑦ Distal tip
⑧ Working channel connector seal 11014 DS
⑨ Leakage tester 13242 XL
⑩ Suction valve, sterile, for single use 091011-20
⑪ Suction valve, reusable 11301 CB1
⑫ Pressure compensation cap 11025 E
⑬ LUER port with stopcocks 11014 L
⑭ Double LUER port 11014 LU
⑮ Cleaning adaptor for working channel 11014 RA
⑯ Cleaning adaptor 11301 CD1
⑰ Right navigation button
⑱ Left navigation button
⑲ Menu button

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3 Intended use

Flexible cysto-urethroscopes (fiber endoscopes and videoscopes) are used for endoscopic imaging of the urethra and the bladder for the purposes of diagnosis and treatment. Flexible cysto-urethroscopes with an integrated working channel allow for the use of flexible auxiliary instruments such as forceps, stone baskets etc., coagulation electrodes and laser fibers for lithotripsy or for the ablation of benign/malignant tissue. With flexible cysto-urethroscopes, the urethra and the entirety of the bladder can be visualized/treated. The special feature of the flexible cysto-urethroscope lies in the fact that it can be used percutaneously in the kidneys or suprapubically in the bladder via an operating sheath.

Flexible cysto-urethroscopes (fiber endoscopes and videoscopes) are designed for temporary or short-term use in endoscopic interventions performed through natural orifices or surgical invasive interventions.

3.1 Indications

Description of all known indications:

The use of cysto-urethroscopes is indicated if urethro-cystoscopy is indicated in the opinion of the attending physician.

This includes, amongst other things:

- Clarification of micro and macro hematuria
- Clarification/biopsy in the event of suspected bladder tumors
- Control urethro-cystoscopy as per TUR-B or P
- Clarification in particular in the event of suspected tumor infiltration due to adjacent organ cancer (e.g., colon carcinoma, cervical carcinoma)
- Clarification or treatment in particular in the event of suspected vesicoenteric or obstetric fistulas
- Urethral or vesical foreign body removal
- Clarification of recurring urinary tract infections
- Clarification of LUTS/BES, subvesical obstructions, BOO, neurogenic micturition disorder

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- Visual assessment and possibly treatment (e.g., injection) of the orifices of the ureters (localization, condition, etc.)
- Visual clarification of the bladder sphincter in the event of incontinence
- Visual assessment and possibly treatment of the urethra (tumor, stricture/stenosis, fistula)
- Visual assistance for indwelling catheter system, e.g., with urethral trauma in the event of a pelvic fracture
- Retrograde pyelography, ureteral splinting

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3.2 Contraindications

The use of the cysto-urethroscopes is contraindicated if, in the opinion of the attending physician, the surgical or examination method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Cysto-urethroscopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

The following contraindications also apply:

- Acute urethritis
- Acute epididymitis
- Acute prostatitis
- Acute/florid cystitis/urinary tract infection
- Possible coagulation disorders
- Restless, agitated patients

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4 User qualification

The flexible cystoscope HD VIEW from KARL STORZ may only be used by persons with an appropriate medical qualification and who are acquainted with the flexible cystoscopy technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the flexible cystoscope HD VIEW. It is not intended as an introduction to the techniques of flexible cystoscopy.

5 Warnings and cautions

⚠ **HAZARD SYMBOL:** *This is a hazard symbol. It warns of the risk of injury. To prevent injuries, all the measures marked with the hazard symbol must be followed.*

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6 Safety instructions

The flexible cystoscope HD VIEW must only be used according to the medical rules and procedures of endoscopy recognized for this method.

⚠ WARNING:

Risk of injury and damage to the product:

➢ *It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other devices at any time.*

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➢ *Never look directly into the light emitted by the flexible videoscope.*

⚠ WARNING:

Risk of explosion: If flammable anesthetics are used in the immediate proximity of the flexible videoscope.

➢ *The flexible videoscope must not be used in the presence of flammable or combustible media.*

⚠ WARNING:

Risk of injury: When in operation, the integrated LED light source heats up the flexible videoscope.

➢ *Switch off the flexible videoscope before putting it down, and do not place it in the vicinity of the patient (as per chapter 16 of the 3rd edition of IEC 60601-1).*

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⚠ WARNING:

Risk of injury: High energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the fiberoptic light cable and at the tip of the endoscope.

➢ *Prolonged direct contact of light-emitting surfaces on tissue can lead to burns; this must therefore be avoided at all costs.*

⚠ WARNING:

Risk of infection: The flexible videoscopes are not sterile when delivered, and must be cleaned and disinfected/sterilized prior to initial use and each subsequent reuse.

⚠ WARNING:

Risk of injection: To avoid infection and maintain functionality, the chapter on reprocessing in the instruction manual must be observed. The cleaning, disinfection and sterilization procedures are explained therein in detail.

⚠ WARNING:

Risk of infection: Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding precleaning.

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⚠ WARNING:

The use of accessories other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased interference immunity of the flexible videoscope. The accessories listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories other than those specified here, it is the responsibility of the operator to ensure that they comply with EN/IEC 60601-1-2.

⚠ WARNING:

Risk of injury: The videoscope does not offer any electric insulation against high frequency (HF) voltages. Protection against HF voltage must be achieved using the appropriate HF instrument set.

➢ *Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.*

⚠ WARNING:

Combinations of medical devices are only assured to be safe if they are listed as such in the respective instruction manuals or the intended use and interface specifications of the products used in the combination permit this.

➢ *The instruction manuals and interface specifications of the medical devices used in combination must be observed precisely.*

⚠ WARNING:

The pressure compensation cap must be removed from the vent port for:

- endoscopic interventions,
- immersing in liquid, e.g., disinfectant solution and for machine reprocessing

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⚠ WARNING:

The pressure compensation cap must be in place on the vent port for:

- gas sterilization (EO, FO, hydrogen peroxide)
- transport and shipping

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⚠ WARNING:

Risk of injury: The flexible videoscope must not be used during the discharge of a cardiac defibrillator.

➢ *Remove the flexible videoscope from the operation area before discharge takes place.*

⚠ WARNING:

Risk of injury: Only KARL STORZ devices may be used in combination with the flexible videoscope. Disregarding this warning can endanger the patient and user.

***ATTENTION:** Risk of damage: Only use accessories (guide probes, cleaning brushes, etc.) from KARL STORZ. Accessories from other manufacturers may damage the videoscope.*

***ATTENTION:** Do not allow any live components or instruments to come into contact with the flexible videoscope.*

***NOTE:** Do not hit the distal tip of the flexible videoscope against hard objects, and never use a sharp object to remove dirt.*

***NOTE:** Never place heavy objects on the flexible videoscope.*

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NOTE: *Never clamp or fold the sheath of the flexible videoscope. This can result in damage to the videoscope.*

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NOTE: *The flexible videoscope must not be carried by the sheath. Do not pull, clamp or twist the sheath of the videoscope.*

NOTE: *The sheath of the flexible videoscope must never be clamped or coiled up too tightly. It is best to store the video videoscope in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39406AS).*

⚠ WARNING:

Risk of injury from clouding of the lens.

➢ *Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!*

⚠ WARNING:

When inserting the flexible videoscope into the patient, the distal tip of the flexible videoscope must be in the straight, non-deflected position.

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⚠ WARNING:

Risk of injury: Before inserting or removing the flexible videoscope, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

⚠ WARNING:

Risk of injury: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible videoscope must not be used.

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➢ *Place the distal tip of the flexible videoscope in the neutral, non-deflected position and remove the flexible videoscope slowly and carefully from the patient.*

NOTE: *The resistance that can be felt when deflecting the distal tip using the deflection lever increases as the angle of deflection increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.*

⚠ WARNING:

Risk of injury: The flexible videoscopes from KARL STORZ are not intended for use with light sources. The use of non-compatible devices can lead to patient injury or device damage.

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⚠ WARNING:

Risk of injury: Whenever the flexible videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged flexible videoscopes or damaged accessories must not be used.

NOTE: *There should be no instruments in the suction channel whilst suctioning is being performed. The instrument channel inlet should be sealed to prevent additional air from being drawn in.*

NOTE: *Do not squeeze the distal tip of the flexible videoscope too tightly.*

⚠ WARNING:

Risk of injury: Actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device.

➢ *Move the deflection lever smoothly and slowly!*

⚠ WARNING:

Risk of injury: Before each use, the functionality of the deflection lever and the connected distal tip must be checked.

➢ *Defective flexible videoscopes must not be used.*

⚠ WARNING:

Risk of injury: If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the device from the patient's body gently and carefully. If the flexible videoscope continues to be used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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⚠ WARNING:

Risk of injury: Only insert suitable accessories, e.g., biopsy forceps, into the instrument channel of the videoscope when the distal tip is in a neutral, non-deflected position.

- ➢ *Once the instrument has been introduced, carefully actuate the lever gently as far as the limit stop.*
- ➢ *The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device.*

⚠ WARNING:

Risk of injury: When HF electrodes or laser probes are used, it must be ensured that power is not applied until the active parts have left the instrument channel and are in the user's field of vision, and the site of the operation is clearly recognizable. Activation in the instrument channel may lead to serious injury to the patient and damage to the videoscope.

NOTE: *Do not apply force to insert accessories, e.g., biopsy forceps, through the working channel, to prevent damaging the device. It should be possible to move the accessory easily, and without any hindrance, through the entire channel.*

NOTE: *When the distal tip is deflected, do not insert instruments into the videoscope. Only insert the instrument when the distal tip is in the neutral position and deflect the deflection lever gently and carefully up to the limit stop.*

NOTE: *Inserting sharp-tipped accessories into the instrument channel may damage the videoscope.*

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⚠ WARNING:

Risk of injury: Foreign bodies and items to be extracted must not be removed via the instrument channel with the aid of instruments. Instead, the entire videoscope must be taken out of the body and the foreign body must then be removed.

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NOTE: *Accessories such as biopsy forceps that are kinked or bent will not function properly, are difficult to guide through the channel and may cause severe damage to it. They must therefore be replaced.*

⚠ WARNING:

Ensure that the working channel is completely dry prior to sterilization to avoid damage to the working channel.

⚠ WARNING:

Always unplug the flexible videoscope before carrying out any maintenance and cleaning procedures on it.

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6.2 Basic equipment

Flexible cystoscope HD VIEW 11272 VHK/VHU with following accessories:

Set 11272VHK

11272 VH	Flexible cystoscope HD VIEW
11301 P	Seal for suction valve seat
11014 L	LUER port with stopcocks
11014 DS	Working channel connector seal
11014 RA	Cleaning adaptor for working channel
11301 CD1	Irrigation adaptor

Set 11272 VHUK

11272 VHU	Flexible cystoscope HD VIEW
11301 P	Seal for suction valve seat
11014 LU	Double LUER port
11014 DS	Working channel connector seal
11014 RA	Cleaning adaptor for working channel
11301 CD1	Irrigation adaptor

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8 Preparation before use

1. Assemble LUER attachment. To do this, insert the lug of the LUER attachment into the recess on the working channel. Push and turn clockwise until it audibly clicks into place.

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
2. Assemble the suction valve/valve seal.

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9 Handling

NOTE: *Any damage to the flexible cystoscope HD VIEW resulting from incorrect operation is not covered by the manufacturer's warranty.*

9.1 Operation

The following instructions outline the correct operation of the flexible cystoscope HD VIEW and the corresponding accessories. This is not an introduction to the techniques of cystoscopy. Instructions on endoscopic technique can be found in the relevant medical literature.

9.2 Holding the device correctly

Grip the handle of the flexible cystoscope HD VIEW in your left or right hand, as shown in the picture. With this kind of control grip, the flexible cystoscope HD VIEW can be operated with one hand.

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9.3 Operating the suction valve

When the suction valve is assembled, a suction pump must be connected to the suction valve via a tube. Suction is activated upon pressing the valve button.

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- ➢ *Do not carry the flexible cystoscope HD VIEW by the sheath.*
- ➢ *Do not pull, clamp or twist the sheath of the flexible cystoscope HD VIEW.*
- ➢ *The sheath of the flexible cystoscope HD VIEW must not be kinked or coiled up too tightly.*
- ➢ *It is best to store the flexible cystoscope HD VIEW in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39403 AS).*

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If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the sheath of the flexible cystoscope HD VIEW out of the patient's body.

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9.5 Handling the deflection mechanism

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Deflection functions as follows with type 11272 VH endoscopes:

If the lever is pushed upwards, the tip is lowered.

If the lever is pushed downwards, the tip is deflected upwards.

This mechanism is inverted with model 11272 VHU.

ATTENTION: *When the flexible videoscope is inserted into the patient, the distal tip must be in a straight, non-deflected position.*

⚠ WARNING:

The deflection lever must never be operated with force or abruptly to avoid injuring the patient or damaging the device!

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⚠ WARNING:

Risk of injury and damage from faulty deflection mechanisms: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible videoscope must not be used.

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⚠ WARNING:

Risk of injury: If the flexible videoscope malfunctions during use on the patient, discontinue use immediately.

➢ *Place the distal tip of the flexible videoscope in the neutral, non-deflected position and remove slowly and carefully from the patient.*

⚠ WARNING:

Risk of damage: The lens at the distal end of the flexible videoscope may become detached, e.g., due to mechanical forces.

➢ *For this reason, it is essential to check that the lens is undamaged and present before and after every decontamination.*

9.6 Handling the working channel

Do not force instruments or cleaning brushes through the channels This can lead to damage to the working channel, especially if the tip of the instrument is deflected.

When connecting and operating suction/ irrigation devices, pay attention to the instruction manuals for the particular devices. The irrigation devices recommended by KARL STORZ are set to a maximum pressure of 54 kpa (400 mmHg/0.54 bar/7.8 psi).

- The LUER adaptor 27014 Y or the seal 27550 C can be used to seal the instrument channel when inserting instruments.
- Instruments can be introduced through the instrument channel of the flexible cystoscope HD VIEW. Make sure that the accessories are of the appropriate length and diameter. Unsuitable accessories can damage the flexible cystoscopes HD VIEW. If you have any questions concerning accessories, please contact KARL STORZ directly or an authorized specialist dealer.

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Handling

- The flexible cystoscope HD VIEW must be in the straight, non-deflected position when inserting an instrument into the working channel.

Forceps, scissors and other accessories whose jaw can be opened or closed must always be inserted into the instrument channel with closed jaws. The jaws may only be opened when the distal tip of the instrument is completely clear of the channel. Opening the jaws in the instrument channel can cause damage to the videoscope.

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ATTENTION: *Risk of damage: Only use accessories (guide probes, cleaning brushes, etc.) from KARL STORZ. Accessories from other manufacturers may damage the videoscope.*

9.7 Handling the accessories (forceps, cleaning brushes, HF electrodes etc.)

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9.8 Inspecting the surfaces

⚠ WARNING:

Risk of injury: The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

***NOTE:** Do not squeeze the distal tip of the flexible videoscope too tightly.*

Before each use, verify that the device is in good condition by checking the surface of the flexible cystoscope HD VIEW, especially the sheath area and the distal tip, for sharp edges. Cystoscopes HD VIEW with sharp edges must not be used.

Otherwise they may cause serious injury to the patient or user.

Check the connecting cable for damage (such as fractures, cracks, twisted or crushed areas, etc.) and the connecting plug for corrosion.

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⚠ WARNING:

Risk from damage: Faulty flexible videoscopes must not be used.

➢ *Before each use, check the functionality of the deflection lever and the distal tip.*

The distal tip is deflected by actuating the deflection lever slowly.

The distal end will deflect in accordance with the direction of actuation (forward or backward).

9.10 Connecting to IMAGE1 S™

The flexible cystoscope HD VIEW is connected to the camera control unit for operation.

For operation, an IMAGE1 S™ system with an IMAGE1 S Connect® and an IMAGE1 X-Link are required. Both camera modules must be connected to each other.

The flexible cystoscope HD VIEW is then connected to the IMAGE 1X-Link.

9.11 Focus

The image focus is adjusted automatically. Manual focusing is not possible/necessary.

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ATTENTION: *Risk from damage to the device: Flexible endoscopes must never be steam sterilized (autoclaved).*

ATTENTION: *Risk from damage to the device: To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.*

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NOTE: *Damage to the device. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be found at www.karlstorz.com*

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Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in 'Leakage test – dry' and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

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ATTENTION: *In this case, fill out the accompanying repairs form (see Appendix or Art. no. 96096004DER), mark the endoscope, which has been cleaned by wiping to remove organic soiling/residue, as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the transport protection sheet 13990 SFN.*

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10.5 Disassembly

The flexible cystoscope HD VIEW must be completely dismantled for cleaning/disinfection.

1. Dismantle the suction valve by pressing and turning.

2. Dismantle the LUER attachment by pressing and turning

3. Dismantle the stopcock (if present) on the LUER attachment.

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
5. Clean the channels and lumina of the LUER attachment with a suitable brush.

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6. Dismantle the reusable suction valve 11301 CB1 by unscrewing the metal cap.

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8. Brush the working channel connector.

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9. Brush the working channel and suction channel.

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
10. Connect the cleaning adaptor 11301 CD1 to the valve seat.

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10.6 Manual disinfection

The endoscope must be completely immersed in a cleaning solution. The suction and working channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

10.6.1 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The endoscope can be disinfected with the following high level disinfection solutions which have been tested by way of example:

- 2.4% glutaraldehyde solution (e.g., CIDEX® Activated Dialdehyde Solution)
- 0.55% ortho-phthalaldehyde solution (e.g., CIDEX® OPA Solution)
- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

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10.7 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’ (Art. no. 96216003DER):

10.7.1 Manual preparation

The endoscope must be completely immersed in a cleaning solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channels are brushed until the brush head emerges and is visibly clean. Before retracting the brush head that emerges at the distal end, clean it. Never brush back and forth within the channel.

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NOTE: *KARL STORZ recommends using the adaptors 11301CD1 and 11014RA for simpler filling of the lumina.*

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10.7.2 Connecting

In order to ensure effective machine cleaning and disinfection, the endoscope's working and suction channel must be connected to the washer and disinfector using the recommended adaptors 11301 CD1 and 11014 RA.

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
10.7.3 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

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11 Reprocessing of LUER port

⚠ WARNING:

Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments (art. no. 96216003DER)’ and the accompanying documentation.

⚠ WARNING:

Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

⚠ WARNING:

When carrying out any work on contaminated medical devices, the guidelines of the Employers’ Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

***ATTENTION:** When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.*

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ATTENTION: *Danger of damage to medical devices: The use of chemicals that have not been approved by KARL STORZ or chemical manufacturers may cause damage to the medical devices. Only use chemicals approved for the medical device for reprocessing. A helpful guide is available on the Internet at www.karlstorz.com.*

NOTE: *National laws and regulations must be observed.*

NOTE: *The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.*

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Brushes: 27652, 27650B

Adaptor: 39107ALK (only for 11014L)

11.2 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

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Brushing the lumina

The working channels, lumina and hollow spaces must be precleaned under cold running water with the aid of a suitable brush.

Rinsing out with a pressure water gun

The working channel, lumina and hollow spaces must be rinsed using a pressure water gun and a suitable rinsing attachment. For rinsing, pulse at least 3 times for 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

11.4 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’ (Art. no. 96216003):

Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C must not be exceeded. The relevant national regulations and the chemical manufacturer’s instructions must be taken into account when using this method.

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

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11.5 Connecting

In order to ensure effective machine cleaning, disinfection and rinsing, the instrument must be connected in the washer and disinfector.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Joints, threads and glide surfaces must be treated locally with instrument oil and/or special lubricant. (only valid for 11014L)
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

NOTE: *The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).*

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Steam sterilization using the fractionated prevacuum procedure

The medical device must be sterilized in its fully assembled state using the fractionated pre-vacuum procedure (DIN EN ISO 17665-1) at 132°C–137°C with a minimum exposure time of 3 minutes and a maximum of 18 minutes. This procedure is only suitable for thermostable instruments.

NOTE: *Greased components must be disassembled for sterilization, so that steam penetration is ensured.*

11.9 Limits of reprocessing

The end of the product’s service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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13 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts only.

14 Maintenance

Preventive maintenance is not essential. However, regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be requested from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

15 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this instrument only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by persons authorized by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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The medical device must always be sent to your local subsidiary (see the section 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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17 Directive compliance

In accordance with the Medical Device Directive (MDD):

11272 VH/VHU:
This medical device belongs to Class II a

This medical device bears a CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

18 Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/devices to the sender.

Repairs, modifications or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

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19 Disposal

This device has been marked in accordance with the European Directive 2012/19/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the device as electronic scrap. Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary, or your authorized dealer about the collection point in your area. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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20 Electromagnetic Compatibility user information

The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment.

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⚠ WARNING:

The use of an accessory or cable with the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ other than those specified in this manual may result in increased emissions or decreased immunity of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™. When using an accessory or cable with the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ other than those specified in this manual, it becomes the responsibility of the user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ to determine compliance with IEC 60601-1-2 when using this item.

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Electromagnetic Compatibility EMC

⚠ WARNING:

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

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The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in an electromagnetic environment specified below. The user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines lines ± 1 kV for input and output lines 100 kH repetition	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	Voltage Dip: Dip to 0% for 1 cycles @ 0°phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100% for 250/300 cycles	Voltage Dip: Dip to 0% for 1 cycles @ 0°phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be me asured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3V/m 80 MHz to 2.7 GH	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	

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The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbances acc. to IEC 61000-4-6	3 V_{rms} 150 kHz to 80 MHz	3 V_{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2\sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey[a] should be less than the $d = 1.2\sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz compliance level in each frequency range[b]. Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

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[a] Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is used exceeds the above compliance levels, the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
[b] Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 Emission class and group		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ uses RF energy only for its internal function. The customer or user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5
Recommended separation distances between portable and mobile RF telecommunication devices and the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™

The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ as recommended below, according to the maximum output power of the communications equipment.

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KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang Intenational Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: phillippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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INSTRUCTION MANUAL
Flexible Video Esophagoscope –
Series 13303 E, Monitor 8403 ZX

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It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and, where required, sterilization of the Video-Esophagoscope. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

WARNING: *KARL STORZ instruments are provided unsterilized, and must be cleaned, disinfected and/or sterilized prior to the initial use and before each subsequent use.*

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be deduced from the information, figures and descriptions in this manual.

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Contents

1 Pictures of the instruments 4
2 Symbol description 6
2.1 Pictograms 6
2.2 Symbols for safety instructions 8
3 Intended use 9
3.1 Indications 9
3.2 Contraindications 9
4 User qualification 10
5 Safety information 11
6 Unpacking the device 15
6.1 Basic equipment 15
6.2 Optional accessories 15
6.3 Operating conditions 15
6.4 Electromagnetic compatibility 15
7 External system components ... 16
7.1 Charger 16
7.2 Country specific configuration 16
8 Handling 17
8.1 Operating 17
8.2 Holding correctly 17
8.3 Handling the shaft 18
8.4 Handling the deflection mechanism 19
8.5 Connection to the monitor 8403 ZX 20
8.6 Use of the HDMI output 21
8.7 Functions of the C-MAC® Monitor 8403 ZX 21
8.8 Focus 21
8.9 Function buttons of the endoscope 22
8.10 Monitor function keys 22
8.11 Energy management 35
8.12 Fixing the monitor to a stand 37
9 Preparation prior to use 38
9.1 Inspecting the surfaces 38
9.2 Checking the deflection 39
9.3 Switching the Video-Esophagoscope on/off 40

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Contents

9.4 Inspecting the image quality....41
9.5 Switching off the C-MAC® Monitor41
10 System reset....42
11 Monitor bag 8403 YD....42
12 Disassembly42
13 Reprocessing43
13.1 Overview of reprocessing procedures 43
13.2 Accessories....45
13.3 Preliminary cleaning directly after use 45
13.4 Leakage test....46
13.4.1 Performing the leakage test – dry....46
13.4.2 Performing the leakage test – in liquid46
13.5 Manual reprocessing....47
13.5.1 Manual cleaning47
13.5.2 Manual disinfection48
13.5.3 Manual High Level Disinfection (HLD) .. 48
13.5.4 Drying48
13.6 Automated reprocessing....49
13.6.1 Manual preparation49
13.6.2 Machine cleaning/chemical-thermal disinfection49
13.7 Assembly, inspection and care....49
13.8 Packaging....50
13.9 Sterilization....51
13.9.1 Hydrogen peroxide (H_2O_2) sterilization - ASP STERRAD®....51
13.9.2 Hydrogen peroxide (H_2O_2) sterilization – STERIS® AMSCO® V-PRO™....52
13.9.3 Ethylene oxide sterilization (EO)....52
13.9.4 Low-temperature steam and formaldehyde process (LTSF)....53
13.10 Limits of reprocessing....53

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Contents

14 Storage 53
15 System description 55
16 Technical description 56
17 Servicing and repair 58
18 Maintenance 58
19 Limitation of liability 58
20 Warranty 59
21 Standard compliance 59
22 Directive compliance 60
23 Infection prevention 60
24 Disposal 60
25 Appendix 61
25.1 Troubleshooting 61
25.2 Accessories 63
25.3 Electromagnetic Compatibility (EMC) Information 64
26 Subsidiaries 84

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

① Connecting cord
② Connection for leakage tester or pressure compensation cap (vent port)
③ Distal tip
④ Shaft
⑤ Handle
⑥ Deflection lever
⑦ Leakage tester 13242 XL
⑧ Pressure compensation cap 11025 E
⑨ Monitor
⑩ Connection 2 for video endoscope
⑪ Connection 1 for video endoscope
⑫ Mains connection
⑬ Reset button
⑭ Holder for SD memory card
⑮ Status bar
⑯ Pilot lamp
⑰ On/Off switch
⑱ Pilot lamp charging (orange)
⑲ Softkeys
⑳ Function symbols
㉑ USB port
㉒ HDMI port
㉓ Video recording Start/Stop
㉔ White balance
㉕ Still image recording

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2 Symbol description

2.1 Pictograms

Refer to instruction manual
White balance
Single image capture
Video recording
This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
Device type BF
Electronic information product pollution control (China RoHS)
Manufacturer
UL approval mark, valid in the USA and Canada
Degree of protection in accordance with International Protection Codes: Protection against contact with tools or similar and small foreign bodies, protection against vertically dripping water
Rated ambient temperature/transducer temperature class
DC voltage connection
Consult instructions for use

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2.2 Symbols for safety instructions

WARNING: *A Warning indicates that the personal safety of the patient or personnel may be endangered. Failure to observe a Warning could result in injury to the patient or personnel.*

CAUTION: *A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument set.*

NOTE: *A Note indicates special information regarding operation or clarifies important information.*

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3 Intended use

The video esophagoscope is used for endoscopic imaging and image transmission, thus enabling an endoscopic diagnosis of body cavities and hollow organs in esophagoscopy. It is also used for inspection of the upper and lower respiratory tracts, endotracheal intubation as well as inspection of the tube position during PCT. The video esophagoscope is designed for short-term use in invasive procedures through natural orifices and surgical invasive interventions.

3.1 Indications

The use of the instrument set for esophagoscopy is indicated if esophagoscopy is indicated in the opinion of the attending physician. This includes, amongst other things, the diagnosis and treatment of:

- Removal of foreign bodies from the esophagus
- Eosinophilic esophagitis
- Stenoses
- Diagnosis of disorders of the esophagus, e.g., esophagitis, tumors, bleeding
- Biopsy sample
- Endoscopic ablation with esophagus tumors/ squamous cell carcinoma of the hypopharynx and the upper third of the esophagus

3.2 Contraindications

The use of the instrument set for esophagoscopy and diverticuloscopy is contraindicated if, in the opinion of the responsible physician, the health of the patient is endangered through such use, due to the general condition of the patient or if the endoscopic method itself is contraindicated.

Esophagoscopes and instrument sets must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

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- General anesthesia
- Coagulation disorders
- ASA > 3
- Multiple diverticula
- Esophageal stricture

4 User qualification

The KARL STORZ Video-Esophagoscope may only be used by persons with an appropriate medical qualification and who are acquainted with the esophagoscopy technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the Video-Esophagoscope. It is not intended as an introduction to the technique of esophagoscopy.

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Safety information

Risks of infection

- Warning: The Video-Esophagoscope is delivered unsterilized and must be cleaned and disinfected/sterilized before initial and subsequent use.
- Warning: In order to prevent infection and protect personnel, it is strictly forbidden to ship contaminated medical devices. To avoid contact and aerogen infections, medical devices must be decontaminated on-site.

Risks due to electric energy

- Caution: Do not allow any live components or instruments to come into contact with the Video-Esophagoscope.
- Warning: Danger of glare and burns! Never look directly into the light of the Video-Esophagoscope. The LED light radiated through the Video-Esophagoscope may cause glare or high temperatures at the distal light outlet and therefore lead to burns.
- Warning: Danger of explosion if flammable gases are used in the immediate proximity of the Video-Esophagoscope. The Video-Esophagoscope must not be used in the presence of flammable or combustible media.
- Warning: Always unplug the Video-Esophagoscope before carrying out any maintenance and cleaning work on it.
- Warning: Danger of electric shock! Refer servicing for this device to authorized personnel only. Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims. Exceptions are those that are described in the instruction manual.

Risks due to damaged parts

- Warning: Every time the Video-Esophagoscope is used, you must first check that it, and any accessories used in combination with it, is in perfect condition. Damaged Video-Esophagoscopes or damaged accessories must not be used.

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Safety information

- Warning: The Video-Esophagoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for Video-Esophagoscopes that are repaired by unauthorized personnel.

Risks due to combination of system components

- Warning: Combinations of medical devices are only then assured to be safe if
 - they are identified as such in the respective instruction manual or
 - the intended purpose and the interface specifications of the devices used in combination permit this.
- Warning: Only use accessories that are specifically suited for use with Video-Esophagoscopes and can be employed in line with generally accepted technical standards and safety guidelines.
- Caution: Only use accessories from KARL STORZ. Accessories from other manufacturers may damage fiberscopes and the flexible videoscopes.
- Warning: Only KARL STORZ devices (e.g. monitor 8403 ZX) may be used. Disregarding this warning can endanger the patient and user.
- Warning: The Video-Esophagoscope must not be used during the discharge of a cardiac defibrillator. Remove the Video-Esophagoscope from the operation area before discharge takes place.
- Warning: When in operation, the integrated LED light source heats up the Video-Esophagoscope. Switch off the Video-Esophagoscope before putting it down and do not place it in the proximity of the patient (as per IEC 60601-1).
- Warning: Damage to the Video-Esophagoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

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Safety information

Damage to the instrument

- Warning: If the Video-Esophagoscope should malfunction during application on a patient, stop the application immediately. Place the distal tip of the Video-Esophagoscope in the neutral, nondeflected position and remove the Video-Esophagoscope slowly and carefully from the patient.
- Caution: Do not steam sterilize (autoclave) the Video-Esophagoscope under any conditions.
- Caution: The maximum permissible temperature for decontamination, sterilization, transport and storage is 65 °C.
- Caution: Do not place the Video-Esophagoscope on devices that become very hot. Do not expose it to direct sunlight or excessive heat, to avoid damage to the device.
- Caution: Do not hit the distal tip of the Video-Esophagoscope against hard objects, and never use a sharp object to remove dirt.
- Caution: Never lay heavy objects on the Video-Esophagoscope.
- Caution: To avoid damaging the monitor, do not remove the SD card or connecting cord while a video is recording.
- Caution: To avoid damaging the monitor, always switch it off before replacing the SD card.

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8 Handling

NOTE: *Damage to the Video-Esophagoscope resulting from incorrect operation is not covered by the manufacturer's warranty.*

WARNING: *The Video-Esophagoscope is delivered unsterilized, so it has to be reprocessed prior to initial use and before each subsequent reuse.*

CAUTION: *Do not store or disinfect the flexible Video-Esophagoscope together with sharp objects in order to avoid causing damage to the outer casing.*

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8.1 Operating

The following instructions outline the correct operation of the Video-Esophagoscope and the corresponding accessories. This is not an introduction to the techniques of esophagoscopy. Instructions on the endoscopic technique can be found in the relevant medical literature.

The Video-Esophagoscope is intended for visualisation and diagnostic purposes.

8.2 Holding correctly

Grip the handle of the Video-Esophagoscope in your left or right hand, as shown in the picture. With this kind of control grip, the Video-Esophagoscope can be operated with one hand.

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8.3 Handling the shaft

CAUTION: *Do not carry the Video-Esophagoscope by the shaft. Do not pull, clamp or twist the shaft of the Video-Esophagoscope.*

CAUTION: *The shaft of the Video-Esophagoscope must never be kinked or coiled up too tightly. It is best to store the Video-Esophagoscope in a suspended position or to keep it in a tray specifically designed for sterilizing and storage (39405 AS).*

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If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the shaft of the Video-Esophagoscope out of the patient's body.

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8.4 Handling the deflection mechanism

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WARNING: *The distal tip of the Video-Esophagoscope must be straight, not angled, when the videoscope is inserted into the patient*

WARNING: *Before removing the Video-Esophagoscope from the patient, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.*

WARNING: *If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the Video-Esophagoscope must not be used.*

WARNING: *If a functional defect occurs with the Video-Esophagoscope during use on the patient, discontinue use immediately. Put the distal tip of the Video-Esophagoscope into the neutral, non-deflected position and remove the Video-Esophagoscope slowly and carefully from the patient.*

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WARNING: *The lens at the distal end of the Video-Esophagoscope may become detached, e.g. due to mechanical forces. For this reason, it is essential to check that the lens is undamaged and present before and after decontamination.*

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8.5 Connection to the monitor 8403 ZX


CAUTION: *The Video-Esophagoscope can be used together with the KARL STORZ monitor 8403 ZX. The instruction manual for the monitor must be observed.*

WARNING: *It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces, or voltage-carrying elements of other units, at any time.*

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Insert the connecting cord ① of the flexible Video-Esophagoscope into either connection socket ⑩ or ⑪ on the monitor.

The connecting cord ① can be connected or disconnected during operation. Insert the cable from the power supply unit into the connection socket ⑫ of the C-MAC® Monitor and connect the power supply unit with the mains socket. The power cord can be fixed to the connection socket using the union nut.

The battery is charged automatically as soon as the C-MAC® Monitor is connected to the power supply. This will be indicated by the lighting up of the orange colored LED ⑱.

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8.6 Use of the HDMI output

The device features a port for an external monitor. An HDMI cable can be used to send the video signal to an external monitor, if necessary.

CAUTION: *The C-MAC® Monitor must be switched off when an external monitor is connected.*

NOTE: *The connected monitor must support a resolution of 1024 × 768 or 1280 × 720.*

NOTE: *To connect the C-MAC® Monitor to a KARL STORZ monitor, the HDMI-DVI adaptor with article number **20**9190 10 is required.*

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8.7 Functions of the C-MAC® Monitor 8403 ZX

A more detailed description of the operation and functions of the monitor 8403 ZX can be found in the corresponding instruction manual supplied with every monitor.

In case of loss of the instruction manual, you can also download it from our homepage at www.karlstorz.com.

WARNING: *The KARL STORZ flexible Video-Esophagoscope can not be used in combination with camera control units and light sources. The use of non-compatible units can lead to patient injury or unit damage.*

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8.11 Energy management

The charge state of the battery is displayed on the monitor by the charge indicator in the menu bar.

The symbols mean:

- Battery empty
- 25 %: Approx. 30 minutes of operating time
- 50 %: Approx. 1 hour of operating time
- 75 %: Approx. 1.5 hours of operating time
- 100 %: Approx. 2 hours of operating time
- Battery is being charged
- Battery must be replaced (Maximum battery capacity is less than 50 % of the original maximum capacity)

NOTE: *The C-MAC® Monitor can be powered using either mains electricity or the integrated battery.*

The battery is charged automatically as soon as the monitor is connected to the power supply. This is displayed by the lighting up of the orange-colored LED.

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
Handling

NOTE: *The monitor should be recharged after each use or, at the latest, once the charge state has dropped below 50 %.*

CAUTION: *Damage to the power cord: To remove the power cord from the socket, pull only on the plug housing, never on the cord itself!*

The lithium-ion battery can be used for approx. 2 hours when charged.

The discharged battery is fully charged after 5 hours (when the device is switched off). During the charging time, the orange-colored LED lights up.

If the system is not used for 12 weeks or longer, you should check the charge state using the progress bar described on page 35 and, if necessary, charge.

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8.12 Fixing the monitor to a stand


Attach the stand holder to the rear of the monitor using 4 screws.

The VESA 75 holder (8401 YC), which is included in the scope of supply, is suitable for securing devices to tubes with a diameter of 25 mm.

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Press the holder together to allow the stand bar to be pushed through the openings.

Release the holder in the desired position. The selected position is fixed by the holder's basic tension.

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The optionally available clamp 8401 YB (see accessories) can also be secured to the VESA 75 holder and used to attach the monitor to tubes with a diameter of less than 25 mm and/or to rails with a cross section of less than 25 mm.

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
WARNING: *Every time the Video-Esophagoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged Video-Esophagoscopes or damaged accessories must not be used.*

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CAUTION: *Do not squeeze the distal tip of the Video-Esophagoscope too tightly.*

Every time the Video-Esophagoscope is used, you must first check that its surface, and above all the surface of the shaft and distal tip, are in perfect condition and that there are no sharp edges. Video-Esophagoscope with sharp edges must not be used. Otherwise
serious injuries may be caused to the patient or the user.

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9.3 Switching the Video-Esophagoscope on/off

Switching on/off: The Video-Esophagoscope with the integrated LED light source is switched on by pressing the power button ⑰ on the monitor. Pressing the key a second time switches the unit off.

WARNING: *Danger of glare: Do not look directly into the light beam.*

WARNING: *When in operation, the integrated LED light source heats up the Video-Esophagoscope. Switch off the Video-Esophagoscope before putting it down and do not place it in the proximity of the patient (as per IEC 60601-1).*

NOTE: *The function keys ⑰ on the monitor are used for switching the unit on/off. The Video-Esophagoscope itself does not have any function keys.*

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Video recording: When the control key ‘Video recording’ on the C-MAC® Monitor or the flexible Video-Esophagoscope is pressed, the flexible Video-Esophagoscope records a video stream on the SD memory card. Pressing the key a second time stops the recording.

NOTE: *The video stream can only be viewed on a PC if an MPEG 4-codec is installed on the PC.*

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Single image capture: When the control key ‘Single image capture’ on the C-MAC® Monitor or the flexible Video-Esophagoscope is pressed, the current image on the monitor is saved on the SD memory card.

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
9.4 Inspecting the image quality

WARNING: *If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the Video-Esophagoscope from the patient's body gently and carefully. If the Video-Esophagoscope is used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.*

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9.5 Switching off the C-MAC® Monitor

Standby: Briefly pressing the power button ⑰ will set the C-MAC® Monitor to standby operation mode. In the case of brief interruptions, this enables the C-MAC® Monitor to be ready for operation again quickly. The standby mode is indicated by a flashing green pilot lamp ⑯ above the power button.

After 5 minutes, the system switches itself from standby mode to the off mode.

Power off: Pressing the power button ⑰ for more than approx. 3 seconds will switch off the C-MAC® Monitor. The green pilot lamp ⑯ will go out.

NOTE: *When the C-MAC® Monitor is switched off, the integrated light source in the flexible Video-Esophagoscope also switches off.*

NOTE: *Information on the operation and functions of the C-MAC® Monitor 8403 ZX can be found in the corresponding instruction manual supplied with every C-MAC® Monitor (Art. No. 96076008DER). In case of loss of the instruction manual, you can also download it from our homepage at www.karlstorz.com.*

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This reprocessing step is part of the validated procedure

❍ This reprocessing step may be carried out additionally in order to aid the cleaning. The sterilization procedure is approved in respect of material compatibility; however, it must be validated by the operator on site in respect of effectiveness.

X On the grounds of potential damage, this reprocessing step must NEVER be performed

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pressure compensation cap, manual leakage tester

Place of application	Art. number brush/ adaptor
Clean surface	27652
Ventilation connector	13242 XA, 11025 E

13.3 Preliminary cleaning directly after use

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Wipe-clean the exterior surfaces of the endoscope with a moistened cloth.

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ⓘ ***NOTE:*** *Do not use disinfectants for preliminary cleaning, these cause fixation of proteins and organic residues.*

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13.4.2 Performing the leakage test – in liquid

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Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in ‘Leakage test – dry’ and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously)

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from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

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ⓘ ***NOTE:*** *Never release the pressure in the liquid.*

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❗ ***CAUTION:*** *In this case, fill out the accompanying repairs form (see Appendix or Art. no. 96096004DER), mark the endoscope, which has been cleaned by wiping to remove organic soiling/ residue, as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the transport protection sheet 13990 SFN.*

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CAUTION: *Never remove dirt with a needle or other sharp objects, since these can cause severe damage.*

CAUTION: *The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.*

13.5.2 Manual disinfection

The endoscope must be completely immersed in a cleaning solution. At the end of the exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. For this, water of the best possible quality must be used, taking into account the national regulations. Sterile filtered water must be used for the final rinsing.

13.5.3 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

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13.5.4 Drying

Finally, the surfaces are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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CAUTION: *The pressure must not exceed 100 kPa (1 bar / 14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.*

13.6 Automated reprocessing

The following procedure for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual “Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments”:

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13.6.2 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C and a maximum pressure of 100 kPa (1 bar / 14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer‘s instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out, must take place in consultation with the manufacturer of the device.

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he procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

The following sterilization procedures have been validated and approved by KARL STORZ for this medical device:

13.9.1 Hydrogen peroxide (H_2O_2) sterilization - ASP STERRAD®

WARNING: *Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.*

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ⓘ ***NOTE:*** *Refer to the ‘STERRAD® Sterility Guide’ to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.*

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- STERRAD® 100S Long Cycle with Boosters
- STERRAD® 100S Short Cycle with Boosters
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

13.9.2 Hydrogen peroxide (H_2O_2) sterilization – STERIS® AMSCO® V-PRO™

ⓘ ***NOTE:*** *The ‘STERIS Low Temperature Reprocessing Device Compatibility Matrix’ can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.*

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13.9.3 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 60 minutes.

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The following procedures are approved in respect of material compatibility, but must be validated by the on-site operator with regard to efficacy:

13.9.4 Low-temperature steam and formaldehyde process (LTSF)

KARL STORZ had only the specified sterilization method validated with regard to a reduction in microorganisms on the order of 10^{-6} SAL. However, due to the large number of processes used in the global market, validation in respect of efficacy was not carried out. Therefore, the user must validate his specific processes individually in collaboration with the device manufacturer.

WARNING: *Risk of infection: If FO sterilization is performed, the method must be validated on site.*

13.10 Limits of reprocessing

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

14 Storage

CAUTION: *Do not expose the Video-Esophagoscope to direct sunlight or excessive heat.*

CAUTION: *The distal tip of the Video-Esophagoscope must be in the neutral position during storage.*

CAUTION: *Do not keep the Video-Esophagoscope in the carrying case.*

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Reprocessing

CAUTION: *The shaft of the Video-Esophagoscope must never be kinked or coiled up too tightly. It is best to suspend the Video-Esophagoscope for storage or to keep it in a large container. Video-Esophagoscope should be packaged and stored in sterile cloths or suitable containers in accordance with EN 868, and must be absolutely dry before being put into storage. It is most suitable to store the Video-Esophagoscope suspended vertically with the distal end pointing down.*
The room where they are stored should be at a normal temperature and humidity.
Store dismantled valves and caps separately.

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15 System description

The flexible Video-Esophagoscope from KARL STORZ has a 750 mm pellethane working shaft, 2.90 mm in diameter. The distal tip made of stainless steel is 2.9 mm in diameter and houses the video module. The distal tip of the flexible Video-Esophagoscope can be deflected 210° up or 130° down. The handle houses the deflection lever, the function keys and the LED light source.

The flexible Video Esophagoscope from KARL STORZ is connected to the recommended monitor 8403 ZX using the integrated 2 m supply cable. Images and videos can be saved using a SD memory card and the control keys on the endoscope or on the monitor, which are assigned, for example, the functions of video recording and image capture.

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21 Standard compliance

According to IEC 60601-1, 60601-2-18, UL 60601-1, CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electric shock:
 - In combination with the monitor 8403 ZX: Internally powered during battery operation, Protection class II during line operation
 - In combination with the C-HUB®: Protection class II
- Degree of protection against electric shock: Applied part type BF 
- Type of protection against moisture: Protection against continuous immersion: IPX 8

According to IEC 60601-1-2:

Please observe the Electromagnetic Compatibility information in the appendix (pp. 64–80).

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13303 E: This medical product belongs to Class IIa
8403 ZX: This medical product belongs to Class I

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24 Disposal

This unit has been marked in accordance with the European Directive 2012/19/EC on waste electrical and electronic equipment (WEEE). At the end of its service life, dispose of the device as electronic waste. Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information about your local collection point.
Within the scope of application of this directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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25.3 Electromagnetic Compatibility (EMC) Information

The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment.

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WARNING: *The use of an accessory or cable with the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX other than those specified in this manual may result in increased emissions or decreased immunity of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX. When using an accessory or cable with the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX other than those specified in this manual, it becomes the responsibility of the user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX to determine compliance with IEC 60601-1-2 when using this item.*

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: *Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.*

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The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is intended for use in an electromagnetic environment specified below. The user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines lines ± 1 kV for input and output lines 100 kH repetition	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	Voltage Dip: Dip to 0% for 1 cycles @ 0°phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100% for 250/300 cycles	Voltage Dip: Dip to 0% for 1 cycles @ 0°phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be me asured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3V/m 80 MHz to 2.7 GH	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	

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The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbances acc. to IEC 61000-4-6	3 V_{rms} 150 kHz to 80 MHz	3 V_{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2\sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey[a] should be less than the $d = 1.2\sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz compliance level in each frequency range[b]. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

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[a] Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is used exceeds the above compliance levels, the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
[b] Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 **Guidelines and manufacturer's declaration – electromagnetic emissions**		
The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is intended for use in the electromagnetic environment specified below. The customer or user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX uses RF energy only for its internal function. The customer or user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5
Recommended separation distances between portable and mobile RF telecommunication devices and the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX

The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1{,}2\sqrt{P}$	80 MHz to 800 MHz $d = 1{,}2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2{,}3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE: *At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.*

NOTE: *These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.*

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment.

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this medical device is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality piece of equipment from KARL STORZ.

This instruction manual is intended to aid you in the proper installation, connection, and operation of this medical device.

- Failure to observe and follow this instruction manual and the instruction manuals of medical devices used in combination can result in injury to patients, users and third parties as well as damage to the medical device.
- Please read all relevant instruction manuals and always follow the instructions given. Check the functioning of the medical devices used in combination.
- Take note of all the instructions on the device.
- Keep this instruction manual available for ready reference in a conspicuous location near the device.

KARL STORZ reserves the right to make engineering modifications, in particular in the context of ongoing development and improvement of our devices.

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Contents

1 Symbols employed ... 3
1.1 Warnings and cautions ... 3
1.2 Pictograms/packaging symbols ... 4
2 Pictures of the instruments ... 5
3 Intended use ... 7
3.1 Intended use ... 7
3.2 Indications ... 7
3.3 Contraindications: ... 7
3.4 User qualification ... 8
4 Safety instructions ... 9
5 Initial operation ... 11
5.1 Unpacking ... 11
5.2 Basic equipment ... 11
5.3 Operating conditions ... 11
6 Handling ... 12
6.1 General Warnings ... 12
6.2 Operation ... 13
6.2.1 Holding correctly ... 13
6.2.2 Handling the sheath ... 13
6.2.3 Handling the deflection mechanism ... 15
6.3 Connection to the camera contol unit ... 17
6.3.1 Compatible camera control units ... 17
6.3.2 Connecting the videoscope ... 17
6.3.3 Inspecting the LED function ... 18
6.4 Focus ... 19
6.5 Function buttons ... 19
7 Reprocessing ... 20
7.1 Overview ... 20
7.2 General information ... 21
7.3 Accessories ... 22
7.4 Precleaning directly after use ... 22
7.5 Leakage test ... 23
7.5.1 Performing the leakage test – dry ... 23
7.5.2 Performing the leakage test – in liquid ... 23
7.6 Manual wipe-down disinfection ... 24
7.7 Manual reprocessing ... 24
7.7.1 Manual cleaning ... 24
7.7.2 Manual disinfection ... 25
7.7.3 Manual High Level Disinfection (HLD) ... 25
7.7.4 Drying ... 25
7.8 Machine reprocessing ... 26

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Contents

7.8.1 Manual preparation ... 26
7.8.2 Machine cleaning/chemical-thermal disinfection ... 26
7.9 Assembly, inspection and care ... 26
7.10 Packaging ... 27
7.11 Sterilization ... 27
7.11.1 Hydrogen peroxide (H_2O_2) sterilization – ASP STERRAD® ... 28
7.11.2 Hydrogen peroxide (H_2O_2) sterilization – STERIS® V-PRO™ ... 29
7.11.3 Chemical low-temperature sterilization with peracetic acid – STERIS System 1/1E® ... 29
7.12 Limits of reprocessing ... 29
7.13 Storage ... 29
8 System description ... 31
9 Technical data ... 32
9.1 Operating/storage conditions ... 32
10 Service, repair and disposal ... 33
10.1 Warranty ... 33
10.2 Maintenance ... 33
10.3 Repairing device parts ... 33
10.4 Repair program ... 34
10.5 Infection prevention ... 34
10.6 Disposal ... 34
10.7 Directive compliance ... 34
11 Appendix ... 35
11.1 Troubleshooting ... 35
12 Electromagnetic Compatibility user information ... 36
13 Subsidiaries ... 55

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1 Symbols employed

1.1 Warnings and cautions

HAZARD SYMBOL: *This is a hazard symbol. It warns of the risk of injury. To prevent injuries, all the measures marked with the hazard symbol must be followed.*

Signal words:

WARNING! designates a possible imminent risk. If this is not avoided, it **could** lead to death or serious injuries.

ATTENTION! designates a possibly harmful situation. If it is not avoided, the medical devices could be damaged.

NOTE! indicates special information on the operation of the device or important information.

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1.2 Pictograms/ packaging symbols

	Manufacturer
	Follow instructions for use
REF	Catalogue number
QTY	Number of products in the product packaging
SN	Serial number
	Fragile, handle with care
	Keep dry
	Temperature limit
H2O2	Suitable for hydrogen peroxide sterilization
RxOnly	Federal (USA) law restricts this device to sale by or on the order of a physician.

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
① Function buttons
② Video connector
③ Connection for leakage tester or pressure compensation cap (not visible)
④ Handpiece
⑤ Sheath
⑥ Distal tip
⑦ Deflection lever
⑧ Leakage tester 13242 XL
⑨ Pressure compensation cap 11025 E

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3 Intended use

3.1 Intended use

Flexible endoscopes are used for endoscopic imaging and image transmission, enabling endoscopic diagnostics of body cavities or hollow organs in ENT medicine. The models with working channel can also be used for therapeutic purposes. Flexible endoscopes are designed for transient use in both surgically and naturally invasive interventions.

3.2 Indications

The use of the video rhino laryngoscope is indicated if, in the opinion of the attending physician, a diagnosis or treatment of the nasal, paranasal sinus, nasopharyngeal and laryngeal structures with the flexible endoscope is indicated.

3.3 Contraindications:

The use of the flexible endoscope is contraindicated if, in the opinion of the attending physician, the surgical method is contraindicated or if the patient's general condition does not allow endoscopic surgery or anaesthesia.
The flexible endoscope must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

In addition, the following contraindications apply:

- Bleeding tendency and coagulation disorder that are not reversible or treatable

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3.4 User qualification

The KARL STORZ video rhino laryngoscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the flexible rhino/laryngoscopic technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the video rhino laryngoscopes. It is not intended as an introduction to the technique of flexible rhino-/laryngoscopy.

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4 Safety instructions

The flexible videoscope must only be used according to the medical rules and procedures of endoscopy recognized for this method.

Risks of infection

The instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties.

- ➢ *Inspect reusable instruments for visible contamination. Reprocess contaminated instruments.*
- ➢ *Reprocess reusable instruments before and after each use. When doing so, use validated procedures (see chapter "Reprocessing").*

Combination of system components

***WARNING:** Combinations of medical devices are only assured to be safe if they are listed as such in the respective instruction manuals or the intended use and interface specifications of the products used in the combination permit this.*

- ➢ *The instruction manuals and interface specifications of the medical devices used in combination must be observed precisely.*

***WARNING:** Risk of injury: Whenever the flexible videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition.*

- ➢ *Damaged flexible videoscopes or damaged accessories must not be used.*

***WARNING:** Risk of injury and damage to the products: The flexible HD video rhino laryngoscope from KARL STORZ are not intended for use with light sources. The use of non-compatible devices can lead to patient injury or device damage.*

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Hazards due to unintentional current flow

WARNING: *Risk of injury and damage to the product:*

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WARNING: *Risk of injury: The flexible videoscope must not be used during the discharge of a cardiac defibrillator.*

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WARNING: *Risk from electromagnetic incompatibility: Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC).*

➢ *Observe the stated EMC instructions during installation and operation.*

Risk of explosion

WARNING: *Risk of explosion: If flammable anesthetics are used in the immediate proximity of the flexible videoscope, there is a risk of explosion.*

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➢ *Do not use damaged videoscopes under any circumstances.*

⚠ ***WARNING:***

Risk of injury: The lens at the distal end of the flexible videoscope may become detached, e.g., due to mechanical forces.

➢ *Check that the lens is present and undamaged before and after every reprocessing.*

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NOTE: *Any damage to the flexible videoscope resulting from incorrect operation is not covered by the manufacturer's warranty.*

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6.2 Operation

The following instructions outline the correct operation of the video rhino laryngoscope and the corresponding accessories. This is not an introduction to the techniques of rhino-/laryngoscopy. Instructions on endoscopic technique can be found in the relevant medical literature.

The video rhino laryngoscope is intended for visualization and diagnostic purposes.

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6.2.1 Holding correctly

Grip the handle of the video rhino laryngoscope in your left or right hand, as shown in the picture. With this kind of control grip, the video rhino laryngoscope can be operated with one hand.

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6.2.2 Handling the sheath

ATTENTION: *Risk of damage to the sheath:*

- ➢ *Do not carry the flexible videoscope by the sheath.*
- ➢ *Do not pull, clamp or twist the sheath of the flexible videoscope. Never place heavy objects on the videoscope.*
- ➢ *The sheath of the flexible videoscope must not be kinked, clamped or coiled up too tightly.*
- ➢ *It is best to store the flexible videoscope in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39406 AS).*

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
Handling

Every time the video rhino laryngoscope is used, you must first check that its surface, and above all the surface of the sheath and distal tip, is in perfect condition and that there are no sharp edges. Video rhino laryngoscopes with sharp edges must not be used. Otherwise serious injuries may be caused to the patient or the user.

***CAUTION:** Do not squeeze the distal tip of the video rhino laryngoscope too tightly.*

Inspect the supply hose for damage (such as breaks, cracks, twisted or squashed sections...) and check the supply connector for loose connections.

Observe the lumen the whole time during use. If you lose your orientation, pull the sheath back slightly into the neutral position until vision is restored.

⚠ ***WARNING:***

Risk of injury: Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

***ATTENTION:** Lubricants that contain silicone may cause damage to the sheath.*

➢ *Do not use lubricants that contain silicone.*

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6.3 Connection to the camera contol unit

6.3.1 Compatible camera control units

The flexible videoscope is connected to the camera control unit for operation.

For operation, an IMAGE1 S™ system with an IMAGE1 S CONNECT® and an IMAGE1 X-LINK® module are required. Both camera control unit modules must be connected to each other.

6.3.2 Connecting the videoscope

1. Insert the camera cable as far as it will go into the camera connection on the TC 301 IMAGE1 X-LINK® module.
2. Turn the device on by pressing the ON button on the IMAGE1 CONNECT® module TC 200 or TC 201.

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NOTE: *For the correct settings for the camera control unit, please see the operating instructions of the devices used.*

⚠ ***WARNING:***

If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the device from the patient's body gently and carefully. If the videoscope is used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

6.3.3 Inspecting the LED function

Prepare the camera control unit (CCU) as described above. Connect the videoscope to the CCU. Aim the distal end at a surface. A light emission should be visible. The light intensity is controlled automatically by the camera control unit.

⚠ ***WARNING:***

Risk of burns and glare: The LED light emitted through the flexible videoscope may cause glare or high temperatures, and thereby burns, at the distal light outlet.

- ➢ *Never look directly into the light emitted by the flexible videoscope.*
- ➢ *Prolonged direct contact of light-emitting surfaces with tissue can lead to burns; this must therefore be avoided at all costs.*
- ➢ *Switch off the flexible videoscope before putting it down and do not place it in the vicinity of the patient (as per chapter 16 of the 3rd edition of IEC 60601-1)*

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6.5 Function buttons

The video rhino laryngoscope has three function buttons which can be assigned specific functions in the camera control unit's menu. The assignment set is retained after switching off and on again. Instructions on programming the function buttons can be found in the instruction manual of the camera control unit.

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Function buttons default function as follows

Large button in the middle:
Press briefly: Live menu
Press for longer: Set-up menu

Right button:
Press briefly: Image capture
Press for longer: Video recording

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7.2 General information

⚠ WARNING:

Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.

⚠ WARNING:

Risk of infection from non-sterile delivery. Reprocess the device before each use. Check the medical device for visible soiling before use.

⚠ WARNING:

Before reusing the device, replace any defective parts or send the flexible endoscope for repair.

***ATTENTION:** The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.*

***ATTENTION:** Never clean flexible endoscopes with metal brushes.*

***ATTENTION:** Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.*

***ATTENTION:** Flexible endoscopes must never be steam sterilized (autoclaved).*

***ATTENTION:** To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.*

***NOTE:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.*

***NOTE:** Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.*

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NOTE: *The manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’ can be downloaded or requested by visiting www.karlstorz.com.*

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Place of application	Art. number brush/ adaptor
Clean surface	27652
Ventilation connector	13242 XA (leakage test) / 11025 E (pressure compensation cap)

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Wipe-clean the exterior surfaces of the endoscope with a damp cloth.

NOTE: *To avoid damaging the sheath, do not use pushing and pulling movements.*

NOTE: *Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.*

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7.5.1 Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector securely to the supply connector vent port and lock by turning the bayonet catch clockwise 25°.

NOTE: *The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.*

Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer on the manometer. The test pressure can drop by up to 10 mmHg until the endoscope is completely filled with air.

After stabilizing the pressure on the manometer, wait at least 2 minutes to see whether the pressure remains constant or drops.

7.5.2 Performing the leakage test – in liquid

NOTE: *To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.*

Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in ‘Leakage test – dry’ and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

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Cleaning can then begin.

If a pressure loss occurs, this means that the endoscope leaks and may no longer be used, immersed in liquid or reprocessed by machine.

ATTENTION: *In this case, fill out the accompanying repairs form and send the endoscope, which has been cleaned by wiping to remove organic soiling/residue, for repair or, alternatively, contact KARL STORZ. The endoscope must be labeled as not disinfected. Line the transport case with the transport protection sheet 13990 SFN.*

7.6 Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer‘s specifications regarding material compatibility must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth.

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ATTENTION: *Never remove dirt from the device with a needle or other sharp objects, since these can cause severe damage.*

ATTENTION: *The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.*

7.7.2 Manual disinfection

The endoscope must be completely immersed in a disinfectant solution. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

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- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

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ATTENTION: *The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.*

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7.8.1 Manual preparation

The endoscope must be completely immersed in a cleaning solution. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush.

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7.8.2 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C and pressure of 100 kPa (1 bar / 14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer‘s instructions must be taken into account when using this method.

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7.9 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

NOTE: *Clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement.*

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- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surfaces (e.g., connections) must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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- Afterwards, a functional check must be carried out.

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7.11 Sterilization

⚠ ***WARNING:***

Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

⚠ ***WARNING:***

Risk of infection: Sterilize the medical device in a disassembled state.

ATTENTION: *Flexible endoscopes must never be steam sterilized (autoclaved).*

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The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’. The method must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

7.11.1 Hydrogen peroxide (H_2O_2) sterilization – ASP STERRAD®

ATTENTION: *Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.*

NOTE: *Refer to the ‘STERRAD® Sterility Guide’ to verify whether the medical device concerned can be sterilized using the different STERRAD® devices.*

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S Long Cycle
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

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7.11.2 Hydrogen peroxide (H_2O_2) sterilization – STERIS® V-PRO™

NOTE: *The ‘STERIS Low Temperature Reprocessing Device Compatibility Matrix’ can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.*

The following STERIS V-PRO® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- V-PRO® maX Flexible Cycle
- V-PRO® 60 Flexible Cycle

7.11.3 Chemical low-temperature sterilization with peracetic acid – STERIS System 1/1E®

Detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

ATTENTION: *The pressure compensation cap must not be in place on the vent port.*

7.12 Limits of reprocessing

The product’s service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

7.13 Storage

The storage method must be selected taking into account the respective applicable national requirements.

KARL STORZ only recommends dust-protected, suspended storage for decontaminated endoscopes.

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ATTENTION: *Do not expose flexible endoscopes to direct sunlight or excessive heat.*

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ATTENTION: *The distal tip of the flexible endoscope must be in the neutral position during storage.*

ATTENTION: *Do not store the flexible endoscope in the case.*

ATTENTION: *The sheath of the flexible endoscope must never be kinked or coiled up too tightly.*

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8 System description

The KARL STORZ flexible HD video rhino laryngoscope has a 300 mm long working shaft, 3.70 mm in diameter. The distal tip is 3.70 mm in diameter and houses the CCD video module. The distal tip of the flexible HD video rhino laryngoscope can be deflected 140° up or 140° down. The handle houses the deflection lever, LED light source and the video function buttons, which can be programmed to perform functions such as image capture, zoom, and provide access to the set up menu.

The KARL STORZ HD video rhino laryngoscope is connected to the recommended camera control unit IMAGE 1 S TC 200 / TC 201 /TC 301 via an integrated 2.5 meter connecting cable.

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9 Technical data

Direction of view	0°
Angle of view	100°
Deflection	up/down 140°/140°
Working length:	30 cm
Sheath outer diameter	3.7 mm
Distal outer diameter	3.7 mm
Illumination	integrated LED light source
Weight	510 g max.
Drip water protection	IPX 8
Classification	IIa

9.1 Operating/storage conditions

	Temperature	Rel. humidity
Storage	0°C … +60°C	15% … 90%
Operation	+15°C … +37°C	15% … 80%
Air pressure	700 hPa … 1080 hPa	

ATTENTION: *The specifications in relation to atmospheric pressure can only be guaranteed if the pressure equalization cap is fitted to the vent port during transportation.*

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10 Service, repair and disposal

10.1 Warranty

The warranty conditions provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local KARL STORZ subsidiary (see the section entitled ‘Subsidiaries’), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the device by unauthorized persons shall relieve KARL STORZ of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall invalidate any warranty.

10.2 Maintenance

No further repair measures are required beyond reprocessing.

However, the legal requirements as per the Medical Devices Operator Ordinance (MPBetreibV) for devices to be used in combination apply (see Section 4.4). The relevant warnings in the instruction manuals must be observed.

10.3 Repairing device parts

Only personnel authorized by the manufacturer may repair device parts. The operations described in this manual are exempt from this rule.

The authorized personnel may only repair device parts which the manufacturer designates as repairable.

Upon request, the manufacturer will provide authorized personnel with information on repairs, e.g., circuit diagrams, spare parts lists, setting instructions.

The fact of having such documents at one‘s disposal does not correspond to authorization.

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10.4 Repair program

For this medical device, individual repair is necessary. To bridge the repair period, you will receive a device on loan, which must then be returned as soon as you receive the repaired medical device.

In the case of repairs, please contact your local KARL STORZ subsidiary (see section ‘Subsidiaries’) or authorized dealer.

10.5 Infection prevention

To prevent infections, contaminated devices and instruments must not be returned.

Decontaminate instruments and devices before returning. KARL STORZ returns contaminated instruments/devices to the sender.

10.6 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

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Please ask KARL STORZ, a KARL STORZ subsidiary (see section ‘Subsidiaries’), or your authorized dealer about the collection point in your area.

Within the scope of application of this directive, KARL STORZ is responsible for the proper disposal of this device.

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11 Appendix

11.1 Troubleshooting

Problem Image unclear **Possible reason** 1. Distal objective lens fogged 2. Dirty lens **Recommended remedy** 1. Rinse lens 2. Withdraw video rhino laryngoscope and clean distal lens with cotton tip applicator and alcohol
Problem Inadequate lighting **Possible reason** Light conducting fibers damaged **Recommended remedy** Notify customer service
Problem Light does not function **Possible reason** LED faulty **Recommended remedy** Notify customer service

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12 Electromagnetic Compatibility user information

⚠ WARNING:

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

⚠ WARNING:

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the flexible videoscope 11101 HD together with the IMAGE1 S™ including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

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The flexible videoscope 11101 HD together with the IMAGE1 S™ is intended for use in an electromagnetic environment specified below. The user should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kH repetition	± 2 kV/1kV power lines ± 1 kV for input and output lines 100 kH repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	Voltage dip: Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100% for 250/300 cycles	Voltage dip: Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the flexible videoendoscope 11101 HD together with the IMAGE1 S™ requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the flexible videoendoscope 11101 HD together with the IMAGE1 S™ further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be me asured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GH	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 Vrms on 150 kHz to 80 MHz 1kHz 80% AM modulation 6 Vrms in ISM bands	3 Vrms on 150 kHz to 80 MHz 1kHz 80% AM modulation 6 Vrms in ISM bands	

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The flexible videoscope 11101 HD together with the IMAGE1 S™ is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
			Portable and mobile RF communications equipment should be used no closer to any part of the flexible videoscope 11101 HD together with the IMAGE1 S™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF disturbances acc. to IEC 61000-4-6	3 V_{rms} 150 kHz to 80 MHz	3 V_{rms}	$d = 1.2\sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey[a] should be less than the compliance level in each frequency range[b].
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz
			Interference may occur in the vicinity of equipment marked with the following symbol: 

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[a] Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the flexible videoscope 11101 HD together with the IMAGE1 S™ is used exceeds the above compliance levels, the flexible videoscope 11101 HD together with the IMAGE1 S™ should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the flexible videoscope 11101 HD together with the IMAGE 1 S™.
[b] Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 Guidelines and manufacturer's declaration – electromagnetic emissions		
The flexible videoscope 11101 HD together with the IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the flexible videoscope 11101 HD together with the IMAGE1 S™should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The flexible videoscope 11101 HD together with the IMAGE1 S™ uses RF energy only for its internal function. The customer or user of the flexible videoscope 11101 HD together with the IMAGE1 S™ should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The flexible videoscope 11101 HD together with the IMAGE1 S™ is suitable for use in all establishments other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5
Recommended separation distances between portable and mobile RF telecommunication devices and the flexible videoendoscope 11101 HD together with the IMAGE1 S™

The videoscope 11101 HD together with the IMAGE1 S™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the flexible videoscope 11101 HD together with the IMAGE1 S™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{3}\right]\sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{3}\right]\sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3}\right]\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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