BEROMED GMBH HOSPITAL PRODUCTS®

Postal Address
PF 060108
10051 Berlin
Germany

Office
Kronenstrasse 19
10117 Berlin
Germany

Phone +49-30-20 45 14 60
Fax +49-30-20 45 14 61
E-mail info@beromed.de
http www.beromed.de

Office: BEROMED GMBH HOSPITAL PRODUCTS · Kronenstrasse 19 · 10117 Berlin · Germany

Beromed GmbH Hospital Products

Instruction For Use

Version: A **Date: 20210225**

Product Name: BEROJECT III Disposable Sterile Syringe (Luer Lock) without needle
Specification: 1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50(60)ml, 60 ml
Manufacturer Name: Beromed GmbH Hospital Products, Germany
Registered Address: Kronenstrasse 19, 10117 Berlin, Germany
Production Site: Jiangsu Kanghua Medical Equipment Co., Ltd.
Production Address: Sanhekou, Changzhou, China

Product Performance, Structure and Composition: The product is 3-part, luer lock, and composed by barrel, plunger and gasket. The main raw materials are polypropylene, polyethylene and natural rubber.

Intended Use: Preparation of solutions and reagents for patients. Only for transient use. (as defined in MDD93/42/EEC Appendix IX).

Indication: Absorb or prepare solution

Contraindications:
1. The patient's allergy to the material used (e.g. natural rubber)
2. Use with pressure equipment (e.g. syringe pumps).
3. The product is sterile and non-pyrogen. For single use only.

Instruction for use
1. Open the individual package, then remove the protector.
2. Absorb solution, and vent out the gas inside, then perform the preparation of solution or reagent.

Caution:
1. Sterilized by EO gas.
2. Do not use when individual package is broken or contaminated, protector removed or there is something inside the package.
3. Sterilized, only for single use. Any re-use or re-sterilisation can be harmful to patients. Discard after use and dispose of waste in accordance with national and institutional guidelines.

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Instruction For Use

Version: B **Date: JSKH/CE-2018-01**

Product name ： BEROJECT III Disposable sterile insulin syringe
Manufacturer name: Beromed GmbH Hospital Products, Germany
Registered Address: Kronenstrasse 19, 10117 Berlin, Germany
Production Site: Jiangsu Kanghua Medical Equipment Co., Ltd.
Production Address: Sanhekou, Changzhou, China

1. Scope of application

This product is mainly suitable for the injection of insulin solution into patients.

2. Main performance and structure, specifications

1)The product is composed by hand cap , barrel, plunger, gasket , needle cap and needle

2)Specification :0.5ml (U-100,U-40,U-20)、 1ml (U-100,U-40,U-20)

3) Material specification:

Component	Material
Syringe	
Barrel:	Polypropylene
Plunger:	Polypropylene
Rubber seal (gasket)	synthetic rubber, latex-free
Hypodermic needle	
Needle tube	Stainless steel
Protective cap	Polyethylene

4)Sterilized by EO, bacteria free, non- pyrogenic.

3. Instruction for use

1) Remove protector, absorb solution, and vent out the gas inside.

2) Perform hypodermic injection.

4. Attention

1） Operators using this product must have relevant medical knowledge;

2） This product is valid for five years. Please pay attention to the expiry date at the sealing of individual package before use. Do not use when expired.

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3）Do not use when individual package is broken or contaminated, protector removed or there is something inside the package.

4）Once the single package is opened, use it immediately;

5）Discard after single use.

6）Only the corresponding concentration of insulin can be injected.

5. The storage and transportation of products

1) While transporting, the products should avoid heavy weight, direct sunlight and rain or snow soak.

2) This product should be stored in the environment where relative humidity is not more than 80%, and there is no corrosive gas, cool, dry and well-ventilation.

Labelling and explanation.

1	Single use only		The product can not be reused.
2	Reading the instruction		Read the instruction before use
3	Lot No.	LOT	Lot No. at a same production condition
4	Manufacturing date		The manufacturing date of the product
5	EO sterilization	STERILE EO	The product is sterilized, by EO gas
6	Telephone		The telephone to be contacted
7	Expiry date		The product has to be used before this date
8	Fragile		The product is fragile, easy to be broken, careful handling
9	Keep away from sunlight		The product shall avoid the hot during transport and storage
10	Keep dry		The product shall avoid the humidity during transport and storage

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11	Weight limit for pile up		The maximum weight to pile up when transport and store
12	Manufacturer		Name and address of the company
13	EC representative		The authorized representative in European countries.
14	Don't use if the packing damaged		Don't use the product if the packing is damaged.
15	Caution		Caution, please see instruction for use
16	Do not resterilize		The device should not be re-sterilized after it once has been sterilized
17	Contains or presence of latex		Contain or presence of latex
18	Latex-free		Do not contain or presence of latex
19	Non-pyrogenic		The product is non-pyrogenic

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3.1 Product name: BEROJECT III Disposable sterile Catheter Tip Syringe

3.1.1 General Description

BEROJECT III Disposable sterile Catheter Tip Syringes are manufactured according to ISO 7886-1: 2017 and GB15810-2001 standards. The main products are 3-part, using medical grade Polypropylene as raw material, strictly manufactured from materials input to products assembling at 10,000 grade cleaning workshop according to the requirement. There are up to 4 types. The products are non-bacterial, non-toxic, and non-pyrogenic after EO sterilization. The products are single-use, which can avoid cross infection. They are easy to use for patients' feeding and irrigation.

3.1.2 Scope of Application and Intended Use

50 (60) ml (3-part, central cone): Accurate measurement of the volume of irrigation solution. It is used for surgical and gynaecological irrigation as well as irrigation of body wound surfaces and hollow organs. Use within five years after sterilization.

100 ml (3-part, central cone): Accurate measurement of the volume of irrigation solution. It is used for surgical and gynaecological irrigation as well as irrigation of body wound surfaces and hollow organs. Use within five years after sterilization.

120ml (3-part, central cone): Accurate measurement of the volume of irrigation solution. It is used for surgical and gynaecological irrigation as well as irrigation of body wound surfaces and hollow organs. Use within five years after sterilization.

3.1.3 The Structure and composition

3-part catheter tip syringe: 50 (60)ml, 100ml, 120ml. All are composed of barrel, plunger, latex gasket or latex-free gasket, protector (cap).

3.1.4 Instruction for Use

1) Open the individual package.

2) Remove the protector, suction food then feed the patent. Or remove the protector, suction fluid, then do surgical, gynaecological irrigation or flush the patent's wound surfaces.

3.1.5 Attention

1) This product is valid for five years. Please pay attention to the expiry date at the sealing of individual package before use. Do not use when expired.

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2) Do not use when individual package is broken or contaminated, protector removed or there is something inside the package.
3) Discard after single use.

3.1.6 The storage and transportation

1) While transporting, the products should avoid heavy weight, direct sunlight and rain or snow soak.
2) This product should be stored in the environment where relative humidity is not more than 80%, and there is no corrosive gas, cool, dry and well-ventilation.

3.1.7 Label and explanation (see Table 1).

1	Single use only		The product can not be reused.
2	Reading the instruction		Read the instruction before use
3	Lot No.		Lot No. at a same production condition
4	Manufacturing date		The manufacturing date of the product
5	EO sterilization		The product is sterilized, by EO gas
6	Telephone		The telephone to be contacted
7	Expiry date		The product has to be used before this date
8	Fragile		The product is fragile, easy to be broken, careful handling
9	Keep away from sunlight		The product shall avoid the hot during transport and storage
10	Keep dry		The product shall avoid the humidity during transport and storage
11	Weight limit for pile up		The maximum weight to pile up when transport and store

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12	Manufacturer		Name and address of the company
13	EC representative	EC REP	The authorized representative in European countries.
14	Don't use if the packing damaged		Don't use the product if the packing is damaged.

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Instruction for Use

Version: C **Preparation date:** 20180909

Product name: BEROJECT III single-use sterile syringe (with needle)

Product specifications: 1ml, 2ml, 2.5 ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50 (60) ml, 100ml

Manufacturer name: Beromed GmbH Hospital Products, Germany

Registered Address: Kronenstrasse 19, 10117 Berlin, Germany

Production Site: Jiangsu Kanghua Medical Equipment Co., Ltd.

Production Address: Sanhekou, Changzhou, China

Product performance structure: this product is composed of a barrel, a plunger and a gasket (rubber seal) and an injection needle (it is optional accessories).

The main materials are polypropylene, polyethylene, natural rubber, and stainless steel (variants with needle).

Expected use: Subcutaneous injections, skin tests, blood draws, and preparation of medicinal solutions and reagents. Syringes can only be used for short periods (as defined in MDD 93/42/EEC Annex IX).

This product is sterile and non-pyrogenic. For single use only.

Indications.

1. Subcutaneous injection, skin test and blood draw.
2. Extracting or disposing liquid medicine.
3. Total parenteral infusion.
4. Infusion of drugs like antibiotics.

Contraindications.

1. The patient is allergic to the materials used (e.g. natural rubber).

Instructions.

1. Remove the sheath after removing the single package.
2. Aspirate the medicinal solution, exhaust the internal gas, and perform operations such as injection and skin

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test.

Precautions.

1. This product is sterilized by ethylene oxide.
2. Please check the integrity of the single package before use. If the single package leaks, and/ or it is damaged, and/or the needle protective cover falls off, it is strictly prohibited to use the product.
3. This product is sterile and is for single use only. Any re-use or re-sterilization will cause harm to the patient. Destroy after use, and use the waste in accordance with national and institutional guidelines.
4. The operation should be implemented and monitored by professional nursing staff.
5. This product may contain natural rubber or no natural rubber. It is forbidden to those who are allergic to natural rubber. It is recommended to use non-natural rubber equivalent products.
6. This product is valid for five years after being sterilized by EO. It is forbidden to use it after it expires.
7. Our products have screw lock joints, which should meet the requirements of ISO80369-7 standard.
8. Our company's injection needle should meet the ISO7864: 2016 standard.

Transportation and storage.

1. The product should be protected from heavy pressure, direct sunlight, rain and snow during transportation.
2. This product should be stored in a relative humidity not exceeding 80%, without corrosive gases, in a cool, dry, well-ventilated and clean environment.

Label description.

1	Single use only		The product can not be reused.
2	Reading the instruction		Read the instruction before use
3	Lot No.		Lot No. at a same production condition
4	Manufacturing date		The manufacturing date of the product

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5	EO sterilization	STERILE EO	The product is sterilized, by EO gas
6	Telephone		The telephone to be contacted
7	Expiry date		The product has to be used before this date
8	Fragile		The product is fragile, easy to be broken, careful handling
9	Keep away from sunlight		The product shall avoid the hot during transport and storage
10	Keep dry		The product shall be avoid the humidity during transport and storage
11	Weight limit for pile up		The maximum weight to pile up when transport and store
12	Manufacturer		Name and address of the company
13	EC representative	EC REP	The authorized representative in European countries.
14	Don't use if the packing damaged		Don't use the product if the packing is damaged.
15	Caution		Caution, pls see the operation instruction
16	Do not resterilize	STERILIZE	Non-Repeatable Sterilization
17	Contains or presence of natural rubber latex	LATEX	Including the latex or natural latex
18	Latex-free	LATEX	Exclusive latex or natural latex
19	Non-pyrogenic		Non Pyrogen

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4. Operation should be carried out and supervised by professional nursing staff.
5. The product may contain natural rubber or not. It is prohibited for those who are allergic to natural rubber, and recommended to use non-natural rubber products.
6. This product is valid for five years after sterilization. Please pay attention to the expiry date at the sealing of individual package before use. Do not use when expired.
7. Luer lock cone should conform to ISO80369-7 standard.

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Labelling and explanation.

1	Single use only		The product can not be reused.
2	Reading the instruction		Read the instruction before use
3	Lot No.		Lot No. at a same production condition
4	Manufacturing date		The manufacturing date of the product
5	EO sterilization		The product is sterilized, by EO gas
6	Telephone		The telephone to be contacted
7	Expiry date		The product has to be used before this date
8	Fragile		The product is fragile, easy to be broken, careful handling
9	Keep away from sunlight		The product shall avoid the hot during transport and storage

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10	Keep dry		The product shall avoid the humidity during transport and storage
11	Weight limit for pile up		The maximum weight to pile up when transport and store
12	Manufacturer		Name and address of the company
13	EC representative	EC REP	The authorized representative in European countries.
14	Don't use if the packing damaged		Don't use the product if the packing is damaged.
15	Caution		Caution, please see instruction for use
16	Do not resterilize		The device should not be re-sterilized after it once has been sterilized.
17	Contains or presence of latex		Contain or presence of latex
18	Latex-free		Do not contain or presence of latex
19	Non-pyrogenic		The product is non-pyrogenic

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