BEROMED GMBH HOSPITAL PRODUCTS®

Postal Address
PF 060108
10051 Berlin
Germany

Office
Kronenstrasse 19
10117 Berlin
Germany

Phone +49-30-20 45 14 60
Fax +49-30-20 45 14 61
E-mail info@beromed.de
http www.beromed.de

Office: BEROMED GMBH HOSPITAL PRODUCTS · Kronenstrasse 19 · 10117 Berlin · Germany

Instruction for use

Product name: BEROFINE single-use sterile hypodermic (injection) needle
Manufacturer name: Beromed GmbH Hospital Products, Germany
Registered Address: Kronenstrasse 19, 10117 Berlin, Germany
Production Site: Jiangsu Kanghua Medical Equipment Co., Ltd.
Production Address: Sanhekou, Changzhou, China

Product performance structure and composition: This product is composed of hub, needle tube and protective cap.
The main materials are:

Component	Material
1 Protective cap	Polypropylene
2 Needle tube	Stainless steel SUS 304 (chromium nickel steel)
3 Adhesive	Epoxy adhesive
4 Needle hub	Polypropylene
Lubricant	Silicone oil (dimethyl silicone oil; in compliance with monograph 3.1.8 of the European Pharmacopoeia)

Expected use: Hypodermic needles are used in connection with disposable syringes or disposable infusion / transfusion sets for medicine injection / infusion or blood collection.

Indications.
The needle tube is made of stainless steel for the puncture of skin, muscle tissue and veins to allow injection of medicine, transfusion or withdrawal of blood.
One end of the hypodermic needle is a luer lock connector suitable for connection to syringes or to infusion / transfusion sets.
This product is sterile and non-pyrogenic. For single use only.

Instructions.
Open the unit packing, take out the needle cap, connect with syringe or infusion set, absorb the solution or make the hypodermic injection.

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Caution for use.

- the product is sterilized by EO, non bacteria, non toxic, non pyrogenic.
- the users have to understand some medical knowledge.
- the shelf life of the product is five years after EO sterilization, check the expiry date before use.
- don't use if the unit packing is damaged, contaminated, or some foreign matters inside the unit packing.
- Promptly use once after the unit packing was opened.
- the product is designed for single use only, discard after use.

Transportation and storage.

- Heavy pressure, direct sunlight, raining and snowing soak shall be avoided during the transportation.
- the product has to be stored at a place with relative humidity not more than 80%, no corrosive gas, cool and dry, clean and good ventilated.

8. Labeling and explanation.

1	Single use only		The produt can not be reused.
2	Reading the instruction		Read the instruction before use
3	Lot No.	LOT	Lot No. at a same production condition
4	Manufacturing date		The manufacturing date of the product
5	EO sterilization	STERILE EO	The product is sterilized, by EO gas
6	Weight limit for pile up		The maximum weight to pile up when transport and store
7	Manufacturer		Name and address of the company
8	Telephone		The telephone to be contacted
9	Expiry date		The product have to be used before this date
10	Fragile		The product is fragile, easy to be broken, careful handling
11	Keep away from sunlight		The product shall avoid the hot during transport and storage
12	Keep dry		The product shall avoid the humidity during transport and storage
13	EC representative	EC REP	The authorized representative in European countries.
14	Don't use if the packing damaged		Don't use the product if the packing is damaged.

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