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Anhui Hongyu Wuzhou Medical Manufacturer Co.,Ltd

Add:No.2 Guanyin Road ,Economic Development Zone,Taihu County 246400,Anqing,Anhui,PEOPLE'S REPUBLIC OF CHINA

INSTRUCTION FOR USE OF MEDICAL DEVICE

The disposable hypodermic needle is for single use, Sterile, Non Toxic Non pyrogenic.
And is constructed by hub , cannula ,cap.

Direction of use:
1.Open the dialysis paper, put the hub into the syringe cone head, rotate and close the join point of hub. .
2.Remove the cap pumped drug liquid, the needle upwards, slowly push the plunger to exclude air
3.Hypodermic injection

Caution:
1.Sterile unless package has been opened or damaged
2.Do not use if tip protectors are not in place
3.Avoid freezing and excessive heat.
4.The product should be used immediately after envelope opened,discard after use.
6 For single use only. Do not resterilize or re-use.

Shelf life and storage temperature:
Production from the date of production is valid for three years , must be used within the period of validity. The batch number and expiry date was printing on the seal of individual packing.
This product should be stored in the ralatibe humidity of not more than 80%, no corrosive gas , cool , dry and clean indoor.

ANHUI HONGYU WUZHOU MEDICAL MANUFACTURER CO.,LTD.

Anhui Hongyu Wuzhou Medical Manufacturer Co.,Ltd.

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