BECKMAN COULTER

CE

COULTER DxH Diff Pack

REF 628020

TABLE OF CONTENTS

English . . . 2
Français (FR) . . . 4
Deutsch (DE) . . . 6
Italiano (IT) . . . 8
Español (ES) . . . 10
Português Portugal (PT-PT) . . . 12
Dansk (DA) . . . 14
Svenska (SV) . . . 16
Norsk (NO) . . . 18
Ελληνικά (EL) . . . 20
日本語 (JP) . . . 22
Lietuviškai (LT) . . . 24
Magyar (HU) . . . 26
Polski (PL) . . . 28
Čeština (CZ) . . . 30
Slovenčina (SK) . . . 32
한국어 (KO) . . . 34
Türkçe (TR) . . . 36
Русский (RU) . . . 38
Hrvatski (HR) . . . 40
Български (BG) . . . 42
中文 (ZH-TW) . . . 44
Română (RO) . . . 46
Србија (SR) . . . 48
Українська (UK) . . . 50
Português Brasil (PT-BR) . . . 52
Nederlands . . . 54
Tiếng Việt (VI) . . . 56

© 2018 Beckman Coulter, Inc. All rights reserved.

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INTENDED USE
For use on UniCel DxH Coulter Cellular Analysis Systems to lyse erythrocytes and preserve leukocytes. Use only with COULTER DxH Diluent. Refer to your instrument product manuals and/or online help, as applicable.

SUMMARY
An erythrolytic reagent and a leukocyte preservative produce a sample for performing a five-part WBC differential analysis using VCSn technology.

PRINCIPLE
WBC differential analysis begins by diluting a whole-blood sample with a lytic solution of formic acid that rapidly lyses the erythrocytes and reduces the cellular debris to an insignificant level without altering the leukocytes. The leukocyte preservative, consisting of an aqueous solution containing soluble salts, is added to stop and/or neutralize the lytic activity and restore the ionic balance to the blood, thus preserving the leukocytes in their near-native state for differentiation into five subpopulations (neutrophils, lymphocytes, monocytes, eosinophils, and basophils) using impedance, radio frequency, and VCSn technology.

REACTIVE INGREDIENTS
Coulter DxH Diff Pack contains the following:

COULTER Erythrolyse II Reagent
Wetting Agent 0.3-1.5g/L

Formic Acid 1.2mL/L

COULTER StabiLyse Reagent
Sodium Carbonate 6.0g/L

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GHS HAZARD CLASSIFICATION		
Erythrolyse II Reagent	WARNING	Causes mild skin irritation.
		If skin irritation occurs: Get medical advice/attention.
		May produce an allergic reaction.
		Formic Acid <1%
		Diazolidinyl Urea <0.1%
		Reaction Mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) <0.0015%

STORAGE, STABILITY, AND DISPOSAL
- Store COULTER DxH Diff Pack at 2-25°C.
- Do not use product past expiration date.
- Use product at temperatures stated in the instrument product manuals and/or online help.
- Opened containers are stable for 60 days when stored as recommended.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.

PREPARATION
Reagents in the COULTER DxH Diff Pack are ready to use. Replace the reagent kit as directed in your instrument product manuals and/or online help.

To avoid cross contamination when replacing the DxH Diff Pack, transfer Erythrolyse II and StabiLyse pickup tubes one at a time.

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COULTER DxH Retic Pack

REF 628021

For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.

INTENDED USE
For use on UniCel DxH Coulter Cellular Analysis Systems for clearing red cells and staining reticulocytes. Use only with DxH Diluent.

SUMMARY
Reticulocytes are immature erythrocytes (red blood cells) that retain a small amount of RNA as the nucleus is extruded during the maturation process. To detect the presence of this RNA, the red blood cells are stained while they are still viable in a process called supravital staining. Once a whole-blood sample is prepared, reticulocytes are differentiated from mature red blood cells and other cell populations using VCSn and reticulocyte counting technologies.

PRINCIPLE
New Methylene Blue (Reagent A), a supravital dye, is incubated with a whole-blood sample. The dye precipitates the basophilic RNA network found in reticulocytes.[1] A clearing reagent (Reagent B) removes hemoglobin and unbound stain from the erythrocytes (red blood cells or RBCs) while keeping the cell and its membranes intact and without removing the precipitated dye - RNA complex. This process results in clear spherical mature RBCs and darkly stained reticulocytes that are ready for analysis. The total RBC population is analyzed on a hematology analyzer with reticulocyte counting technology using impedance, radio frequency, and VCSn technology then reticulocytes are differentiated from mature red blood cells and other cell populations.

REACTIVE INGREDIENTS
The COULTER DxH Retic Pack contains the following:

REAGENT A - Retic Stain
New Methylene Blue in Buffered Solution…….0.06% (w/v)

REAGENT B - Retic Clear
Sulfuric Acid with Stabilizers……0.15% (w/v)

WARNING AND PRECAUTIONS
- Do not inhale and/or ingest.
- Avoid eye and skin contact. In case of eye or skin contact, flush affected area with copious amounts of water for at least 15 minutes.
- DO NOT REUSE CONTAINERS.

GHS HAZARD CLASSIFICATION		
REAGENT B	DANGER 	Causes severe skin burns and eye damage.
		Wear protective gloves, protective clothing and eye/face protection.
		IF ON SKIN (or hair): Rinse skin with water.
		IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
		Immediately call a POISON CENTER or doctor/physician.
		Sulfuric Acid 0.1-0.5%

STORAGE, STABILITY, AND DISPOSAL
- Store COULTER DxH Retic Pack at 2-30°C.
- Do not use product past expiration date.
- Use product at temperatures stated in the instrument product manuals and/or online help.
- Opened containers are stable for 60 days when stored as recommended.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.

PREPARATION
Reagents in the COULTER DxH Retic Pack are ready to use. Replace the reagent kit as directed in your instrument product manuals and/or online help.

To avoid cross contamination when replacing the DxH Retic Pack, transfer pickup tubes one at a time beginning with Reagent A (Retic Stain) and then Reagent B (Retic Clear).

IMPORTANT: If product has been partially or completely frozen, allow product to warm to room temperature. Mix product by gentle inversion prior to placement on the instrument. Install and prime the reagent kit as directed in your instrument product manuals and/or online help. Verify background counts are acceptable before analyzing patient samples.

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COULTER DxH Retic Pack

REF 628021

(Contains 1 x 380 mL Reagent A and 1 x 1,900 mL Reagent B)

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

ADDITIONAL INFORMATION

Glossary of Symbols is available at techdocs.beckmancoulter.com (PN C07163)

May be covered by one or more pat. - see www.beckmancoulter.com/patents

SDS Safety Data Sheet is available at techdocs.beckmancoulter.com

Revision History

Revision AH, 08/2017

Changes were made to:

- Added new Bulgarian language

Revision AJ, 01/2018

Changes were made to:

- Added new Serbian language
- Added new Vietnamese language

Revision AK, 09/2018

Changes were made to:

- Added new Patent Information Statement
- Added new Dutch, Traditional Chinese, and Ukrainian languages

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REFERENCES

1. Brecher G. 1949. New methylene blue as a reticulocyte stain. *Am J Clin Path* 19:895.

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INTENDED USE
A cleaning agent for use on UniCel DxH Coulter Cellular Analysis Systems for components that come in contact with blood samples.

SUMMARY
This azide-free, formaldehyde-free, biodegradable cleaner contains a proteolytic enzyme that aids in the removal of protein buildup.

PRINCIPLE
Proteins that flow through apertures may accumulate, and narrow or block aperture openings. The cleaning reagent dissolves proteins at apertures as well as at aspiration and dispensing components of the systems.[1]

REACTIVE INGREDIENTS
A solution of a proteolytic enzyme

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GHS HAZARD CLASSIFICATION		
COULTER DxH Cleaner	DANGER 	May cause an allergic skin reaction.
		May cause allergy or asthma symptoms or breathing difficulties if inhaled.
		Avoid breathing vapours.
		Wear protective gloves, protective clothing and eye/face protection.
		In case of inadequate ventilation, wear respiratory protection.
		IF INHALED: Remove person to fresh air and keep at rest in a position comfortable for breathing.
		If skin irritation or rash occurs: Get medical advice/attention.
		If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.
		Subtilisin 0.1-0.5%
		Reaction Mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) <0.05%

STORAGE, STABILITY, AND DISPOSAL
- Store COULTER DxH Cleaner at 2-25°C.
- Do not use product past expiration date.
- Use product at temperatures stated in the instrument online help and/or product manuals.
- Opened containers are stable for 3 months when stored as recommended.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.

PREPARATION
COULTER DxH Cleaner is ready to use. Replace the reagent container as directed in your instrument online help and/or product manuals.

IMPORTANT: If product has been partially or completely frozen, allow product to warm to room temperature. Mix product by gentle inversion prior to placement on the instrument. Install and prime the reagent as directed in your instrument online help and/or product manuals.

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Revision AL, 01/2018

Changes were made to:

- INTENDED USE
- Added new Serbian language
- Added new Vietnamese language

Revision AM, 09/2018

Changes were made to:

- Removed European Hazard Classification
- Added new Patent Information Statement
- Added new Dutch, Traditional Chinese, and Ukrainian languages

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REFERENCES

1. Lin Y, Means GE, and Feeney RE. 1969. Action of proteolytic enzymes on N,N-dimethyl proteins. Basis for a microassay for proteolytic enzymes. *J Biol Chem,* 244(3):789-793.

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INTENDED USE

For use as an isotonic buffered diluent in conjunction with a cyanide-free lytic agent for counting and sizing blood cells on all UniCel DxH Coulter Cellular Analysis Systems. Use only with COULTER DxH Cell Lyse. Refer to your instrument product manuals and/or online help, as applicable. For use as a rinsing agent on the UniCel DxH Slidemaker Stainers.

SUMMARY

Blood cell analysis comprises diluting a whole-blood sample with a solution that functions as a diluent. The diluent provides the ability to analyze portions of the diluted blood sample for different blood cell types and the measurement of hemoglobin. COULTER DxH Diluent also acts as a rinsing agent on all UniCel DxH Coulter Cellular Analysis Systems.

PRINCIPLE

The chemical composition of COULTER DxH Diluent provides an isotonic solution with buffering capacity for cell component stabilization and measurement of hemoglobin. Biocides are included for product preservation against microbial growth. Alkaline salts provide a suitable osmolality so as not to adversely affect cell volume (iso-osmotic).

When used on the Unicel DxH Coulter Cellular Analysis Systems in conjunction with COULTER DxH Cell Lyse, COULTER DxH Diff Pack, and COULTER DxH Retic Pack, the diluent is used for erythrocyte, hemoglobin, platelet, NRBC and reticulocyte measurements, as well as, for the measurement and differentiation of leukocytes into five subpopulations – neutrophils, lymphocytes, monocytes, eosinophils, and basophils.

REACTIVE INGREDIENTS

Sodium Sulfate13.73g/L

Sodium Chloride1.04g/L

Tetracaine HCL0.02g/L

Imidazole2.85g/L

WARNING AND PRECAUTIONS

- Do not inhale and/or ingest.
- Avoid eye and skin contact. In case of eye or skin contact, flush affected area with copious amounts of water for at least 15 minutes.
- DO NOT REUSE CONTAINER.

GHS HAZARD CLASSIFICATION		
COULTER DxH Diluent	DANGER 	Causes mild skin irritation.
		Obtain special instructions before use.
		Wear protective gloves, protective clothing and eye/face protection.
		IF exposed or concerned: Get medical advice/attention.
		May produce an allergic reaction.
		Imidazole <1%
		Reaction Mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) <0.05%
		May damage fertility or the unborn child.

STORAGE, STABILITY, AND DISPOSAL

- Store Coulter DxH Diluent at 2-40°C.
- Do not use product past expiration date.
- Use product at temperatures stated in the instrument product manuals and/or online help.
- Opened containers are stable for 60 days when stored as recommended.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.

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PREPARATION
Coulter DxH Diluent is ready to use. Replace the reagent as directed in your instrument product manuals and/or online help.

IMPORTANT: If product has been partially or completely frozen, allow product to warm to room temperature. Invert the container 16 times to ensure complete mixing prior to placement on the instrument. Install and prime the diluent as directed in your instrument product manuals and/or online help. Verify background counts are acceptable before analyzing patient samples.

PRODUCT AVAILABILITY
COULTER DxH Diluent - 1 x 10 L

REF 628017

TRADEMARKS
Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

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ADDITIONAL INFORMATION
Glossary of Symbols is available at techdocs.beckmancoulter.com (PN C07163)

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Revision History
Revision AM, 06/2017

Changes were made to:

- INTENDED USE
- Added new Bulgarian language

Revision AN, 01/2018

Changes were made to:

- SUMMARY
- PRINCIPLE
- Added new Serbian language
- Added new Vietnamese language

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.