SWARTZ™ BRAIDED TRANSSEPTAL GUIDING INTRODUCER

English - Instructions for use 1
Deutsch - Gebrauchsanleitung 4
Français - Mode d'emploi 7
Italiano - Istruzioni per l'uso. 10
Español - Instrucciones de uso 13
中文 – 使用说明 16
Ελληνικά - Οδηγίες χρήσης 19
Русский - Инструкция по применению. 22
Magyar - Használati utasítás 25
Português - Instruções de Utilização 28
Nederlands - Gebruiksaanwijzing 31
Suomi - Käyttöohjeet 34
Svenska - Bruksanvisning 37
Dansk - Brugsanvisning. 40
Norsk - Bruksanvisning 43
Polski - Instrukcja użytkowania 46
Česky - Návod k použití 49
Slovensky - Návod na použitie. 52
Eesti - Kasutusjuhend 55
Lietuvių - Naudojimo instrukcijos 58
Latviski - Lietošanas pamācība 61
Türkçe - Kullanım talimatları 64
한국어 – 사용안내 67
Български - Инструкции за употреба. 70
Română - Instrucţiuni de utilizare 73
Srpski – uputstvo za upotrebu 76
العربية - تعليمات الاستخدام 79
SYMBOL 82

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Abbott

LIMITED WARRANTY AND DISCLAIMER

Abbott Medical hereby warrants that if any Abbott product fails to perform within normal tolerances for a patient due to a defect in materials or workmanship, Abbott will provide, at no charge, a replacement Abbott product for the patient's use. This limited warranty applies only if each of the following conditions are met.

1. The product was packaged and labeled by Abbott.
2. The failed product must be returned to Abbott and becomes the property of Abbott.
3. The product has not been mishandled, reprocessed or altered in any way.
4. The product was used before the "USE BY" date marked on the packaging of the product.

No representation or warranty is made that a Abbott product will not fail. Abbott disclaims responsibility for any medical complications, including death, resulting from the use of its products. Except as expressly provided by this limited warranty, Abbott IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE, OR MALFUNCTION OF ITS PRODUCTS, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE. Some states do not allow the exclusion or limitation of incidental or consequential damages however, so the above limitation or exclusion may not apply to you.

Except as expressly provided by the limited warranty, Abbott MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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⚠ WARNING: This product can expose you to chemicals including ethylene oxide, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

™ Indicates a trademark of the Abbott group of companies.
‡ Indicates a third party trademark, which is property of its respective owner.
Pat. http://www.abbott.com/patents
© 2020 Abbott. All Rights Reserved.

SWARTZ™ BRAIDED TRANSSEPTAL GUIDING INTRODUCER Instructions for Use

See product catalog for reorder numbers.

Read instructions for use prior to use of this device.

See individual sterile package label for contents.

Single-use disposable medical device. Contents are sterile if package is unopened and undamaged. Do not resterilize.

The instructions for use is recyclable. Dispose of used product and packaging following standard solid biohazard waste procedures.

The French size specified represents the inner diameter of the introducer sheath.

DEVICE DESCRIPTION

The Abbott Medical Transseptal Catheter Introducer Set consists of a radiopaque sheath and dilator, each with a specially curved distal portion to accommodate positioning against the atrial septum and to accommodate an .032" guidewire and a BRK™ type curved puncture needle. The introducer sheath is fitted with a valve to provide hemostasis during catheter introduction and/or exchange. Located on the valve housing, a sideport with three-way stopcock is provided for aspiration, fluid infusion, blood sampling and pressure monitoring. The dilator is tapered at the distal tip with an internal lumen designed to accept ancillary devices (e.g., needles or guidewires) that have a maximum diameter of 0.032". The dilator inner lumen is also designed with a special geometry at its distal end to limit the exposure of the BRK™ needle.

INDICATIONS

The Abbott Medical Transseptal Catheter Introducer Set is used for introducing various cardiovascular catheters into the left side of the heart through the interatrial septum.

CONTRAINDICATIONS

- Previous interatrial septal patch or prosthetic ASD closure device.
- Any previous thromboembolic event.
- Known or suspected left atrial myxoma.
- Known or suspected myocardial infarction within the last two weeks.
- Unstable angina.
- Recent cerebral vascular accident (CVA).
- Patients who do not tolerate anticoagulation therapy.
- Patients with an active infection.

WARNINGS

- Do not alter this device in any way.
- Federal law (U.S.) restricts this device to sale by or on the order of a physician.
- Do not reuse this device. Thorough cleaning of biological and foreign material is not possible. Adverse patient reactions may result from reuse of this device.
- Store in a cool, dark, dry place.

PRECAUTIONS

Special Patient Populations

Prior to the procedure, the patient must be hemodynamically stable. Certain conditions may require special consideration when using this product. These may include, but are not limited to:

- Rotated heart
- Enlarged aortic root
- Marked right atrial enlargement
- Scoliosis/kyphosis
- Abnormal left atrial geometry
- Congenital malformations
- Vascular malformations
- Inability to access the right atrium through the IVC

Procedural Considerations

- Carefully reading the Instructions before use of this device will help to reduce the potential dangers associated with the transseptal technique such as air emboli or perforation of the aorta or left atrium. Only those physicians who are specially trained in transseptal procedures should use this device. Fluoroscopy should be used to confirm positioning throughout the procedure. Transseptal procedures should be performed only in facilities appropriately equipped and staffed to perform such procedures. Lab capabilities should include, but are not limited to:
 - Intracardiac pressure monitoring capabilities;
 - Systemic pressure monitoring;
 - Contrast media injection, and management of untoward reactions to contrast media;
 - Pericardiocentesis;
 - Surgical backup;
 - Anticoagulation therapy and monitoring
- Maintain monitoring of vital signs throughout the procedure.
- Inspect all components before use.
- Use only a BRK™ type curved needle with stylet.
- During insertion, check for excessive resistance as the tip of the needle advances through the curvature of the sheath/dilator assembly. If there is any resistance to needle advancement, retract the needle, and assess components.
- During insertion, use caution not to create excessive bends in this device. This may inhibit advancement of the needle and may result in inadvertent needle puncture of the dilator/sheath assembly.
- During insertion, always use the stylet to facilitate needle passage through the dilator/sheath assembly. Failure to use the stylet may inhibit advancement of the needle and may result in inadvertent needle puncture of the dilator/sheath assembly or skiving of material from the inner surface of the dilator.
- To minimize the potential for creating a vacuum in the sheath, remove components and make catheter exchanges slowly.
- Once the sheath is inserted into the vasculature and the dilator is removed, aspirate until steady blood return is achieved prior to flushing or infusion.
- All fluid infusion should be through the sideport.
- Thrombus may accumulate in or on the sheath tip during the procedure. Aspirate when removing dilator or catheter.
- In order to minimize embolic risk, either provide a continuous infusion of heparinized solution or periodically aspirate and flush through the sideport while the sheath is positioned in the vasculature.
- Do not remove dilator or catheter rapidly. Damage to the valve may occur, potentially compromising hemostasis.
- If resistance is met when advancing or withdrawing guidewire or introducer, determine the cause and [illegible] this procedure.
- Indwelling intracardiac introducer sheaths should always be supported with a catheter or an obturator.
- Do not manipulate the sheath in the heart without a device extending from its distal tip.
- This device is intended for one time use only. Do not reprocess or reuse. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another.

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POTENTIAL COMPLICATIONS

Complications that may occur during use of this device include, but are not limited to:

- Air embolism
- Cardiac perforation
- Cardiac tamponade
- Conduction system disturbances such as SA node, AV node, or His-Purkinje system block
- Hematoma or excessive bleeding at the vascular access site
- Reaction to contrast medium
- Stroke
- Thromboembolism
- Valvular damage
- Cardiac arrhythmias

Suggested Transseptal Procedure

There are 8 major steps in the transseptal technique:

1. Prepare and assemble equipment.
2. Advance sheath/dilator assembly into superior vena cava.
3. Position BRK™ needle inside assembly.
4. Drag assembly and engage fossa ovalis.
5. Puncture the fossa ovalis with the BRK™ needle.
6. Advance sheath/dilator assembly over fixed needle.
7. Advance sheath over fixed dilator and needle into left atrium.
8. Remove the dilator and needle from the sheath.

1. PREPARE AND ASSEMBLE EQUIPMENT.

- Remove contents from packaging.

 NOTE: Do not assemble the sheath and dilator.
- Flush the dilator and sheath with heparinized saline.
- Close the sheath stopcock.
- Remove the stylet from the needle and flush the needle with heparinized saline.
- Insert the stylet into the needle and lock it onto the hub.
- Insert the dilator fully into the transseptal sheath.

2. ADVANCE SHEATH/DILATOR ASSEMBLY INTO SUPERIOR VENA CAVA.

- Obtain femoral venous access (right femoral preferred). OPT: A larger (≥ 2.5 French sizes greater than the transseptal sheath) standard length sheath may be used to obtain and maintain venous access for device exchange and hemostasis.
- Introduce a 0.032", 150 -180 cm, 3 mm "J" tip guidewire into the superior vena cava (SVC).

 NOTE: 0.032" is the maximum guidewire diameter that can be used with the dilator.
- Insert the sheath/dilator into the vein over the guidewire until the sheath tip is in the SVC. Orient the dilator tip medially.

3. POSITION THE BRK™ NEEDLE AND STYLET ASSEMBLY INSIDE THE SHEATH/DILATOR ASSEMBLY.

- Remove the guidewire from the dilator.
- Fully aspirate and then flush the dilator. Flush the dilator with heparinized saline.
- Separate the sheath and dilator by withdrawing the dilator approximately 5 cm (See Fig. 1). This will facilitate passage of the needle curve through the rigid hubs of the dilator and sheath.



Fig. 1


- Lock the stylet onto the hub.
- Insert the needle/stylet into the dilator. Advance the needle/stylet under fluoroscopy until it is approximately 2 cm from the dilator tip inside the sheath (See Fig. 2).

1. needle
2. stylet
3. dilator tip
4. sheath

Fig. 2

- Reconnect the sheath and dilator by advancing the dilator and needle/stylet into the sheath simultaneously under fluoroscopy. Ensure the needle's position within the dilator does not change as both components are advanced together. Make certain to maintain the sheath tip position in the SVC.
- Under fluoroscopy, advance the needle/stylet until the stylet is positioned just proximal of the dilator tip (See Fig. 3).


1. sheath tip
2. dilator
3. needle
4. stylet

Fig. 3

- Remove the stylet. Do not discard.
- Aspirate the needle until blood return is observed. Discard the syringe.

 NOTE: The use of a slip-tip (non Luer-Lok‡) syringe may prevent aspirating air.
- Flush the needle with heparinized saline. Close the stopcock.
- Advance the needle under fluoroscopy until the needle tip is just proximal of the dilator tip (See Fig. 4).

1. dilator tip
2. needle tip

Fig. 4


- OPT: Attach a 3-way rotating stopcock to the needle hub.
- OPT: Attach a syringe with radiopaque contrast media to the stopcock. Aspirate the BRK™ Transseptal Needle until blood return is observed. Then load the needle with the contrast media under fluoroscopic guidance.
- OPT: Connect a pressure monitoring line to the stopcock.
- OPT: Use a standard 3-port manifold setup to connect contrast, pressure, and flush lines.

4. ENGAGE THE FOSSA OVALIS.

- Visualize and identify anatomic landmarks.
- Set the fluoro unit to an appropriate angle, parallel to the plane of the mitral valve and orthogonal to the plane of the septum. This will typically be LAO, approximately 30 to 40 degrees.
- OPT: During electrophysiology procedures, the coronary sinus and His bundle catheter positions can serve as useful anatomic landmarks. In the appropriate LAO view, the coronary sinus catheter will be seen ight anterior oblique (RAO) view the His bundle catheter will appear in profile. The fossa ovalis is located at or slightly below the level of His bundle catheter and superior and posterior to the coronary sinus ostium.

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- OPT: Placing a pigtail angiographic/hemodynamic monitoring catheter in the non-coronary cusp of the aortic valve can serve as a useful anatomic landmark.
- OPT: Observe the pressure waveform being recorded through the BRK™ Transseptal Needle.
- Adjust the pointer flange so that the needle is perpendicular to the fossa ovalis (typically between 3:00 and 5:00 o'clock, as viewed from the foot end of the patient). (See Fig 5.)

1. needle pointer flange

Fig. 5

- Confirm that the needle tip is inside the dilator by fluoroscopy.
- After confirming that the tip of the needle is within the dilator, drag the entire sheath/dilator/needle assembly slowly. Prevent any movement of the assembly parts relative to each other. It is critical to maintain the previous orientation of the pointer flange while dragging assembly.
- In the LAO view (orthogonal to the interatrial septum) observe the tip of the dilator during the drag for abrupt medial (or rightward) movement, indicating the tip has engaged the fossa ovalis. (See Figs. 6a, 6b, & 6c.)


 NOTE: If the fossa ovalis is probe patent, the dilator tip will now move into the left atrium with ease.

a. Starting position in SVC

b. Initial medial movement in RA

c. Abrupt medial movement onto fossa ovalis

Fig. 6

- OPT: If pressure is being monitored through the needle, note that the pressure through the needle will not be accurate at this point, since the tip is against the fossa ovalis.

5. PUNCTURE THE FOSSA OVALIS WITH THE BRK™ NEEDLE.

- Confirm the correct location of the sheath/dilator/needle assembly on the fossa ovalis before advancing the needle.
- Once the correct location is confirmed, extend the needle to full engagement within the sheath/dilator assembly and advance across the interatrial septum.
- OPT: Under pressure monitoring, entry into the left atrium is confirmed when the pressure tracing shows a left atrial pressure waveform.
- OPT: Left atrial access can be confirmed via fluoroscopy with contrast injections.
- If there is any resistance to needle advancement, retract the needle and re-evaluate the anatomic landmarks.

CAUTION: If pericardial or aortic entry occurs, do not advance the dilator over the needle. If the needle has penetrated the pericardium or aorta, it must be withdrawn. Monitor vital signs closely.

6. ADVANCE THE SHEATH/DILATOR ASSEMBLY.

- While maintaining a fixed needle position within the left atrium, advance the sheath/dilator assembly fully over the needle into the left atrial cavity.

7. ADVANCE THE SHEATH OVER THE FIXED DILATOR AND NEEDLE INTO THE LEFT ATRIUM.

- Maintain the position of the needle and dilator across the septum.
- While maintaining the dilator in a fixed location, advance the sheath fully over the dilator into the left atrial cavity.

8. WITHDRAW THE BRK™ NEEDLE AND THE DILATOR.

CAUTION: There is a risk of air infiltration when withdrawing objects from the hemostasis valve of the sheath. Take precautions to prevent air infiltration by withdrawing objects slowly to prevent vacuum buildup in the sheath and fluoroscopically monitor the sheath during ensuing device insertion for the presence of air.

- Turn the needle stopcock to the off position and disconnect any attachments to the needle hub.
- Remove the BRK™ Transseptal Needle from the dilator. The needle may be cleaned and set aside for repeat use in this procedure. Otherwise, discard by appropriate means for contaminated sharp objects.
- Immediately attach a syringe to the dilator and aspirate. Continue aspirating blood while holding the sheath in position and withdrawing the dilator. The blood should be arterial blood.
- Once the dilator is removed, aspirate blood through the side arm of the sheath, and then flush it with heparinized saline, taking care to avoid air bubbles.
- The sheath is now in place in the left atrium.

NOTE: The symbols section contains all the symbols that may be used on product labels. Product is labeled as required.

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State: Released **Date:** 2020.09.08 14:25 GMT **Effectivity:** Upon Release

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Abbott Medical
5050 Nathan Lane North
Plymouth, MN 55442 USA
+1 651 756 5833

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600108988

CE
2797

Abbott Medical Peel-Away Introducer

English - Instructions for use 2
Deutsch - Gebrauchsanleitung 4
Français - Mode d'emploi 5
Italiano - Istruzioni per l'uso 6
Español - Instrucciones de uso 7
中文 – 使用说明 8
Ελληνικά - Οδηγίες χρήσης 9
Русский - Инструкция по применению 10
Magyar - Használati utasítás 12
Português - Instruções de Utilização 13
Nederlands - Gebruiksaanwijzing 14
Suomi - Käyttöohjeet 15
Svenska - Bruksanvisning 16
Dansk - Brugsanvisning 17
Norsk - Bruksanvisning 18
Polski - Instrukcja obsługi 19
Česky - Návod k použití 20
Slovensky - Návod na použitie 21
Eesti - Kasutusjuhend 22
Lietuvių - Naudojimo instrukcijos 23
Latviski - Lietošanas pamācība 24
Türkçe - Kullanım talimatları 25
한국어 – 사용안내 26
Български - Инструкции за употреба 27
Română - Instrucţiuni de utilizare 28
Srpski - uputstvo za upotrebu 29
العربية - تعليمات الاستخدام 31
SYMBOL 32

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Abbott

English

Abbott Medical Peel-Away Introducer Instructions for Use

See Product Catalog for Reorder Numbers.

THE FRENCH SIZE SPECIFIED REPRESENTS THE INNER DIAMETER OF THE INTRODUCER SHEATH.

DO NOT ATTEMPT TO INSERT A CATHETER HAVING A DISTAL TIP OR BODY SIZE LARGER THAN THE INTRODUCER SIZE INDICATED.

SEE INDIVIDUAL STERILE PACKAGE LABEL FOR CONTENTS.

SINGLE-USE DISPOSABLE MEDICAL DEVICE. CONTENTS ARE STERILE IF PACKAGE IS UNOPENED AND UNDAMAGED. DO NOT RESTERILIZE.

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CAUTION

Do Not alter this device.

Federal law (U.S.) restricts this device to sale by or on the order of a physician.

Do Not use forceps to break handle and/or peel the sheath and cause its withdrawal from the body.

The peel-away sheath should withdraw freely from the body. If resistance is encountered, do not force the peeling action of the sheath or its withdrawal. Make procedural corrections, then continue.

Do Not attempt to use a guidewire over maximum diameter specified on package label.

This device should only be used by physicians thoroughly trained in the technique of angiography.

Individual patient anatomy and physician technique may require procedural variations.

Insertion into artery may cause excessive bleeding and/or outer complications.

Damage to the guidewire may result if withdrawn through the needle.

The need for symmetrical peeling of the sheath cannot be overemphasized. If catheter withdraw is a problem when the sheath is peeled, completely withdraw the sheath from the introduction site while holding the catheter in position. Then peel the sheath symmetrically from the catheter by grasping both handle halves.

STORE IN A COOL, DARK, DRY PLACE.

PRECAUTIONS

This device is intended for one time use only. Do not reprocess or reuse. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another.

POSSIBLE COMPLICATIONS

Complications may occur during certain procedures but may not be limited to pneumothorax, hematoma formation, brachial plexus injury, sepsis/infection, excessive bleeding, vessel damage.

APPLICATION

The introduction of catheters into a vessel.

SUGGESTED DIRECTIONS FOR USE

Introducer

1. Follow normal accepted practice for vessel puncture and guidewire insertion.

NOTE: A syringe may be attached to the needle for added support during insertion and as an aid in observing venous return.

NOTE: Leave an appropriate amount of guidewire exposed. Exposed portions of the guidewire should be more than the length of the percutaneous catheter introducer assembly.

2. Advance the dilator sheath introducer assembly over the exposed guidewire into the vessel using a twisting motion.

 CAUTION: DO NOT ALLOW GUIDEWIRE TO INADVERTENTLY ADVANCE TOTALLY INTO PATIENT.

3. Withdraw the vessel dilator and guidewire, leaving the introducer sheath in position. Place thumb over the exposed orifice of the sheath to reduce air aspiration and/or blood loss.
4. Insert and advance catheter.
5. With catheter advanced, break handle of introducer sheath and peel sheath apart while withdrawing from the vessel. Care must be exercised **NOT** to withdraw the catheter during this process. (See Cautions.)

SUGGESTED DIRECTIONS FOR USE

Introducer Needle

1. Use proper aseptic technique.
2. Introduce needle slowly into vessel and observe venous return.
3. Remove the needle, leaving the introducer sheath in position. Place thumb over exposed orifice of the introducer sheath to reduce air aspiration and/or blood loss.

 CAUTION: DO NOT RE-INSERT NEEDLE WHILE SHEATH IS IN VESSEL.

4. Insert and advance catheter.
5. With catheter advanced, break handle of introducer sheath and peel sheath apart while withdrawing from the vessel. Care must be exercised **NOT** to withdraw the catheter during this process. (See Cautions.)

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USA ONLY

LIMITED WARRANTY AND DISCLAIMER

Abbott Medical hereby warrants that if any Abbott Medical product fails to perform within normal tolerances for a patient due to a defect in materials or workmanship, Abbott Medical will provide, at no charge, a replacement Abbott Medical product for the patient's use. This limited warranty applies only if each of the following conditions are met.

1. The product was packaged and labeled by Abbott Medical.
2. The failed product must be returned to Abbott Medical and becomes the property of Abbott Medical.
3. The product has not been mishandled, reprocessed or altered in any way.
4. The product was used before the "USE BY" date marked on the packaging of the product.

No representation or warranty is made that an Abbott Medical product will not fail. Abbott Medical disclaims responsibility for any medical complications, including death, resulting from the use of its products. Except as expressly provided by this limited warranty, ABBOTT MEDICAL IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE, OR MALFUNCTION OF ITS PRODUCTS, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE. Some states do not allow the exclusion or limitation of incidental or consequential damages however, so the above limitation or exclusion may not apply to you.

Except as expressly provided by the limited warranty, ABBOTT MEDICAL MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE.

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™ Indicates a trademark of the Abbott group of companies.
‡ Indicates a third-party trademark, which is property of its respective owner.
Pat. http://www.abbott.com/patents
© 2020 Abbott. All Rights Reserved.

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CE
2797

FAST-CATH™ DUO and FAST-CATH™ TRIO
HEMOSTASIS INTRODUCER

English - Instructions for use 2
Deutsch - Gebrauchsanleitung 4
Français - Mode d'emploi 5
Italiano - Istruzioni per l'uso 6
Español - Instrucciones de uso 7
中文 – 使用说明 8
Ελληνικά - Οδηγίες χρήσης 9
Русский - Инструкция по применению 10
Magyar - Használati utasítás 11
Português - Instruções de Utilização 12
Nederlands - Gebruiksaanwijzing 13
Suomi - Käyttöohjeet 14
Svenska - Bruksanvisning 15
Dansk - Brugsanvisning 16
Norsk - Bruksanvisning 17
Polski - Instrukcja obsługi 18
Česky - Návod k použití 19
Slovensky - Návod na použitie 20
Eesti - Kasutusjuhend 21
Lietuvių - Naudojimo instrukcijos 22
Latviski - Lietošanas pamācība 23
Türkçe - Kullanım talimatları 24
한국어 – 사용안내 25
Български - Инструкции за употреба 26
Română - Instrucţiuni de utilizare 27
Srpski – uputstvo za upotrebu 28
العربية - تعليمات الاستخدام 29
SYMBOL 30

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Abbott

English

FAST-CATH™ DUO and FAST-CATH™ TRIO HEMOSTASIS INTRODUCER

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Refer to the individual package label for the appropriate catheter French sizes and catheter combinations the FAST-CATH™ DUO or FAST-CATH™ TRIO Hemostasis Introducer will accommodate.

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The closing force of the hemostasis valve may alter or impair the function of some catheters.

Damage to guidewire may result if withdrawn through a metal needle cannula. The cannula must be removed first.

Do Not attempt to advance or withdraw guidewire if resistance is felt. Use flouroscopy to determine the cause.

Damage to the valve assembly may occur if inner catheter is withdrawn rapidly.

Before proceeding, verify that the desired catheter/medical device is of proper size. Refer to the individual package label of the appropriate FAST-CATH™ DUO or FAST-CATH™ TRIO Hemostasis Introducer.

Advance dilator/sheath assembly with a twisting motion to avoid damage to the sheath or vessel.

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APPLICATIONS

The introduction of angiographic catheters, closed end catheters, balloon catheters, and electrodes into a vessel where minimizing blood loss is essential.

SUGGESTED DIRECTIONS FOR USE

INSPECT SHEATH INTRODUCER, FAST-CATH™ DUO OR FAST-CATH™ TRIO HEMOSTASIS ADAPTER ASSEMBLY AND ACCESSORIES FOR DEFECTS. DO NOT USE ANY DEFECTIVE DEVICE.

FLUSH THROUGH ALL SIDEPORTS OF THE ABBOTT MEDICAL FAST-CATH™ DUO OR ABBOTT MEDICAL FAST-CATH™ TRIO HEMOSTASIS ADAPTER WITH HEPARINIZED NORMAL SALINE.

1. Follow normal accepted practice for vessel puncture and guidewire insertion, vessel dilator and introducer use.
2. Advance the dilator sheath introducer assembly over the exposed guidewire into the vessel, using a twisting motion.

 CAUTION: DO NOT ATTEMPT ADVANCEMENT WITHOUT GUIDEWIRE. SEVERE VASCULAR DAMAGE AND/OR INJURY MAY OCCUR.

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 CAUTION: INSERTION INTO ARTERY MAY CAUSE EXCESSIVE BLEEDING AND/OR OTHER COMPLICATIONS.
3. Withdraw the vessel dilator and guidewire, leaving the introducer sheath in position. Place thumb over the exposed orifice of the sheath to reduce air aspiration and/or blood loss.
4. Attach Abbott Medical FAST-CATH™ DUO or Abbott Medical FAST-CATH™ TRIO hemostasis adapter to introducer sheath.

 CAUTION: ASPIRATE THROUGH ALL SIDEPORTS OF THE ABBOTT MEDICAL FAST-CATH™ DUO OR ABBOTT MEDICAL FAST-CATH™ TRIO HEMOSTASIS ADAPTER, INTRODUCER SHEATH ASSEMBLY TO REMOVE ALL AIR PRIOR TO INSERTION OF A DEVICE.
5. Follow manufacturer's recommendations for the catheter or device being introduced via the hemostasis introducer.
6. Follow normal practice of using a continuous drip of anticoagulant fluid under pressure through the sideports when the hemostasis adapter/introducer assembly is in the vessel.

 CAUTION: UPON WITHDRAWAL OF CATHETER OR OTHER INSERTED DEVICE, ASPIRATION FROM EACH PORT OF THE ADAPTER IS RECOMMENDED TO REMOVE ANY FIBRIN ACCUMULATION.

 CAUTION: UPON REMOVAL OF HEMOSTASIS ADAPTER AND/OR INTRODUCER ASSEMBLY, PROPER PRECAUTIONS SHOULD BE TAKEN TO PREVENT BLEEDING.

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600108996

CE
2797

FAST-CATH™
TRANSSEPTAL GUIDING INTRODUCER

English - Instructions for use 1
Deutsch - Gebrauchsanleitung 4
Français - Mode d'emploi 7
Italiano - Istruzioni per l'uso. 10
Español - Instrucciones de uso 13
中文 – 使用说明 16
Ελληνικά - Οδηγίες χρήσης 19
Русский - Инструкция по применению. 22
Magyar - Használati utasítás 25
Português - Instruções de Utilização 28
Nederlands - Gebruiksaanwijzing 31
Suomi - Käyttöohjeet 34
Svenska - Bruksanvisning 37
Dansk - Brugsanvisning. 40
Norsk - Bruksanvisning 43
Polski - Instrukcja obsługi 46
Česky - Návod k použití 49
Slovensky - Návod na použitie. 52
Eesti - Kasutusjuhend 55
Lietuvių - Naudojimo instrukcijos 58
Latviski - Lietošanas pamācība 61
Türkçe - Kullanım talimatları 64
한국어 – 사용안내 67
Български - Инструкции за употреба. 70
Română - Instrucţiuni de utilizare 73
Srpski - uputstvo za upotrebu 76
العربية - تعليمات الاستخدام 79
SYMBOL 82

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FAST-CATH™ TRANSSEPTAL GUIDING INTRODUCER Instructions for Use

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Read Instructions for Use prior to use of this device.

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Single-use disposable medical device. Contents are sterile if package is unopened and undamaged. Do not re-sterilize.

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PACKAGE CONTENTS

All FAST-CATH™ Transseptal Guiding Introducer kits consist of an introducer and dilator. Many kits also contain a guidewire.

The FAST-CATH™ introducer and accessories are supplied sterile in a single sterile barrier. The device is packaged in a tray which is placed inside a Mylar‡ and Tyvek‡ pouch. The pouch is heat-sealed to create a sterile barrier.

DESCRIPTION

The FAST-CATH™ Transseptal Guiding Introducer set consists of a radiopaque introducer and dilator, each with a specially curved distal portion to accommodate positioning against the atrial septum. The introducer is also compatible with a .032" guidewire and BRK™ type curved puncture needle. The introducer is fitted with a valve to provide hemostasis during catheter introduction and/or exchange over a guidewire. Located on the valve housing, a sideport with three-way stopcock is provided for aspiration, fluid infusion, blood sampling and pressure monitoring. The dilator is tapered at the distal tip with an internal lumen designed to accept ancillary devices (e.g., needles or guidewires) that have a maximum diameter of 0.032 inches. The dilator inner lumen is also designed with a special geometry at its distal end to limit the exposure of the BRK™ needle.

INDICATIONS FOR USE

The FAST-CATH™ Transseptal Guiding Introducer set is indicated for introducing various cardiovascular catheters into the left side of the heart through the interatrial septum.

CONTRAINDICATIONS

- Previous interatrial septal patch or prosthetic atrial septal defect closure device
- Any previous thromboembolic event
- Known or suspected left atrial myxoma
- Known or suspected myocardial infarction within the last two weeks
- Unstable angina
- Recent cerebral vascular accident (CVA)
- Patients who do not tolerate anticoagulation therapy
- Patients with an active infection

WARNINGS

- Do not alter this device in any way.
- Do not reuse this device. After use, thorough cleaning of biological and foreign material is not possible. Adverse patient reactions may result from reuse of this device.

PRECAUTIONS

- CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician.
- Store in a cool, dark, dry place.
- This device is intended for one time use only. Do not reprocess or reuse. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another.

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- Rotated heart
- Enlarged aortic root
- Marked right atrial enlargement
- Scoliosis/kyphosis
- Abnormal left atrial geometry
- Congenital malformations
- Vascular malformations
- Inability to access the right atrium through the inferior vena cava

Procedural Considerations

- Carefully reading the Instructions before use of this device will help to reduce the potential risks associated with the transseptal catheterization technique such as air embolism or perforation of the aorta or left atrium.
- Only those physicians who are specially trained in transseptal procedures should use this device.
- Fluoroscopy should be used to confirm positioning throughout the procedure.
- Transseptal procedures should be performed only in facilities appropriately equipped and staffed to perform such procedures. Lab capabilities should include, but are not limited to:
 - Intracardiac pressure monitoring capabilities
 - Systemic pressure monitoring
 - Contrast media injection and management of untoward reactions to contrast media
 - Pericardiocentesis
 - Surgical backup
 - Anticoagulation therapy and monitoring
- Maintain monitoring of vital signs throughout the procedure.
- Inspect all components before use.
- During insertion, check for excessive resistance as the tip of the needle advances through the curvature of the sheath/dilator assembly. If there is any resistance to needle advancement, retract the needle, and assess components.
- During insertion, use caution not to create excessive bends in this device. This may inhibit advancement of the needle and may result in inadvertent needle puncture of the dilator/sheath assembly.
- During insertion, always use the stylet to facilitate needle passage through the dilator/sheath assembly. Failure to use the stylet may inhibit advancement of the needle and may result in inadvertent needle puncture of the dilator/sheath assembly or skiving of material from the inner surface of the dilator.

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- Cardiac perforation
- Cardiac tamponade
- Conduction system disturbances such as SA node, AV node or His-Purkinje system block
- Hematoma or excessive bleeding at the vascular access site
- Reaction to contrast medium
- Stroke
- Thromboembolism
- Valvular damage
- Cardiac arrhythmias

DIRECTIONS FOR USE - TRANSSEPTAL PROCEDURE

1. PREPARE AND ASSEMBLE EQUIPMENT.

- Remove contents from packaging.
 NOTE: Do not assemble the sheath and dilator.
- Flush the dilator and sheath with heparinized saline.
- Close the sheath stopcock.
- Remove the stylet from the needle and flush the needle with heparinized saline.
- Insert the stylet into the needle and lock it onto the hub.
- Insert the dilator fully into the transseptal sheath.

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- Obtain femoral venous access (right femoral preferred). OPT: A larger (≥ 2.5 French sizes greater than the transseptal sheath) standard length sheath may be used to obtain and maintain venous access for device exchange and hemostasis.
- Introduce a 0.032", 150 -180 cm, 3 mm "J" tip guidewire into the superior vena cava (SVC).
 NOTE: 0.032" is the maximum guidewire diameter that can be used with the dilator.
- Insert the sheath/dilator into the vein over the guidewire until the sheath tip is in the SVC. Orient the dilator tip medially.

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


4. ENGAGE THE FOSSA OVALIS.

- Visualize and identify anatomic landmarks.
- Set the fluoro unit to an appropriate angle, parallel to the plane of the mitral valve and orthogonal to the plane of the septum. This will typically be LAO, approximately 30 to 40 degrees.
- OPT: During electrophysiology procedures, the coronary sinus and His bundle catheter positions can serve as useful anatomic landmarks. In the appropriate LAO view, the coronary sinus catheter will be seen in profile. In the appropriate right anterior oblique (RAO) view the His bundle catheter will appear in profile. The fossa ovalis is located at or slightly below the level of His bundle catheter and superior and posterior to the coronary sinus ostium.
- OPT: Placing a pigtail angiographic/hemodynamic monitoring catheter in the non-coronary cusp of the aortic valve can serve as a useful anatomic landmark.
- OPT: Observe the pressure waveform being recorded through the BRK™ Transseptal Needle.
- Adjust the pointer flange so that the needle is perpendicular to the fossa ovalis (typically between 3:00 and 5:00 o'clock, as viewed from the foot end of the patient). (See Fig 5.)

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