Carl Zeiss Vision International GmbH P.O. Box 1865 73408 Aalen

to be presented to the competend authorities in Republic of Belarus

Carl Zeiss Vision International GmbH

Turnstr. 27
73430 Aalen
Germany

Phone: +49 73 61 591 - 310

Fax: +49 73 61 591 - 734

E-Mail: herbert.krug@zeiss.de

Division/Dept.: VIS-QR-DE
Your contact: H. Krug

Yours ref.:
Yours of:
Our ref:
Date: 11.11.2020

Instructions for Use

We

Carl Zeiss Vision International GmbH, Turnstrasse 27, 73430 Aalen, Germany

HEREBY CERTIFY that all

Product:	Spectacle lens Brand: ZEISS
Product type:	Single Vision lenses, Bi- and Multi-focal lenses, Progressive lenses (Uncoated, hard-coated or anti-reflex-coated)
Commercial Name/Variants	see Declaration of Conformity
Manufacturer:	Carl Zeiss Vision International GmbH
Manufacturing sites:	see Declaration of Conformity

follow the same characteristics in terms of instructions for use

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The product is to be used as per the instruction of Eye Care Professional based on the patients' prescription.

The intended use of ZEISS ophthalmic spectacle lenses is vision correction. Unless specifically ordered as safety glasses, they do not claim to protect the eye from physical impact because they are neither shatterproof nor unbreakable and are not intended for use as personal protective equipment.

Do not clean spectacle lenses with clothing, towels, paper or household cleaners. Use ZEISS cleaning products or just water, a pH-neutral cleaning agent and a microfiber or soft cotton cloth. Dry cleaning with an ordinary cloth or household paper can cause damage to the lens surface. Only use dedicated lens cleaning cloth.

Organic lenses should not be exposed to short-time high temperature fluctuations; this can cause stress within the material and coating damage may result.

Short exposure time to temperatures above 80 °C (e.g. sauna visit), or 60 °C of long exposure time (e.g. on the dashboard of a car) can cause damage to the lens / coating.

Lenses with less than 8% transmission are not suitable for driving.

Lenses with less than 25% transmission are not suitable for driving at night.

If treated carefully avoiding damage of the lens, the service life (life-time, durability) of a lens is only limited by a potential change of the prescription of the wearer resulting in the need for a new prescription (lens).

Authorised Signatory:

Name: Sabine Bock

Job Title: Director Quality EMEA, QP&CI

Signature: _____

Date: 11.11.2020

Carl Zeiss Vision GmbH ZEISS Regulatory Affairs

Name: Herbert Krug

Job Title: Director Regulatory Affairs

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