2014/NOV/11 at 10:20 a.m. Doc number: M053018T001 Rev. 1A [Uni+]
Printspec a5-55 - 148 x 210 mm

Medtronic

IN.PACT™ Falcon™

Paclitaxel-eluting PTCA Balloon Catheter 4
Излъчващ паклитаксел балонен катетър за перкутанна транслуминална коронарна ангиопластика (PTCA) 6
Balónkový katétr pro PTCA s uvolňováním paclitaxelu 8
Paclitaxel-eluerende PTCA-ballonkateter 11
Paclitaxel-freisetzender PTCA-Ballonkatheter 13
Καθετήρας μπαλονιού PTCA (διαδερμικής διαυλικής στεφανιαίας αγγειοπλαστικής) έκλυσης πακλιταξέλης 16
Catéter de balón liberador de paclitaxel para ACTP 18
Paklitakseeli eritav PTKA balloonkateeter 21
Paklitakselia vapauttava PTCA-pallokatetri 23
Cathéter d'ACTP à ballonnet à élution de paclitaxel 25
PTCA balonski kateter koji ispušta paklitaksel 28
Paclitaxelt kibocsátó PTCA ballonkatéter 30
Catetere a palloncino per PTCA a eluizione di paclitaxel 32
파클리탁셀-방출 PTCA 풍선 카테터 35
paklitakselį išskiriantis PTCA balioninis kateteris 37
PTCA balona katetrs ar paklitaksela pārklājumu 39
Балон катетер што ослободува паклитаксел за PTCA (перкутана транслуминална коронарна ангиопластика) што ослободува паклитаксел 41
Paclitaxel-afgevende PTCA-ballonkatheter 44
Paclitaxel-avgivende PTCA-ballongkateter 46
Cewnik balonowy do PTCA uwalniający paklitaksel 48
Cateter de balão com eluição de paclitaxel para PTCA 51
Cateter de balão com eluição de paclitaxel para PTCA 53
Cateter cu balon pentru ACTP (angioplastie coronariană transluminală percutanată) cu eluţie de paclitaxel 56
Баллонный катетер для ЧТКА с покрытием, содержащим паклитаксел 58
Balónový katéter uvoľňujúci paclitaxel na perkutánnu transluminálnu koronárnu angioplastiku (PTCA) 63
Balonski kateter za perkutano transluminalno koronarno angioplastiko, ki izloča paklitaksel 61
PTCA balon-kateter koji otpušta paklitaksel 66
Paklitaxel-avgivande PTCA-ballongkateter 68
Paklitaksel salınımlı PTCA Balon Kateteri 70
紫杉醇洗脱 PTCA 球囊导管 72

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning
Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet
Mode d'emploi • Upute za upotrebu • Használati útmutató • Istruzioni per l'uso • 사용 지침
Naudojimo instrukcijos • Lietošanas pamācība • Упатства за употреба • Gebruiksaanwijzing
Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instruções de utilização
Instrucţiuni de utilizare • Инструкция по эксплуатации • Pokyny na používanie
Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • 使用说明

2014-07-14
M053018T001 Rev. 1A

CE 0086

multilang2col28:20140514, alldocMod:20140905, alldocStyle:20130225, master_pageset:20130225, languageLookup:20140514, table:20140620, characterJoin:20110211, 3of9:20110211, 3of9_svg:20110211

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Paclitaxel-eluting PTCA Balloon Catheter

1. Description

IN.PACT™ Falcon™ is a paclitaxel-eluting, percutaneous transluminal coronary angioplasty (PTCA) balloon catheter manufactured by Medtronic, Inc.

The balloon is coated with FreePac™, a paclitaxel eluting formulation that is intended to reduce the proliferative response that is associated with restenosis. Paclitaxel stabilizes microtubules to reduce cell proliferation.

The IN.PACT Falcon features a rapid exchange (RX) design. It is composed of a single-lumen proximal shaft, a dual-lumen distal shaft, and a balloon close to the catheter tip. The proximal shaft consists of a stainless steel hypotube with a PTFE coating, and a proximal end luer-lock connector hub for balloon inflation. At the opposite side, a special transition construction provides an optimal push-torque transmission through the full catheter length. The first lumen of the distal shaft is dedicated to guidewire passage while the other, which continues through the proximal shaft up to the hub, is dedicated to the inflation of the balloon.

The guidewire lumen permits the use of guidewires to facilitate advancement of the catheter to and through the stenosis to be dilated. The guidewire access port is located at the distal tip of the catheter; the guidewire exit port is located 23 cm from the distal tip. Maximum guidewire diameter is 0.014 in (0.36 mm). The balloon is designed to reach specific diameters at specific pressures (see compliance chart included in the package).

Two radiopaque markers are available in order to correctly position the balloon under fluoroscopy. The markers are placed on the shaft and under the balloon, defining its cylindrical area. A flushing needle with a luer port is provided in the sterile packaging to facilitate the flushing of the guidewire lumen prior to usage. The IN.PACT Falcon does not contain natural rubber latex; however, during the manufacturing process, it may have incidental contact with latex.

2. Indications

The IN.PACT Falcon is intended to be used in the native coronary arteries.

It is indicated for the dilatation of stenotic portions, including total occlusions and acute myocardial infarction. It is also indicated for the postdilatation of stents.

3. Contraindications

IN.PACT Falcon is contraindicated for use in:

- The presence of artery spasms in absence of a significant stenosis
- Peripheral, renal, and cerebrovascular arteries
- Lesions of the left main coronary artery, internal mammary artery, aortic ostium, and saphenous vein grafts
- Pregnant or breast-feeding women
- Patients with known allergies or hypersensitivities to paclitaxel

4. Warnings

- This device is designed and intended for single use only. DO NOT RESTERILIZE AND/OR REUSE. The drug coated surface is effective for a single use only. Reuse or resterilization may create a risk of contamination of the device and/or cause patient infection or cross-infection, including the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. Reuse or resterilization may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness, and death. Medtronic will not be responsible for any direct, incidental or consequential damages resulting from resterilization or reuse.
- **Drug-eluting stents must not be implanted into the vessel segment that has been treated with an IN.PACT Falcon balloon.**
- Inspect the device prior to procedure to verify that the product is intact and functional. Do not use the device if the outer or the inner packaging is damaged or opened.
- When the catheter is in the body, it should be manipulated while under sufficient and/or high quality fluoroscopy. Prior to withdrawing the balloon catheter from the lesion, the balloon must be completely deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Do not use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium (equal parts contrast medium and saline solution).
- Do not expose the device to organic solvents, eg, alcohol.
- Do not exceed the Rated Burst Pressure (RBP). Use a pressure monitoring device to prevent over pressurization. Use of pressures higher than those specified on the product label may result in a ruptured balloon, causing possible intimal damage and dissection.
- Only physicians thoroughly trained and educated in the performance of PTCA / PTA should use this device. Physicians should keep themselves updated regarding recent publications about PTCA / PTA techniques.
- To reduce the risk of vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- Use by the labeled Use By date.
- Do not use with Lipiodol™ or Ethiodol™ contrast media (or other such contrast media that incorporate components of these agents).
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration. This includes possible hemodynamic support during PTCA for this patient population with special risk.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.

5. Precautions

- Administer appropriate drug therapy (anticoagulant, vasodilator, etc.) to the patient according to standard protocols for PTCA before inserting the IN.PACT Falcon.
- Exercise care during handling in order to avoid possible damage to the catheter. Avoid acute bending or kinking of the catheter. Do not use a catheter that has been damaged.
- To minimize the introduction of air, aspirate and flush the system and keep a tight catheter connection throughout the procedure.
- Take precautions to prevent or reduce clotting when any catheter is used. Flush or rinse all products entering the vascular system with sterile isotonic saline or a similar solution via the guidewire access port prior to use. Consider the use of systemic heparinization.
- Use with caution for procedures involving calcified lesions due to the abrasive nature of these lesions.
- Identify allergic reactions to contrast media, antiplatelet therapy, balloon catheters, and FreePac coating before treatment.
- Catheter applications vary. Select the technique on the basis of the patient's condition and the experience of the interventionalist.
- Never advance the IN.PACT Falcon without the guidewire extending from the tip.
- Store at controlled room temperature in a dry place. Keep away from sunlight.
- Use a pressure-monitoring device to prevent overpressurization.

 Caution: Larger sizes of the IN.PACT Falcon may exhibit slower deflation times, particularly on long catheter shafts.
- After use, this product may be a biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable hospital, administrative, and government regulations.

6. Potential Complications/Adverse Effects

Complications associated with the use of IN.PACT Falcon are similar to the ones associated with standard PTCA / PTA procedures. Possible complications may include, but are not limited to:

- Death
- Stroke and / or TIA
- Emergency CABG
- Acute myocardial infarction
- Allergic reaction to drug and/or to contrast medium
- Hemorrhage or hematoma
- Restenosis of the dilated vessel
- Vessel dissection, perforation, rupture or injury
- Arrhythmias, including ventricular fibrillation
- Total occlusion of the artery or bypass graft
- Unstable angina
- Pain and tenderness
- Endocarditis
- Hypo/hypertension
- Infection
- Arteriovenous fistula
- Pseudo-aneurysm
- Artery spasm
- Thrombosis and / or embolization

multilang2col28:20140514, alldocMod:20140905, alldocStyle:20130225, master_pageset:20130225, languageLookup:20140514, table:20140620, characterJoin:20110211, 3of9:20110211, 3of9_svg:20110211

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CE 0086

INVATEC Innovative Technologies
We make ideas come alive

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
EN

DEVICE NAME
IN.PACT™ ADMIRAL is a paclitaxel-eluting Balloon Dilatation Catheter manufactured by Invatec Technology Center GmbH.

DESCRIPTION
The IN.PACT™ ADMIRAL is an Over the Wire (OTW) peripheral balloon catheter, specifically designed for percutaneous transluminal angioplasty in atherosclerotic obstructed vessels. The catheter has a dual-lumen shaft. The double-lumen shaft is branched at the proximal end so that one tube forms the entrance to the central lumen for the guide wire, while the other tube is used to inflate and deflate the dilatation balloon with a mixture of contrast medium and saline solution. The catheter construction and the special balloon material are designed so that a specific balloon diameter can be reached, depending on the balloon size and defined pressure. The length of each balloon is specified. Maximum guide wire diameter is 0,035" (0,89 mm).
The IN.PACT™ ADMIRAL is available in different balloon sizes. Nominal balloon diameter and lengths are printed on the hub.

INDICATIONS
The IN.PACT™ ADMIRAL is indicated for percutaneous transluminal angioplasty (PTA) in patients with obstructive disease of peripheral arteries.

CONTRAINDICATIONS
- The IN.PACT™ ADMIRAL PTA catheter is contraindicated for the use in coronary arteries and supra-aortic/cerebrovascular vessels.
- Inability to cross lesion with a guide wire.
- The IN.PACT™ ADMIRAL should not be used in pregnant or breast feeding women or in patients with known hypersensitivity to paclitaxel.

WARNINGS
- This device is designed and intended for single use only. DO NOT RESTERILIZE AND/OR REUSE. Reuse or resterilization may create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Reuse or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness and death. INVATEC will not be responsible for any direct, incidental or consequential damages resulting from resterilization or reuse. The coating excludes any reuse of the device.
- Inspect the device, prior to procedure, to verify functionality and lack of damaged parts. Do not use the device if the outer or the inner package is damaged or opened.
- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just distal to the stenosis.
- When the catheter is in the body, it should be manipulated while under sufficient and/or high quality fluoroscopy. Prior to withdrawing the balloon catheter from the lesion, the balloon must be completely deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Never use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium.
- Do not expose the device to organic solvents, e.g. alcohol.
- Do not manipulate the PTA balloon in inflated state. The position of the PTA balloon may only be changed with the guidewire in place.
- If resistance occurs during manipulation, the cause must first be ascertained by fluoroscopy, road mapping or DSA before the balloon catheter is moved backwards or forwards.
- The guidewire may under no circumstances be moved during inflation of the PTA balloon.
- The balloon must be completely deflated before retrieving the catheter from the vascular system.
- Do not exceed the Rated Burst Pressure (RBP). The RBP is based on the results of in vitro testing. At least 99.9% of the balloon (with 95% confidence) will not burst at or below their RBP. Use of a pressure monitoring device is recommended to prevent over pressurization. Inflation in excess of the rated burst pressure may cause the balloon to rupture.
- Only physicians thoroughly trained and educated in the performance of PTA should use this device. Physicians should keep themselves informed and updated on recent publications about PTA techniques.
- To reduce the risk of vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just distal to the stenosis.
- Use only a mixture of contrast medium and saline solution to fill the balloon (1:1). Never use air or any gaseous medium to inflate the balloon of the catheter.
- Use by the labeled date.
- Do not use with Lipiodol™ or Ethiodol™ contrast media (or other such contrast media which incorporates the components of these agents).
- Never apply positive pressure to the balloon during preparation.

PRECAUTIONS
- Appropriate drug therapy (anticoagulant, vasodilator, etc.) should be administered to the patient according to standard protocols for PTA before insertion of the dilatation catheter.
- To minimize the possible introduction of air into the system, it is imperative that before proceeding careful attention is paid to the maintenance of tight catheter connections and through aspiration and flushing of the system.
- Precautions to prevent or reduce clotting should be taken when any catheter is used. Flush or rinse all products entering the vascular system with sterile isotonic saline or a similar solution via the guidewire access port prior to use. Consider the use of systemic heparinization.
- The IN.PACT™ ADMIRAL PTA balloon should be used with caution for procedures involving calcified lesions due to the abrasive nature of these lesions.
- Allergic reactions to contrast medium should be identified before treatment.
- Catheter applications vary and the technique must be selected on the basis of the patient's condition and the experience of the interventionalist.
- Never advance the angioplasty catheter without the guidewire extending from the tip.
- Never attempt to move the guidewire when the balloon is inflated.
- Do not advance the catheter against significant resistance. The cause of resistance should be determined via fluoroscopy and remedial action taken.
- Store at controlled room temperature, in a dry place. Keep away from sunlight.
- Don't exceed dilatation pressures to "rated burst pressure". In any case, always use a manometer controlled inflation.

CAUTION: Larger sizes of the IN.PACT™ ADMIRAL PTA balloon catheter may exhibit slower deflation times particularly on long catheter shafts.

POTENTIAL COMPLICATIONS/ADVERSE EFFECTS
Complications associated with the use of the IN.PACT™ ADMIRAL balloon catheter are similar to the ones associated with standard PTA procedures. Possible complications may include, but are not limited to:
Puncture related
- Local hematoma
- Local hemorrhage
- Local or distal thromboembolic episodes
- Thrombosis
- Arterio-venous fistula
- Pseudoaneurysm
- Local infections

Dilatation related
- Dissection in the dilated artery wall
- Perforation of the artery wall
- Prolonged spasms
- Acute re-occlusion necessitating surgical intervention
- Restenosis of the dilated artery
- Total occlusion of the peripheral artery

Angiography related
- Hypotension
- Pain and tenderness
- Arrhythmias
- Sepsis/infection
- Systemic embolization
- Endocarditis
- Short-term hemodynamic deterioration
- Death
- Drug reactions
- Allergic reaction to contrast medium
- Pyrogenic reaction

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PHARMACOLOGICAL INTERACTION
Due to the low dosage and local administration pharmacological interactions are not to be expected. Pharmacological interactions have not been reported.

INSTRUCTION FOR USE
Prior to angioplasty, carefully examine all equipment to be used during the procedure, including the dilatation catheter, to verify proper function. Verify that the catheter and sterile packaging have not been damaged in shipment and that the catheter size is suitable for the specific procedure for which it is intended.
Handle the device with extreme caution in order to avoid any damage to the folded balloon. Moistening of the device with saline, insertion through the sheath and contact with blood will not reduce the efficacy of the FreePac™ coating as long as the balloon remains uninflated.
1. Inflation Device Preparation
a. Prepare the inflation device according to the manufacturer's instructions.
b. Carefully remove the balloon protective sheath. Prior to use, carefully examine the unit to verify that the catheter or sterile package has not been damaged in shipment. In case of resistance, twist the protection with one hand holding the shaft with the other.
2. Dilatation Catheter Selection
a. The nominal balloon size must be chosen equal to the inner diameter of the artery distal to the lesion. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter (non drug coated) catheter to pre-dilate the lesion to facilitate the passage of a more appropriate-sized dilatation catheter.
b. In cases of total occlusion or subocclusive lesions a predilatation with a (non drug coated) balloon is recommended.
3. Dilatation Catheter Preparation
a. The catheter is packaged in a protective hoop; carefully remove the catheter from the package.
b. The folded balloon catheter contains tiny air bubbles that should be purged prior use. To do this, keep the balloon catheter in a downward vertical position. Keep the protective sheath in place during the purging procedure. Connect a Luer-lock syringe partially filled with saline-contrast mixture solution to the inflation port of the catheter hub. Apply negative pressure for 15 seconds until air is completely evacuated and release the plunger. Repeat this operation a couple of times until migration of air bubbles towards the syringe stops. **Never use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium. Never apply positive pressure to the balloon during preparation.**
c. Remove the stylet and the protective sheath from the balloon.
d. Flush the wire lumen properly through the connector luer hub.
4. Inflation Device connection to Catheter
a. To remove any air lodged in the distal Luer fitting of the inflation device, purge approximately 1 ml (cc) of contrast medium.
b. With the stopcock in the closed position, disconnect the syringe used in preparation applying a slight positive pressure. A meniscus of contrast medium will appear in the balloon port when the syringe is removed. Verify that a meniscus of contrast medium is evident in both the dilatation catheter balloon port (hub) and the inflation device connection. Securely couple the inflation device to the balloon port of the balloon dilatation catheter.
5. Use of the IN.PACT™ ADMIRAL PTA Balloon Angioplasty Catheter
a. Insert a guidewire through the hemostatic valve following the manufacturer's instructions or standard practice. Advance the guidewire carefully into the guiding catheter/ introducer sheath. When complete, withdraw the guidewire/ introducer, if used.
b. Attach a torque device to the wire, if desired. Under fluoroscopy, advance the guidewire to the desired vessel, then across the stenosis.
c. Thoroughly aspirate and flush the guiding catheter in preparation (when guiding catheter is used) for introduction of the dilatation catheter. Back load the distal tip of the dilatation catheter onto the guidewire.
NOTE: To avoid kinking, advance the dilatation catheter slowly, in small increments until the proximal end of the guidewire emerges from the catheter.
d. Advance the catheter through the hemostatic valve slowly, while the balloon is fully deflated. It should be observed that the hemostatic valve is only closed as much as to prevent blood return yet permitting easy movements of the dilatation catheter. If resistance is encountered, do not advance the catheter through the adapter.
e. Under fluoroscopy, use the balloon radiopaque markers to position the balloon within the lesion to be dilated and inflate the balloon to the appropriate pressure (refer to balloon compliance table).
f. The majority of the drug is released within the first 30 seconds of balloon inflation. The duration of the inflation should therefore be between 30 seconds and 1 minute for optimal drug release. In order to optimize lesion dilatation, longer inflation times are possible at the discretion of the operator.
g. Maintain negative pressure on the balloon between inflations.
h. Completely deflate the balloon catheter. Withdraw the deflated dilatation catheter and guidewire from the guiding catheter/ introducer sheath, through the hemostatic valve. Tighten the knurled knob on the hemostatic valve.
i. If necessary the balloon catheter can be exchanged over the guidewire which remains in the vessel, for different balloon types or sizes.
j. WARNING: The same vessel segment must not be dilated with more than one IN.PACT™ ADMIRAL. Should a postdilatation be necessary after the utilization of an IN.PACT™ ADMIRAL balloon, so should this be done with a standard (non drug coated) dilatation balloon. Drug Eluting Stents must not be implanted into the vessel segment that has been treated with an IN.PACT™ ADMIRAL balloon.
k. If more than one balloon is required to treat a single lesion balloons must overlap by at least 1 cm. A second new balloon catheter has to be used because the drug is almost completely released during the first expansion. Additional overlapping should be avoided.
6. Pre- and postprocedure antiplatelet regimen
a. Dual antiplatelet therapy (ASS and clopidrogel or ticlopidine) should be administered pre-procedure and for a minimum of 4 weeks after the intervention. Prolonged antiplatelet therapy can be given at the discretion of the physician and should be considered after placement of stents.

HOW SUPPLIED
The IN.PACT™ ADMIRAL PTA balloon catheter is supplied sterile and intended for single use only. The IN.PACT™ ADMIRAL is sterilized by ethylene oxide gas. It will remain sterile as long as the packaging remains unopened and undamaged. Use product prior to labeled Expiration Date.
CAUTION: Do not use if the inner package is open or damaged.

STORAGE
Store at controlled room temperature, in a dry place. Keep away from sunlight. Do not expose to organic solvents (e.g. alcohol), ionizing radiation or ultraviolet light. Rotate inventory so that catheters are used prior to the expiration date on package label.

WARRANTY/LIABILITY
The product and each component of its system have been designed, manufactured, tested and packaged with all reasonable care. The warnings contained in Invatec's instructions for use are expressly to be considered as an integral part of this provision. Invatec guarantees the product until the expiry date indicated on the same. The warranty is valid provided that the use of the product was consistent with the instructions for use. Invatec disclaims any warranty of merchantability or fitness for a particular purpose of the product. Invatec is not liable for any direct, indirect, incidental or consequential damages caused by the product. Except in the case of fraud or grave fault on Invatec's part, compensation of any damage to the buyer will not, in any event, be greater than the invoice price of the disputed products. The guarantee contained in this provision incorporates and substitutes the legal guarantees for defects and compliance, and excludes any other possible liability of Invatec, however originating, from the product supplied. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects Invatec's legitimate interest in limiting its liability or warranty. No person has any authority to bind Invatec to any warranty or liability regarding the product.
Due to biological variability between different individuals no product can be 100% effective under all conditions. For this reason and as Invatec has no influence on the diagnosis of the patients, the methods applied, and the handling of the device after it has left the Invatec warehouse, Invatec does not guarantee the effectiveness nor the absence of complications associated with its use.

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© 2014 Medtronic, Inc.