© DePuy Synthes 2018-2021

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CE
2797

DePuy Synthes

PART OF THE Johnson & Johnson FAMILY OF COMPANIES

SCREWS

Revised 2021-02
© DePuy Synthes 2018-2021. All rights reserved.

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en INSTRUCTIONS FOR USE

SCREWS

IMPORTANT INFORMATION – Please read before use

IMPORTANT MEDICAL INFORMATION

The use of orthopaedic screws provides the orthopaedic surgeon a means of ancillary bone fixation in reconstructive surgeries. These implants are intended to provide additional stability to implanted orthopaedic devices, and are NOT intended to replace normal body structure or bear the weight of the body in the presence of incomplete bone healing.

DePuy's Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use. Reuse can potentially compromise device performance and patient safety.

INDICATIONS:

DePuy Bone Screws are indicated for use as an ancillary fixation device to be used with tibial and acetabular components in total knee and hip arthroplasty.

Refer to the Instructions for Use packaged with total knee and hip arthroplasty components for indications specific to those systems.

CONTRAINDICATIONS:

Refer to the Instructions for Use packaged with total knee and hip arthroplasty components for contraindications specific to those systems.

WARNINGS:

- This device is not approved for screw attachment of fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.
- Use of other manufacturer's components with this implant is not advised. Use of components other than those recommended could lead to loosening, wear, fracture during assembly, and premature failure.

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POSSIBLE ADVERSE EFFECTS:

- Loosening, bending, cracking or fracture of the screw or loss of fixation in bone.
- Infection, both deep and superficial.
- Allergies and other reaction to the device material.

PRECAUTIONS:

1. Preoperative and operative procedures, including knowledge of surgical techniques, good reduction and proper selection and placement of the implants are important considerations in the successful utilization of ancillary fixation devices. See the surgical technique for specific surgical procedures.
2. In evaluating patients for orthopaedic appliance application the patient's weight, occupation, activity level and any degenerative diseases are of extreme importance to the eventual success of the procedure. These conditions must be evaluated as part of the preoperative planning.
3. CORRECT HANDLING OF IMPLANTS IS EXTREMELY IMPORTANT. The device should not be bent, notched or scratched. All of these operations can produce defects in the surface finish and internal stress concentrations, which may become the focal point for eventual failure of the appliance.
4. METALLIC SURGICAL IMPLANT SHOULD NOT BE REUSED. Any metal implant, once used, should be discarded. Even though it appears undamaged, it may already have small defects and internal stress patterns which may lead to fatigue failure.
5. Postoperative care is extremely important. The patient must be warned that noncompliance with surgeon's instruction may delay postoperative recovery and/or increase the risk of adverse effects including implant or implant fixation failure.
6. Strict adherence by the patient to the surgeon's instructions and warnings is extremely important. Accepted practices should be followed in postoperative care. Postoperatively, the patient should be cautioned to govern his/her activity level to protect the replaced joint from unreasonable stresses. The patient should be released from the hospital with complete written instructions and warnings regarding exercises

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and therapies and any limitations on his/her activities. Periodic follow-up including x-rays are recommended for close comparison with previous post-op conditions to detect long-term evidence of changes in position, loosening, bending, or cracking of components. If evidence of one or more of these conditions is discovered, patients should be closely observed, the possibilities of further deterioration evaluated, and the benefits of early revision considered.

7. Store at room temperature.

MRI SAFETY INFORMATION:

The Screws have not been tested for safety and compatibility in the MR environment. The devices have not been tested for heating or migration in the MR environment. Scanning a patient who has this device may result in patient injury.

Risks associated with other passive implants in an MR environment have been evaluated and are known to include heating, migration, and image artifacts at or near the implant site.

HOW SUPPLIED:

Refer to product label to determine if the product is supplied **STERILE** or **NONSTERILE**.

When labeled **STERILE**, the integrity of the packaging should be checked to ensure that the sterility of the package content is not compromised. Contents are sterile unless package is damaged or opened.

DO NOT USE if the package is damaged or broken and sterility may be compromised.

Caution: Federal Law (USA) restricts this device to the sale by or on the order of a physician.

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QTY
QUANTITY

MATL
MATERIAL

MFG and
MFG/MANUFACTURED BY:
MANUFACTURED BY:

MADE IN
MADE IN

DIST
DISTRIBUTED BY

Size and **SZ**
SIZE

℞ Only
CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

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DePuy International Ltd.