0902-00-252
Rev. P

CE 2797

DePuy Synthes

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0902-00-785
Rev. G

CE
2797

DePuy Synthes
PART OF THE Johnson & Johnson FAMILY OF COMPANIES

P.F.C.® Modular Plus Tibial Wedges and P.F.C.® SIGMA® Femoral Augments

Revised 2021-02
© DePuy Synthes 2018-2021. All rights reserved.

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en **INSTRUCTIONS FOR USE**

P.F.C.® Modular Plus Tibial Wedges and P.F.C.® SIGMA® Femoral Augments

IMPORTANT INFORMATION – Please read before use

STERILE **See product label for method**

P.F.C. MODULAR PLUS TIBIAL WEDGES:

The P.F.C. Modular Plus Tibial Wedges are provided for optional use with the P.F.C. Modular Plus and Modular Plus Offset Tibial Trays to augment a tibial plateau with bone deficiencies. The wedges are attached to the P.F.C. Modular Plus or Modular Plus Offset Tibial Tray with the use of two or four Co-Cr-Mo attachment screws. The reversible wedges are provided in assorted hemispherical stepped (STEP), hemispherical angled (HEMI), and full angled (FULL) configurations to facilitate surgeon selection for addressing typical bony defects. These wedges have a recessed waffle pattern and an anterior pocket for cement fixation.

Co-Cr-Mo FEMORAL AUGMENTATION COMPONENTS DISTAL AND POSTERIOR

The Femoral Augmentation Components provide the surgeon the option of building up the distal and/or posterior surfaces of the Femoral Component in cases of bony insufficiency. The Femoral Augmentation Components utilize a Co-Cr-Mo collet and set screw that locks into a cavity on the inner surface of the Femoral Component. Applying the specified torque to the set screw expands the collet, securing the Augmentation Component in place. The Augmentation Components are available in distal or posterior configurations in various thicknesses, for use with the Cruciate Substituting or TC3 Femoral Components.

INDICATIONS AND USAGE:

The P.F.C. Modular Plus Tibial Wedges and P.F.C. SIGMA® Femoral Augments are indicated for use as a total knee replacement for patients suffering from severe pain and disability due to permanent structural damage resulting from rheumatoid arthritis, osteoarthritis, post-traumatic arthritis, collagen disorders or pseudogout. This damage may also be the result of trauma or failed prior surgical intervention.

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Choice of specific femoral and tibial components will depend on whether the posterior cruciate ligament is retained or excised, the extent and nature of anterior/posterior and medial/lateral stabilization, the femoral and/or tibial augmentation required, and the need for enhanced mechanical fixation provided by femoral and tibial stems or fluted rods.

Components are indicated for use only with Polymethylmethacrylate (PMMA) bone cement.

The P.F.C. Tibial Wedge is indicated to provide augmentation of a tibial plateau with bony defects.

INFORMATION FOR USE:

An instrumentation system is available to assure proper fit and alignment of the prosthesis. Correct fit and alignment will reduce stresses at interface surfaces to enhance implant fixation. The surgeon should refer to the appropriate surgical technique manual for details on the use of the instrument system and implantation of the prosthesis. This manual is available from your local DePuy sales representative or distributor.

CONTRAINDICATIONS:

Use is contraindicated in cases with active or recent joint sepsis, insufficient bone stock, skeletal immaturity, inadequate soft tissue coverage or support, Charcot arthropathy, significant muscular or neuromuscular compromise or psychosocial disorders that would limit rehabilitation.

WARNINGS:

Improper prosthesis selection or alignment, inadequate fixation, use where contraindicated or in patients where medical, physical, mental or occupational conditions will likely result in extreme stresses to the implant, may result in premature failure due to loosening, fracture, or wear. Postoperative care is extremely important. The patient should be instructed on the limitations of the device and should be cautioned regarding load bearing, ranges of motion, and activity levels permissible. Early motion and load bearing should be carefully controlled.

Refer to “Precautions” for important component compatibility information. Failure to adhere to the compatibility recommendations will result in poor function, loosening, wear, fracture, or premature failure. Use of other manufacturer’s components with this system is not advised.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

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DePuy's Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use. Reuse can potentially compromise device performance and patient safety.

MRI SAFETY INFORMATION:

The P.F.C. Modular Plus Tibial Wedges and P.F.C. SIGMA Femoral Augments have not been tested for safety and compatibility in the MR environment. The devices have not been tested for heating or migration in the MR environment. Scanning a patient who has this device may result in patient injury.

Risks associated with other passive implants in an MR environment have been evaluated and are known to include heating, migration, and image artifacts at or near the implant site.

PRECAUTIONS:

COMPONENT COMPATIBILITY - P.F.C. TIBIAL WEDGES:

P.F.C. Modular Plus Tibial Wedges can only be attached, with the use of screws, to the P.F.C. Modular Plus or Modular Plus Offset Tibial Tray of the corresponding size. The full P.F.C. Modular Plus Tibial Wedge must be fixed to the tibial tray with two long screws on the thick side and two short screws on the thin side of the wedge.

COMPONENT COMPATIBILITY - FEMORAL AUGMENTATION COMPONENTS:

Femoral Augmentation Components, available in either Distal or Posterior configurations, may be used only with either Cruciate Substituting or TC3 Femoral Components and have a one-to-one size relationship with these femoral components. Certain individual Posterior and Distal Augmentation Components can be combined. When a combination of individual components are utilized, the compatibilities in Table 1 must be observed.

An implant should never be reused. Any implant, once used, should be discarded. Even though it appears undamaged, it may have small defects and internal stress patterns that may lead to failure. Likewise, a new implant should be handled carefully to avoid damage that could compromise the mechanical integrity of the device and cause early failure or loosening.

The wear rate of prosthesis contact surfaces is greatly accelerated if loose fragments of bone cement become detached and act as an abrasive in the bearing surfaces. When using bone cement, care should be taken to remove all excess cement from the periphery of the implant.

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ADVERSE EFFECTS:

Fracture may occur due to improper preparation of the implant site or if excessive force is used during seating of the implant.

Transient peroneal palsy has been reported following total knee arthroplasty especially after correction or severe flexion or valgus deformities.

Patients have complained of persistent pain and stiffness following total knee arthroplasty. In addition, patellar tendon rupture, femoral-tibial subluxation or dislocation and persistent ligamentory laxity have been reported with use of total knee implants. Infection and loosening have been reported following total knee arthroplasty, as have wear and failure due to fracture of the knee prosthesis components.

Histological reactions have been reported as an apparent response to exposure to a foreign material.

Serious adverse effects may necessitate surgical intervention.

STERILITY AND HANDLING:

The P.F.C. Modular Plus Tibial Wedges and P.F.C. SIGMA Femoral Augments are provided STERILE using gamma irradiation. Remove from the package using accepted aseptic technique only after the correct size has been determined.

The operating room staff should be conservative in component selection since opened packages will not be accepted for credit.

HOW SUPPLIED:

The P.F.C. Modular Plus Tibial Wedges and P.F.C. SIGMA Femoral Augments are supplied in a range of sizes. Please refer to the current price list for product codes and sizes available.

® P.F.C. and SIGMA are trademarks of DePuy Synthes

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"X" denotes compatibility (R)

* The 12mm and 16mm Distal Augmentation Components are indicated for use in conjunction with a stem for additional implant fixation and stability. (S)

\+ The 16mm Distal Augmentation Components are indicated for use with a TC3 Femoral Component. (T)

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SIZE

Rx Only

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

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© DePuy Synthes 2015-2021. All rights reserved.
DePuy International Ltd.