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DENTIN CONDITIONER

POLYACRYLIC ACID CAVITY CONDITIONER

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RECOMMENDED INDICATIONS
DENTIN CONDITIONER is a mild polyacrylic acid solution designed to remove the dentinal smear layer and to condition dentin, thus enhancing the bond between glass ionomer / glass hybrid and the dentin / enamel.

CONTRAINDICATIONS
Avoid use of this product in patients with known allergies to polyacrylic acid.

DIRECTIONS FOR USE
1. After tooth preparation, apply DENTIN CONDITIONER to the bonding surfaces for 20 seconds using a cotton pellet or sponge.
2. Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe. DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).
3. Apply GC Fuji BOND LC, GC Fuji II LC or other glass ionomer / glass hybrid as a bonding agent, base / liner or restorative material.
4. Close the bottle immediately after use.

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PACKAGE
Bottle of 25 g (23.8 mL)

CAUTION
1. In case of contact with oral tissue or skin, flush immediately with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
4. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

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Last revised: 07/2019

MANUFACTURED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
DISTRIBUTED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU : GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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G-COAT PLUS

NANOFILLED SELF ADHESIVE LIGHT CURED PROTECTIVE COATING

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RECOMMENDED INDICATIONS
1. To seal and protect the surface of glass ionomer, direct/indirect composite and provisional restorations.
2. To seal and protect the adhesive interface between restorations and tooth structure.

CONTRAINDICATIONS
Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

DO NOT TO BE USED
In combination with desensitizers and eugenol containing materials as these may hinder G-COAT PLUS from setting or bonding properly.

DIRECTIONS FOR USE
Preparation
Preparation of newly placed restorations
1) Finish the restoration with a superfine diamond bur. Alternatively, polish the surface according to the manufacturer's instructions.
2) Spray preparation dust away with water. Dry by gently blowing with oil free air.
Note:
Surfaces to be coated should appear dry.
Preparation of existing restorations
1) Roughen the surfaces to be coated with a superfine diamond bur.
2) Spray preparation dust away with water. Dry by gently blowing with oil free air.
Note:
Surfaces to be coated should appear dry.
Preparation of marginal area
Non cut enamel should be treated with a 35-40% phosphoric acid according to the manufacturer's instructions. Rinse and dry with oil free air.
Note:
Surfaces to be coated should appear dry.

Application
1. Dispense a few drops into a disposable dispensing dish. Replace bottle cap immediately after use.
2. IMMEDIATELY apply (within 1 minute after dispensing) to the surfaces to be coated using a micro-tip applicator. Make sure that a disposable micro-tip applicator is firmly fixed on the applicator holder. DO NOT AIR BLOW.
Note:
Should the coated surface be contaminated with water, blood or saliva prior to light curing, wash and dry the surface and repeat the coating procedure.

Light curing
1. IMMEDIATELY light cure all coated surfaces with a visible light curing unit (> 500 mW/cm² : Halogen, GC G-Light, LED) for 20 seconds.
Note:
1) Place light source as close as possible to the coated surface.
2) If surface is tacky or yellowish, repeat light curing.
3) Use a protective light shield or similar protective eye wear during light curing.

STORAGE
Recommended for optimal performance, store at room temperature (4-25°C)(39.2-77.0°F).

PACKAGE
G-COAT PLUS 4 ml, (1), Disposable Dispensing Dish (20), Micro Tip Applicator (50), Applicator Holder (1)

CAUTION
1. G-COAT PLUS is flammable. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
2. G-COAT PLUS is volatile. Use in a well ventilated place. Replace cap immediately.
3. In case of contact with eyes, flush immediately with water and seek medical attention.
4. In case of contact with oral tissue or skin, immediately remove with a sponge or cotton pellet. After the restorative treatment is finished, rinse thoroughly with water.
5. If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks. To avoid contact, it is recommended to apply COCOA BUTTER to the area where rubber dam cannot be placed.
6. Avoid inhalation or ingestion of material.
7. Avoid getting material on clothing.
8. Do not mix with other products.
9. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
10. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: http://www.gceurope.com
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

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CONTRAINDICATIONS
1. Direct pulp capping.
2. Avoid use on patient with known allergies to methacrylate monomer, methacrylate polymer, or acetone.

NOT TO BE USED
1. In combination with chemically cured composite resin.
2. In combination with eugenol containing materials as eugenol may hinder G-BOND from setting or bonding properly.
3. In combination with dual cured luting and core build up composites in case these materials will not be light cured.

DIRECTIONS FOR USE
1. CAVITY PREPARATION
Prepare tooth using standard techniques. Use rubber dam for protection. Dry the prepared tooth surfaces by gently blowing with an air syringe.
Note:
For direct pulp capping, use calcium hydroxide.
2. APPLICATION
a) Prior to dispensing, shake the bottle of G-BOND thoroughly (Fig. 1).
b) Dispense a few drops into a clean dispensing dish (Fig. 2).
c) Dispense bottle cap immediately after use (Fig. 3).
d) IMMEDIATELY apply to the prepared enamel and dentin surfaces using the DISPOSABLE APPLICATOR II (Fig. 4).
e) Leave undisturbed for 5-10 seconds after the end of application (Fig.5).
f) Then, dry thoroughly for 5 seconds with oil free air under MAXIMUM air pressure. Use vacuum suction to prevent splatter of the adhesive (Fig. 6). The final result should be a thin, rough, adhesive film with the appearance of frosted glass and which does not visibly move under further air pressure.
Note:
1) When removing from refrigerator after long-term storage, leave the bottle of G-BOND to stand at room temperature for several minutes prior to use.
2) Non cut enamel should be treated with a 35-40% phosphoric acid gel, rinsed with water and gently driedbefore using G-BOND.
3) Apply G-BOND immediately as the material contains a volatile solvent.
4) Remove excess material remaining on the tooth (other than the surfaces to be bonded) using a sponge or cotton pellet as residue material is hard to remove after light curing.

Fig.2 / Фиг. 2 / obr. 2 / slika 2 / 2. ábra / ryc. 2 / рис. 2

Fig.3 / Фиг. 3 / obr. 3 / slika 3 / 3. ábra / ryc. 3 / рис. 3

Fig.4 / Фиг. 4 / obr. 4 / slika 4 / 4. ábra / ryc. 4 / рис. 4

Fig.5 / Фиг. 5 / obr. 5 / slika 5 / 5. ábra / ryc. 5 / рис. 5

In case where the light guide tip is more than 10mm away from the surface to be irradiated, light cure for the following time.

Halogen / LED (700 mW/cm²)	10 seconds
High power LED (more than 1200 mW/cm²)	5 seconds

In case where the light guide tip is more than 10mm away from the surface to be irradiated, light cure for the following time.

Halogen / LED (700 mW/cm²)	20 seconds
High power LED (more than 1200 mW/cm²)	10 seconds

Note:
1) Light cure completely for an effective bonding strength.
2) Use a protective light shield or similar protective eye wear during light curing.

Fig.6 / Фиг. 6 / obr. 6 / slika 6 / 6. ábra / ryc. 6 / рис. 6

4A. PLACEMENT OF LIGHT CURED COMPOSITES AND COMPOMERS
After having light cured the adhesive, follow manufacturer's instructions for use in order to place, contour and light cure the composite resin material of choice.
4B. PLACEMENT OF DUAL CURED COMPOSITES
After having light cured the adhesive, make sure to light cure the dual cured composite material separately. Self curing only will result in improper adhesion.
5. FINISHING
Adjust restoration and polish using standard techniques.

Fig.7 / Фиг. 7 / obr. 7 / slika 7 / 7. ábra / ryc. 7 / рис. 7

STORAGE
If not in use for a prolonged period of time, store in refrigerator, otherwise store at room temperature (1-28°C) (33.8-82.4°F).

PACKAGE
1. Kit:
5 mL liquid (1), dispensing dish (No.1) (1), DISPOSABLE APPLICATOR II (50)
2. Refill
5 mL liquid (1)
Note:
All packaging types may not be available in all countries.

CAUTION
1. G-BOND is flammable. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
2. G-BOND is volatile. Use in a well ventilated place.
3. In case of contact with eyes, flush immediately with water and seek medical attention.
4. In case of contact with oral tissue or skin, immediately remove with a sponge or cotton pellet. After the restorative treatment is finished, rinse thoroughly with water.
5. If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks. To avoid contact, it is recommended to apply COCOA BUTTER to the area where rubber dam cannot cover.
6. Avoid inhalation or digestion of material.
7. When put on the table or floor, immediately wipe off with a dry cloth.
8. Do not mix with other products.
9. Dispose of all waste according to local regulations.
10. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
11. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
12. Do not place under operating light.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at
http://www.gceurope.com
They can also be obtained from your supplier.

CLEANING AND DISINFECTING
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

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Last revised: 06/2019

CE 2797

MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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Prior to use, carefully read the instructions for use. EN

G-ænial™ Bond

ONE COMPONENT SELF-ETCHING LIGHT-CURED ADHESIVE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS
1. Bonding of light cured composites and acid modified composites (compomers) to tooth structure.
2. Bonding of dual cured luting and core build up composites to tooth structure as long as these materials are light cured.

CONTRAINDICATIONS
1. Pulp capping.
2. Avoid use of this product in patients with known allergies to methacrylate monomer, methacrylate polymer, or acetone.

NOT TO BE USED
1. In combination with chemically cured composite resin.
2. In combination with eugenol containing materials as eugenol may hinder G-ænial Bond from setting or bonding properly.
3. In combination with desensitizers as desensitizers may hinder G-ænial Bond from setting or bonding properly.
4. In combination with dual cured luting and core build up composites unless both G-ænial Bond and these materials are individually light cured.

DIRECTIONS FOR USE
1. CAVITY PREPARATION
 Prepare tooth using standard techniques. Make sure that the temporary material is thoroughly removed from the prepared cavity. Use rubber dam for protection. Dry the prepared tooth surfaces by gently blowing with an air syringe.
 Note:
 For pulp capping, use calcium hydroxide.
2. TECHNIQUE SELECTION
 Select out of the following 2 technique options
 a) Self-etch technique – apply G-ænial Bond to bur cut enamel and dentin, without a separate etching step.
 b) Selective etching of enamel – prior to applying G-ænial Bond to enamel and dentin, etch (bur cut) enamel with 35-40% phosphoric acid gel for 10 seconds, rinse for 5 seconds and gently dry.
 Note:
 Non bur cut enamel should always be treated with a 35-40% phosphoric acid gel for 10 seconds, rinsed with water for 5 seconds and gently dried.
3. APPLICATION
 a) Prior to dispensing, shake the bottle of G-ænial Bond thoroughly (Fig.1). Dispense a few drops into a clean dispensing dish (Fig. 2). Replace bottle cap immediately after use (Fig. 3).
 b) IMMEDIATELY apply to the prepared enamel and dentin surfaces using the disposable applicator (Fig. 4).
 c) Leave undisturbed for 10 seconds after the end of application (Fig.5).
 d) Then, dry thoroughly for 5 seconds with oil free air under MAXIMUM air pressure. Use sufficient strength to prevent splatter of the adhesive (Fig. 6). The final result should be a thin adhesive film with the appearance of frosted glass and which does not visibly move under further air pressure.
 Note:
 1) When removing from refrigerator after long-term storage, leave G-ænial Bond to stand at room temperature for several minutes prior to use.
 2) Apply G-ænial Bond immediately as the material contains a volatile solvent.
 3) Remove excess material remaining on the tooth (other than the surfaces to be bonded) using a sponge or cotton pellet as residue material is hard to remove after light curing.
 4) Should the applied material be contaminated with water, blood or saliva prior to light curing, wash and dry the tooth and repeat the procedure by re-applying material.
4. LIGHT CURING
 Light cure using a visible light curing unit (Fig. 7).

Irradiation time	
Halogen / LED (700 mW/cm²)	: 10 seconds
High power LED (more than 1200 mW/cm²)	: 5 seconds

 In case where the light guide tip is more than 10mm away from the surface to be irradiated, light cure for the following time.

Halogen / LED (700 mW/cm²)	: 20 seconds
High power LED (more than 1200 mW/cm²)	: 10 seconds

 Note:
 1) Light cure completely for an effective bonding strength. Lower light intensity may cause insufficient adhesion.
 2) Use a protective light shield or similar protective eye wear during light curing.

5A. PLACEMENT OF LIGHT CURED COMPOSITES AND COMPOMERS
After having light cured the adhesive, follow manufacturer's instructions for use in order to place, contour and light cure the composite resin material of choice.

5B. PLACEMENT OF DUAL CURED COMPOSITES
After having light cured the adhesive, make sure to light cure the dual cured composite material separately. Self curing only will result in insufficient bonding strength.

6. FINISHING
 Adjust restoration, finish and polish using standard techniques.

STORAGE
Store at room temperature (1-28°C) (33.8-82.4°F).
If not in use for a prolonged period of time, store in refrigerator.

PACKAGE
1. G-ænial Bond Bottle Kit: 5 mL liquid (1). DISPOSABLE DISPENSING DISH (20). DISPOSABLE APPLICATOR II (50)
2. G-ænial Bond Bottle Refill: 5 mL liquid (1)
3. G-ænial Bond 3-bottle Pack: 5 mL liquid (3)

CAUTION
1. G-ænial Bond is flammable. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
2. G-ænial Bond is volatile. Use in a well ventilated place. Advise patients to breath through nose.
3. In case of contact with eyes, flush immediately with water and seek medical attention.
4. In case of contact with oral tissue or skin, immediately remove with a sponge or cotton pellet. After the restorative treatment is finished, rinse thoroughly with water.
5. If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks. If contact, it is recommended to apply cocoa butter to the area where rubber dam cannot cover.
6. Avoid inhalation or digestion of material.
7. When spill on the table or floor, immediately wipe off with a dry cloth. Do not mix with other products.
8. Dispose of all waste according to local regulations.
10. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
11. In rare cases this may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
http://www.gceurope.com
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires high level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by user of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

Last revised: 06/2019

CE 2797

MANUFACTURED by:
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED by:
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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RECOMMENDED INDICATIONS
GC Fuji PLUS CONDITIONER is a tooth conditioning agent designed for use with GC Fuji PLUS glass ionomer luting cement. It helps to enhance the bond between GC Fuji PLUS and tooth structure.

DIRECTIONS FOR USE
1. After tooth preparation, apply GC Fuji PLUS CONDITIONER to the bonding surface for 20 seconds using a cotton pellet or sponge.
2. Rinse thoroughly with water and dry. DO NOT DESICCATE.
3. Proceed according to GC Fuji PLUS directions.
4. Close the bottle immediately after use.

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PACKAGES
Bottle of 7 g (6.5 mL)

CAUTION
1. In case of contact with oral tissue or skin, flush immediately with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. Not recommended for use with products that contain Epinephrine, which is known to react with iron containing products and may stain oral tissue.
4. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
5. In rare cases the product may cause sensitivity to some persons. If such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS.
Always familiarize yourself with the safety data sheets available at:
http://www.gceurope.com
or for The Americas
http://www.gcamerica.com
They can also be obtained from your supplier.

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OPTIGLAZE™ color

LIGHT-CURED CHARACTERIZATION MATERIAL

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RECOMMENDED INDICATIONS

1. Characterization of direct & indirect composite restorations, acrylic denture base and artificial acrylic teeth.
2. For obtaining surface smoothness and wear resistance of restorations made of composite resin, acrylic denture base and artificial acrylic teeth.
 NOTE:
 Not intended for characterization of natural teeth.

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to acrylate / methacrylate monomers or acrylate / methacrylate polymers.

DIRECTIONS FOR USE

Surface characterization

1) After contouring the resin surface, finish with a carbide bur or coarse silicone point. The oxygen inhibition layer must be removed as it inhibits the polymerization of OPTIGLAZE color.
 Laboratory use:
 Sandblast with 25-50μm alumina (0.15MPa/1.5bar.), clean and dry.
 NOTE:
 1) If surface still contains oil residue, the surface can be cleaned further with ethanol (alcohol.) The use of ethanol (alcohol) on acrylic resin may cause micro cracks. A mild detergent can be used to clean the acrylic resin.
 2) When applying OPTIGLAZE color on hybrid ceramic blocks such as CERASMART, treat the surface using CERAMIC PRIMER II.
2) Shake the bottle of OPTIGLAZE color well. Dispense a few drops in a dispensing dish. Dilute by adding color to the desired shade. Apply a thin layer to the resin surface using the brush provided. Do not air blow.
3) Light cure with a suitable light curing device. (See table 1) Place the light tip as close as possible to the surface.
 Laboratory use:
 Set a 2cm platform in the curing device and place the prosthesis on it with the coated surface facing the light. Reverse the prosthesis and repeat the light curing to fully cure the undercut or area in shadow.

Internal characterization

1) Build up the dentin layer with composite and light cure. Follow the instructions of your composite system.
 Note:
 Do not remove the oxygen inhibition layer. Apply OPTIGLAZE color on the oxygen inhibition layer. If the oxygen inhibition layer is removed, treat the surface with CERAMIC PRIMER II (blow dry) and COMPOSITE PRIMER (light cure) before applying OPTIGLAZE color.
2) Shake the bottle of OPTIGLAZE color well. Dispense a few drops in a disposable dispensing dish. Apply a thin layer of OPTIGLAZE color of choice to the dentin surface with an instrument or brush (e.g. the convex and concave of dentin manekin in the picture). Mix multiple colors to obtain the desired shade if needed. Do not air blow.
3) Light cure with a suitable light curing device. (See table 1)
4) Continue building up with enamel shade on the oxygen inhibition layer. Follow the instructions of your composite system.
5) Finish and polish using standard technique or apply OPTIGLAZE color CLEAR or CLEAR HV on the entire surface.

Repair of restoration (e.g. when the gloss is lost)

1) Roughen the surface of restoration using a carbide bur and paper cone
2) Clean the surface with a steam cleaner or ultrasonic cleaner and then dry with oil free air.
 NOTE:
 If surface is contaminated by oil residue, the surface can be cleaned further with ethanol (alcohol). The use of ethanol (alcohol) on acrylic resin may cause micro cracks. A mild detergent can be used to clean the acrylic resin.
3) Apply OPTIGLAZE color CLEAR or CLEAR HV shade to the surface. Do not air blow.
4) Light cure with a suitable light curing device. (See table 1)

Table 1: Light Curing Time

GC LABOLIGHT DUO		90 sec.
GC LABOLIGHT LV-III, II		5 min.
GC STEPLIGHT SL-1	Distance from light source <3cm	10 sec.
	Distance from light source >3cm	20 sec.
Halogen Light		40 sec.
Plasma Arc		8 sec.
LED Light (with light source having wavelength lower than 430nm)		40 sec.

Table 2: Depth of Cure

Orange, Yellow, Redbrown, Red	0.07 mm
A Plus, B Plus, C Plus, White, Ivory White, Pink Orange, Pink, Olive, Lavender, Grey, Blue	0.15 mm
Clear, Clear HV	0.5 mm

STORAGE

Recommended for optimal performance, store at room temperature of (4-25°C) (39.2-77.0°F).

PACKAGE:

Complete set:
A Plus 2.6mL, B Plus 2.6mL, C Plus 2.6mL, White 2.6mL, Ivory White 2.6mL, Yellow 2.6mL, Orange 2.6mL, Pink Orange 2.6mL, Pink 2.6mL, Red Brown 2.6mL, Olive 2.6mL, Lavender 2.6mL, Grey 2.6mL, Blue 2.6mL, Red 2.6mL, Clear 5mL, Clear HV 5mL, Disposable Dispensing Dish, Flat Brush, Round Brush

CAUTION

1. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
2. In particular avoid contact of uncured material with skin.
3. Avoid getting material on clothing.
4. OPTIGLAZE color is highly flammable liquid and vapor. DO not use near open flame. Keep away from source of ignition.
5. OPTIGLAZE color is volatile and may cause respiratory irritation. Use in a well ventilated space. Replace cap immediately after dispensing.
6. It is not recommended to use OPTIGLAZE color on high pressure bearing area.
7. In case of contact with eyes, flush immediately with water and seek medical attention.
8. In case of contact with skin, immediately remove with a sponge or cotton pellet. Flush with water.
9. The oxygen inhibited layer of the resin surface to be coated should be removed before applying OPTIGLAZE color as it may hinder the proper curing of OPTIGLAZE color.
10. Dispense a few drops of OPTIGLAZE color into a disposable dispensing dish immediately before applying as the material contains a volatile solvent.
11. Store locked up
12. Dispose all waste according to local regulation.
13. In rare cases the product may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

CLEANING AND DISINFECTING

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

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Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

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OPTIGLAZE™

NANOFILLED LIGHT CURED PROTECTIVE COATING FOR COMPOSITE AND ACRYLIC INDIRECT RESTORATIONS

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RECOMMENDED INDICATIONS
1. Coating of indirect composite surfaces (fissures and interproximal areas)
2. Coating of denture base surfaces
3. Coating of provisional materials

CONTRAINDICATIONS
Avoid use of this product in patients with known allergies to acrylate or methacrylate monomer.

PRECAUTIONS
1. It is not recommended to use OPTIGLAZE on high stress bearing areas (occlusal surfaces).
2. The use of OPTIGLAZE does not replace the mechanical polishing of aesthetic indirect composite restorations such as GRADIA.

DIRECTIONS FOR USE
Preparation
After contouring the resin surface, finish with carbide bur or coarse silicone point. Clean the surface with a steam cleaner or ultrasonic cleaner and then dry. If surface still contains oil residue, the surface can be further cleaned with ethanol (alcohol).
NOTE:
1) OPTIGLAZE does not adhere to a polished surface. Make sure surfaces to be coated are rough.
2) The air-inhibited layer of the resin surface to be coated should be removed before applying OPTIGLAZE as it may inhibit the proper curing of OPTIGLAZE.
3) The use of ethanol (alcohol) on acrylic resin may cause micro-cracks. A mild detergent can be used to clean the acrylic resin.

Application
1. Dispense a few drops into a dispensing dish. Replace bottle cap immediately after dispensing.
2. IMMEDIATELY apply (within 1 minute after dispensing) to the surface with a brush. DO NOT AIR BLOW.

Light Curing
Light cure with a suitable light curing device.

LABOLIGHT DUO		90 sec.
LABOLIGHT LV-III, II		3 min.
STEPLIGHT SL-I	Distance from light source <3 cm	10 sec.
	Distance from light source >3 cm	20 sec.
Halogen Light		40 sec.
LED Light (with light source having wavelength lower than 430 nm)		40 sec.

NOTE:
1) If surface is tacky, repeat light curing.
2) If necessary, rotate the prosthesis a half-turn or turn-over and light cure again to ensure all surfaces are completely cured.
3) Additional build-up cannot be performed on surfaces coated with OPTIGLAZE. Remove the coated layer with a bur and then build -up.
4) OPTIGLAZE cannot be cured by a conventional LED light curing device with blue LED only.
5) When using other manufacturers' light curing devices, use final curing mode. Refer to the respective manufacturers' instructions for use.

STORAGE
Recommended for optimal performance, store at temperature of 4-25°C (39.2-77.0°F).

PACKAGE
OPTIGLAZE 15 mL (1)

CAUTION
1. OPTIGLAZE is flammable. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
2. OPTIGLAZE is volatile. Use in a well ventilated area. Replace cap immediately after dispensing.
3. In case of contact with eyes, flush immediately with water and seek medical attention.
4. In case of contact with skin, immediately remove with a sponge or cotton pellet. Flush with water.
5. Wear protective gloves. In particular, avoid contact of uncured material with skin.
6. Avoid getting material on clothing.
7. Avoid inhalation or ingestion of material.
8. Do not mix with other products.
9. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
10. In rare cases the product may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

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Undesired effects- Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant competent authority, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

Last revised: 06/2019

MANUFACTURED by:
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan
DISTRIBUTED by:
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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Windows® Mac® Chromebook™ www.gcamerica.com/products/manuals

EN

Prior to use, carefully read the instructions for use.
For use only by a dental professional in the recommended indications.

INDICATIONS FOR USE

1) Class I restorations
2) Stress bearing Class II restorations (see separate technique guide for suggested cavity preparations)
3) Non-stress bearing Class II restorations
4) Intermediate restorative
5) Class V and root surface restorations
6) Core build up

EQUIA Forte Coat is used to seal, strengthen and protect the surface of EQUIA Forte HT Fil restorations.

CONTRAINDICATIONS

1) Pulp capping.
2) Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer.
3) EQUIA Forte Coat is not suitable for patients with a history of hypersensitivity to acrylate/methacrylate monomers.

NOT TO BE USED

EQUIA Forte Coat should not be used in combination with desensitizers and eugenol containing materials as these may hinder EQUIA Forte Coat from setting or bonding properly.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	0.40 / 0.13
Mixing Time (sec.)	10"
Time of initial extrusion after mixing	10 seconds maximum
Working Time (23°C/73 °F) (min., sec.)	1'30"
Net Setting Time (37 °C/99 °F) (min., sec.)	2'30"
Finishing Time (37 °C/99 °F, from start of mixing) (min.,sec.)	2'30"
VLC curing time of EQUIA Forte Coat	20"

1. CAVITY PREPARATION

a) Prepare tooth using standard techniques. Extensive mechanical retention is unnecessary. For pulp capping use calcium hydroxide. In class V or II situations, we advise the use of an anatomical matrix system. Metal bands should be lubricated with petroleum jelly.
Note:
Optionally, CAVITY CONDITIONER or DENTIN CONDITIONER can be applied to the bonding surfaces.

b) Rinse thoroughly with water. Blot away excess water with a cotton pellet or dry by gently blowing with an air syringe.
DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).

2. CAPSULE ACTIVATION AND MIXING

a) Before activation, shake the capsule or tap its side on a hard surface to loosen the powder.

b) To activate the capsule, push the plunger until it is fully depressed with the main body and hold it down for 2 seconds.
Note: 1) Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.
2) The capsule should be activated just before mixing and used immediately.

c) Immediately set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000 RPM).
Note: When using a rotational mixing device, it is recommended not to use the centrifuge option.

Riva Luxo MAX Polyacrylic
Version 1.0 | 2021-10-01 | DFU-0000-EN-01

3. RESTORATIVE TECHNIQUE

a) Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.
b) Make two clicks to prime the capsule then syringe. The working time is 1 minute 30 seconds from start of mixing at 23°C (73°F). Higher temperatures will shorten working time.
c) Within 10 seconds maximum after mixing, start to extrude the mixture directly into the preparation.
d) Form the preliminary contour and cover with a matrix if required.
 Note:
 it is recommended to wait until the material is not glossy anymore (approximately 10 seconds) to start contouring.
e) During the first 2 minutes 30 seconds from start of mix extra care should be taken to avoid moisture contamination or drying-out. In case this cannot be guaranteed, immediately apply EQUIA Forte Coat and light cure.
 Note:
 1) To adjust the direction of the nozzle, hold the applier with the capsule towards you and turn the capsule body.
 2) To remove the used capsule, push the applier release button. Twist the capsule and pull upwards.

4. FINISHING

Finish under water spray using superfine diamond burs after 2 minutes 30 seconds from start of mixing.

Note:

When applying EQUIA Forte Coat on the existing EQUIA Forte HT Fil restorations, roughen the surface to be coated with superfine diamond burs.

5. COATING

a) Spray preparation dust away with water. Dry by gently blowing with oil free air. Surfaces to be coated should be dry. Do not desiccate.

b) Bottle type: Prior to dispensing, shake the bottle of EQUIA Forte COAT thoroughly. To open the flip-top cap, pinch both cap release buttons on the neck of the bottle at the same time, as illustrated in Figure 1. Dispense a few drops into a clean dispensing dish. Close the cap immediately.

Note:

1) If the cap detaches from the bottle, it should be reattached by reseating the cap HORIZONTALLY on the plastic hinge pegs, and firmly pushing it back into place. The metal spring should be INSIDE the cap, as shown in Figure 2.
2) Wipe any excess liquid from the nozzle before closing the cap.

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

Fig. 1

Fig. 2

Unitdose type: Before opening unit dose capsule, ensure top of the capsule is facing up. Unscrew cap counterclockwise.
Note:
If the cap of the unit dose is too tight to open, use another unit dose. Do not try to open the unit dose capsule with excessive force as this may cause the liquid inside to spill and cause accidental contact with skin and eyes.

c) IMMEDIATELY apply (within 1 minute after dispensing) to the surfaces to be coated using the disposable micro-tip applicator. Make sure that a disposable micro-tip applicator is firmly fixed on the applicator holder. Use floss to apply to approximal surface. DO NOT AIR BLOW.

Note:

1) Non bur cut enamel should be pre-treated with a 35-40% phosphoric acid according to the manufacturer's instructions. Rinse and dry with oil free air.
2) Should the coated surface be contaminated with water, blood or saliva prior to light curing, wash and dry the surface and repeat the coating procedure.
3) Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the application procedure.
4) If the cap of the unit dose is too tight to open, do not try to open it with force to avoid accidental contact of the liquid with skin and eyes. Use another unit dose.
5) Unit dose is disposable. Do not recap or reuse to avoid cross contamination between patients.

6. LIGHT CURING

IMMEDIATELY light cure all coated surfaces with a visible light curing unit (> 500mW/cm2: Halogen/LED) for 20 seconds.

Note:

1) Place light source as close as possible to the coated surface.
2) If surface is tacky or yellowish, repeat light curing.

Ritm Lipase Mill & Pro Lipase manual
Ritm Lipase Mill & Pro Lipase manual

3) Use a protective light shield or similar protective eye wear during light curing.
4) When curing, EQUIA Forte Coat may produce a brief exothermic reaction which can be minimized by applying in a thin layer.

Do not cure on soft tissue.

STORAGE

Recommended for optimal performance, store the EQUIA Forte HT Fil capsules in the original aluminium foil in a cool and dark place (4-25°C) (39-77°F).
If not in use for a prolonged period of time, store EQUIA Forte Coat in refrigerator, otherwise store at room temperature
(4-25°C) (39-77°F).

SHADE

Available in Refills: A1, A2, A3, A3.5, B1, B2, B3, C4
Shade numbers according to Vita® shade guide.
Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGE

EQUIA Forte HT Fil

1) Single shade package: 20 capsules (A2, A3, B2)
2) Single shade package: 50 capsules (A1, A2, A3, A3.5, B1, B2, B3, C4)
3) Assortment package : 50 capsules (10 each of A2, A3, A3.5, B1, B3)
4) Single shade package: 100 capsules (A2, A3)
5) Double shade package: 100 capsules (A2-A3, A3-B2)
6) Single shade package: 200 capsules (A2, A3, B2)

Average contents per capsule: 0.40g powder and 0.13g (0.10mL) liquid
Minimum net volume of mixed cement per capsule: 0.14mL

EQUIA Forte Coat

1) Bottle 4mL (1), MICRO TIP APPLICATOR (25)
2) Unitdose (0.1mL) (50), MICRO TIP APPLICATOR (25)

Option: GC CAPSULE APPLIER (1 piece)

CAUTION

1. In case of contact with oral tissue or skin, immediately remove with a sponge or cotton pellet. After the restorative treatment is finished, rinse thoroughly with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. EQUIA Forte Coat is a highly flammable liquid and vapor. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
4. EQUIA Forte Coat is volatile. Use in well ventilated place. Replace cap immediately.
5. If the tissue contacted by EQUIA Forte Coat turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks. To avoid contact, it is recommended to apply COCOA BUTTER to the area where rubber dam cannot cover.
6) Avoid inhalation or ingestion of material.
7) Avoid getting material on clothing
8) Do not mix with other products.
9) Instruct the patient not to apply pressure for 1 hour.
10) Use EQUIA Forte HT Fil with EQUIA Forte Coat.
11) Dispose of all waste according to local regulations.

12) In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
13) Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: http://www.gceurope.com; or for The Americas: http://www.gcamerica.com
They can also be obtained from your supplier.

Undesired effects - Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

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RECOMMENDED INDICATIONS

1. Class I restorations
2. Non-stress bearing Class II restorations
3. Stress bearing Class II restorations when isthmus is less than half of intercuspal distance
4. Intermediate restorative
5. Class V and root surface restorations
6. Core build up

EQUIA Coat is used to seal and protect the surface of EQUIA Fil restorations.

CONTRAINDICATIONS

1. Pulp capping.
2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

NOT TO BE USED

EQUIA Coat should not be used in combination with desensitizers and eugenol containing materials as these may hinder EQUIA Coat from setting or bonding properly.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	0.40 / 0.12
Mixing Time (sec.)	10"
Time of initial extrusion after mixing	10 seconds maximum
Working Time (min., sec.)	1' 15"
Net Setting Time (min., sec.)	2' 00"
Final Finishing Commencing Time (37°C, after starting mix) (min.,sec.)	2' 30"
VLC curing time of EQUIA Coat	20"

Test conditions : Temperature (23 +/-1°C), Relative humidity (50 +/-5%)

1. CAVITY PREPARATION

a) Prepare tooth using standard techniques. Extensive mechanical retention is unnecessary. For pulp capping use calcium hydroxide. In class V or II situations, make sure to use a right matrix system. Metal bands should be lubricated with petroleum jelly.
b) Apply GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) to the bonding surfaces using a cotton pellet or sponge.
c) Rinse thoroughly with water. Blot away excess water with a cotton pellet or dry by gently blowing with an air syringe.
DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).

2. CAPSULE ACTIVATION AND MIXING

a) Before activation, shake the capsule or tap its side on a hard surface to loosen the powder.
b) To activate the capsule, push the plunger until it is fully depressed with the main body and hold it down for 2 seconds.
Note: 1) Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.
2) The capsule should be activated just before mixing and used immediately.
c) Immediately set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4,000 RPM).
Note: When using a rotational mixing device, it is recommended not to use the centrifuge option.

Reinforced Metal-Free Glass-ionomer
Luting and Restorative System

3. RESTORATIVE TECHNIQUE

a) Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.
b Make two clicks to prime the capsule then syringe. The working time is 1 minute 15 seconds from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
c) Within 10 seconds maximum after mixing, start to extrude the mixture directly into the preparation.
d) Form the preliminary contour, and cover with a matrix if required.
e) During the first 2 minutes 30 seconds from start of mix extra care should be taken to avoid moisture contamination or drying-out. In case this cannot be guaranteed, immediately apply EQUIA Coat and light cure.
Note : 1) To adjust the direction of the nozzle, hold the applier with the capsule towards you and turn the capsule body.
2) To remove the used capsule, push the applier release button. Twist the capsule and pull upwards.

4. FINISHING

Finish under water spray using superfine diamond burs after 2 minutes 30 seconds from start of mixing.
Note: When applying EQUIA Coat on the existing EQUIA restorations, roughen the surface to be coated with superfine diamond burs.

5. COATING

a) Spray preparation dust away with water. Dry by gently blowing with oil free air. Surfaces to be coated should be dry. Do not desiccate.
b) Dispense a few drops of EQUIA Coat into a dispensing dish. Replace bottle cap immediately after use.
c) IMMEDIATELY apply (within 1 minute after dispensing) to the surfaces to be coated using the disposable applicator. Use floss to apply to approximal surface. DO NOT AIR BLOW.

Note : 1) Non bur cut enamel should be treated with a 35-40% phosphoric acid according to the manufacturer's instructions. Rinse and dry with oil free air.
2) Should the coated surface be contaminated with water, blood or saliva prior to light curing, wash and dry the surface and repeat the coating procedure.

6. LIGHT CURING

IMMEDIATELY light cure all coated surfaces with a visible light curing unit (> 500mW/cm^2 : Halogen, GC G-Light, LED) for 20 seconds.

Note: 1) Place light source as close as possible to the coated surface.
2) If surface is tacky or yellowish, repeat light curing.

3) Use a protective light shield or similar protective eye wear during light curing.

STORAGE

Recommended for optimal performance, store the EQUIA Fil capsules in the original aluminium foil in a cool and dark place (4-25°C) (39.2-77.0°F).
If not in use for a prolonged period of time, store EQUIA Coat in refrigerator, otherwise store at room temperature (4-25°C)(39.2-77.0°F).
Shelf life: EQUIA Fil - 2 years from date of manufacture
EQUIA Coat - 3 years from date of manufacture

SHADES

Available in Refills: A1, A2, A3, A3.5, B1, B2, B3, C4, SW (Standard White)*
Shade numbers according to Vita® shade guide.
Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.
*Only available in limited countries and packages.

PACKAGES

1. Single shade package : 50 capsules.
2. Assortment package : 50 capsules (10 each of A2, A3, A3.5, B1, B3).
 Average contents per capsule : 0.40g powder and 0.12g (0.10mL) liquid.
 Minimum net volume of mixed cement per capsule : 0.14mL.
3. EQUIA Coat 4mL (1) / 6mL* (1), disposable applicator (25).

Option : GC CAPSULE APPLIER (1 piece).
*Only available in limited countries and packages.

CAUTION

1. In case of contact with oral tissue or skin, immediately remove with a sponge or cotton pellet. After the restorative treatment is finished, rinse thoroughly with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.

Reproduction prohibited without permission

3. EQUIA Coat is flammable. Do not use near naked flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.
4. EQUIA Coat is volatile. Use in a well ventilated place. Replace cap immediately.
5. If the tissue contacted by EQUIA Coat turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks. To avoid contact, it is recommended to apply cocoa butter to the area where rubber dam cannot cover.
6. Do not use EQUIA Coat in combination with desensitizers and eugenol containing materials as these may hinder EQUIA Coat from setting or bonding properly.
7. Avoid inhalation or ingestion of material.
8. Avoid getting material on clothing
9. Do not mix with other products.
10. Instruct the patient not to apply pressure for 1 hour.
11. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Manufactured by GC Europe N.V. Interleuvenlaan 33, B-3001 Leuven, Belgium

Some products referenced in the present IFU may be classified as hazardous according to GHS.
Always familiarize yourself with the safety data sheets available at:
http://www.gceurope.com
They can also be obtained from your supplier.

Undesired effects - Reporting:
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

Last revised : 08/2019

Reproduction and/or distribution of this document and communication of its contents to others without express authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent or the registration of a utility model or design.

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