myra™ BMS — Instructions for use - English

Balloon Expandable Peripheral Stent System

1. Device Description:

The Myra™ BMS - Balloon Expandable Peripheral Stent System comprises following components -

- A balloon-expandable L605 Cobalt Chromium Peripheral Stent.
- An over-the-wire stent delivery PTA Balloon Dilatation Catheter.

The stent is pre mounted on Balloon Catheter & placed between two platinum-iridium radio opaque markers bands.

1.1 Device Components Description:

Table - 1

1.1.1	Stent	Electropolished L605 Cobalt Chromium alloy, laser-cut from seamless tubing in a hybrid design pattern.
1.1.2	Stent delivery balloon catheter system	Semi-compliant Polyamide balloon, nominally 1.50 mm longer than the stent length. Mounted stent length & location is defined by two platinum-iridium swaged radiopaque markers.
1.1.3	Delivery system usable length	0.035" OTW PTA Balloon Dilatation Catheter available with usable catheter length viz. 80 cm and 135 cm.
1.1.4	Guide wire lumen	Guide wire lumen covers entire length of the delivery system from distal tip to the Y Luer Hub.
1.1.5	Shaft outer profile	5F
1.1.6	Stent dilatation/Balloon inflation pressure	Refer Box Label
1.1.7	Recommended Minimum Sheath / Introducer	6/7F***
1.1.8	Guide wire compatibility	OTW 0.035" (0.89 mm)

***See individual manufacturer specifications for (F) equivalent on box label and pouch label

2. How Supplied:

Sterile: This device is sterilized with Ethylene Oxide (EO) gas and is non-pyrogenic. It is intended for single use only. Do not resterilize. Do not use the device if the package is opened or damaged.

Contents: One (1) Myra™ BMS - Balloon Expandable Peripheral Stent System housed in a protective tray, one (1) Instructions for Use, two (2) Stent Implant Card.

Storage: Store between 08-30°C (46.4-86°F) temperature in a dry, dark, cool place.

3. Indications:

The Myra™ BMS - Balloon Expandable Peripheral Stent System is intended for use in the treatment of atherosclerotic diseases of peripheral arteries below the aortic arch, with reference vessel diameter of 5.00 mm to 10.00 mm for Iliac / Sub-Clavian and for other protected peripheral arteries in patients eligible for Percutaneous Transluminal Angioplasty (PTA) and Stenting procedures.

4. Contraindications:

- In general, balloon expandable stents are contraindicated for the implantation in unprotected vascular districts, exposed to a potential risk of external or other mechanical impacts which may result in stent crush.
- Highly calcified lesions, resistant to PTA.

5. Warnings:

- Do not use if the inner package is opened or damaged.
- Myra™ BMS - Balloon Expandable Peripheral Stent System intended for single use only. Do not resterilise or reuse the device as this can potentially result in compromised device performance and increased risk of inappropriate resterilisation and cross contamination.
- Inflation at a high rate may damage the balloon.
- Persons with allergic reactions to a Cobalt-Chromium alloy, whose main constituents are Cobalt, Chromium, Tungsten, and Nickel may suffer an allergic response to this implant.
- Judicious patient selection is necessary during use of this device since it carries the associated risks of subacute thrombosis, vascular complications and/or bleeding events.

5.1. Use in special population:

Safety and effectiveness of Myra™ BMS in patients below 18 years of age has not been established.

6. Precautions:

6.1. General Precautions:

- Only physicians who have received adequate training should perform implantation of the stent.
- Subsequent blockage may require repeat dilatation of the arterial segment containing the stent. The long term outcome following repeat dilatation of the endothelialized stents is not well characterized.

6.2. Stent Handling Precautions:

- Use the device before the "Use by" date as specified on the product label.
- Do not remove the stent from the delivery system as removal may damage the stent and / or lead to stent embolization. The Myra™ BMS - Balloon Expandable Peripheral Stent System is intended to perform as a system.
- The stent should not be removed for use in conjunction with other dilatation catheters.
- Special care must be taken not to handle or in any way disrupt the stent position on the delivery device. This is especially important during catheter removal from the packaging, placement of the guide wire, advancement through a guiding catheter or introducer sheath.
- Do not manipulate, touch or handle the stent with fingers or contact with liquids prior to the preparation and delivery as this may result in contamination or dislodgement of the stent from the balloon catheter.
- Do not expose or wipe the device with organic solvents such as alcohol or detergents.
- Use only the appropriate balloon inflation media. Do not use any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in deployment of the stent.
- When back loading catheter on the guide wire, provide adequate support to shaft segments.

6.3. Stent Placement Precautions:

- Do not prepare or pre-inflate the balloon prior to stent deployment, other than as directed.
- Do not induce vacuum (negative pressure) on the delivery balloon catheter before reaching the target lesion.
- Implantation of a stent may lead to dissection of the vessel distal and / or proximal to the stented portion and may cause acute closure of the vessel requiring additional intervention (e.g. further dilatation or placement of additional stents).
- Do not expand the stent if it is not properly positioned in the vessel.
- Long term outcome following repeat dilatation of endothelialized peripheral stents is unknown at present.
- Do not exceed Rated Burst Pressure as indicated on labelling. Use of pressures higher than those specified on the product label may result in a ruptured balloon and potential intimal damage and dissection.
- Guiding catheters/introducer sheath used must have lumen sizes that are suitable to accommodate the introduction Myra™ BMS stent.
- Stent retrieval methods (use of additional wires, snares or forceps) may result in additional trauma to vessel / duct and or the access site. Complications may include bleeding, hematoma or pseudoaneurysm.
- To avoid the possibility of dissimilar metal corrosion, do not implant the stents of different material in tandem overlap or contact if possible.
- If reinserting the catheter, flush the guide wire lumen using the flushing needle before insertion.

6.4. Stent/System Removal Precautions:

- Should any unusual resistance be felt at any time during either lesion access or removal of the stent delivery system pre-stent implantation, the entire system must be removed as a single unit.
- When removing the delivery system as a single unit, do not retract the delivery system into the guiding catheter.
- Failure to follow these steps and/or applying excessive force to the stent delivery system can potentially result in loss or damage to the stent and/or stent delivery system components.

6.5. Post Implant Precautions:

Great care must be exercised when crossing a newly deployed stent with a guide wire or balloon catheter to avoid disrupting the stent geometry.

7. Adverse Effects:

Undesirable effects/adverse events that may be associated with the implantation of a peripheral stent in native peripheral arteries include but are not limited to:

- Allergic reactions to contrast medium
- Amputation
- Aneurysm
- Arrhythmias
- Arterio-venous fistula
- Death
- Dissection
- Embolisation
- Hematoma
- Hemorrahage
- Hypotension/Hypertension
- Infection
- Pain and tenderness
- Perforation
- Pseudoaneurysm
- Pyrogenic reactions
- Reactions to antiplatelet/anticoagulation agents
- Restenosis
- Spasms
- Stent migration
- Thrombosis
- Total occulsion of Peripheral artery

There may be other potential adverse events that are unforeseen at this time.

8. Clinical Use Information:

8.1. Inspection prior to use:

- Open the product carton to remove the inner Tyvek pouch containing blister tray with sterile Myra™ BMS- Balloon Expandable Peripheral Stent System.
- Carefully inspect the sterile package before opening.
- Do not use if the package has been damaged or opened.
- The product should not be used after the "Use by" Date.
- If the sterile package appears intact, tear open the sterile pouch to carefully remove the product and pass on or drop the contents into the sterile field using aseptic technique.
- Carefully inspect for kinks, bends or other damages.
- Verify that the stent is located between the radiopaque markers.
- Do not use if any defects are noted.

8.2. Materials Required:

- Introducer sheath
- 2-3 syringes (10-20 cc)
- 1000 u/500 cc Normal heparinised saline (HepNS)
- 0.035" (0.89mm) diameter guide wire
- Contrast diluted 1:1 with normal saline
- Inflation device
- Three-way stopcock
- Torque device (If applicable)
- Guide wire introducer

8.3. Preparation:

8.3.1 Guidewire Lumen Flush:

- Remove the protective stylet from the guide wire lumen and discard.
- Flush the guide wire with sterile heparinised solution until the fluid exits the guide wire exit port.

Caution: Avoid manipulation of stent during flushing of guide wire lumen, as this may disrupt the placement of the stent on the balloon.

8.3.2 Delivery system Preparation:

- Prepare an inflation device with diluted contrast medium
- Attach inflation device to stopcock; attach to hub (balloon inflation port)

Caution: Do not apply negative or positive pressure to balloon at this time.

- With the tip down orient the delivery system vertically.
- Open stopcock to stent delivery system; pull negative for 30 sec; release to neutral for contrast fill.
- Close the stop cock to the delivery system; purge the inflation device / syringe of all air.
- Repeat previous three steps until air is expelled.
- If a syringe was used, attach a prepared inflation device to stopcock.
- Open the stopcock to the delivery system, leave on neutral.

8.3.3 Delivery Procedure:

- Prepare vascular access site according to standard practice.
- Wipe the exposed guide wire with Heparinized saline to remove residual blood or contrast medium.
- Maintain neutral pressure on inflation device.
- Backload delivery system onto proximal portion of guide wire while maintaining guide wire position across target lesion.
- Advance the stent delivery system over guide wire to target lesion. Use radiopaque balloon markers to position stent across lesion; perform angiography to confirm stent position.

Note: Should unusual resistance be felt at any time during the either lesion access or removal of the stent delivery system before stent implantation, the entire system should be removed. See Stent/System Removal Precautions, section 6.4 for specific stent delivery removal instructions.

- Final internal stent diameter should match reference vessel. ASSURE THAT THE STENT IS NOT UNDERDILATED.

8.3.6 Disposal:

After usage, dispose the system as per accepted regulations for medical waste management.

Note: Updated guidelines for the procedure have to be consulted through new publications.

9. Disclaimer of Warranty and Limitation of Remedy:

There is no express or implied warranty, including without limitation any implied warranty of merchantability of fitness for a particular purpose, on the Meril Life Sciences Pvt. Ltd. product(s) described in this publication. Under no circumstances shall Meril Life Sciences Pvt. Ltd. be liable for any direct, indirect, incidental or consequential damages. Some jurisdictions do not allow the exclusion or limitations on implied warranty. For jurisdictions which do not allow the exclusion or limitations of incidental or consequential damages, some of the above exclusions may not apply. The patient may also have other rights, which vary from jurisdiction to jurisdiction.

Descriptions or specifications in Meril Life Sciences Pvt. Ltd. printed matter, including this publications, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Meril Life Sciences Pvt. Ltd. will not be responsible for any direct, indirect, incidental, or consequential damages resulting from reuse of the product.

Symbols used in labeling

Stent inner diameter — Stent length
Contains one unit — Keep dry
Consult instruction for use — Non-Pyrogenic
Reference number — LOT
Max. guide wire diameter
Sterilized using ethylene oxide — Manufacturer
For single use only do not reuse — Manufacturing date
Do not use if box opened or damaged — Do not resterilize
EC REP Authorized representative in the European community — Temperature Limitation

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