DOCUMENT VERIFICATION

Prior to using this document, the user is responsible for verifying that the revision and effective date are current.

REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MCZH	
Rev.	**Effective Date**	**Change Order No.**	**Changes Made to Document**
100	28-MAR-2019	CO-10073976	**PWR:** 5477 This artwork replaces 10001117-001 (LLS10001117-001 REV 100). Updates include the following changes; • Legal Manufacturer address changed from Czech Republic to Nathan Lane, Minneapolis. • EU Representative Czech Republic address added. • New Notified Body, changed from CE0473 to CE2797. • Copyright date updated to 2019 in Trademark Statement. • EU Rep symbol and text added to symbol glossary.

LANGUAGES

en English	fr French	de German	it Italian	es Spanish	pt Portuguese	nl Dutch	sv Swedish	da Danish	no Norwegian
Y	Y	N	N	Y	Y	N	N	N	N
fi Finnish	**el** Greek	**ja** Japanese	**cs** Czech	**sk** Slovakian	**hu** Hungarian	**pl** Polish	**sl** Slovenian	**et** Estonian	**lv** Latvian
N	N	N	N	N	N	N	N	N	N
lt Lithuanian	**tr** Turkish	**bg** Bulgarian	**ro** Romanian	**ru** Russian	**hr** Croatian	**zh** Chinese	**ko** Korean	**ar** Arabic	
N	N	N	N	N	N	N	N	N	

SECTION 1 - Folding Diagram
SECTION 2 - IFU Layout
SECTION 3 - Czech Language IFU, Reference Text

ARTWORK DESIGN NOTES


Part Number: 10001117-002
Size: 297mm (W) x 210mm (H). **Finished Folded Size:** 148mm (W) x 210mm (H).
IFU Format: Printed Instruction Booklet, 20 pages.


A) For detailed IFU specifications refer to the purchase specification.
B) Artwork Content and Specification Approval via electronic approval.
C) Approval of this document is for final format and layout only. All translations have been previously approved on LLS10001117-001 Rev. 100.

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
SECTION 1 - Folding Diagram

= Stitch locations.
Internal pages stitched twice to cover.

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SECTION 2 - IFU Layout **(To scale)**

PORTEX®

Percutaneous Dilation Tracheostomy Kit

With Single Stage Dilator and Soft Introducers for 7, 8 and 9mm tubes only

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
en

Portex® Percutaneous Dilation Tracheostomy Kit

With Single Stage Dilator and Soft Introducers for 7, 8 and 9mm tubes only Without Tracheostomy Tube

INSTRUCTIONS FOR USE

Fig. 1

These instructions are intended for use with the following Portex® products:

Catalogue No:

100/572/000 *Percutaneous Dilation Tracheostomy Kit with Single Stage Dilator, and soft Introducers for 7, 8 and 9mm tubes* ***only****.* ***Without*** *Tracheostomy Tube.*

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THESE KITS MUST ONLY BE USED BY CLINICIANS TRAINED TO PERFORM PERCUTANEOUS DILATIONAL TRACHEOSTOMY.

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Single Stage Dilator with soft tip and ergonomically designed handle to provide a better grip and control during the procedure.

KIT CONTENTS ***(see Fig 1)***

1	Scalpel
2	10ml Syringe
3	14G Needle & Cannula
4	Guidewire and Introducer
5	14FR Short Dilator
6	Guiding Catheter
7	Pre-lubricated Single Stage Dilator
8	Soft Tracheostomy Tube Introducers (21, 24, 27FR) for size 7, 8 and 9mm Tracheostomy Tubes only
9	Swabs (not shown)

2 INDICATIONS FOR USE

To create a percutaneous dilational tracheostomy using guidewire, single stage dilator and components of this product which allow for tracheal access for airway management.

3 ABSOLUTE CONTRA-INDICATIONS

3.1 Emergency surgical airway management (in emergency, cricothyrotomy is the procedure of-choice).

3.2 Children.

3.3 Pre-existing infection at the tracheostomy site.

3.4 Pre-existing malignancy at the tracheostomy site.

3.5 Uncertainty in identifying the anatomical landmarks.

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4.1 Enlarged thyroid gland.

4.2 Previous surgery at the tracheostomy site (e.g. thyroidectomy).

4.3 Bleeding diathesis (e.g. due to anticoagulant therapy).

4.4 Pre-existing tracheomalacia.

4.5 Morbid obesity (skin to trachea distance may render tracheostomy tube too-short).

5 WARNING

This product contains DEHP – data derived from animal and epidemiological studies outline the following risks, but are not limited to: Adverse effects on fetal development, kidney and liver toxicity, negative impact on fertility. Special attention on the application of phthalate products should be given for the following patient populations: pediatrics, pregnant and nursing women. For further information on phthalates please refer to the Smiths Medical website: www.smiths-medical.com/phthalates

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6 CAUTION

Follow universal precautions as specified by the Center for Disease Control and Prevention, USA, or relevant advisory bodies in the country where this device is being-used.

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8 PREPARATION OF COMPONENTS

Caution: This kit does not contain a Portex® tracheostomy tube. Read and follow the complete instructions for use supplied with your tracheostomy tube before use.

8.1 A small amount of sterile lubricating gel should be applied to the outside of the tracheostomy tube to reduce insertion forces and to the introducer to ease its movement in and out of the tracheostomy tube.

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8.3 Select an appropriately sized purple Introducer ***(Fig 1.8)*** to ensure a close fi t in the tracheostomy tube tip. Ensure the tracheostomy tube introducer moves freely in and can be easily removed from the tube.

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8.5 Check that the guiding catheter ***(Fig 1.6)*** moves freely through the relevant dilators, introducers and over the guidewire.

9 PREPARATION OF PATIENT

9.1 The patient should be positioned supine, with a pillow to support the neck and shoulders in order to hyper extend the neck ***(Fig. 2)***.

9.2 Locate the thyroid cartilage between the thumb and forefinger, and identify and mark the anatomical landmarks ***(Fig. 3)***.

(a) Thyroid cartilage
(b) Cricoid cartilage
(c) Subcricoid Space
(d) Tracheal rings
(e) Sternal notch
(f) Possible insertion sites

Suggested sites for tube insertion are marked **(x)** and are either between the first and second or second and third tracheal rings.

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9.4 Suction pharynx and, deflate the tracheal tube cuff . With bronchoscopic guidance, position the tracheal tube cuff just below the level of the vocal cords and re-inflate the cuff to re-establish a seal. This is necessary to avoid the tube interfering with the procedure. Hold the tube in position ensuring that the head and neck are in the mid-line and that the airway is being maintained.

10 PREPARATION OF THE SITE FOR SURGERY

10.1 Disinfect the skin and apply sterile drapes.

10.2 Palpate cricoid cartilage and infi ltrate local anaesthetic. Injection of local anaesthetic solutions containing adrenaline may help reduce bleeding at the incision site.

WARNING:

- **Performance of this procedure under bronchoscopic guidance is recommended to confirm midline placement of the 14G needle and cannula, guidewire and subsequent tracheostomy tube. Confirming tracheostomy tube placement following the procedure will prevent paratracheal placement that can lead to hypoxia, subcutaneous emphysema and/or pneumothorax.**
- **Care must be taken not to cause trauma to the carina and posterior-tracheal wall with dilators/introducers in small/short/low weight patients.**

11 THE PROCEDURE

WARNING:

- **The patient's airway must be established with a tracheal tube before tracheostomy is attempted to secure a patent airway, ensure adequate oxygenation during the procedure and to allow rapid re-intubation in cases of difficulty. If the patient is intubated, this is not an issue.**

11.1 Make a horizontal or vertical incision ***(Fig. 4)*** at the chosen insertion site of sufficient size to accept a tracheostomy tube (1.5-2.0cm). At this stage it may be of benefit to perform an exploratory blunt dissection **(keeping in the midline)** to identify the anatomical landmarks.

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WARNING:
- **The presence of anomalous blood vessels may cause bleeding during the procedure. If bleeding occurs, direct pressure will control most cases. If bleeding remains uncontrolled, the wound must be explored and the bleeding controlled using an open surgical technique.**

Caution: Confirm the tracheal tube has been withdrawn away from the insertion site to avoid puncturing the cuff or impaling the tube tip with the needle and cannula.

11.2 Insert the 14G needle and cannula (with syringe attached) in the midline into the insertion site ***(Fig. 5)*** in a caudal direction, to avoid the likelihood of the guidewire being passed into the pharynx at **stage 10.5**. Advance the needle (if necessary stepping off a tracheal ring) until free withdrawal of air into the syringe confirms the entry of the needle and cannula into the trachea. This may be facilitated by using a fluid filled syringe, allowing visualisation of withdrawn air bubbles. Some secretions and mucus may be aspirated during this process, which is normal.

WARNING:
- **Once the needle and cannula has entered the trachea, it must be held immobile so that the tip does not drift out of the trachea which will result in incorrect placement of the guidewire.**
- **Care must be exercised to ensure the needle, guidewire, guiding catheter and dilators do not perforate the posterior tracheal wall resulting in potentially life threatening injury.**

11.3 Gentle in and out movement of the tracheal tube will determine if it has been-impaled by the needle as it will swing up and down as the tube is moved in-and out. In such cases the needle and cannula should be removed and-re-inserted after the tracheal tube has been withdrawn at least 2cm.

11.4 Leaving the cannula in place, withdraw the needle and syringe. Connect the syringe to the-cannula and withdraw air to re-confirm that the cannula is in the trachea. Remove the-syringe.

WARNING:
- **Once the needle has been partially or completely withdrawn from the cannula, never re-insert it as this could result in the catheter being cut and left in the trachea.**

11.5 Ease the guidewire introducer out from its sheath and straighten the 'J' tip, leaving a sufficient length of exposed guidewire (2-3cm) to enable its dispensing with the forefinger and thumb ***(Fig. 6)***. Insert the introducer into the cannula and feed the guidewire out of the sheath until the first mark on the guidewire is positioned at the skin ***(Fig. 7)***.

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Caution: As the guidewire will be used to guide components into the trachea, it is important to ensure that it is not damaged (e.g. kinked) during the procedure. To help avoid this, always anchor the free end of the wire with forefinger and thumb whilst advancing components into the patient. During the procedure, check the guidewire is not kinked by ensuring its free movement through the components and the trachea.

WARNING:
- **If the guidewire is damaged during the procedure, it may be difficult to continue resulting in trauma to the trachea and misplacement of the tracheostomy tube. In these cases it may be possible to proceed only if:**
 - **(a) the damaged section can be advanced into the trachea, and**
 - **(b) there is a sufficient length of undamaged guidewire remaining.**

 If the guidewire cannot be salvaged, it will be necessary to use a new kit or a new guidewire.

11.7 Pass the 14FR short dilator over the guidewire through the soft tissues, until resistance is felt on the tracheal wall. With a gentle twisting movement of the dilator ***(Fig. 8)*** push the dilator forward to penetrate the anterior tracheal wall, dilating both the tissues and tracheal wall. Remove the 14FR short dilator, leaving the guidewire in place.

Caution: Confirm the free movement of the guidewire in the trachea and short dilator before removing the dilator, leaving the guidewire in place.

11.8 Pass the long guiding catheter over the guidewire into the trachea in the direction of the arrow marked on the catheter (safety stop end first) until the safety stop on the guiding catheter is located at the skin. Align the proximal end of the guiding catheter with the proximal band mark on the guidewire to determine the depth of insertion ***(Fig.9)***.

11.9 Single Stage Dilation
Immediately prior to insertion, immerse the distal end of the "single stage" dilator in sterile water or saline to activate the lubricious coating on the dilator. Pass the dilator over the guiding catheter until it reaches the "safety stop". In this position, the proximal mark on the guiding catheter will just be visible at the handle end of the dilator (***Fig. 10***).

Whilst stabilising the guidewire and guiding catheter to ensure they remain stable and in position in the trachea, insert and partially remove the dilator several times in order to slightly over-dilate the trachea to a size appropriate for the tracheostomy tube to be inserted. The dilator is marked, for guidance, with 38FR and maximum insertion depth ***(Fig. 10)***.

Caution: Do not pull the guiding catheter out of the trachea when removing the dilator. An assistant must stabilize the guidewire and guiding catheter during the dilation procedure to ensure that they remain correctly positioned within the trachea.

WARNING:
- **It is essential to ensure the correct alignment of the proximal end of the guiding catheter with the proximal mark on the guidewire. This will ensure the distal end of the guiding catheter does not extend beyond the tip of the guidewire where it could cause trauma to the posterior tracheal wall.**
- **Do not insert the thick black "Maximum Insertion Depth" mark on the "single stage" dilator beyond the level of the skin as this may cause damage to the trachea and/or carina.**
- **Do not use excessive twisting motion when dilating the stoma because the angled tip of the dilator may cause trauma or perforate the trachea.**
- **Portex® Purple Introducers with soft tips must only be used for introduction of Tracheostomy tubes, and not for dilation of the stoma. Their flexibility could cause them to buckle under the forces required for dilation resulting in trauma to the trachea and neck.**

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11.10 Tracheostomy Tube Insertion
Insert the lubricated tracheostomy tube loaded on its lubricated introducer over the guiding catheter through the stoma with a slight twisting motion (Fig. 11).

WARNING:
- **Do not use excessive twisting motion when inserting the tracheostomy tube because the tip of the tube may cause trauma or perforate the trachea.**

11.11 Remove the introducer, guiding catheter and guidewire leaving the tracheostomy tube in place.

11.12 Suction the trachea and tracheostomy tube to establish a clear airway.

WARNING:
- **Blood, which is not suctioned from the airway readily, clots and may produce serious airway obstruction.**

11.13 Transfer the breathing system to the tracheostomy tube.

11.14 Inflate the tracheostomy tube cuff with the minimum volume of air to form a seal.

11.15 Confirm successful tube placement using currently accepted medical techniques such as auscultation, capnography, bronchoscopic visualisation and/or chest X-ray and oximetry.

WARNING:
- **Verify the tracheostomy tube position by bronchoscopic visualization and/or chest X-ray to ensure correct placement. Incorrect placement could result in trauma to the trachea or respiratory obstruction.**

11.16 Secure the tracheostomy tube.

11.17 Deflate the tracheal tube cuff and remove the tube.

11.18 A tracheostomy dressing may be used at the stoma for improved patient comfort and hygiene.

Caution: Dispose of this product in a safe manner according to local guidelines for disposal of contaminated waste.


Portex, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2019 Smiths Medical. All rights reserved.

Fig. 2

Fig. 3

(a) Thyroid cartilage
(b) Cricoid cartilage
(c) Subcricoid Space
(d) Tracheal rings
(e) Sternal notch
(f) Possible insertion sites

Fig. 4

Fig. 5


Fig. 6


Fig. 7


Fig. 8


Fig. 9


Fig. 10



Fig. 11

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Title: **INSTRUCTIONS FOR USE - PERCUTANEOUS DILATION TRACHEOSTOMY KIT WITHOUT BLU (10018928-001)**

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REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MCZH	
Rev.	**Effective Date**	**Change Order No.**	**Changes Made to Document**
100	28-MAR-2019	CO-10073976	**PWR:** 5477 This artwork replaces 004/005/245A (L6394 Revision 103). Updates include the following changes; • Legal Manufacturer address changed from Czech Republic to Nathan Lane, Minneapolis. • EU Representative Czech Republic address added. • New Notified Body, changed from CE0473 to CE2797. • Copyright date updated to 2019 in Trademark Statement. • EU Rep symbol and text added to symbol glossary.

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ARTWORK DESIGN NOTES

Part Number: 10018928-001
Finished Folded Size: 210mm (W) x 260mm (H).
IFU Format: Printed Instruction Booklet, 52 pages.

A) For detailed IFU specifications refer to the purchase specification.
B) Artwork Content and Specification Approval via electronic approval.
C) Approval of this document is for final format and layout only. All translations have been previously approved on L6394 ISSUE - 100.

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Portex® Percutaneous Dilation Tracheostomy Kit

INSTRUCTIONS FOR USE

These instructions apply to the use of the following Portex® products:

Catalogue No. 100/562/000 Internal REF: 100/562/000CZ — Percutaneous Dilatational Tracheostomy Kit with Single Stage Dilator, and soft Introducers for 7, 8 and 9mm tubes **only**. **Without** Tracheostomy Tube.

Catalogue No. 100/566/000 Internal REF: 100/566/000CZ — Percutaneous Dilation Tracheostomy Kit with Serial Dilators/Introducers for 7, 8, and 9mm tubes **only**. **Without** Tracheostomy Tube.

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

Portex® single use Percutaneous Dilation Tracheostomy Kits allow the percutaneous insertion of a tracheostomy tube using either a serial (sequential) or single (one stage) circumferencial dilatational Seldinger guidewire technique.

The kits are intended for use in a controlled setting such as an Intensive Care Unit or operating room with the assistance of trained personnel. A **minimum** of 2 operators are required – one to maintain the patient's airway, anaesthesia, breathing and circulation, and the other to perform the procedure.

Dilators have a lubricious hydrophilic coating when wetted, to improve ease of insertion.

WARNING

- ***READ THE COMPLETE INSTRUCTIONS FOR USE BEFORE USE.***
- ***THESE KITS MUST ONLY BE USED BY CLINICIANS TRAINED TO PERFORM PERCUTANEOUS DILATIONAL TRACHEOSTOMY.***

1 Description

The Portex® Percutaneous Dilation Tracheostomy Kit incorporates the following features:

- Single Stage Dilator with soft tip and ergonomically designed handle to provide a better grip and control during the procedure.

2 Indications for Use

Controlled, elective subcricoid percutaneous insertion of a tracheostomy tube for airway management using a Seldinger guidewire dilation technique.

3 Absolute Contra-indications

3.1 Emergency surgical airway management (in emergency, cricothyrotomy is the procedure of choice).
3.2 Children.
3.3 Pre-existing infection at the tracheostomy site.
3.4 Pre-existing malignancy at the tracheostomy site.
3.5 Uncertainty in identifying the anatomical landmarks.
3.6 Unstable cervical spine fracture.

4 Relative Contra-indications

Risks in these situations, must be balanced against the benefits of the procedure.

4.1 Enlarged thyroid gland.
4.2 Previous surgery at the tracheostomy site (e.g. thyroidectomy).
4.3 Bleeding diathesis (e.g. due to anticoagulant therapy).
4.4 Pre-existing tracheomalacia.
4.5 Morbid obesity (skin to trachea distance may render tracheostomy tube too short).

Kit Contents

Contents *(see Fig. 1)*	**100/562**	**100/566**
1 Scalpel	✓	✓
2 10ml Syringe	✓	✓
3 14G Needle & Cannula	✓	✓
4 Guidewire and Introducer	✓	✓
5 14FR Short Dilator	✓	✓
6 Guiding Catheter	✓	✓
7 Pre-lubricated Single Stage Dilator	✓	
8 Pre-lubricated Serial Dilator Set (18, 21, 24, 27, 30, 34, 38 FR)		✓
9 Soft Tracheostomy Tube Introducers (21, 24, 27 FR) for size 7, 8 and 9mm Tracheostomy Tubes **only**	✓	
10 Sterile Lubricating Jelly 5g *(not shown)*	✓	✓
11 Swabs *(not shown)*	✓	✓

5 Warnings

5.1 The patient's airway must be established with a tracheal tube before tracheostomy is attempted.
5.2 Performance of this procedure under bronchoscopic guidance is recommended to confirm midline placement of the 14G needle and cannula, guidewire and subsequent tracheostomy tube to prevent paratracheal placement that can lead to hypoxia, subcutaneous emphysema and/or pneumothorax.
5.3 The presence of anomalous blood vessels may cause bleeding during the procedure. If bleeding occurs, direct pressure will control most cases. If bleeding remains uncontrolled, the wound must be explored using an open surgical technique and the bleeding stopped.
5.4 Care must be taken not to cause trauma to the carina with dilators/introducers in small/short/low weight patients.
5.5 If the guidewire is damaged during the procedure, it may be difficult to continue. In these cases it may be possible to proceed only if:
a) the damaged section can be advanced into the trachea, and
b) there is a sufficient length of undamaged guidewire remaining.
If the guidewire cannot be salvaged, it will be necessary to use a new kit or a new guidewire.
5.6 Once the needle and cannula has entered the trachea, it must be held absolutely immobile so that the tip does not drift out of the trachea.
5.7 Care must be exercised to ensure the needle, guidewire, guiding catheter and dilators do not perforate the posterior tracheal wall.
5.8 Once the needle has been partially or completely withdrawn from the cannula, never re-insert it as this could result in the catheter being cut and left in the trachea.
5.9 It is essential to ensure the correct alignment of the proximal mark on the guiding catheter with the proximal mark on the guidewire to ensure the distal end of the guiding catheter does not extend beyond the tip of the guidewire and cause trauma to the posterior tracheal wall.

This document and all information contained herein is confidential, proprietary and may not be disclosed, disseminated or reproduced without written approval from Smiths Medical.

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5.10 Do not insert the thick black "Maximum Insertion Depth" mark on the "single stage" dilator beyond the level of the skin to prevent damage to the trachea and/or carina.

5.11 Do not pull the guiding catheter out of the trachea when removing the dilators. An assistant must stabilize the guidewire and guiding catheter during the dilation procedure to ensure that they remain correctly positioned within the trachea.

5.12 Portex® Purple Introducers with soft tips must only be used for introduction of Tracheostomy tubes, and not dilation of the stoma due to their flexibility.

5.13 Do not use excessive twisting motion when dilating the stoma to avoid trauma to the trachea.

5.14 Do not use excessive twisting motion when inserting the tracheostomy tube to avoid trauma to the trachea.

5.15 Blood, which is not suctioned from the airway readily, clots and may produce serious airway obstruction.

5.16 Verify the tracheostomy tube position by bronchoscopic visualization and/or chest X-ray to ensure correct placement.

5.17 Prior to removal of cuffed Tracheostomy tubes, all air must be completely removed from the cuff to prevent damage to the trachea and stoma.

5.18 Ensure the sterile lubricating jelly does not occlude the lumen of the tube preventing ventilation of the patient.

5.19 In order to avoid tracheal damage from over inflation of the cuff, the cuff should not be inflated with a measured volume of air, or by feel of pressure from the syringe since little resistance should be felt during cuff inflation.

5.20 During anaesthesia, nitrous oxide may diffuse into the cuff causing an increase or decrease in cuff pressure.

5.21 Cuff pressures should be monitored and adjusted routinely. Overinflation of the cuff may result in permanent damage to the trachea.

5.22 Contact with electrosurgery electrodes or laser surgery beams must be avoided because PVC will produce toxic fumes in air or ignite in an enriched oxygen environment (e.g. anaesthesia).

5.23 Do not use Portex® tracheostomy tubes in patients with abnormal upper airway anatomy or pathology as this may result in partial or total airway obstruction.

5.24 During and after attachment of the breathing system to the tracheostomy tube connector, avoid application of excessive rotational or linear forces on the tube to prevent accidental disconnection or occlusion.

6 Cautions

6.1 The security of all breathing system connectors should be checked when the circuit is established and frequently thereafter. Disconnection may be facilitated with the use of a disconnection wedge (100/555/000).

6.2 Patients should be adequately humidified e.g. using a Portex® Thermovent® HEPA+ (100/586/000) or Portex® Thermovent® T (100/570/015) to minimise encrustation of the Tracheostomy tube and/or inner cannula lumen and prevent mucosal damage.

6.3 Ensure the tracheostomy tube inner cannula, if supplied, is removed before inserting the tracheostomy tube introducer.

6.4 Confirm the tracheal tube has been withdrawn away from the insertion site to avoid puncturing the cuff or impaling the tube tip with the needle and cannula.

6.5 As the guidewire will be used to guide components into the trachea, it is important to ensure that it is not damaged (e.g. kinked) during the procedure. To help avoid this, always anchor the free end of the wire with forefinger and thumb whilst advancing components into the patient. During the procedure, check the guidewire is not kinked by ensuring its free movement through the components and the trachea.

6.6 Confirm the free movement of the guidewire in the trachea and shortdilator before removing the dilator, leaving the guidewire and guiding catheter in place.

6.7 The patency of the tracheostomy tube lumen must be assured by regular suctioning and regular cleaning of the inner cannula (if one is used). Check routinely and replace as required to maintain a patent airway.

6.8 Devices used in or during inflation of the cuff must be clean and free from all foreign matter. The inflation device should be removed from the inflation valve immediately after use, and the dust cap fitted. Test inflate the cuff before insertion.

6.9 Guard against cuff damage by avoiding contact with sharp edges.

6.10 The inflation line valve may interfere with Magnetic Resonance Imaging (MRI) picture clarity. Ensure the valve is positioned away from the area being scanned.

6.11 Tracheostomy tubes must be changed regularly to suit individual patient's needs.

6.12 Only sterile saline should be used to clean accessible parts of the tube whilst in-situ. If the tracheostomy tube is removed from the patient, it must be discarded.

6.13 Repositioning of the in-situ tracheostomy tube while the cuff is inflated should be avoided.

6.14 **Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

6.15 Dispose of this product in a safe manner according to local guidelines for disposal of contaminated waste.

6.16 Follow universal precautions as specified by the Centre for Disease Control and Prevention, USA, or relevant advisory bodies in the country where this device is being used.

7 Adverse Reactions

Reported adverse reactions associated with tracheostomy tubes are many and diverse. Consult standard textbooks and medical literature for specific adverse reaction information.

8 Preparation of Components

8.1 A small amount of the sterile lubricating jelly (supplied) should be applied to the outside of the tracheostomy tube to reduce insertion forces and to the introducer to ease its movement in and out of the tracheostomy tube.

WARNING:

Ensure the sterile lubricating jelly does not occlude the lumen of the tube preventing ventilation of the patient.

8.2 Check the tracheostomy tube cuff by test inflation.

8.3 For Kit 100/562/000 select an appropriately sized purple Introducer *(Fig. 1.9)* or, for Kit 100/566/000 a white dilator *(Fig. 1.8)* to ensure a close fit in the tracheostomy tube tip. Ensure the tracheostomy tube introducer moves freely in and can be easily removed from the tube.

8.4 If satisfactory, **fully deflate** the cuff (smoothing it back towards the flange) in order to avoid rupture during tube insertion.

8.5 Check that the guiding catheter *(Fig. 1.6)* moves freely through the relevant dilators, introducers and over the guidewire.

Caution:

- Follow universal precautions as specified by the Centre for Disease Control and Prevention, USA, or relevant advisory bodies in the country where this device is being used.
- Ensure the tracheostomy tube inner cannula, if supplied, is removed before inserting the tracheostomy tube introducer.

This document and all information contained herein is confidential, proprietary and may not be disclosed, disseminated or reproduced without written approval from Smiths Medical.

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9 Preparation of Patient

9.1 The patient should be positioned supine, with a pillow to support the neck and shoulders in order to hyper extend the neck *(Fig. 2)*.

9.2 Locate the thyroid cartilage between the thumb and forefinger, and identify and mark the anatomical landmarks *(Fig. 3)*.

(a) Thyroid cartilage

(b) Cricoid cartilage

(c) Tracheal rings

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(e) Possible insertion sites

Suggested sites for tube insertion are marked (x) and are either between the first and second or second and third tracheal rings.

9.3 Prior to the procedure, increase FiO_2 to 100% and monitor the patient using oximetry, capnography, ECG and blood pressure. Adjust the ventilator settings to compensate for leaks during the procedure.

9.4 Suction pharynx and, deflate the tracheal tube cuff. With bronchoscopic guidance, position the tracheal tube cuff just below the level of the vocal cords and re-inflate the cuff to re-establish a seal. This is necessary to avoid the tube interfering with the procedure. Hold the tube in position ensuring that the head and neck are in the mid-line and that the airway is being maintained.

10 Preparation of the site for surgery

10.1 Disinfect the skin and apply sterile drapes.

10.2 Palpate cricoid cartilage and infiltrate local anaesthetic. Injection of local anaesthetic solutions containing adrenaline may help reduce bleeding at the incision site.

WARNING:

- ***Performance of this procedure under bronchoscopic guidance is recommended to confirm midline placement of the 14G needle and cannula, guidewire and subsequent tracheostomy tube to prevent paratracheal placement that can lead to hypoxia, subcutaneous emphysema and/or pneumothorax.***
- ***Care must be taken not to cause trauma to the carina with dilators/ introducers in small/short/low weight patients.***

11 The Procedure

WARNING:

The patient's airway must be established with a tracheal tube before tracheostomy is attempted.

11.1 Make a horizontal or vertical incision *(Fig. 4)* at the chosen insertion site of sufficient size to accept a tracheostomy tube (1.5-2.0cm). At this stage it may be of benefit to perform an exploratory blunt dissection (keeping in the midline) to identify the anatomical landmarks.

WARNING:

The presence of anomalous blood vessels may cause bleeding during the procedure. If bleeding occurs, direct pressure will control most cases. If bleeding remains uncontrolled, the wound must be explored using an open surgical technique and the bleeding stopped.

Caution:

Confirm the tracheal tube has been withdrawn away from the insertion site to avoid puncturing the cuff or impaling the tube tip with the needle and cannula.

11.2 Insert the 14G needle and cannula (with syringe attached) in the midline into the insertion site *(Fig. 5)* in a caudal direction, to avoid the likelihood of the guidewire being passed into the pharynx at **stage 11.5**. Advance the needle (if necessary stepping off a tracheal ring) until free withdrawal of air into the syringe confirms the entry of the needle and cannula into the trachea. This may be facilitated by using a fluid filled syringe, allowing visualisation of withdrawn air bubbles. Some secretions and mucus may be aspirated during this process, which is normal.

WARNING:

- ***Once the needle and cannula has entered the trachea, it must be held absolutely immobile so that the tip does not drift out of the trachea.***
- ***Care must be exercised to ensure the needle, guidewire, guiding catheter and dilators do not perforate the posterior tracheal wall.***

11.3 Gentle in and out movement of the tracheal tube will determine if it has been impaled by the needle as it will swing up and down as the tube is moved in and out. In such cases the needle and cannula should be removed and re-inserted after the tracheal tube has been withdrawn at least 2cm.

11.4 Leaving the cannula in place, withdraw the needle and syringe. Connect the syringe to the cannula and withdraw air to re-confirm that the cannula is in the trachea. Remove the syringe.

WARNING:

Once the needle has been partially or completely withdrawn from the cannula, never re-insert it as this could result in the catheter being cut and left in the trachea.

11.5 Ease the guidewire introducer out from its sheath and straighten the 'J' tip, leaving a sufficient length of exposed guidewire (2-3cm) to enable its dispensing with the forefinger and thumb *(Fig. 6)*. Insert the introducer into the cannula and feed the guidewire out of the sheath until the first mark on the guidewire is positioned at the skin *(Fig. 7)*.

11.6 Confirm free movement of the guidewire in the cannula before removing the cannula, leaving the guidewire in place.

Caution:

As the guidewire will be used to guide components into the trachea, it is important to ensure that it is not damaged (e.g. kinked) during the procedure. To help avoid this, always anchor the free end of the wire with forefinger and thumb whilst advancing components into the patient. During the procedure, check the guidewire is not kinked by ensuring its free movement through the components and the trachea.

WARNING:

If the guidewire is damaged during the procedure, it may be difficult to continue. In these cases it may be possible to proceed only if:

(a) the damaged section can be advanced into the trachea, and

(b) there is a sufficient length of undamaged guidewire remaining.

If the guidewire cannot be salvaged, it will be necessary to use a new kit or a new guidewire.

11.7 Pass the 14FR short dilator over the guidewire through the soft tissues, until resistance is felt on the tracheal wall. With a gentle twisting movement of the dilator *(Fig. 8)* push the dilator forward to penetrate the anterior tracheal wall, dilating both the tissues and tracheal wall. Remove the 14FR short dilator, leaving the guidewire in place.

Caution:

Confirm the free movement of the guidewire in the trachea and short dilator before removing the dilator, leaving the guidewire in place.

11.8 Pass the long guiding catheter over the guidewire into the trachea in the direction of the arrow marked on the catheter (safety stop end first) until the safety stop on the guiding catheter is located at the skin. Align the proximal end of the guiding catheter with the proximal band mark on the guidewire to determine the depth of insertion *(Fig.9)*.

11.9 Single Stage Dilation

(For serial dilation go to step 11.10)

Immediately prior to insertion, immerse the distal end of the "single stage" dilator in sterile water or saline to activate the lubricious coating on the dilator. Pass the dilator over the guiding catheter until it reaches the "safety stop". In this position, the proximal mark on the guiding catheter will just be visible at the handle end of the dilator *(Fig. 10)*.

This document and all information contained herein is confidential, proprietary and may not be disclosed, disseminated or reproduced without written approval from Smiths Medical.

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Whilst stabilising the guidewire and guiding catheter to ensure they remain stable and in position in the trachea, insert and partially remove the dilator several times in order to slightly over-dilate the trachea to a size appropriate for the tracheostomy tube to be inserted. The dilator is marked, for guidance, with 38FR and maximum insertion depth *(Fig. 10)*.

WARNING:

- ***It is essential to ensure the correct alignment of the proximal mark on the guiding catheter with the proximal mark on the guidewire to ensure the distal end of the guiding catheter does not extend beyond the tip of the guidewire and cause trauma to the posterior tracheal wall.***
- ***Do not insert the thick black "Maximum Insertion Depth" mark on the "single stage" dilator beyond the level of the skin to prevent damage to the trachea and/or carina.***
- ***Do not use excessive twisting motion when dilating the stoma to avoid trauma to the trachea.***
- ***Portex® Purple Introducers with soft tips must only be used for introduction of Tracheostomy tubes, and not dilation of the stoma due to their flexibility.***
- ***Do not pull the guiding catheter out of the trachea when removing the dilators. An assistant must stabilize the guidewire and guiding catheter during the dilation procedure to ensure that they remain correctly positioned within the trachea.***

Remove the dilator leaving the guidewire and guiding catheter in the trachea. Air will exit through the tracheostomy stoma confirming the correct placement of the dilator, guiding catheter and guidewire in the trachea.

Proceed to stage 11.12.

11.10 Serial dilation

Immediately prior to insertion, immerse the distal end of each dilator in sterile water or saline to activate the lubricious coating. Pass the dilator over the guiding catheter until it reaches the safety stop. In this position, the proximal mark on the guiding catheter will be just visible at the handle end of the dilator *(Fig. 11)*.

Dilate the stoma by passing each dilator, in turn (beginning with the smallest), over the guiding catheter into the trachea ensuring the correct alignment of the dilator, guiding catheter and guidewire *(Fig. 11)*. The final dilator should be selected in order to slightly over dilate the trachea to a size appropriate for the size of tracheostomy tube to be inserted.

WARNING:

- ***It is essential to ensure the correct alignment of the proximal mark on the guiding catheter with the proximal mark on the guidewire to ensure the distal end of the guiding catheter does not extend beyond the tip of the guidewire and cause trauma to the posterior tracheal wall.***
- ***Do not insert the thick black mark on the "single stage" dilator beyond the level of the skin in order to prevent damage to the trachea.***
- ***Portex® Purple Introducers with soft tips must only be used for introduction of Tracheostomy tubes, and not dilation of the stoma due to their flexibility.***
- ***Do not pull the guiding catheter out of the trachea when removing the dilators. An assistant must stabilize the guidewire and guiding catheter during the dilation procedure to ensure that they remain correctly positioned within the trachea.***
- ***Do not use excessive twisting motion when dilating the stoma to avoid trauma to the trachea.***

11.11 Maintain manual pressure with a sponge over the insertion site between dilators to control bleeding. Air should bubble freely from the stoma between dilator insertions confirming correct placement of the dilators, guiding catheter and guidewire in the trachea.

Tracheostomy Tube Insertion

11.12 Insert the lubricated tracheostomy tube loaded on its lubricated introducer over the guiding catheter through the stoma with a slight twisting motion *(Fig. 12)*.

WARNING:

- **Do not use excessive twisting motion when inserting the tracheostomy tube to avoid trauma to the trachea.**
- **Do not use Portex® Tracheostomy tubes in patients with abnormal upper airway anatomy or pathology as this may result in partial or total airway obstruction.**

11.13 Remove the introducer, guiding catheter and guidewire leaving the tracheostomy tube in place.

11.14 Suction the trachea and tracheostomy tube to establish a clear airway.

WARNING:

Blood, which is not suctioned from the airway readily, clots and may produce serious airway obstruction.

11.15 Transfer the breathing system to the tracheostomy tube.

WARNING:

During and after attachment of the breathing system to the tracheostomy tube connector, avoid application of excessive rotational or linear forces on the tube to prevent accidental disconnection or occlusion.

11.16 Inflate the tracheostomy tube cuff with the minimum volume of air to form a seal.

WARNING:

In order to avoid tracheal damage from over inflation of the cuff,

the cuff should not be inflated with a measured volume of air, or by feel of pressure from the syringe since little resistance should be felt during cuff inflation.

11.17 Confirm successful tube placement using currently accepted medical techniques such as auscultation, capnography, bronchoscopic visualisation and/or chest X-ray and oximetry.

WARNING:

Verify the tracheostomy tube position by bronchoscopic visualisation and/or chest X-ray to ensure correct placement.

11.18 Secure the tracheostomy tube.

11.19 Deflate the tracheal tube cuff and remove the tube.

11.20 A tracheostomy dressing may be used at the stoma for improved patient comfort and hygiene.

Caution:

- Tracheostomy tubes must be changed regularly to suit individual patient's needs.
- Dispose of this product in a safe manner according to local guidelines for disposal of contaminated waste.

Thermovent, Portex, Portex design mark and Smiths Medical design mark are trademarks Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.
© 2019 Smiths Medical. All rights reserved.

This document and all information contained herein is confidential, proprietary and may not be disclosed, disseminated or reproduced without written approval from Smiths Medical.

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REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MCZH	
Rev.	**Effective Date**	**Change Order No.**	**Changes Made to Document**
100	18-MAR-2019	CO-10074820	**PWR:** 5477 This artwork replaces 004/005/225A (L6402 Issue 100). Updates include the following changes; • Legal Manufacturer address changed from Czech Republic to Nathan Lane, Minneapolis. • EU Representative Czech Republic address added. • New Notified Body, changed from CE0473 to CE2797. • Copyright date updated to 2019 in Trademark Statement. • EU Rep symbol and text added to symbol glossary.

LANGUAGES

en English	fr French	de German	it Italian	es Spanish	pt Portuguese	nl Dutch	sv Swedish	da Danish	no Norwegian
Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
fi Finnish	**el** Greek	**ja** Japanese	**cs** Czech	**sk** Slovakian	**hu** Hungarian	**pl** Polish	**sl** Slovenian	**et** Estonian	**lv** Latvian
Y	Y	N	Y	N	Y	Y	N	N	N
lt Lithuanian	**tr** Turkish	**bg** Bulgarian	**ro** Romanian	**ru** Russian	**hr** Croatian	**zh** Chinese	**ko** Korean	**ar** Arabic	
N	N	N	N	N	N	N	N	N	

ARTWORK DESIGN NOTES

Part Number: 10019334-001
Size: 210mm (W) x 160mm (H). **Finished Folded Size:** 105mm (W) x 160mm (H).
IFU Format: Printed Instruction Booklet, 24 pages.


A) For detailed IFU specifications refer to the purchase specification.
B) Artwork Content and Specification Approval via electronic approval.
C) Approval of this document is for final format and layout only. All translations have been previously approved on L6402 Issue 100.

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

IFU ARTWORK (To scale)

PORTEX®

Thermovent® T Heat and Moisture Exchanger

en INSTRUCTIONS FOR USE
fr MODE D'EMPLOI
de GEBRAUCHSANWEISUNG
it ISTRUZIONI PER L'USO
es INSTRUCCIONES DE USO
pt INSTRUÇÕES DE UTILIZAÇÃO
nl GEBRUIKSAANWIJZING
sv BRUKSANVISNING
da BRUGSANVISNING
no BRUKSANVISNING
fi KÄYTTÖOHJEET
el ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
cs NÁVOD K POUŽITÍ
hu HASZNÁLATI UTASÍTÁS
pl INSTRUKCJA STOSOWANIA

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English

en

Portex® Thermovent® T Heat and Moisture Exchanger

INSTRUCTIONS FOR USE

These instructions are intended for use with the following Portex® product:

Catalogue No: 100/570/015 *Thermovent® T Heat and Moisture Exchanger*
Internal Ref: 100/570/015CZ

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Description

Thermovent® T is a sterile, single use heat and moisture exchanger which provides effective humidification for spontaneously breathing adult and paediatric patients whose upper airways are bypassed by a tracheal or tracheostomy tube.

Indications for use

The Portex® Thermovent® T is indicated for:

- Preserving humidity and minimising damage to tracheal epithelial cells
- Reduction of heat loss through the bronco pulmonary tree
- Helping to prevent thickened secretions and changes in lung function during long term anaesthesia or ventilatory support.

Instructions for use

1 The device should be connected to the 15mm termination of the tracheal or tracheostomy tube.
2 Where tenacious secretions or mucous plugging develop suitable therapy must be instituted.
3 The patient should be closely monitored for any complications resulting from blockage to the artificial airway or Thermovent® T. Should any complications develop, the device should be removed or replaced.
4 Change the device every 24 hours or more frequently as determined by expert clinical assessment to prevent accumulation of secretions in the device.

Precautions

1 As with all airway management products, the security of all connections and patency of the system and device should be checked prior to use.
2 This device should only be used with spontaneously breathing patients.
3 The effect of the additional deadspace volume of this device should be evaluated on an individual patient basis.
4 The device must not be used on conjunction with heated or nebulized water humidification systems or nebulizers.

Warnings

1 **The device is intended for single use only and must not be re-used or cleaned. Wetting of the heat and moisture exchange elements with any cleaning solution may dangerously increase the resistance to breathing or may result in the retention of harmful residues or transmission of infectious agents.**
2 **Expert clinical judgement must be used in assessing the patient's humidification requirements. This device must be used as described in these instructions and in accordance with currently accepted medical techniques.**
3 **If a heat and moisture exchange element is displaced from the housing (for example by patient coughing), the device should be discarded and a new device used. Under no circumstances must replacement of the element be attempted.**
4 **The device may block if the patient develops copious secretions, pulmonary oedema or bleeding when the device is in use. If the device becomes blocked then it must be replaced. Under no circumstances should any attempt be made to clear or clean out the blockage as this may damage the heat and moisture exchange elements and could lead to further blockage.**
5 **Caution should be taken if this device is used on children as it contains components which may be removed and could be harmful if ingested or inhaled.**

Technical Specification

1 **Moisture Output**
25mg/L H_2O @ 10 breath/min
Tidal Volume 1,000 ml
2 **Resistance to Flow**
0.9 hPa (cmH_2O) at gas flow rate 0.5L/sec (30L/min)
2.0 hPa (cmH_2O) at gas flow rate 1.0L/sec (60L/min)
3.5 hPa (cmH_2O) at gas flow rate 1.5L/sec (90L/min)

Thermovent, Portex, Portex design mark and Smiths Medical design mark are trademarks Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. © 2019 Smiths Medical. All rights reserved.

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en

Single use.

Do Not Reuse: Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only. Any attempt to re-process the device for subsequent re-use may adversely affect the integrity of the device or lead to deterioration in performance.

Sterile unless unit container is opened or damaged.

Destroy after single use. Do not resterilize.

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REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MCZH	
Rev.	**Effective Date**	**Change Order No.**	**Changes Made to Document**
100	19-MAR-2019	CO-10074966	**PWR:** 5477 This artwork replaces 004/005/241A (L6389 Issue 100). Updates include the following changes; • Legal Manufacturer address changed from Czech Republic to Nathan Lane, Minneapolis. • EU Representative Czech Republic address added. • Copyright date updated to 2019 in Trademark Statement. • EU Rep symbol and text added to symbol glossary.

LANGUAGES

en English	fr French	de German	it Italian	es Spanish	pt Portuguese	nl Dutch	sv Swedish	da Danish	no Norwegian
Y	Y	Y	Y	Y	Y	Y	Y	Y	N
fi Finnish	**el** Greek	**ja** Japanese	**cs** Czech	**sk** Slovakian	**hu** Hungarian	**pl** Polish	**sl** Slovenian	**et** Estonian	**lv** Latvian
Y	Y	N	Y	Y	N	Y	Y	N	N
lt Lithuanian	**tr** Turkish	**bg** Bulgarian	**ro** Romanian	**ru** Russian	**hr** Croatian	**zh** Chinese	**ko** Korean	**ar** Arabic	
Y	N	N	N	N	N	N	N	N	

ARTWORK DESIGN NOTES

Part Number: 10019361-001
Size: 297mm (W) x 210mm (H). **Finished Folded Size:** 148mm (W) x 105mm (H).
IFU Format: Printed Instruction Booklet, 16 pages.

A) For detailed IFU specifications refer to the purchase specification.
B) Artwork Content and Specification Approval via electronic approval.
C) Approval of this document is for final format and layout only. All translations have been previously approved on L6389 Issue 100.

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



Decannulation Cap

INSTRUCTIONS FOR USE

These instructions apply to the use of the following Portex® product:
Catalogue Number: 100/534/010 Internal REF: 100/534/010CZ

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1 Description:
A plastic molded cap for fitment to the 15 mm connection of Portex® Fenestrated Tracheostomy Tubes for the purpose of occluding the tracheostomy tube termination. The Decannulation Cap is colored red for easy identification.

2 Indications for Use:
The Portex® Decannulation Cap can be used to aid weaning or vocalization of patients that are fitted with a Portex® Fenestrated Tracheostomy Tube.

3 Contra-Indication:
There are no known contraindications.

4 Warnings:
4.1 The Decannulation Cap is to be used with Portex® Fenestrated Tracheostomy Tubes ONLY. Failure to do this may cause occlusion of the patients' airway; resulting in patient death orpermanent injury.

4.2 Do not force the Decannulation Cap onto the 15 mm termination as this may cause difficulty in removal, resulting in unintentional prolonged occlusion of the tracheostomy tube termination and may result in patient death or permanent injury. If difficulty occurs a Portex® Tracheostomy Wedge can be used to ease disconnection.

4.3 Caution must be exercised when using a Decannulation Cap in conjunction with a cuffed fenestrated tracheostomy tube. The user must ensure that the cuff is completely deflated and that a fenestrated inner cannula is used, or that the inner cannula has been removed. Failure to do so may result in total or partial obstruction of the airway. Observe the patient for signs of respiratory distress and *IMMEDIATELY* remove the decannulation cap if required, replace the inner cannula and reconnect the ventilator circuit.

5 Instructions for Use
5.1 Before fitting the Decannulation Cap, the patient's airway should be cleared by coughing and/or suctioning.

NOTE: *A Disposable Inner Cannula should always be used during suctioning to prevent the suction catheter protruding through the fenestrations and damaging the tracheal wall.*

5.2 Remove the Disposable Inner Cannula (if used).
5.3 Deflate the Tracheostomy Tube cuff (if fitted).
5.4 Attach the Decannulation Cap and closely monitor the patient's condition.
5.5 Remove the Decannulation Cap after the desired time. Check for airway patency. A Portex® Tracheostomy Wedge can be used to ease disconnection if required. Four wedges are supplied in each carton.
5.6 Insert a new Disposable Inner Cannula and connect to the breathing system if required.
5.7 If the single use Disposable Inner Cannula or Decannulation Cap become contaminated with accumulated secretions remove and dispose of according to hospital recommendations and replace with a new device.

CAUTION: Follow standard infection control procedures as specified by the Centers for Disease Control and Prevention (USA) or local equivalent.

Portex, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

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REVISION HISTORY

Owning Site: MASH		Additional Sites Using Document: MCZH	
Rev.	**Effective Date**	**Change Order No.**	**Changes Made to Document**
100	22-MAR-2019	CO-10073623	**PWR:** 5477 This artwork replaces 10001116-001 (LLS10001116-001 REV 100). Updates include the following changes; • Legal Manufacturer address changed from Czech Republic to Nathan Lane, Minneapolis. • EU Representative Czech Republic address added. • New Notified Body, changed from CE0473 to CE2797. • Copyright date updated to 2019 in Trademark Statement. • EU Rep symbol and text added to symbol glossary.

LANGUAGES

en English	fr French	de German	it Italian	es Spanish	pt Portuguese	nl Dutch	sv Swedish	da Danish	no Norwegian
Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
fi Finnish	**el** Greek	**ja** Japanese	**cs** Czech	**sk** Slovakian	**hu** Hungarian	**pl** Polish	**sl** Slovenian	**et** Estonian	**lv** Latvian
Y	Y	N	N	N	N	N	N	N	N
lt Lithuanian	**tr** Turkish	**bg** Bulgarian	**ro** Romanian	**ru** Russian	**hr** Croatian	**zh** Chinese	**ko** Korean	**ar** Arabic	
N	N	N	N	N	N	N	N	N	

SECTION 1 - Folding Diagram
SECTION 2 - IFU Artwork
SECTION 3 - Czech Language IFU, Reference Text

ARTWORK DESIGN NOTES


Part Number: 10001116-002
Size: 592mm (W) x 210mm (H). **Finished Folded Size:** 148mm (W) x 210mm (H).
IFU Format: Printed Instruction Booklet, 52 pages.


A) For detailed IFU specifications refer to the purchase specification.
B) Artwork Content and Specification Approval via electronic approval.
C) Approval of this document is for final format and layout only. All translations have been previously approved on LLS10001116-001 Rev. 100.

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

SECTION 1 - FOLDING DIAGRAM

Finished Booklet

◌ = Stitch locations.
Internal pages stitched twice to cover.

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PORTEX®

Percutaneous Suctionaid® Tracheostomy Kit

With Guidewire Dilating Forceps – 100/547/070-090
Without Guidewire Dilating Forceps – 100/548/070-090

en INSTRUCTIONS FOR USE
fr MODE D'EMPLOI
de GEBRAUCHSANLEITUNG
it ISTRUZIONI PER L'USO
es INSTRUCCIONES DE USO
pt INSTRUÇÕES DE UTILIZAÇÃO
nl GEBRUIKSAANWIJZING
sv BRUKSANVISNING
da BRUGSANVISNING
no BRUKSANVISNING
fi KÄYTTÖOHJEET
el ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ

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en

Single use.

Do Not Reuse: Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only. Any attempt to re-process the device for subsequent re-use may adversely affect the integrity of the device or lead to deterioration in performance.
Sterile unless unit container is opened or damaged.
Destroy after single use.
Do not resterilize.

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English

en

Portex® Percutaneous Suctionaid® Tracheostomy Kit
INSTRUCTIONS FOR USE

These instructions apply to the use of the following Portex® products:
Catalogue No. 100/547/070-090 With Guidewire Dilating Forceps Internal REF: 100/547/070CZ-090CZ
Catalogue No. 100/548/070-090 Without Guidewire Dilating Forceps Internal REF: 100/548/070CZ-090CZ

It is essential that only the Portex® Guidewire Dilating Forceps are used in the procedure.
For details call your local Smiths Medical distributor or Smiths Medical International direct on +44 (0)1303)260551.

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Introduction

The Percutaneous Tracheostomy technique is a new, minimally invasive procedure which offers an easier, more rapid method of insertion of a tracheostomy tube than conventional open surgical techniques. It also has the advantage in that it can be performed at the bedside.

In a comparative clinical trial involving over 150 patients[1] there was no mortality associated with this new technique and the overall complication rate was 3.9%. The group who had undergone the traditional surgical technique had one fatality and a significantly greater incidence of complications of 18.9% ($p<0.001$).

This new procedure uses the Seldinger technique to guide the specially designed Guidewire Dilating Forceps into the trachea, which then are used to dilate the trachea. The design of the forceps allows the outward force to be distributed evenly along the large surface area of its jaws, thereby minimising the force at any given point of the trachea, and so minimising any risk of tracheal trauma. The guidewire is then finally used to position the tracheostomy tube.

N.B. As the Percutaneous Tracheostomy technique is new, we strongly advise that doctors are fully trained in the procedure and familiar with the Seldinger technique. We therefore recommend that before attempting the technique, the following be carried out:

(1) The instructions for use are read and understood.

(2) The training video is viewed.

(3) The wallchart illustrating the technique is displayed for easy reference.

(4) The technique is observed being performed first hand.

For advice on obtaining the video and the wallchart, please call your local Smiths Medical distributor or Smiths Medical International direct on +44 (0)1303 260551.

A **minimum** of 2 operators is needed – one to maintain the patient's airway, anaesthesia, breathing and circulation, and the other to perform the procedure.

The Suctionaid® tracheostomy tube incorporates a clear suction line terminating above the cuff with a Luer port to facilitate connection to either a syringe or the vacuum control valve included in the pack which allows any pooled secretions to be aspirated.

If vocalisation is required we recommend the use of Portex® Vocalaid® (100/517/060-100).

Kit Contents *[Fig. 1]*

(1) Scalpel
(2) 14G IV Cannula assembly
(3) 10cc Syringe
(4) Guidewire with Introducer
(5) Dilator
(6) Guidewire Dilating Forceps (GWDF)
(7) Tracheostomy tube and obturator with lumen
(8) Suction Control Valve

Cuff Resting Diameter (mm)

Tracheostomy Tube I.D.	7.0	8.0	9.0
Profile® Cuff	22	27	30

Preparation of Components

1 Check the tracheostomy tube cuff for leaks by test inflation. Ensure that the obturator moves freely in the tracheostomy tube, having first released the retaining clips by means of a slight twist. Re-insert the obturator ensuring that it 'clicks' into position. If satisfactory, **fully deflate** the cuff in order to avoid cuff rupture during tube insertion, and clip flange wings back onto the obturator.

2 Check that the guidewire moves freely through the forceps and the tracheostomy tube obturator.

This docu confidential, proprietary and may not be disclosed, disseminated or reproduced without written approval from Smiths Medical.

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Preparation of the Patient

3 The patient should be positioned supine, with a pillow to support the neck and shoulders in order to hyper extend the neck *[Fig. 2]*.

4 Locate the thyroid cartilage between thumb and forefinger and identify and mark the anatomical landmarks *[Fig. 3]*.

(a) Thyroid cartilage

(b) Cricoid cartilage

(c) Tracheal rings

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(e) Possible insertion sites

5 Suggested sites for tube insertion are marked (x) and are either between the first and second or second and third tracheal rings.

6 Prior to the procedure, increase FiO_2 to 100% and monitor the patient using oximetry, capnography, ECG and blood pressure.

7 Suction pharynx and, if present, deflate the tracheal tube cuff. Withdraw the tube to laryngeal inlet and reinflate with the cuff entirely above the vocal cords to re-establish a seal. This is necessary to avoid the tube interfering with the procedure. Hold the tube in position ensuring that the head and neck are in the mid-line and that the airway is being maintained.

N.B. It is recommended that a bronchoscope is used to confirm midline placement of the guidewire and subsequent tracheostomy tube which may prevent paratracheal placement, that can lead to subcutaneous emphysema and/or pneumothorax.

Preparation of the site for surgery

8 Disinfect the skin and apply sterile drapes.

9 Palpate cricoid cartilage and infiltrate local anaesthetic if appropriate. Injection of solutions containing adrenaline may help reduce bleeding at the incision site.

The Procedure

10 Make a horizontal or vertical incision *[Fig. 4]* at the chosen insertion site, of sufficient size to accept the tracheostomy tube (1.5-2 cm). At this stage it may be of benefit to perform an exploratory blunt dissection **(keeping in the midline)** to identify the anatomical landmarks.

N.B. WHERE PRESENT, CONFIRM THAT THE TRACHEAL TUBE HAS BEEN WITHDRAWN AWAY FROM THE INSERTION SITE TO AVOID PUNCTURING THE CUFF.

11 Insert the needle and cannula (with syringe attached) in the midline, into the insertion site *[Fig. 5]* in a caudal direction, to avoid the likelihood of the guidewire being passed into the pharynx at step **13**. Advance the needle (if necessary, stepping off a tracheal ring), until free withdrawal of air into the syringe confirms the entry of the needle and cannula into the trachea. This may be facilitated by using a fluid filled syringe, allowing visualisation of withdrawn air. Some secretions of mucus may be aspirated during this process, which is normal.

12 Leaving the cannula in place, withdraw the needle and syringe. Connect the syringe to the cannula and withdraw air to re-confirm that the cannula is in the trachea. Remove the syringe.

13 Ease the guidewire introducer out from its sheath and straighten the 'J' tip, leaving a sufficient length of exposed guidewire (2-3 cm) to enable its dispensing with the forefinger and thumb *[Fig. 6]*.

Insert the Introducer into the cannula and feed the guidewire out of the sheath until at least 10 cm of its lies in the trachea *[Fig. 7]*, leaving about 30 cm free externally.

14 Confirm free movement of the guidewire in the cannula before removing the cannula, leaving the guidewire in place.

N.B. As the guidewire will be used to guide components into the trachea, it is important to ensure that it is not damaged (e.g. kinked) during the procedure. To help avoid this, always anchor the free end of the wire with forefinger and thumb whilst advancing components into the patient. (In case of damage to the wire at any stage, please refer to warnings at the back of this leaflet). During the procedure, you will need to check that the guidewire is not kinked by ensuring its free movement through the components and the trachea.

15 Pass the dilator over the guidewire through the soft tissues, until resistance is felt on the tracheal wall. With a gentle sideways movement of the dilator *[Fig. 8]*, push the dilator forward to penetrate the anterior tracheal wall, dilating both the tissues and tracheal wall.

N.B. Confirm free movement of the guidewire in the trachea and dilator before removing dilator, leaving guidewire in place.

16 Clamp the Guidewire Dilating Forceps (GWDF) shut and thread the guidewire through the bore at the tip, which then re-appears through the hole on the curved arms *[Fig. 9]*.

17 Hold the free end of the wire and with the wire and forceps tip held at the same angle as the initial insertion of the cannula, advance the forceps until resistance is felt against the anterior tracheal wall. Occasionally, this first pass may result in the forceps entering the trachea. If this has **CLEARLY** occurred, move to step **20**.

18 Gradually open the forceps with both hands, dilating tissues sufficiently to accept the tracheostomy tube. **Withdraw the forceps in the opened position.** *[Fig. 10]*.

N.B. Confirm free movement of the guidewire in the trachea.

19 Repeat stages **16** and **17.** Keeping the forceps closed, advance the forceps through the tracheal wall *[Fig. 11]*. A loss of resistance is felt as the tracheal wall is penetrated.

CAUTION: Before proceeding to the next stage, confirm the successful penetration of the tracheal wall to avoid severe guidewire damage. Repeat* 19 *if necessary. If a loss of resistance during wall penetration was not evident, a pre-tracheal false passage may have been formed** *[Fig. 12]*. ***Continuing the next procedure after this will severely damage the guidewire (see warnings).

20 With the forceps handles along the midline, raise the handles into the vertical position so that the forceps tip penetrates the tracheal wall further and lies longitudinally in the trachea *[Fig. 13]*.

N.B. Confirm free movement of the guidewire in the forceps and trachea.

21 Using 2 hands, open the forceps to dilate the tracheal wall sufficiently, to accept the tracheostomy tube *[Fig. 14]*. **Withdraw the forceps in the open position.**

22 Thread the guidewire through the bore at the tip of the obturator and advance the obturator and tracheostomy tube into the trachea *[Fig. 15]*. Remove the obturator and guidewire leaving the tracheostomy tube in place.

23 Suction trachea and tracheostomy tube to establish a clear airway. Blood which is not suctioned from the airway readily clots and may produce severe airway obstruction.

24 Transfer the breathing system to the tracheostomy tube.

25 Inflate the Soft-Seal® Profile® cuff of the tracheostomy tube with the minimum of air to form a seal. The Soft-Seal® Profile® cuff is designed to seal at very low pressures, to minimise tracheal wall damage.

26 Confirm successful tube placement by currently accepted medical techniques such as chest movement, auscultation, end-tidal CO_2, chest X-ray and SaO_2.

27 Secure the tracheostomy tube with the supplied cotton tapes or a Portex® Tracheostomy Tube Holder (100/503/200).

28 If present, deflate the tracheal tube cuff and remove tube.

29 A tracheostomy dressing may be used at the tracheostomy site for improved patient comfort and hygiene.

30 With the cuff inflated attach the suction device, either syringe or low level suction to the suction line using the suction control valve, and slowly aspirate the secretions up to a maximum suction of 300mm Hg.

31 If resistance is met on aspiration, flush the suction line with either air or sterile saline and aspirate again. If resistance to instillation of air or saline is still met and the instilled fluid cannot be withdrawn, the suction line may be blocked and the tracheostomy tube should be changed or an alternative means of suction instituted.

32 Once suctioning is complete remove the suction device including the control valve and cap the Luer port at the proximal end of the suction tube. The suction control valve is intended for single patient use and can be cleaned using sterile saline solution. Replacement suction control valves are also available (100/515/200).

33 If using continuous low level suction, monitor the suction pressure regularly. If the suction pressure increases disconnect the suction and follow the above procedure for resistance to suction.

Absolute Contra-Indications

1 Emergency surgical airway management.
(In emergency, cricothyrotomy is the procedure of choice).

2 Children.

3 Pre-existing infection at the tracheostomy site.

4 Pre-existing malignancy at the tracheostomy site.

5 Uncertainty in identifying the anatomical landmarks.

Relative Contra-Indications

Risks in these situations must be balanced against the benefits of the procedure.

1 Enlarged thyroid gland.

2 Previous surgery at the tracheostomy site (e.g. thyroidectomy).

3 Bleeding diathesis (e.g. due to anticoagulant therapy).

Precautions

1 The security of all breathing system connectors should be checked when the circuit is established and frequently thereafter. Disconnection may be facilitated by the use of a disconnection wedge (100/555/100).

2 Patients should be adequately humidified to minimise encrustation of the tracheostomy tube lumen and prevent tracheal mucosal damage.

3 The patency of the tracheostomy tube lumen must be assured by regular suctioning. Check routinely and replace as required to maintain a patent airway. Recommend period of use 30 days.

4 Cuff pressure should be monitored and adjusted routinely.

5 Devices used in or during inflation of the cuff must be clean and free from all foreign matter. The inflation device should be removed from the inflation valve immediately after use, and the dust cap fitted.

6 Guard against cuff damage by avoiding contact with sharp edges.

7 The inflation line valve may interfere with Magnetic Resonance Imaging (MRI) picture clarity. Ensure the valve is positioned away from the area being scanned.

8 If Portex® tracheostomy tubes are used outside the hospital, the patient must be instructed by a health care professional in the safe use and handling of this product.

9 Do not use solution other than sterile saline to clean any part of this tracheostomy tube.

10 When suctioning via the suction line has been completed the suction device including the control valve should be removed, and the suction line capped to minimise the risk of infection.

11 Dispose of per local procedure.

This docu confidential, proprietary and may not be disclosed, disseminated or reproduced without written approval from Smiths Medical.

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WARNINGS

1 **If bleeding occurs, direct pressure will control most cases. If bleeding remains uncontrolled, the wound must be explored using an open surgical technique and the bleeding stopped.**

2 **If the guidewire is damaged during the procedure, it may be difficult to continue. In these cases it may be possible to proceed only if:**

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If the guidewire cannot be salvaged, it will be necessary to use a new kit or a new guidewire *(see below for details).*

3 **Prior to removal of cuffed tracheostomy tubes all air must be completely removed from the cuff to prevent damage to the trachea and stoma.**

4 **If the tracheostomy tube is lubricated prior to insertion ensure the lubricant does not occlude the lumen of the tube preventing ventilation of the patient.**

5 **Do not inflate the cuff with a measured volume of air or by feel of pressure from the syringe, since little resistance should be felt during cuff inflation.**

6 **During anaesthesia nitrous oxide may diffuse into the cuff causing an increase or decrease in cuff pressure.**

7 **Contact with electrosurgery electrodes or laser surgery beams must be avoided because PVC will produce toxic fumes in air or ignite in an enriched oxygen environment (e.g. anaesthesia).**

8 **Do not use this tracheostomy tube in patients with abnormal upper airway anatomy or pathology as this may result in partial or total airway obstruction.**

9 **Never insert the guidewire into the Tracheostomy tube with the red J-tip straightener in place.**

10 **This product contains DEHP – data derived from animal and epidemiological studies outline the following risks, but are not limited to: Adverse effects on fetal development, kidney and liver toxicity, negative impact on fertility. Special attention on the application of phthalate products should be given for the following patient populations: pediatrics, pregnant and nursing women. For further information on phthalates please refer to the Smiths Medical website: www.smiths-medical.com/phthalates**

FOR FRANCE ONLY: As of July 1st 2015 the Use of products containing DEHP is banned in paediatric-, neonatal- and maternity departments. Please refer to French Law "LOI no 2012-1442 du 24 décembre 2012".

Further Information

The Guidewire Dilating Forceps are subjected to high stress during the procedure and may be damaged in use as evidenced by the forceps arms having been bent or the tips not meeting correctly.

If this is the case, the forceps must be discarded.

For further information please call your local Smiths Medical distributor or Smiths Medical International direct on: +44 (0)1303 260551.

Ordering Details

Product	Code	Boxed in
Portex® Percutaneous Suctionaid® Tracheostomy Kit with Guidewire Dilating Forceps	100/547/070-090	1
Portex® Percutaneous Suctionaid® Tracheostomy Kit without Guidewire Dilating Forceps	100/548/070-090	1
Spare Guidewires	100/544/000	5
Suction Control Valves	100/515/200	10

References

[1] GRIGGS WM, MYBURGH JA, WORTHLEY LIG. A prospective comparison of a percutaneous tracheostomy technique with standard surgical tracheostomy. *Intensive care Med* (1991); 17:261-3.

Profile, Soft-Seal, Vocalaid, Suctonaid, Portex, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2019 Smiths Medical . All rights reserved.

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
IFU ARTWORK

(Not to scale)

smiths medical bringing technology to life	WORK INSTRUCTION FM-DVSOP2002-01 Rev. 001 Eff. Date: 25-OCT-2016 CO No.: CO-10029546	Doc. Number AS10007787-005	Revision 100
WI Type: ARTWORK SPECIFICATION	Effective Date: 06-NOV-2019	Page 29 of 29	
Title: INSTRUCTIONS FOR USE - PORTEX® BLUGRIGGS® PERCUTANEOUS DILATION TRACHEOSTOMY PROCEDURAL KIT OR TRAY WITH BLUSELECT® TRACHEOSTOMY TUBE WITH OR WITHOUT FORCEPS (10007787-005)			

This document and all information contained herein is confidential, proprietary and may not be disclosed, disseminated or reproduced without written approval from Smiths Medical.
- INTERNAL -

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Portex® Percutaneous Tracheostomy Kit

INSTRUCTIONS FOR USE

These instructions apply to the use of the following Portex® products:
Catalogue No. 100/540/070-090 With Guidewire Dilating Forceps Internal REF: 100/540/070CZ-090CZ
Catalogue No. 100/542/070-090 Without Guidewire Dilating Forceps Internal REF: 100/542/070CZ-090CZ

It is essential that only the Portex® Guidewire Dilating Forceps are used in the procedure. For details call your local Smiths Medical distributor or Smiths Medical International direct on +44 (0)1303)260551.

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A minimum of 2 operators is needed – one to maintain the patient's airway, anaesthesia, breathing and circulation, and the other to perform the procedure.

Kit Contents *(Fig. 1)*

(1) Scalpel
(2) 14G IV Cannula assembly
(3) 10cc Syringe
(4) Guidewire with Introducer
(5) Dilator
(6) Guidewire Dilating Forceps (GWDF)
(7) Tracheostomy tube and obturator with lumen

Cuff Resting Diameter (mm)

Tracheostomy Tube I.D.	7.0	8.0	9.0
Traditional Cuff	14	16	18
Profile® Cuff	22	27	30

Preparation of Components

1 Check the tracheostomy tube cuff for leaks by test inflation. Ensure that the obturator moves freely in the tracheostomy tube. If satisfactory, **fully deflate** the cuff in order to avoid cuff rupture during tube insertion, and clip flange wings back onto the obturator.

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(c) Tracheal rings

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(e) Possible insertion sites

5 Suggested sites for tube insertion are marked (e) and are either between the first and second or second and third tracheal rings.

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The Procedure

10 Make a horizontal or vertical incision *[Fig. 4]* at the chosen insertion site, of sufficient size to accept the tracheostomy tube (1.5–2 cm). At this stage it may be of benefit to perform an exploratory blunt dissection **(keeping in the midline)** to identify the anatomical landmarks.

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11 Insert the needle and cannula (with syringe attached) in the midline, into the insertion site *[Fig. 5]* in a caudal direction, to avoid the likelihood of the guidewire being passed into the pharynx at step **13**. Advance the needle (if necessary, clamping off a tracheal ring), until free withdrawal of air into the syringe confirms the entry of the needle and cannula into the trachea. This may be facilitated by using a fluid filled syringe, allowing visualisation of withdrawn air. Some secretions of mucus may be aspirated during this process, which is normal.

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Insert the introducer into the cannula and feed the guidewire out of the sheath until at least 10 cm of its lies in the trachea *[Fig. 7]*, leaving about 30 cm free externally.

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N.B. As the guidewire will be used to guide components into the trachea, it is important to ensure that it is not damaged (e.g. kinked) during the procedure. To help avoid this, always anchor the free end of the wire with forefinger and thumb whilst advancing components into the patient. (In case of damage to the wire at any stage, please refer to warnings at the back of this leaflet). During the procedure, you will need to check that the guidewire is not kinked by ensuring is free movement through the components and the trachea.

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19 Repeat stages **16** and **17**, keeping the forceps closed, advance the forceps through the tracheal wall *[Fig. 11]*. A loss of resistance is felt as the tracheal wall is penetrated.

CAUTION: Before proceeding to the next stage, confirm the successful penetration of the tracheal wall to avoid severe guidewire damage. Repeat 19 if necessary. If a loss of resistance during wall penetration was not evident, a pre-tracheal false passage may have been formed** [Fig. 12]. **Continuing the next procedure after this will severely damage the guidewire (see warnings).

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22 Thread the guidewire through the bore at the tip of the obturator and advance the obturator and tracheostomy tube into the trachea *[Fig. 15]*. Remove the obturator and guidewire, leaving the tracheostomy tube in place.

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25 Inflate the Profile* cuff of the tracheostomy tube with the minimum of air to form a seal. The Profile* cuff is designed to seal at very low pressures, to minimise tracheal wall damage.

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27 Secure the tracheostomy tube with the supplied cotton tapes or a Portex* Tracheostomy Tube Holder (100/503/200).

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8 If Portex* tracheostomy tubes are used outside the hospital, the patient must be instructed by a health care professional in the safe use and handling of this product.

9 Do not use solutions other than sterile saline to clean any part of this tracheostomy tube.

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WARNINGS

1 If bleeding occurs, direct pressure will control most cases. If bleeding remains uncontrolled, the wound must be explored using an open surgical technique and the bleeding stopped.

2 If the guidewire is damaged during the procedure, it may be difficult to continue. In these cases it may be possible to proceed only if:

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If the guidewire cannot be salvaged, it will be necessary to use a new kit or a new guidewire *(see below for details)*.

3 Prior to removal of cuffed tracheostomy tubes all air must be completely removed from the cuff to prevent damage to the trachea and stoma.

4 If the tracheostomy tube is lubricated prior to insertion ensure the lubricant does not occlude the lumen of the tube preventing ventilation of the patient.

5 Do not inflate the cuff with a measured volume of air or by feel of pressure from the syringe, since little resistance should be felt during cuff inflation.

6 During anaesthesia nitrous oxide may diffuse into the cuff causing an increase or decrease in cuff pressure.

7 Contact with electrosurgery electrodes or laser surgery beams must be avoided because PVC will produce toxic fumes in air or ignite in an enriched oxygen environment (e.g. anaesthesia).

8 Do not use this tracheostomy tube in patients with abnormal upper airway anatomy or pathology as this may result in partial or total airway obstruction.

9 Never insert the guidewire into the Tracheostomy tube with the red J-tip straightener in place.

10 This product contains DEHP – data derived from animal and epidemiological studies outline the following risks, but are not limited to: Adverse effects on foetal development, kidney and liver toxicity, negative impact on fertility. Special attention on the application of phthalate products should be given for the following patient populations: pediatrics, pregnant and nursing women. For further information on phthalates please refer to the Smiths Medical website: www.smiths-medical.com/phthalates

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Further Information

The Guidewire Dilating Forceps are subjected to high stress during the procedure and may be damaged in use as evidenced by the forceps arms having been bent or the tips not meeting correctly. If this is the case, the forceps must be discarded.

For further information please call your local Smiths Medical distributor or Smiths Medical International direct on: +44 (0)1303 266551.

Ordering Details

Product	Code	Boxed in
Portex® Percutaneous Tracheostomy Kit with Guidewire Dilating Forceps	100/540/070-090CZ	1
Portex® Percutaneous Tracheostomy Kit without Guidewire Dilating Forceps	100/542/070-090CZ	1
Spare Guidewires	100/544/000CZ	5

References

[1] GRIGGS WM, MYBURGH JA, WORTHLEY LIG. A prospective comparison of a percutaneous tracheostomy technique with standard surgical tracheostomy. Intensive care Med (1991); 17:261-3.

Profile, Portex, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.
© 2019 Smiths Medical. All rights reserved.

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