VITEK 2™ — technology

AST-GN90

IVD

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S. pneumoniae Susceptibility

SUMMARY AND EXPLANATION

The VITEK® 2 Gram Positive Susceptibility Card is intended for use with the VITEK 2 Systems in clinical laboratories as an in vitro test to determine the susceptibility of ***Streptococcus pneumoniae*** to antimicrobial agents when used as instructed in the Product Information manual.

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**SP(P&M) = *S. pneumoniae* (pneumonia and meningitis)
**Spn except P(R) = *S. pneumoniae* (excluding penicillin-resistant strains)
[c] = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.

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Numerical values are expressed in µg/ml.
[c] = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.

LIMITATIONS

The ability of the AST card to detect resistance with the following combination(s) is unknown because resistant strains were not available at the time of comparative testing:
- Amoxicillin: *Streptococcus pneumoniae*
- Benzylpenicillin: *S. pneumoniae* (based on the breakpoints for pneumonia)
- Linezolid: *Streptococcus pneumoniae*
- Telithromycin: *Streptococcus pneumoniae*

The ability of the VITEK 2 System to detect **intermediate isolates of *Streptococcus pneumoniae*** with **erythromycin** is unknown because intermediate organisms were not available at the time of comparative testing. Erythromycin-intermediate strains of *Streptococcus pneumoniae* are very rare.

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⇩ MUST enter the following barcodes into "Flex Panel Entry" program before first use of this Susceptibility Card.

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Staphylococcus spp., *Enterococcus* spp. and *S. agalactiae* Susceptibility

SUMMARY AND EXPLANATION

The VITEK® 2 Gram Positive Susceptibility Card is intended for use with the VITEK 2 Systems in clinical laboratories as an in vitro test to determine the susceptibility of ***Staphylococcus* spp., *Enterococcus* spp., and *S. agalactiae*** to antimicrobial agents when used as instructed in the Product Information manual.

INSTRUCTIONS FOR USE

A package insert is provided in the kit or downloadable from **www.biomerieux.com/techlib**

See the Product Information manual for additional Instructions for Use.

STORAGE CONDITIONS

Store at 2° to 8° C.

CONTENTS OF THE CARD

Antimicrobic	Code	Concentration §	Calling Range ≤	Calling Range ≥	FDA Indications for Use
Benzylpenicillin	P	0.125, 0.25, 1	0.03	0.5	*Staphylococcus* spp.
Cefoxitin Screen	OXSF	6	NEG	POS	*Staphylococcus* spp.
Clindamycin	CM	0.06, 0.25, 1	0.125	4	MSSA**, MSSE**
Daptomycin [NS]	DAP	0.5, 1, 2, 4, 16	0.12	8	*S. aureus*, VSEfaeca**
Erythromycin	E❷	1, 2, 4, 8	0.25	8	
Fosfomycin	FOS	8, 32	8	128	N/A**
Fusidic Acid	FA	0.5, 1, 4	0.5	32	N/A**
Gentamicin	GM	8, 16, 64	0.5	16	*Staphylococcus* spp.
Inducible Clindamycin Resistance	ICR❷	CM 0.5, CM/E 0.25/0.5	NEG	POS	*Staphylococcus* spp.
Levofloxacin	LEV	0.25, 2, 8	0.12	8	*Staphylococcus* spp., *Enterococcus* spp., *S. agalactiae*
Linezolid	LNZ	0.5, 1, 2	0.5	8	*S. agalactiae*, *E. faecalis*, *E. faecium*, *S. aureus*, *S. epidermidis*, *S. haemolyticus*
Mupirocin	MUP❷	2, 4	2	512	
Oxacillin	OX1	0.5, 1, 2	0.25	4	*Staphylococcus* spp.
Rifampicin	RA	0.015, 0.03, 0.125, 0.5	0.03	4	N/A**
Teicoplanin	TEC	1, 4, 8, 16	0.5	32	N/A**
Tetracycline	TE	0.5, 1, 2	1	16	*Staphylococcus* spp., *Enterococcus* spp., *S. agalactiae*
Tigecycline [NS]	TGC	0.25, 0.5, 1	0.12	2	*E. faecalis*, *E. faecium*, *S. aureus*, *S. epidermidis*, *S. haemolyticus*, *S. agalactiae*, *E. casseliflavus*
Trimethoprim/Sulfamethoxazole	SXT❷	8/152, 16/304, 32/608	10 (0.5/9.5)	320 (16/304)	N/A**
Vancomycin	VA❷	1, 2, 4, 8, 16	0.5	32	*Enterococcus* spp., *Staphylococcus* spp., *S. agalactiae*

Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
NEG = Negative
POS = Positive
**MSSA = Methicillin-susceptible *S. aureus*
**MSSE = Methicillin-susceptible *S. epidermidis*
[NS] = The current absence of resistant isolates precludes defining any results other than susceptible. Isolates yielding MIC results suggestive of Nonsusceptible category should be submitted to a reference laboratory for further testing.
**VSEfaeca = *E. faecalis* (vancomycin susceptible strains)

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❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**N/A = No specific FDA Indications for Use available

QUALITY CONTROL

		CLSI® Quality Control Organisms VITEK 2 Results				
Antimicrobic	Code	*E. faecalis* ATCC® 29212™	*S. aureus* ATCC® 29213™	*S. aureus* ATCC® BAA-1026™	*S. aureus* ATCC® BAA-976™	*S. aureus* ATCC® BAA-977™
Benzylpenicillin	P	-	0.25 – ≥0.5	-	-	-
Cefoxitin Screen	OXSF	-	NEG	POS	-	-
Clindamycin	CM	≥4	≤0.12 – 0.25	-	-	-
Daptomycin[NS]	DAP	1 – 4	0.25 – 1	-	-	-
Erythromycin	E	1 – 4	≤0.25 – 1	-	-	-
Fosfomycin	FOS	-	≤8	-	-	-
Fusidic Acid	FA	-	≤0.5	-	-	-
Gentamicin	GM	-	≤0.5 – 1	-	-	-
Inducible Clindamycin Resistance	ICR	-	-	-	NEG	POS
Levofloxacin	LEV	0.25 – 2	≤0.12 – 0.5	-	-	-
Linezolid	LNZ	1 – 4	1 – 4	-	-	-
Mupirocin	MUP	-	≤2	-	-	-
Oxacillin	OX1	-	≤0.25 – 0.5	-	-	-
Rifampicin	RA	0.5 – ≥4	≤0.03	-	-	-
Teicoplanin	TEC	≤0.5	≤0.5 – 1	-	-	-
Tetracycline	TE	8 – ≥16	≤1	-	-	-
Tigecycline[NS]	TGC	≤0.12	≤0.12 – 0.25	-	-	-
Trimethoprim/Sulfamethoxazole	SXT	≤10 (0.5/9.5)	≤10 (0.5/9.5)	-	-	-
Vancomycin	VA	1 – 4	≤0.5 – 2	-	-	-

Numerical values are expressed in µg/ml.
NEG = Negative
POS = Positive
[NS] = The current absence of resistant isolates precludes defining any results other than susceptible. Isolates yielding MIC results suggestive of Nonsusceptible category should be submitted to a reference laboratory for further testing.

LIMITATIONS

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):
- Daptomycin: *Streptococcus agalactiae*
- Rifampicin: *Streptococcus agalactiae*

Perform an alternative method of testing prior to reporting of results when a Positive (+) result is obtained with the following antibiotic/organism combination(s):
- Cefoxitin Screen: *Staphylococcus saprophyticus*

The ability of the AST card to detect resistance with the following combination(s) is unknown because resistant strains were not available at the time of comparative testing:
- Linezolid: *Enterococcus* spp., *Staphylococcus* spp., *Streptococcus agalactiae*
- Tigecycline: *Enterococcus* spp., *Staphylococcus* spp., *Streptococcus* spp.

NOTE: A result for an antibiotic/organism combination which may have a limitation may be suppressed from reporting. Refer to the software user manual for instructions.

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⇩ MUST enter the following barcodes into “Flex Panel Entry” program before first use of this Susceptibility Card.

01 A320JT0U---K

02 BAST-P63201R

03 C---ZY20G37I

04 D473X300T0UY

05 E193I132N3Q8

06 F20480Y1H3KL

07 G0K3L-V-U-TD

08 H-P-K-J-I-65

INDEX OF SYMBOLS

Symbol	Meaning
REF	Catalogue number
IVD	In Vitro Diagnostic Medical Device
	Manufacturer
	Date of Manufacture
	Temperature limitation
	Use by
LOT	Batch code
	Consult Instructions for Use
Σ	Contains sufficient for <n> tests
EC REP	Authorized representative in the European Community

PATENTS

Product covered by one or more of U.S. Patent Nos. D414,272; D437,797; 5,609,828; 5,746,980; 5,804,437; 5,869,005; 5,932,177; 5,951,952; 6,267,929; 6,309,890 and 6,340,573; and Foreign Counterparts. Other Patents Pending.

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABLITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

BIOMERIEUX, the blue logo, VITEK and bioLiaison are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries, or one of its companies. The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection. Any other name or trademark is the property of its respective owner.

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli
**Kox(-ESBL) = *K. oxytoca* (excluding ESBL producing strains)
Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid. For ESBL positive strains, the test interpretation should be reported as resistant for all penicillins, cephalosporins, and aztreonam.
NEG = Negative
POS = Positive
A negative ESBL test result does not rule out the presence of an ESBL masked by an AmpC beta-lactamase.
**N/A = No specific FDA Indications for Use available
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

QUALITY CONTROL

		CLSI® Quality Control Organisms VITEK 2 Results			
Antimicrobic	Code	*E. coli* ATCC® 25922™	*P. aeruginosa* ATCC® 27853™	*E. coli* ATCC® 35218™	*K. pneumoniae* ssp *pneumoniae* ATCC® 700603™
Amikacin	AN	≤2 – 4	≤2 – 4	-	-
Amoxicillin/Clavulanic Acid	AMC	≤2/1 – 8/4	-	4/2 – 16/8	-
Ampicillin	AM	≤2 – 8	-	-	-
Cefepime	FEP	≤1	≤1 – 8	-	-
Cefotaxime	CTX	≤1	8 – 32	-	-
Ceftazidime	CAZ	≤1	≤1 – 4	-	-
Ciprofloxacin	CIP	≤0.25	≤0.25 – 1	-	-
Ertapenem	ETP	≤0.5	2 – ≥8	-	-
ESBL'	ESB	NEG	-	-	POS
Fosfomycin	FOS	≤16	≤16	-	-
Gentamicin	GM	≤1	≤1 – 2	-	-
Imipenem	IPM	≤0.25	1 – 4	-	-
Meropenem	MEM	≤0.25	≤0.25 – 1	-	-
Nitrofurantoin	FT	≤16	-	-	-
Norfloxacin	NOR	≤0.5	1 – 4	-	-
Piperacillin/Tazobactam	TZP	≤4/4	≤4/4 – 8/4	≤4/4	-
Trimethoprim/Sulfamethoxazole	SXT	≤20 (1/19)	160 (8/152) – ≥320 (16/304)	-	-

Numerical values are expressed in µg/ml.
Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid. For ESBL positive strains, the test interpretation should be reported as resistant for all penicillins, cephalosporins, and aztreonam.
NEG = Negative
POS = Positive
A negative ESBL test result does not rule out the presence of an ESBL masked by an AmpC beta-lactamase.

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℞ only IVD

INTENDED USE

These Instructions for Use correspond to the VITEK® 2 Systems 7.01 or higher software. If you are not using VITEK® 2 Systems 7.01 or higher software, please refer to the VITEK® 2 Systems Product Information that you received with your current software version.

The VITEK® 2 Gram-Positive identification card (GP) is intended for use with VITEK® 2 Systems for the automated identification of most significant Gram-positive organisms. The VITEK® 2 GP identification card is a single-use disposable. For a list of claimed species, see the Organisms Identified section.

DESCRIPTION

The GP identification card is based on established biochemical methods[2,3,7,8,9,10,11,14,20,21,22,23,27,32,37,39] and newly developed substrates. There are 43 biochemical tests measuring carbon source utilization, enzymatic activities and resistance. Final identification results are available in approximately eight hours or less.

For a list of well contents, see the GP Well Contents table.

Table 1: GP Well Contents

Well	Test	Mnemonic	Amount/Well
2	D-AMYGDALIN	AMY	0.1875 mg
4	PHOSPHATIDYLINOSITOL PHOSPHOLIPASE C	PIPLC	0.015 mg
5	D-XYLOSE	dXYL	0.3 mg
8	ARGININE DIHYDROLASE 1	ADH1	0.111 mg
9	BETA-GALACTOSIDASE	BGAL	0.036 mg
11	ALPHA-GLUCOSIDASE	AGLU	0.036 mg
13	Ala-Phe-Pro ARYLAMIDASE	APPA	0.0384 mg
14	CYCLODEXTRIN	CDEX	0.3 mg
15	L-Aspartate ARYLAMIDASE	AspA	0.024 mg
16	BETA GALACTOPYRANOSIDASE	BGAR	0.00204 mg
17	ALPHA-MANNOSIDASE	AMAN	0.036 mg
19	PHOSPHATASE	PHOS	0.0504 mg
20	Leucine ARYLAMIDASE	LeuA	0.0234 mg
23	L-Proline ARYLAMIDASE	ProA	0.0234 mg
24	BETA GLUCURONIDASE	BGURr	0.0018 mg
25	ALPHA-GALACTOSIDASE	AGAL	0.036 mg
26	L-Pyrrolydonyl-ARYLAMIDASE	PyrA	0.018 mg
27	BETA-GLUCURONIDASE	BGUR	0.0378 mg
28	Alanine ARYLAMIDASE	AlaA	0.0216 mg
29	Tyrosine ARYLAMIDASE	TyrA	0.0276 mg
30	D-SORBITOL	dSOR	0.1875 mg
31	UREASE	URE	0.15 mg
32	POLYMIXIN B RESISTANCE	POLYB	0.00093 mg
37	D-GALACTOSE	dGAL	0.3 mg

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Note: Other well numbers between 1 and 64 not designated in this table are empty.

PRECAUTIONS

Note: For industry customers that need assistance on selecting the correct VITEK® 2 identification card, please refer to the VITEK® 2 Compact Instrument User Manual chapter, "Guidance to Select a VITEK® 2 Identification Card."

- For *In Vitro* Diagnostic Use Only.
- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For professional use only.
- Suspensions not within the appropriate zone on the VITEK® 2 DENSICHEK™ Plus or the VITEK® 2 DENSICHEK™ may compromise card performance.
- Do not use the card after the expiration date shown on the package liner.
- Store the card unopened in the package liner. Do not use the card if the protective package liner is damaged or if no desiccant is present.
- Allow the card to come to room temperature before opening the package liner.
- Do not use powdered gloves. Powder may interfere with the optics.
- Use of culture media other than the recommended types must be validated by the customer laboratory for acceptable performance.
- A Gram stain should be performed to determine an organism’s Gram reaction and morphology prior to selecting the identification card to inoculate.
- The card performs as intended only when used in conjunction with VITEK® 2 Systems, following the instructions contained in these Instructions for Use.
- **Do not use glass test tubes**. Use clear plastic (polystyrene) tubes only. Variation exists among test tubes of standard diameter. Carefully place the tube into the cassette. If resistance is encountered, discard and try another tube that does not require pressure to insert.
- Prior to inoculation, inspect cards for tape tears or damage to the tape and discard any that are suspect. Check the saline level in the tubes after the cassette has been processed to ensure proper filling of card.
 - VITEK® 2 60 or VITEK® 2 XL: Eject improperly filled cards.
 - VITEK® 2 Compact: Do not load improperly filled cards.
- Give special consideration to specimen source and patient drug or antimicrobic regimen.

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- Interpretation of test results requires the judgment and skill of a person knowledgeable in microbial identification testing. Additional testing may be required. (See the Supplemental Tests section.)
- Do not clean saline dispenser with chemical agents. The use of chemical agents may impact card performance.

Warning: All patient specimens, microbial cultures, and inoculated VITEK® 2 cards, along with associated materials, are potentially infectious and should be treated with universal precautions.[30,35]

Warning: All hazardous waste must be disposed of by following your local inspecting agency's guidelines.

STORAGE CONDITIONS

Upon receipt, store VITEK® 2 GP cards unopened in their original package liner at 2°C to 8°C.

SPECIMEN PREPARATION

For specimen preparation information, see the Culture Requirements Table.

Table 2: Culture Requirements Table

VITEK® 2 Card	Media	Age of Culture[1]	Incubation Conditions	Inoculum Density	Dilution for AST	Age of Suspension Before Loading Instrument
GP	TSAB[2,3] CBA[2,3] TSA[2,3] BP CHBA CHOC CHOC PVX CNT CPS ID MRSA ID MSA SAID TSAHB TSAL VRE	12 to 48 hours	35°C to 37°C 5% to 10% CO_2 or aerobic, non CO_2	0.50 to 0.63 McFarland Standard	N/A[4]	≤30 minutes
GP and AST GP pair	TSAB CBA CPS ID	18 to 24 hours	35°C to 37°C 5% to 10% CO_2 or aerobic, non-CO_2	0.50 to 0.63 McFarland Standard	280 µL in 3.0 mL saline	< 30 minutes
GP and AST ST pair	TSAB CBA	18 to 24 hours	35°C to 37°C 5% to 10% CO_2	0.50 to 0.63 McFarland Standard	280 µL in 3.0 mL saline	< 30 minutes

[1]Cultures with scant or poor growth may give unidentified or incorrect results even when the Age of Culture requirements are met.

[2]These media were used in the identification product database development and will give optimal performance.

[3]OMA Official Methods of Analysis validated medium.

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Culture Requirements Table — Media Abbreviations

BP = Baird Parker

CBA = Columbia Blood Agar with 5% Sheep Blood

CHBA = Columbia Horse Blood Agar

CHOC = Chocolate Agar

CHOC PVX = Chocolate Polyvitex

CNT = Count-TACT® (irradiated) Trypticase Soy Agar

CPS ID = chromID™ CPS (CPS ID agar)

MRSA ID = chromID™ (MRSA ID Agar)

MSA = Mannitol Salt Agar

SAID = chromID™ S. aureus (S. aureus ID Agar)

TSA = Trypticase Soy Agar

TSAB = Trypticase Soy Agar with 5% Sheep Blood

TSAHB = Trypticase Soy Agar with 5% Horse Blood

TSAL = TSA with Lecithin and P80

VRE = chromID™ VRE

TEST PROCEDURE

Materials

When used with VITEK® 2 instrumentation, the GP card is a complete system for routine identification testing of most clinically significant Gram-positive organisms.

Required materials are:

- VITEK® 2 GP Card
- DENSICHEK™ Plus Kit or VITEK® DENSICHEK® Kit
- DENSICHEK™ Plus Standards Kit or DENSICHEK® Standards Kit
- VITEK® 2 Cassette
- Sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- 12 mm x 75 mm clear plastic (polystyrene) disposable test tubes
- Sterile sticks or swabs
- Appropriate agar medium (see Culture Requirements table).

Optional accessories:

- Adjustable volume saline dispenser
- Loops
- Pre-dispensed saline test tubes (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- Test tube caps
- Vortex

Procedure

Warning: Failure to follow instructions and recommendations provided in this section for performing laboratory tasks may cause erroneous or delayed results.

For product-specific information, see the Culture Requirements table.

Note: Prepare the inoculum from a pure culture, according to good laboratory practices. In case of mixed cultures, a re-isolation step is required. It is recommended that a purity check plate be done to ensure that a pure culture was used for testing.

1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.

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2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).
3. Use a sterile stick or swab to transfer a sufficient number of morphologically similar colonies to the saline tube prepared in step 2. Prepare a homogenous organism suspension with a density equivalent to a McFarland No. 0.50 to 0.63 using a calibrated VITEK® 2 DENSICHEK™ Plus or VITEK® 2 DENSICHEK™.
 Note: Age of suspension must not exceed 30 minutes before inoculating card.
4. Place the suspension tube and GP card in the cassette.
5. Refer to the appropriate Instrument User Manual for instructions on data entry and how to load the cassette into the instrument.
6. Follow your local inspecting agency's guidelines for disposal of hazardous waste.

RESULTS

Identification Analytical Techniques

VITEK® 2 Systems identify an organism by using a methodology based on the characteristics of the data and knowledge about the organism and reactions being analyzed. Sufficient data have been collected from known strains to estimate the typical reactions of the claimed species to a set of discriminating biochemicals. If a unique identification pattern is not recognized, a list of possible organisms is given, or the strain is determined to be outside the scope of the database.

The printed lab report contains suggestions for any supplemental tests necessary to complete the identification. If the tests are not sufficient to complete the identification, then standard microbiology references and literature should be consulted.

Certain species may belong to slashline (mixed) taxa identification. This occurs when the biopattern is the same for the taxa listed. Supplemental tests may be used to separate slashline taxa. The species in the GP Slashline Taxa table belong to the GP slashline taxa.

Table 3: GP Slashline Taxa

Slashline Name	Species Belonging to the Slashline
For 7.01 Software Users	
Micrococcus luteus/lylae	*Micrococcus luteus* *Micrococcus lylae*
For 7.01 or Higher Software Users	
Dermacoccus nishinomiyaensis/ *Kytococcus sedentarius*	*Dermacoccus nishinomiyaensis* *Kytococcus sedentarius*
Listeria ivanovii	*Listeria ivanovii* ssp. *ivanovii* *Listeria ivanovii* ssp. *londoniensis*
Streptococcus mitis/ *Streptococcus oralis*	*Streptococcus mitis* *Streptococcus oralis*

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Table 4: Identification Card Qualifying Messages

ID Message Confidence Level	Choices	% Probability	Comments
Excellent	1	96 to 99	N/A
Very Good	1	93 to 95	N/A
Good	1	89 to 92	N/A
Acceptable	1	85 to 88	N/A
Low Discrimination	2 to 3	Sum of choices = 100; after resolution to one choice, percent probability reflects the number associated with selected choice.	Two to three taxa exhibit same biopattern. Separate by supplemental testing.
Inconclusive or Unidentified Organism	> 3 or 0	N/A	Either > 3 taxa exhibit same biopattern or Very atypical biopattern. Does not correspond to any taxon in the database. Check Gram stain and purity.

PERCENT PROBABILITY

As part of the identification process, the software compares the test set of reactions to the expected set of reactions of each organism, or organism group, that can be identified by the product. A quantitative value, the percent probability, is calculated and relates to how well the observed reactions compare to the typical reactions of each organism. A perfect match between the test reaction pattern and the unique reaction pattern of a single organism, or organism group, would provide a percent probability of 99. When a perfect match is not obtained, it is still possible for the reaction pattern to be sufficiently close to that of an expected reaction pattern such that a clear decision can be provided about the organism identification. The range of percent probabilities in the one-choice case is 85 to 99. Values closer to 99 indicate a closer match to the typical pattern for the given organism.

When the reaction pattern is not sufficient to discriminate between two to three organisms, the percent probabilities reflect this ambiguity. The reported probability values indicate, relatively, the order in which the reaction pattern best corresponds to the listed possibilities. The order does not, however, suggest that the pattern match to one of the possible identifications is clearly superior to another. The probability characteristic of an overall sum of 100 is retained through the calculation process. After resolution to one choice, the probability characteristic of the single choice is retained.

ADDITIONAL INFORMATION ON LAB REPORT

Supplemental test — External (offline) test that allows the user to resolve a slashline or Low Discrimination identification. Numbers in parentheses indicate percent positive reaction for the species/test listed.

Contraindicating test — Test result that is unusual for a reported taxon.

Table 5: Notes Associated with Certain Taxa

Taxa	Note
For 7.01 or Higher Software Users	
Enterococcus durans	Possibility of *Enterococcus villorum* if veterinary.
Listeria monocytogenes	Critical pathogen, check CAMP test and beta hemolysis. The species identified may have significance to patient or sample outcome and can be stopped for review.
Staphylococcus warneri	Possibility of *Staphylococcus pasteuri* if yellow pigmented.
For 8.01 or Higher Software Users	

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Taxa	Note
Listeria innocua	Possibility of *Listeria monocytogenes*. Check for beta hemolysis. *Listeria innocua* strains are non-hemolytic.

Notes Associated with an Improperly Filled Card or with a Negative Profile (Biopattern)

- For the case where the time between two readings is greater than 40 minutes: “CARD ERROR — Missing data.”
- For the case where there is a negative profile: “Organism with low reactivity biopattern — please check viability.”
- When a biopattern is calculated for an unknown organism that is completely negative or consists of both negative tests and tests that fall within the uncertainty zone, the identification call will be “Non or low reactive biopattern.”

The following non-reactive species could potentially trigger this note if a test was atypical or fell within the uncertainty zone:

- *Alloiococcus otitis*
- *Dermacoccus nishinomiyaensis*
- *Gemella bergeri*
- *Kocuria rosea*
- *Kocuria varians*
- *Kytococcus sedentarius*
- *Leuconostoc mesenteroides* ssp. *cremoris*
- *Micrococcus lylae*
- *Staphylococcus auricularis*
- *Streptococcus pluranimalium*

QUALITY CONTROL

Quality control organisms and their expected results are listed in the VITEK® 2 GP Quality Control Tables. Process these according to the procedure for test isolates outlined in this document.

Certification Statement

This is to certify that bioMérieux complies with ISO 13485 and FDA Quality System Regulation (QSR) requirements for design, development, and manufacture of microbial identification systems.

Frequency of Testing

Currently, it is recommended that you use your most stringent inspecting agency’s guidelines for frequency of identification product testing.

Common practice is to perform QC upon receipt of shipment of the test kits. Reactions must follow Instructions for Use results.

If the results do not meet the criteria, subculture for purity and repeat the test. If discrepant results are repeated, perform an alternate identification method and contact bioMérieux.

Testing and Storage of QC Organisms

1. Rehydrate the organism according to the manufacturer's instructions.
2. Use Trypticase Soy with 5% sheep blood agar (TSAB) and incubate at 35°C to 37°C in 5% to 10% CO_2 for approximately 18 to 24 hours.
3. Check for purity. Perform second subculture for testing.
4. Use Trypticase Soy with 5% sheep blood agar (TSAB) and incubate at 35°C to 37°C in 5% to 10% CO_2 for approximately 18 to 24 hours.

Short-Term Storage Conditions

1. Subculture to a TSAB plate or slant.
2. Incubate for 24 hours at 35°C to 37°C in 5% to 10% CO_2.
3. Refrigerate at 2°C to 8°C for up to two weeks.
4. Subculture once as described above and use for QC.

Long-Term Storage Conditions

1. Make a heavy suspension in Tryptic Soy Broth (TSB) with 15% glycerol.

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2. Freeze at -70°C.
3. Subculture to TSAB twice before running QC.

Note: Avoid repeated thawing and refreezing by either freezing in single-use aliquots or removing a small portion of frozen organism preparation with a sterile applicator stick.

STREAMLINED QUALITY CONTROL

Note: Industrial Use Only laboratories should perform quality control following the Streamlined Quality Control section. No additional testing is required for these users.

As there are no substrates that are consistently sensitive to degradation during shipping conditions, streamlined quality control may be conducted by testing two strains: one that is mostly positive and the other, which is mostly negative for reactions on GP. (See GP Quality Control tables for more details.)

COMPREHENSIVE QUALITY CONTROL

Customers who do not qualify for streamlined quality control testing are required to perform comprehensive quality control testing, which entails demonstration of a positive and negative reaction for each substrate of an identification product.[6]

In order to qualify initially for streamlined quality control testing, the CLSI® M50-A standard requires that the user perform and document either of the following:[5]

- Verification testing to show that performance is equivalent to the manufacturer's claims.
- Comprehensive quality control testing of at least three lots over at least three different seasons.

Refer to the complete CLSI® M50-A standard for information regarding continued qualification and further details of requirements and responsibilities for both the user and the manufacturer related to streamlined quality control testing.

GP Quality Control Tables:

***Enterococcus casseliflavus* ATCC® 700327™** (for streamlined or comprehensive quality control)

***Streptococcus salivarius* ssp. *thermophilus* ATCC® 19258™** (for comprehensive quality control)

***Kocuria kristinae* ATCC® BAA-752™** (for comprehensive quality control)

***Listeria monocytogenes* ATCC® BAA-751™** (for comprehensive quality control)

***Streptococcus pneumoniae* ATCC® 49619™** (for comprehensive quality control)

***Staphylococcus saprophyticus* ATCC® BAA-750™** (for streamlined or comprehensive quality control)

***Staphylococcus sciuri* ATCC® 29061™** (for comprehensive quality control)

***Streptococcus equi* ssp. *zooepidemicus* ATCC® 43079™** (for comprehensive quality control)

***Enterococcus saccharolyticus* ATCC® 43076™** (for comprehensive quality control)

The GP card typically identifies the quality control organisms as one-choice or within a low discrimination or slashline identification. However, strains are chosen for reaction performance over identification performance. Therefore, an unidentified or misidentified result may occur when all expected quality control reactions are correct.

Table 6: QC Organism: *Enterococcus casseliflavus* ATCC® 700327™ (for streamlined or comprehensive quality control)

AMY	+	CDEX	-	BGURr	-	URE	-	dMAL	+	PUL	-
PIPLC	-	AspA	v[1]	AGAL	+	POLYB	+	BACI	+	dRAF	+
dXYL	+	BGAR	+	PyrA	+	dGAL	+	NOVO	+	O129R	+
ADH1	v[2]	AMAN	v	BGUR	-	dRIB	+	NC6.5	+	SAL	+
BGAL	+	PHOS	-	AlaA	v	ILATk	-	dMAN	+	SAC	+
AGLU	v	LeuA	v	TyrA	+	LAC	+	dMNE	+	dTRE	+
APPA	v	ProA	-	dSOR	v	NAG	+	MBdG	+	ADH2s	v
										OPTO	+

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

[1]Reaction is mostly positive although occasional negative reaction may occur.

[2]Reaction updated to variable but will not appear in QC software until 9.02.

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+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 8: QC Organism: *Kocuria kristinae* ATCC® BAA-752™ (for comprehensive quality control)

AMY	v	CDEX	v	BGURr	v	URE	v	dMAL	v	PUL	v
PIPLC	v	AspA	–	AGAL	–	POLYB	v	BACI	–	dRAF	v
dXYL	v	BGAR	v	PyrA	v	dGAL	v	NOVO	–	O129R	v
ADH1	v	AMAN	v	BGUR	v	dRIB	v	NC6.5	v	SAL	v
BGAL	–	PHOS	v	AlaA	v	lLATk	+	dMAN	v	SAC	v
AGLU	+	LeuA	+	TyrA	v	LAC	–	dMNE	v	dTRE	v
APPA	–	ProA	+	dSOR	v	NAG	v	MBdG	v	ADH2s	v
										OPTO	v

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 9: QC Organism: *Listeria monocytogenes* ATCC® BAA-751™ (for comprehensive quality control)

AMY	+	CDEX	+	BGURr	v	URE	v	dMAL	v	PUL	v
PIPLC	+	AspA	v	AGAL	v	POLYB	+	BACI	v	dRAF	–
dXYL	v	BGAR	–	PyrA	v	dGAL	–	NOVO	v	O129R	v
ADH1	–	AMAN	+	BGUR	v	dRIB	v	NC6.5	+	SAL	v
BGAL	–	PHOS	v	AlaA	v	lLATk	v	dMAN	–	SAC	–
AGLU	+	LeuA	v	TyrA	v	LAC	v	dMNE	v	dTRE	v
APPA	v	ProA	v	dSOR	v	NAG	+	MBdG	v	ADH2s	v
										OPTO	v

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 10: QC Organism: *Streptococcus pneumoniae* ATCC® 49619™ (for comprehensive quality control)

AMY	v	CDEX	v	BGURr	v	URE	v	dMAL	v	PUL	v
PIPLC	v	AspA	v	AGAL	v	POLYB	v	BACI	–	dRAF	v[1]
dXYL	v	BGAR	v	PyrA	v	dGAL	v	NOVO	v	O129R	–
ADH1	v	AMAN	–	BGUR	v	dRIB	–	NC6.5	–	SAL	v[1]
BGAL	v	PHOS	v	AlaA	+	lLATk	v	dMAN	v	SAC	v
AGLU	v	LeuA	v	TyrA	v	LAC	v	dMNE	v	dTRE	v
APPA	v[1]	ProA	v	dSOR	v	NAG	v	MBdG	v	ADH2s	v
										OPTO	–

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+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

[1]Reaction updated to variable but will not appear in QC software until 9.02

Table 11: QC Organism: *Staphylococcus saprophyticus* ATCC® BAA-750™ (for streamlined or comprehensive quality control)

AMY	-	CDEX	-	BGURr	-	URE	+	dMAL	+	PUL	-
PIPLC	-	AspA	-	AGAL	-	POLYB	-	BACI	v	dRAF	-
dXYL	-	BGAR	-	PyrA	v	dGAL	v	NOVO	+	O129R	v
ADH1	v	AMAN	-	BGUR	-	dRIB	v	NC6.5	+	SAL	-
BGAL	+	PHOS	v	AlaA	-	lLATk	v	dMAN	+	SAC	+
AGLU	v	LeuA	-	TyrA	-	LAC	+	dMNE	v[1]	dTRE	+
APPA	v	ProA	-	dSOR	-	NAG	v	MBdG	-	ADH2s	-
										OPTO	+

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

[1]Reaction is mostly negative although occasional positive reaction may occur.

Table 12: QC Organism: *Staphylococcus sciuri* ATCC® 29061™ (for comprehensive quality control)

AMY	v	CDEX	v	BGURr	+	URE	–	dMAL	v	PUL	v
PIPLC	v	AspA	v	AGAL	v	POLYB	–	BACI	v	dRAF	–
dXYL	–	BGAR	v	PyrA	v	dGAL	v	NOVO	v	O129R	v
ADH1	+	AMAN	v	BGUR	+	dRIB	v	NC6.5	v	SAL	v
BGAL	v	PHOS	+	AlaA	–	lLATk	v	dMAN	v	SAC	v
AGLU	v	LeuA	–	TyrA	–	LAC	–	dMNE	v	dTRE	+
APPA	–	ProA	v	dSOR	v	NAG	v	MBdG	+	ADH2s	–
										OPTO	v

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 13: QC Organism: *Streptococcus equi* ssp. *zooepidemicus* ATCC® 43079™ (for comprehensive quality control)

AMY	v	CDEX	v	BGURr	v	URE	v	dMAL	v	PUL	v[1]
PIPLC	v	AspA	v	AGAL	v	POLYB	v	BACI	v	dRAF	v
dXYL	v	BGAR	v	PyrA	v	dGAL	+	NOVO	v	O129R	v
ADH1	+[2]	AMAN	v	BGUR	v	dRIB	v[3]	NC6.5	v	SAL	v
BGAL	v	PHOS	+	AlaA	v	lLATk	v	dMAN	v	SAC	v
AGLU	v	LeuA	v	TyrA	v	LAC	v	dMNE	v	dTRE	–
APPA	+[2]	ProA	v	dSOR	v	NAG	v	MBdG	v	ADH2s	+
										OPTO	v

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

[1]Reaction is mostly positive although occasional negative reaction may occur.

[2]Reaction updated to positive but will not appear in QC software until 9.02.

[3]Reaction updated to variable but will not appear in QC software until 9.02.

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+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

LIMITATIONS

The VITEK® 2 GP card cannot be used with a direct clinical specimen or sample or other sources containing mixed flora. Any change or modification in the procedure may affect the results.

Newly described or rare species may not be included in the GP database. Selected species will be added as strains become available.

Warning: Testing of unclaimed species may result in an unidentified result or a misidentification.

PERFORMANCE CHARACTERISTICS

For 7.01, 8.01, and 9.01 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 GP identification card was evaluated using 457 clinical and stock isolates of both commonly and rarely observed species of gram-positive cocci. The reference identification was determined with API® STAPH and API® 20 STREP identification kits. Overall, the VITEK® 2 GP correctly identified 96.1% of the isolates, including 3.9% low discrimination with the correct species listed. Misidentifications occurred at 3.5% and no identifications occurred at 0.4%.

For 9.02 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 GP identification card was evaluated using 457 clinical and stock isolates of both commonly and rarely observed species of gram-positive cocci. The reference identification was determined with API® STAPH and API® 20 STREP identification kits. Overall, the VITEK® 2 GP correctly identified 95.8% of the isolates, including 3.5% low discrimination with the correct species listed. Misidentifications occurred at 3.5% and no identifications occurred at 0.7%.

*Data on file at bioMérieux, Inc.

ORGANISMS IDENTIFIED

Claims are for all software users unless otherwise stated.

- *Abiotrophia defectiva*
- *Aerococcus urinae*
- *Aerococcus viridans*
- *Alloiococcus otitis*
- *Dermacoccus nishinomiyaensis/Kytococcus sedentarius*
- *Enterococcus avium*
- *Enterococcus casseliflavus*
- *Enterococcus cecorum*
- *Enterococcus columbae*
- *Enterococcus durans*
- *Enterococcus faecalis*
- *Enterococcus faecium*
- *Enterococcus gallinarum*
- *Enterococcus hirae*
- *Enterococcus raffinosus*
- *Enterococcus saccharolyticus*

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- *Vagococcus fluvialis*

Additional Claims For 8.01 or Higher Software Users

- *Listeria fleischmannii*
- *Listeria rocourtiae*
- *Streptococcus iniae*

Taxonomy Changes For 9.02 Software Users

- *Streptococcus suis* (formerly *Streptococcus suis* I and II)

**Staphylococcus aureus* claim contains only the subspecies *aureus.*

+ OMA Official Methods of Analysis validated claim.

SUPPLEMENTAL TESTS

Table 15: GP Supplemental Tests

Abbreviation	Test Name	Description	Comments	Reference
For 7.01 or Higher Software Users				
A-HEM	ALPHA HEMOLYSIS	Certain species produce incomplete hemolysis resulting in a green zone around colonies on blood based media.	N/A	21, 22, 23, 26, 27
	Acidification of:	Acidification of carbon source observed with pH indicators (e.g., phenol red, bromcresol purple).	Some tests also appear on the GP card but are recommended as supplemental tests since results of conventional macromethods may differ from rapid commercial micromethods.	2, 3, 4, 8, 10, 14, 16, 17, 21, 22, 23, 26 , 27, 28, 31, 32, 33, 34, 36, 40, 42
AMD/STARCH	AMIDON/STARCH			
GLYCOGENac	GLYCOGEN			
lARABINOSE	L-ARABINOSE acid.			
INULIN	INULIN			
MdG	METHYL-A-DGLUCOPYRANOSIDE			
MdM	METHYL-A-DMANNOPYRANOSIDE			
PULLULAN	PULLULAN			
SACCHAROSE	SACCHAROSE (SUCROSE)			
dGALACTOSE	D-GALACTOSE			
dMALTOSE	D-MALTOSE			
dMANNITOL	D-MANNITOL			
dMANNOSE	D-MANNOSE			
dMELEZIT	D-MELEZITOSE			
dMELIBIOSE	D-MELIBIOSE			
dRAFFINOSE	D-RAFFINOSE			
dRIBOSE	D-RIBOSE			
dSORBITOL	D-SORBITOL			
dTREHALOSE	D-TREHALOSE			
dXYLOSE	D-XYLOSE			
lRHAMNOSE	L-RHAMNOSE			
ANANE	ALPHA-D-N-ACETYLNEURAMINIDASE	Presence of respective enzyme cleaves substrate generating detectable leaving group (e.g., p-nitrophenol, methyl umbelliferone, beta-naphthylamide, beta-naphthol, p-nitroaniline, 7-amidomethyl-coumarin).	Presence of enzyme is indicated by generation of a colored or fluorescent product, or a noncolored product that forms color upon addition of a specific reagent.	7, 11, 21, 22, 23, 27, 28, 29, 31, 34, 39
AIFUC	ALPHA-L-FUCOSIDASE			
BGLU	BETA-GLUCOSIDASE			
BGURase	BETA-GLUCURONIDASE			
BNAG	BETA-N-ACETYLGLUCOSAMINIDASE			
BNAGA	BETA-N-ACETYLGALACTOSAMINIDASE			
BdFUC	BETA-D-FUCOSIDASE			
PAL	ALKALINE PHOSPHATASE			
Pyrro. Ary.	Pyrrolidonyl ARYLAMIDASE			

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Abbreviation	Test Name	Description	Comments	Reference
Adherence	Adherence to agar	Sticking of colonies to the agar surface	Characteristic of *Rothia mucilaginosa*	27
AER.GROWTH	AEROBIC GROWTH	Growth in air	N/A	22
Arg.hydr.	ARGININE dihydrolase	Hydrolysis of arginine releases an amine resulting in alkalinization of the medium observed with a pH indicator (e.g., purple color formation in the presence of bromcresol purple).	N/A	2, 21, 22, 27, 36
B-HEM	BETA HEMOLYSIS	Certain species possess hemolysins that give a transparent zone around colonies on blood-based agars.	N/A	21, 22, 27, 37
BILE SOL	BILE SOLUBILITY	Pneumococcal colonies completely lyse and disappear when exposed to a 10% solution of deoxycholate.	Rapid test for *Streptococcus pneumoniae*	27
CAMP (S.au)	CAMP TEST (*Staph. aureus*)	Synergistic hemolysis of *Listeria monocytogenes* colonies by beta-toxin producing colonies of *Staphylococcus aureus*.	N/A	27
CAT	CATALASE	Colony placed on a drop of 3% hydrogen peroxide produces gas bubbles. The bacteria that contain cytochrome enzyme are catalase positive.	Differentiation of *Micrococcaceae* (+) from *Streptococcaceae* (–)	21, 22, 27, 38
CLINDA.S	Clindamycin susceptible	Zone of inhibition around the clindamycin disk > 20mm	Used to differentiate *Lactococcus lactis* and *Lactococcus garvieae*.	15
ESCULIN	ESCULIN hydrolysis	Hydrolysis of esculin forms esculetin which produces a black pigment in the presence of iron salts.	N/A	3, 18, 21, 22, 24, 27
Gas prod.	Gas production	Production of CO_2 from degradation carbohydrate (e.g., glucose) metabolism.	N/A	27
HIP	HIPPURATE hydrolysis	Hydrolysis of sodium hippurate releases glycine that produces a blue colored product after addition of ninhydrin.	N/A	12, 21, 22, 26, 27
LAP	LEUCINE AMINOPEPTIDASE	The substrate leucine-beta-naphthylamide is hydrolyzed by the enzyme leucine aminopeptidase and released beta-naphthylamide combines with the cinnamaldehyde reagent to form a bright pink to cherry red pigment.	N/A	19
LitmusMILK	Litmus Milk Medium	Acid production in Litmus Milk	N/A	16

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Abbreviation	Test Name	Description	Comments	Reference
NaCl 6.5%	GROWTH IN 6.5% NaCl	Growth in 6.5% NaCl broth	N/A	16, 19
NO3	NITRATE REDUCTION	Test for the ability to reduce nitrate to nitrite or nitrogen gas.	N/A	21, 22, 38
NOVO_R OPTO_R VANCO_R	NOVOBIOCIN_RESISTANCE OPTOCHIN_RESISTANCE VANCOMYCIN_RESISTANCE	Ability of certain species to grow in the presence of specific antibacterial compounds	N/A	20, 21, 22, 27, 40, 41
NaCl 7.5%	GROWTH IN 7.5% NACL	Ability of certain species to grow in the presence of a high concentration of NaCl	N/A	22
PI/OR/RED	PINK/ORANGE/RED PIGMENT	Ability of certain species to produce pink, orange, or red colonies on nondifferential media	Characteristic of *Kocuria rosea*	22, 27, 28
PVATE	PYRUVATE	Ability to use pyruvate as a sole carbon source	N/A	28
SATELLITE	SATELLITE behavior	Appearance of satellite colonies of nutritionally deficient *Streptococcaceae* around colonies of *Staphylococcus epidermidis*.	Nutritionally deficient *Streptococcaceae* require nutritional factors supplied by metabolism of colonies of *Staphylococcus epidermidis*.	9, 27
Str.sero.A Str.sero.B Str.sero.C Str.sero.D Str.sero.G	Strepto Serology A Strepto Serology B Strepto Serology C Strepto Serology D Strepto Serology G	Agglutination tests for *Streptococcus* groups A,B,C,D, and G	N/A	1, 13, 18, 21, 22, 24, 27, 28, 36
UREASE	Urease	Hydrolysis of urea releases ammonia, resulting in alkalinization of the medium observed with a pH indicator (e.g., red color formation in the presence of phenol red).	N/A	21, 22, 25, 27
VP	VOGES PROSKAUER	Ability of some species to produce acetoin from glucose fermentation.	N/A	20, 21, 27, 34
YELLOW	YELLOW PIGMENT	Ability of certain species to produce yellow pigmented colonies on nondifferential media.	For example, used to differentiate *E. casseliflavus* (+) from *E. gallinarum* (–).	21, 22, 27, 28
For 7.01 Software Users Only				
BILE ESC	BILE ESCULIN	Bile-esculin positive organisms are able to grow in the presence of 40% bile and to hydrolyze esculin.	N/A	18, 24
CAROTENOID	CAROTENOID PIGMENT	Presence of red, pink, or orange pigment	N/A	28

REFERENCES

1. Balows A, Hausler Jr. WJ, Herrmann KL, Isenberg HD, Shadomy HJ. *Manual of Clinical Microbiology* 5th edition. American Society of Microbiology, Washington, D.C.1991.

2. Barros RR, Carvalho GS, Peralta JM, Facklam RR, Teixeira LM. Phenotypic and Genotypic Characterization of *Pediococcus* Strains Isolated from Human Clinical Sources. J. Clin. Microbiol. 2001. 39:1241- 1246.

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30. National Committee for Clinical Laboratory Standards, M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue* — Approved Guideline, 1997.
31. Poyart C, Quesne G, Trieu-Cuot P. Taxonomic dissection of the *Streptococcus bovis* group by analysis of manganese-dependent superoxide dismutase gene (sodA) sequences: reclassification of *Streptococcus infantarius* subsp. *coli* as *Streptococcus lutetiensis* sp.nov. and of *Streptococcus bovis* biotype II.2 as *Streptococcus pasteurianus* sp nov. Int. J. Syst. Evol. Microbiol. 2002. 52:1247-1255.
32. Schlegel L, Grimont F, Collins MD, Regnault B, Grimont PAD, Bouvet A. *Streptococcus infantarius* sp. nov., *Streptococcus infantarius* subsp *infantarius* subsp. nov. and *Streptococcus infantarius* subsp *coli* subsp. nov., isolated from humans and food. Int. J. Syst. Evol. Microbiol. 2000. 50:1425-1434.
33. Schlegel, L., F. Grimont, E. Ageron, P. A. D. Grimont, and A. Bouvet. Reappraisal of the taxonomy of the *Streptococcus bovis/Streptococcus equinus* complex and related species: description of *Streptococcus gallolyticus* subsp. *gallolyticus* subsp. nov., *S. gallolyticus* subsp. *macedonicus* subsp. nov. and *S. gallolyticus* subsp. *pasteurianus* subsp. nov. Int. J. Syst. Evol. Microbiol. 2003. 53:631-645.
34. Takashi, S., K. Kikuchi, Y. Tanaka, N. Takahashi, S. Kamata and K. Hiramatsu. Reclassification of Phenotypically Identified *Staphylococcus intermedius* Strains. J.Clin.Microbiol. 2007. 45(9): 2770-2778.
35. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.
36. Versalovic J, Carroll KC, Funke G, Jorgensen JH, Landry ML, Warnock DW, editors. Manual Of Clinical Microbiology, Volume 1, 10th edition. American Society for Microbiology, Washington, D.C. 2011.
37. Viera VV, Teixeira LM, Zahner V, Momen H, Facklam RR, Steigerwalt AG, Brenner DJ, Castro ACD. Genetic relationships among the different phenotypes of *Streptococcus dysgalactiae* strains. Int. J. Syst. Bacteriol. 1998. 48:1231-1243.
38. Von Graevenitz, A. *Rothia dentocariosa*: taxonomy and differential diagnosis. Clin.Microbiol. and Infection, 2004. 10:399-402.
39. Whiley RA, Hall LMC, Hardie JM, Beighton D. A study of small colony beta hemolytic, Lancefield group C streptococci within the anginosus group: description of *Streptococcus constellatus* subsp. *pharyngis* subsp.nov., associated with the human throat and pharyngitis. Int. J. Syst. Bacteriol. 1999. 49:1443-1449.
40. Whiley, R.A. & Hardie ,J.M. 2009. *Streptococcaceae*. In *Bergey's Manual of Systematic Bacteriology,* 2nd ed., Volume 3. DeVos, P., Garrity, G., Jones, D., Krieg, N.R., Ludwig, W., Rainey, F.A., Schleifer, K.-H., and Whitman, W.B. (eds), Springer, New York, pp. 707-708.
41. Jorgensen, JH., Pfaller, M., Carroll, K., Funke, G., Landry, ML., Richter, S., Warnock, D. *Manual of Clinical Bicrobiology*, Volume 1, 11th edition. American Society for Microbiology, Washington, D.C., 2015.
42. Winn, W.C., S. Allen, W. Janda, E. Koneman, G. Procop, P. Schreckenberger, G. Woods. *Color Atlas and Textbook of Diagnostic Microbiology*, 6th ed., 2006, Lippincott Williams & Wilkins, Philadelphia, Baltimore.

Use this Instructions for Use with VITEK® 2 Product No. 21342.

INDEX OF SYMBOLS

Symbol	Meaning
REF	Catalog number
IVD	In Vitro Diagnostic Medical Device
	Legal Manufacturer
	Temperature limitation
	Use by date
LOT	Batch code
	Consult Instructions for Use

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Symbol	Meaning
	Date of manufacture
	Contains sufficient for <n> tests
EC REP	Authorized representative in the European Community
℞ only	For US Only : Caution : US Federal Law restricts this device to sale by or on the order of a licensed practitioner

Instructions for Use provided in the kit or downloadable from www.biomerieux.com/techlib

LIMITED WARRANTY

bioMérieux warrants the performance of the product for its stated intended use provided that all procedures for usage, storage and handling, shelf life (when applicable), and precautions are strictly followed as detailed in the instructions for use (IFU).

Except as expressly set forth above, bioMérieux hereby disclaims all warranties, including any implied warranties of merchantability and fitness for a particular purpose or use, and disclaims all liability, whether direct, indirect or consequential, for any use of the reagent, software, instrument and disposables (the “System”) other than as set forth in the IFU.

WASTE DISPOSAL

All hazardous waste must be disposed of by following your local inspecting agency's guidelines.

REVISION HISTORY TABLE

Change type categories	
N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user
Note :	Minor typographical, grammar, and formatting changes are not included in the revision history.

Release Date	Part Number	Change Type	Change Summary
2019-03	043900-03	Technical change	Updated for 9.02 software release. Updated sections: • Intended Use • Precautions • Additional Information on Lab Report • Testing of QC Organisms • Performance Characteristics • Organisms Identified • References
2016-10	043900-02	Technical change	• Updated content to reflect the 8.01 Product Information Manual

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Release Date	Part Number	Change Type	Change Summary
2016-05	043900-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• New IFU derived from product chapter in the Product Information Manual • Updated Limited Warranty section • Updated with RX only information

BIOMÉRIEUX, the BIOMÉRIEUX logo, VITEK, API, Count-TACT, chromID, DensiCHEK and bioLiaison are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

This product may be protected by one or more patents, see: http://www.biomerieux-usa.com/patents.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

Any other name or trademark is the property of its respective owner.

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℞ only IVD

INTENDED USE

These Instructions for Use correspond to the VITEK® 2 Systems 7.01 or higher software. If you are not using VITEK® 2 Systems 7.01 or higher software, please refer to the VITEK® 2 Systems Product Information that you received with your current software version.

The VITEK® 2 Neisseria-Haemophilus identification card (NH) is intended for use with the VITEK® 2 Systems for the automated identification of most clinically significant fastidious organisms. The VITEK® 2 NH identification card is a single-use disposable. For a list of claimed species, see the Organisms Identified section.

DESCRIPTION

The NH Card is based on established biochemical methods and newly developed substrates measuring carbon source utilization and enzymatic activities. There are 30 biochemical tests. Final identification results are available in approximately six hours.

For a list of well contents, see the NH Well Contents table.

Table 1: NH Well Contents

Well	Test	Mnemonic	Amount/Well
1	Arginine ARYLAMIDASE	ArgA	0.0324 mg
2	GAMMA-GLUTAMYL-TRANSFERASE	GGT	0.0228 mg
3	L-Lysine-ARYLAMIDASE	LysA	0.0228 mg
4	D-GALACTOSE	dGAL	0.3 mg
5	Leucine ARYLAMIDASE	LeuA	0.023 mg
6	ELLMAN	ELLM	0.03 mg
7	Phenylalanine ARYLAMIDASE	PheA	0.026 mg
8	L-Proline ARYLAMIDASE	ProA	0.023 mg
10	L-Pyrrolidonyl-ARYLAMIDASE	PyrA	0.018 mg
13	Tyrosine ARYLAMIDASE	TyrA	0.0279 mg
15	Ala-Phe-Pro-ARYLAMIDASE	APPA	0.038 mg
18	D-GLUCOSE	dGLU	0.3 mg
19	GLYCOGENE	GLYG	0.18 mg
20	D-MANNOSE	dMNE	0.3 mg
22	D-MALTOSE	dMAL	0.3 mg
28	SACCHAROSE/SUCROSE	SAC	0.3 mg
33	N-ACETYL-D-GLUCOSAMINE	NAG	0.3 mg
36	UREASE	URE	0.15 mg
39	BETA-GALACTOPYRANOSIDASE Indoxyl	BGALi	0.006 mg
40	ORNITHINE DECARBOXYLASE	ODC	0.15 mg
41	ALPHA-ARABINOSIDASE	AARA	0.0324 mg
45	PYRUVATE	PVATE	0.15 mg
46	PHOSPHORYL CHOLINE	PHC	0.0366 mg
47	D-MALATE	dMLT	0.15 mg

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Note: Other well numbers between 1 and 64 not designated in this table are empty.

PRECAUTIONS

Note: For industry customers that need assistance on selecting the correct VITEK® 2 identification card, please refer to the VITEK® 2 Compact Instrument User Manual chapter, "Guidance to Select a VITEK® 2 Identification Card."

- For *In Vitro* Diagnostic Use Only.
- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For professional use only.
- Suspensions not within the appropriate zone on the VITEK® 2 DENSICHEK™ Plus or the VITEK® 2 DENSICHEK™ may compromise card performance.
- Do not use the card after the expiration date shown on the package liner.
- Store the card unopened in the package liner. Do not use the card if the protective package liner is damaged or if no desiccant is present.
- Allow the card to come to room temperature before opening the package liner.
- Do not use powdered gloves. Powder may interfere with the optics.
- Use of culture media other than the recommended types must be validated by the customer laboratory for acceptable performance.
- A Gram stain should be performed to determine an organism's Gram reaction and morphology prior to selecting the identification card to inoculate.
- The card performs as intended only when used in conjunction with VITEK® 2 Systems, following the instructions contained in these Instructions for Use.
- **Do not use glass test tubes**. Use clear plastic (polystyrene) tubes only. Variation exists among test tubes of standard diameter. Carefully place the tube into the cassette. If resistance is encountered, discard and try another tube that does not require pressure to insert.
- Prior to inoculation, inspect cards for tape tears or damage to the tape and discard any that are suspect. Check the saline level in the tubes after the cassette has been processed to ensure proper filling of card.
 - VITEK® 2 60 or VITEK® 2 XL: Eject improperly filled cards.
 - VITEK® 2 Compact: Do not load improperly filled cards.
- Give special consideration to specimen source and patient drug or antimicrobic regimen.
- Give special consideration to sample source.
- Interpretation of test results requires the judgment and skill of a person knowledgeable in microbial identification testing. Additional testing may be required. (See the Supplemental Tests section.)
- Do not clean saline dispenser with chemical agents. The use of chemical agents may impact card performance.

Warning: All patient specimens, microbial cultures, and inoculated VITEK® 2 cards, along with associated materials, are potentially infectious and should be treated with universal precautions.[18,20]

Warning: All hazardous waste must be disposed of by following your local inspecting agency's guidelines.

STORAGE CONDITIONS

Upon receipt, store VITEK® 2 NH cards unopened in their original package liner at 2°C to 8°C.

SPECIMEN PREPARATION

For specimen preparation information, see the Culture Requirements Table.

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Table 2: Culture Requirements Table

VITEK® 2 Card	Media	Age of Culture[1]	Incubation Conditions	Inoculum Density	Dilution for AST	Age of Suspension Before Loading Instrument
NH	*Campylobacter*: TSAB[2] CBA CHBA TSAHB	*Campylobacter*: 18 to 24 hours	*Campylobacter*: 35°C to 37°C or 40°C to 42°C microaerobic conditions	2.70 – 3.30 McFarland Standard	N/A [5]	≤ 30 minutes
	Haemophilus: CHOC[2] CHOC PVX[2] CBA CHOC + B	Fastidious: 18 to 24 hours	Fastidious: 35°C to 37°C with 5% to 10% CO_2			
	Neisseria: CHOC[2] CHOC PVX[2] CHOC VCAT CHBA ML[3] NYC[4] TM[3] TSAB					
	Other Fastidious: CHOC[2] CHOC PVX[2] CBA CHBA ML[3] TM[3] TSAB TSAHB					

[1]Cultures with scant or poor growth may give unidentified or incorrect results even when the Age of Culture requirements are met.

[2]These media were used in the identification product database development and will give optimal performance.

[3]These media were validated for *Neisseria gonorrhoeae*, *Neisseria meningitidis* and *Moraxella catarrhalis*.

[4]This medium was validated for *Neisseria gonorrhoeae*.

[5]N/A = not applicable

Culture Requirements Table — Media Abbreviations

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CBA = Columbia Blood Agar with 5% Sheep Blood

CHBA = Columbia Horse Blood Agar

CHOC = Chocolate Agar

CHOC + B = Chocolate Agar with Bacitracin

CHOC PVX = Chocolate Polyvitex Agar

CHOC VCAT = Chocolate Polvitex Agar with VCAT

ML = Martin-Lewis Agar

NYC = New York City Medium

TM = Thayer-Martin Agar

TSAB = Trypticase Soy Agar with 5% Sheep Blood

TSAHB = Trypticase Soy Agar with 5% Horse Blood

TEST PROCEDURE

Materials

When used with VITEK® 2 instrumentation, the NH Card is a complete system for routine identification testing of most significant fastidious organisms.

Required materials are:

- VITEK® 2 NH Card
- DENSICHEK™ Plus Kit or VITEK® DENSICHEK® Kit
- DENSICHEK™ Plus Standards Kit or DENSICHEK® Standards Kit
- VITEK® 2 Cassette
- Sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- 12 mm x 75 mm clear plastic (polystyrene) disposable test tubes
- Sterile sticks or swabs
- Appropriate agar medium (see Culture Requirements table).

Optional accessories:

- Adjustable volume saline dispenser
- Loops
- Pre-dispensed saline test tubes (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- Test tube caps
- Vortex

Procedure

Warning: Failure to follow instructions and recommendations provided in this section for performing laboratory tasks may cause erroneous or delayed results.

For product-specific information, see the Culture Requirements table.

Note: Prepare the inoculum from a pure culture, according to good laboratory practices. In case of mixed cultures, a re-isolation step is required. It is recommended that a purity check plate be done to ensure that a pure culture was used for testing.

1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.
2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).
3. Use a sterile stick or swab to transfer a sufficient number of morphologically similar colonies to the saline tube prepared in step 2. Prepare a homogenous organism suspension with a density equivalent to a McFarland No. 2.70 to 3.30 using a calibrated VITEK® 2 DENSICHEK™ Plus or VITEK® 2 DENSICHEK™.
 Note: Age of suspension must not exceed 30 minutes before inoculating card.
4. Place the suspension tube and NH card in the cassette.

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5. Refer to the appropriate Instrument User Manual for instructions on data entry and how to load the cassette into the instrument.
6. Follow your local inspecting agency's guidelines for disposal of hazardous waste.

RESULTS

Identification Analytical Techniques

VITEK® 2 Systems identify an organism by using a methodology based on the characteristics of the data and knowledge about the organism and reactions being analyzed. Sufficient data have been collected from known strains to estimate the typical reactions of the claimed species to a set of discriminating biochemicals. If a unique identification pattern is not recognized, a list of possible organisms is given, or the strain is determined to be outside the scope of the database.

The printed lab report contains suggestions for any supplemental tests necessary to complete the identification. If the tests are not sufficient to complete the identification, then standard microbiology references and literature should be consulted.

Table 3: Identification Card Qualifying Messages

ID Message Confidence Level	Choices	% Probability	Comments
Excellent	1	96 to 99	N/A
Very Good	1	93 to 95	N/A
Good	1	89 to 92	N/A
Acceptable	1	85 to 88	N/A
Low Discrimination	2 to 3	Sum of choices = 100; after resolution to one choice, percent probability reflects the number associated with selected choice.	Two to three taxa exhibit same biopattern. Separate by supplemental testing.
Inconclusive or Unidentified Organism	> 3 or 0	N/A	Either > 3 taxa exhibit same biopattern or Very atypical biopattern. Does not correspond to any taxon in the database. Check Gram stain and purity.

PERCENT PROBABILITY

As part of the identification process, the software compares the test set of reactions to the expected set of reactions of each organism, or organism group, that can be identified by the product. A quantitative value, the percent probability, is calculated and relates to how well the observed reactions compare to the typical reactions of each organism. A perfect match between the test reaction pattern and the unique reaction pattern of a single organism, or organism group, would provide a percent probability of 99. When a perfect match is not obtained, it is still possible for the reaction pattern to be sufficiently close to that of an expected reaction pattern such that a clear decision can be provided about the organism identification. The range of percent probabilities in the one-choice case is 85 to 99. Values closer to 99 indicate a closer match to the typical pattern for the given organism.

When the reaction pattern is not sufficient to discriminate between two to three organisms, the percent probabilities reflect this ambiguity. The reported probability values indicate, relatively, the order in which the reaction pattern best corresponds to the listed possibilities. The order does not, however, suggest that the pattern match to one of the possible identifications is clearly superior to another. The probability characteristic of an overall sum of 100 is retained through the calculation process. After resolution to one choice, the probability characteristic of the single choice is retained.

ADDITIONAL INFORMATION ON LAB REPORT

Supplemental test — External (offline) test that allows the user to resolve a slashline or Low Discrimination identification. Numbers in parentheses indicate percent positive reaction for the species/test listed.

Contraindicating test — Test result that is unusual for a reported taxon.

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Table 4: Notes Associated with Certain Taxa

Taxa	Note			
For 7.01 or Higher Software Users				
Haemophilus influenzae	*Haemophilus aegyptius* is a recognized species, but there is controversy as to whether it should be retained as a legitimate species. *Haemophilus aegyptius* is indistinguishable from *H. influenzae* by either DNA/DNA hybridization or by any single phenotypic test. Isolates of *H. aegyptius* exhibit distinct pathogenicity and are associated with cases of acute purulent conjunctivitis. *Haemophilus influenzae* biogroup *Aegyptius* is also indistinguishable from *H. aegyptius* and *H. influenzae*, but is considered the etiologic agent of Brazilian purpuric fever, which is a systemic pediatric infection that is typically preceded by purulent conjunctivitis that resolves before the onset of the systemic infection. Consequently, isolates of *H. aegyptius*, *H. influenzae* biogroup *Aegyptius*, and other biogroups of *H. influenzae* will all identify as *H. influenzae* when tested with the NH card.			
Neisseria gonorrhoeae *Neisseria meningitidis*	Critical pathogen The species identified may have significance to patient or sample outcome and can be stopped for review.			
Neisseria sicca	Possibility of *N. flavescens* or *N. mucosa*. Isolates of these species may be misidentified as *N. sicca*. In order to rule them out, perform the following tests:			
		YELLOW	GLU	NO3
	N. flavescens	+	-	-
	N. mucosa	+	+	+
	N. sicca	-	+	-
For 9.02 Software Users				
Neisseria cinerea	Possibility of *Neisseria gonorrhoeae*			

Notes Associated with an Improperly Filled Card or with a Negative Profile (Biopattern)

- For the case where the time between two readings is greater than 40 minutes: “CARD ERROR — Missing data.”
- For the case where there is a negative profile: “Organism with low reactivity biopattern — please check viability.”
- When a biopattern is calculated for an unknown organism that is completely negative or consists of both negative tests and tests that fall within the uncertainty zone, the identification call will be “Non or low reactive biopattern.”

Campylobacter jejuni ssp. *jejuni* could potentially trigger “Non or low reactive biopattern” if a test was atypical or fell within the uncertainty zone.

QUALITY CONTROL

Quality control organisms and their expected results are listed in the VITEK® 2 NH Quality Control Tables. Process these according to the procedure for test isolates outlined in this document.

Note: *Staphylococcus epidermidis* ATTC® 12228™ needs to be tested at a McFarland Standard No. 0.5 to 0.63. All other QC strains are tested at a McFarland Standard No. 2.70 to 3.30.

Certification Statement

This is to certify that bioMérieux complies with ISO 13485 and FDA Quality System Regulation (QSR) requirements for design, development, and manufacture of microbial identification systems.

Frequency of Testing

Currently, it is recommended that you use your most stringent inspecting agency’s guidelines for frequency of identification product testing.

Common practice is to perform QC upon receipt of shipment of the test kits. Reactions must follow Instructions for Use results.

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If the results do not meet the criteria, subculture for purity and repeat the test. If discrepant results are repeated, perform an alternate identification method and contact bioMérieux.

Testing and Storage of QC Organisms

1. Rehydrate the organism according to the manufacturer's instructions.
2. Use Chocolate agar and incubate at 35°C to 37°C in 5% to 10% CO_2. Incubate for 18 to 24 hours or until sufficient growth is obtained.
3. Check for purity. Perform second subculture for testing.
4. Use Chocolate agar and incubate at 35°C to 37°C in 5% to 10% CO_2. Incubate for 18 to 24 hours.

Short-Term Storage Conditions

Short-term storage is not recommended. Cultures maintained by other methods, specifically those maintained on agar plates or slants for long periods of time at room temperature or at 2°C to 8°C are subject to loss or change of important biochemical characteristics.

Long-Term Storage Conditions

1. Make a heavy suspension in Tryptic Soy Broth (TSB) with 15% glycerol.
2. Freeze at -70°C.
3. Subculture twice to Chocolate agar before running QC.

Note: Avoid repeated thawing and refreezing by either freezing in single-use aliquots or removing a small portion of frozen organism preparation with a sterile applicator stick.

STREAMLINED QUALITY CONTROL

Note: Industrial Use Only laboratories should perform quality control following the Streamlined Quality Control section. No additional testing is required for these users.

Streamlined quality control may be used to confirm acceptable performance of the NH card after shipping/storage. This methodology may be performed with the NH card by following the instructions for quality control testing as described in the NH Instructions for Use and meeting the criteria stated in the CLSI® M50-A Quality Control for Commercial Microbial Identification Systems.

Testing may be conducted using *Eikenella corrodens* ATCC® BAA-1152™ and evaluating the performance of the PHOS well. Testing at bioMérieux, Inc. has shown that the PHOS well is the most labile well on the NH card and *E. corrodens* ATCC® BAA-1152™ is the most sensitive strain for detecting degradation of this well with a false positive reaction. (See NH Quality Control table for more details).

COMPREHENSIVE QUALITY CONTROL

Customers who do not qualify for streamlined quality control testing are required to perform comprehensive quality control testing, which entails demonstration of a positive and negative reaction for each substrate of an identification product.[4]

In order to qualify initially for streamlined quality control testing, the CLSI® M50-A standard requires that the user perform and document either of the following:[3]

- Verification testing to show that performance is equivalent to the manufacturer's claims.
- Comprehensive quality control testing of at least three lots over at least three different seasons.

Refer to the complete CLSI® M50-A standard for information regarding continued qualification and further details of requirements and responsibilities for both the user and the manufacturer related to streamlined quality control testing.

NH Quality Control Tables:

***Eikenella corrodens* ATCC® BAA-1152™** (for streamlined or comprehensive quality control)

***Aggregatibacter aphrophilus* ATCC® 33389™** (for comprehensive quality control)

***Haemophilus influenzae* ATCC® 9007™** (for comprehensive quality control)

***Neisseria gonorrhoeae* ATCC® 19424™** (for comprehensive quality control)

***Neisseria lactamica* ATCC® 23970™** (for comprehensive quality control)

***Oligella urethralis* ATCC® 17960™** (for comprehensive quality control)

***Enterobacter aerogenes* ATCC® 13048™** (for comprehensive quality control)

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The NH card typically identifies the quality control organisms as one-choice or within a low discrimination or slashline identification. However, strains are chosen for reaction performance over identification performance. Therefore, an unidentified or a misidentified result may occur when all expected quality control reactions are correct.

Note: The NH card uses unclaimed taxa for quality control testing. These strains will give an unidentified or a misidentified result.

Table 5: QC Organism: *Eikenella corrodens* ATCC® BAA-1152™ (for streamlined or comprehensive quality control)

ArgA	-	PheA	-	GLYG	-	BGALi	-	MTE	-
GGT	-	ProA	+	dMNE	-	ODC	+	IGLM	v
LysA	-	PyrA	-	dMAL	-	AARA	-	*PHOS**	-
dGAL	-	TyrA	-	SAC	-	PVATE	-	dRIB2	-
LeuA	+	APPA	+	NAG	-	PHC	-	OPS	-
ELLM	+	dGLU	-	URE	-	dMLT	v	dXYL	-

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive.

*Key well for streamlined quality control.

Table 6: QC Organism *Aggregatibacter aphrophilus* ATCC® 33389™ (for comprehensive quality control)

ArgA	v	PheA	v	GLYG	v	BGALi	+	MTE	+
GGT	+	ProA	–	dMNE	+	ODC	–	IGLM	–
LysA	v	PyrA	v	dMAL	+	AARA	v	PHOS	+
dGAL	v	TyrA	v	SAC	+	PVATE	–	dRIB2	v
LeuA	v	APPA	–	NAG	v	PHC	v	OPS	v
ELLM	v	dGLU	+	URE	–	dMLT	v	dXYL	v

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive.

Table 7: QC Organism: *Haemophilus influenzae* ATCC® 9007™ (for comprehensive quality control)

ArgA	v	PheA	+	GLYG	v	BGALi	–	MTE	v
GGT	–	ProA	–	dMNE	v	ODC	v	IGLM	v
LysA	v	PyrA	–	dMAL	–	AARA	v	PHOS	+
dGAL	+	TyrA	v	SAC	v	PVATE	v	dRIB2	+
LeuA	+	APPA	–	NAG	v	PHC	+	OPS	+
ELLM	v	dGLU	+	URE	+	dMLT	+	dXYL	+

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive.

Table 8: QC Organism: *Neisseria gonorrhoeae* ATCC® 19424™ (for comprehensive quality control)

ArgA	+	PheA	v	GLYG	v	BGALi	v	MTE	v
GGT	v	ProA	v	dMNE	v	ODC	–	IGLM	–
LysA	v	PyrA	v	dMAL	v	AARA	–	PHOS	v
dGAL	v	TyrA	v	SAC	v	PVATE	v	dRIB2	v
LeuA	v	APPA	+	NAG	v	PHC	–	OPS	v
ELLM	–	dGLU	v	URE	v	dMLT	–	dXYL	v

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive.

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+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive.

Table 10: QC Organism: *Oligella urethralis* ATCC® 17960™ (for comprehensive quality control)

ArgA	–	PheA	+	GLYG	–	BGALi	v	MTE	–
GGT	+	ProA	+	dMNE	–	ODC	v	IGLM	+
LysA	v	PyrA	v	dMAL	v	AARA	v	PHOS	–
dGAL	–	TyrA	+	SAC	–	PVATE	+	dRIB2	–
LeuA	v	APPA	v	NAG	–	PHC	v	OPS	v
ELLM	+	dGLU	–	URE	v	dMLT	+	dXYL	–

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive.

Table 11: QC Organism: *Enterobacter aerogenes* ATCC® 13048™ (for comprehensive quality control)

ArgA	v	PheA	v	GLYG	v	BGALi	v	MTE	v
GGT	v	ProA	v	dMNE	v	ODC	v	IGLM	v
LysA	+	PyrA	+	dMAL	v	AARA	v	PHOS	v
dGAL	v	TyrA	v	SAC	v	PVATE	v	dRIB2	v
LeuA	v	APPA	v	NAG	+	PHC	v	OPS	v
ELLM	v	dGLU	v	URE	v	dMLT	v	dXYL	v

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive.

Note: *Enterobacter aerogenes* is an unclaimed taxa for the NH card.

Table 12: QC Organism: *Paenibacillus polymyxa* ATCC® 7070™ (for comprehensive quality control)

ArgA	v	PheA	v	GLYG	+	BGALi	v	MTE	v
GGT	v	ProA	v	dMNE	v	ODC	v	IGLM	v
LysA	v	PyrA	v	dMAL	v	AARA	v	PHOS	v
dGAL	v	TyrA	v	SAC	v	PVATE	v	dRIB2	v
LeuA	v	APPA	v	NAG	v	PHC	v	OPS	v
ELLM	v	dGLU	v	URE	v	dMLT	v	dXYL	v

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive.

Note: *Paenibacillus polymyxa* is an unclaimed taxa for the NH card.

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+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive.

Note: *Staphylococcus epidermidis* is an unclaimed taxa for the NH card.

LIMITATIONS

The VITEK® 2 NH card cannot be used with a direct clinical specimen or sample or other sources containing mixed flora. Any change or modification in the procedure may affect the results.

Newly described or rare species may not be included in the NH database. Selected species will be added as strains become available.

Warning: Testing of unclaimed species may result in an unidentified result or a misidentification.

PERFORMANCE CHARACTERISTICS

For 7.01 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 NH identification card was evaluated using 371 clinical and stock isolates of both commonly and rarely observed species of fastidious organisms. The reference identification was determined by 16S rRNA gene sequencing. Overall, the VITEK® 2 NH correctly identified 96.5% of these isolates, including 10.2% low discrimination with the correct species listed. Misidentifications occurred at 2.7% and no identifications occurred at 0.8%.

For 8.01, 9.01, and 9.02 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 NH identification card was evaluated using 371 clinical and stock isolates of both commonly and rarely observed species of fastidious organisms. The reference identification was determined by 16S rRNA gene sequencing. Overall, the VITEK® 2 NH correctly identified 95.7% of these isolates, including 10.5% low discrimination with the correct species listed. Misidentifications occurred at 3.2% and no identifications occurred at 1.1%.

*Data on file at bioMérieux, Inc.

ORGANISMS IDENTIFIED

Claims are for all software users unless otherwise stated.

- *Actinobacillus ureae*
- *Aggregatibacter actinomycetemcomitans*
- *Aggregatibacter aphrophilus*
- *Aggregatibacter segnis*
- *Campylobacter coli*
- *Campylobacter fetus* ssp. *fetus*
- *Campylobacter jejuni* ssp. *jejuni*
- *Capnocytophaga* spp.
- *Cardiobacterium hominis*
- *Eikenella corrodens*
- *Gardnerella vaginalis*
- *Haemophilus haemolyticus*
- *Haemophilus influenzae*
- *Haemophilus parahaemolyticus*
- *Haemophilus parainfluenzae*
- *Kingella denitrificans*
- *Kingella kingae*

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Additional Claims For 8.01 or Higher Software Users

- *Actinobacillus pleuropneumoniae*
- *Actinobacillus suis*
- *Histophilus somni*
- *Moraxella (Neisseria) ovis*
- *Neisseria weaveri*
- *Riemerella anatipestifer*

SUPPLEMENTAL TESTS

Table 14: NH Supplemental Tests

Abbreviation	Test Name	Description	Comments	Reference
For 7.01 or Higher Software Users				
25C	GROWTH AT 25degC	Ability of certain species to grow at 25°C.	N/A	14, 17
42C	GROWTH AT 42degC	Ability of certain species to grow at 42°C.	N/A	17
AGAR 35	GROWTH AT 35C (NUTRIENT AGAR)	Ability of certain species to grow at 35°C on nutrient agar.	N/A	16, 17
CAT	CATALASE	Colony placed on a drop of hydrogen peroxide produces gas bubbles. The bacteria that contain cytochrome enzyme are catalase positive.	N/A	11, 12, 14, 17
COCCI	COCCUS SHAPE	Coccus (round) shape of the bacterial cell examined with gram stain.	N/A	11, 14, 17
DNAse	DNAse	Ability of certain species to produce DNAse resulting in the degradation of DNA.	N/A	12, 16, 17
ESCULIN	ESCULIN hydrolysis	Hydrolysis of esculin forms esculetin that produces a black pigment in the presence of iron salts.	N/A	11, 17, 21
HEMO-horse	Horse blood hemolysis	Certain species possess hemolysins that give a transparent zone around the colonies on blood-based agars.	Hemolysis on horse blood is used as a differential test in the identification of *Haemophilus* spp.	17
HIP	HIPPURATE hydrolysis	Hydrolysis of sodium hippurate releases glycine that produces a blue colored product after addition of ninhydrin.	Of the *Campylobacter* species, only *Campylobacter jejuni* gives a positive hippurate reaction.	17

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Abbreviation	Test Name	Description	Comments	Reference
IND	INDOLE	Ability of certain species to split indole from tryptophan detected by a colored product revealed with a specific reagent (e.g., Kovacs, Ehrlichs, DMAC reagents).	N/A	11, 17, 22
MOB	MOTILITY	Test for motility using hanging drop procedure or wet mount.	Bacterial motility can be observed by placing a drop of bacterial suspension on a slide and viewing it under a microscope.	17, 22
NO3	NITRATE REDUCTION	Test for the ability to reduce nitrate to nitrite or nitrogen gas.	N/A	11, 12, 17, 22
ONPG	BETA-GALACTOSIDASE	Presence of beta-galactosidase cleaves o-nitrophenol-beta-D-galactopyranoside to produce a yellow colored product.	N/A	10, 11, 16, 17
OX	OXIDASE	Detection of the presence of cytochrome C.	N/A	11, 17, 19
THAYER M.	Thayer Martin Agar	Growth on selective medium used in the differentiation of *Neisseria* spp.	Thayer Martin Agar, New York City Agar, or Chocolate Polyvitex Agar with VCAT can be used for this test.	11, 17
UREASE	Urease	Hydrolysis of urea releases ammonia resulting in alkalinization of the medium observed with a pH indicator (e.g., red color formation in the presence of phenol red).	Some tests also appear on the NH card but are recommended as supplemental tests since results of conventional macromethods often differ from rapid commercial micromethods.	11, 17
V FACTOR	V FACTOR (NAD) REQUIREMENT	NAD required for growth.	N/A	16, 17, 19
X FACTOR	X FACTOR (HEMIN) REQUIREMENT	Hemin required for growth.	N/A	11, 17, 19
dFRUCTOSE	D-FRUCTOSE acidification	Acidification of carbon source observed with pH indicator (e.g., phenol red, bromcresol purple).	Some tests also appear on the NH card but are recommended as supplemental tests since results of conventional macromethods often differ from rapid commercial micromethods.	8, 9, 11, 12, 14, 17, 19, 21, 22
dGALACTOSE	D-GALACTOSE acidification			
dGLUCOSE	D-GLUCOSE acidification			
LACTOSE	LACTOSE acidification			
dMALTOSE	D-MALTOSE acidification			
dMANNITOL	D-MANNITOL acidification			
dMANNOSE	D-MANNOSE acidification			
dRAFFINOSE	D-RAFFINOSE acidification			
SACCHAROSE	SACCHAROSE/SUCROSE acidification			
dTREHALOSE	D-TREHALOSE acidification			
For 8.01 or Higher Software Users				
22C	GROWTH AT 22degC	Ability of certain species to grow at 22 °C.	N/A	12
A-HEM	ALPHA HEMOLYSIS	Certain species produce incomplete hemolysis resulting in a green zone around colonies on blood based media.	N/A	3

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Abbreviation	Test Name	Description	Comments	Reference
MICROAER.	MICROAEROBIC GROWTH	Growth requiring a lower concentration of oxygen than is present in the atmosphere.	N/A	14
NO2	NITRITE REDUCTION	Test for the ability to reduce nitrite to nitrogen gas (NO2), nitrate to nitrite, and/or nitrogen gas from nitrate (NO3 to N2).	N/A	12
YELLOW	YELLOW PIGMENT	Ability of certain species to produce yellow pigmented colonies on non-differential media.	N/A	3
dXYLOSE	D-XYLOSE acidifcation	Acidification of carbon source observed with pH indicator (e.g., phenol red, bromcresol purple).	Some tests also appear on the NH card but are recommended as supplemental tests since results of conventional macromethods often differ from rapid commercial micromethods.	8, 9, 11, 12, 14, 17, 19, 21, 22

REFERENCES

1. Balows A, Hausler WJ, Herrmann KL, Isenberg HD, Shadomy HJ, editors. Manual of Clinical Microbiology, 5th ed. American Society for Microbiology, Washington, D.C. 1991.
2. Balows A, Truper HG, Dworkin M, Harder W, and Schleifer K-H, editors. The Prokaryotes - a Handbook on the Biology of Bacteria: Exophysiologoy, Isolation, Identification, Applications, 2nd ed., Volume II. Springer-Verlag, New York 1992.
3. Brenner DJ, Krieg NR, Staley JT, Garrity GM, editors. *Bergey's Manual of Systematic Bacteriology*, 2nd Edition. Springer, New York, NY. 2005.
4. Clinical and Laboratory Standards Institute, M50-A, Quality Control for Commercial Microbial Identification Systems; Approved Guideline, Vol. 28 No. 23.
5. Clinical Laboratory Improvement Amendments of 1988. 42 U.S.C 263a. PL 100-578. 1988.
6. Difco Manual Dehydrated Culture Media and Reagents for Microbiology. 10th ed. 1984.
7. Holmes B, Costas M, On SL, Vandamme P, Falsen E, Kersters K. Neisseria weaveri sp. nov. (formerly CDC group M-5), from dog bite wounds of humans. Int J Syst Bacteriol. 1993 Oct; 43(4):687-93.
8. Holt JG, Krieg NR, Sneath PHA, Staley JT, Williams ST, editors. Bergey's Manual of Determinative Bacteriology, 9th ed. Williams and Wilkins, Baltimore 1994.
9. Kilian M, Nicolet J, Biberstein EL, Biochemical and Serological Characterization of Haemophilus pleuropneumoniae (Matthews and Pattison 1961) Shope 1964 and Proposal of a Neotype Strain. Int J Syst Bacteriol. 1978. 28:20-26.
10. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC, editors. Color Atlas and Textbook of Diagnostic Microbiology 4th ed. Lippincott, Philadelphia, PA 1992.
11. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC, editors. Color Atlas and Textbook of Diagnostic Microbiology 5th ed. Lippincott, Philadelphia, PA. 1997.
12. Koneman EW, Allen SD, Janda WM, Procop GW, Schreckenberger PC, Winn WC, Woods GL, editors. Color Atlas and Textbook of Diagnostic Microbiology 6th ed. Lippincott, Philadelphia, PA. 2006.
13. Krieg NR, Holt JG, editors. *Bergey's Manual of Systematic Bacteriology*, 9th ed. Williams and Wilkins, Baltimore 1984.
14. Lindqvist, K., A Neisseria Species Associated with Infectious Keratoconjunctivitis of Sheep: Neisseria ovis nov. spec. *The Journal of Infectious Diseases*. 1960. 106:162-165.
15. Manafi M, Kneifel W, Bascomb S. Fluorogenic and chromogenic substrates used in bacterial diagnostics. Microbiol. 1991; Rev. 55:335- 348.
16. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Yolken RH, editors. Manual of Clinical Microbiology, 7th ed. American Society for Microbiology, Washington, D.C. 1999.
17. Murray PR, Baron EJ, Jorgensen JH, Pfaller MA and Yolken RH, editors. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue — Approved Guideline, 1997.

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20. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.

21. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry and D.W. Warnock. 2011. *Manual of Clinical Microbiology*, 10th edition. American Society for Microbiology, Washington, D.C.

22. Weyant RS, Moss CW, Weaver RE, Hollis DG, Jordon JG, Cook EC, and Daneshvar MI. Identification of Unusual Pathogenic and Gram-Negative Aerobic and Facultatively Anaerobic Bacteria 2nd ed. Williams & Wilkins, Philadelphia, PA 1996.

Use this Instructions for Use with VITEK® 2 Product No. 21346.

INDEX OF SYMBOLS

Symbol	Meaning
REF	Catalog number
IVD	In Vitro Diagnostic Medical Device
	Legal Manufacturer
	Temperature limitation
	Use by date
LOT	Batch code
	Consult Instructions for Use
	Date of manufacture
	Contains sufficient for <n> tests
EC REP	Authorized representative in the European Community
Rx only	For US Only : Caution : US Federal Law restricts this device to sale by or on the order of a licensed practitioner

Instructions for Use provided in the kit or downloadable from www.biomerieux.com/techlib

LIMITED WARRANTY

bioMérieux warrants the performance of the product for its stated intended use provided that all procedures for usage, storage and handling, shelf life (when applicable), and precautions are strictly followed as detailed in the instructions for use (IFU).

Except as expressly set forth above, bioMérieux hereby disclaims all warranties, including any implied warranties of merchantability and fitness for a particular purpose or use, and disclaims all liability, whether direct, indirect or consequential, for any use of the reagent, software, instrument and disposables (the “System”) other than as set forth in the IFU.

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WASTE DISPOSAL

All hazardous waste must be disposed of by following your local inspecting agency's guidelines.

REVISION HISTORY TABLE

Change type categories	
N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user
Note :	Minor typographical, grammar, and formatting changes are not included in the revision history.

Release Date	Part Number	Change Type	Change Summary
2019-03	043902-03	Technical change	Updated for 9.02 software release. Updated sections: • Intended Use • Precautions • Testing of QC Organisms • Organisms Identified
2016-10	043902-02	Technical change	• Updated content to reflect the 8.01 Product Information Manual
2016-05	043902-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• New IFU derived from product chapter in the Product Information Manual • Updated Limited Warranty section • Updated with RX only information

BIOMÉRIEUX, the BIOMÉRIEUX logo, VITEK, API, Count-TACT, chromID, DensiCHEK and bioLiaison are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

This product may be protected by one or more patents, see: http://www.biomerieux-usa.com/patents.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

Any other name or trademark is the property of its respective owner.

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℞ only IVD

INTENDED USE

These Instructions for Use correspond to the VITEK® 2 Systems 7.01 or higher software. If you are not using VITEK® 2 Systems 7.01 or higher software, please refer to the VITEK® 2 Systems Product Information that you received with your current software version.

The VITEK® 2 Anaerobic and Corynebacteria identification card (ANC) is intended for use with VITEK® 2 Systems for the automated identification of most clinically significant anaerobic organisms and *Corynebacterium* species. The VITEK® 2 ANC identification card is a single-use disposable. For a list of claimed species, see the Organisms Identified section.

DESCRIPTION

The ANC card is based on established biochemical methods and newly developed substrates. There are 36 biochemical tests measuring carbon source utilization and enzymatic activities. Final results are available in approximately six hours.

For a list of well contents, see the ANC Well Contents table.

Table 1: ANC Well Contents

Well	Test	Mnemonic	Amount/Well
4	D-GALACTOSE	dGAL	0.3 mg
5	Leucine ARYLAMIDASE	LeuA	0.023 mg
6	ELLMAN	ELLM	0.03 mg
7	Phenylalanine ARYLAMIDASE	PheA	0.026 mg
8	L-Proline ARYLAMIDASE	ProA	0.023 mg
10	L-Pyrrolidonyl-ARYLAMIDASE	PyrA	0.018 mg
11	D-CELLOBIOSE	dCEL	0.3 mg
13	Tyrosine ARYLAMIDASE	TyrA	0.0279 mg
15	Ala-Phe-Pro-ARYLAMIDASE	APPA	0.038 mg
18	D-GLUCOSE	dGLU	0.3 mg
20	D-MANNOSE	dMNE	0.3 mg
22	D-MALTOSE	dMAL	0.3 mg
28	SACCHAROSE/SUCROSE	SAC	0.3 mg
30	ARBUTIN	ARB	0.1875 mg
33	*N*-ACETYL-D-GLUCOSAMINE	NAG	0.3 mg
34	5-Bromo-4-chloro-3-indoxyl-beta-glucoside	BGLUi	0.006 mg
36	UREASE	URE	0.15 mg
37	5-Bromo-4-chloro-3-indoxyl-beta-glucuronide	BGURi	0.006 mg
39	BETA-GALACTOPYRANOSIDASE Indoxyl	BGALi	0.006 mg
41	ALPHA-ARABINOSIDASE	AARA	0.0324 mg
42	5-Bromo-4-chloro-3-indoxyl-alpha-galactoside	AGALi	0.006 mg
43	BETA-MANNOSIDASE	BMAN	0.036 mg
44	ARGININE GP	ARG	0.15 mg
45	PYRUVATE	PVATE	0.15 mg
51	MALTOTRIOSE	MTE	0.3 mg

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Note: Other well numbers between 1 and 64 not designated in this table are empty.

PRECAUTIONS

Note: For industry customers that need assistance on selecting the correct VITEK® 2 identification card, please refer to the VITEK® 2 Compact Instrument User Manual chapter, "Guidance to Select a VITEK® 2 Identification Card."

- For *In Vitro* Diagnostic Use Only.
- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For professional use only.
- Suspensions not within the appropriate zone on the VITEK® 2 DENSICHEK™ Plus or the VITEK® 2 DENSICHEK™ may compromise card performance.
- Do not use the card after the expiration date shown on the package liner.
- Store the card unopened in the package liner. Do not use the card if the protective package liner is damaged or if no desiccant is present.
- Allow the card to come to room temperature before opening the package liner.
- Do not use powdered gloves. Powder may interfere with the optics.
- Use of culture media other than the recommended types must be validated by the customer laboratory for acceptable performance.
- A Gram stain should be performed to determine an organism's Gram reaction and morphology prior to selecting the identification card to inoculate.
- The card performs as intended only when used in conjunction with VITEK® 2 Systems, following the instructions contained in these Instructions for Use.
- **Do not use glass test tubes**. Use clear plastic (polystyrene) tubes only. Variation exists among test tubes of standard diameter. Carefully place the tube into the cassette. If resistance is encountered, discard and try another tube that does not require pressure to insert.
- Prior to inoculation, inspect cards for tape tears or damage to the tape and discard any that are suspect. Check the saline level in the tubes after the cassette has been processed to ensure proper filling of card.
 - VITEK® 2 60 or VITEK® 2 XL: Eject improperly filled cards.
 - VITEK® 2 Compact: Do not load improperly filled cards.
- Give special consideration to specimen source and patient drug or antimicrobic regimen.
- Interpretation of test results requires the judgment and skill of a person knowledgeable in microbial identification testing. Additional testing may be required. (See the Supplemental Tests section.)
- Do not clean saline dispenser with chemical agents. The use of chemical agents may impact card performance.

Warning: All patient specimens, microbial cultures, and inoculated VITEK® 2 cards, along with associated materials, are potentially infectious and should be treated with universal precautions.[17,18]

Warning: All hazardous waste must be disposed of by following your local inspecting agency's guidelines.

STORAGE CONDITIONS

Upon receipt, store VITEK® 2 ANC cards unopened in their original package liner at 2°C to 8°C.

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SPECIMEN PREPARATION

For specimen preparation information, see the Culture Requirements Table.

Table 2: Culture Requirements

VITEK® 2 Card	Media	Age of Culture[1]	Incubation Conditions	Inoculum Density	Dilution for AST	Age of Suspension Before Loading Instrument
ANC	Corynebacteria: CBA[2] CNA TSAB TSAHB	Corynebacteria: 18 to 24 hours	Corynebacteria: 35°C to 37°C CO_2 or non-CO_2	2.70 to 3.30 McFarland Standard	N/A[3]	≤ 30 minutes
	Anaerobes: CBA[2] CDC[2] BRU CHBA TSAB TSAHB	Anaerobes: 18 to 72 hours	Anaerobes: 35°C to 37ºC anaerobic environment			
	Gram-Positive Anaerobes ONLY: CNA CDC PEA PEA					

[1]Cultures with scant or poor growth may give unidentified or incorrect results even when the Age of Culture requirements are met.

[2]These media were used in the identification product database development and will give optimal performance.

[3]N/A = not applicable

Culture Requirements Table — Media Abbreviations

BRU = Brucella Agar with 5% Sheep Blood, Hemin, and Vitamin K

CBA = Columbia Blood Agar with 5% Sheep Blood

CDC = CDC Anaerobe Agar with 5% Sheep Blood

CDC PEA = CDC Blood Agar with PEA

CHBA = Columbia Horse Blood Agar

CNA = Columbia CNA Agar with 5% Sheep Blood

PEA = Phenylethyl Alcohol Agar with 5% Sheep Blood

TSAB = Trypticase Soy Agar with 5% Sheep Blood

TSAHB = Trypticase Soy Agar with 5% Horse Blood

Materials

When used with VITEK® 2 instrumentation, the ANC card is a complete system for routine identification testing of most clinically significant anaerobic organisms and *Corynebacterium* species.

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Optional accessories:

- Adjustable volume saline dispenser
- Loops
- Pre-dispensed saline test tubes (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- Test tube caps
- Vortex

Procedure

Warning: Failure to follow instructions and recommendations provided in this section for performing laboratory tasks may cause erroneous or delayed results.

For product-specific information, see the Culture Requirements table.

Note: Prepare the inoculum from a pure culture, according to good laboratory practices. In case of mixed cultures, a re-isolation step is required. It is recommended that a purity check plate be done to ensure that a pure culture was used for testing.

1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.
2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).
3. Use a sterile stick or swab to transfer a sufficient number of morphologically similar colonies to the saline tube prepared in step 2. Prepare a homogenous organism suspension with a density equivalent to a McFarland No. 2.70 to 3.30 using a calibrated VITEK® 2 DENSICHEK™ Plus or VITEK® 2 DENSICHEK™.
 Note: Age of suspension must not exceed 30 minutes before inoculating card.
4. Place the suspension tube and ANC card in the cassette.
5. Refer to the appropriate Instrument User Manual for instructions on data entry and how to load the cassette into the instrument.
6. In addition to the internal tests included on the card, three offline tests are required in the ANC ID algorithm. The offline tests selected for use in the ANC ID product are Gram stain, morphology, and aerotolerance. The ANC Offline Test results can be entered at the Smart Carrier Station (VITEK® 2 60 or VITEK® 2 XL only) or the workstation.

Table 3: ANC Offline Tests

Test Name	Test	Result	Definition
AERO	Aerotolerance	–	Anaerobe
		+	Aerobe
		?	Facultative
GRAM	Gram stain results	–	Gram Negative
		+	Gram Positive
		?	Gram Variable

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7. Follow your local inspecting agency's guidelines for disposal of hazardous waste.

RESULTS

Identification Analytical Techniques

VITEK® 2 Systems identify an organism by using a methodology based on the characteristics of the data and knowledge about the organism and reactions being analyzed. Sufficient data have been collected from known strains to estimate the typical reactions of the claimed species to a set of discriminating biochemicals. If a unique identification pattern is not recognized, a list of possible organisms is given, or the strain is determined to be outside the scope of the database.

The printed lab report contains suggestions for any supplemental tests necessary to complete the identification. If the tests are not sufficient to complete the identification, then standard microbiology references and literature should be consulted.

Certain species may belong to a slashline (mixed) taxa identification. This occurs when the biopattern is the same for the taxa listed. Supplemental tests may be used to separate slashline taxa. The species in the ANC Slashline Taxa table belong to the ANC slashline taxa.

Table 4: ANC Slashline Taxa

Slashline Name	Species Belonging to the Slashline
For 7.01, 8.01, and 9.01 Software Users	
Clostridium group	*Clostridium innocuum* *Clostridium limosum* *Clostridium novyi*
For 9.02 Software Users	
Clostridium group	*Clostridium innocuum* *Hathewaya limosa* (formerly known as *Clostridium limosum*) *Clostridium novyi* *Clostridium sporogenes* *Colstridium subterminale*

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Table 5: Identification Card Qualifying Messages

ID Message Confidence Level	Choices	% Probability	Comments
Excellent	1	96 to 99	N/A
Very Good	1	93 to 95	N/A
Good	1	89 to 92	N/A
Acceptable	1	85 to 88	N/A
Low Discrimination	2 to 3	Sum of choices = 100; after resolution to one choice, percent probability reflects the number associated with selected choice.	Two to three taxa exhibit same biopattern. Separate by supplemental testing.
Inconclusive or Unidentified Organism	> 3 or 0	N/A	Either > 3 taxa exhibit same biopattern or Very atypical biopattern. Does not correspond to any taxon in the database. Check Gram stain and purity.

PERCENT PROBABILITY

As part of the identification process, the software compares the test set of reactions to the expected set of reactions of each organism, or organism group, that can be identified by the product. A quantitative value, the percent probability, is calculated and relates to how well the observed reactions compare to the typical reactions of each organism. A perfect match between the test reaction pattern and the unique reaction pattern of a single organism, or organism group, would provide a percent probability of 99. When a perfect match is not obtained, it is still possible for the reaction pattern to be sufficiently close to that of an expected reaction pattern such that a clear decision can be provided about the organism identification. The range of percent probabilities in the one-choice case is 85 to 99. Values closer to 99 indicate a closer match to the typical pattern for the given organism.

When the reaction pattern is not sufficient to discriminate between two to three organisms, the percent probabilities reflect this ambiguity. The reported probability values indicate, relatively, the order in which the reaction pattern best corresponds to the listed possibilities. The order does not, however, suggest that the pattern match to one of the possible identifications is clearly superior to another. The probability characteristic of an overall sum of 100 is retained through the calculation process. After resolution to one choice, the probability characteristic of the single choice is retained.

ADDITIONAL INFORMATION ON LAB REPORT

Supplemental test — External (offline) test that allows the user to resolve a slashline or Low Discrimination identification. Numbers in parentheses indicate percent positive reaction for the species/test listed.

Contraindicating test — Test result that is unusual for a reported taxon.

Table 6: Notes Associated with Certain Taxa

Taxa	Note
For 7.01 or Higher Software Users	
Actinomyces israelii	*Actinomyces israelii* represents a complex of two closely related species, *A. israelii* and *A. gerencseriae* (formerly known as *A. israelii serotype II*).
Corynebacterium diphtheriae	Critical pathogen. The species identified may have significance to patient or sample outcome and can be stopped for review.
For 9.02 Software Users	

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Taxa	Note
Rhodococcus hoagii *Turicella otitidis*	*Rhodococcus hoagii* and *Turicella otitidis* can be differentiated from each other by pigment and body site. *R. hoagii* produces a salmon color and *T. otitidis* is only found in the region of the ear. In addition, *R. hoagii* is most commonly found as a veterinary pathogen while *T. otitidis* is found as human pathogen.

Notes Associated with an Improperly Filled Card or with a Negative Profile (Biopattern)

- For the case where the time between two readings is greater than 40 minutes: “CARD ERROR — Missing data.”
- For the case where there is a negative profile: “Organism with low reactivity biopattern — please check viability.”
- When a biopattern is calculated for an unknown organism that is completely negative or consists of both negative tests and tests that fall within the uncertainty zone, the identification call will be “Non or low reactive biopattern.”

The following non-reactive species could potentially trigger this note if a test was atypical or fell within the uncertainty zone:

- *Clostridium clostridioforme*
- *Fusobacterium nucleatum*
- *Fusobacterium mortiferum*

QUALITY CONTROL

Quality control organisms and their expected results are listed in the VITEK® 2 ANC Quality Control Tables. Process these according to the procedure for test isolates outlined in this document.

Certification Statement

This is to certify that bioMérieux complies with ISO 13485 and FDA Quality System Regulation (QSR) requirements for design, development, and manufacture of microbial identification systems.

Frequency of Testing

Currently, it is recommended that you use your most stringent inspecting agency’s guidelines for frequency of identification product testing.

Common practice is to perform QC upon receipt of shipment of the test kits. Reactions must follow Instructions for Use results.

If the results do not meet the criteria, subculture for purity and repeat the test. If discrepant results are repeated, perform an alternate identification method and contact bioMérieux.

Testing and Storage of QC Organisms

1. Rehydrate the organism according to the manufacturer's instructions.
2. *Corynebacterium*: Use Columbia Blood agar with 5% sheep blood (CBA) and incubate at 35°C to 37°C in non-CO_2 aerobic conditions. Incubate for 18 to 24 hours or until sufficient growth is obtained.
3. Anaerobes: Use Columbia Blood agar with 5% sheep blood and incubate at 35°C to 37°C in anaerobic conditions for 18 to 24 hours or until sufficient growth is obtained.
4. Check for purity. Perform second subculture for testing.
5. *Corynebacterium*: Use Columbia Blood agar with 5% sheep blood and incubate at 35°C to 37°C in non-CO_2 aerobic conditions. Incubate for 18 to 24 hours.
6. Anaerobes: Use Columbia Blood agar with 5% sheep blood and incubate at 35°C to 37°C in anaerobic conditions for 18 to 24 hours.

Short-Term Storage Conditions - *Corynebacterium*

1. Streak to a CBA plate or slant.
2. Incubate at 35°C to 37°C in non-CO_2 aerobic conditions. Incubate for 18 to 24 hours.
3. Refrigerate at 2°C to 8°C for up to five days.
4. Subculture to CBA. Incubate at 35°C to 37°C in non-CO_2 aerobic conditions for 18 to 24 hours. Use for QC.

Short-Term Storage Conditions - Anaerobes

1. Streak to a CBA plate or slant.
2. Incubate at 35°C to 37°C in anaerobic conditions for 18 to 24 hours or until sufficient growth is obtained.

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3. Store at room temperature in anaerobic conditions for up to five days.
4. Subculture to CBA. Incubate at 35°C to 37°C in anaerobic conditions for 18 to 24 hours. Use for QC.

Long-Term Storage Conditions

1. Make a heavy suspension in Tryptic Soy Broth (TSB) with 15% glycerol.
2. Freeze at -70°C.
3. Subculture to CBA twice before running QC.

Note: Avoid repeated thawing and refreezing by either freezing in single-use aliquots or removing a small portion of frozen organism preparation with a sterile applicator stick.

STREAMLINED QUALITY CONTROL

Note: Industrial Use Only laboratories should perform quality control following the Streamlined Quality Control section. No additional testing is required for these users.

As there are no substrates that are consistently sensitive to degradation during shipping conditions, streamlined quality control may be conducted by testing two strains: one that is mostly positive and the other, which is mostly negative for reactions on ANC. (See ANC Quality Control tables for more details).

COMPREHENSIVE QUALITY CONTROL

Customers who do not qualify for streamlined quality control testing are required to perform comprehensive quality control testing, which entails demonstration of a positive and negative reaction for each substrate of an identification product.[4]

In order to qualify initially for streamlined quality control testing, the CLSI® M50-A standard requires that the user perform and document either of the following:[3]

- Verification testing to show that performance is equivalent to the manufacturer's claims.
- Comprehensive quality control testing of at least three lots over at least three different seasons.

Refer to the complete CLSI® M50-A standard for information regarding continued qualification and further details of requirements and responsibilities for both the user and the manufacturer related to streamlined quality control testing.

ANC Quality Control Tables:

***Clostridium septicum* ATCC® 12464™** (for streamlined or comprehensive quality control)

***Bacteroides ovatus* ATCC® BAA-1296™** (for streamlined or comprehensive quality control)

***Bacteroides vulgatus* ATCC® 8482™** (for comprehensive quality control)

***Clostridium perfringens* ATCC® 13124™** (for comprehensive quality control)

***Corynebacterium striatum* ATCC® BAA-1293™** (for comprehensive quality control)

***Parabacteroides distasonis* ATCC® BAA-1295™** (for comprehensive quality control)

For 7.01, 8.01, and 9.01 Software Users

***Clostridium sordellii* ATCC® 9714™** (for comprehensive quality control)

For 9.02 Software Users

***Paeniclostridium sordellii* ATCC® 9714™** (formerly known as ***Clostridium sordellii* ATCC® 9714™**) for comprehensive quality control

The ANC card typically identifies the quality control organisms as one-choice or within a low discrimination or slashline identification. However, strains are chosen for reaction performance over identification performance. Therefore, an unidentified or misidentified result may occur when all expected quality control reactions are correct.

Table 7: QC Organism: *Clostridium septicum* ATCC® 12464™ (for streamlined or comprehensive quality control)

dGAL	-	dCEL	-	SAC	-	BGALi	+	MTE	-	PHOS	-	GRAM	+
LeuA	-	TyrA	-	ARB	-	AARA	v	ESC	-	IARA	-	MORPH	-
ELLM	-	APPA	-	NAG	-	AGALi	-	BdFUC	+	dRIB2	-	AERO	-
PheA	-	dGLU	-	BGLUi	-	BMAN	-	BNAGi	-	OPS	+		
ProA	-	dMNE	-	URE	-	ARG	-	AMANi	v	AARAF	-		
PyrA	v	dMAL	-	BGURi	-	PVATE	-	AIFUC	-	dXYL	-		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

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Table 8: QC Organism: *Bacteroides ovatus* ATCC® BAA-1296™ (for streamlined or comprehensive quality control)

dGAL	+	dCEL	+	SAC	v	BGALi	+	MTE	+	PHOS	v	GRAM	-
LeuA	-	TyrA	-	ARB	v	AARA	+	ESC	+	IARA	+	MORPH	-
ELLM	+	APPA	+	NAG	+	AGALi	+	BdFUC	v	dRIB2	+	AERO	-
PheA	-	dGLU	+	BGLUi	v	BMAN	v	BNAGi	-	OPS	v		
ProA	-	dMNE	+	URE	-	ARG	-	AMANi	v[1]	AARAF	+		
PyrA	-	dMAL	+	BGURi	v	PVATE	v	AIFUC	v	dXYL	v		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive.

[1]Reaction is mostly positive although occasional negative reaction may occur.

Table 9: QC Organism: *Bacteroides vulgatus* ATCC® 8482™ (for comprehensive quality control)

dGAL	v	dCEL	v	SAC	v	BGALi	v	MTE	v	PHOS	+	GRAM	–
LeuA	v	TyrA	v	ARB	v	AARA	+	ESC	v	IARA	v	MORPH	–
ELLM	+	APPA	+	NAG	v	AGALi	v	BdFUC	+	dRIB2	v	AERO	–
PheA	v	dGLU	v	BGLUi	v	BMAN	v[1]	BNAGi	v[1]	OPS	v		
ProA	v	dMNE	v	URE	v	ARG	v	AMANi	-	AARAF	+		
PyrA	v	dMAL	v	BGURi	v[1]	PVATE	v	AIFUC	+	dXYL	+		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

[1]Reaction is mostly positive although occasional negative reaction may occur.

Table 10: QC Organism: *Clostridium perfringens* ATCC® 13124™ (for comprehensive quality control)

dGAL	v	dCEL	v	SAC	+	BGALi	v	MTE	+	PHOS	+	GRAM	+
LeuA	v	TyrA	v	ARB	v	AARA	v	ESC	v	IARA	v	MORPH	–
ELLM	–	APPA	–	NAG	v	AGALi	+	BdFUC	v	dRIB2	+	AERO	–
PheA	v	dGLU	v	BGLUi	v	BMAN	v	BNAGi	v	OPS	+		
ProA	v	dMNE	v	URE	v	ARG	+	AMANi	v	AARAF	–		
PyrA	+	dMAL	+	BGURi	v	PVATE	–	AIFUC	v	dXYL	v		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 11: QC Organism: *Clostridium sordellii ATCC® 9714™ (for comprehensive quality control)**

dGAL	–	dCEL	v	SAC	–	BGALi	–	MTE	v	PHOS	v	GRAM	+
LeuA	v	TyrA	v	ARB	v	AARA	–	ESC	–	IARA	v	MORPH	–
ELLM	v	APPA	v	NAG	v	AGALi	–	BdFUC	–	dRIB2	v	AERO	–
PheA	v	dGLU	v	BGLUi	–	BMAN	–	BNAGi	v	OPS	–		
ProA	+	dMNE	-	URE	+	ARG	v	AMANi	–	AARAF	v		
PyrA	v	dMAL	v	BGURi	v	PVATE	v	AIFUC	v	dXYL	–		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

*For 7.01, 8.01, and 9.01 software users, *Clostridium sordellii*

*For 9.02 software users, *Paeniclostridium sordellii*, formerly known as *Clostridium sordellii*

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Table 12: QC Organism: *Corynebacterium striatum* ATCC® BAA-1293™ (for comprehensive quality control)

dGAL	+	dCEL	–	SAC	+	BGALi	–	MTE	–	PHOS	–	GRAM	+
LeuA	+	TyrA	+	ARB	–	AARA	v	ESC	v	IARA	–	MORPH	–
ELLM	v	APPA	v	NAG	–	AGALi	v	BdFUC	–	dRIB2	–	AERO	+
PheA	v	dGLU	+	BGLUi	v	BMAN	v	BNAGi	v	OPS	v		
ProA	+	dMNE	+	URE	v	ARG	v	AMANi	v	AARAF	v		
PyrA	–	dMAL	–	BGURi	v	PVATE	+	AIFUC	v	dXYL	v		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 13: QC Organism: *Parabacteroides distasonis* ATCC® BAA-1295™ (for comprehensive quality control)

dGAL	v	dCEL	v	SAC	v	BGALi	v	MTE	v	PHOS	v	GRAM	–
LeuA	v	TyrA	v	ARB	+	AARA	v	ESC	+	IARA	v	MORPH	–
ELLM	v	APPA	v	NAG	+	AGALi	v	BdFUC	v	dRIB2	v	AERO	–
PheA	v[1]	dGLU	v	BGLUi	+	BMAN	v	BNAGi	v	OPS	v		
ProA	v	dMNE	v	URE	v	ARG	v	AMANi	v	AARAF	v		
PyrA	+	dMAL	v	BGURi	–	PVATE	v	AIFUC	–	dXYL	v		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

[1]Reaction is mostly positive although occasional negative reaction may occur.

LIMITATIONS

The VITEK® 2 ANC card cannot be used with a direct clinical specimen or sample or other sources containing mixed flora.

Newly described or rare species may not be included in the ANC database. Selected species will be added as strains become available.

Warning: Testing of unclaimed species may result in an unidentified result or a misidentification.

PERFORMANCE CHARACTERISTICS

For 7.01 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 ANC identification card was evaluated using 365 clinical and stock isolates of both commonly and rarely observed species. The reference identification was determined using 16S rRNA gene sequencing. Overall, the VITEK® 2 ANC correctly identified 93.9% of these isolates, including 9.0% low discrimination with the correct species listed. Misidentifications occurred at 5.8% and no identifications occurred at 0.3%.

For 8.01 and 9.01 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 ANC identification card was evaluated using 365 clinical and stock isolates of both commonly and rarely observed species. The reference identification was determined using 16S rRNA gene sequencing. Overall, the VITEK® 2 ANC correctly identified 94.0% of these isolates, including 9.0% low discrimination with the correct species listed. Misidentifications occurred at 5.8% and no identifications occurred at 0.3%.

For 9.02 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 ANC identification card was evaluated using 365 clinical and stock isolates of both commonly and rarely observed species. The reference identification was determined using 16S rRNA gene sequencing. Overall, the VITEK® 2 ANC correctly identified 94.3% of these isolates, including 10.4% low discrimination with the correct species listed. Misidentifications occurred at 5.5% and no identifications occurred at 0.3%.

*Data on file at bioMérieux, Inc.

ORGANISMS IDENTIFIED

Claims are for all software users unless otherwise stated.

- *Actinobaculum schaalii*
- *Actinomyces bovis*
- *Actinomyces israelii*
- *Actinomyces meyeri*

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Taxonomy Changes For 8.01 or Higher Software Users

- *Terrisporobacter glycolicus* (formerly known as *Clostridium glycolicum*)

Additional Claims and Taxonomy Changes For 9.02 Software Users

- *Actinotignum schaalii* (formerly known as*Actinobaculum schaalii*)
- *Actinomyces canis*
- *Bacteroides pyogenes*
- *Bacteroides xylanisolvens*
- *Corynebacterium glucuronolyticum*
- *Hathewaya histolytica* (formerly known as *Clostridium histolyticum*)
- *Hathewaya limosa* (formerly known as *Clostridium limosum*)
- *Paeniclostridium sordellii* (formerly known as *Clostridium sordellii*)
- *Paraclostridium bifermentans* (formerly known as *Clostridium bifermentans*)
- *Porphyromonoas asaccharolytica*
- *Porphyromonas gulae*
- *Porphyromonas uenonus*
- *Prevotella nanceiensis*
- *Prevotella salivae*

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SUPPLEMENTAL TESTS

Table 14: ANC Supplemental Tests

Abbreviation	Test Name	Description	Comments	Reference
For 7.01 or Higher Software Users				
AFUC	Alpha-Fucosidase	Presence of enzyme cleaves substrate generating detectable leaving group (e.g., p-nitrophenol, methyl umbelliferone, beta-naphthylamide, p-nitroaniline, 7-amino-methyl-coumarin).	Presence of enzyme is indicated by generation of a colored or fluorescent product, or a non-colored product that forms color upon addition of a specific reagent.	9, 16, 25, 28
BNAG	BETA-N-ACETYL-GLUCOSAMINIDASE	Presence of respective enzyme cleaves substrate generating detectable leaving group (e.g., p-nitrophenol, methyl umbelliferone, beta-naphthylamide, beta-naphthol, p-nitroaniline, 7-amidomethyl-coumarin).	Presence of enzyme is indicated by generation of a colored or fluorescent product, or a noncolored product that forms color upon addition of a specific reagent.	16, 27
Branch.flt	BRANCHING FILAMENTS	Appearance of branching filaments on microscopic examination.	N/A	8, 9, 15
CAT	CATALASE	Colony placed on a drop of hydrogen peroxide produces gas bubbles. The bacteria that contain cytochrome enzyme are catalase positive.	N/A	7, 9, 15, 16, 23, 25
ESCULIN	ESCULIN hydrolysis	Hydrolysis of esculin forms esculetin, which produces a black pigment in the presence of iron salts.	Some tests also appear on the ANC card but are recommended as supplemental tests since results of conventional macromethods may differ from rapid commercial micromethods.	5, 8, 9, 15, 16, 28, 30
GELATIN	Gelatin hydrolysis	Mediated by a gelatinase enzyme. A positive reaction is evidenced by liquefaction of the gelatin substrate.	N/A	7, 9, 15, 16,20, 28, 30
IND	INDOLE	Ability of certain species to split indole from tryptophan detected by a colored product revealed with a specific reagent (e.g., Kovacs, Ehrlichs, DMAC reagents).	N/A	7, 15, 16, 24
LECITHIN.	LECITHINASE	A precipitate surrounding the colony on egg yolk agar indicates lecithinase activity of alphatoxin produced by the organism.	N/A	7, 9, 16
LIP	LIPASE	An iridescent pearly sheen on the surface of the colony on egg yolk agar indicates lipase activity.	N/A	7, 9, 16

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Abbreviation	Test Name	Description	Comments	Reference
LIPOPHILY	LIPOPHILY	Lipophilic	Enhanced growth in the presence of lipids in the culture medium.	16
NO3	NITRATE REDUCTION	Test for the ability to reduce nitrate to nitrite or nitrogen gas.	N/A	7, 9, 15, 16
PAL	ALKALINE PHOSPHATASE	Presence of enzyme cleaves substrate generating detectable leaving group.	Presence of enzyme is indicated by generation of a colored or fluorescent product, or a non-colored product that forms color upon addition of a specific reagent.	16
PIGMENT	Pigment	Ability of certain species to produce pigmented colonies on nondifferential media.	N/A	16, 28, 29, 30
Point.ends	POINTED ENDS	Appearance of thin gram-negative rods with pointed ends is a microscopic characteristic for *Fuscobacterium nucleatum*.	N/A	9
PYRAZINAM.	PYRAZINAMIDASE	Test for the activity of the enzyme pyrazinamidase, which hydrolyses pyrazinamide to pyrazinoic acid.	N/A	16
SPOR	SPORE	Microscopic examination for spores. Phase-contrast microscopy is recommended.	N/A	9, 16
UREASE	Urease	Hydrolysis of urea releases ammonia resulting in alkalinization of the medium observed with a pH indicator (e.g., red color formation in the presence of phenol red).	Some tests also appear on the ANC card but are recommended as supplemental tests since results of conventional macromethods often differ from rapid commercial micromethods.	7, 8, 9, 15, 16
lARABINOSE	L-ARABINOSE acidification	Anaerobe: PRAS* method *Corynebacterium*: acid from carbohydrate. Acidification of carbon source observed with pH indicator (e.g., phenol red, bromcresol purple).	Some tests also appear on the ANC card but are recommended as supplemental tests since results of conventional macromethods often differ from rapid commercial micromethods.	7, 8, 9, 15, 16, 19, 20, 23, 25, 27, 28, 30
dCELLOB	DCELLOBIOSE acidification			
dFRUCTOSE	D-FRUCTOSE acidification			
dGALACTOSE	DGALACTOSE acidification			
dGLUCOSE	D-GLUCOSE acidification			
LACTOSE	Lactose acidification			
dMALTOSE	D-MALTOSE acidification			
dMANNITOL	D-MANNITOL acidification			
dMANNOSE	D-MANNOSE acidification			
dRAFFINOSE	dRAFFINOSE acidification			
lRHAMNOSE	L-RHAMNOSE acidification			
dRIBOSE	dRIBOSE acidification			
SACCHAROSE	SACCHAROSE/SUCROSE acidification			
SALICIN	SALICIN			
STARCHac	STARCH acidification			
dTREHALOSE	dTREHALOSE acidification			
XYL	XYLOSE			
XYLAN	XYLAN acidification			
For 9.02 Software Users				

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Abbreviation	Test Name	Description	Comments	Reference
CASEIN	CASEIN hydrolysis	Mediated by an exoenzyme casease. A positive reaction is evidenced by liquefaction of the casein substrate.	N/A	21, 26

*PRAS: Pre-reduced anaerobically sterilized media

REFERENCES

1. Balows, A., W.J. Hausler, K.L. Herrmann, H.D. Isenberg, H.J. Shadomy (ed.). 1991. *Manual of Clinical Microbiology*, 5th ed. American Society for Microbiology, Washington, D.C.
2. Balows, A., H.G. Truper, M. Dworkin, W. Harder, and K-H. Schleifer (ed.). 1992. *The Prokaryotes- a Handbook on the Biology of Bacteria: Exophysiology, Isolation, Identification, Applications*, 2nd ed., Volume II. Springer-Verlag, New York.
3. Clinical and Laboratory Standards Institute, M50-A, Quality Control for Commercial Microbial Identification Systems; Approved Guideline, Vol. 28 No. 23.
4. Clinical Laboratory Improvement Amendments of 1988. 42 U.S.C 263a. PL 100-578. 1988.
5. De Vos, P., Garrity, G., Jones, D., Krieg, N., Ludwig, W., Rainey, F., Schleifer, K., Whitman, W. Bergey's' Manual of Systematic Bacteriology, Second Edition, Volume Three the Firmicutes, Springer Publishing Dordrecht, 2009.
6. Difco Manual Dehydrated Culture Media and Reagents for Microbiology. 10th ed. 1984
7. Holdeman, L.V., Cato, E.P., Moore, W.E.C., Anaerobe Laboratory Manual 4th ed. Blacksburg VA.: VPI Anaerobe Laboratory, 1977.
8. Holt, J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, S.T. Williams (ed.). 1994. *Bergey's Manual of Determinative Bacteriology*, 9th ed. Williams and Wilkins, Baltimore.
9. Jousimies-Somer, H., P. Summanen, D. M. Citron, E.J. Baron, H.M. Wexler, S.M. Finegold (ed) 2002 *Wadsworth - KTL Anaerobic Bacteriology Manual*, 6th ed. Star Publishing Belmont California.
10. Koneman, E.W., S.D. Allen, W. M. Janda, P.C. Schreckenberger, W.C. Winn (ed.). 1992. *Color Atlas and Textbook of Diagnostic Microbiology*, 4th ed. Lippincott Publishing Company, Philadelphia,PA.
11. Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, W.C. Winn (ed.). *Color Atlas and Textbook of Diagnostic Microbiology*, 5th ed. 1997 Lippincott, Philadelphia, New York.
12. Krieg, N.R., and J.G. Holt. (ed.) 1984. *Bergey's Manual of Systematic Bacteriology*, 9th ed. Williams and Wilkins, Baltimore.
13. Manafi, M., W. Kneifel, and S. Bascomb. 1991. Fluorogenic and chromogenic substrates used in bacterial diagnostics. Microbiol. Rev. 55:335-348.
14. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Yolken (ed.) 1999. *Manual of Clinical Microbiology*, 7th ed. American Society for Microbiology, Washington, D.C.
15. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Yolken (ed.) 2003. *Manual of Clinical Microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.
16. Murray, P.R., Baron, E., Jorgensen, J.H., Landry, M.L., Pfaller, M.A. *Manual of Clinical Microbiology*. 9th ed. Washington D.C.: ASM Press, 2007.
17. National Committee for Clinical Laboratory Standards, M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue* - Approved Guideline, 1997.
18. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, *Biosafety in Microbiological and Biomedical Laboratories*, 1988.
19. Winn, W.C. Jr, Allen, S.D., Janda, W.M., Koneman, E.W., Propcop, G.W., Schreckenberger, P.C., Woods, G.L. *Koneman's Color Atlas and Textbook of Diagnostic Microbiology 6th* ed. Philadelphia: Lippincott Williams & Wilkins, 2006.
20. Alauzet, C., Mory, F., Carlier, J.-P., Marchandin, H., Jumas-Bilak, E., & Lozniewski, A., 2007. *Prevotella nanceiensi*s sp. nov., isolated from human clinical samples. Int J Syst Evol Microbiol 57, 2216-2220.
21. Benno, Y., Watabe, J. & Mitsuoka, T. 1983. *Bacteroides pyogenes* sp. nov., *Bacteroides suis* sp. nov., and *Bacteroides helcogenes* sp. nov., New Species from Abscesses and Feces of Pigs. System. Appl. Microbiol. 4, p. 396-407.
22. Chassard, C., Delmas, E., Lawson, P.A. & Bernalier-Donadille, A. 2008. *Bacteroides xylanisolvens* sp. nov., a Xylan Degrading Bacterium Isolated from Human Faeces. Int J Syst Evol Microbiol 58, 1008-1013.
23. Devriese, L.A., Riegel, P., Hommez, J., Vaneechoutte, M., De Baere, T. & Haesebrouck, F. 2000. Identification of *Corynebacterium glucuronolyticum* Strains from the Urogenital Tract of Humans and Pigs. J. Clin. Microbiol. 38, p. 4657-4659.
24. Finegold, S.M., Vaisanen, M.-L., Rautio, M., Eerola, E., Summanen, P., Molitoris, D., Song, Y., Liu, C., & Jousimies-Somer, H. 2004. *Porphyromonas uenonis* sp. nov., a Pathogen for Humans Distinct from *P. asaccharolytica* and *P. endodontalis*. J. Clin. Microb, 42, p. 5298-5301.

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25. Fournier, D., Mouton, C., Lapierre, P., Kato, T., Okuda, K., & Ménard, C. 2001. Porphyromonas gulae sp. nov., an anaerobic, Gram-negative coccobacillus from the gingival sulcus of various animal hosts. Int J Syst Evol Microbiol 51, 1179-1189.
26. Holdeman, L.V. and Johnson, J.L. 1977. *Bacteroides disiens* sp. nov. and *Bacteroides bivius* sp. nov. from Human Clinical Infections. Int. J. Syst. Bacteriol. 27, p. 337-345.
27. Hoyles, L., Falsen, E., Foster, G., Pascual, C., Greko, C. & Collins, M.D. 2000. *Actinomyces canis* sp. nov., isolated from dogs. Int J Syst Evol Microbiol 50, 1547-1551.
28. Jorgensen, J.H., M. A Phaller, K.C. Carroll, G. Funke, M.L. Landry, S.S. Richter, and D.W. Warnock. 2015. Manual of Clinical Microbiology, 11th edition. American Society for Microbiology, Washington, D.C.
29. Kedlaya, I., Ing, M.B. & Wong, S.S. 2001. *Rhodococcus equi* Infections and Immunocompetent Hosts: Case Report and Review. Clinical Infectious Disease, 32, p. e39-e47.
30. Love, D.N., Johnson, J.L., Jones, R.F., Bailey, M. & Calverley, A. 1986. *Bacteroides tectum* sp. nov. and Characteristics of Other Nonpigmented *Bacteroides* Isolates from Soft-Tissue Infections from Cats and Dogs. Int J Syst Bacteriol. 36, 123-128.

Use this Instructions for Use with VITEK® 2 Product No. 21347.

INDEX OF SYMBOLS

Symbol	Meaning
REF	Catalog number
IVD	In Vitro Diagnostic Medical Device
	Legal Manufacturer
	Temperature limitation
	Use by date
LOT	Batch code
	Consult Instructions for Use
	Date of manufacture
	Contains sufficient for <n> tests
EC REP	Authorized representative in the European Community
Rx only	For US Only : Caution : US Federal Law restricts this device to sale by or on the order of a licensed practitioner

Instructions for Use provided in the kit or downloadable from www.biomerieux.com/techlib

LIMITED WARRANTY

bioMérieux warrants the performance of the product for its stated intended use provided that all procedures for usage, storage and handling, shelf life (when applicable), and precautions are strictly followed as detailed in the instructions for use (IFU).

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Except as expressly set forth above, bioMérieux hereby disclaims all warranties, including any implied warranties of merchantability and fitness for a particular purpose or use, and disclaims all liability, whether direct, indirect or consequential, for any use of the reagent, software, instrument and disposables (the “System”) other than as set forth in the IFU.

REVISION HISTORY TABLE

Change type categories	
N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user
Note :	Minor typographical, grammar, and formatting changes are not included in the revision history.

Release Date	Part Number	Change Type	Change Summary
2019-03	043907-03	Technical change	Updated for 9.02 software release. Updated sections: • Intended Use • Precautions • Culture Requirements • Additional Information on Lab Report • Testing of QC Organisms • Performance Characteristics • Organisms Identified • References
2016-10	043907-02	Technical change	• Updated content to reflect the 8.01 Product Information Manual
2016-05	043907-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• New IFU derived from product chapter in the Product Information Manual • Updated Limited Warranty section • Updated with RX only information

BIOMÉRIEUX, the BIOMÉRIEUX logo, VITEK, API, Count-TACT, chromID, DensiCHEK and bioLiaison are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

This product may be protected by one or more patents, see: http://www.biomerieux-usa.com/patents.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

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Rx only IVD

INTENDED USE

These Instructions for Use correspond to the VITEK® 2 Systems 7.01 or higher software. If you are not using VITEK® 2 Systems 7.01 or higher software, please refer to the VITEK® 2 Systems Product Information that you received with your current software version.

The VITEK® 2 Yeast identification card (YST) is intended for use with VITEK® 2 Systems for the automated identification of most clinically significant yeasts and yeast-like organisms. The VITEK® 2 YST identification card is a single-use disposable. For a list of claimed species, see the Organisms Identified section.

DESCRIPTION

The YST card is based on established biochemical methods[1,6,8,10,11] and newly developed substrates. There are 46 biochemical tests measuring carbon source utilization, nitrogen source utilization, and enzymatic activities. Final results are available in approximately 18 hours.

For a list of well contents, see the YST Well Contents table.

Table 1: YST Well Contents

Well	Test	Mnemonic	Amount/Well
3	L-Lysine-ARYLAMIDASE	LysA	0.0228 mg
4	L-MALATE assimilation	IMLTa	0.15 mg
5	Leucine-ARYLAMIDASE	LeuA	0.0234 mg
7	ARGININE	ARG	0.15 mg
10	ERYTHRITOL assimilation	ERYa	0.3 mg
12	GLYCEROL assimilation	GLYLa	0.16 µL
13	Tyrosine ARYLAMIDASE	TyrA	0.0276 mg
14	BETA-N-ACETYL-GLUCOSAMINIDASE	BNAG	0.0408 mg
15	ARBUTIN assimilation	ARBa	0.3 mg
18	AMYGDALIN assimilation	AMYa	0.3 mg
19	D-GALACTOSE assimilation	dGALa	0.3 mg
20	GENTOBIOSE assimilation	GENa	0.3 mg
21	D-GLUCOSE assimilation	dGLUa	0.3 mg
23	LACTOSE assimilation	LACa	0.96 mg
24	METHYL-A-D-GLUCOPYRANOSIDE assimilation	MAdGa	0.3 mg
26	D-CELLOBIOSE assimilation	dCELa	0.3 mg
27	GAMMA-GLUTAMYL-TRANSFERASE	GGT	0.0228 mg
28	D-MALTOSE assimilation	dMALa	0.3 mg
29	D-RAFFINOSE assimilation	dRAFa	0.3 mg
30	PNP-N-acetyl-BD-galactosaminidase 1	NAGA1	0.0306 mg
32	D-MANNOSE assimilation	dMNEa	0.3 mg
33	D-MELIBIOSE assimilation	dMELa	0.3 mg
34	D-MELEZITOSE assimilation	dMLZa	0.3 mg
38	L-SORBOSE assimilation	ISBEa	0.3 mg

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Note: Other well numbers between 1 and 64 not designated in this table are empty.

PRECAUTIONS

Note: For industry customers that need assistance on selecting the correct VITEK® 2 identification card, please refer to the VITEK® 2 Compact Instrument User Manual chapter, "Guidance to Select a VITEK® 2 Identification Card."

- For *In Vitro* Diagnostic Use Only.
- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For professional use only.
- Suspensions not within the appropriate zone on the VITEK® 2 DENSICHEK™ Plus or the VITEK® 2 DENSICHEK™ may compromise card performance.
- Do not use the card after the expiration date shown on the package liner.
- Store the card unopened in the package liner. Do not use the card if the protective package liner is damaged or if no desiccant is present.
- Allow the card to come to room temperature before opening the package liner.
- Do not use powdered gloves. Powder may interfere with the optics.
- Use of culture media other than the recommended types must be validated by the customer laboratory for acceptable performance.
- A Gram stain should be performed to determine an organism’s Gram reaction and morphology prior to selecting the identification card to inoculate.
- The card performs as intended only when used in conjunction with VITEK® 2 Systems, following the instructions contained in these Instructions for Use.
- **Do not use glass test tubes**. Use clear plastic (polystyrene) tubes only. Variation exists among test tubes of standard diameter. Carefully place the tube into the cassette. If resistance is encountered, discard and try another tube that does not require pressure to insert.

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- Prior to inoculation, inspect cards for tape tears or damage to the tape and discard any that are suspect. Check the saline level in the tubes after the cassette has been processed to ensure proper filling of card.
 - VITEK® 2 60 or VITEK® 2 XL: Eject improperly filled cards.
 - VITEK® 2 Compact: Do not load improperly filled cards.
- Give special consideration to specimen source and patient drug or antimicrobic regimen.
- Interpretation of test results requires the judgment and skill of a person knowledgeable in microbial identification testing. Additional testing may be required. (See the Supplemental Tests section.)
- Do not clean saline dispenser with chemical agents. The use of chemical agents may impact card performance.

Warning: All patient specimens, microbial cultures, and inoculated VITEK® 2 cards, along with associated materials, are potentially infectious and should be treated with universal precautions.[13,17]

Warning: All hazardous waste must be disposed of by following your local inspecting agency's guidelines.

STORAGE CONDITIONS

Upon receipt, store VITEK® 2 YST cards unopened in their original package liner at 2°C to 8°C.

SPECIMEN PREPARATION

For specimen preparation information, see the Culture Requirements Table.

Table 2: Culture Requirements Table

VITEK® 2 Card	Media	Age of Culture[1]	Incubation Conditions	Inoculum Density	Dilution for AST	Age of Suspension Before Loading Instrument
YST	SDA[2] SDA-E[2] TSAB[2] CBA IMA TSA TSAHB CHBA CID CPS ID	18 to 72 hours	30°C to 37°C aerobic, non-CO_2 (or 25°C to 30°C for species that do not tolerate 30°C to 37°C)	1.80 to 2.20 McFarland Standard	N/A[3]	≤ 30 minutes
YST and AST-YST Pair	SDA SDA-E TSAB CBA TSA CHBA CID CPS ID	18 to 72 hours	35°C to 37°C aerobic, non-CO_2	1.80 to 2.20 McFarland Standard	280 µL in 3.0 mL saline	≤ 30 minutes

[1]Cultures with scant or poor growth may give unidentified or incorrect results even when the Age of Culture requirements are met.

[2]These media were used in the identification product database development and will give optimal performance.

[3]N/A = not applicable

Culture Requirements Table — Media Abbreviations

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CBA = Columbia Blood Agar with 5% Sheep Blood

CHBA = Columbia Horse Blood Agar

CID = chromID™Candida (Candida ID2 agar)

CPS ID = chromID™ CPS (CPS ID agar)

IMA = Inhibitory Mould Agar

SDA = Sabouraud Dextrose Agar

SDA-E = Sabouraud Dextrose Agar (Emmons)

TSA = Trypticase Soy Agar

TSAB = Trypticase Soy Agar with 5% Sheep Blood

TSAHB = Trypticase Soy Agar with 5% Horse Blood

Materials

When used with VITEK® 2 instrumentation, the YST card is a complete system for routine identification testing of most clinically significant yeasts and yeast-like organisms.

Required materials are:

- VITEK® 2 YST Card
- DENSICHEK™ Plus Kit or VITEK® DENSICHEK® Kit
- DENSICHEK™ Plus Standards Kit or DENSICHEK® Standards Kit
- VITEK® 2 Cassette
- Sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- 12 mm x 75 mm clear plastic (polystyrene) disposable test tubes
- Sterile sticks or swabs
- Appropriate agar medium (see Culture Requirements table).

Optional accessories:

- Adjustable volume saline dispenser
- Loops
- Pre-dispensed saline test tubes (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- Test tube caps
- Vortex

Procedure

Warning: Failure to follow instructions and recommendations provided in this section for performing laboratory tasks may cause erroneous or delayed results.

For product-specific information, see the Culture Requirements table.

Note: Prepare the inoculum from a pure culture, according to good laboratory practices. In case of mixed cultures, a re-isolation step is required. It is recommended that a purity check plate be done to ensure that a pure culture was used for testing.

1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.
2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).
3. Use a sterile stick or swab to transfer a sufficient number of morphologically similar colonies to the saline tube prepared in step 2. Prepare the homogenous organism suspension with a density equivalent to a McFarland No. 1.80 to 2.20 using a calibrated VITEK® 2 DENSICHEK™ Plus or VITEK® 2 DENSICHEK™.
 Note: Filamentous species may pick up small amounts of glucose from the isolation media. This has the potential for causing false positive reactions. Avoid scraping or rubbing the agar when preparing the organism suspension. For strains

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that do not readily form a smooth suspension in saline, it is recommended that a wet swab be used to make the suspension. Do not rub the agar surface when preparing a suspension using a wet swab.

Note: Age of suspension must not exceed 30 minutes before inoculating card.

4. Place the suspension tube and YST card in the cassette.
5. Refer to the appropriate Instrument User Manual for instructions on data entry and how to load the cassette into the instrument.
6. Follow your local inspecting agency's guidelines for disposal of hazardous waste.

RESULTS

Identification Analytical Techniques

VITEK® 2 Systems identify an organism by using a methodology based on the characteristics of the data and knowledge about the organism and reactions being analyzed. Sufficient data have been collected from known strains to estimate the typical reactions of the claimed species to a set of discriminating biochemicals. If a unique identification pattern is not recognized, a list of possible organisms is given, or the strain is determined to be outside the scope of the database.

The printed lab report contains suggestions for any supplemental tests necessary to complete the identification. If the tests are not sufficient to complete the identification, then standard microbiology references and literature should be consulted.

Certain species may belong to a slashline (mixed) taxa identification. This occurs when the biopattern is the same for the taxa listed. Supplemental tests may be used to separate slashline taxa. The species in the Slashline (Mixed) Taxa Identification table belong to the YST slashline taxa.

Table 3: Slashline (Mixed) Taxa Identification

Slashline Name	Species Belonging to the Slashline
C. inconspicua/C. lambica	*Candida inconspicua* *Candida lambica*
Kloeckera spp.	*Kloeckera apiculata* *Kloeckera apis* *Kloeckera japonica*
*Rhodotorula glutinis/mucilaginosa/(Crypto. laurentii)**	*Rhodotorula glutinis* *Rhodotorula mucilaginosa*

*This is also a pseudoslashline.

Certain species may belong to a pseudoslashline (mixed) taxa identification. A pseudoslashline indicates a rare isolate or rare occurrence of the same biopattern. Supplemental tests may be used to separate pseudoslashline taxa. The species in Pseudoslashline Taxa Table belong to the pseudoslashline taxa.

Table 4: Pseudoslashline Taxa

Pseudoslashline name	Species Belonging to Pseudoslashline
Candida sake/(C. famata/C. lipolytica)	*Candida famata* *Candida lipolytica*
*Rhodotorula glutinis/mucilaginosa/(Crypto. laurentii)**	*Cryptococcus laurentii*

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Table 5: Identification Card Qualifying Messages

ID Message Confidence Level	Choices	% Probability	Comments
Excellent	1	96 to 99	N/A
Very Good	1	93 to 95	N/A
Good	1	89 to 92	N/A
Acceptable	1	85 to 88	N/A
Low Discrimination	2 to 3	Sum of choices = 100; after resolution to one choice, percent probability reflects the number associated with selected choice.	Two to three taxa exhibit same biopattern. Separate by supplemental testing.
Inconclusive or Unidentified Organism	> 3 or 0	N/A	Either > 3 taxa exhibit same biopattern or Very atypical biopattern. Does not correspond to any taxon in the database. Check Gram stain and purity.

PERCENT PROBABILITY

As part of the identification process, the software compares the test set of reactions to the expected set of reactions of each organism, or organism group, that can be identified by the product. A quantitative value, the percent probability, is calculated and relates to how well the observed reactions compare to the typical reactions of each organism. A perfect match between the test reaction pattern and the unique reaction pattern of a single organism, or organism group, would provide a percent probability of 99. When a perfect match is not obtained, it is still possible for the reaction pattern to be sufficiently close to that of an expected reaction pattern such that a clear decision can be provided about the organism identification. The range of percent probabilities in the one-choice case is 85 to 99. Values closer to 99 indicate a closer match to the typical pattern for the given organism.

When the reaction pattern is not sufficient to discriminate between two to three organisms, the percent probabilities reflect this ambiguity. The reported probability values indicate, relatively, the order in which the reaction pattern best corresponds to the listed possibilities. The order does not, however, suggest that the pattern match to one of the possible identifications is clearly superior to another. The probability characteristic of an overall sum of 100 is retained through the calculation process. After resolution to one choice, the probability characteristic of the single choice is retained.

ADDITIONAL INFORMATION ON LAB REPORT

Supplemental test — External (offline) test that allows the user to resolve a slashline or Low Discrimination identification. Numbers in parentheses indicate percent positive reaction for the species/test listed.

Contraindicating test — Test result that is unusual for a reported taxon.

Table 6: Notes Associated with Certain Taxa

Taxa	Note			
For 7.01 or Higher Software Users				
Candida krusei	Possibility of *C. inconspicua* or *C. lambica*. Isolates of these infrequent species may be misidentified as *C. krusei*; in order to rule them out, perform the following tests:			
		HYPH/PH	dGLUf	dXYLOSEa
	C. inconspicua	-	-	-
	C. krusei	+	+	-
	C. lambica	+	+	+

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Taxa	Note
Rhodotorula glutinis/ mucilaginosa *Cryptococcus laurentii*	Possibility of *Cryptococcus albidus*
Geotrichum klebahnii	Possibility of *Geotrichum candidum*
Cryptococcus neoformans	Critical pathogen The species identified may have significance to patient or sample outcome and can be stopped for review.
For 8.01 or Higher Software Users	
Candida glabrata	*Candida nivariensis* and *Candida bracarensis* have similar morphology and biochemical utilization patterns compared to *Candida glabrata*, and they exhibit similar multi-resistance to anti-fungals. These three species may be differentiated from each other by molecular methods, such as MALDI-TOF, since phenotypic tests do not differentiate them. Retrospective studies of culture collections showed that up to 0.1% of isolates previously identified as *C. glabrata* were strains of *C. nivariensis* and that 0.2-2.2% of isolates previously identified as *C. glabrata* were strains of *C. bracarensis*.[3,9]
For 9.02 Software Users	
Candida parapsilosis	*Candida metapsilosis* and *Candida orthopsilosis* have similar morphology and biochemical utilization patterns compared to *Candida parapsilosis*. Although data are limited, antifungal susceptibility testing of both species show a similar profile to that of *C. parapsilosis*. These three species may be differentiated from each other by molecular methods, such as MALDI-TOF, since phenotypic tests do not differentiate them. *C.orthopsilosis* and *C. metapsilosis* are rare clinical isolates but are reported emerging pathogens associated with candidemia.

Notes Associated with an Improperly Filled Card or with a Negative Profile (Biopattern)

- For the case where the time between two readings is greater than 40 minutes: “CARD ERROR — Missing data.”
- For the case where there is a negative profile: “Organism with low reactivity biopattern — please check viability.”
- When a biopattern is calculated for an unknown organism that is completely negative or consists of both negative tests and tests that fall within the uncertainty zone, the identification call will be “Non or low reactive biopattern.”

The following non-reactive species could potentially trigger this note if a test was atypical or fell within the uncertainty zone:

- *Candida sake*
- *Candida zeylanoides*
- *Malassezia furfur*
- *Malassezia pachydermatis*

For 7.01 Software Users

- *Zygosaccharomyces bailii*

For 8.01 or Higher Software Users

- *Zygosaccharomyces* species

QUALITY CONTROL

Quality control organisms and their expected results are listed in the VITEK® 2 YST Quality Control Tables. Process these according to the procedure for test isolates outlined in this document.

Note: *Staphylococcus epidermidis* ATCC® 12228™ needs to be tested at a McFarland Standard No. 0.5 to 0.63. All other QC strains are tested at a McFarland Standard No. 1.80 to 2.20.

Certification Statement

This is to certify that bioMérieux complies with ISO 13485 and FDA Quality System Regulation (QSR) requirements for design, development, and manufacture of microbial identification systems.

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Frequency of Testing

Currently, it is recommended that you use your most stringent inspecting agency's guidelines for frequency of identification product testing.

Common practice is to perform QC upon receipt of shipment of the test kits. Reactions must follow Instructions for Use results.

If the results do not meet the criteria, subculture for purity and repeat the test. If discrepant results are repeated, perform an alternate identification method and contact bioMérieux.

Testing and Storage of QC Organisms

1. Rehydrate the organism according to the manufacturer's instructions.
2. Yeast: Streak onto Sabouraud dextrose agar (SDA) or SDA (Emmons) and incubate aerobically at 35°C to 37°C for 18 to 24 hours, or until sufficient growth is achieved, except for:
 - *Prototheca wickerhamii* ATCC® 16529™, *Zygosaccharomyces parabailii* ATCC® MYA-4549™, and *Kloeckera japonica* ATCC® 58370™,which are incubated at 28°C to 30°C.
 - *Sporobolomyces salmonicolor* ATCC® MYA-4550™, which is incubated at 25°C to 27°C.
3. Bacteria: Use trypticase soy agar with 5% sheep blood (TSAB). Incubate aerobically at 35°C to 37°C for 18 to 24 hours.
4. Check for purity. Perform second subculture for testing.
5. Yeast: Streak onto SDA or SDA (Emmons) and incubate aerobically at 35°C to 37°C for 18 to 24 hours, or until sufficient growth is achieved, except for:
 - *Prototheca wickerhamii* ATCC® 16529™, *Zygosaccharomyces parabailii* ATCC® MYA-4549™, and *Kloeckera japonica* ATCC® 58370™, which are incubated at 28°C to 30°C.
 - *Sporobolomyces salmonicolor* ATCC® MYA-4550™, which is incubated at 25°C to 27°C.
6. Bacteria: Use trypticase soy agar with 5% sheep blood (TSAB). Incubate aerobically at 35°C to 37°C for 18 to 24 hours.

Short-Term Storage Conditions

1. Streak to a SDA or SDA (Emmons) plate or slant for yeast QC organisms and TSAB for bacteria QC organisms.
2. Incubate for 24 hours at the appropriate temperature.
3. Refrigerate at 2°C to 8°C for up to one week.
4. Subculture once as described above and use for QC.

Long-Term Storage Conditions

1. Make a heavy suspension in Tryptic Soy Broth (TSB) with 15% glycerol.
2. Freeze at -70°C.
3. Subculture to appropriate medium twice before running QC.

Note: Avoid repeated thawing and refreezing by either freezing in single-use aliquots or removing a small portion of frozen organism preparation with a sterile applicator stick.

STREAMLINED QUALITY CONTROL

Note: Industrial Use Only laboratories should perform quality control following the Streamlined Quality Control section. No additional testing is required for these users.

Streamlined quality control may be used to confirm acceptable performance of the YST card after shipping/storage. This methodology may be performed with the YST card by following the instructions for quality control testing as described in the YST Instructions for Use and meeting the criteria stated in CLSI® M50-A Quality Control for Commercial Microbial Identification Systems.

Testing may be conducted using *Candida albicans* ATCC® 14053™ and evaluating the performance of the NAGA1 well. Testing at bioMérieux, Inc. has shown that the NAGA1 well is the most labile well on the YST card and *Candida albicans* ATCC® 14053™ is the most sensitive strain for detecting degradation of this well with a false negative reaction. (See YST Quality Control table for more details).

COMPREHENSIVE QUALITY CONTROL

Customers who do not qualify for streamlined quality control testing are required to perform comprehensive quality control testing, which entails demonstration of a positive and negative reaction for each substrate of an identification product.[5]

In order to qualify initially for streamlined quality control testing, the CLSI® M50-A standard requires that the user perform and document either of the following:[4]

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- Verification testing to show that performance is equivalent to the manufacturer's claims.
- Comprehensive quality control testing of at least three lots over at least three different seasons.

Refer to the complete CLSI® M50-A standard for information regarding continued qualification and further details of requirements and responsibilities for both the user and the manufacturer related to streamlined quality control testing.

YST Quality Control Tables:

***Candida albicans* ATCC® 14053™** (for streamlined or comprehensive quality control)

***Candida glabrata* ATCC® MYA-2950™** (for comprehensive quality control)

***Candida lusitaniae* ATCC® 34449™** (for comprehensive quality control)

***Candida utilis* ATCC® 9950™** (for comprehensive quality control)

***Kloeckera japonica* ATCC® 58370™** (for comprehensive quality control)

***Prototheca wickerhamii* ATCC® 16529™** (for comprehensive quality control)

***Sporobolomyces salmonicolor* ATCC® MYA-4550™** (for comprehensive quality control)

***Trichosporon mucoides* ATCC® 204094™** (for comprehensive quality control)

***Oligella ureolytica* ATCC® 43534™** (for comprehensive quality control)

***Staphylococcus epidermidis* ATCC® 12228™** (for comprehensive quality control)

***Zygosaccharomyces parabailii*ATCC® MYA-4549™** (for comprehensive quality control)

Note: The comprehensive quality control strain *Zygosaccharomyces bailii* ATCC® MYA-4549™ had a taxonomic update to *Zygosaccharomyces parabailii* ATCC® MYA-4549™. This strain identifies as *Zygosaccharomyces bailii* in the 7.01 software and as *Zygosaccharomyces* species in the 8.01 software.

The YST card typically identifies the quality control organisms as one-choice or within a low discrimination or slashline identification. However, strains are chosen for reaction performance over identification performance. Therefore, an unidentified or misidentified result may occur when all expected quality control reactions are correct.

Note: The YST card uses unclaimed taxa for quality control testing. These strains will give an unidentified or a misidentified result.

Table 7: QC Organism: *Candida albicans* ATCC® 14053™ (for streamlined or comprehensive quality control)

LysA	-	ARBa	-	GGT	v	IRHAa	-	NO3a	-	CITa	+
IMLTa	+	AMYa	v	dMALa	+	XLTa	+	IARAa	v	GRTas	v
LeuA	+	dGALa	+	dRAFa	-	dSORa	+	dGATa	v	IPROa	+
ARG	+	GENa	-	NAGA1	+	SACa	+	ESC	-	2KGa	+
ERYa	-	dGLUa	+	dMNEa	+	URE	-	IGLTa	+	NAGa	+
GLYLa	v	LACa	-	dMELa	-	AGLU	+	dXYLa	+	dGNTa	+
TyrA	v	MAdGa	+	dMLZa	-	dTURa	+	LATa	+		
BNAG	-	dCELa	-	ISBEa	-	dTREa	+	ACEa	+		

+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Note: The NAGA1 well is used for streamlined quality control.

Table 8: QC Organism: *Candida glabrata* ATCC® MYA-2950™ (for comprehensive quality control)

LysA	–	ARBa	–	GGT	–	IRHAa	–	NO3a	v	CITa	–
IMLTa	–	AMYa	v	dMALa	–	XLTa	v	IARAa	–	GRTas	–
LeuA	v	dGALa	–	dRAFa	v	dSORa	–	dGATa	–	IPROa	v
ARG	–	GENa	v	NAGA1	–	SACa	–	ESC	–	2KGa	–
ERYa	–	dGLUa	v	dMNEa	v	URE	–	IGLTa	v	NAGa	–
GLYLa	–	LACa	–	dMELa	–	AGLU	–	dXYLa	–	dGNTa	v
TyrA	–	MAdGa	–	dMLZa	v	dTURa	–	LATa	–		

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+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 9: QC Organism: *Candida lusitaniae* ATCC® 34449™ (for comprehensive quality control)

LysA	v	ARBa	+	GGT	v	lRHAa	+	NO3a	v	CITa	+
lMLTa	+	AMYa	+	dMALa	v	XLTa	v	lARAa	–	GRTas	v
LeuA	+	dGALa	v	dRAFa	–	dSORa	+	dGATa	v	lPROa	+
ARG	v	GENa	+	NAGA1	v	SACa	v	ESC	+	2KGa	v
ERYa	v	dGLUa	v	dMNEa	v	URE	v	lGLTa	+	NAGa	+
GLYLa	v	LACa	v	dMELa	v	AGLU	v	dXYLa	v	dGNTa	v
TyrA	v	MAdGa	v	dMLZa	v	dTURa	+	LATa	v		
BNAG	v	dCELa	+	ISBEa	+	dTREa	v	ACEa	+		

+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 10: QC Organism: QC Organism *Candida utilis* ATCC® 9950™ (for comprehensive quality control)

LysA	v	ARBa	v	GGT	v	lRHAa	v	NO3a	+	CITa	v
lMLTa	v	AMYa	+	dMALa	v	XLTa	–	lARAa	v	GRTas	v
LeuA	v	dGALa	v	dRAFa	+	dSORa	–	dGATa	v	lPROa	v
ARG	v	GENa	v	NAGA1	–	SACa	+	ESC	v	2KGa	–
ERYa	v	dGLUa	+	dMNEa	+	URE	v	lGLTa	v	NAGa	–
GLYLa	+	LACa	v	dMELa	v	AGLU	v	dXYLa	v	dGNTa	v
TyrA	v	MAdGa	v	dMLZa	+	dTURa	v	LATa	v		
BNAG	v	dCELa	v	ISBEa	v	dTREa	v	ACEa	v		

+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 11: QC Organism: *Kloeckera japonica* ATCC® 58370™ (for comprehensive quality control)

LysA	v	ARBa	v	GGT	v	lRHAa	v	NO3a	v	CITa	v
lMLTa	v	AMYa	v	dMALa	–	XLTa	v	lARAa	v	GRTas	v
LeuA	v	dGALa	v	dRAFa	v	dSORa	v	dGATa	v	lPROa	v
ARG	v	GENa	v	NAGA1	v	SACa	v	ESC	v	2KGa	v
ERYa	v	dGLUa	v	dMNEa	v	URE	v	lGLTa	v	NAGa	v
GLYLa	v	LACa	v	dMELa	v	AGLU	v	dXYLa	v	dGNTa	v
TyrA	v	MAdGa	v	dMLZa	v	dTURa	v	LATa	v		
BNAG	v	dCELa	v	ISBEa	v	dTREa	v	ACEa	–		

+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 12: QC Organism: *Prototheca wickerhamii* ATCC® 16529™ (for comprehensive quality control)

LysA	v	ARBa	v	GGT	–	lRHAa	v	NO3a	–	CITa	–
lMLTa	–	AMYa	–	dMALa	v	XLTa	v	lARAa	v	GRTas	v
LeuA	v	dGALa	v	dRAFa	v	dSORa	v	dGATa	v	lPROa	v
ARG	v	GENa	–	NAGA1	v	SACa	–	ESC	v	2KGa	v
ERYa	v	dGLUa	v	dMNEa	v	URE	v	lGLTa	v	NAGa	v
GLYLa	+	LACa	v	dMELa	v	AGLU	v	dXYLa	v	dGNTa	v

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+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 13: QC Organism: *Sporobolomyces salmonicolor* ATCC® MYA-4550™ (for comprehensive quality control)

LysA	+	ARBa	v	GGT	v	lRHAa	v	NO3a	v	CITa	v
lMLTa	v	AMYa	v	dMALa	v	XLTa	v	lARAa	v	GRTas	v
LeuA	v	dGALa	v	dRAFa	v	dSORa	v	dGATa	v	lPROa	v
ARG	v	GENa	v	NAGA1	v	SACa	v	ESC	v	2KGa	v
ERYa	v	dGLUa	v	dMNEa	v	URE	v	lGLTa	v	NAGa	v
GLYLa	v	LACa	v	dMELa	v	AGLU	v	dXYLa	v	dGNTa	v
TyrA	v	MAdGa	v	dMLZa	v	dTURa	v	LATa	v		
BNAG	v	dCELa	v	ISBEa	v	dTREa	v	ACEa	v		

+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 14: QC Organism: *Trichosporon mucoides* ATCC® 204094™ (for comprehensive quality control)

LysA	v	ARBa	+	GGT	+	lRHAa	+	NO3a	v	CITa	v
lMLTa	v	AMYa	–	dMALa	+	XLTa	+	lARAa	+	GRTas	+
LeuA	v	dGALa	+	dRAFa	+	dSORa	v	dGATa	+	lPROa	v
ARG	+	GENa	+	NAGA1	+	SACa	v	ESC	+	2KGa	+
ERYa	+	dGLUa	v	dMNEa	v	URE	+	lGLTa	v	NAGa	v
GLYLa	v	LACa	+	dMELa	+	AGLU	+	dXYLa	+	dGNTa	+
TyrA	+	MAdGa	+	dMLZa	+	dTURa	v	LATa	+		
BNAG	+	dCELa	+	ISBEa	v	dTREa	v	ACEa	v		

+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 15: QC Organism: *Oligella ureolytica* ATCC® 43534™ (for comprehensive quality control)

LysA	v	ARBa	v	GGT	v	lRHAa	v	NO3a	v	CITa	v
lMLTa	v	AMYa	v	dMALa	v	XLTa	v	lARAa	v	GRTas	v
LeuA	v	dGALa	v	dRAFa	v	dSORa	v	dGATa	v	lPROa	v
ARG	v	GENa	v	NAGA1	v	SACa	v	ESC	v	2KGa	v
ERYa	v	dGLUa	–	dMNEa	–	URE	v	lGLTa	v	NAGa	v
GLYLa	v	LACa	v	dMELa	v	AGLU	v	dXYLa	v	dGNTa	v
TyrA	v	MAdGa	v	dMLZa	v	dTURa	v	LATa	v		
BNAG	v	dCELa	v	ISBEa	v	dTREa	v	ACEa	v		

+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Note: *Oligella ureolytica* is an unclaimed taxa for the YST card.

Table 16: QC Organism: *Staphylococcus epidermidis* ATCC® 12228™ (for comprehensive quality control)

LysA	v	ARBa	v	GGT	v	lRHAa	v	NO3a	v	CITa	v
lMLTa	v	AMYa	v	dMALa	v	XLTa	v	lARAa	v	GRTas	v
LeuA	–	dGALa	v	dRAFa	v	dSORa	v	dGATa	v	lPROa	v
ARG	v	GENa	v	NAGA1	v	SACa	v	ESC	v	2KGa	v

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+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Note: *Staphylococcus epidermidis* is an unclaimed taxa for the YST card.

Table 17: QC Organism: *Zygosaccharomyces parabailii* ATCC® MYA-4549™ (for comprehensive quality control)

LysA	v	ARBa	v	GGT	v	lRHAa	v	NO3a	v	CITa	v
lMLTa	v	AMYa	v	dMALa	v	XLTa	v	lARAa	v	GRTas	v
LeuA	v	dGALa	v	dRAFa	v	dSORa	v	dGATa	v	lPROa	–
ARG	v	GENa	v	NAGA1	v	SACa	v	ESC	v	2KGa	v
ERYa	v	dGLUa	v	dMNEa	v	URE	v	lGLTa	–	NAGa	v
GLYLa	v	LACa	v	dMELa	v	AGLU	v	dXYLa	v	dGNTa	–
TyrA	v	MAdGa	v	dMLZa	v	dTURa	v	LATa	v		
BNAG	v	dCELa	v	ISBEa	v	dTREa	–	ACEa	v		

+ = 95% to100% positive; v = 6% to 94% positive; – = 0% to 5% positive

For 7.01 Software Users

Zygosaccharomyces parabailii ATCC® MYA-4549™ identifies as *Zygosaccharomyces bailii*.

For 8.01 or Higher Software Users
Zygosaccharomyces parabailii ATCC® MYA-4549™ identifies as *Zygosaccharomyces* species.

LIMITATIONS
The VITEK® 2 YST card cannot be used with a direct clinical specimen, sample, or another source containing mixed flora. Any change or modification in the procedure may affect the results.

Newly described or rare species may not be included in the YST database. Selected species will be added as strains become available.

Warning: Testing of unclaimed species may result in an unidentified result or a misidentification.

PERFORMANCE CHARACTERISTICS
For 7.01 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 YST identification card was evaluated using 623 clinical and stock isolates of both commonly and rarely observed species of yeast and yeast-like organisms. The reference identification was determined with API® 20C AUX identification kits. Overall, the VITEK® 2 YST correctly identified 98.2% of the isolates, including 8.1% low discrimination with the correct species listed. Misidentifications occurred at 1.5% and no identifications occurred at 0.3%.

For 8.01 and 9.01 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 YST identification card was evaluated using 621 clinical and stock isolates of both commonly and rarely observed species of yeast and yeast-like organisms. The reference identification was determined with API® 20C AUX identification kits. Overall, the VITEK® 2 YST correctly identified 97.9% of the isolates, including 7.2% low discrimination with the correct species listed. Misidentifications occurred at 1.8% and no identifications occurred at 0.5%.

For 9.02 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 YST identification card was evaluated using 621 clinical and stock isolates of both commonly and rarely observed species of yeast and yeast-like organisms. The reference identification was determined with API® 20C AUX identification kits. Overall, the VITEK® 2 YST correctly identified 97.6% of the isolates, including 6.0% low discrimination with the correct species listed. Misidentifications occurred at 1.9% and no identifications occurred at 0.5%.

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Additional Claims For 8.01 or Higher Software Users

- *Candida auris*
- *Candida ciferrii* (formerly known as *Stephanoascus ciferrii*)
- *Candida duobushaemulonii*
- *Candida haemulonii* var *vulnera*
- *Cryptococcus gattii*
- *Zygosaccharomyces* species (includes *Zygosaccharomyces bailii*; *Zygosaccharomyces bailii* is no longer a single species claim)

SUPPLEMENTAL TESTS

Table 18: YST Supplemental Tests

Abbreviation	Test Name	Description	Comment	Reference
For 7.01 or Higher Software Users				
2KG	2-KETO-D-GLUCONATE	Ability to use 2-keto-D-gluconate as a sole carbon source.	N/A	2, 6, 10
4ASCOSPOR.	4 Ascospores	Microscopic examination for the presence of four ascospores per ascus.	N/A	2, 6, 10
Apic.CELLS	APICULATE CELLS	Microscopic examination for the presence of apiculate (lemon-drop shaped) cells.	N/A	2, 6, 10
Arthro.	Arthroconidia	Microscopic examination for the presence of arthroconidia (fragmentation of hyphae into rectangular cells) on morphology agar (e.g., cornmeal agar).	N/A	2, 6, 10
CAROTENOID	CAROTENOID PIGMENT	Presence of red, pink, or orange pigment on Sabouraud dextrose agar.	N/A	2, 6, 7, 10, 15, 16
dCELLOB.a	D-CELLOBIOSE Assimilation	Ability to use cellobiose as the sole carbon source.	N/A	2, 6, 7, 10, 15, 16
CHLS	Chlamydospores	Microscopic examination for the presence of chlamydospores on morphology agar (e.g., cornmeal agar).	N/A	2, 6, 10
DULCITOLa	DULCITOL Assimilation	Ability to use dulcitol (galactitol) as the sole carbon source.	N/A	2, 6, 10

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Abbreviation	Test Name	Description	Comment	Reference
ERYTHRIT.a	ERYTHRITOL Assimilation	Ability to use erythritol as the sole carbon source.	N/A	2, 6, 10
dGALACT.a	D-GALACTOSE Assimilation	Ability to use galactose as the sole carbon source.	N/A	2, 6, 10, 15
dGALf	D-GALACTOSE Fermentation	Production of gas from fermentation of galactose.	N/A	2, 6, 10
dGLUf	D-GLUCOSE Fermentation	Production of gas from fermentation of glucose.	N/A	2, 6, 10
w/o OIL	GROWTH WITHOUT OIL	Ability to grow on Sabouraud dextrose agar without the addition of a fatty acid source (e.g., olive oil).	N/A	2, 6, 7, 10, 16
HYPH/PH	HYPHAE/ PSEUDOHYHPAE	Microscopic examination for the presence of filaments on morphology agar (e.g., cornmeal agar).	N/A	2, 6, 7, 10, 16
INOSITOLa	myo-INOSITOL Assimilation	Ability to use inositol as the sole carbon source.	N/A	2, 6, 10, 15
NITRATEa	NITRATE Assimilation	Ability to use potassium nitrate as the sole nitrogen source.	N/A	2, 6, 7, 10, 14, 16
LACTOSEa	LACTOSE Assimilation	Ability to use lactose as the sole carbon source.	N/A	2, 6, 10
lARABIN.a	L-ARABINOSE Assimilation	Ability to use arabinose as the sole carbon source.	N/A	2, 6, 10
dMALTOSEa	D-MALTOSE Assimilation	Ability to use maltose as the sole carbon source.	N/A	2, 6, 10
dMALf	D-MALTOSE Fermentation	Production of gas from fermentation of maltose.	N/A	2, 6, 10
dMELIBIO.a	D-MELIBIOSE Assimilation	Ability to use melibiose as the sole carbon source.	N/A	2, 6, 10
OX_Phe	Phenol Oxidase	Ability to produce brown to black pigment from phenol oxidase activity on phenolic substrates (e.g., caffeic acid or birdseed agar).	N/A	12

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Abbreviation	Test Name	Description	Comment	Reference
dRAFFIN.a	D-RAFFINOSE Assimilation	Ability to use raffinose as the sole carbon source.	N/A	2, 6, 7, 10, 16
lRHAMNOSEa	L-RHAMNOSE Assimilation	Ability to use rhamnose as the sole carbon source.	N/A	2, 6, 10
SACCHAR.a	SACCHAROSE/ SUCROSE Assimilation	Ability to use sucrose as the sole carbon source.	N/A	2, 6, 10
SACf	SACCAROSE/ SUCROSE Fermentation	Production of gas from fermentation of sucrose.	N/A	2, 6, 10
SATELLITE	SATELLITE behavior	Formation of satellite colonies on Sabouraud dextrose agar.	N/A	2, 6, 10
Sphe.CELLS	Spherical CELLS	Microscopic examination for the presence of spherical cells.	*Candida famata* can be differentiated from *Candida guilliermondii* by the shape of the cells. *Candida famata* has mostly spherical cells while *Candida guilliermondii* has mostly ovoid cells.	2, 6, 10
SPORANGE	SPORANGE	Microscopic examination for the presence of sporangia.	N/A	11
dTREHAL.a	D-TREHALOSE Assimilation	Ability to use trehalose as the sole carbon source.	N/A	2, 6, 10
dTREf	D-TREHALOSE Fermentation	Production of gas from fermentation of trehalose.	N/A	2, 6, 10
UREASE	Urease	Hydrolysis of urea releases ammonia resulting in alkalinization of the medium observed with a pH indicator (e.g., red color formation in the presence of phenol red).	N/A	2, 6, 10
dXYLOSEa	D-XYLOSE Assimilation	Ability to use xylose as the sole carbon source.	N/A	2, 6, 10
For 8.01 or Higher Software Users				
37C	GROWTH AT 37C	Ability to grow at 37°C.	N/A	15
42C	Growth at 42degC	Ability to grow at 42°C.	N/A	2, 15
Cser.AorD	Capsular Serotype A or D	Agglutination tests for capsular serotype A , D or AD.	N/A	18

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Abbreviation	Test Name	Description	Comment	Reference
Cser.BorC	Capsular Serotype B or C	Agglutination tests for capsular serotype B or C.	N/A	18
GLYCEROLa	Glycerol assim.	Ability to use glycerol as the sole carbon source.	N/A	2, 15
INUa	INULIN assimilation	Ability to use inulin as the sole carbon source.	N/A	15
RAFf	RAFFINOSE fermentation	Production of gas from fermentation of raffinose.	N/A	2, 15

REFERENCES

1. Atlas RA. *Handbook of Microbiological Media.* CRC Press, Ann Arbor. 1993.
2. Barnett JA, Payne RW, Yarrow D, editors. *Yeasts: Characteristics and Identification,* 3rd ed. Cambridge University Press, New York. 2000.
3. Bishop JA, Chase N, Magill SS, Kurtzman CP, Fiandaca MJ, Merz WG. *Candida bracarensis* detected among isolates of *Candida glabrata* by peptide nucleic acid fluorescence in situ hybridization: susceptibility data and documentation of presumed infection. J Clin Microbial. 2008; 46:443-446.
4. Clinical and Laboratory Standards Institute, M50-A, Quality Control for Commercial Microbial Identification Systems; Approved Guideline, Vol. 28 No. 23.
5. Clinical Laboratory Improvement Amendments of 1988. 42 U.S.C 263a. PL 100-578. 1988.
6. Kreger-van Rij NJW, editor. *The yeasts — a taxonomic study,* 3rd ed. Elsevier Science Publishers B.V. Amsterdam. 1984.
7. Kurtzman, CP, JW Fell, T Boekhout, editors. *The Yeasts, a Taxonomic Study,* 5th ed. Elsevier, San Diego, CA. 2011.
8. Larone DH. *Medically Important Fungi — a guide to identification.* 3rd ed. ASM Press. American Society for Microbiology. Washington, D.C. 1995.
9. Lockhart SR, Messer SA, Gherna M, Bishop JA, Merz WG, Pfaller MA, Diekema DJ. Identification of *Candida nivariensis* and *Candida bracarensis* in a large global collection of *Candida glabrata* isolates: comparison to the literature. J Clin Microbial. 2009; 47: 1216-1217.
10. Lodder J. *The Yeasts,* Second Edition. North Holland Publishing Company, Netherlands. 1971.
11. McGinnis MR. *Laboratory Handbook of Medical Mycology,* Academic Press, New York. 1980.
12. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Yolken RH, editors. *Manual of Clinical Microbiology,* 7th Edition. American Society for Microbiology, Washington, D.C. 1999.
13. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue—Approved Guideline. 1997.
14. Pincus DH, Salkin IF, Hurd NH, Levy IL, Kemna MA. Modification of Potassium Nitrate Assimilation Test for Identification of Clinically Important Yeasts. J. Clin. Microbiol. 1988; 26:366-368.
15. Satoh K., K. Makimura, Y.Hasumi, Y.Nishiyama, K. Uchida and H. Yamaguchi. *Candida auris* sp.nov., a novel ascomycetous yeast isolated from the external ear canal of an inpatient in a Japanese hospital. Microbiol Immunol 2009; 53: 41-44.
16. Suh S-O., P. Gujjari, C. Beres, B.Beck and J. Zhou. Proposal of Zygosaccharomyces parabailii sp. nov. and Zygosaccharomyces pseudobailii sp. nov., novel species closely related to Zygosaccharomyces bailii. Int. J. Syst. Evol. Microbial., May 2013; 63: 1922-1929.
17. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories. 1988.
18. Versalovic, J., G. Funke, K.C. Carroll, J.H. Jorgensen, M.L. Landry and D.W. Warnock. 2011. *Manual of Clinical Microbiology*, 10th edition. American Society for Microbiology, Washington, D.C.
19. Kurtzman, C.P., Fell, J.W., Boekhout, T. (eds.) *The Yeasts - A Taxonomic Study*, 5th edn. Amsterdam: Elsevier Science B.V., 2011.
20. Lockhart, S.R., Messer, S.A., Pfaller, M.A., Diekema, D.J. Geographic distribution and antifungal susceptibility of the newly described specied *Candida orthopsilosis* and *Candida metapsilosis* in comparison to the closely related species *Candida parapsilosis*. J Clin Microbiol. 2008; 46:2659-2664.

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Use this Instructions for Use with VITEK® 2 Product No. 21343.

INDEX OF SYMBOLS

Symbol	Meaning
	Catalog number
	In Vitro Diagnostic Medical Device
	Legal Manufacturer
	Temperature limitation
	Use by date
	Batch code
	Consult Instructions for Use
	Date of manufacture
	Contains sufficient for <n> tests
	Authorized representative in the European Community
	For US Only : Caution : US Federal Law restricts this device to sale by or on the order of a licensed practitioner

Instructions for Use provided in the kit or downloadable from www.biomerieux.com/techlib

LIMITED WARRANTY

bioMérieux warrants the performance of the product for its stated intended use provided that all procedures for usage, storage and handling, shelf life (when applicable), and precautions are strictly followed as detailed in the instructions for use (IFU).

Except as expressly set forth above, bioMérieux hereby disclaims all warranties, including any implied warranties of merchantability and fitness for a particular purpose or use, and disclaims all liability, whether direct, indirect or consequential, for any use of the reagent, software, instrument and disposables (the “System”) other than as set forth in the IFU.

WASTE DISPOSAL

All hazardous waste must be disposed of by following your local inspecting agency's guidelines.

REVISION HISTORY TABLE

Change type categories	
N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user

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Note : Minor typographical, grammar, and formatting changes are not included in the revision history.

Release Date	Part Number	Change Type	Change Summary
2019-03	043908-03	Technical change	Updated for 9.02 software release. Updated sections: • Intended Use • Precautions • Culture Requirements • Additional Information on Lab Report • Testing of QC Organisms • Performance Characteristics • Organisms Identified • References
2016-10	043908-02	Technical change	• Updated content to reflect the 8.01 Product Information Manual
2016-05	043908-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• New IFU derived from product chapter in the Product Information Manual • Updated Limited Warranty section • Updated with RX only information

BIOMÉRIEUX, the BIOMÉRIEUX logo, VITEK, API, Count-TACT, chromID, DensiCHEK and bioLiaison are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

This product may be protected by one or more patents, see: http://www.biomerieux-usa.com/patents.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

Any other name or trademark is the property of its respective owner.

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Rx only IVD

INTENDED USE

These Instructions for Use correspond to the VITEK® 2 Systems 7.01 or higher software. If you are not using VITEK® 2 Systems 7.01 or higher software, please refer to the VITEK® 2 Systems Product Information that you received with your current software version.

The VITEK® 2 Gram-Negative identification card (GN) is intended for use with VITEK® 2 Systems for the automated identification of most clinically significant fermenting and non-fermenting Gram-negative bacilli. The VITEK® 2 GN identification card is a single-use disposable. For a list of claimed species, see the Organisms Identified section.

DESCRIPTION

The GN card is based on established biochemical methods[1,2,4,8,9,10,11,12,17,18,20,21,24,25,27] and newly developed substrates measuring carbon source utilization, enzymatic activities, and resistance. There are 47 biochemical tests and one negative control well. The Decarboxylase Negative Control Well (well 52) is used as a baseline reference for the Decarboxylase test wells. Final results are available in approximately 10 hours or less.

For a list of well contents, see the GN Well Contents table.

Table 1: GN Well Contents

Well	Test	Mnemonic	Amount/Well
2	Ala-Phe-Pro-ARYLAMIDASE	APPA	0.0384 mg
3	ADONITOL	ADO	0.1875 mg
4	L-Pyrrolydonyl-ARYLAMIDASE	PyrA	0.018 mg
5	L-ARABITOL	lARL	0.3 mg
7	D-CELLOBIOSE	dCEL	0.3 mg
9	BETA-GALACTOSIDASE	BGAL	0.036 mg
10	H2S PRODUCTION	H2S	0.0024 mg
11	BETA-N-ACETYL-GLUCOSAMINIDASE	BNAG	0.0408 mg
12	Glutamyl Arylamidase pNA	AGLTp	0.0324 mg
13	D-GLUCOSE	dGLU	0.3 mg
14	GAMMA-GLUTAMYL-TRANSFERASE	GGT	0.0228 mg
15	FERMENTATION/ GLUCOSE	OFF	0.45 mg
17	BETA-GLUCOSIDASE	BGLU	0.036 mg
18	D-MALTOSE	dMAL	0.3 mg
19	D-MANNITOL	dMAN	0.1875 mg
20	D-MANNOSE	dMNE	0.3 mg
21	BETA-XYLOSIDASE	BXYL	0.0324 mg
22	BETA-Alanine arylamidase pNA	BAlap	0.0174 mg
23	L-Proline ARYLAMIDASE	ProA	0.0234 mg
26	LIPASE	LIP	0.0192 mg
27	PALATINOSE	PLE	0.3 mg
29	Tyrosine ARYLAMIDASE	TyrA	0.0276 mg
31	UREASE	URE	0.15 mg

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Note: Other well numbers between 1 and 64 not designated in this table are empty.

PRECAUTIONS

Note: For industry customers that need assistance on selecting the correct VITEK® 2 identification card, please refer to the VITEK® 2 Compact Instrument User Manual chapter, "Guidance to Select a VITEK® 2 Identification Card."

- For *In Vitro* Diagnostic Use Only.
- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For professional use only.
- Suspensions not within the appropriate zone on the VITEK® 2 DENSICHEK™ Plus or the VITEK® 2 DENSICHEK™ may compromise card performance.
- Do not use the card after the expiration date shown on the package liner.
- Store the card unopened in the package liner. Do not use the card if the protective package liner is damaged or if no desiccant is present.
- Allow the card to come to room temperature before opening the package liner.
- Do not use powdered gloves. Powder may interfere with the optics.
- Use of culture media other than the recommended types must be validated by the customer laboratory for acceptable performance.
- A Gram stain should be performed to determine an organism’s Gram reaction and morphology prior to selecting the identification card to inoculate.
- The card performs as intended only when used in conjunction with VITEK® 2 Systems, following the instructions contained in these Instructions for Use.

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- **Do not use glass test tubes**. Use clear plastic (polystyrene) tubes only. Variation exists among test tubes of standard diameter. Carefully place the tube into the cassette. If resistance is encountered, discard and try another tube that does not require pressure to insert.
- Prior to inoculation, inspect cards for tape tears or damage to the tape and discard any that are suspect. Check the saline level in the tubes after the cassette has been processed to ensure proper filling of card.
 - VITEK® 2 60 or VITEK® 2 XL: Eject improperly filled cards.
 - VITEK® 2 Compact: Do not load improperly filled cards.
- Give special consideration to specimen source and patient drug or antimicrobic regimen.
- Interpretation of test results requires the judgment and skill of a person knowledgeable in microbial identification testing. Additional testing may be required. (See the Supplemental Tests section.)
- Do not clean saline dispenser with chemical agents. The use of chemical agents may impact card performance.

Warning: All patient specimens, microbial cultures, and inoculated VITEK® 2 cards, along with associated materials, are potentially infectious and should be treated with universal precautions.[23,26] It is suggested that highly pathogenic species such as *Brucella melitensis*, *Burkholderia mallei*, *Burkholderia pseudomallei*, *Escherichia coli O157*, *Francisella tularensis*, and *Yersinia pestis* be sent to your state health laboratory or other suitable reference laboratory for confirmation.

Warning: All hazardous waste must be disposed of by following your local inspecting agency's guidelines.

STORAGE CONDITIONS

Upon receipt, store VITEK® 2 GN cards unopened in their original package liner at 2°C to 8°C.

SPECIMEN PREPARATION

For specimen preparation information, see the Culture Requirements Table.

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Table 2: Culture Requirements Table

VITEK® 2 Card	Media	Age of Culture[1]	Incubation Conditions	McFarland Standards	Dilution for AST	Age of Suspension Before Loading Instrument
GN	TSA[2,3] CBA[2,3] MAC[2,3] BCP CET CLED CHOC CHOC PVX CHBA CNT CPS ID DENA DRIG HEK SM ID TSAHB TSAB TSAL VRBG XLD	18 to 24 hours	35°C to 37°C Aerobic, non-CO_2	0.50 to 0.63 McFarland Standard	N/A[4]	≤ 30 minutes
GN and AST GN pair	CBA MAC TSAB CPS ID	18 to 24 hours	35°C to 37°C Aerobic, non-CO_2	0.50 to 0.63 McFarland Standard	145 µL in 3.0 mL saline	< 30 minutes

[1]Cultures with scant or poor growth may give unidentified or incorrect results even when the Age of Culture requirements are met.

[2]These media were used in the identification product database development and will give optimal performance.

[3]OMA Official Methods of Analysis validated medium.

[4]N/A = not applicable

Culture Requirements Table — Media Abbreviations

BCP = Bromcresol Purple Agar

CBA = Columbia Blood Agar with 5% Sheep Blood

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CLED = Cystine Lactose Electrolyte Deficient Agar

CNT = Count-TACT® (irradiated) Trypticase Soy Agar

CPS ID = chromID™ CPS (CPS ID agar)

DENA = DE Neutralizing Agar

DRIG = Drigalski Agar

HEK = Hektoen Agar

MAC = MacConkey Agar

SM ID = chromID™Salmonella (SM ID2 Agar)

TSA = Trypticase Soy Agar

TSAB = Trypticase Soy Agar with 5% Sheep Blood

TSAHB = Trypticase Soy Agar with 5% Horse Blood

TSAL = TSA with Lecithin and P80

VRBG = Violet Red Bile Glucose Agar

XLD = Xylose Lysine Desoxycholate

Materials

When used with VITEK® 2 instrumentation, the GN card is a complete system for routine identification testing of most significant fermenting and non-fermenting Gram-negative bacilli.

Required materials are:

- VITEK® 2 GN Card
- DENSICHEK™ Plus Kit or VITEK® DENSICHEK® Kit
- DENSICHEK™ Plus Standards Kit or DENSICHEK® Standards Kit
- VITEK® 2 Cassette
- Sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- 12 mm x 75 mm clear plastic (polystyrene) disposable test tubes
- Sterile sticks or swabs
- Appropriate agar medium (see Culture Requirements table).

Optional accessories:

- Adjustable volume saline dispenser
- Loops
- Pre-dispensed saline test tubes (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- Test tube caps
- Vortex

Procedure

Warning: Failure to follow instructions and recommendations provided in this section for performing laboratory tasks may cause erroneous or delayed results.

For product-specific information, see the Culture Requirements table.

Note: Prepare the inoculum from a pure culture, according to good laboratory practices. In case of mixed cultures, a re-isolation step is required. It is recommended that a purity check plate be done to ensure that a pure culture was used for testing.

1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.
2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).

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3. Use a sterile stick or swab to transfer a sufficient number of morphologically similar colonies to the saline tube prepared in step 2. Prepare a homogenous organism suspension with a density equivalent to a McFarland No. 0.50 to 0.63 using a calibrated VITEK® 2 DENSICHEK™ Plus or VITEK® 2 DENSICHEK™.
 Note: Age of suspension must not exceed 30 minutes before inoculating card.
4. Place the suspension tube and GN card in the cassette.
5. Refer to the appropriate Instrument User Manual for instructions on data entry and how to load the cassette into the instrument.
6. Follow your local inspecting agency's guidelines for disposal of hazardous waste.

RESULTS

Identification Analytical Techniques

VITEK® 2 Systems identify an organism by using a methodology based on the characteristics of the data and knowledge about the organism and reactions being analyzed. Sufficient data have been collected from known strains to estimate the typical reactions of the claimed species to a set of discriminating biochemicals. If a unique identification pattern is not recognized, a list of possible organisms is given, or the strain is determined to be outside the scope of the database.

The printed lab report contains suggestions for any supplemental tests necessary to complete the identification. If the tests are not sufficient to complete the identification, then standard microbiology references and literature should be consulted.

Certain species may belong to a slashline (mixed) taxa identification. This occurs when the biopattern is the same for the taxa listed. Supplemental tests may be used to separate slashline taxa. The species in the GN Slashline Taxa table belong to the GN slashline taxa.

Table 3: GN Slashline Taxa

Slashline Name	Species Belonging to the Slashline
For 7.01 or Higher Software Users	
Acinetobacter baumannii complex	*Acinetobacter baumannii* *Acinetobacter calcoaceticus* *Acinetobacter pittii* (*Acinetobacter* genomospecies 3) *Acinetobacter nosocomialis* (*Acinetobacter* genomospecies TU13)
Brevundimonas diminuta/vesicularis	*Brevundimonas diminuta* *Brevundimonas vesicularis*
Burkholderia cepacia group	*Burkholderia cepacia* *Burkholderia multivorans* *Burkholderia stabilis* *Burkholderia vietnamiensis*
Enterobacter cloacae complex	*Enterobacter cloacae* ssp. *cloacae* *Enterobacter hormaechei* *Enterobacter kobei* *Enterobacter ludwigii* *Enterobacter cloacae* ssp. *dissolvens*
Moraxella group	*Moraxella lacunata* *Moraxella nonliquefaciens* *Moraxella osloensis*
Neisseria animaloris/zoodegmatis	*Neisseria animaloris* *Neisseria zoodegmatis*

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Table 4: Identification Card Qualifying Messages

ID Message Confidence Level	Choices	% Probability	Comments
Excellent	1	96 to 99	N/A
Very Good	1	93 to 95	N/A
Good	1	89 to 92	N/A
Acceptable	1	85 to 88	N/A
Low Discrimination	2 to 3	Sum of choices = 100; after resolution to one choice, percent probability reflects the number associated with selected choice.	Two to three taxa exhibit same biopattern. Separate by supplemental testing. Must resolve to mate with susceptibility card.
Inconclusive or Unidentified Organism	>3 or 0	N/A	Either > 3 taxa exhibit same biopattern or Very atypical biopattern. Does not correspond to any taxon in the database. Check Gram stain and purity.

PERCENT PROBABILITY

As part of the identification process, the software compares the test set of reactions to the expected set of reactions of each organism, or organism group, that can be identified by the product. A quantitative value, the percent probability, is calculated and relates to how well the observed reactions compare to the typical reactions of each organism. A perfect match between the test reaction pattern and the unique reaction pattern of a single organism, or organism group, would provide a percent probability of 99. When a perfect match is not obtained, it is still possible for the reaction pattern to be sufficiently close to that of an expected reaction pattern such that a clear decision can be provided about the organism identification. The range of percent probabilities in the one-choice case is 85 to 99. Values closer to 99 indicate a closer match to the typical pattern for the given organism.

When the reaction pattern is not sufficient to discriminate between two to three organisms, the percent probabilities reflect this ambiguity. The reported probability values indicate, relatively, the order in which the reaction pattern best corresponds to the listed possibilities. The order does not, however, suggest that the pattern match to one of the possible identifications is clearly superior to another. The probability characteristic of an overall sum of 100 is retained through the calculation process. After resolution to one choice, the probability characteristic of the single choice is retained.

ADDITIONAL INFORMATION ON LAB REPORT

Supplemental test — External (offline) test that allows the user to resolve a slashline or Low Discrimination identification. Numbers in parentheses indicate percent positive reaction for the species/test listed.

Contraindicating test — Test result that is unusual for a reported taxon.

Table 5: Notes Associated with Certain Taxa

Taxa	Note
For 7.01 or Higher Software Users	

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Taxa	Note
Brucella melitensis	**Important!** Presumptive Identification Highly pathogenic organism. The following are included in an identification of Brucella melitensis: *Brucella melitensis abortus* *Brucella melitensis canis* *Brucella melitensis melitensis* *Brucella melitensis neotamae* *Brucella melitensis ovis* *Brucella melitensis suis*
Burkholderia mallei	**Important!** Presumptive Identification Highly pathogenic organism.
Burkholderia pseudomallei	Highly pathogenic organism. Isolates of *Burkholderia thailandensis* are biochemically similar to *Burkholderia pseudomallei.* Since the possibility of *Burkholderia thailandensis* exists, the user should send the isolate to their state lab or other suitable reference laboratory for confirmation.
Escherichia coli O157	Confirm by serological tests. Highly pathogenic organism.
Francisella tularensis	Confirm by serological tests. Highly pathogenic organism.
Salmonella enterica ssp. *arizonae* *Salmonella enterica* ssp. *diarizonae* *Salmonella* group *Salmonella* ser. Gallinarum *Salmonella* ser. Paratyphi A *Salmonella* ser. Typhi	Confirm by serological tests.
Shigella group *Shigella sonnei*	Confirm by serological tests.
Vibrio cholerae	Critical pathogen. The species identified may have significance to patient or sample outcome and can be stopped for review.
Yersinia pestis	**Important!** Presumptive Identification Highly pathogenic organism.
For 9.02 Software Users	
Ochrobactrum anthropi	Possibility of *Brucella* spp.

Notes Associated with an Improperly Filled Card or with a Negative Profile (Biopattern)

- For the case where the time between two readings is greater than 40 minutes: “CARD ERROR — Missing data.”
- For the case where there is a negative profile: “Organism with low reactivity biopattern — please check viability.”
- When a biopattern is calculated for an unknown organism that is completely negative or consists of both negative tests and tests that fall within the uncertainty zone, the identification call will be “Non or low reactive biopattern.”

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The following species could potentially trigger this note if a test was atypical or fell within the uncertainty zone:

- *Acinetobacter haemolyticus*
- *Acinetobacter lwoffii*
- *Actinobacillus ureae*
- *Aeromonas salmonicida*
- *Brucella melitensis*
- *Francisella tularensis*
- *Methylobacterium* spp.
- *Moraxella lacunata*
- *Moraxella nonliquefaciens*
- *Moraxella osloensis*
- *Pasteurella multocida*
- *Pseudomonas alcaligenes*
- *Pseudomonas fluorescens*
- *Pseudomonas stutzeri*

QUALITY CONTROL

Quality control organisms and their expected results are listed in the VITEK® 2 GN Quality Control Tables. Process these according to the procedure for test isolates outlined in this document.

Certification Statement

This is to certify that bioMérieux complies with ISO 13485 and FDA Quality System Regulation (QSR) requirements for design, development, and manufacture of microbial identification systems.

Frequency of Testing

Currently, it is recommended that you use your most stringent inspecting agency's guidelines for frequency of identification product testing.

Common practice is to perform QC upon receipt of shipment of the test kits. Reactions must follow Instructions for Use results.

If the results do not meet the criteria, subculture for purity and repeat the test. If discrepant results are repeated, perform an alternate identification method and contact bioMérieux.

Testing and Storage of QC Organisms

1. Rehydrate the organism according to the manufacturer's instructions.
2. Use Trypticase Soy agar with 5% sheep blood agar (TSAB). Incubate aerobically at 35°C to 37°C for approximately 18 to 24 hours.
3. Check for purity. Perform second subculture for testing.
4. Use Trypticase Soy agar with 5% sheep blood agar (TSAB). Incubate aerobically at 35°C to 37°C for approximately 18 to 24 hours.

Short-Term Storage Conditions

1. Streak to a TSAB plate or slant.
2. Incubate for 24 hours at 35°C to 37°C.
3. Refrigerate at 2°C to 8°C for up to two weeks.
4. Subculture once as described above and use for QC.

Long-Term Storage Conditions

1. Make a heavy suspension in Tryptic Soy Broth (TSB) with 15% glycerol.
2. Freeze at -70°C.
3. Subculture to TSAB twice before running QC.

Note: Avoid repeated thawing and refreezing by either freezing in single-use aliquots or removing a small portion of frozen organism preparation with a sterile applicator stick.

STREAMLINED QUALITY CONTROL

Note: Industrial Use Only laboratories should perform quality control following the Streamlined Quality Control section. No additional testing is required for these users.

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As there are no substrates that are consistently sensitive to degradation during shipping conditions, streamlined quality control may be conducted by testing two strains: one that is mostly positive and the other which is mostly negative for reactions on GN. See the GN Quality Control tables.

COMPREHENSIVE QUALITY CONTROL

Customers who do not qualify for streamlined quality control testing are required to perform comprehensive quality control testing, which entails demonstration of a positive and negative reaction for each substrate of an identification product.[6]

In order to qualify initially for streamlined quality control testing, the CLSI® M50-A standard requires that the user perform and document either of the following:[5]

- Verification testing to show that performance is equivalent to the manufacturer's claims.
- Comprehensive quality control testing of at least three lots over at least three different seasons.

Refer to the complete CLSI® M50-A standard for information regarding continued qualification and further details of requirements and responsibilities for both the user and the manufacturer related to streamlined quality control testing.

GN Quality Control Tables:

***Enterobacter hormaechei* ATCC® 700323™** (for streamlined or comprehensive quality control)

***Stenotrophomonas maltophilia* ATCC® 17666™** (for streamlined or comprehensive quality control)

***Acinetobacter baumannii* ATCC® BAA-747™** (for comprehensive quality control)

***Elizabethkingia meningoseptica* ATCC® 13253™** (for comprehensive quality control)

***Klebsiella oxytoca* ATCC® 700324™** (for comprehensive quality control)

***Ochrobactrum anthropi* ATCC® BAA-749™** (for comprehensive quality control)

***Proteus vulgaris* ATCC® 6380™** (for comprehensive quality control)

***Pseudomonas aeruginosa* ATCC® 9721™** (for comprehensive quality control)

***Pseudomonas aeruginosa* ATCC® BAA-1744™** (for comprehensive quality control)

Note: *Pseudomonas aeruginosa* ATCC® BAA-1744™ may contain two morphologically distinct colony types; however, either will provide proper expected reactions when tested for quality control.

For 7.01 Software Users

***Shigella sonnei* ATCC® 25931™** (for comprehensive quality control)

For 8.01 or Higher Software Users

***Escherichia coli* ATCC® 25922™** (for comprehensive quality control)

The GN card typically identifies the quality control organisms as one-choice or within a low discrimination or slashline identification. However, strains are chosen for reaction performance over identification performance. Therefore, an unidentified or misidentified result may occur when all expected quality control reactions are correct.

Table 6: QC Organism: *Enterobacter hormaechei* ATCC® 700323™ (for streamlined or comprehensive quality control)

APPA	-	AGLTp	-	BXYL	+	SAC	+	SUCT	v	CMT	-
ADO	+	dGLU	+	BAlap	-	dTAG	-	NAGA	+	BGUR	v
PyrA	-	GGT	+	ProA	v	dTRE	+	AGAL	+	O129R	+
IARL	-	OFF	+	LIP	v	CIT	+	PHOS	v	GGAA	-
dCEL	+	BGLU	-	PLE	+	MNT	+	GlyA	v	IMLTa	-
BGAL	+	dMAL	+	TyrA	v	5KG	-	ODC	+	ELLM	-
H2S	-	dMAN	+	URE	-	ILATk	v	LDC	-	ILATa	-
BNAG	+	dMNE	+	dSOR	+	AGLU	-	IHISa	-		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 7: QC Organism: *Stenotrophomonas maltophilia* ATCC® 17666™ (for streamlined or comprehensive quality control)

APPA	+	AGLTp	-	BXYL	-	SAC	-	SUCT	v	CMT	-

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+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 8: QC Organism: *Acinetobacter baumannii* ATCC® BAA-747™ (for comprehensive quality control)

APPA	v	AGLTp	v	BXYL	v	SAC	v	SUCT	+	CMT	v
ADO	v	dGLU	v	BAlap	v	dTAG	v	NAGA	v	BGUR	v
PyrA	v	GGT	v	ProA	v	dTRE	v	AGAL	v	O129R	v
lARL	v	OFF	v	LIP	v	CIT	+	PHOS	–	GGAA	v
dCEL	v	BGLU	v	PLE	v	MNT	+	GlyA	v	IMLTa	v
BGAL	v	dMAL	v	TyrA	+	5KG	v	ODC	v	ELLM	v
H2S	v	dMAN	v	URE	v	ILATk	+	LDC	v	ILATa	v
BNAG	v	dMNE	v	dSOR	v	AGLU	v	lHISa	+		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 9: QC Organism: *Elizabethkingia meningoseptica* ATCC® 13253™ (for comprehensive quality control)

APPA	+	AGLTp	+	BXYL	v	SAC	v	SUCT	–	CMT	v
ADO	v	dGLU	–	BAlap	v	dTAG	v	NAGA	+	BGUR	v
PyrA	+	GGT	v	ProA	v	dTRE	v	AGAL	v	O129R	v
lARL	v	OFF	–	LIP	v	CIT	v	PHOS	v	GGAA	+
dCEL	v	BGLU	v	PLE	v	MNT	v	GlyA	+	IMLTa	v
BGAL	v	dMAL	v	TyrA	v	5KG	v	ODC	v	ELLM	v
H2S	v	dMAN	v	URE	v	ILATk	–	LDC	v	ILATa	v
BNAG	+	dMNE	v	dSOR	v	AGLU	+	lHISa	v		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

Table 10: QC Organism: *Klebsiella oxytoca* ATCC® 700324™ (for comprehensive quality control)

APPA	–	AGLTp	v	BXYL	v	SAC	v	SUCT	v	CMT	v
ADO	+	dGLU	+	BAlap	v	dTAG	+	NAGA	v	BGUR	–
PyrA	v	GGT	–	ProA	–	dTRE	+	AGAL	+	O129R	v
lARL	+	OFF	+	LIP	–	CIT	v	PHOS	v	GGAA	–
dCEL	+	BGLU	+	PLE	+	MNT	v	GlyA	–	IMLTa	v
BGAL	+	dMAL	v	TyrA	v[2]	5KG	v[1]	ODC	–	ELLM	v
H2S	v	dMAN	+	URE	+	ILATk	v	LDC	+	ILATa	v
BNAG	–	dMNE	+	dSOR	v	AGLU	–	lHISa	v		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

[1]Reaction is mostly positive although occasional negative reaction may occur.

[2]Reaction is mostly negative although occasional positive reaction may occur.

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[1]Reaction is mostly positive although occasional negative reaction may occur.

Note: Culture may contain two morphologically distinct colony types; however, either will provide proper expected reactions when tested for quality control.

For 7.01 Software Users

Table 15: QC Organism: *Shigella sonnei* ATCC® 25931™ (for comprehensive quality control)

APPA	v	AGLTp	v	BXYL	–	SAC	–	SUCT	v	CMT	+
ADO	v	dGLU	v	BAlap	v	dTAG	v	NAGA	–	BGUR	+
PyrA	v	GGT	–	ProA	v	dTRE	v	AGAL	v	O129R	v
lARL	v	OFF	v	LIP	v	CIT	–	PHOS	v	GGAA	v
dCEL	v	BGLU	–	PLE	–	MNT	–	GlyA	v	IMLTa	v
BGAL	v	dMAL	+	TyrA	+	5KG	v	ODC	+	ELLM	v
H2S	v	dMAN	v	URE	v	lLATk	v	LDC	v	lLATa	v
BNAG	–	dMNE	v	dSOR	v	AGLU	v	lHISa	v		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

For 8.01 or Higher Software Users

Table 16: QC Organism: *Escherichia coli* ATCC® 25922™ (for comprehensive quality control)

APPA	v	AGLTp	v	BXYL	–	SAC	–	SUCT	v	CMT	+
ADO	v	dGLU	v	BAlap	v	dTAG	v	NAGA	–	BGUR	+
PyrA	v	GGT	–	ProA	v	dTRE	v	AGAL	v	O129R	v
lARL	v	OFF	v	LIP	v	CIT	–	PHOS	v	GGAA	v
dCEL	v	BGLU	–	PLE	–	MNT	–	GlyA	v	IMLTa	v
BGAL	v	dMAL	+	TyrA	+	5KG	v	ODC	+	ELLM	v
H2S	v	dMAN	v	URE	v	lLATk	v	LDC	v	lLATa	v
BNAG	–	dMNE	v	dSOR	v	AGLU	v	lHISa	v		

+ = 95% to 100% positive; v = 6% to 94% positive; – = 0% to 5% positive

LIMITATIONS

The VITEK® 2 GN card cannot be used with direct clinical samples or other sources containing mixed flora. Any change or modification in the procedure may affect the results.

Newly described or rare species may not be included in the GN database. Selected species will be added as strains become available.

Warning: Testing of unclaimed species may result in an unidentified result or a misidentification.

PERFORMANCE CHARACTERISTICS

For 7.01 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 GN identification card was evaluated using 562 clinical and stock isolates of both commonly and rarely observed species of gram-negative bacilli, including 153 non-fermentative strains. The reference identification was determined with API® 20 E and API® 20 NE identification kits. Overall, the VITEK® 2 GN correctly identified 96.2% of the isolates, including 6.8% low discrimination with the correct species listed. Misidentifications occurred at 3.4% and no identifications occurred at 0.4%.

For 8.01 and 9.01 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 GN identification card was evaluated using 562 clinical and stock isolates of both commonly and rarely observed species of gram-negative bacilli, including 153 non-fermentative strains. The reference identification was determined with API® 20 E and API® 20 NE identification kits. Overall, the VITEK® 2 GN correctly identified 95.4% of the isolates, including 6.6% low discrimination with the correct species listed. Misidentifications occurred at 4.1% and no identifications occurred at 0.5%.

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For 9.02 Software Users

In a multi-site clinical study*, the performance of the VITEK® 2 GN identification card was evaluated using 562 clinical and stock isolates of both commonly and rarely observed species of gram-negative bacilli, including 153 non-fermentative strains. The reference identification was determined with API® 20 E and API® 20 NE identification kits. Overall, the VITEK® 2 GN correctly identified 95.2% of the isolates, including 6.4% low discrimination with the correct species listed. Misidentifications occurred at 4.3% and no identifications occurred at 0.5%.

*Data on file at bioMérieux, Inc.

ORGANISMS IDENTIFIED

Claims are for all software users unless otherwise stated.

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Additional Claims and Taxonomy Changes For 8.01 or Higher Software Users

- *Hafnia paralvei*
- *Lelliottia amnigena* 1* (formerly known as *Enterobacter amnigenus* 1)
- *Lelliottia amnigena* 2* (formerly known as *Enterobacter amnigenus* 2)
- *Pandoraea* spp.
- *Pluralibacter gergoviae** (formerly known as *Enterobacter gergoviae*)
- *Ralstonia insidiosa*
- *Tatumella ptyseos*

Additional Claims For 9.02 Software Users

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Additional Claims For 8.01 or Higher Software Users

- *Pandoraea* species
- *Ralstonia insidiosa*

Additional Claims and Taxonomy Changes For 9.02 Software Users

- *Aeromonas hydrophila/Aeromonas punctata* (formerly known as *Aeromonas caviae*)
- *Bergeyella zoohelcum*

Highly Pathogenic Organisms

- *Brucella melitensis**
- *Burkholderia mallei**
- *Burkholderia pseudomallei**
- *Escherichia coli* O157*
- *Francisella tularensis**
- *Yersinia pestis**

* OMA Official Methods of Analysis validated claim.

+ Species within this group or complex that are OMA Official Methods of Analysis validated claims are *Burkholderia cepacia*, *Cronobacter sakazakii*, and *Enterobacter cloacae*.

SUPPLEMENTAL TESTS

Table 17: GN Supplemental Tests

Abbreviation	Test Name	Description	Comments	Reference
For 7.01 or Higher Software Users				
41C	GROWTH AT 41°C	Ability of certain species to grow at 41°C.	N/A	18, 20
42C	GROWTH AT 42°C	Ability of certain species to grow at 42°C.	N/A	20, 22
44C	GROWTH AT 44°C	Ability of certain species to grow at 44°C.	N/A	21

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Abbreviation	Test Name	Description	Comments	Reference
ADONITOL dCELLOB dMALTOSE dMANNITOL dMELIBIOSE dSORBITOL dTREHALOSE dTURANOSE DUL INOSITOL LACTOSE IRHAMNOSE SACCHAROSE SALICIN	ADONITOL acidification D-CELLOBIOSE acidification D-MALTOSE acidification D-MANNITOL acidification D-MELIBIOSE acidification SORBITOL acidification D-TREHALOSE acidification TURANOSE acidification DULCITOL acidification INOSITOL acidification LACTOSE acidification L-RHAMNOSE acidification SACCHAROSE/SUCROSE acidification SALICIN acidification	Acidification of carbon source observed with pH indicator (e.g., phenol red, bromcresol purple, etc.).	Some tests also appear on the GN card but are recommended as supplemental tests since results of conventional macromethods may differ from rapid commercial micromethods.	2, 8, 10, 12, 13, 14, 16, 17, 19, 21, 22, 27, 28
Arg.hydr.	ARGININE dihydrolase	Hydrolysis of arginine releases an amine resulting in alkalinization of the medium observed with a pH indicator (e.g., red color formation in the presence of phenol red).	N/A	7, 10, 12, 17, 18, 19, 20, 22, 25, 27
B-HEM	BETA HEMOLYSIS	Certain species possess hemolysins that give a transparent zone around colonies on blood-based agar.	N/A	3, 9, 20, 27
DNAse	DNAse test	Ability of certain species to produce DNAse resulting in the degradation of DNA.	N/A	17, 20, 27
ESCULIN	ESCULIN hydrolysis	Hydrolysis of esculin forms esculetin that produces a black pigment in the presence of iron salts.	N/A	12, 17, 19, 20, 27
GELATIN	GELATIN hydrolysis	Mediated by a gelatinase enzyme, a positive reaction is observed by liquefaction of the gelatin substrate.	N/A	3, 9, 18, 19, 20, 22, 24
dGLUf	Glucose fermentation	Fermentation of glucose observed with pH indicators (e.g. phenol red, bromcresol purple, etc.).	Some tests also appear on the GN card but are recommended as supplemental tests since results of conventional macromethods may differ from rapid commercial micromethods.	29
IND	INDOLE	Ability of certain species to split indole from tryptophan detected by a colored product revealed with a specific reagent (e.g., Kovacs, Ehrlich's, DMAC reagents, etc.).	N/A	10, 12, 16, 17, 19, 20, 27
JordanTART	Jordan_Tartrate	Fermentation of tartrate results in acidification of the medium observed with a pH indicator (e.g., yellow color formation in the presence of phenol red).	N/A	19
Lysine dec.	Lysine decarboxylase	Hydrolysis of lysine releases an amine resulting in alkalinization of the medium observed with a pH indicator (e.g., purple color formation in the presence of bromcresol purple).	Some tests also appear on the GN card but are recommended as supplemental tests since results of conventional macromethods may differ from rapid commercial micromethods.	21, 22

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Abbreviation	Test Name	Description	Comments	Reference
MNTka	MALONATE alkalinization	Utilization of malonate as sole carbon source.	N/A	15, 16, 30
MOB	MOTILITY	Test for motility using hanging drop procedure or wet mount.	Bacterial motility can be observed by placing a drop of bacterial suspension on a slide and viewing it under a microscope.	4, 12, 17, 19, 20, 25, 27, 28, 30
NAT	SODIUM-ACETATE alkalinization	Ability of certain species to utilize acetate as a sole source of carbon.	N/A	29
NO2 NO3 NO3→N2	NITRITE REDUCTION NITRATE REDUCTION NITROGEN PRODUCTION FROM NO3	Test for the ability to reduce nitrite to nitrogen gas (NO2) nitrate to nitrite and/ or nitrogen gas from nitrate (NO3→N2)	N/A	10, 20, 22, 29, 30
NaCl 0% NaCl 6%	GROWTH IN 0% NaCl GROWTH IN 6% NaCl	Ability of certain species to grow in the presence or absence of 6.0% NaCl.	N/A	7, 8, 20, 21, 22
O/129 R	O/129 RESISTANCE	Ability of certain species to grow in the presence of the vibriostatic compound O/129.	Some tests also appear on the GN card but are recommended as supplemental tests since results of conventional macromethods may differ from rapid commercial micromethods.	8, 11
ONPG	BETA_ GALACTOSIDASE	Presence of beta-galactosidase cleaves o-nitrophenol-beta-D-galactopyranoside to produce a yellow colored product.	N/A	8, 12, 17, 19, 20
Ornith.dec	Ornithine decarboxylase	Hydrolysis of ornithine releases an amine resulting in alkalinization of the medium observed with a pH indicator (e.g., purple color formation in the presence of bromcresol purple).	Some tests also appear on the GN card but are recommended as supplemental tests since results of conventional macromethods may differ from rapid commercial micromethods.	8, 10, 17, 19, 20, 27
OX	OXIDASE	Detection of the presence of cytochrome C.	Characteristic useful in identifying many species of non-fermenters. All members of *Enterobacteriaceae* are oxidase negative.	10, 12, 17, 18, 19, 20, 21, 22, 25, 27, 28
PURPLE	PURPLE PIGMENT	Ability of certain species to produce purple colonies on non-differential media.	Characteristic of *Chromobacterium violaceum*.	19, 20
PYOCYANIN	PYOCYANIN pigment	Ability of species to produce blue pigment (pyocyanin) or fluorescent pigment (pyoverdin).	Presence of both pyocyanin and pyoverdin is characteristic of *Pseudomonas aeruginosa* producing greenish fluorescent colonies.	1, 20
PYOVERDIN	PYOVERDIN pigment			
RM	Methyl Red	Test for acid production, requiring positive organisms to produce acid from glucose.	N/A	21
UREASE	Urease	Hydrolysis of urea releases ammonia resulting in alkalinization of the medium observed with a pH indicator (e.g., red color formation in the presence of phenol red).	N/A	10, 12, 17, 19, 20, 25, 27
VP	VOGES PROSKAUER	Ability of some species to produce acetoin from glucose fermentation.	N/A	12, 17, 19, 20, 25, 30

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Abbreviation	Test Name	Description	Comments	Reference
YELLOW	YELLOW PIGMENT	Ability of certain species to produce yellow pigmented colonies on nondifferential media.	N/A	12, 17, 19, 20, 29
For 7.01 Software Users Only				
dFRUCTOSEa dGLUCOSEa dMANNITOLa dMELa lSORBOSEa	D-FRUCTOSE assimilation D-GLUCOSE assimilation D-MANNITOL assimilation D-MELIBIOSE assimilation L-SORBOSE assimilation	Capability of organisms to grow using a specific sole carbon source.	N/A	2, 4, 17, 18
dMLZ	MELEZITOSE acidification	Acidification of carbon source observed with pH indicator (e.g., phenol red, bromcresol purple, etc.).	Some tests also appear on the GN card but are recommended as supplemental tests since results of conventional macromethods may differ from rapid commercial micromethods.	8, 10, 12, 13, 14, 16, 17, 19, 21, 22, 27
For 8.01 or Higher Software Users				
dGLUCOSE dMELEZIT. dXYLOSE lSORBOSE	D-GLUCOSE acidification MELEZITOSE acidification D-XYLOSE acidification L-SORBOSE acidification	Acidification of carbon source observed with pH indicator (e.g., phenol red, bromcresol purple, etc.).	Some tests also appear on the GN card but are recommended as supplemental tests since results of conventional macromethods may differ from rapid commercial micromethods.	2, 8, 10, 12, 13, 14, 16, 17, 19, 21, 22, 27, 28
COL R	COLISTIN RESISTANCE	Ability of certain species to grow in the presence of the colistin.	N/A	28

REFERENCES

1. American Society for Microbiology. 98th General Meeting Workshop Program. Practical Approach to the Identification of the Medically Important Glucose Non-Fermenting Gram-Negative Bacilli. American Society for Microbiology, Washington, D.C. 1998.
2. Brenner DJ, Grimont PAD, Steigerwalt AG, Fanning GR, Ageron E, Riddle CF. Classification of Citrobacteria by DNA Hybridization: Designation of *Citrobacter farmeri* sp.nov., *Citrobacter youngae* sp.nov., *Citrobacter braakii* sp.nov., *Citrobacter werkmanii* sp.nov., *Citrobacter sedlakii* sp.nov., and Three Unnamed Citrobacter Genomospecies. Int. J. Syst. Bacteriol. 1993;43:645-658.
3. Brenner DJ, Krieg NR, Staley JT, Garrity GM, editors. Bergey's Manual of Systematic Bacteriology, 2nd Edition. Springer, New York, NY. 2005
4. Chang YH, Han J, Chun J, Lee KC, Rhee MS, Kim YB, Bae KS. *Comamonas koreensis* sp.nov., a non-motile species from wetland in Woopo, Korea. Int. J. Syst. Evol. Microbiol. 2002;52:377-381.
5. Clinical and Laboratory Standards Institute, M50-A, Quality Control for Commercial Microbial Identification Systems; Approved Guideline, Vol. 28 No. 23.
6. Clinical Laboratory Improvement Amendments of 1988. 42 U.S.C. 263a. PL 100-578.1988.
7. Coenye, T., Falsen, E., Hoste, B., Ohlen, M., Goris, J., Govan, J.R.W., Gillis, M. and Vandamme, P. Description of *Pandoraea* gen. nov. with *Pandoraea apista* sp. nov., *Pandoraea pulmonicola* sp. nov., *Pandoraea pnomenusa* sp. nov., *Pandoraea sputorum* sp. nov. and *Pandoraea norimbergensis* comb. nov. Int. J . Syst. Evol. Microbial. 2000; 50:887-889.
8. Coenye T, Mahenthiralingam E, Henry D, Lipuma JJ, Laevens S, Gillis M, Speert DP, Vandamme P. Burkholderia ambifaria sp nov., a novel member of the *Burkholderia cepacia* complex including biocontrol and cystic fibrosis-related isolates. Int. J. Syst. Evol. Microbiol. 2001; 51:1481-1490.
9. Coenye T, Vandamme P, Gowan JRW, Lipuma JJ. Taxonomy and Identification of the *Burkholderia cepacia* Complex. J. Clin. Microbiol. 2001;39:3427-3436.
10. De Baere T, Steyaert, Wauters G, De Vos P, Goris J, Coenye T, Suyama T, Verschraegen G, Vaneechoutte M. Classification of *Ralstonia pickettii* biovar 3/ 'thomasii' strains (Pickett 1994) and of new isolates related to nosocomial recurrent meningitis as Ralstonia mannitolytica sp.nov. Int. J. Syst. Evol. Microbiol. 2001;51:547-558.
11. Freney J, Renaud F, Hansen W, Bollet C. *Précis de bactériologie clinique*. ESKA, Paris, France. 2000.

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23. National Committee for Clinical Laboratory Standards, M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue* — Approved Guideline, 1997.

24. Richard C, Kiredjian M. *Laboratory methods for the Identification of the Medically Important Glucose Nonfermenting Gram-Negative Bacilli*. Institut Pasteur, Paris, France. 1992.

25. Smith S.K., Sutton D.C., Fuerst J.A., Reichelt J.L.. Evaluation of the Genus *Listonella* and the reassignment of *Listonella damsela* (Love et al.) MacDonell and Colwell to the Genus *Photobacterium* as *Photobacterium damsela* comb. nov. with an Emended Description. Int. J. Syst. Bacteriol. 1991;41:529-534.

26. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, *Biosafety in Microbiological and Biomedical Laboratories*, 1988.

27. Vandamme P, Goris J, Coenye T, Hoste B, Janssens D, Kersters K, DeVos P, Falsen E. Assignment of Centers for Disease Control group Ivc-2 to the genus *Ralstonia* as *Ralstonia paucula* sp.nov. Int. J. Syst. Bacteriol. 1999;49:663-669.

28. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry and D.W. Warnock. 2011. *Manual of Clinical Microbiology*, 10th edition. American Society for Microbiology, Washington, D.C.

29. Weyant R.S., Moss C.W., Weaver R.E., Hollis D.G., Jordan J.G., Cook E.C., and Daneshvar M.I. *Identification of Unusual Pathogenic Gram-Negative Aerobic and Facultatively Anaerobic Bacteria*. 2nd Edition. Williams & Wilkins. Baltimore, Maryland. 1996.

30. J.H. Jorgensen, M.A. Phaller, K.C. Carroll, G. Funke, M.L. Landry, S.S. Richter, and D.W. Warnock. 2015. *Manual of Clinical Microbiology*, 11th edition. American Society for Microbiology, Washington, D.C.

Use this Instructions for Use with VITEK® 2 Product No. 21341.

INDEX OF SYMBOLS

Symbol	Meaning
REF	Catalog number
IVD	In Vitro Diagnostic Medical Device

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Symbol	Meaning
	Legal Manufacturer
	Temperature limitation
	Use by date
	Batch code
	Consult Instructions for Use
	Date of manufacture
	Contains sufficient for <n> tests
	Authorized representative in the European Community
	For US Only : Caution : US Federal Law restricts this device to sale by or on the order of a licensed practitioner

Instructions for Use provided in the kit or downloadable from www.biomerieux.com/techlib

LIMITED WARRANTY

bioMérieux warrants the performance of the product for its stated intended use provided that all procedures for usage, storage and handling, shelf life (when applicable), and precautions are strictly followed as detailed in the instructions for use (IFU).

Except as expressly set forth above, bioMérieux hereby disclaims all warranties, including any implied warranties of merchantability and fitness for a particular purpose or use, and disclaims all liability, whether direct, indirect or consequential, for any use of the reagent, software, instrument and disposables (the “System”) other than as set forth in the IFU.

WASTE DISPOSAL

All hazardous waste must be disposed of by following your local inspecting agency's guidelines.

REVISION HISTORY TABLE

Change type categories	
N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user
Note :	Minor typographical, grammar, and formatting changes are not included in the revision history.

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Release Date	Part Number	Change Type	Change Summary
2019-03	044066-03	Technical change	Updated for 9.02 software release. Updated sections: • Intended Use • Precautions • Culture Requirements • Additional Information on Lab Report • Testing of QC Organisms • Performance Characteristics • Organisms Identified • References
2016-10	044066-02	Technical change	• Updated content to reflect the 8.01 Product Information Manual
		Correction	• Performance Characteristics
2016-05	044066-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• New IFU derived from product chapter in the Product Information Manual • Updated Limited Warranty section • Updated with RX only information

BIOMÉRIEUX, the BIOMÉRIEUX logo, VITEK, API, Count-TACT, chromID, DensiCHEK and bioLiaison are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

This product may be protected by one or more patents, see: http://www.biomerieux-usa.com/patents.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

Any other name or trademark is the property of its respective owner.

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Rx only IVD

INTENDED USE

The VITEK® 2 Gram-negative Susceptibility Card is intended for use with the VITEK® 2 Systems in clinical laboratories as an *in vitro* test to determine the susceptibility of clinically significant aerobic Gram-negative bacilli to antimicrobial agents when used as instructed.

SUMMARY AND EXPLANATION

Susceptibility testing is indicated for any organism that contributes to an infectious process warranting antimicrobial chemotherapy. Susceptibility tests are most often indicated when the causative organism is thought to belong to a species capable of exhibiting resistance to commonly used agents. Isolated colonies of each type of organism that may play a pathogenic role are selected from an agar plate and tested for susceptibility. These tests are then examined and the Minimum Inhibitory Concentration (MIC) is determined. The MIC obtained using a dilution test may tell the physician the concentration of an antimicrobial agent needed at the site of infection to inhibit the infecting organism.

MICs have traditionally been determined using antimicrobial concentrations derived from serial twofold dilutions.[2] The MIC is then determined from the lowest concentration that exhibits inhibition of growth. An interpretive criterion (Susceptible, Intermediate, or Resistant) can then be assigned to MIC results to aid in the direction of therapy.

For some antimicrobials (e.g., high-level gentamicin, high-level streptomycin) a qualitative result is generated.

The standard and reference procedures are based on susceptibility tests requiring 16 to 24 hours of incubation for bacteria. Various manufacturers have now developed automated procedures designed to generate results more rapidly by using shortened incubation times. Laboratories worldwide use either variations of the standard reference procedure or a commercially available product to determine the MICs of infectious organisms.

Extended-Spectrum Beta Lactamases (ESBLs)

ESBLs are enzymes that arise by mutations in genes for common plasmid-mediated beta-lactamases. Strains of *Klebsiella* spp. and *E. coli* that produce ESBLs may be clinically resistant to therapy with penicillins, cephalosporins, or aztreonam despite apparent *in vitro* susceptibility to some of these agents. Some of these strains will show MICs above the normal susceptible population but below the standard breakpoints for certain extended-spectrum cephalosporins or aztreonam. Before reporting a strain as an ESBL producer, a confirmatory method may be performed. The VITEK® 2 ESBL Test is a confirmatory test for those ESBLs inhibited by clavulanic acid, and it utilizes cefepime, cefotaxime, and ceftazidime, with and without clavulanic acid, to determine a positive or negative result.

STORAGE CONDITIONS

Upon receipt, store VITEK® 2 AST cards unopened in their original package liner at 2°C to 8°C.

PRINCIPLE OF THE TEST

The AST card for VITEK® 2 Systems is an automated test methodology based on the MIC technique reported by MacLowry and Marsh and Gerlach.[15,16] The AST card is essentially a miniaturized and abbreviated version of the doubling dilution technique for MICs determined by the microdilution method.[1]

Each AST card contains a control well containing only microbiological culture medium. The remaining microwells contain premeasured amounts of specific antimicrobials combined with culture medium.

The organism suspension to be tested must be diluted to a standardized concentration in 0.45% saline before being used to rehydrate the antimicrobial medium within the card. The card is then filled, sealed, and placed into the instrument incubator/reader, either automatically (as with VITEK® 2 60 or VITEK® 2 XL) or manually (as with VITEK® 2 Compact). The instrument monitors the growth of each well in the card over a defined period of time (up to 18 hours for Gram-negatives). At the completion of the incubation cycle, MIC values (or test results, as appropriate) are determined for each antimicrobial contained on the card.

PRECAUTIONS

- For *In Vitro* Diagnostic Use Only.
- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For professional use only.

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- Suspensions not within the appropriate zone on the VITEK® 2 DensiCHEK™ Plus or VITEK® 2 DensiCHEK™ may compromise card performance.
- Do not use the card after the expiration date shown on the package liner.
- Store the card unopened in the package liner. Do not use the card if the protective package liner is damaged or if no desiccant is present.
- Allow the card to come to room temperature before opening the package liner.
- Do not use powdered gloves. Powder may interfere with the optics.
- Use of culture media other than the recommended types must be validated by the customer laboratory for acceptable performance.
- The card performs as intended only when used in conjunction with VITEK® 2 Systems, following the instructions contained in the Instructions for Use.
- It is highly recommended that Good Laboratory Practices (e.g., FDA, CLSI, ISO, etc.) also be followed, according to local guidelines or requirements.
- **Do not use glass test tubes**. Use clear plastic (polystyrene) tubes only. Variation exists among test tubes of standard diameter. Carefully place the tube into the cassette. If resistance is encountered, discard and try another tube that does not require pressure to insert.
- Prior to inoculation, inspect cards for tape tears or damage to the tape and discard any that are suspect. Check the saline level in the tubes after the cassette has been processed to ensure proper filling of card.
 - VITEK® 2 60 or VITEK® 2 XL: Eject improperly filled cards.
 - VITEK® 2 Compact: Do not load improperly filled cards.
- Give special consideration to specimen source and patient therapy regimen. AST cards may contain some antimicrobials that are not proven to be effective for treatment of infections due to all organisms that may be tested. For interpreting and reporting of antimicrobial results that have been shown to be active against organism groups both *in vitro* and in clinical infections, refer to the individual pharmaceutical antimicrobial labeling or local therapy guidelines.
- Interpretation of test results requires the judgment and skill of a person knowledgeable in AST. Additional testing may be required.[17]

Warning: All patient specimens, microbial cultures, and inoculated VITEK® 2 cards, along with associated materials, are potentially infectious and should be treated with universal precautions [18,20]

REAGENTS

When used with VITEK® 2 instrumentation, the AST card is a complete system for routine susceptibility testing. Each AST card contains selected antimicrobials in varying concentrations, dried with a microbiological culture medium.

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1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.
2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).
3. Use a sterile stick or swab to transfer a sufficient number of morphologically similar colonies to the saline tube prepared in step 2. Prepare a homogenous organism suspension with a density equivalent to the appropriate McFarland standard using the VITEK® 2 DensiCHEK™ Plus (see the Culture Requirements Table). The inoculum density required for the GN, GP, ST, or YST can be diluted automatically for AST cards by the VITEK® 2 instrument (VITEK® 2 60 or VITEK® 2 XL). **Note:** The age of the suspension before loading the instrument for AST testing must be less than one hour when using VITEK® 2 60 or VITEK® 2 XL, and less than 30 minutes when using VITEK® 2 Compact.
4. Choose one of the following:
 - **For an automatic dilution (VITEK® 2 60 or VITEK® 2 XL only):** Place the suspension tube prepared in step 3 into the cassette with or without an identification card. In the next cassette slot, place an empty tube and an AST card. The instrument will automatically dilute the bacterial suspension to prepare an inoculum suitable for the susceptibility card.
 - **For a manual dilution (VITEK® 2 Compact, VITEK® 2 60 or VITEK® 2 XL):** In a second tube containing 3.0 mL of saline, transfer 145 µL of the suspension prepared in step 3. Then place this tube in the cassette with a susceptibility card. The tube with the initial bacterial suspension can also be used for inoculation of an identification card.

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- amoxicillin/clavulanic acid (amc01n) (µg/mL) 2/1, 4/2, 8/4, 16/8, 32/16
- piperacillin/tazobactam (µg/mL) 4/4, 8/4, 16/4, 32/4, 64/4, 128/4
- trimethoprim/sulfamethoxazole: **Note exception** — This drug is listed on the laboratory and patient reports as the sum of the two antimicrobial concentrations: 20 µg/mL = 1/19, 40 µg/mL = 2/38, 80 µg/mL = 4/76, 160 µg/mL = 8/152, 320 µg/mL = 16/304

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Indications for use specific to the FDA are found in each VITEK® 2 AST card package insert in the column labeled, "FDA Indications for Use". This list contains the antimicrobial/organism combinations that have been cleared by the FDA for testing and reporting on the VITEK® 2 System. FDA clearance has been granted according to the FDA approved pharmaceutical labeling and the VITEK® 2 AST clinical trial data. To report only organisms listed in the FDA Indications for Use section of the package insert, enable bioART Indications for Use suppression rules.

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Perform an alternative method of testing prior to reporting of results when an intermediate or resistant result is obtained with the following antibiotic/organism combination(s):

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Category agreement (CA) occurs when the VITEK® 2 and the reference interpretative results agree (Susceptible, Intermediate, and Resistant). There are instances when the category agreement for an antimicrobial falls below the essential

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agreement. This can occur when a significant number of MICs cluster around a category breakpoint during clinical trial testing, resulting in interpretative errors. For a description of interpretive errors, refer to the footnotes below the table that follows (Performance Characteristics). When the majority of the errors are of the minor type, a high corresponding essential agreement percentage demonstrates that the antimicrobial retains an acceptable overall performance.

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- *Leclercia adecarboxylata**
- *Lelliottia amnigena* (formerly known as Enterobacter amnigenus)*
- *Lelliottia amnigena 1* (formerly known as Enterobacter amnigenus 1)*
- *Lelliottia amnigena 2* (formerly known as Enterobacter amnigenus 2)
- *Mannheimia haemolytica*
- *Moraxella* group
- *Moraxella lacunata*
- *Moraxella nonliquefaciens*
- *Moraxella osloensis*
- *Morganella morganii**
- *Morganella morganii* ssp. *morganii**
- *Morganella morganii* ssp. *sibonii**
- *Myroides* ssp.
- *Pantoea agglomerans*
- *Pantoea dispersa*
- *Pasteurella aerogenes*
- *Pasteurella multocida**
- *Pasteurella pneumotropica*
- *Plesiomonas shigelloides*
- *Pluralibacter gergoviae* (formerly known as *Enterobacter gergoviae*)*
- *Proteus hauseri**
- *Proteus mirabilis**
- *Proteus penneri**
- *Proteus vulgaris**
- *Providencia alcalifaciens**
- *Providencia rettgeri**
- *Providencia rustigianii**
- *Providencia stuartii**
- *Pseudomonas aeruginosa**
- *Pseudomonas aeruginosa* ATCC® 27853™
- *Pseudomonas alcaligenes*
- *Pseudomonas fluorescens*
- *Pseudomonas luteola*
- *Pseudomonas mendocina*
- *Pseudomonas oleovorans*
- *Pseudomonas oryzihabitans*
- *Pseudomonas putida*
- *Pseudomonas* spp.
- *Pseudomonas stutzeri*
- *Ralstonia pickettii*
- *Raoultella ornithinolytica**
- *Raoultella planticola**
- *Raoultella terrigena**
- *Salmonella enterica* ssp. *arizonae*
- *Salmonella enterica* ssp. *enterica**
- *Salmonella* group*
- *Salmonella* ser. Enteritidis*
- *Salmonella* ser. Paratyphi A*
- *Salmonella* ser. Paratyphi B*
- *Salmonella* ser. Paratyphi C*
- *Salmonella* ser. Typhi*
- *Salmonella* ser. Typhimurium*

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REFERENCES

1. Barry, AL The Antimicrobic Susceptibility Test, Principles and Practices, Lea and Febiger, Philadelphia, PA. 1976.
2. Clinical Laboratory Standards Institute (CLSI®), Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically, M7- A7, Wayne, Pennsylvania, January 2006.
3. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Eighteenth Informational Supplement, M100-S18, Vol. 27, No. 1, January 2008.
4. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-third Informational Supplement, M100-S22, January 2012.
5. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fourth Informational Supplement; M100-S24, January 2014.
6. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fifth Informational Supplement, M100-S25, January 2015.
7. Comité de l'Antibiogramme de la Société Française de Microbiologie. Communiqué 1996. Path Biol, 1996, 44, n° 8, I-VIII.
8. Comité de l'Antibiogramme de la Société Française de Microbiologie, Communiqué 2007.

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9. Comite de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM), Recommendations 2012.
10. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2014.
11. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2015.
12. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 2.0, January 2012.
13. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 4.0, January 2014.
14. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 5.0, January 2015.
15. Gerlach, EH Microdilution 1: A Comparative Study, p. 63-76, In: Balows, A. (ed.), Current Techniques for Antibiotic Susceptibility Testing, Charles C. Thomas, Springfield, IL. 1974.
16. MacLowry, JD, and HH Marsh. 1968. Semi-automatic microtechnique for serial dilution antibiotic sensitivity testing in the clinical laboratory. J. Lab. Clin. Med. 1968;72:685-687.
17. Murray, PR, Baron EJ, Pfaller MA, Tenover FC, and Yolken RH, editors. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue – Approved Guideline (1997).
19. National Committee for Clinical Laboratory Standards, Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard — Third Edition, M27-A3, Vol. 22, No. 15, 2008.
20. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.

Permission to incorporate portions of M100 (Performance Standards for Antimicrobial Susceptibility Testing: Informational Supplement) in the bioMérieux clinical microbiology instrumentation and System has been granted by CLSI®. The current standard and supplements to it may be obtained from CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087, USA.

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Release Date	Part Number	Change Type	Change Summary
2018-01	045516-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• Combined product package insert content with VITEK® 2 Product Information Manual AST content • Updated Limited Warranty section • Updated with RX only information • Updated with EUCAST limitations

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Rx only IVD

INTENDED USE

The VITEK® 2 Gram-positive Susceptibility Card is intended for use with the VITEK® 2 Systems in clinical laboratories as an *in vitro* test to determine the susceptibility of *Staphylococcus* spp., *Enterococcus* spp., and *S. agalactiae* to antimicrobial agents when used as instructed.

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Methicillin-Resistant Staphylococci (MRS)

Resistance to oxacillin is used to detect the presence of MRS. Most MRS are also usually resistant to multiple antimicrobials, including other beta-lactams, aminoglycosides, macrolides, clindamycin, and tetracycline. However, a “community-acquired methicillin-resistant *S. aureus*” phenotype has been described. These are not multi-drug resistant strains, but are typically resistant to penicillin and oxacillin, susceptible or resistant to erythromycin, and susceptible to gentamicin, clindamycin, and tetracycline.

Oxacillin (OX1)

This test provides an MIC determination and category interpretation for oxacillin. Interpretive results for MRS are reported as resistant to all beta-lactams (except anti-MRSA beta-lactams) when the instrument is operating in the CLSI® mode or the User-Defined (based on CLSI®) mode. The results from this test correlate to results that would be obtained from standard dilution testing of oxacillin. The calling range is 0.25 µg/mL to 4.0 µg/mL.

Cefoxitin Screen

This test may be used to predict *mec*A-mediated oxacillin resistance, and it is based on the cefoxitin disk screen test. The cefoxitin screen and oxacillin work in combination to determine the final interpretation reported for oxacillin.

Synergy Screen

Since the use of penicillin or ampicillin alone often results in treatment failure of serious enterococcal endocarditis, combination therapy is usually indicated to enhance bactericidal activity. The synergy between a cell wall active agent (such as penicillin, ampicillin, or vancomycin) and an aminoglycoside (such as gentamicin, kanamycin or streptomycin) is best predicted for enterococci by screening for high-level resistance to the aminoglycoside. When enterococci are susceptible *in vitro* to the high-level aminoglycoside and a cell wall active agent, this is predictive of the effectiveness of this combination therapy. Results are reported as SYN-S (high-level synergy screen is susceptible) and SYN-R (high-level synergy screen is resistant).

Inducible Clindamycin Resistance Test (ICR)

A positive ICR test is indicative of inducible MLS_b resistance, which confers resistance to macrolides, lincosamides, and type B streptogramin. An isolate with a positive ICR test should be reported as resistant to clindamycin; however, clindamycin may

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still be effective in some patients. If the ICR test is positive and the clindamycin result is susceptible or intermediate, the clindamycin result will be forced resistant by the ICR test (in CLSI® mode or User-defined based on CLSI®).

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The organism suspension to be tested must be diluted to a standardized concentration in 0.45% saline before being used to rehydrate the antimicrobial medium within the card. The card is then filled, sealed, and placed into the instrument incubator/reader, either automatically (as with VITEK® 2 60 or VITEK® 2 XL) or manually (as with VITEK® 2 Compact). The instrument monitors the growth of each well in the card over a defined period of time (up to 18 hours for *Staphylococcus, Enterococcus* and *S. agalactiae*). At the completion of the incubation cycle, MIC values (or test results, as appropriate) are determined for each antimicrobial contained on the card.

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Numerical values are expressed in µg/mL.

§ Equivalent standard method concentration by efficacy.

NEG = Negative

POS = Positive

❶, ❷ etc. = See performance characteristics identified by the drug code with this symbol.

**N/A = No specific FDA Indications for Use available

**MSSA = Methicillin-susceptible *S. aureus*

[NS] = The current absence of resistant isolates precludes defining any results other than susceptible. Isolates yielding MIC results suggestive of Nonsusceptible category should be submitted to a reference laboratory for further testing.

PRECAUTIONS

- For *In Vitro* Diagnostic Use Only.
- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For professional use only.
- Suspensions not within the appropriate zone on the VITEK® 2 DensiCHEK™ Plus or VITEK® 2 DensiCHEK™ may compromise card performance.
- Do not use the card after the expiration date shown on the package liner.
- Store the card unopened in the package liner. Do not use the card if the protective package liner is damaged or if no desiccant is present.
- Allow the card to come to room temperature before opening the package liner.
- Do not use powdered gloves. Powder may interfere with the optics.
- Use of culture media other than the recommended types must be validated by the customer laboratory for acceptable performance.
- The card performs as intended only when used in conjunction with VITEK® 2 Systems, following the instructions contained in the Instructions for Use.
- It is highly recommended that Good Laboratory Practices (e.g., FDA, CLSI, ISO, etc.) also be followed, according to local guidelines or requirements.
- **Do not use glass test tubes**. Use clear plastic (polystyrene) tubes only. Variation exists among test tubes of standard diameter. Carefully place the tube into the cassette. If resistance is encountered, discard and try another tube that does not require pressure to insert.
- Prior to inoculation, inspect cards for tape tears or damage to the tape and discard any that are suspect. Check the saline level in the tubes after the cassette has been processed to ensure proper filling of card.
 - VITEK® 2 60 or VITEK® 2 XL: Eject improperly filled cards.
 - VITEK® 2 Compact: Do not load improperly filled cards.

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- Give special consideration to specimen source and patient therapy regimen. AST cards may contain some antimicrobials that are not proven to be effective for treatment of infections due to all organisms that may be tested. For interpreting and reporting of antimicrobial results that have been shown to be active against organism groups both *in vitro* and in clinical infections, refer to the individual pharmaceutical antimicrobial labeling or local therapy guidelines.
- Interpretation of test results requires the judgment and skill of a person knowledgeable in AST. Additional testing may be required.[17]

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INSTRUMENT

The VITEK® 2 instruments are a family of *in vitro* diagnostic devices intended to rapidly assess the antimicrobial susceptibility of bacterial and yeast pathogens to available antimicrobial agents. For detailed information on the use and operation of these devices, refer to the appropriate Instrument User Manual.

SPECIMEN PREPARATION

Culture Requirements Table

VITEK® 2 Card	Media	Age of Culture	Incubation Conditions	McFarland Standards	Dilution for AST	Age of Suspension Before Loading Instrument
AST Gram-Positive	TSAB CBA CPS ID	18 to 24 hours	35°C to 37°C 5% to 10% CO_2 or aerobic, non-CO_2	0.50 to 0.63	280 µL in 3.0 mL 0.45% saline	VITEK® 2 Compact: ≤30 minutes VITEK® 2: ≤ 1 hour
GP and AST GP pair	TSAB[1] CBA[1] CPS ID	18 to 24 hours	35°C to 37°C 5% to 10% CO_2 or aerobic, non-CO_2	0.50 to 0.63	280 µL in 3.0 mL 0.45% saline	≤ 30 minutes

[1] These media were used in the identification product database developments and will give optimal performance.

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TEST PROCEDURE

Warning: Failure to properly follow instructions and recommendations provided in this section for performing laboratory tasks may cause erroneous or delayed results.

Materials

Materials provided are:

- VITEK® 2 DensiCHEK™ Plus Kit
- DensiCHEK™ Plus Standards Kit
- VITEK® 2 Cassette
- Adjustable volume saline dispenser
- 12 mm x 75 mm clear plastic (polystyrene) disposable test tubes
- VITEK® 2 60 or VITEK® 2 XL only: VITEK® 2 Pipettor/Diluter Accessory Kit (containing instrument pipette tips and saline hookup) and 0.45% saline bag

Materials required, but not provided are:

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- Sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- Loops, sterile sticks or swabs
- Appropriate agar medium (See the Culture Requirements Table.)
- QC isolates
- VITEK® 2 AST Cards
- Micropipettors to deliver 280 µL
- Disposable pipette tips

Optional Accessories:

- Pre-dispensed saline test tubes (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- Test tube caps
- Vortex

Test Card Setup Procedure

The following procedure contains general information applicable to all susceptibility products. (See the Culture Requirements Table for product-specific information.)

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1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.
2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).
3. Use a sterile stick or swab to transfer a sufficient number of morphologically similar colonies to the saline tube prepared in step 2. Prepare a homogenous organism suspension with a density equivalent to the appropriate McFarland standard using the VITEK® 2 DensiCHEK™ Plus (see the Culture Requirements Table). The inoculum density required for the GN, GP, ST, or YST can be diluted automatically for AST cards by the VITEK® 2 instrument (VITEK® 2 60 or VITEK® 2 XL).
 Note: The age of the suspension before loading the instrument for AST testing must be less than one hour when using VITEK® 2 60 or VITEK® 2 XL, and less than 30 minutes when using VITEK® 2 Compact.
4. Choose one of the following:
 - **For an automatic dilution (VITEK® 2 60 or VITEK® 2 XL only):** Place the suspension tube prepared in step 3 into the cassette with or without an identification card. In the next cassette slot, place an empty tube and an AST card. The instrument will automatically dilute the bacterial suspension to prepare an inoculum suitable for the susceptibility card.
 - **For a manual dilution (VITEK® 2 Compact, VITEK® 2 60 or VITEK® 2 XL):** In a second tube containing 3.0 mL of saline, transfer 280 µL of the suspension prepared in step 3. Then place this tube in the cassette with a susceptibility card. The tube with the initial bacterial suspension can also be used for inoculation of an identification card.

 Note: Check the saline level in the tubes after filling. When it is evident by the saline level in the tube that a card has been improperly filled, do not load the card if using VITEK® 2 Compact; **or**, eject the card if using VITEK® 2 60 or VITEK® 2 XL.

 Note: Refer to the appropriate Instrument User Manual for detailed information regarding data entry, processing, etc.
5. Follow your local inspecting agency's guidelines for disposal of hazardous waste.

QUALITY CONTROL

Quality Control organisms should be processed according to the Test Card Setup Procedure.

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Numerical values are expressed in µg/mL.

NEG = Negative

POS = Positive

[NS] = The current absence of resistant isolates precludes defining any results other than susceptible. Isolates yielding MIC results suggestive of Nonsusceptible category should be submitted to a reference laboratory for further testing.

Certification Statement

This is to certify that bioMérieux complies with ISO 13485 and FDA Quality System Regulation (QSR) requirements for design, development, and manufacture of antimicrobial susceptibility systems.

Frequency of QC Testing

Refer to *Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically,* CLSI® and/or your local guidelines.[2]

Preparation of QC Organisms

1. Rehydrate the organism according to the manufacturer's instructions.
2. Subculture to Trypticase Soy Agar with 5% sheep blood (TSAB).
3. Incubate at 35°C for 24 hours.
 Note: Gram-positive organisms may require an atmosphere of CO_2. (Refer to the Culture Requirements Table).
4. Check for purity.

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5. Subculture to a TSAB plate.
6. Incubate 16 to 18 hours at 35°C.

Short-Term Storage Conditions

1. Streak to a TSAB plate or slant.
2. Incubate for 24 hours.
3. Refrigerate at 2°C to 8°C for up to two weeks.
4. Subculture once as described above and use for QC.

Long-Term Storage Conditions

1. Make a heavy suspension in Tryptic Soy Broth (TSB) with 15% glycerol.
2. Freeze at -70°C.
3. Subculture to TSAB twice before running QC.
 Note: Avoid repeated thawing and refreezing by either freezing in single-use aliquots or removing a small portion of frozen organism preparation with a sterile applicator stick.

RESULTS

Susceptibility Analytical Techniques

The system evaluates each organism's growth pattern in the presence of the antimicrobial in relation to the growth control well. Several parameters based on the growth characteristics are used to determine the MIC or qualitative result (for example, ESBL POS/NEG). The MIC result must be linked to an organism identification to determine a category interpretation. Accurate identification is critical, especially with certain organism/antimicrobial combinations (e.g., *Staphylococcus aureus*/oxacillin).

In cases where the identification of an organism is in question, confirmatory testing is necessary to ensure correct interpretation of susceptibility results.

A category interpretation will be reported along with a MIC, according to the interpretations defined by the Food and Drug Administration (FDA), CLSI®, Comité de l'Antibiogramme de la Société Française de Microbiologie (CASFM), European Committee for Antimicrobial Susceptibility Testing (EUCAST), or to an adaptation of the global settings according to other local guidelines.

Note: When FDA and CLSI® breakpoints differ, VITEK® 2 Systems AST tests are cleared for use with FDA breakpoints applied.

Combination Antimicrobials

The MICs for the combination antimicrobials are listed on the laboratory and patient reports as the first concentration (e.g., ampicillin/sulbactam ≤ 8/4 µg/mL is reported as ≤ 8 µg/mL.). The actual concentrations for each value in the antimicrobial calling range are as follows:

- trimethoprim/sulfamethoxazole: **Note exception** — This drug is listed on the laboratory and patient reports as the sum of the two antimicrobial concentrations: 20 µg/mL = 1/19, 40 µg/mL = 2/38, 80 µg/mL = 4/76, 160 µg/mL = 8/152, 320 µg/mL = 16/304

Antimicrobial Deduction

Antimicrobials that have been deduced will only report an interpretive result and will be noted with a **+**.

Clinical Efficacy and Indications for Use

AST cards may contain some antimicrobials that are not proven to be effective for treatment of infections caused by all organisms that may be tested. For interpreting and reporting of antimicrobial results that have shown to be active against organism groups both *in vitro* and in clinical infections, refer to the individual pharmaceutical antimicrobial labeling or the local therapy guidelines.

Urinary Use Only Antimicrobials

Certain antimicrobial agents are limited to use in treating urinary tract infections. Accordingly, these agents should not be reported against pathogens recovered from infection sites other than the urinary tract (except as noted). Please refer to the current *CLSI Performance Standards for Susceptibility Testing, M100* and/or local guidelines for additional information.[3]

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LIMITATIONS

A VITEK® 2 AST card cannot be used with a direct clinical specimen or sample or other sources containing mixed flora. Any change or modification in the procedure may affect the results.

A result for an antibiotic/organism combination, which may have a limitation, may be suppressed from reporting. This can be accomplished through the use of bioART rules in the VITEK® 2 Systems software. Refer to the software user manual for instructions.

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):

- Erythromycin (e02n): *Streptococcus agalactiae*
- Cefoxitin Screen (oxsf01n): *Staphylococcus pseudintermedius* (specific to VITEK® 2 Systems 8.01 software or higher)

Perform an alternative method of testing prior to reporting of results when a Positive (+) result is obtained with the following antibiotic/organism combination(s):

- Cefoxitin Screen (oxsf01n): *Staphylococcus saprophyticus*

The ability of the AST card to detect resistance with the following combination(s) is unknown because resistant strains were not available at the time of comparative testing:

- Linezolid (lnz02n): *Enterococcus* spp., *Staphylococcus* spp., *Streptococcus agalactiae*
- Tigecycline (tgc02n): *Enterococcus* spp., *Staphylococcus* spp., *Streptococcus* spp.

EUCAST Limitations

It is recommended that the existing bioART suppression rules be enabled, or new rules created and enabled, for these limitations if the EUCAST breakpoints are applied.

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):

- Erythromycin (e02n): *Streptococcus agalactiae, Staphylococcus* spp.
- Teicoplanin (tec01n): *Staphylococcus epidermidis, Staphylococcus haemolyticus*
- Cefoxitin Screen (oxsf01n): *Staphylococcus pseudintermedius* (specific to VITEK® 2 Systems 8.01 software or higher)

Perform an alternative method of testing prior to reporting of results when a Positive (+) result is obtained with the following antibiotic/organism combination(s):

- Cefoxitin Screen (oxsf01n): *Staphylococcus saprophyticus*

The ability of the AST card to detect resistance with the following combination(s) is unknown because resistant strains were not available at the time of comparative testing:

- Linezolid (lnz02n): *Streptococcus agalactiae*
- Nitrofurantoin (ft01n): *Staphylococcus saprophyticus*
- Tigecycline (tgc02n): *Streptococcus agalactiae, Staphylococcus* spp.

EXPECTED VALUES

Expected results for susceptibility tests will vary based on location and institution. VITEK® 2 Systems were tested at several geographically diverse locations to ensure that trends that occurred by location were integrated into the performance characteristics of the system. Organism resistance patterns will differ by institution; therefore, expected values will be directly related to the population of organisms at each site.

PERFORMANCE CHARACTERISTICS

The performance characteristics of the antimicrobial agents included in VITEK® 2 AST cards were established using the manual and autodilution modes (on a VITEK® 2 System) at multiple clinical laboratories. The VITEK® 2 AST card results were compared to results from a CLSI® reference method. Essential agreement (EA) represents VITEK® 2 results which agree exactly or are within a ± twofold dilution (± two doubling dilutions for antifungal) of the reference result.

Category agreement (CA) occurs when the VITEK® 2 and the reference interpretative results agree (Susceptible, Intermediate, and Resistant). There are instances when the category agreement for an antimicrobial falls below the essential agreement. This can occur when a significant number of MICs cluster around a category breakpoint during clinical trial testing, resulting in interpretative errors. For a description of interpretive errors, refer to the footnotes below the table that

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follows (Performance Characteristics). When the majority of the errors are of the minor type, a high corresponding essential agreement percentage demonstrates that the antimicrobial retains an acceptable overall performance.

There are instances when the performance is based solely on category agreement (CA) because at the time performance was established, less than five discrete twofold dilutions were evaluated. A minimum of five dilutions is necessary to calculate essential agreement (EA) based on ± one twofold dilution. These instances are designated by a "c" footnote in the Contents of the Card table. The following performance tables have values for CA only when EA was not established at the time of FDA clearance.

The reproducibility of the VITEK® 2 system was established by testing a set of on-scale organisms.*

*Data on file at bioMérieux, Inc

Performance Characteristics for Gram-Positive Antimicrobial Susceptibility Testing

Antimicrobial	Antimicrobial Code	Antibiotic Version	Bp[1]	Comment[2]	Essential Agreement				Category Agreement				% Reproducibility
					% Error				% Error				
					% EA	VME	ME	mE	% CA	VME	ME	mE	
Benzylpenicillin (Penicillin)	P	p04n	CLSI	#, E	97.4	0.4	0.3	N/A	98.5	0.9	1.0	N/A	99.3
Cefoxitin Screen	OXSF	oxsf01n	CLSI	E, Ref. = CLSI Disk Diffusion	–	–	–	–	98.3	2.0	1.4	N/A	100
				E, Ref. = mecA PCR	–	–	–	–	97.2	2.3	3.2	N/A	
Oxacillin + Cefoxitin Screen			CLSI	#, Ref. = CLSI Cefoxitin Disk Diffusion	–	–	–	–	98.5	0.4	2.8	N/A	–
Oxacillin + Cefoxitin Screen			CLSI	E, Ref. = mecA PCR	–	–	–	–	97.8	1.2	3.2	N/A	–
Clindamycin	CM	cm01n	CLSI	#, E	98.7	2.3	0.0	0.7	99.1	2.3	0.0	0.9	98.5
			CA-SFM	I	96.1	2.5	4.3	0.0	-	-	-	-	
Erythromycin	E	e02n	CLSI	#, E ❶	95.9	0.4	0.4	5.1	92.8	0.4	0.4	7.5	95.2
			CA-SFM	I ❶	94.7	9.4	0.0	0.3	95.0	9.4	0.0	0.3	
Fosfomycin	FOS	fos01n	CA-SFM	I, Staph	97.1	6.5	0.3	N/A	95.8	11.8	1.7	N/A	100
Fusidic Acid	FA	fa01n	CA-SFM	I, Staph	100	0.0	0.0	0.0	99.2	0.0	0.0	0.8	100
Gentamicin	GM	gm01n	CLSI	#, E, Staph	99.2	0.0	0.0	0.8	95.1	0.0	0.0	4.9	100
			CA-SFM	E, Staph	99.2	0.0	0.0	0.8	95.1	0.0	0.0	4.9	
Inducible Clindamycin Resistance	ICR	icr02n	CLSI	#, E	N/A	N/A	N/A	N/A	99.5	2.0	0.4	N/A	100
Levofloxacin	LEV	lev01n	CLSI	#, E	99.3	0.0	0.0	0.2	95.8	0.0	0.0	4.2	100
Linezolid	LNZ	lnz02n	CLSI	#, E	98.7	0.0	0.1	0.1	98.9	0.0	0.1	1.0	100
			CA-SFM	E	98.7	0.0	0.1	0.7	92.6	0.0	0.1	7.3	
Moxifloxacin	MXF	mxf02n	CLSI	#, E	99.2	0.0	0.0	0.0	83.3	0.0	0.0	16.7	100
			CA-SFM	E	98.5	0.0	0.0	0.1	92.4	0.0	0.0	7.6	
Mupirocin	MUP	mup01n	CA-SFM	E ❶	99.7	0.0	0.0	0.0	100	0.0	0.0	0.0	100
			Global	E ❶	99.7	0.0	0.3	0.0	99.2	0.0	0.9	0.0	
Nitrofurantoin	FT	ft01n	CLSI	#, E	98.3	0.0	1.4	0.0	91.7	0.0	1.4	7	100
Oxacillin	OX1	ox101n	CLSI	#, E, Staph	97.4	1.6	1.1	0.0	97.1	2.2	3.4	0.0	98.1
Rifampicin (Rifampin)	RA	ra01n	CLSI	#, E	100	0.0	0.0	0.0	100	0.0	0.0	0.0	100
			CA-SFM	I, Staph	99.5	0.0	0.0	0.3	–	–	–	–	

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[1] Abbreviations — Bp = breakpoint committee; EA = essential agreement; CA = category agreement; VME = Very Major Error (susceptible result with resistant reference result); ME = Major Error (resistant result with susceptible reference result); mE = minor Error (susceptible or resistant result with an intermediate reference result, or an intermediate result with a susceptible or resistant reference result).

[2] Comment — Specific organism groups are designated as Staph for *staphylococci*, Enc for *enterococci*, S aga for Group B *streptococci*, Sau for *Staphylococcus aureus*.

[3] FDA breakpoints are used in the CLSI Interpretation Standard (breakpoint committee) in the VITEK® 2 Systems software.

For all gram-positive Beta-lactams with indications for *Staphylococcus* spp., with the exception of Benzylpenicillin (Penicillin), performance in this table represents forcing-resistance based on oxacillin resistance.

† Non-defined error; susceptible (S) breakpoints only, no I or R breakpoints exist.

Key:

= US Food and Drug Administration 510(k) cleared
CLSI® = Clinical and Laboratory Standards Institute
CA-SFM = Comité de l'Antibiogramme de la Société Française de Microbiologie
E = External performance data
I = Internal performance data
– = Not available
N/A = Not applicable
Ref. = Reference method for clinical performance study.
❶ ❷ = Symbol identifies the performance characteristics for a specific antimicrobial version.

EUCAST Performance Characteristics for Gram-Positive Antimicrobial Susceptibility Testing

Antimicrobial	Antimicrobial Code	Antibiotic Version	Comment[1]	Essential Agreement				Category Agreement			
				% Error				% Error			
				%EA[2]	VME	ME	mE	%CA	VME	ME	mE
Benzylpenicillin	P	p04n	*Staphylococcus, S. agalactiae*	98.3	0.8	1.9	0.0	98.3	0.8	3.9	0.0
Clindamycin	CM	cm01n	*Staphylococcus, S. agalactiae*	98.1	3.3	0.3	0.0	98.1	3.3	0.3	0.7
Erythromycin	E	e02n	*Staphylococcus* ❶	94.7	9.4	0.0	0.0	95.3	9.4	0.0	0.0
Fosfomycin	FOS	fos01n	*Staphylococcus*	97.1	6.5	0.3	N/A	95.8	11.8	1.7	N/A
Fusidic Acid	FA	fa01n	*Staphylococcus*	100.0	0.0	0.0	0.0	100.0	0.0	0.0	0.0
Gentamicin	GM	gm01n	*Staphylococcus*	99.2	0.0	0.0	0.0	99.7	0.0	0.4	0.0
Levofloxacin	LEV	lev01n	N/A	99.7	1.4	0.0	0.0	98.4	1.4	0.0	1.3

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[1] Comment — Unless stated otherwise, performance is for *Staphylococcus*, *Enterococcus*, and *S. agalactiae*.

[2] Abbreviations — Bp = breakpoint committee; EA = essential agreement; CA = category agreement; VME = Very Major Error (susceptible result with resistant reference result); ME = Major Error (resistant result with susceptible reference result); mE = minor Error (susceptible or resistant result with an intermediate reference result, or an intermediate result with a susceptible or resistant reference result).

Key:

= EUCAST = European Committee on Antimicrobial Susceptibility
❶ ❷ = Symbol identifies the performance characteristics for a specific antimicrobial version.

LIST OF CLAIMS

Note: If the organism is not in the VITEK® 2 susceptibility database, results will not be reported.

Note: Organisms listed with an asterisk (*) indicate an AES claimed organism. No asterisk is displayed for a group; however, when an individual species (with an asterisk) is contained within a group, it is expertised.

Gram-Positive Organisms Claimed for AST-GP (keyID)

- *Coagulase-negative Staphylococcus**
- *Coagulase-positive Staphylococcus**
- *Enterococcus avium*
- *Enterococcus casseliflavus**
- *Enterococcus durans*
- *Enterococcus faecalis**
- *Enterococcus faecalis* ATCC® 29212™
- *Enterococcus faecalis* ATCC® 51299™
- *Enterococcus faecium**

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REFERENCES

1. Barry, AL The Antimicrobic Susceptibility Test, Principles and Practices, Lea and Febiger, Philadelphia, PA. 1976.

2. Clinical Laboratory Standards Institute (CLSI®), Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically, M7- A7, Wayne, Pennsylvania, January 2006.

3. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Eighteenth Informational Supplement, M100-S18, Vol. 27, No. 1, January 2008.

4. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-third Informational Supplement, M100-S22, January 2012.

5. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fourth Informational Supplement; M100-S24, January 2014.

6. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fifth Informational Supplement, M100-S25, January 2015.

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11. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2015.
12. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 2.0, January 2012.
13. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 4.0, January 2014.
14. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 5.0, January 2015.
15. Gerlach, EH Microdilution 1: A Comparative Study, p. 63-76, In: Balows, A. (ed.), Current Techniques for Antibiotic Susceptibility Testing, Charles C. Thomas, Springfield, IL. 1974.
16. MacLowry, JD, and HH Marsh. 1968. Semi-automatic microtechnique for serial dilution antibiotic sensitivity testing in the clinical laboratory. J. Lab. Clin. Med. 1968;72:685-687.
17. Murray, PR, Baron EJ, Pfaller MA, Tenover FC, and Yolken RH, editors. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue – Approved Guideline (1997).
19. National Committee for Clinical Laboratory Standards, Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard — Third Edition, M27-A3, Vol. 22, No. 15, 2008.
20. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.

Permission to incorporate portions of M100 (Performance Standards for Antimicrobial Susceptibility Testing: Informational Supplement) in the bioMérieux clinical microbiology instrumentation and System has been granted by CLSI®. The current standard and supplements to it may be obtained from CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087, USA.

BARCODES

User MUST enter the following barcodes into "Flex Panel Entry" program before first use of this Susceptibility Card.

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Release Date	Part Number	Change Type	Change Summary
2018-01	045573-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• Combined product package insert content with VITEK® 2 Product Information Manual AST content • Updated Limited Warranty section • Updated with RX only information • Updated with EUCAST limitations

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PRECAUTIONS

- For *In Vitro* Diagnostic Use Only.
- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For professional use only.
- Suspensions not within the appropriate zone on the VITEK® 2 DensiCHEK™ Plus or VITEK® 2 DensiCHEK™ may compromise card performance.
- Do not use the card after the expiration date shown on the package liner.
- Store the card unopened in the package liner. Do not use the card if the protective package liner is damaged or if no desiccant is present.
- Allow the card to come to room temperature before opening the package liner.
- Do not use powdered gloves. Powder may interfere with the optics.
- Use of culture media other than the recommended types must be validated by the customer laboratory for acceptable performance.
- The card performs as intended only when used in conjunction with VITEK® 2 Systems, following the instructions contained in the Instructions for Use.
- It is highly recommended that Good Laboratory Practices (e.g., FDA, CLSI, ISO, etc.) also be followed, according to local guidelines or requirements.

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- **Do not use glass test tubes**. Use clear plastic (polystyrene) tubes only. Variation exists among test tubes of standard diameter. Carefully place the tube into the cassette. If resistance is encountered, discard and try another tube that does not require pressure to insert.
- Prior to inoculation, inspect cards for tape tears or damage to the tape and discard any that are suspect. Check the saline level in the tubes after the cassette has been processed to ensure proper filling of card.
 - VITEK® 2 60 or VITEK® 2 XL: Eject improperly filled cards.
 - VITEK® 2 Compact: Do not load improperly filled cards.
- Give special consideration to specimen source and patient therapy regimen. AST cards may contain some antimicrobials that are not proven to be effective for treatment of infections due to all organisms that may be tested. For interpreting and reporting of antimicrobial results that have been shown to be active against organism groups both *in vitro* and in clinical infections, refer to the individual pharmaceutical antimicrobial labeling or local therapy guidelines.
- Interpretation of test results requires the judgment and skill of a person knowledgeable in AST. Additional testing may be required.[17]

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1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.
2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).
3. Use a sterile stick or swab to transfer a sufficient number of morphologically similar colonies to the saline tube prepared in step 2. Prepare a homogenous organism suspension with a density equivalent to the appropriate McFarland standard using the VITEK® 2 DensiCHEK™ Plus (see the Culture Requirements Table). The inoculum density required for the GN, GP, ST, or YST can be diluted automatically for AST cards by the VITEK® 2 instrument (VITEK® 2 60 or VITEK® 2 XL). **Note:** The age of the suspension before loading the instrument for AST testing must be less than one hour when using VITEK® 2 60 or VITEK® 2 XL, and less than 30 minutes when using VITEK® 2 Compact.
4. Choose one of the following:
 - **For an automatic dilution (VITEK® 2 60 or VITEK® 2 XL only):** Place the suspension tube prepared in step 3 into the cassette with or without an identification card. In the next cassette slot, place an empty tube and an AST card. The instrument will automatically dilute the bacterial suspension to prepare an inoculum suitable for the susceptibility card.
 - **For a manual dilution (VITEK® 2 Compact, VITEK® 2 60 or VITEK® 2 XL):** In a second tube containing 3.0 mL of saline, transfer 280 µL of the suspension prepared in step 3. Then place this tube in the cassette with a susceptibility card. The tube with the initial bacterial suspension can also be used for inoculation of an identification card.

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Preparation of QC Organisms

1. Rehydrate the organism according to the manufacturer's instructions.
2. Subculture to Trypticase Soy Agar with 5% sheep blood (TSAB).
3. Incubate at 35°C for 24 hours.
 Note: Gram-positive organisms may require an atmosphere of CO_2. (Refer to the Culture Requirements Table).
4. Check for purity.
5. Subculture to a TSAB plate.
6. Incubate 16 to 18 hours at 35°C.

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EXPECTED VALUES

Expected results for susceptibility tests will vary based on location and institution. VITEK® 2 Systems were tested at several geographically diverse locations to ensure that trends that occurred by location were integrated into the performance

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characteristics of the system. Organism resistance patterns will differ by institution; therefore, expected values will be directly related to the population of organisms at each site.

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Category agreement (CA) occurs when the VITEK® 2 and the reference interpretative results agree (Susceptible, Intermediate, and Resistant). There are instances when the category agreement for an antimicrobial falls below the essential agreement. This can occur when a significant number of MICs cluster around a category breakpoint during clinical trial testing, resulting in interpretative errors. For a description of interpretive errors, refer to the footnotes below the table that follows (Performance Characteristics). When the majority of the errors are of the minor type, a high corresponding essential agreement percentage demonstrates that the antimicrobial retains an acceptable overall performance.

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Ref. = Reference method for clinical performance study.

❶ ❷ = Symbol identifies the performance characteristics for a specific antimicrobial version.

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Key:

= EUCAST = European Committee on Antimicrobial Susceptibility

❶ ❷ = Symbol identifies the performance characteristics for a specific antimicrobial version.

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2. Clinical Laboratory Standards Institute (CLSI®), Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically, M7- A7, Wayne, Pennsylvania, January 2006.
3. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Eighteenth Informational Supplement, M100-S18, Vol. 27, No. 1, January 2008.
4. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-third Informational Supplement, M100-S22, January 2012.
5. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fourth Informational Supplement; M100-S24, January 2014.
6. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fifth Informational Supplement, M100-S25, January 2015.
7. Comité de l'Antibiogramme de la Société Française de Microbiologie. Communiqué 1996. Path Biol, 1996, 44, n° 8, I-VIII.
8. Comité de l'Antibiogramme de la Société Française de Microbiologie, Communiqué 2007.
9. Comite de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM), Recommendations 2012.
10. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2014.
11. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2015.
12. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 2.0, January 2012.
13. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 4.0, January 2014.
14. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 5.0, January 2015.
15. Gerlach, EH Microdilution 1: A Comparative Study, p. 63-76, In: Balows, A. (ed.), Current Techniques for Antibiotic Susceptibility Testing, Charles C. Thomas, Springfield, IL. 1974.
16. MacLowry, JD, and HH Marsh. 1968. Semi-automatic microtechnique for serial dilution antibiotic sensitivity testing in the clinical laboratory. J. Lab. Clin. Med. 1968;72:685-687.
17. Murray, PR, Baron EJ, Pfaller MA, Tenover FC, and Yolken RH, editors. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue – Approved Guideline (1997).
19. National Committee for Clinical Laboratory Standards, Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard — Third Edition, M27-A3, Vol. 22, No. 15, 2008.
20. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.

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Symbol	Meaning
	Contains sufficient for <n> tests
EC REP	Authorized representative in the European Community
Rx only	For US Only : Caution : US Federal Law restricts this device to sale by or on the order of a licensed practitioner

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REVISION HISTORY TABLE

Change type categories	
N/A	Not applicable (First publication)
Correction	Correction of documentation anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user
Note:	Minor typographical, grammar, and formatting changes are not included in the revision history.

Release Date	Part Number	Change Type	Change Summary
2018-01	045610-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• Combined product package insert content with VITEK® 2 Product Information Manual AST content • Updated Limited Warranty section • Updated with RX only information • Updated with EUCAST limitations

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- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For professional use only.
- Suspensions not within the appropriate zone on the VITEK® 2 DensiCHEK™ Plus or VITEK® 2 DensiCHEK™ may compromise card performance.
- Do not use the card after the expiration date shown on the package liner.
- Store the card unopened in the package liner. Do not use the card if the protective package liner is damaged or if no desiccant is present.
- Allow the card to come to room temperature before opening the package liner.
- Do not use powdered gloves. Powder may interfere with the optics.
- Use of culture media other than the recommended types must be validated by the customer laboratory for acceptable performance.
- The card performs as intended only when used in conjunction with VITEK® 2 Systems, following the instructions contained in the Instructions for Use.
- It is highly recommended that Good Laboratory Practices (e.g., FDA, CLSI, ISO, etc.) also be followed, according to local guidelines or requirements.
- **Do not use glass test tubes**. Use clear plastic (polystyrene) tubes only. Variation exists among test tubes of standard diameter. Carefully place the tube into the cassette. If resistance is encountered, discard and try another tube that does not require pressure to insert.
- Prior to inoculation, inspect cards for tape tears or damage to the tape and discard any that are suspect. Check the saline level in the tubes after the cassette has been processed to ensure proper filling of card.
 - VITEK® 2 60 or VITEK® 2 XL: Eject improperly filled cards.
 - VITEK® 2 Compact: Do not load improperly filled cards.
- Give special consideration to specimen source and patient therapy regimen. AST cards may contain some antimicrobials that are not proven to be effective for treatment of infections due to all organisms that may be tested. For interpreting and reporting of antimicrobial results that have been shown to be active against organism groups both *in vitro* and in clinical infections, refer to the individual pharmaceutical antimicrobial labeling or local therapy guidelines.
- Interpretation of test results requires the judgment and skill of a person knowledgeable in AST. Additional testing may be required.[17]

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[1]FDA/CLSI Broth Microdilution expected QC range = 1/4 - 4/4 µg/mL

[2]FDA/CLSI Broth Microdilution expected QC range = 1/4 - 8/4 µg/mL

[3]FDA/CLSI Broth Microdilution expected QC range = 0.5/4 - 2/4 µg/mL

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Preparation of QC Organisms

1. Rehydrate the organism according to the manufacturer's instructions.
2. Subculture to Trypticase Soy Agar with 5% sheep blood (TSAB).
3. Incubate at 35°C for 24 hours.
4. Check for purity.
5. Subculture to a TSAB plate.
6. Incubate 16 to 18 hours at 35°C.

Short-Term Storage Conditions

1. Streak to a TSAB plate or slant.
2. Incubate for 24 hours.

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3. Refrigerate at 2°C to 8°C for up to two weeks.
4. Subculture once as described above and use for QC.

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- piperacillin/tazobactam (µg/mL) 4/4, 8/4, 16/4, 32/4, 64/4, 128/4
- trimethoprim/sulfamethoxazole: **Note exception** — This drug is listed on the laboratory and patient reports as the sum of the two antimicrobial concentrations: 20 µg/mL = 1/19, 40 µg/mL = 2/38, 80 µg/mL = 4/76, 160 µg/mL = 8/152, 320 µg/mL = 16/304

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There are instances when the performance is based solely on category agreement (CA) because at the time performance was established, less than five discrete twofold dilutions were evaluated. A minimum of five dilutions is necessary to calculate

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essential agreement (EA) based on ± one twofold dilution. These instances are designated by a "c" footnote in the Contents of the Card table. The following performance tables have values for CA only when EA was not established at the time of FDA clearance.

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*Data on file at bioMérieux, Inc

Table 4: Performance Characteristics for Gram-Negative Antimicrobial Susceptibility Testing

Antimicrobial	Antimicrobial Code	Antibiotic Version	Bp[1]	Comment[2]	Essential Agreement				Category Agreement				% Reproducibility
					% Error				% Error				
					% EA	VME	ME	mE	% CA	VME	ME	mE	
Ampicillin	AM	am01n	CLSI	#, E	93.3	1.3	0.0	0.8	90.9	1.3	0.0	7.8	99.3
			CA-SFM	I	94.9	0.6	2.2	4.3	–	–	–	–	
Cefotaxime	CTX	ctx01n	CLSI	#, E	95.8	1.1	0.2	2.6	92.9	1.1	0.2	6.7	98.1
			CA-SFM	I	94.0	0.4	0.1	3.5	–	–	–	–	
Cefpodoxime	CPD	cpd01n	CLSI	#, E	96.0	0.0	0.4	0.8	96.3	0.0	0.4	3.5	100
			CA-SFM	I	94.1	1.2	2.9	1.3	–	–	–	–	
Ceftazidime	CAZ	caz01n	CLSI	#, E	97.5	1.1	0.2	1.0	97.3	1.1	0.2	3.5	95.4
			CA-SFM	I	92.3	5.1	0.0	4.1	–	–	–	–	
Cefuroxime	CXM	cxm01n	CLSI	#, E	95.3	1.5	1.0	1.6	93.9	1.5	1.0	4.8	99.6
Ciprofloxacin	CIP	cip01n	CLSI	#, E	97.9	0.0	0.0	0.5	96.6	0.0	0.0	3.8	98.9
			CA-SFM	E	97.9	0.0	0.0	0.5	96.6	0.0	0.0	3.8	
Ertapenem	ETP	etp01n	CLSI	#, E	–	–	–	–	98.8	0.0	0.0	1.2	97.8
Extended-Spectrum Beta-Lactamase (ESBL)	ESBL	esb01n	CLSI	#, E	N/A	N/A	N/A	N/A	97.7	2.7	1.9	N/A	100
Gentamicin	GM	gm01n	CLSI	#, E	97.1	1.9	0.0	0.6	96.4	1.9	0.0	2.9	100
			CA-SFM	E	97.1	1.9	0.0	0.6	96.4	1.9	0.0	2.9	
Imipenem	IPM	ipm04n	CLSI	#, E ❷	95.7	2.4	0.3	0.4	94.8	2.4	0.3	4.4	99.6
Meropenem	MEM	mem02n	CLSI	#, E ❷	97.6	1.4	0.0	0.3	96.8	1.4	0.0	2.8	99.6
			Global	#, E ❷	97.6	0.0	0.0	0.5	97.7	0.0	0.0	2.3	
Moxifloxacin	MXF	mxf01n	CLSI	#, E	98.7	0.0	0.0	0.0	97.8	2.2	0.0	0.0	100
Piperacillin	PIP	pip02n	CLSI	#, E ❷	94.6	2.0	0.2	3.7	92.4	2.7	0.2	6.8	99.3
			CA-SFM	E ❷	94.6	2.0	0.3	2.5	91.0	2.0	0.3	8.1	
			Global	#, E ❷	94.9	1.2	0.5	3.8	91.9	1.2	0.5	7.4	
Piperacillin/ Tazobactam*	TZP	tzp03n	CLSI	#, E	N/A	N/A	N/A	N/A	93.2	0.9	2.2	5.0	96.7
Tigecycline	TGC	tgc01n	Global	#, E	94.4	0.0	2.0	0.0	90.1	0.0	2.0	8.5	100
Trimethoprim/ Sulfamethoxazole	SXT	sxt02n	CLSI	#, E ❷	–	–	–	–	100	0.0	0.0	N/A	100
			CA-SFM	E ❷	99.5	0.0	0.0	0.0	100	0.0	0.0	0.0	
Trimethoprim/ Sulfamethoxazole	SXT	sxt02n	CLSI[3]	#, E ❷ Applies to V2S SW PC version 5.02 and above	-	-	-	-	100	0.0	0.0	N/A	97.8

[1] Abbreviations — Bp = breakpoint committee; EA = essential agreement; CA = category agreement; VME = Very Major Error (susceptible result with resistant reference result); ME = Major Error (resistant result with susceptible reference result); mE =

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minor Error (susceptible or resistant result with an intermediate reference result, or an intermediate result with a susceptible or resistant reference result).

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* The overall essential agreement for VITEK® 2 piperacillin/tazobactam was not calculated since the test contains less than five discrete concentrations. However, of the 2453 isolates tested, 348 (14.2%) yielded on-scale VITEK® 2 MIC values (i.e., results of 8, 16, 32, and 64 µg/ml). Those results were evaluated for MIC trends and the following results were obtained when VITEK® 2 was compared to CLSI broth microdilution method; 74.4% (259/348) of the VITEK® MICs were within +/- one doubling dilution; 18.4% (64/348) of the VITEK® MICs were higher by more than one doubling dilution; and 7.2% (25/348) of the VITEK® MICs were lower by more than one doubling dilution.

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Table 5: EUCAST Performance Characteristics for Gram-Negative Antimicrobial Susceptibility Testing

Antimicrobial	Antimicrobial Code	Antibiotic Version	Comment[1]	Essential Agreement				Category Agreement			
				% Error				% Error			
				%EA[2]	VME	ME	mE	%CA	VME	ME	mE
Ampicillin	AM	am01n	*Enterobacteriaceae* (excludes: *Citrobacter, Enterobacter, Pantoea, Serratia* spp.)	93.2	2.7	0.9	0.0	96.6	4.3	1.8	0.0
Cefotaxime	CTX	ctx01n	*Enterobacteriaceae*	95.9	7.4	0.8	0.9	94.1	7.4	0.8	3.8
Cefpodoxime	CPD	cpd01n	*Enterobacteriaceae* excluding *Serratia* spp.	93.8	3.8	2.2	0.0	94.1	6.3	5.8	0.0
Ceftazidime	CAZ	caz01n	*Enterobacteriaceae, Pseudomonas*	96.7	3.0	0.3	0.8	97.9	3.0	1.0	1.3
Cefuroxime	CXM	cxm01n	*Enterobacteriaceae*	94.0	2.8	2.6	0.0	93.7	3.7	7.8	0.0
Ciprofloxacin	CIP	cip01n	N/A	98.1	0.0	0.3	0.2	95.3	0.0	0.3	4.5
Ertapenem	ETP	etp01n	*Enterobacteriaceae*	99.5	0.0	0.0	0.5	98.9	0.0	0.0	1.1
Gentamicin	GM	gm01n	N/A	98.3	0.4	0.0	0.8	96.8	0.0	0.0	3.2
Imipenem	IPM	ipm04n	N/A	95.8	1.0	0.4	0.8	92.2	1.0	0.4	7.3
Meropenem	MEM	mem02n	❷	97.6	0.5	0.0	0.6	95.9	0.5	0.0	3.9
Moxifloxacin	MXF	mxf01n	*Enterobacteriaceae*	98.7	0.0	0.6	0.2	95.4	0.0	0.6	4.1
Piperacillin	PIP	pip02n	*Enterobacteriaceae, Pseudomonas* ❷	93.9	2.4	1.6	1.9	88.5	4.3	2.9	8.5
Piperacillin/Tazobactam	TZP	tzp03n	*Enterobacteriaceae, Pseudomonas*	96.2	0.8	0.9	1.0	94.8	1.2	1.5	3.8
Tigecycline	TGC	tgc01n	(excludes: *Serratia, Morganella/Proteus/ Providencia*, and *Enterobacter*)	96.9	0.0	0.4	0.8	94.6	0.0	0.4	5.0
Trimethoprim/ Sulfamethoxazole	SXT	sxt02n	(also includes *S. maltophilia*) ❷	98.4	0.0	0.0	0.4	99.5	0.0	0.0	0.5

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- *Edwardsiella hoshinae**
- *Edwardsiella tarda**
- *Elizabethkingia meningoseptica* (formerly known as *Chryseobacterium meningosepticum*)
- *Enterobacter aerogenes**
- *Enterobacter asburiae**
- *Enterobacter cancerogenus**
- *Enterobacter cloacae**
- *Enterobacter cloacae* complex*
- *Enterobacter cloacae* ssp. *cloacae**
- *Enterobacter cloacae* ssp. *dissolvens* *
- *Escherichia coli**
- *Escherichia coli* ATCC® 25922™
- *Escherichia coli* ATCC® 35218™
- *Escherichia fergusonii*
- *Escherichia hermannii**
- *Escherichia vulneris**
- *Ewingella americana*
- *Hafnia alvei**
- *Hafnia paralvei**
- *Klebsiella oxytoca**
- *Klebsiella pneumoniae**
- *Klebsiella pneumoniae* ssp. *ozaenae**
- *Klebsiella pneumoniae* ssp. *pneumoniae**
- *Klebsiella pneumoniae* ssp. *pneumoniae* ATCC® 700603 ™
- *Klebsiella pneumoniae* ssp. *rhinoscleromatis**
- *Klebsiella* spp.*
- *Kluyvera ascorbata**
- *Kluyvera cryocrescens**
- *Kluyvera intermedia* (formerly known as *Enterobacter intermedius*)*
- *Leclercia adecarboxylata**
- *Lelliottia amnigena* (formerly known as Enterobacter amnigenus)*
- *Lelliottia amnigena 1* (formerly known as Enterobacter amnigenus 1)*
- *Lelliottia amnigena 2* (formerly known as Enterobacter amnigenus 2)
- *Mannheimia haemolytica*
- *Moraxella* group
- *Moraxella lacunata*
- *Moraxella nonliquefaciens*
- *Moraxella osloensis*
- *Morganella morganii**
- *Morganella morganii* ssp. *morganii**
- *Morganella morganii* ssp. *sibonii**
- *Myroides* ssp.
- *Pantoea agglomerans*
- *Pantoea dispersa*
- *Pasteurella aerogenes*
- *Pasteurella multocida**
- *Pasteurella pneumotropica*
- *Plesiomonas shigelloides*
- *Pluralibacter gergoviae* (formerly known as *Enterobacter gergoviae*)*
- *Proteus hauseri**
- *Proteus mirabilis**
- *Proteus penneri**

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REFERENCES

1. Barry, AL The Antimicrobic Susceptibility Test, Principles and Practices, Lea and Febiger, Philadelphia, PA. 1976.
2. Clinical Laboratory Standards Institute (CLSI®), Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically, M7- A7, Wayne, Pennsylvania, January 2006.
3. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Eighteenth Informational Supplement, M100-S18, Vol. 27, No. 1, January 2008.
4. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-third Informational Supplement, M100-S22, January 2012.
5. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fourth Informational Supplement; M100-S24, January 2014.
6. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fifth Informational Supplement, M100-S25, January 2015.
7. Comité de l'Antibiogramme de la Société Française de Microbiologie. Communiqué 1996. Path Biol, 1996, 44, n° 8, I-VIII.
8. Comité de l'Antibiogramme de la Société Française de Microbiologie, Communiqué 2007.
9. Comite de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM), Recommendations 2012.
10. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2014.
11. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2015.
12. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 2.0, January 2012.
13. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 4.0, January 2014.
14. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 5.0, January 2015.
15. Gerlach, EH Microdilution 1: A Comparative Study, p. 63-76, In: Balows, A. (ed.), Current Techniques for Antibiotic Susceptibility Testing, Charles C. Thomas, Springfield, IL. 1974.
16. MacLowry, JD, and HH Marsh. 1968. Semi-automatic microtechnique for serial dilution antibiotic sensitivity testing in the clinical laboratory. J. Lab. Clin. Med. 1968;72:685-687.
17. Murray, PR, Baron EJ, Pfaller MA, Tenover FC, and Yolken RH, editors. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue – Approved Guideline (1997).
19. National Committee for Clinical Laboratory Standards, Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard — Third Edition, M27-A3, Vol. 22, No. 15, 2008.
20. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.

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**MSSA = Methicillin-susceptible *S. aureus*

PRECAUTIONS

- For *In Vitro* Diagnostic Use Only.
- For US Only: Caution: US Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- For professional use only.
- Suspensions not within the appropriate zone on the VITEK® 2 DensiCHEK™ Plus or VITEK® 2 DensiCHEK™ may compromise card performance.
- Do not use the card after the expiration date shown on the package liner.
- Store the card unopened in the package liner. Do not use the card if the protective package liner is damaged or if no desiccant is present.
- Allow the card to come to room temperature before opening the package liner.
- Do not use powdered gloves. Powder may interfere with the optics.
- Use of culture media other than the recommended types must be validated by the customer laboratory for acceptable performance.
- The card performs as intended only when used in conjunction with VITEK® 2 Systems, following the instructions contained in the Instructions for Use.
- It is highly recommended that Good Laboratory Practices (e.g., FDA, CLSI, ISO, etc.) also be followed, according to local guidelines or requirements.
- **Do not use glass test tubes**. Use clear plastic (polystyrene) tubes only. Variation exists among test tubes of standard diameter. Carefully place the tube into the cassette. If resistance is encountered, discard and try another tube that does not require pressure to insert.
- Prior to inoculation, inspect cards for tape tears or damage to the tape and discard any that are suspect. Check the saline level in the tubes after the cassette has been processed to ensure proper filling of card.
 - VITEK® 2 60 or VITEK® 2 XL: Eject improperly filled cards.
 - VITEK® 2 Compact: Do not load improperly filled cards.
- Give special consideration to specimen source and patient therapy regimen. AST cards may contain some antimicrobials that are not proven to be effective for treatment of infections due to all organisms that may be tested. For interpreting and reporting of antimicrobial results that have been shown to be active against organism groups both *in vitro* and in clinical infections, refer to the individual pharmaceutical antimicrobial labeling or local therapy guidelines.

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- Interpretation of test results requires the judgment and skill of a person knowledgeable in AST. Additional testing may be required.[17]

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- Sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- Loops, sterile sticks or swabs
- Appropriate agar medium (See the Culture Requirements Table.)
- QC isolates

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REFERENCES

1. Barry, AL The Antimicrobic Susceptibility Test, Principles and Practices, Lea and Febiger, Philadelphia, PA. 1976.
2. Clinical Laboratory Standards Institute (CLSI®), Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically, M7- A7, Wayne, Pennsylvania, January 2006.
3. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Eighteenth Informational Supplement, M100-S18, Vol. 27, No. 1, January 2008.
4. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-third Informational Supplement, M100-S22, January 2012.

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12. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 2.0, January 2012.

13. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 4.0, January 2014.

14. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 5.0, January 2015.

15. Gerlach, EH Microdilution 1: A Comparative Study, p. 63-76, In: Balows, A. (ed.), Current Techniques for Antibiotic Susceptibility Testing, Charles C. Thomas, Springfield, IL. 1974.

16. MacLowry, JD, and HH Marsh. 1968. Semi-automatic microtechnique for serial dilution antibiotic sensitivity testing in the clinical laboratory. J. Lab. Clin. Med. 1968;72:685-687.

17. Murray, PR, Baron EJ, Pfaller MA, Tenover FC, and Yolken RH, editors. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C. 2003.

18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue – Approved Guideline (1997).

19. National Committee for Clinical Laboratory Standards, Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard — Third Edition, M27-A3, Vol. 22, No. 15, 2008.

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08 

INDEX OF SYMBOLS

Symbol	Meaning
REF	Catalog number
IVD	In Vitro Diagnostic Medical Device
	Legal Manufacturer
	Temperature limitation
	Use by date
LOT	Batch code
	Consult Instructions for Use
	Date of manufacture
Σ	Contains sufficient for <n> tests

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Symbol	Meaning
EC REP	Authorized representative in the European Community
Rx only	For US Only : Caution : US Federal Law restricts this device to sale by or on the order of a licensed practitioner

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Release Date	Part Number	Change Type	Change Summary
2018-01	045616-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• Combined product package insert content with VITEK® 2 Product Information Manual AST content • Updated Limited Warranty section • Updated with RX only information • Updated with EUCAST limitations

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CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.
Any other name or trademark is the property of its respective owner.
©BIOMÉRIEUX 2018

bioMérieux, Inc.
100 Rodolphe Street
Durham, North Carolina 27712 USA
www.biomerieux.com

EC REP **bioMérieux SA**
376 Chemin de l'Orme
69280 Marcy-l'Etoile - France
673 620 399 RCS LYON
Tel. 33 (0)4 78 87 20 00
Fax 33 (0)4 78 87 20 90

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[c] = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.

**Kox(-ESBL) = *K. oxytoca* (excluding ESBL producing strains)

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**Eco(+ETEC) = E. coli (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

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1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.
2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).
3. Using sterile technique, prepare a homogenous organism suspension with a density equivalent to the appropriate McFarland standard or McFarland Reference using a compatible benchtop densitometer (see the Culture Requirements Table).
 Note: The age of the suspension before loading the instrument for AST testing must be less than one hour when using VITEK® 2 60 or VITEK® 2 XL, and less than 30 minutes when using VITEK® 2 Compact.
4. Choose one of the following:
 - **For an automatic dilution (VITEK® 2 60 or VITEK® 2 XL only):** Place the suspension tube prepared in step 3 into the cassette with or without an identification card. In the next cassette slot, place an empty tube and an AST card. The instrument will automatically dilute the bacterial suspension to prepare an inoculum suitable for the susceptibility card.
 - **For a manual dilution (VITEK® 2 Compact, VITEK® 2 60 or VITEK® 2 XL):** In a second tube containing 3.0 mL of saline, transfer 145 µL of the suspension prepared in step 3. Then place this tube in the cassette with a susceptibility card. The tube with the initial bacterial suspension can also be used for inoculation of an identification card.

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- *Klebsiella pneumoniae* ssp. *pneumoniae* ATCC® 700603 ™
- *Klebsiella pneumoniae* ssp. *rhinoscleromatis**
- *Klebsiella* spp.*
- *Kluyvera ascorbata**
- *Kluyvera cryocrescens**
- *Kluyvera intermedia* (formerly known as *Enterobacter intermedius*)*
- *Leclercia adecarboxylata**
- *Lelliottia amnigena* (formerly known as Enterobacter amnigenus)*
- *Lelliottia amnigena 1* (formerly known as Enterobacter amnigenus 1)*
- *Lelliottia amnigena 2* (formerly known as Enterobacter amnigenus 2)
- *Mannheimia haemolytica*
- *Moraxella* group
- *Moraxella lacunata*
- *Moraxella nonliquefaciens*
- *Moraxella osloensis*
- *Morganella morganii**
- *Morganella morganii* ssp. *morganii**
- *Morganella morganii* ssp. *sibonii**
- *Myroides* ssp.
- *Pantoea agglomerans*
- *Pantoea dispersa*
- *Pasteurella aerogenes*
- *Pasteurella multocida**
- *Pasteurella pneumotropica*
- *Plesiomonas shigelloides*
- *Pluralibacter gergoviae* (formerly known as *Enterobacter gergoviae*)*
- *Proteus hauseri**
- *Proteus mirabilis**
- *Proteus penneri**
- *Proteus vulgaris**
- *Providencia alcalifaciens**
- *Providencia rettgeri**
- *Providencia rustigianii**
- *Providencia stuartii**
- *Pseudomonas aeruginosa**
- *Pseudomonas aeruginosa* ATCC® 27853™
- *Pseudomonas alcaligenes*
- *Pseudomonas fluorescens*
- *Pseudomonas luteola*
- *Pseudomonas mendocina*
- *Pseudomonas oleovorans*
- *Pseudomonas oryzihabitans*
- *Pseudomonas putida*
- *Pseudomonas* spp.
- *Pseudomonas stutzeri*
- *Ralstonia pickettii*
- *Raoultella ornithinolytica**
- *Raoultella planticola**
- *Raoultella terrigena**
- *Salmonella enterica* ssp. *arizonae*
- *Salmonella enterica* ssp. *enterica**
- *Salmonella* group*

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5. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fourth Informational Supplement; M100-S24, January 2014.
6. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fifth Informational Supplement, M100-S25, January 2015.
7. Comité de l'Antibiogramme de la Société Française de Microbiologie. Communiqué 1996. Path Biol, 1996, 44, n° 8, I-VIII.
8. Comité de l'Antibiogramme de la Société Française de Microbiologie, Communiqué 2007.
9. Comite de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM), Recommendations 2012.
10. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2014.
11. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2015.
12. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 2.0, January 2012.
13. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 4.0, January 2014.
14. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 5.0, January 2015.
15. Gerlach, EH Microdilution 1: A Comparative Study, p. 63-76, In: Balows, A. (ed.), Current Techniques for Antibiotic Susceptibility Testing, Charles C. Thomas, Springfield, IL. 1974.
16. MacLowry, JD, and HH Marsh. 1968. Semi-automatic microtechnique for serial dilution antibiotic sensitivity testing in the clinical laboratory. J. Lab. Clin. Med. 1968;72:685-687.
17. Murray, PR, Baron EJ, Pfaller MA, Tenover FC, and Yolken RH, editors. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue – Approved Guideline (1997).
19. National Committee for Clinical Laboratory Standards, Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard — Third Edition, M27-A3, Vol. 22, No. 15, 2008.
20. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.

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INDEX OF SYMBOLS

Symbol	Meaning
	Catalog number
	In Vitro Diagnostic Medical Device
	Legal Manufacturer
	Temperature limitation
	Use by date
	Batch code
	Consult Instructions for Use
	Date of manufacture
	Contains sufficient for <n> tests

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- *Klebsiella pneumoniae* ssp. *pneumoniae**
- *Klebsiella pneumoniae* ssp. *pneumoniae* ATCC® 700603 ™
- *Klebsiella pneumoniae* ssp. *rhinoscleromatis**
- *Klebsiella* spp.*
- *Kluyvera ascorbata**
- *Kluyvera cryocrescens**
- *Kluyvera intermedia* (formerly known as *Enterobacter intermedius*)*
- *Leclercia adecarboxylata**
- *Lelliottia amnigena* (formerly known as Enterobacter amnigenus)*
- *Lelliottia amnigena 1* (formerly known as Enterobacter amnigenus 1)*
- *Lelliottia amnigena 2* (formerly known as Enterobacter amnigenus 2)
- *Mannheimia haemolytica*
- *Moraxella* group
- *Moraxella lacunata*
- *Moraxella nonliquefaciens*
- *Moraxella osloensis*
- *Morganella morganii**
- *Morganella morganii* ssp. *morganii**
- *Morganella morganii* ssp. *sibonii**
- *Myroides* ssp.
- *Pantoea agglomerans*
- *Pantoea dispersa*
- *Pasteurella aerogenes*
- *Pasteurella multocida**
- *Pasteurella pneumotropica*
- *Plesiomonas shigelloides*
- *Pluralibacter gergoviae* (formerly known as *Enterobacter gergoviae*)*
- *Proteus hauseri**
- *Proteus mirabilis**
- *Proteus penneri**
- *Proteus vulgaris**
- *Providencia alcalifaciens**
- *Providencia rettgeri**
- *Providencia rustigianii**
- *Providencia stuartii**
- *Pseudomonas aeruginosa**
- *Pseudomonas aeruginosa* ATCC® 27853™
- *Pseudomonas alcaligenes*
- *Pseudomonas fluorescens*
- *Pseudomonas luteola*
- *Pseudomonas mendocina*
- *Pseudomonas oleovorans*
- *Pseudomonas oryzihabitans*
- *Pseudomonas putida*
- *Pseudomonas* spp.
- *Pseudomonas stutzeri*
- *Ralstonia pickettii*
- *Raoultella ornithinolytica**
- *Raoultella planticola**
- *Raoultella terrigena**
- *Salmonella enterica* ssp. *arizonae*
- *Salmonella enterica* ssp. *enterica**

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4. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-third Informational Supplement, M100-S22, January 2012.
5. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fourth Informational Supplement; M100-S24, January 2014.
6. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fifth Informational Supplement, M100-S25, January 2015.
7. Comité de l'Antibiogramme de la Société Française de Microbiologie. Communiqué 1996. Path Biol, 1996, 44, n° 8, I-VIII.
8. Comité de l'Antibiogramme de la Société Française de Microbiologie, Communiqué 2007.
9. Comite de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM), Recommendations 2012.
10. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2014.
11. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2015.
12. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 2.0, January 2012.
13. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 4.0, January 2014.
14. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 5.0, January 2015.
15. Gerlach, EH Microdilution 1: A Comparative Study, p. 63-76, In: Balows, A. (ed.), Current Techniques for Antibiotic Susceptibility Testing, Charles C. Thomas, Springfield, IL. 1974.
16. MacLowry, JD, and HH Marsh. 1968. Semi-automatic microtechnique for serial dilution antibiotic sensitivity testing in the clinical laboratory. J. Lab. Clin. Med. 1968;72:685-687.
17. Murray, PR, Baron EJ, Pfaller MA, Tenover FC, and Yolken RH, editors. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue – Approved Guideline (1997).
19. National Committee for Clinical Laboratory Standards, Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard — Third Edition, M27-A3, Vol. 22, No. 15, 2008.
20. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.

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INDEX OF SYMBOLS

Symbol	Meaning
	Catalog number
	In Vitro Diagnostic Medical Device
	Legal Manufacturer
	Temperature limitation
	Use by date
	Batch code
	Consult Instructions for Use
	Date of manufacture

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Symbol	Meaning
	Contains sufficient for <n> tests
	Authorized representative in the European Community
	For US Only : Caution : US Federal Law restricts this device to sale by or on the order of a licensed practitioner

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Release Date	Part Number	Change Type	Change Summary
2018-01	046017-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• Combined product package insert content with VITEK® 2 Product Information Manual AST content • Updated Limited Warranty section • Updated with RX only information • Updated with EUCAST limitations

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**CSAGNB = Clinically significant aerobic gram-negative bacilli

Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid.

NEG = Negative

POS = Positive

A negative ESBL test result does not rule out the presence of an ESBL masked by an AmpC beta-lactamase.

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c = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.

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SPECIMEN PREPARATION

Culture Requirements Table

VITEK® 2 Card	Media	Age of Culture	Incubation Conditions	McFarland Standards	Dilution for AST	Age of Suspension Before Loading Instrument
AST Gram-Negative	TSAB CBA MAC CPS ID	8 to 24 hours	35°C to 37°C aerobic, non-CO_2	0.50 to 0.63	145 µL in 3.0 mL 0.45% saline	VITEK® 2 Compact: ≤ 30 minutes VITEK® 2: ≤ 1 hour
GN and AST GN pair	CBA[1] MAC[1] TSAB CPS ID	18 to 24 hours	35°C to 37°C aerobic, non-CO_2	0.50 to 0.63	145 µL in 3.0 mL 0.45% saline	≤ 30 minutes

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- Sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- Loops, sterile sticks or swabs
- Appropriate agar medium (See the Culture Requirements Table.)
- QC isolates
- VITEK® 2 AST Cards
- Micropipettors to deliver 145 µL

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1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.
2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).
3. Use a sterile stick or swab to transfer a sufficient number of morphologically similar colonies to the saline tube prepared in step 2. Prepare a homogenous organism suspension with a density equivalent to the appropriate McFarland standard using the VITEK® 2 DensiCHEK™ Plus (see the Culture Requirements Table). The inoculum density required for the GN, GP, ST, or YST can be diluted automatically for AST cards by the VITEK® 2 instrument (VITEK® 2 60 or VITEK® 2 XL).
 Note: The age of the suspension before loading the instrument for AST testing must be less than one hour when using VITEK® 2 60 or VITEK® 2 XL, and less than 30 minutes when using VITEK® 2 Compact.
4. Choose one of the following:
 - **For an automatic dilution (VITEK® 2 60 or VITEK® 2 XL only):** Place the suspension tube prepared in step 3 into the cassette with or without an identification card. In the next cassette slot, place an empty tube and an AST card. The instrument will automatically dilute the bacterial suspension to prepare an inoculum suitable for the susceptibility card.
 - **For a manual dilution (VITEK® 2 Compact, VITEK® 2 60 or VITEK® 2 XL):** In a second tube containing 3.0 mL of saline, transfer 145 µL of the suspension prepared in step 3. Then place this tube in the cassette with a susceptibility card. The tube with the initial bacterial suspension can also be used for inoculation of an identification card.

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Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid.

NEG = Negative

POS = Positive

A negative ESBL test result does not rule out the presence of an ESBL masked by an AmpC beta-lactamase.

† *Escherichia coli* ATCC® 35218™ is not intended for routine QC testing of ampicillin, piperacillin, or ticarcillin. This organism is recommended only for routine QC of ß-lactamase inhibitor combinations. However, since this strain contains a plasmid-encoded ß-lactamase (non-ESBL), it is resistant to many penicillinase-labile antimicrobials, but susceptible to ß-lactamase/ß-lactamase inhibitor combinations. To ensure that the plasmid is present, testing of the ß-lactam alone (e.g., AM, PIP, TIC), in addition to the ß-lactamase/ß-lactamase inhibitor combination (e.g., AMC, SAM, TZP, TCC), can be performed. If the strain has lost its plasmid, it will be susceptible to the ß-lactam alone (i.e., AM, PIP, TIC), indicating that the QC test is invalid, and a new culture of *Escherichia coli* ATCC® 35218™ must be used.

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Preparation of QC Organisms

1. Rehydrate the organism according to the manufacturer's instructions.

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2. Subculture to Trypticase Soy Agar with 5% sheep blood (TSAB).
3. Incubate at 35°C for 24 hours.
4. Check for purity.
5. Subculture to a TSAB plate.
6. Incubate 16 to 18 hours at 35°C.

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Ticarcillin/Clavulanic Acid (tcc01n) has been approved for use only with *Pseudomonas aeruginosa*. An alternative method of testing must be performed prior to reporting of results for all other organisms.

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Perform an alternative method of testing prior to reporting results when a resistant result is obtained with the following antibiotic/organism combination(s):

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Perform an alternative method of testing prior to reporting of susceptible results for the following antibiotic/organism combination(s):

- Colistin (cs01n): Clinically-significant, aerobic, Gram-negative bacilli

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Perform an alternative method of testing prior to reporting of susceptible results for the following antibiotic/organism combination(s):

- Colistin (cs01n): Clinically-significant, aerobic, Gram-negative bacilli

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Category agreement (CA) occurs when the VITEK® 2 and the reference interpretative results agree (Susceptible, Intermediate, and Resistant). There are instances when the category agreement for an antimicrobial falls below the essential agreement. This can occur when a significant number of MICs cluster around a category breakpoint during clinical trial

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testing, resulting in interpretative errors. For a description of interpretive errors, refer to the footnotes below the table that follows (Performance Characteristics). When the majority of the errors are of the minor type, a high corresponding essential agreement percentage demonstrates that the antimicrobial retains an acceptable overall performance.

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*Data on file at bioMérieux, Inc

Performance Characteristics for Gram-Negative Antimicrobial Susceptibility Testing

Antimicrobial	Antimicrobial Code	Antibiotic Version	Bp[1]	Comment[2]	Essential Agreement				Category Agreement				% Reproducibility
					% Error				% Error				
					% EA	VME	ME	mE	% CA	VME	ME	mE	
Aztreonam	ATM	atm01n	CLSI	#, E	98.9	1.6	0.3	0.3	96.6	1.6	0.3	0.3	100
			CA-SFM	I	93.7	0.8	0.2	5.2	–	–	–	–	
Cefepime	FEP	fep02n	CLSI	#, E	94.2	4.0	1.4	6.3	92.1	–	–	–	98.6
Cefixime	CFM	cfm01n	CLSI	I	96.3	0.3	1.5	1.1	92.3	0.4	1.5	6.8	99.8
			CA-SFM	I	96.3	0.3	1.5	1.1	92.3	0.4	1.5	6.8	
Ceftriaxone	CRO	cro01n	CLSI	#, E	95.8	1.1	0.5	3.0	90.2	1.1	0.5	9.3	98.5
			CA-SFM	I	93.0	0.4	0.6	4.2	–	–	–	–	
Cefuroxime	CXM	cxm01n	CLSI	#, E	95.3	1.5	1.0	1.6	93.9	1.5	1.0	4.8	99.6
Chloramphenicol	C	c01n	CLSI	I	95.3	0.8	0.0	3.0	86.2	0.8	0.0	13.4	96.1
			CA-SFM	I	95.3	0.8	0.0	3.0	86.2	0.8	0.0	13.4	
Colistin	CS	cs01n	CA-SFM	I	96.1	3.4	0.7	0.0	97.3	5.4	1.4	N/A	99.2
Extended-Spectrum Beta-Lactamase (ESBL)	ESBL	esb01n	CLSI	#, E	N/A	N/A	N/A	N/A	97.7	2.7	1.9	N/A	100
Levofloxacin	LEV	lev02n	CLSI	#, E ❷	98.6	0.5	0.2	0.3	95.8	0.5	0.2	3.9	100
			CA-SFM	E ❷	98.6	0.0	0.7	0.1	96.5	0.0	0.7	3.0	
Meropenem	MEM	mem02n	CLSI	#, E ❷	97.6	1.4	0.0	0.3	96.8	1.4	0.0	2.8	99.6
			Global	#, E ❷	97.6	0.0	0.0	0.5	97.7	0.0	0.0	2.3	
Minocycline	MNO	mno01n	CLSI	I	94.7	0.3	1.3	2.3	84.8	0.3	1.3	14.3	98
			CA-SFM	I	94.7	0.3	1.3	2.3	84.8	0.3	1.3	14.3	
Moxifloxacin	MXF	mxf01n	CLSI	#, E	98.7	0.0	0.0	0.0	97.8	2.2	0.0	0.0	100
Piperacillin	PIP	pip02n	CLSI	#, E ❷	94.6	2.0	0.2	3.7	92.4	2.7	0.2	6.8	99.3
			CA-SFM	E ❷	94.6	2.0	0.3	2.5	91.0	2.0	0.3	8.1	
			Global	#, E ❷	94.9	1.2	0.5	3.8	91.9	1.2	0.5	7.4	
Tetracycline	TE	te01n	CLSI	#, E	95.5	1.0	0.0	1.4	93.0	1.0	0.0	6.6	100
Ticarcillin/Clavulanic Acid	TCC	tcc01n	CLSI	#, E	98.9	0.0	1.9	0.6	92.7	0.0	1.9	19.6	100
			CA-SFM	I	94.8	2.2	0.3	4.1	–	–	–	–	
Tigecycline	TGC	tgc01n	Global	#, E	94.4	0.0	2.0	0.0	90.1	0.0	2.0	8.5	100

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[2] Comment — Specific organism groups are designated as P. aeurg. for *Pseudomonas aeruginosa*, others for species other than *Pseudomonas aeruginosa*, and Acineto. for *Acinetobacter*.

Key:

= US Food and Drug Administration 510(k) cleared
CLSI® = Clinical and Laboratory Standards Institute
CA-SFM = Comité de l'Antibiogramme de la Société Française de Microbiologie
E = External performance data
I = Internal performance data
– = Not available
❶ ❷ = Symbol identifies the performance characteristics for a specific antimicrobial version.
Ref. = Reference method for clinical performance study.

EUCAST Performance Characteristics for Gram-Negative Antimicrobial Susceptibility Testing

Antimicrobial	Antimicrobial Code	Antibiotic Version	Comment[1]	Essential Agreement				Category Agreement			
				% Error				% Error			
				%EA[2]	VME	ME	mE	%CA	VME	ME	mE
Aztreonam	ATM	atm01n	*Enterobacteriaceae, Pseudomonas*	98.1	0.0	0.9	1.0	97.4	0.0	0.9	1.9
Cefepime	FEP	fep02n	*Enterobacteriaceae, Pseudomonas*	97.0	2.0	1.3	0.8	93.9	9.0	2.6	2.4
Cefixime	CFM	cfm01n	*Enterobacteriaceae*	95.6	2.5	1.6	0.0	96.4	8.6	2.3	0.0
Ceftriaxone	CRO	cro01n	*Enterobacteriaceae*	98.0	0.0	0.4	1.3	98.0	0.0	0.4	1.6
Cefuroxime	CXM	cxm01n	*Enterobacteriaceae*	94.0	2.8	2.6	0.0	93.7	3.7	7.8	0.0
Chloramphenicol	C	c01n	*Enterobacteriaceae* (excludes: *Serratia*)	97.2	2.3	3.1	0.0	89.1	4.1	14.8	0.0
Colistin	CS	cs01n	*Pseudomonas, Acinetobacter, Enterobacteriaceae*	96.3	4.7	0.6	0.0	97.8	5.0	0.7	0.0
Levofloxacin	LEV	lev02n	❷	99.4	0.0	1.3	0.2	96.7	0.0	0.7	2.9
Meropenem	MEM	mem02n	❷	97.6	0.5	0.0	0.6	95.9	0.5	0.0	3.9
Moxifloxacin	MXF	mxf01n	*Enterobacteriaceae*	98.7	0.0	0.6	0.2	95.4	0.0	0.6	4.1
Piperacillin	PIP	pip02n	*Enterobacteriaceae, Pseudomonas* ❷	93.9	2.4	1.6	1.9	88.5	4.3	2.9	8.5
Ticarcillin/Clavulanic Acid	TCC	tcc01n	*Enterobacteriaceae* (excludes: *Serratia*), *Pseudomonas*	94.7	2.9	0.8	1.7	87.0	6.9	1.9	8.9
Tigecycline	TGC	tgc01n	(excludes: *Serratia, Morganella/Proteus/Providencia*, and *Enterobacter*)	96.9	0.0	0.4	0.8	94.6	0.0	0.4	5.0
Trimethoprim	TMP	tmp01n	*Enterobacteriaceae*	94.2	3.8	0.5	1.3	93.5	3.8	0.5	5.0

[1] Comment — Unless stated otherwise, performance is for *Enterobacteriaceae*, *Pseudomonas*, and *Acinetobacter*.

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- *Edwardsiella tarda**
- *Elizabethkingia meningoseptica* (formerly known as *Chryseobacterium meningosepticum*)
- *Enterobacter aerogenes**
- *Enterobacter asburiae**
- *Enterobacter cancerogenus**
- *Enterobacter cloacae**
- *Enterobacter cloacae* complex*
- *Enterobacter cloacae* ssp. *cloacae**
- *Enterobacter cloacae* ssp. *dissolvens* *
- *Escherichia coli**
- *Escherichia coli* ATCC® 25922™
- *Escherichia coli* ATCC® 35218™
- *Escherichia fergusonii*
- *Escherichia hermannii**
- *Escherichia vulneris**
- *Ewingella americana*
- *Hafnia alvei**
- *Hafnia paralvei**
- *Klebsiella oxytoca**
- *Klebsiella pneumoniae**
- *Klebsiella pneumoniae* ssp. *ozaenae**
- *Klebsiella pneumoniae* ssp. *pneumoniae**
- *Klebsiella pneumoniae* ssp. *pneumoniae* ATCC® 700603 ™
- *Klebsiella pneumoniae* ssp. *rhinoscleromatis**
- *Klebsiella* spp.*
- *Kluyvera ascorbata**
- *Kluyvera cryocrescens**
- *Kluyvera intermedia* (formerly known as *Enterobacter intermedius*)*
- *Leclercia adecarboxylata**
- *Lelliottia amnigena* (formerly known as Enterobacter amnigenus)*
- *Lelliottia amnigena 1* (formerly known as Enterobacter amnigenus 1)*
- *Lelliottia amnigena 2* (formerly known as Enterobacter amnigenus 2)
- *Mannheimia haemolytica*
- *Moraxella* group
- *Moraxella lacunata*
- *Moraxella nonliquefaciens*
- *Moraxella osloensis*
- *Morganella morganii**
- *Morganella morganii* ssp. *morganii**
- *Morganella morganii* ssp. *sibonii**
- *Myroides* ssp.
- *Pantoea agglomerans*
- *Pantoea dispersa*
- *Pasteurella aerogenes*
- *Pasteurella multocida**
- *Pasteurella pneumotropica*
- *Plesiomonas shigelloides*
- *Pluralibacter gergoviae* (formerly known as *Enterobacter gergoviae*)*
- *Proteus hauseri**
- *Proteus mirabilis**
- *Proteus penneri**
- *Proteus vulgaris**

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- Sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0)
- Loops, sterile sticks or swabs
- Appropriate agar medium (See the Culture Requirements Table.)
- QC isolates
- VITEK® 2 AST Cards
- Micropipettors to deliver 145 µL
- Micropipettors to deliver 280 µL
- Disposable pipette tips

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1. Do one of the following:
 - Select isolated colonies from a primary plate if culture requirements are met.
 - Subculture the organism to be tested to appropriate agar medium and incubate accordingly.
2. Aseptically transfer 3.0 mL of sterile saline (aqueous 0.45% to 0.50% NaCl, pH 4.5 to 7.0) into a clear plastic (polystyrene) test tube (12 mm x 75 mm).
3. Using sterile technique, prepare a homogenous organism suspension with a density equivalent to the appropriate McFarland standard or McFarland Reference using a compatible benchtop densitometer (see the Culture Requirements Table).
 Note: The age of the suspension before loading the instrument for AST testing must be less than one hour when using VITEK® 2 60 or VITEK® 2 XL, and less than 30 minutes when using VITEK® 2 Compact.
4. Choose one of the following:
 - **For an automatic dilution (VITEK® 2 60 or VITEK® 2 XL only):** Place the suspension tube prepared in step 3 into the cassette with or without an identification card. In the next cassette slot, place an empty tube and an AST card.

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The instrument will automatically dilute the bacterial suspension to prepare an inoculum suitable for the susceptibility card.

- **For a manual dilution (VITEK® 2 Compact, VITEK® 2 60 or VITEK® 2 XL):** In a second tube containing 3.0 mL of saline, transfer 145 µL of the suspension prepared in step 3. Then place this tube in the cassette with a susceptibility card. The tube with the initial bacterial suspension can also be used for inoculation of an identification card.

Note: Check the saline level in the tubes after filling. When it is evident by the saline level in the tube that a card has been improperly filled, do not load the card if using VITEK® 2 Compact; **or**, eject the card if using VITEK® 2 60 or VITEK® 2 XL.

Note: Refer to the appropriate Instrument User Manual for detailed information regarding data entry, processing, etc.

5. Follow your local inspecting agency's guidelines for disposal of hazardous waste.

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Gram-Negative Organisms Claimed for AST-GN (keyID)

- *Achromobacter denitrificans*
- *Achromobacter xylosoxidans*
- *Acinetobacter baumannii**
- *Acinetobacter baumannii complex**
- *Acinetobacter calcoaceticus**

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- *Klebsiella oxytoca**
- *Klebsiella pneumoniae**
- *Klebsiella pneumoniae* ssp. *ozaenae**
- *Klebsiella pneumoniae* ssp. *pneumoniae**
- *Klebsiella pneumoniae* ssp. *pneumoniae* ATCC® 700603 ™
- *Klebsiella pneumoniae* ssp. *rhinoscleromatis**
- *Klebsiella* spp.*
- *Kluyvera ascorbata**
- *Kluyvera cryocrescens**
- *Kluyvera intermedia* (formerly known as *Enterobacter intermedius*)*
- *Leclercia adecarboxylata**
- *Lelliottia amnigena* (formerly known as Enterobacter amnigenus)*
- *Lelliottia amnigena 1* (formerly known as Enterobacter amnigenus 1)*
- *Lelliottia amnigena 2* (formerly known as Enterobacter amnigenus 2)
- *Mannheimia haemolytica*
- *Moraxella* group
- *Moraxella lacunata*
- *Moraxella nonliquefaciens*
- *Moraxella osloensis*
- *Morganella morganii**
- *Morganella morganii* ssp. *morganii**
- *Morganella morganii* ssp. *sibonii**
- *Myroides* ssp.
- *Pantoea agglomerans*
- *Pantoea dispersa*
- *Pasteurella aerogenes*
- *Pasteurella multocida**
- *Pasteurella pneumotropica*
- *Plesiomonas shigelloides*
- *Pluralibacter gergoviae* (formerly known as *Enterobacter gergoviae*)*
- *Proteus hauseri**
- *Proteus mirabilis**
- *Proteus penneri**
- *Proteus vulgaris**
- *Providencia alcalifaciens**
- *Providencia rettgeri**
- *Providencia rustigianii**
- *Providencia stuartii**
- *Pseudomonas aeruginosa**
- *Pseudomonas aeruginosa* ATCC® 27853™
- *Pseudomonas alcaligenes*
- *Pseudomonas fluorescens*
- *Pseudomonas luteola*
- *Pseudomonas mendocina*
- *Pseudomonas oleovorans*
- *Pseudomonas oryzihabitans*
- *Pseudomonas putida*
- *Pseudomonas* spp.
- *Pseudomonas stutzeri*
- *Ralstonia pickettii*
- *Raoultella ornithinolytica**
- *Raoultella planticola**

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REFERENCES

1. Barry, AL The Antimicrobic Susceptibility Test, Principles and Practices, Lea and Febiger, Philadelphia, PA. 1976.
2. Clinical Laboratory Standards Institute (CLSI®), Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically, M7- A7, Wayne, Pennsylvania, January 2006.

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3. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Eighteenth Informational Supplement, M100-S18, Vol. 27, No. 1, January 2008.
4. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-third Informational Supplement, M100-S22, January 2012.
5. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fourth Informational Supplement; M100-S24, January 2014.
6. Clinical and Laboratory Standards Institute (CLSI®), Performance Standards for Antimicrobial Susceptibility Testing; Twenty-fifth Informational Supplement, M100-S25, January 2015.
7. Comité de l'Antibiogramme de la Société Française de Microbiologie. Communiqué 1996. Path Biol, 1996, 44, n° 8, I-VIII.
8. Comité de l'Antibiogramme de la Société Française de Microbiologie, Communiqué 2007.
9. Comite de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM), Recommendations 2012.
10. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2014.
11. Comité de l'Antibiogramme de la Société Française de Microbiologie (CA-SFM). Communiqué 2015.
12. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 2.0, January 2012.
13. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 4.0, January 2014.
14. European Committee on Antimicrobial Susceptibility Testing (EUCAST), version 5.0, January 2015.
15. Gerlach, EH Microdilution 1: A Comparative Study, p. 63-76, In: Balows, A. (ed.), Current Techniques for Antibiotic Susceptibility Testing, Charles C. Thomas, Springfield, IL. 1974.
16. MacLowry, JD, and HH Marsh. 1968. Semi-automatic microtechnique for serial dilution antibiotic sensitivity testing in the clinical laboratory. J. Lab. Clin. Med. 1968;72:685-687.
17. Murray, PR, Baron EJ, Pfaller MA, Tenover FC, and Yolken RH, editors. Manual of Clinical Microbiology, 8th ed. American Society for Microbiology, Washington, D.C. 2003.
18. National Committee for Clinical Laboratory Standards, M29-A, Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissue – Approved Guideline (1997).
19. National Committee for Clinical Laboratory Standards, Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; Approved Standard — Third Edition, M27-A3, Vol. 22, No. 15, 2008.
20. U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health, Office of Health and Safety, Biosafety in Microbiological and Biomedical Laboratories, 1988.

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Release Date	Part Number	Change Type	Change Summary
2019-01	046128-01	Administrative	• Formatting changes do not affect the fit, form, or function of the product
		Technical change	• Combined product package insert content with VITEK® 2 Product Information Manual AST content • Updated Limited Warranty section • Updated with RX only information • Updated with EUCAST limitations

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli
**N/A = No specific FDA Indications for Use available
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

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LIMITATIONS

Ticarcillin/Clavulanic Acid has been approved for use only with *Pseudomonas aeruginosa*. An alternative method of testing must be performed prior to reporting of results for all other organisms.

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):
- Amikacin: *Acinetobacter baumannii*
- Aztreonam: *Pseudomonas* spp.
- Imipenem: *Serratia marcescens*
- Piperacillin/Tazobactam: *Serratia marcescens*

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- Colistin: Clinically significant aerobic gram negative bacilli

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
NEG = Negative
POS = Positive
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**MSSA = Methicillin-susceptible *S. aureus*
**VREF = Vancomycin-resistant *E. faecium*
[NS] = The current absence of resistant isolates precludes defining any results other than susceptible. Isolates yielding MIC results suggestive of Nonsusceptible category should be submitted to a reference laboratory for further testing.
[c] = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.

QUALITY CONTROL

		CLSI® Quality Control Organisms VITEK 2 Results					
Antimicrobic	Code	*E. faecalis* ATCC® 29212™	*S. aureus* ATCC® 29213™	*E. faecalis* ATCC® 51299™	*S. aureus* ATCC® BAA-1026™	*S. aureus* ATCC® BAA-976™	*S. aureus* ATCC® BAA-977™
Ampicillin *Enterococcus* spp. *S. agalactiae*	AM	≤2	0.5 – 2	-	-	-	-
Benzylpenicillin *Enterococcus* spp. *Staphylococcus* spp. *S. agalactiae*	P	1 – 4	0.25 – ≥0.5	-	-	-	-
Cefoxitin Screen	OXSF	-	NEG	-	POS	-	-
Ciprofloxacin	CIP	≤0.5 – 2	≤0.5	-	-	-	-
Clindamycin	CM	4 – ≥8	≤0.25	-	-	-	-
Erythromycin	E	1 – 4	≤0.25 – 1	-	-	-	-
Gentamicin	GM	-	≤0.5 – 1	-	-	-	-
Gentamicin High Level (synergy)	HLG	S	-	R	-	-	-
Inducible Clindamycin Resistance	ICR	-	-	-	-	NEG	POS
Levofloxacin	LEV	0.25 – 2	≤0.12 – 0.5	-	-	-	-
Linezolid	LNZ	1 – 4	1 – 4	-	-	-	-
Moxifloxacin	MXF	≤0.25 – 0.5	≤0.25	-	-	-	-
Nitrofurantoin	FT	≤16	≤16 – 32	-	-	-	-
Oxacillin	OX1	-	≤0.25 – 0.5	-	-	-	-
Quinupristin/Dalfopristin	QDA	2 – 8	≤0.25 – 1	-	-	-	-
Rifampicin	RA	-	≤0.5	-	-	-	-
Streptomycin High Level (synergy)	HLS	S	-	R	-	-	-
Tetracycline	TE	8 – ≥16	≤1	-	-	-	-
Tigecycline[NS]	TGC	≤0.12	≤0.12 – 0.25	-	-	-	-
Trimethoprim/Sulfamethoxazole[c]	SXT	-	≤10 (0.5/9.5)	-	-	-	-
Vancomycin	VA	1 – 4	≤0.5 – 2	-	-	-	-

Numerical values are expressed in µg/ml.
NEG = Negative
POS = Positive
[NS] = The current absence of resistant isolates precludes defining any results other than susceptible. Isolates yielding MIC results suggestive of Nonsusceptible category should be submitted to a reference laboratory for further testing.
[c] = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.

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The ability of the AST card to detect resistance with the following combination(s) is unknown because resistant strains were not available at the time of comparative testing:
- Ampicillin: *Streptococcus agalactiae*
- Benzylpenicillin: *Streptococcus agalactiae*
- Linezolid: *Enterococcus* spp., *Staphylococcus* spp., *Streptococcus agalactiae*
- Tigecycline: *Enterococcus* spp., *Staphylococcus* spp., *Streptococcus* spp.

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli
[c] = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.
**Kox(-ESBL) = *K. oxytoca* (excluding ESBL producing strains)
Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid. For ESBL positive strains, the test interpretation should be reported as resistant for all penicillins, cephalosporins, and aztreonam.
NEG = Negative
POS = Positive
A negative ESBL test result does not rule out the presence of an ESBL masked by an AmpC beta-lactamase.
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

QUALITY CONTROL

		CLSI® Quality Control Organisms VITEK 2 Results			
Antimicrobic	Code	*E. coli* ATCC® 25922™	*P. aeruginosa* ATCC® 27853™	*E. coli* ATCC® 35218™	*K. pneumoniae* ssp *pneumoniae* ATCC® 700603™
Amoxicillin/Clavulanic Acid	AMC	≤2/1 – 8/4	-	4/2 – 16/8	-
Ampicillin	AM	≤2 – 8	-	-	-
Ampicillin/Sulbactam	SAM	≤2/1 – 8/4	-	8/4 – ≥32/16	-
Cefazolin	CZ	≤4	-	-	-
Cefepime	FEP	≤1	≤1 – 8	-	-
Ceftazidime	CAZ	≤1	≤1 – 4	-	-
Ceftriaxone	CRO	≤1	8 – ≥64	-	-
Ciprofloxacin	CIP	≤0.25	≤0.25 – 1	-	-
Ertapenem[c]	ETP	≤0.5	2 – ≥8	-	-
ESBL'	ESB	NEG	-	-	POS
Gentamicin	GM	≤1	≤1 – 2	-	-
Imipenem	IPM	≤0.25	1 – 4	-	-
Levofloxacin	LEV	≤0.12	0.5 – 4	-	-
Nitrofurantoin	FT	≤16	-	-	-
Piperacillin/Tazobactam[c]	TZP	≤4/4* *(*FDA/CLSI Broth Microdilution expected QC range = 1/4 – 4/4 µg/ml)*	≤4/4 – 8/4* *(*FDA/CLSI Broth Microdilution expected QC range = 1/4 – 8/4 µg/ml)*	≤4/4* *(*FDA/CLSI Broth Microdilution expected QC range = 0.5/4 – 2/4 µg/ml)*	-
Tobramycin	TM	≤1	≤1	-	-
Trimethoprim/Sulfamethoxazole	SXT	≤20 (1/19)	160 (8/152) – ≥320 (16/304)	-	-

Numerical values are expressed in µg/ml.
[c] = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.
Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid. For ESBL positive strains, the test interpretation should be reported as resistant for all penicillins, cephalosporins, and aztreonam.
NEG = Negative
POS = Positive
A negative ESBL test result does not rule out the presence of an ESBL masked by an AmpC beta-lactamase.

LIMITATIONS

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):
- Amoxicillin/Clavulanic Acid: *Providencia* spp.

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Perform an alternative method of testing prior to reporting results when a resistant result is obtained with the following antibiotic/organism combination(s):
- Imipenem: *Aeromonas* spp.
- Piperacillin/Tazobactam: *Pseudomonas aeruginosa*

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli
**Kox(-ESBL) = *K. oxytoca* (excluding ESBL producing strains)
**N/A = No specific FDA Indications for Use available
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli
**N/A = No specific FDA Indications for Use available
**Kox(-ESBL) = *K. oxytoca* (excluding ESBL producing strains)
Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid. For ESBL positive strains, the test interpretation should be reported as resistant for all penicillins, cephalosporins, and aztreonam.
NEG = Negative
POS = Positive
A negative ESBL test result does not rule out the presence of an ESBL masked by an AmpC beta-lactamase.
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

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AST-N215

IVD

Gram Negative Susceptibility Card

SUMMARY AND EXPLANATION

The VITEK® 2 Gram Negative Susceptibility Card is intended for use with the VITEK 2 Systems in clinical laboratories as an in vitro test to determine the susceptibility of **clinically significant aerobic gram-negative bacilli** to antimicrobial agents when used as instructed in the Product Information manual.

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli

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❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

QUALITY CONTROL

		CLSI® Quality Control Organisms VITEK 2 Results		
Antimicrobic	Code	*E. coli* ATCC® 25922™	*P. aeruginosa* ATCC® 27853™	*E. coli* ATCC® 35218™
Amoxicillin/Clavulanic Acid	AMC	≤2/1 – 8/4	-	4/2 – 16/8
Ampicillin	AM	≤2 – 8	-	-
Cefalotin	CF	4 – 16	-	-
Cefotaxime	CTX	≤1	8 – 32	-
Cefoxitin	FOX	≤4 – 8	-	-
Cefpodoxime	CPD	≤0.25 – 1	-	-
Ceftazidime	CAZ	≤1	≤1 – 4	-
Cefuroxime	CXM	2 – 8	-	-
Ciprofloxacin	CIP	≤0.25	≤0.25 – 1	-
Gentamicin	GM	≤1	≤1 – 2	-
Imipenem	IPM	≤0.25	1 – 4	-
Meropenem	MEM	≤0.25	≤0.25 – 1	-
Nitrofurantoin	FT	≤16	-	-
Norfloxacin	NOR	≤0.5	1 – 4	-
Piperacillin/Tazobactam	TZP	≤4/4	≤4/4 – 8/4	≤4/4
Tetracycline	TE	≤1 – 2	8 – ≥16	-
Tobramycin	TM	≤1	≤1	-
Trimethoprim/Sulfamethoxazole	SXT	≤20 (1/19)	160 (8/152) – ≥320 (16/304)	-

Numerical values are expressed in µg/ml.

LIMITATIONS

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):
- Amoxicillin/Clavulanic Acid: *Providencia* spp.
- Ampicillin: *Citrobacter* spp., *Enterobacter* spp., *Pantoea* spp., *Serratia* spp., *Cronobacter sakazakii*
- Cefpodoxime: *Morganella morganii, Serratia* spp.
- Imipenem: *Serratia marcescens*
- Piperacillin/Tazobactam: *Serratia marcescens*

Perform an alternative method of testing prior to reporting of results when an intermediate or resistant result is obtained with the following antibiotic/organism combination(s):
- Amoxicillin/Clavulanic Acid: *Proteus* spp.

Perform an alternative method of testing prior to reporting results when a resistant result is obtained with the following antibiotic/organism combination(s):
- Imipenem: *Aeromonas* spp.
- Meropenem: *Aeromonas* spp.
- Piperacillin/Tazobactam: *Pseudomonas aeruginosa*

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AST-N340

IVD

Gram Negative Susceptibility Card

SUMMARY AND EXPLANATION

The VITEK® 2 Gram Negative Susceptibility Card is intended for use with the VITEK 2 Systems in clinical laboratories as an in vitro test to determine the susceptibility of **clinically significant aerobic gram-negative bacilli** to antimicrobial agents when used as instructed in the Product Information manual.

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli
**N/A = No specific FDA Indications for Use available
**Kox(-ESBL) = *K. oxytoca* (excluding ESBL producing strains)
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

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Numerical values are expressed in µg/ml.

LIMITATIONS

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):
- Amikacin: *Acinetobacter baumannii*
- Amoxicillin/Clavulanic Acid: *Providencia* spp.
- Ampicillin: *Citrobacter* spp., *Enterobacter* spp., *Pantoea* spp., *Serratia* spp., *Cronobacter sakazakii*
- Piperacillin/Tazobactam: *Serratia marcescens*

Perform an alternative method of testing prior to reporting of results when an intermediate or resistant result is obtained with the following antibiotic/organism combination(s):
- Amoxicillin/Clavulanic Acid: *Proteus* spp.

Perform an alternative method of testing prior to reporting results when a resistant result is obtained with the following antibiotic/organism combination(s):
- Piperacillin/Tazobactam: *Pseudomonas aeruginosa*

The ability of the AST card to detect resistance with the following combination(s) is unknown because resistant strains were not available at the time of comparative testing:
- Ertapenem: *E. coli, K. pneumoniae, C. freundii, C. koseri, E. aerogenes, E. cloacae, K. oxytoca* (excluding ESBL producing strains), *M. morganii, P. mirabilis, P. vulgaris, S. marcescens*

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli
NS = The current absence of resistant isolates precludes defining any results other than susceptible. Isolates yielding MIC results suggestive of Nonsusceptible category should be submitted to a reference laboratory for further testing.
c = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.
**Kox(-ESBL) = *K. oxytoca* (excluding ESBL producing strains)
Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid. For ESBL positive strains, the test interpretation should be reported as resistant for all penicillins, cephalosporins, and aztreonam.
NEG = Negative
POS = Positive
A negative ESBL test result does not rule out the presence of an ESBL masked by an AmpC beta-lactamase.
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

QUALITY CONTROL

		CLSI® Quality Control Organisms VITEK 2 Results			
Antimicrobic	Code	*E. coli* ATCC® 25922™	*P. aeruginosa* ATCC® 27853™	*E. coli* ATCC® 35218™	*K. pneumoniae* ssp *pneumoniae* ATCC® 700603™
Amikacin	AN	≤2 – 4	≤2 – 4	-	-
Ampicillin	AM	≤2 – 8	-	-	-
Ampicillin/Sulbactam	SAM	≤2/1 – 8/4	-	8/4 – ≥32/16	-
Cefazolin	CZ	≤4	-	-	-
Cefepime	FEP	≤1	≤1 – 8	-	-
Ceftazidime	CAZ	≤1	≤1 – 4	-	-
Ceftriaxone	CRO	≤1	8 – ≥64	-	-
Ciprofloxacin	CIP	≤0.25	≤0.25 – 1	-	-
Doripenem[NS]	DOR	≤0.12	≤0.12 – 0.5	-	-
Ertapenem[c]	ETP	≤0.5	2 – ≥8	-	-
ESBL[·]	ESB	NEG	-	-	POS
Gentamicin	GM	≤1	≤1 – 2	-	-
Levofloxacin	LEV	≤0.12	0.5 – 4	-	-
Nitrofurantoin	FT	≤16	-	-	-
Piperacillin/Tazobactam[c]	TZP	≤4/4* *(*FDA/CLSI Broth Microdilution expected QC range = 1/4 – 4/4 µg/ml)*	≤4/4 – 8/4* *(*FDA/CLSI Broth Microdilution expected QC range = 1/4 – 8/4 µg/ml)*	≤4/4* *(*FDA/CLSI Broth Microdilution expected QC range = 0.5/4 – 2/4 µg/ml)*	-
Tobramycin	TM	≤1	≤1	-	-
Trimethoprim/Sulfamethoxazole	SXT	≤20 (1/19)	160 (8/152) – ≥320 (16/304)	-	-

Numerical values are expressed in µg/ml.
NS = The current absence of resistant isolates precludes defining any results other than susceptible. Isolates yielding MIC results suggestive of Nonsusceptible category should be submitted to a reference laboratory for further testing.
c = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.
Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid. For ESBL positive strains, the test interpretation should be reported as resistant for all penicillins, cephalosporins, and aztreonam.
NEG = Negative
POS = Positive
A negative ESBL test result does not rule out the presence of an ESBL masked by an AmpC beta-lactamase.

LIMITATIONS

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):
- Amikacin: *Acinetobacter baumannii*
- Ampicillin: *Citrobacter* spp., *Enterobacter* spp., *Pantoea* spp., *Serratia* spp., *Cronobacter sakazakii*

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LIMITATIONS

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):
- Amikacin: *Acinetobacter baumannii*
- Ampicillin: *Citrobacter* spp., *Enterobacter* spp., *Pantoea* spp., *Serratia* spp., *Cronobacter sakazakii*
- Ampicillin/Sulbactam: *Citrobacter* spp., *Enterobacter* spp., *Pantoea* spp., *Serratia* spp., *Cronobacter sakazakii*

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- Aztreonam: *Pseudomonas* spp.
- Imipenem: *Serratia marcescens*

Perform an alternative method of testing prior to reporting results when a resistant result is obtained with the following antibiotic/organism combination(s):
- Imipenem: *Aeromonas* spp.
- Meropenem: *Aeromonas* spp.

The ability of the AST card to detect resistance with the following combination(s) is unknown because resistant strains were not available at the time of comparative testing:
- Ertapenem: *E. coli, K. pneumoniae, C. freundii, C. koseri, E. aerogenes, E. cloacae, K. oxytoca* (excluding ESBL producing strains), *M. morganii, P. mirabilis, P. vulgaris, S. marcescens*
- Tigecycline: *C. freundii, E. cloacae, E. coli, K. oxytoca, K. pneumoniae, C. koseri, E. aerogenes, S. marcescens*

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
[c] = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.

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**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

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LIMITATIONS

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):
- Amoxicillin/Clavulanic Acid: *Providencia* spp.
- Ampicillin: *Citrobacter* spp., *Enterobacter* spp., *Pantoea* spp., *Serratia* spp., *Cronobacter sakazakii*
- Cefpodoxime: *Morganella morganii, Serratia* spp.
- Piperacillin/Tazobactam: *Serratia marcescens*

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Numerical values are expressed in µg/ml.
§ Equivalent standard method concentration by efficacy.
**CSAGNB = Clinically significant aerobic gram negative bacilli
[c] = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.
**Kox(-ESBL) = *K. oxytoca* (excluding ESBL producing strains)
❶,❷ etc. = See performance characteristics identified by the drug code with this symbol in the Comment column in Systems Product Information.
**Eco(+ETEC) = *E. coli* (including susceptible enterotoxigenic strains implicated in traveler's diarrhea)

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LIMITATIONS

Perform an alternative method of testing prior to reporting of results for the following antibiotic/organism combination(s):
- Amikacin: *Acinetobacter baumannii*
- Ampicillin/Sulbactam: *Citrobacter* spp., *Enterobacter* spp., *Pantoea* spp., *Serratia* spp., *Cronobacter sakazakii*
- Aztreonam: *Pseudomonas* spp.
- Piperacillin/Tazobactam: *Serratia marcescens*

Perform an alternative method of testing prior to reporting results when a resistant result is obtained with the following antibiotic/organism combination(s):
- Meropenem: *Aeromonas* spp.
- Piperacillin/Tazobactam: *Pseudomonas aeruginosa*

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Numerical values are expressed in µg/ml.
[c] = Category agreement was established at the time of FDA clearance. Essential agreement was not established since test contains less than five discrete dilutions.
Note: For ESBL, FEP is Cefepime, CTX is Cefotaxime, CAZ is Ceftazidime and CA is Clavulanic Acid. For ESBL positive strains, the test interpretation should be reported as resistant for all penicillins, cephalosporins, and aztreonam.
NEG = Negative
POS = Positive

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1. Thorpe TC, Wilson ML, Turner JE, DiGuiseppi JL, Willert M, Mirrett S, et al: BacT/ALERT: an Automated Colorimetric Microbial Detection System. *J Clin Micro* 1990;28(7), 1608-1612.
2. *Biosafety in Microbiological and Biomedical Laboratories (BMBL)* 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
3. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC. *Color Atlas and Textbook of Diagnostic Microbiology*, 6th ed. 2006, pp. 446,590.
4. Baron EJ, Weinstein MP, Dunne Jr. WM, Yagupsky P, Welch DF, Wilson DM. 2005. Cumitech 1C, Blood Cultures IV. Coordinating ed. Baron EJ. ASM Press, Washington, DC.
5. CLSI. *Principles and Procedures for Blood Cultures; Approved Guideline.* CLSI document M47-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
6. CLSI/NCCLS. *Quality Control for Commercially Prepared Microbiological Culture Media*; Approved Standard—Third Edition. CLSI/NCCLS document M22-A3. Wayne, PA: NCCLS; 2004.

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1. Thorpe TC, Wilson ML, Turner JE, DiGuiseppi JL, Willert M, Mirrett S, et al: BacT/ALERT: an Automated Colorimetric Microbial Detection System. *J Clin Micro* 1990;28(7), 1608-1612.
2. *Biosafety in Microbiological and Biomedical Laboratories (BMBL)* 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
3. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC. *Color Atlas and Textbook of Diagnostic Microbiology*, 6th ed. 2006, pp. 446,590.
4. Baron EJ, Weinstein MP, Dunne Jr. WM, Yagupsky P, Welch DF, Wilson DM. 2005. Cumitech 1C, Blood Cultures IV. Coordinating ed. Baron EJ. ASM Press, Washington, DC.
5. CLSI. *Principles and Procedures for Blood Cultures; Approved Guideline*. CLSI document M47-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
6. CLSI. *Quality Control for Commercially Prepared Microbiological Culture Media*; Approved Standard—Third Edition. CLSI/NCCLS document M22-A3. Wayne, PA: NCCLS; 2004.
7. Kondratovich M. Comparing Two Medical Tests When Results of Reference Standard Are Unavailable for Those Negative via Both Tests. *J. Biopharm Stat* 2008;18:1,145-166.

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1. Thorpe TC, Wilson ML, Turner JE, DiGuiseppi JL, Willert M, Mirrett S, et al: BacT/ALERT: an Automated Colorimetric Microbial Detection System. *J Clin Micro* 1990;28(7), 1608-1612.
2. *Biosafety in Microbiological and Biomedical Laboratories (BMBL)* 5th Edition. U.S. Department of Health and Human Services Centers for Disease Control and Prevention and National Institutes of Health. Fifth Edition. US Government Printing Office. Washington: Feb 2007.
3. Koneman EW, Allen SD, Janda WM, Schreckenberger PC, Winn WC. *Color Atlas and Textbook of Diagnostic Microbiology*, 6th ed. 2006, pp. 446,590.
4. Baron EJ, Weinstein MP, Dunne Jr. WM, Yagupsky P, Welch DF, Wilson DM. 2005. Cumitech 1C, Blood Cultures IV. Coordinating ed. Baron EJ. ASM Press, Washington, DC.
5. CLSI. *Principles and Procedures for Blood Cultures; Approved Guideline.* CLSI document M47-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
6. CLSI. *Quality Control for Commercially Prepared Microbiological Culture Media*; Approved Standard—Third Edition. CLSI/NCCLS document M22-A3. Wayne, PA: NCCLS; 2004.
7. Kondratovich M. Comparing Two Medical Tests When Results of Reference Standard Are Unavailable for Those Negative via Both Tests. *J. Biopharm Stat* 2008;18:1,145-166.

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