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5 ºC 35 ºC
Temperature limitation
Temperatura de conservación

LATEX
Does not contain latex
No contiene látex

Keep away from sunlight
No exponer a la luz solar

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ENGLISH

DESCRIPTION

RRS® HA Long Lasting is a sterile medical device, Class III, made in conformity with Directive 93/42/EEC, dermal implant that contains cross-linked, resorbable Hyaluronic acid (HA). HA is solubilized in a protective buffer solution composed by amino acids with hydration and antioxidant properties that stabilize the HA. Hyaluronic acid is from biotechnological non-animal origin. This medical device only can be injected by a legally authorized physician in a legally authorized clinic.

PRESENTATION

Box of 1 or more syringe(s) of 3 mL (0,1 Fl Oz), traceability labels, instructions for use.
The number of syringes is printed on the box.

COMPOSITION

BDDE cross-linked hyaluronic acid	21 mg / 3 mL
Protective buffer	14,52 mg / 3 mL

INDICATIONS

RRS® HA Long Lasting is a resorbable dermal implant for the treatment of skin photoaging and its consequences, treatment of dermal atrophy.

HOW TO USE

RRS® HA Long Lasting is recommended for deep dermis injections. If the needle gets obstructed, it has to be changed for a new one.
Any injection may leave small ecchymoses points or even haematomas, mainly on patients taking acetylsalicylic acid, anti-inflammatory of taking medication delaying coagulation. It does not often appear when the product is injected at a depth of more than 1mm. Sensitive skins may benefit from the application of an anaesthetic cream prior to the treatment; said cream may cause redness or skin hypersensitivity.
Before implantation, physician should check if patient is healthy and should inform him/her about the possibilities of side effects (pain, redness, ecchymoses, stinging sensation and swelling, local inflammation, usually disappearing in 24-48 hours).

RECOMMENDATIONS FOR PATIENTS

After the treatment and during 12/24 hours: avoid extreme temperatures, saunas and Hammam, direct exposure to sun or UV. Avoid to make-up or apply products other than those advised by the doctor. The physician must to inform the patient of possible local reactions related to the implantation of this resorbable device: redness and/or local inflammation which disappears within 12-24hours. Slight oedema or small ecchymoses may occur, which usually disappear within 24 or 48 hours. Should any inflammatory reaction or any other side effect not disappeared within a week, the patient must inform the clinician about it immediately. The practitioner will treat the patient via appropriate treatment adapted to the problem. Rare cases of abscesses or granulomas have been reported in medical literature after injection of cross-linked hyaluronic acid.
The doctor is requested to report any notable side effect related to the medical devices **RRS® HA Long Lasting** to **Skin Tech Pharma Group SL.** through the website **www.skintechpharmagroup.com**, by e-mail at **info@skintechpharmagroup.com** or by post.

CONTRAINDICATIONS

Do not inject **RRS® HA Long Lasting** in the following cases:
- Allergy to any of the ingredients.
- Patients presenting any skin alteration, skin disease, infections, sequelae of streptococcal infections.
- Patients taking immunosuppressants, undergoing cortical therapy, with autoimmune disease history, patients with uncompensated diabetes, acute joint rheumatics, repetitive angina, endocarditis.
- **RRS® HA Long Lasting** must not be injected in anatomic hazard zones.
- No studies are available for use during pregnancy and breastfeeding or children or minors under 18.

PRECAUTIONS

In case of allergic risk, a test should be done before the first injection. The test results must be attached to the patient's medical record. Do not re-use the syringe once opened. Do not use if the external protecting package has been seriously damaged or if you see some liquid inside of the box. Do not re-use or re-sterilize the unused leftover product; dispose the leftover products away; dispose the needle and the syringe into specific containers.
Inform the patients under anticoagulant treatment about the risk of ecchymoses or haematomas due to injections. It is recommended that the patient should avoid taking acetylsalicylic acid or anti-inflammatory drug during the week preceding injections. Do not inject into the blood vessels, bones, tendons, muscles, ligaments or any organ different of skin.
Prior to injecting, check the expiry date of syringes and needles; also check the integrity of the syringe seal. Do not use a product which is expired or looks damaged. Caution: the outside of the syringe is not sterile. Once the syringe has been opened, use immediately. Do not heat, do not frost the presentation.
Inform the patient of product characteristics and possible side effects, precautions and contraindications. In the absence of available clinical data, it is not recommended to use **RRS® HA Long Lasting** in anatomic hazard zones. There is no study available regarding cross tolerance or possible reaction between the medical device **RRS® HA Long Lasting** and other filler products, even if it is not injected in the same plane or area. Also, no clinical data is available about associating injections of **RRS® HA Long Lasting** with lasers, dermabrasion or any energy-based device etc.
Hyaluronic acid precipitates when coming into contact with quaternary ammonium cations (such as Benzalkonium chloride). Therefore, do not recommend the medical device **RRS® HA Long Lasting** coming into contact with said product.

STORAGE AND DISPOSAL

Keep away from the light, under the temperature between 5 and 35 ºC. Do not freeze, do not heat.
Dispose of each component of the used syringe and packaging of product in accordance with the laws of the country where it is used.

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