Exadrop®

GB Instructions for use 3

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



GB Precision flow control device with/ for I.V. administration set. Only for gravity infusion.

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GB Instructions for use

Features
- Precision flow control device for gravity infusions
- Figures on the scale (= app. ml/h) facilitate the setting of the required flow rate
- with integrated I. V. administration set Intrafix® Air acc. to DIN EN ISO 8536-4, 15 µm fluid filter
- constant flow rate for prolonged administration
- one-handed use
- no resetting of flow rate
 flow rate once set is not influenced by cold flow of tubing
- slide clamp for short interruption of administration
- Lock connector
- pressure resistant to 50 kPa (0.5 bar)
- loop on handle is used as tube holding device
- Sterile, non-pyrogenic

Material used
PE, PS, SB, PP, PA, silicone rubber, borosilicate glass fibrous web, ABS, PC, PVC (in REF 4062264: PUR instead of PVC)

Indications
Gravity infusions in general

Contra-indications
Pressure > 50 kPa

Side-effects
not known

Precautions
Do not resterilize!
Do not use with pressure infusion devices
(infusion pumps)!
Not suitable for conserved blood.

Warning
Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

Handling
1. Rotate wheel of control device towards max. in order to release initial position.
 Set control device to „0" (Fig. 1).
2. Slide clamp open.
3. Prepare infusion as usual. Insert Intrafix® Air administration set with closed snap vent in fluid container (Fig. 2).
4. Adjust fluid level in drop chamber.
5. Open snap vent.
6. Set control device to "max" and fill I. V. administration set after having expelled air from tubing (Fig. 3).
7. Remove protective cap of Luer Lock connector.
8. Establish connection to the patient by means of an I. V. cannula for example Vasofix® Braunule® (Fig. 4).
9. Set desired flow rate on the markings of the handle by working back from "max" (Fig. 5). **The figures on the scale (= app. ml/h) are only approximate ones.** They only refer to the infusion of 0.9 % NaCl solution administered by means of an Vasofix® Braunule® G 18 being approximately 76 cm statical difference of level without CVP.
 Therefore the figures on the scale must always be controlled by verification.
 Otherwise there is a risk of over- or underdosage for the patient with all subsequent consequences.
10. The **loop** on the Exadrop®-handle is used as **tube holding device** when the connection of the I. V. administration set to the I. V. cannula cannot be established immediately.
11. With the **slide clamp** above the control device the infusion can shortly be interrupted without setting back the control device to "0". After having opened the clamp the flow rate remains the same if the pressure conditions in the infusion system have not changed in the meantime.

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B. Braun Melsungen AG,