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CARTO™ 3 System External Reference Patches

- **NON-STERILE**
- **For single use only**
- **Do not use if package is opened or damaged**

DEVICE DESCRIPTION
The Biosense Webster CARTO™ 3 System External Reference Patches allow non-fluoroscopic catheter navigation when used with the CARTO™ 3 System. The patches provide a means for attachment of the patch sensor cables to the patient in order to convey catheter location. Please refer to the CARTO™ 3 System's User Manual for additional connectivity information.

INDICATIONS
The Biosense Webster CARTO™ 3 System External Reference Patches are indicated for use with the CARTO™ 3 System to provide a means of attachment of the patch sensor cable to the patient in order to convey catheter location.

CONTRAINDICATIONS
There are no known contraindications for the use of the external reference patches.

PATIENT TARGET GROUP
This device is targeted at patients who have been diagnosed with cardiac arrythmias and are undergoing an electrophysiology procedure.

INTENDED USERS
The intended users of this device are appropriately trained personnel in a fully equipped electrophysiological laboratory.

CLINICAL BENEFIT
The external reference patch provides a means for attaching the patch sensor cables to the patient in order to convey catheter location and information on the patient's respiratory cycle.

WARNINGS

- Discontinue use of this product if skin sensitivity is observed.

PRECAUTIONS
- Do not overlap the External Reference Patches with any other patches on the body including but not limited to EKG patches, defibrillator patches, indifferent electrode patches, etc. or any other External Reference Patch.
- Do not place the External Reference Patches over devices such as pacemakers or ICD's.
- Avoid large bony areas such as the scapula and sternum when placing the External Reference Patches.
- Do not expose the external reference patch to organic solvents such as alcohol. If an organic solvent is used on patient, allow to dry for 1 minute prior to applying patches.
- The External Reference Patches should only be connected to the CARTO™ 3 System.
- Do not attempt to use the External Reference Patches prior to completely reading and understanding these Instructions For Use. This device is for external use only.
- Refer to the product label for recommended storage conditions. Inspect the packaging and external reference patches prior to use.
- If all External Reference Patches are not used during a procedure, ensure that unused patches are stored in the original, sealed packaging. If an External Reference Patch has not been stored in its original, sealed packaging, the patches should not be used.
- The External Reference Patches are intended for single patient use only.
- Do not sterilize or reuse.
- This device is packaged for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- If the External Reference Patch does not adhere well to the patient's skin, replace with a new patch.

"USE BY" DATE
Use the device prior to the "Use By" date on the package label.

REPORTING INCIDENTS
Per regulation 2017/745/EU on medical devices, if this device is used in the European Union (or in a country with an identical regime) and a serious incident occurs during the use of this device or as a result of use of this device, report the incident to the manufacturer and/or its authorized representative and to your national authority.

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
DISPOSAL

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

INSTRUCTIONS FOR USE

Please refer to this Instructions for Use and the CARTO™ 3 System's User Manual when using the CARTO™ 3 System External Reference Patches in conjunction with the CARTO™ 3 System.

1. Ensure that the skin surface is completely shaved and dry.
2. Remove each External Reference Patch from the package.
3. Inspect the External Reference Patch to ensure that it is in good physical condition. Care should be taken when opening the packaging to ensure that the sealing mechanism is not damaged.
4. The patches should be placed surrounding the heart and within the CARTO™ 3 mapping zone. The following is the recommended scheme for patch placement.
5. Back Patch Placement:
 - Attach three (3) patches to the patch sensor cables, and then attach to the patient's back while the patient is sitting.
 - It is recommended to identify the center of the heart, slightly left of the spine, approximately at the level of the 7th thoracic vertebra (T7).
 - Place the center of the first External Reference Patch approximately 12 cm (4.725") superior to the center of the heart.
 - Place the center of the second External Reference Patch approximately 12 cm (4.725") inferior and to the right of the center of the heart.
 - Place the center of the third External Reference Patch approximately 12 cm (4.725") inferior and to the left of the center of the heart.
 - The back patches should be positioned in a triangle as demonstrated in the Figure 1 below. After the patches have been attached to the patient's back, the patient should lie down.
6. Slide the location pad under the patient's table, such that the three back patches are in the center of the location pad. Verify the position of the patches using the "Location Set Up Phase" on the CARTO™ 3 System. Adjust the location pad location if required.
7. Chest Patch Placement - The chest patches should be positioned anatomically according to the instructions below. Refer to Figure 2.
 - While the patient is lying down, attach three (3) patches to the patch sensor cables.
 - Place the first patch on the right lower rib cage (5th / 6th intercostals space) at the right mid-clavicular line. Use the Location Setup phase to verify the patch is within the accurate mapping zone.
 - Place second patch on the left lower rib cage (5th / 6th intercostals space) at the anterior axillary line (under V5 ECG patch). Use the Location Setup phase to verify the patch is within the accurate mapping zone.
 - Place the third patch under the clavicle, left or right to the sternum (depending of the Defibrillation Pad position and presence of a device such as ICD or pace maker, etc.). Use the Location Setup phase to verify the patch is within the accurate mapping zone.
8. Using the Location Setup screen, verify that all chest patches are located within the boundaries of the accurate mapping zone.

Patch Placement


Figure 2: Chest Patch Placement

For additional education and materials for patch placement, please contact your local Biosense Webster representative.

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ADVERSE REACTIONS

A number of adverse reactions have been documented for electrophysiology procedures including:

heart block, pulmonary vein stenosis, esophageal fistula and/or injury, stroke (cerebrovascular accident), other arrhythmias (outside diagnosis), life threatening arrhythmias, myocardial infarction, cardiac perforation, pericardial effusion, cardiac tamponade, thrombosis, embolism, pulmonary embolism, air embolism, valvular damage, phrenic nerve injury, vagal nerve injury, pericarditis, coronary artery stenosis, vessel perforation (peripheral and/or central), soft tissue injury, persistent atrial communication, device related infection, embolization of components, device entrapment, surgical intervention (additional), skin burns, and localized skin reaction.

The following complications associated with cardiac catheterization have also been reported in the literature:

major bleed, hematoma, reaction to medications, allergic reaction, vascular access complication, damage to vasculature, implanted device interactions, renal artery stenosis, pneumothorax, ST segment changes, fluid overload, urinary catheter complications, hypotension, sepsis, wound infection, respiratory failure, heart injury, renal injury, heart failure, cardiac arrest and death.

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

ELECTRONIC INSTRUCTIONS FOR USE

This document is available at www.e-ifu.com.

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© Biosense Webster, Inc. 2006-2021
are trademarks of Biosense Webster, Inc.
The third-party trademarks used herein are trademarks of their respective owners.