DePuy Synthes™ – CMW™
GENERAL CEMENTING INSTRUMENTS

(en) Instructions for Use
(de) Gebrauchsanweisung
(es) Instrucciones De Uso
(fr) Mode D'emploi
(it) Indicazioni D'uso
(pt) Instruções De Utilização
(nl) Gebruiksaanwijzing
(da) Brugsanvisning
(sv) Bruksanvisning
(no) Bruksanvisning
(fi) Käyttöohjeet
(bg) Инструкции За Употреба
(hr) Upute Za Uporabu
(cz) Návod K Použití
(et) Kasutusjuhend
(el) Οδηγιες Χρησης
(hu) Használati Útmutató
(lv) Lietošanas Instrukcija
(lt) Naudojimo Instrukcijos
(pl) Instrukcja Użycia
(ro) Instrucțiuni De Utilizare
(sk) Návod Na Použitie
(sl) Navodila Za Uporabo
(tu) Kullanma Talimatlari
(id) Petunjuk Penggunaan
(ko) 사용 지침
(ru) Инструкция по применению

IFU-103587318
Rev A
Date: 2019/09

DePuy Synthes
PART OF THE Johnson & Johnson FAMILY OF COMPANIES

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
103587318 | Rev:1
Released: 05 Mar 2020
CO: 103622056
Release Level: 4. Production

(en) VACUMIX MIXING STAND, SYRINGE HOLDER

(de) Vacu-Mix Mischständer, Spritzenhalter
(es) Soporte de pie de mezclado Vacu-Mix, Soporte de jeringa
(fr) Support de mélange Vacu-Mix, Support seringue
(it) Supporto di miscelazione Vacu-Mix, Porta siringa
(pt) Suporte de Mistura Vacu-Mix, Porta-seringas
(nl) Vacu-Mix-mengstandaard, Injectiespuithouder
(da) Vacu-Mix blandestativ, Sprøjteholder
(sv) Vacu-Mix-blandningsställ, Spruthållare
(no) Vacu-Mix blandingsstativ, Sprøyteholder
(fi) Vacu-Mix-sekoitusteline, Ruiskupidike
(bg) Статив за смесване Vacu-Mix, Държач за шприц
(hr) Stalak za miješanje Vacu-Mix, Držač štrcaljke
(cz) Stojan na míchání Vacu-Mix, Držák stříkačky
(et) Vacu-Mix segamisalus, Süstla hoidik
(el) Βάση ανάμιξης Vacu-Mix, Στήριγμα σύριγγας
(hu) Vacu-Mix keverőállvány, Fecskendőtartó
(lv) Vacu-Mix jaukšanas statīvs, Šļirces turētājs
(lt) „Vacu-Mix“ maišymo stovas, Švirkšto laikiklis
(pl) Stojak do mieszania Vacu-Mix, Uchwyt strzykawki
(ro) Suport pentru amestecare Vacu-Mix, Suport pentru seringă
(sk) Stojan na miešanie Vacu-Mix, Držiak striekačky
(sl) Stojalo za mešanje Vacu-Mix, Držalo za brizgo
(tu) Vacu-Mix Mikser Standı, Şırınga Tutacağı
(id) Tiang Pencampur Vacu-Mix, Tempat Alat Suntik
(ko) Vacu-Mix 혼합 스탠드, 주사기 홀더
(ru) Смесительная стойка Vacu-Mix, Держатель шприца

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

This instruction leaflet provides guidance on the use of DePuy Synthes™ – CMW™ – General Cementing Instruments and includes: Nozzle Cutter, Vacu-Mix™ Mixing Stand, Acetabular T-Handle and Syringe Holder.

These instruments can be used during joint replacement surgery in which bone cement and mixing systems are used.

The Nozzle Cutter consists of a cropping blade attached to a spring-loaded assembly.

Acetabular T-Handle is a rod with flange tip and handle.

The Vacu-Mix Mixing Stand and Syringe Holder consists of a tube attached to a circular base.

All General Cementing Instrumentation are supplied in a non-sterile state and must be cleaned and sterilised according to the relevant Instructions For Use (IFU) prior to clinical use.

Materials & Restricted Substances

- Nozzle Cutter: Polyacetal and Stainless Steel
- Acetabular T-Handle: Stainless Steel
- Vacu-Mix Mixing Stand: Aluminium
- Syringe Holder: Stainless Steel

Intended Use

- Nozzle Cutter - The Nozzle Cutter is a reusable device that is used to cut the nozzle of the delivery syringe towards the end of cement extrusion to allow continued pressurisation of the cement in the femoral cavity.
- Acetabular T-Handle - The Acetabular T-Handle is reusable and forms part of the SmartSeal™ Pressuriser System, fitting securely into the locating hole in the back of the SmartSeal Acetabular Pressuriser to hold it in place during cement pressurisation to allow for optimal cement penetration and therefore enhance bone cement interlock.
- Vacu-Mix Mixing Stand - The Mixing Stand is reusable and is used to firmly hold and support the Vacu-Mix syringe barrel and the ampoules of bone cement liquid component during the preparation and mixing of the bone cement.
- Syringe Holder - The Syringe Holder is reusable and is used to firmly hold and support the Standard Cemvac™ Syringe Barrel during the preparation and mixing of the bone cement.

Indications

There are no indications defined specific to these device as indications are based upon the implant devices rather than the instruments themselves.

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Warnings and Cautions

The General Cement Instrumentation must be cleaned and sterilised following the appropriate instructions prior to clinical use.

The General Cement Instrumentation should only be used with the compatible device.

The General Cement Instrumentation should only be operated by or under the direct supervision of a suitably qualified orthopaedic surgeon within a clinical setting.

Care should be taken to not pierce surgical gloves on the General Cement Instrumentation.

The General Cement Instrumentation are only intended to be used within a controlled orthopaedic clinical setting.

Information for Use

All General Cementing Instrumentation are supplied in a non-sterile state and must be cleaned and sterilised according to the relevant IFU prior to clinical use.

Please refer to the diagrammatic which shows the use of the General Cementing Instrumentation.

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Adverse Events & Complications

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

Performance of Instruments

T-Handle – Acetabular pressurisers:

- The size of the pressuriser head is selected according to the size of the acetabulum to be pressurised.
- Secure the chosen pressuriser head to the end of the SmartSeal™ Acetabular 'T' Handle by holding the rod in one hand and the pressuriser head in the other, then place the 'T' Handle in the locating hole in the SmartSeal Acetabular Pressuriser.
- Apply pressure from the end of the 'T' Handle for at least 30 seconds to ensure good pressurisation of the cement into cancellous bone.
- After completion of the pressurisation process, the SmartSeal Acetabular Pressuriser and SmartSeal Acetabular 'T' Handle should be removed together.

Device Processing

Prior to sterilisation, store all non-sterile devices in a suitably clean and dry location.

All General Cementing Instrumentation must be cleaned and sterilised prior to clinical use.

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After sterilisation, store the General Cementing Instrumentation in-line with local hospital procedures.

This device should be cleaned and sterilised prior to disposal, removing any source of biological contamination. The device can then be disposed in-line with local regulations.

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Compatible Device

Class I Medical Device (Non-Sterile)	Medical Devices used in conjunction with Class I Medical Devices (Non-Sterile)					
	Vacu-Mix Plus (inc. Prefilled)	Vacu-Mix Syringe Kit	Syringe Kit	SmartSeal Pressurisers	SmartMix Cemvac Syringe (inc. Prefilled)	Standard Cemvac Syringe
Nozzle Cutter	✓	✓	✓			✓
VacuMix Mixing Stand	✓	✓				
Acetabular T-Handle				✓		
Syringe Holder					✓	✓

✓ = Compatible Device

Symbols Guide

REF — Catalogue Number

LOT — Batch Code

SN — Serial Number

Caution

Manufacturer

NON STERILE — Non-sterile

Consult Instructions for Use

www.e-IFU.com
US / Canada 800-255-2500
EU: +800 8888 2020
EU: +32 2 4037222
e-IFU website

Keep Dry

Packaging Quantity

EC | MD — Medical Device

MADE IN — Made In

Do Not Disassemble

Recyclable

℞ Only — Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

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DePuy Synthes

PART OF THE Johnson & Johnson FAMILY OF COMPANIES

* DePuy Ireland UC
Loughbeg, Ringaskiddy,
Co. Cork, Ireland
Tel: +353 21 4914 000
Fax: +353 21 4914 199
DePuySynthes.com

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

(en) CEMVAC 1 PIECE GUN, SMARTMIX CEMVAC GUN

(de) 1-teilige Cemvac Pistole, SmartMix Cemvac Pistole
(es) Pistola de una pieza Cemvac, Pistola SmartMix Cemvac
(fr) Pistolet monopièce Cemvac, Pistolet SmartMix Cemvac
(it) Pistola monopezzo Cemvac, Siringa SmartMix Cemvac
(pt) Pistola de Peça Única Cemvac, Pistola SmartMix Cemvac
(nl) Cemvac 1-delig pistool, SmartMix Cemvac-pistool
(da) Cemvac pistol i ét stykke, SmartMix Cemvac pistol
(sv) Cemvac-pistol i ett stycke, SmartMix Cemvac-pistol
(no) Cemvac I-delspistol, SmartMix Cemvac-pistol
(fi) Yksiosainen Cemvac-pistooli, SmartMix Cemvac -pistooli
(bg) Единичен пистолет Cemvac, Пистолет SmartMix Cemvac
(hr) Jednodijelni pištolj Cemvac, Pištolj SmartMix Cemvac
(cz) Jednodílná pistole Cemvac, Pistole GunMix Cemvac
(et) Cemvac üheosaline püstol, SmartMix Cemvac püstol
(el) Μονοκόμματο πιστόλι Cemvac, Πιστόλι SmartMix Cemvac
(hu) Cemvac egyrészes pisztoly, SmartMix Cemvac pisztoly
(lv) Cemvac viengabala pistole, SmartMix Cemvac cementa pistole
(lt) Neardomas „Cemvac“ pistoletas, „SmartMix Cemvac“ pistoletas
(pl) Jednoczęściowy pistolet Cemvac, Pistolet SmartMix Cemvac
(ro) Pistol aplicator Cemvac cu o componentă, Pistol aplicator SmartMix Cemvac
(sk) Jednodielny aplikátor Cemvac, Aplikátor SmartMix Cemvac
(sl) Enodelna pištola Cemvac, Pištola SmartMix Cemvac
(tu) Cemvac Tek Parça Tabanca, SmartMix Cemvac Tabanca
(id) Injektor Bagian Cemvac I, Injektor SmartMix Cemvac
(ko) Cemvac 원피스 건, SmartMix Cemvac 건
(ru) Однокомпонентный пистолет Cemvac, Пистолет SmartMix Cemvac

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IMPORTANT INFORMATION - Please Read Before Use

en **DEPUY SYNTHES™ – CMW™ – CEMENTING GUNS**

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.
For access to the eIFU and any additional languages visit www.e-IFU.com

This instruction leaflet provides guidance on the use of DePuy Synthes™ – CMW™ – Cementing Guns. This includes; CMW MK III Cement Gun, Cemvac™ 1 Piece Gun and SmartMix™ Cemvac Gun.

All CMW Cementing guns are supplied in a non-sterile state and must be cleaned and sterilised according to the respective Instructions For Use (IFU) prior to clinical use.

Materials & Restricted Substances

The components of the DePuy Synthes – CMW – Cementing Guns are constructed from Stainless Steel, Delrin, Aluminium and Phosphorus Bronze.

Intended Use

- CMW MK III Cement Gun - The CMW MK III Cement Gun is reusable and allows the delivery of high, medium and low viscosity bone cements to the operative site via the compatible cement syringe.
- Cemvac I Piece Gun - The Cemvac I Piece Gun allows the delivery of bone cement to the operative site via the Cemvac syringe system.
- SmartMix Cemvac Gun - The SmartMix Cemvac Gun is reusable and allows the delivery of bone cement to the operative site via the Cemvac and SmartMix Cemvac syringe systems.

The DePuy Synthes – CMW – Cementing Guns within this instruction leaflet are intended for use with patients undergoing orthopaedic surgery.

Indications

There are no indications defined specific to these device as indications are based upon the implant devices rather than the instrument themselves.

Contraindications

There are no contraindications defined specific to these device as contraindications are based upon the implant devices rather than the instruments themselves.

Warnings and Cautions

All CMW Cementing guns must be cleaned and sterilised according to the respective IFU prior to clinical use.

DePuy Synthes – CMW – Cementing Guns should only be used with the compatible mixing or cement delivery systems.

DePuy Synthes – CMW – Cementing Guns should only be operated by or under the direct supervision of a suitably qualified orthopaedic surgeon within an orthopaedic clinical setting.

Care should be taken to not pierce surgical gloves on the cementing guns.

The cementing guns are only intended to be used within an orthopaedic clinical setting.

Information for Use

All CMW Cementing guns must be cleaned prior to clinical use.

The CMW Cementing guns are intended for use during cement application/extrusion from a compatible mixing or cement delivery system. Please consult the compatible device chart at the end of this instruction leaflet for more details.

Detailed diagrammatic showing the correct use of the CMW Cementing guns is within this instruction leaflet.

Adverse Events & Complications

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your local authority.

Serious incidents associated with the use of these devices includes:

- Mechanical failure of the gun resulting in a delay to surgery.

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Device Processing

Prior to sterilisation, store all non-sterile device in a suitably clean and dry location.

All CMW Cementing guns must be cleaned and sterilised prior to clinical use.

Detailed instructions for processing of reusable devices, instrument trays and cases are described in the DePuy Synthes 'Instructions for Processing Medical Devices' IFU.

After sterilisation, store the guns in-line with local hospital procedures to prevent possible contamination.

The device should be cleaned and sterilised prior to disposal, removing any sources of biological contamination. The device can then be disposed in-line with local regulations.

Where further information is desired, please contact your local DePuy sales representative. For instruments produced by another manufacturer, reference the manufacturer's instructions for use.

CMW MK III Cement Gun

1. Rotate the piston drive rod through 90° to disengage the ratchet.
2. Pull back the piston drive rod until it hits the stop.
3. Pull back the economiser drive rod.
4. Rotate the piston drive rod through 90° to re-engage the ratchet, ensuring that the piston driving plate locates fully within the gun housing.
5. Place the economiser rod supplied with the syringe system into the hole in the piston drive rod. IF THE NOZZLE IS TO BE REDUCED IN LENGTH AT ANY STAGE OF CEMENT DELIVERY OR PRESSURISATION THE ECONOMISER ROD MUST NOT BE USED.
6. Push the syringe into the cage ensuring it is pushed fully home at each end.
7. Extrude the cement by pulling the trigger. THE ECONOMISER DRIVE ROD SHOULD ONLY BE ADVANCED USING THE TRIGGER. NEVER USE FORCE ON THE ECONOMISER DRIVE ROD.
8. To remove the syringe, repeat steps 1, 2, and 4, and withdraw empty barrel.

Cemvac 1 Piece Gun and SmartMix Cemvac Gun

1. Press and hold the piston drive rod lever and fully pull back the drive rod.
2. Maintain the syringe in an upright position. Insert the syringe into the metal sleeve.
 * Apply the locking device
3. Extrude the cement using the trigger handle. Inspect cement in the nozzle prior to extrusion.
4. To remove the empty syringe from the gun, remove the nozzle locking device, press and hold the piston drive rod lever and fully pull back the drive rod. Withdraw the empty barrel.

(*locking device only applicable to the Cemvac One Piece Gun)

Compatible Device

	Medical Devices used in conjunction with DePuy Synthes - CMW - Cement Guns			
	Vacu-Mix Syringe Kit (inc. Prefilled)	Syringe Kit	SmartMix Cemvac Syringe (inc. Prefilled)	Standard Cemvac Syringe
CMW MK III Cement Gun	✓	✓		
SmartMix Cemvac Gun			✓	✓
Cemvac One Piece Gun				✓

✓ = Compatible Device

Symbols Guide

Catalogue Number

Batch Code

Serial Number

Caution

Manufacturer

Non-sterile

Consult Instructions for Use

e-IFU website

Keep Dry

Packaging Quantity

Medical Device

Made In

Do Not Disassemble

Recyclable

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

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